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References

1. Crookes RL, Hillier CD. Fresh frozen plasma and related products. In: Hillier CD, Silberstein LE, Ness PM, et al. eds. *Blood Banking and Transfusion Medicine Basic Principles and Practice*. 2nd ed. Philadelphia, PA: Churchill Livingstone Elsevier; 2007:259–269.
2. Cardigan R, Lawrie AS, Mackie IJ, Williamson LM. The quality of fresh-frozen plasma produced from whole blood stored at 4°C overnight. *Transfusion*. 2005;45:1342–1348.
3. Solheim BG, Seghatchian J. Update on pathogen reduction technology for therapeutic plasma: an overview. *Transfus Apher Sci*. 2006;35:83–90.
4. Council of Europe. *Guide to the preparation, use and quality assurance of blood components - 13th edition (2007)*. ISBN 978-92-871-6137-6. Available at: http://book.coe.int/EN/ficheouvrage.php?PAGEID=36&lang=EN&produit_aliasid=2151 Accessed November 2007.
5. UK Blood Transfusion and Tissue Transplantation Services. *Guidelines for the Blood Transfusion Service*. 7th ed. Norwich, England: TSO (The Stationery Office); 2005.
6. EU Directive. Commission Directive 2004/33/EC of 22 March 2004 implementing Directive 2002/98/EC of the European Parliament and of the Council as regards certain technical requirements for blood and blood components.
7. Data on File. Navigant Biotechnologies LLC.
8. Heiden M, Seitz R. Quality of therapeutic plasma – requirements for marketing authorization. *Thromb Res*. 2002;107:S47–S51.
9. European Pharmacopoeia 5th Edition. 01/2007:1646 Human Plasma.
10. Buytaert-Hoelen KA, Hendrix BK, Camp D, et al. Protein quality in Mirasol-treated fresh frozen plasma after one year in storage. Presented at the AABB Annual Meeting, Anaheim, USA, 20–23 October, 2007.
11. Ramasamy I, Farrugia A, Tran E, et al. Biological activity of von Willebrand Factor during the manufacture of therapeutic factor VIII concentrates as determined by the collagen binding assay. *Biologicals*. 1998;26:155–166.
12. Larrea L, Calabuig M, Rivera J, et al. Influence of riboflavin photoinactivation treatment on coagulation factors from fresh frozen plasma. Presented at the ISBT XVIIIth Regional Congress Asia, Hanoi, Vietnam, 10–13 November, 2007.
13. Kumar V, Motheral T, Luzniak G, et al. Mirasol™ Pathogen Reduction Technology: Riboflavin-based process retains functional activities of immunoglobulins. *Transfusion*. 2003;43(Suppl.):60A.
14. Fast LD, DiLeone G, Li J, et al. Functional inactivation of white blood cells by Mirasol treatment. *Transfusion*. 2006;46:642–648.

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References

1. Goodrich RP, Li J, Pieters H, et al. Correlation of in vitro platelet quality measurements with in vivo platelet viability in human subjects. *Vox Sang.* 2006;90:279-285.
2. AuBuchon JP, Herschel L, Roger J, et al. Efficacy of apheresis platelets treated with riboflavin and ultraviolet light for pathogen reduction. *Transfusion.* 2005;45:1335-1341.
3. Li J, Lockerbie O, de Korte D, et al. Evaluation of platelet mitochondria integrity after treatment with Mirasol pathogen reduction technology. *Transfusion.* 2005;45:920-926.
4. Fast LD, DiLeone G, Li J, Goodrich RP. Functional inactivation of white blood cells by Mirasol treatment. *Transfusion.* 2006;46:642-648.

† The Mirasol PRT System for Platelets is not available for sale in the United States.

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