MIRASOL®

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A Novel Approach to Blood Safety

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References

1. Brookmeyer R, Goedert JJ. Censoring in an epidemic with an application to hemophilia-associated AIDS. *Biometrics*. 1989;45:325–335.
2. Janatpour KA, Holland PV: A brief history of blood transfusion. In: Hillyer CD, Silberstein LE, Ness PM, et al., eds. *Blood Banking and Transfusion Medicine Basic Principles and Practice*. 2nd ed, Philadelphia, PA: Churchill Livingstone Elsevier, 2007:3–11.
3. Hervig T. Where will pathogen inactivation have the greatest impact? *ISBT Science Series*. 2007;2:25–29.
4. Chamberland ME. Emerging infectious agents: do they pose a risk to the safety of transfused blood and blood products? *Clin Infect Dis*. 2002;34:797–805.
5. Soldan K, Barbara J. The risks of infection transmission by blood transfusions in England. *J Clin Pathol*. 1999;52:405–408.
6. Blajchman MA, Beckers EA, Dickmeiss E, et al. Bacterial detection of platelets: current problems and possible resolutions. *Transfus Med Rev*. 2005;19:259–272.
7. Chiu ERW, Yuen KY, Lie AKW, et al. A prospective study of symptomatic bacteremia following platelet transfusion and of its management. *Transfusion*. 1994;34:960–964.
8. Sullivan MT, McCullough J, Schreiber GB, Wallace EL. Blood collection and transfusion in the United States in 1997. *Transfusion*. 2002;42:1253–1260.
9. Andreu G, Morel P, Forestier F, et al. Haemovigilance network in France: organization and analysis of immediate transfusion incident reports from 1984 to 1998. *Transfusion*. 2002;42:1356–1364.
10. The Trial to Reduce Alloimmunization to Platelets Study Group. Leukocyte reduction and ultraviolet B irradiation of platelets to prevent alloimmunization and refractoriness to platelet transfusions. *N Engl J Med*. 1997;337:1861–1869.
11. Andreu G, Dewailly J. Prevention of HLA alloimmunization by using leukocyte-depleted components. *Curr Stud Hematol Blood Transfus*. 1994;60:29–40.
12. Anderson NA, Gray S, Copplesstone JA, et al. A prospective randomized study of three types of platelet concentrates in patients with haematological malignancy: corrected platelet increments and frequency of non haemolytic febrile transfusion reactions. *Transfus Med*. 1997;7:33–39.
13. Enright H, Davis K, Gernsheimer T, et al. Factors influencing moderate to severe reactions to PLT transfusions: experience of the TRAP multicenter clinical trial. *Transfusion*. 2003;43:1545–1552.
14. Thaler M, Shamiss A, Orgad S, et al. The role of the blood from HLA-homozygous donors in fatal transfusion associated graft-versus-host disease after open heart surgery. *N Engl J Med*. 1989;321:25–28.
15. Wagner FF, Flegel WA. Transfusion-associated graft-versus host disease: risk due to homozygous HLA haplotypes. *Transfusion*. 1995;35:284–291.
16. Seghatchian J, de Sousa G. Pathogen-reduction systems for blood components: the current position and future trends. *Transfus Apher Sci*. 2006;35:189–196.
17. Pelletier JPR, Transue S, Snyder EL. Pathogen inactivation techniques. *Best Pract Res Clin Haematol*. 2006;19:205–242.
18. Moeremans K, Warie H, Annemans L. Assessment of the economic value of the INTERCEPT blood system in Belgium. *Transfus Med*. 2006;16:17–30.
19. Weusten JJ, Van Drimmelen HA, Lelie PN. Mathematic modeling of the risk of HBV, HCV, and HIV transmission by window-phase donations not detected by NAT. *Transfusion*. 2002;42:537–548.
20. Alter HJ, Stramer SL, Dodd RY. Emerging infectious diseases that threaten the blood supply. *Semin Hematol*. 2007;44:32–41.
21. Watson R. Europe witnesses first local transmission of chikungunya fever in Italy. *BMJ*. 2007;335:532–533.
22. Eder AF, Kennedy JM, Dy BA, et al. American Red Cross Regional Blood Centers. Bacterial screening of apheresis platelets and the residual risk of septic transfusion reactions: the American Red Cross experience (2004–2006). *Transfusion*. 2007;47:1134–1142.
23. Akahoshi M, Takahashi M, Masuda M, et al. A case of transfusion-associated graft-versus-host disease not prevented by white cell-reduction filters. *Transfusion*. 1992;32:169–172.
24. Visconti MR, Pennington J, Garner SF, et al. Assessment of removal of human cytomegalovirus from blood components by leukocyte depletion filters using real-time quantitative PCR. *Blood*. 2004;103:1137–1139.
25. Pennington J, Taylor GP, Sutherland J, et al. Persistence of HTLV-I in blood components after leukocyte depletion. *Blood*. 2002;100:677–681.
26. Roberts P. Resistance of vaccinia virus to inactivation by solvent detergent treatment of blood products. *Biologicals*. 2000;28:29–32.

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References

1. Kumar V, Lockerbie O, Keil SD, et al. Riboflavin and UV-light based pathogen reduction: extent and consequence of DNA damage at the molecular level. *Photochem Photobiol.* 2004;80:15–21.
2. World Health Organization. Thiamine, riboflavin, niacin, vitamin B6, pantothenic acid, and biotin. In: *Vitamin and Mineral Requirements in Human Nutrition*. 2nd ed. 2004:Section 9.3.5;172.
3. Scientific Committee for Food (SCF) (2000). Tolerable upper intake level of vitamin B2. (Opinion expressed on 22 November 2000.)
4. Rivlin SR. Riboflavin metabolism. *N Engl J Med.* 1970;283:463–472.
5. Kale H, Hanikumar P, Nair PM, Netrawali MS. Assessment of the genotoxic potential of riboflavin and lumiflavin. *Mutat Res.* 1992;298:9–16.
6. United States Food and Drug Administration (FDA) Center for Food Safety and Applied Nutrition (CFSAN) Priority-based Assessment of Food Additives (PAFA). EAFUS: A Food Additive Database. Available at: <http://vm.cfsan.fda.gov/~7Edms/efalus.html> Accessed November 2007.
7. Warburg VO, Christian W. Über das neue Oxydationsfenment. *Naturwissenschaften.* 1932;20:980–981.
8. Speck WT, Rosenkranz S, Rosenkranz HS. Further observations on the photo-oxidation of DNA in the presence of riboflavin. *Biochim Biophys Acta.* 1976;435:39–44.
9. Oster G, Bellin JS, Holmstrom B, et al. Photochemistry of riboflavin. *Experientia.* 1962;18:249–253.
10. Tsugita A, Okada Y, Uehara K. Photosensitized inactivation of ribonucleic acids in the presence of riboflavin. *Biochim Biophys Acta.* 1965;103:360–363.
11. Speck WT, Chen CC, Rosenkranz HS. In vitro studies of effects of light and riboflavin on DNA and HeLa cells. *Pediatr Res.* 1975;9:150–153.
12. Kuratomi K, Kobayashi Y. Studies on the interactions between DNA and flavins. *Biochim Biophys Acta.* 1977;476:207–217.
13. Korycka-Dahl M, Richardson T. Photodegradation of DNA with fluorescent light in the presence of riboflavin and photoprotection by flavin triplet-state quenchers. *Biochim Biophys Acta.* 1980;610:229–230.
14. Cairns WL, Metzler DE. Photochemical degradation of flavins. VI. A new photoproduct and its use in studying the photolytic mechanism. *J Am Chem Soc.* 1971;93:2772–2777.
15. Hardwick CC, Herivel TR, Hernandez SC, et al. Separation, identification and quantification of riboflavin and its photoproducts in blood products using high-performance liquid chromatography with fluorescence detection: a method to support pathogen reduction technology. *Photochem Photobiol.* 2004;80:609–615.
16. Sisson TR. Photodegradation of riboflavin in neonates. *Federation Proc.* 1987;46:1883–1885.
17. Olsen JH, Hertz H, Kjaer SK, et al. Childhood leukemia following phototherapy for neonatal hyperbilirubinemia (Denmark). *Cancer Causes Control.* 1996;7:411–414.
18. Goodrich RP, Edrich RA, Goodrich LL, et al. The antiviral and antibacterial properties of riboflavin and light: applications to blood safety and transfusion medicine. In: *Comprehensive Series in Photochemical & Photobiological Sciences.* 2006;6:83–113.
19. Martin CB, Willong E, Ruane P, et al. An action spectrum of the riboflavin-photosensitized inactivation of lambda phage. *Photochem Photobiol.* 2005;81:474–480.
20. Dardare N, Platz MS. Binding affinities of commonly employed sensitizers of viral inactivation. *Photochem Photobiol.* 2002;75:561–564.
21. Joshi PC. Ultraviolet radiation-induced photodegradation and 1O_2 , O_2^- production by riboflavin, lumichrome and lumiflavin. *Indian J Biochem Biophys.* 1989;26:186–189.
22. Data on file. Navigant Biotechnologies LLC.
23. Goodrich RP, Edrich RA, Li J, et al. The Mirasol™ PRT system for pathogen reduction of platelets and plasma: an overview of current status and future trends. *Transfus Apher Sci.* 2006;35:5–17.

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References

1. Scientific Committee for Food (SCF) (2000). Tolerable upper intake level of vitamin B2. (Opinion expressed on 22 November 2000.)
2. US Food and Drug Administration. Department of Health and Human Substances. Direct Food Substances Affirmed as Generally Recognized as Safe: Riboflavin. 21 CFR 184.1695, 2001.
3. World Health Organization. Thiamine, riboflavin, niacin, vitamin B6, pantothenic acid, and biotin. In: *Vitamin and Mineral Requirements in Human Nutrition*, 2nd ed. 2004;Section 9.3.5:172.
4. Martindale W, Sweetman SC, eds. *The Complete Drug Reference*. 33rd ed. London: Pharmaceutical Press; 2002:1386.
5. Unna K, Greslin JG. Studies on the toxicity and pharmacology of riboflavin. *J Pharmacol*. 1942;76:75–80.
6. Studer A, Zbinden G, Uehlinger E. Die pathologie der avitaminosen and hypervitaminosen. In: Buchner F, Lettner E, Roulet F, eds. *Handbuch der Allgemeinen Pathologie*. Berlin: Springer Verlag; 1962 (Band II):734–987.
7. Muñoz N, Hayashi M, Bang LJ, et al. Effect of riboflavin, retinol, and zinc on micronuclei of buccal mucosa and of esophagus: a randomized double-blind intervention study in China. *J Natl Cancer Inst*. 1987;79:687–691.
8. Hayashi M, Kishi M, Sotani S, et al. Micronucleus tests in mice on 39 food additives and eight miscellaneous chemicals. *Food Chem Toxicol*. 1988;26:487–500.
9. Sisson TR. Photodegradation of riboflavin in neonates. *Fed Proc*. 1987;46:1863–1865.
10. Gromisch DS, Lopez R, Cole HS, et al. Light (phototherapy)-induced riboflavin deficiency in the neonate. *J Pediatr*. 1977;90:118–122.
11. Speck WT, Rosenkranz HS. Phototherapy for neonatal hyperbilirubinemia – a potential environmental health hazard to newborn infants: a review. *Environmental Mutagenesis*. 1979;1:321–336.
12. Olsen JH, Hertz H, Kjaer SK, et al. Childhood leukemia following phototherapy for neonatal hyperbilirubinemia (Denmark). *Cancer Causes Control*. 1996;7:411–414.
13. Treadwell GE, Cairns WL, Metzler DE. Photochemical degradation of flavins. V. Chromatographic studies of the products of photolysis of riboflavin. *J Chromatogr*. 1968;35:376–388.
14. Cairns WL, Metzler DE. Photochemical degradation of flavins. VI. A new photoproduct and its use in studying the photolytic mechanism. *J Am Chem Soc*. 1971;93:2772–2777.
15. Oka M, McCormick DB. Urinary lumichrome-level catabolites of riboflavin are due to microbial and photochemical events and not rat tissue enzymatic cleavage of the ribityl chain. *J Nutr*. 1985;115:496–499.
16. Hardwick CC, Herivel TR, Hernandez SC, et al. Separation, identification and quantification of riboflavin and its photoproducts in blood products using high-performance liquid chromatography with fluorescence detection: a method to support pathogen reduction technology. *Photochem Photobiol*. 2004;80:609–615.
17. Kale H, Hanikumar P, Nair PM, et al. Assessment of the genotoxic potential of riboflavin and lumiflavin. A. Effect of metabolic enzymes. *Mutat Res*. 1992;293:9–16.
18. Huang R, Choe E, Min DB. Kinetics for singlet oxygen formation by riboflavin photosensitization and the reaction between riboflavin and singlet oxygen. *J Food Sci*. 2004;69:C726–C732.
19. Reddy et al. *Transfusion*. Manuscript in press, April 2008.
20. BIBRA. Riboflavin and its Derivatives: Toxicity Profile. 1990:1–7. Available at: <http://www.bibra-information.co.uk>. Accessed 28 November 2007.
21. Fouty B, Freyman F, Reyes R. Riboflavin to treat nucleoside analogue-induced lactic acidosis. *Lancet*. 1998;352:291–292.
22. Schoenen J, Jacquy J, Lenaerts M. Effectiveness of high-dose riboflavin in migraine prophylaxis: a randomized controlled trial. *Neurology*. 1998;50:466–470.
23. Luzzati R, Del Bravo P, Di Perti G, et al. Riboflavine and severe lactic acidosis. *Lancet*. 1999;353:901–902.
24. Expert Group on Vitamins and Minerals. Safe Upper Levels for Vitamins and Minerals. May 2003. Available at: <http://www.food.gov.uk/multimedia/pdfs/vitmin2003.pdf>. Accessed 23 November 2007.
25. Zempleni J, Galloway JR, McCormick DB. Pharmacokinetics of orally and intravenously administered riboflavin in healthy humans. *Am J Clin Nutr*. 1996;63:54–66.
26. International Organization for Standardization. ISO 10993-1:2003.

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References

1. FDA Guidance for Industry: OSA Viral safety evaluation of biotechnology products derived from cell lines of human or animal origin, September 1998.
2. Data on file at Navigant Biotechnologies LLC.
3. Goodrich RP, Keil S, Busch MP. Pathogen inactivation in blood components and the definition of adequacy in performance. Manuscript in preparation.
4. Committee for Proprietary Medicinal Products CPMP/ICH/295/95 – Guidance for Industry: OSA Viral safety evaluation of biotechnology products derived from cell lines of human or animal origin, March 4, 1997.
5. Slamer S. Project lead: Emerging Infectious Disease Project. American Association of Blood Banks. Personal communication, September 2006.
6. World Health Organization guidelines, WHO technical report, Series No. 924, 2004:61–162.
7. Ruane P, Edrich R, Gampp D, et al. Photochemical inactivation of selected viruses and bacteria in platelet concentrates using riboflavin and light. *Transfusion*. 2004;44:877–885.
8. Goodrich RP, Edrich RA, Li J, Seghatchian J. The Mirasol™ PRT system for pathogen reduction of platelets and plasma: an overview of current status and future trends. *Transfus Apher Sci*. 2006;35:5–17.
9. Piatak M Jr, Saag MS, Yang LC, et al. High levels of HIV-1 in plasma during all stages of infection determined by competitive PCR. *Science*. 1993;259:1749–1754.
10. Brown KE, Young NS, Alving BM, Barbosa LH. Parvovirus B19: implications for transfusion medicine. Summary of a workshop. *Transfusion*. 2001;41:130–135.
11. Syridou G, Skevaki C, Kafetzis DA. Intrauterine infection with parvovirus B19 and CMV: implications in early and late gestation fetal demise. *Expert Rev Anti Infect Ther*. 2005;3:661–661.
12. West Nile virus activity – United States, October 10–16, 2002, and update on West Nile virus infections in recipients of blood transfusions. *MMWR Morb Mortal Wkly Rep*. 2002;51:929–931.
13. Pealer LN, Marfin AA, Petersen LR, et al. Transmission of West Nile virus through blood transfusion in the United States in 2002. *N Engl J Med*. 2003;349:1236–1245.
14. Stramer SL, Fang CT, Foster GA, et al. West Nile virus among blood donors in the United States, 2003 and 2004. *N Engl J Med*. 2006;353:451–459.
15. Busch MP, Caglioti S, Robertson EF, et al. Screening the blood supply for West Nile virus RNA by nucleic acid amplification testing. *N Engl J Med*. 2005;353:460–467.
16. Macedo De Oliveira A, Beecham BD, Montgomery SP, et al. West Nile virus blood transfusion-related infection despite nucleic acid testing. *Transfusion*. 2004;44:1695–1699.
17. Kusov Y, Kanda T, Palmenberg A, et al. Silencing of hepatitis A virus infection by small interfering RNAs. *J Virol*. 2006;80:5599–5610.
18. Strauss RG. Leukocyte-reduction to prevent transfusion-transmitted cytomegalovirus infections. *Pediatr Transplant*. 1999;3(Suppl 1):19–22.
19. Varmakas EC. Is white blood cell reduction equivalent to antibody screening in preventing transmission of cytomegalovirus by transfusion? A review of the literature and meta-analysis. *Transfus Med Rev*. 2005;19:181–199.
20. Ness P, Braine H, King K, et al. Single-donor platelets reduce the risk of septic platelet transfusion reactions. *Transfusion*. 2001;41:857–861.
21. Engelriet CP, Reesink HW, Blajchman MA, et al. Bacterial contamination of blood components. *Vox Sang*. 2000;78:59–67.
22. Brecher M, Holland PV, Pineda AA, et al. Growth of bacteria in inoculated platelets: implications for bacterial detection and the extension of platelet storage. *Transfusion*. 2000;40:1308–1312.
23. Goodrich RP, Gilmour D, Howenga N, Keil SD. A laboratory comparison of Mirasol PRT treatment and culture of platelet products for addressing bacterial contamination concerns. Unpublished manuscript.
24. Foley MJ, Doherty C, Salamit A, et al. Re-testing platelet concentrations for bacterial contamination at day 4 prior to extended storage for 7 days. *Transfusion*. 2007;47(Suppl):36A(S95–040E).
25. Eder AF, Kennedy JM, Oy B, et al. Bacterial screening of apheresis platelets and the residual risk of septic transfusion reactions: the American Red Cross experience (2004–2006). *Transfusion*. 2007;47:1134–1142.
26. Blajchman MA, Beckers EAM, Dickmeiss E, et al. Bacterial detection of platelets: current problems and possible resolutions. *Transfus Med News*. 2005;19:259–272.
27. Leiby DA, Rentas FJ, Nelson KE, et al. Evidence of *Trypanosoma cruzi* infection (Chagas' disease) among patients undergoing cardiac surgery. *Circulation*. 2000;102:2978–2982.
28. Kelly DJ, Richards AL, Temenak J, et al. The past and present threat of rickettsial diseases to military medicine and international public health. *Clin Infect Dis*. 2002;34:5145–5169.
29. Young C, Losikoff P, Chawla A, et al. Transfusion-acquired *Trypanosoma cruzi* infection. *Transfusion*. 2007;47:540–544.
30. Dey A, Singh S. Transfusion transmitted leishmaniasis: a case report and review of literature. *Indian J Med Microbiol*. 2006;24:165–170.
31. Kitchen AD, Chiodini PL. Malaria and blood transfusion. *Vox Sang*. 2006;90:77–84.
32. Cardo LJ, Rentas FJ, Kelchum L, et al. Pathogen inactivation of *Leishmania donovani infantum* in plasma and platelet concentrates using riboflavin and ultraviolet light. *Vox Sang*. 2006;90:85–91.
33. Lippert L, Watson R, Doane S, et al. Inactivation of *P. falciparum* by riboflavin and light. *Vox Sang*. 2002;83:163.
34. Rentas F, Harman R, Gomez C, et al. Inactivation of *Orientia tsutsugamushi* in red blood cells, plasma, and platelets with riboflavin and light as demonstrated in an animal model. *Transfusion*. 2007;47:240–247.

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References

1. Schroeder ML. Transfusion-associated graft-versus-host disease. *Br J Haematol*. 2002;117:275-287.
2. Lee TH, Paglieroni T, Utter GH, et al. High-level long-term white blood cell microchimerism after transfusion of leukoreduced blood components to patients resuscitated after severe traumatic injury. *Transfusion*. 2006;45:1280-1290.
3. Utter GH, Reed WF, Lee TH, Busch MP. Transfusion-associated microchimerism. *Vox sang*. 2007;93:188-195.
4. Toy P, Popovsky MA, Abraham E, et al. Transfusion-related acute lung injury: definition and review. *Crit Care Med*. 2005;33:721-726.
5. Siliman CC. The two-event model of transfusion-related acute lung injury. *Crit Care Med*. 2006;34(5 Suppl):S124-S131.
6. Heddle NM, Klama LN, Griffith L, et al. A prospective study to identify the risk factors associated with acute reactions to platelet and red cell transfusions. *Transfusion*. 1993;33:794-797.
7. Yazer MH, Podlosky L, Clarke G, et al. The effect of prestorage WBC reduction on the rates of febrile nonhemolytic transfusion reactions to platelet concentrates and RBC. *Transfusion*. 2004;44:10-15.
8. Rogers MA, Blumberg N, Saint SK, et al. Allogeneic blood transfusions explain increased mortality in women after coronary artery bypass graft surgery. *Am Heart J*. 2006;152:1028-1034.
9. Council of Europe Expert Committee in Blood Transfusion Study Group on Pathogen Inactivation of Labile Blood Components. Pathogen inactivation of labile blood components. *Transfusion Med*. 2001;11:149-175.
10. Ibojic J, Greiss MA, Urbaniak SJ. Limited efficacy of universal leukodepletion in reducing the incidence of febrile nonhemolytic reactions in red cell transfusions. *Transfusion Med*. 2002;12:181-185.
11. Vamvakas EC. Is white blood cell reduction equivalent to antibody screening in preventing transmission of cytomegalovirus by transfusion? A review of the literature and meta-analysis. *Transfus Med Rev*. 2005;19:181-199.
12. Akahoshi M, Takahashi M, Masuda M, et al. A case of transfusion-associated graft-versus-host disease not prevented by white cell-reduction filters. *Transfusion*. 1992;32:169-172.
13. Fast LD, DiLeone G, Li J, et al. Functional inactivation of white blood cells by Mirasol treatment. *Transfusion*. 2006;46:642-648.
14. Fast LD, Marschner S, DiLeone G, et al. Inactivation of human WBCs in RBC products using the Mirasol® system for whole blood. In: Abstracts of the American Society of Hematology 49th Annual Meeting and Exposition; December 8-11, 2007; Atlanta, GA, USA.
15. Fast LD, DiLeone G, Cardarelli G, et al. Mirasol PRT treatment of donor white blood cells prevents the development of xenogeneic graft-versus-host disease in Rag2-/-gcl-/- double knockout mice. *Transfusion*. 2006;46:1553-1560.
16. Data on file. Navigant Biotechnologies LLC.
17. Sachs UJ, Hattar K, Weissmann N, et al. Antibody-induced neutrophil activation as a trigger for transfusion-related acute lung injury in an ex vivo rat lung model. *Blood*. 2006;107:1217-1219.
18. Ambruso DR, Thurman G, Marschner S, et al. Mirasol™ pathogen reduction of packed red blood cells (PRBCs) and apheresis platelet (Apl) concentrates decreases neutrophil priming activity generated during storage. In: Abstracts of the American Society of Hematology 49th Annual Meeting and Exposition; December 8-11, 2007; Atlanta, GA, USA.
19. Pelszynski MM, Moroff G, Luban NL, et al. Effect of gamma irradiation of red blood cell units on T-cell inactivation as assessed by limiting dilution analysis: implications for preventing transfusion-associated graft-versus-host disease. *Blood*. 1994;83:1683-1689.
20. Luban NL, Drothler D, Moroff G, Quinones R. Irradiation of platelet components: inhibition of lymphocyte proliferation assessed by limiting-dilution analysis. *Transfusion*. 2000;40:346-352.

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
References

1. AuBuchon JP, Herschel L, Roger J, et al. Efficacy of apheresis platelets treated with riboflavin and ultraviolet light for pathogen reduction. *Transfusion*. 2005;45:1335–1341.
2. Goodrich RP, Edrich PA, Li J, et al. The Mirasol™ PRT system for pathogen reduction of platelets and plasma: an overview of current status and future trends. *Transfus Apher Sci*. 2006;35:5–17.
3. Li J, Goodrich L, Hansen E, et al. Platelet glycolytic flux increases stimulated by ultraviolet-induced stress is not the direct cause of platelet morphology and activation changes: possible implications for the role of glucose in platelet storage. *Transfusion*. 2005;45:1750–1758.
4. Holme S, Moroff G, Murphy S. A multi-laboratory evaluation of in vitro platelet assays: the tests for extent of shape change and response to hypotonic shock. *Transfusion*. 1998;38:31–40.
5. Pieterz RNI, Engellriet CP, Reesink HW. Evaluation of stored platelets. *Vox Sang*. 2004;86:203–223.
6. Goodrich RP, Li J, Pieters H, et al. Correlation of in vitro platelet quality measurements with in vivo platelet viability in human subjects. *Vox Sang*. 2006;90:279–286.
7. Li J, De Korte D, Woolum MD, et al. Pathogen reduction of buffy coat platelet concentrates using riboflavin and light: comparisons with pathogen-reduction technology-treated apheresis platelet products. *Vox Sang*. 2004;87:82–90.
8. Perez-Pujol S, Tonda R, Lozano M, et al. Effects of a new pathogen-reduction technology (Mirasol PRT) on functional aspects of platelet concentrates. *Transfusion*. 2006;45:911–919.
9. Verhoeven AJ, Verhaar R, Gouwerok EG, et al. The mitochondrial membrane potential in human platelets: a sensitive parameter for platelet quality. *Transfusion*. 2005;45:82–89.
10. Li J, Lockerbie D, de Korte D, et al. Evaluation of platelet mitochondria integrity after treatment with Mirasol pathogen reduction technology. *Transfusion*. 2006;45:920–926.
11. Data on file. Navigant Biotechnologies LLC.

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References

1. Crookes RL, Hillier CD. Fresh frozen plasma and related products. In: Hillier CD, Silberstein LE, Ness PM, et al. eds. *Blood Banking and Transfusion Medicine Basic Principles and Practice*. 2nd ed. Philadelphia, PA: Churchill Livingstone Elsevier; 2007:259–269.
2. Cardigan R, Lawrie AS, Mackie IJ, Williamson LM. The quality of fresh-frozen plasma produced from whole blood stored at 4°C overnight. *Transfusion*. 2005;45:1342–1348.
3. Solheim BG, Seghatchian J. Update on pathogen reduction technology for therapeutic plasma: an overview. *Transfus Apher Sci*. 2006;35:83–90.
4. Council of Europe. *Guide to the preparation, use and quality assurance of blood components - 13th edition (2007)*. ISBN 978-92-871-6137-6. Available at: http://book.coe.int/EN/ficheouvrage.php?PAGEID=36&lang=EN&produit_aliasid=2151 Accessed November 2007.
5. UK Blood Transfusion and Tissue Transplantation Services. *Guidelines for the Blood Transfusion Service*. 7th ed. Norwich, England: TSO (The Stationery Office); 2005.
6. EU Directive. Commission Directive 2004/33/EC of 22 March 2004 implementing Directive 2002/96/EC of the European Parliament and of the Council as regards certain technical requirements for blood and blood components.
7. Data on File. Navigant Biotechnologies LLC.
8. Heiden M, Seitz R. Quality of therapeutic plasma – requirements for marketing authorization. *Thromb Res*. 2002;107:S47–S51.
9. European Pharmacopoeia 5th Edition. 01/2007:1646 Human Plasma.
10. Buytaert-Hoelen KA, Hendrix BK, Camp D, et al. Protein quality in Mirasol-treated fresh frozen plasma after one year in storage. Presented at the AABB Annual Meeting, Anaheim, USA, 20–23 October, 2007.
11. Ramasamy I, Farrugia A, Tran E, et al. Biological activity of von Willebrand Factor during the manufacture of therapeutic factor VIII concentrates as determined by the collagen binding assay. *Biologicals*. 1998;26:155–166.
12. Larrea L, Calabuig M, Rivera J, et al. Influence of riboflavin photoinactivation treatment on coagulation factors from fresh frozen plasma. Presented at the ISBT XVIIIth Regional Congress Asia, Hanoi, Vietnam, 10–13 November, 2007.
13. Kumar V, Motheral T, Luzniak G, et al. Mirasol™ Pathogen Reduction Technology: Riboflavin-based process retains functional activities of immunoglobulins. *Transfusion*. 2003;43(Suppl.):60A.
14. Fast LD, DiLeone G, Li J, et al. Functional inactivation of white blood cells by Mirasol treatment. *Transfusion*. 2006;46:642–648.

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References

1. Goodrich RP, Li J, Pieters H, et al. Correlation of in vitro platelet quality measurements with in vivo platelet viability in human subjects. *Vox Sang.* 2006;90:279-285.
2. AuBuchon JP, Herschel L, Roger J, et al. Efficacy of apheresis platelets treated with riboflavin and ultraviolet light for pathogen reduction. *Transfusion.* 2005;45:1335-1341.
3. Li J, Lockerbie O, de Korte D, et al. Evaluation of platelet mitochondria integrity after treatment with Mirasol pathogen reduction technology. *Transfusion.* 2005;45:920-926.
4. Fast LD, DiLeone G, Li J, Goodrich RP. Functional inactivation of white blood cells by Mirasol treatment. *Transfusion.* 2006;46:642-648.

† The Mirasol PRT System for Platelets is not available for sale in the United States.

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References

1. Goodman C, Chan S, Collins P, et al. Ensuring blood safety and availability in the US: technological advances, costs, and challenges to payment—final report. *Transfusion*. 2003;43:35–46S.
2. AuBuchon JP, Birkmeyer JD, Busch MP. Safety of the blood supply in the United States: opportunities and controversies. *Ann Intern Med*. 1997;127:904–909.
3. Postma MJ, van Hulst M, De Wolf JT, et al. Cost-effectiveness of pathogen inactivation for platelet transfusions in the Netherlands. *Transfus Med*. 2005;15:379–387.
4. Pathogen inactivation: making decisions about new technologies, consensus conference. Toronto, Canada, 29–30 March 2007.
5. Owens DK. Interpretation of cost-effectiveness analyses. *J Gen Intern Med*. 1998;13:716–717.
6. Jackson BR, Busch MP, Stramer SL, AuBuchon JP. The cost-effectiveness of NAT for HIV, HCV, and HBV in whole-blood donations. *Transfusion*. 2003;43:721–729.
7. Pereira A. Cost-effectiveness of transfusing virus-inactivated plasma instead of standard plasma. *Transfusion*. 1999;39:479–487.
8. Jackson BR, AuBuchon JP, Birkmeyer JD. Update of cost-effectiveness analysis for solvent-detergent-treated plasma. *JAMA*. 1999;282:329.
9. Riedler GF, Haycox AR, Duggan AK, Dakin HA. Cost-effectiveness of solvent/detergent-treated fresh-frozen plasma. *Vox Sang*. 2003;85:88–95.
10. Bell CF, Botteman MF, Gao X, et al. Cost-effectiveness of transfusion of platelet components prepared with pathogen inactivation treatment in the United States. *Clin Ther*. 2003;25:2464–2486.
11. Moeremans K, Warie H, Annemans L. Assessment of the economic value of the INTERCEPT blood system in Belgium. *Transfus Med*. 2006;16:17–30.