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References

1. FDA Guidance for Industry: OSA Viral safety evaluation of biotechnology products derived from cell lines of human or animal origin, September 1998.
2. Data on file at Navigant Biotechnologies LLC.
3. Goodrich RP, Keil S, Busch MP. Pathogen inactivation in blood components and the definition of adequacy in performance. Manuscript in preparation.
4. Committee for Proprietary Medicinal Products Guidelines CPMP/CH/295/95 – Guidance for Industry: OSA Viral safety evaluation of biotechnology products derived from cell lines of human or animal origin, March 4, 1997.
5. Stamer S. Project lead: Emerging Infectious Disease Project. American Association of Blood Banks. Personal communication, September 2006.
6. World Health Organization guidelines, WHO technical report, Series No. 924, 2004:61–162.
7. Ruane P, Edrich R, Gampp D, et al. Photochemical inactivation of selected viruses and bacteria in platelet concentrates using riboflavin and light. *Transfusion*. 2004;44:877–885.
8. Goodrich RP, Edrich RA, Li J, Seghatchian J. The Mirasol™ PRT system for pathogen reduction of platelets and plasma: an overview of current status and future trends. *Transfus Apher Sci*. 2006;35:5–17.
9. Piatak M Jr, Saag MS, Yang LC, et al. High levels of HIV-1 in plasma during all stages of infection determined by competitive PCR. *Science*. 1993;259:1749–1754.
10. Brown KE, Young NS, Alving BM, Barbosa LH. Parvovirus B19: implications for transfusion medicine. Summary of a workshop. *Transfusion*. 2001;41:130–135.
11. Syridou G, Skevaki C, Kafetzis DA. Intrauterine infection with parvovirus B19 and CMV: implications in early and late gestation fetal demise. *Expert Rev Anti Infect Ther*. 2005;3:651–661.
12. West Nile virus activity – United States, October 10–16, 2002, and update on West Nile virus infections in recipients of blood transfusions. *MMWR Morb Mortal Wkly Rep*. 2002;51:929–931.
13. Pealer LN, Marfin AA, Petersen LR, et al. Transmission of West Nile virus through blood transfusion in the United States in 2002. *N Engl J Med*. 2003;349:1236–1245.
14. Stramer SL, Fang CT, Foster GA, et al. West Nile virus among blood donors in the United States, 2003 and 2004. *N Engl J Med*. 2006;353:451–459.
15. Busch MP, Caglioti S, Robertson EF, et al. Screening the blood supply for West Nile virus RNA by nucleic acid amplification testing. *N Engl J Med*. 2005;353:460–467.
16. Macedo De Oliveira A, Beecham BD, Montgomery SP, et al. West Nile virus blood transfusion-related infection despite nucleic acid testing. *Transfusion*. 2004;44:1695–1699.
17. Kusov Y, Kanda T, Palmenberg A, et al. Silencing of hepatitis A virus infection by small interfering RNAs. *J Virol*. 2006;80:5599–5610.
18. Strauss RG. Leukocyte-reduction to prevent transfusion-transmitted cytomegalovirus infections. *Pediatr Transplant*. 1999;3(Suppl 1):19–22.
19. Vamvakas EC. Is white blood cell reduction equivalent to antibody screening in preventing transmission of cytomegalovirus by transfusion? A review of the literature and meta-analysis. *Transfus Med Rev*. 2005;19:181–199.
20. Ness P, Braine H, King K, et al. Single-donor platelets reduce the risk of septic platelet transfusion reactions. *Transfusion*. 2001;41:857–861.
21. Engelfriet CP, Reesink HW, Blajchman MA, et al. Bacterial contamination of blood components. *Vox Sang*. 2000;78:59–67.
22. Brecher M, Holland PV, Pineda AA, et al. Growth of bacteria in inoculated platelets: implications for bacterial detection and the extension of platelet storage. *Transfusion*. 2000;40:1308–1312.
23. Goodrich RP, Gilmour D, Howenga N, Keil SD. A laboratory comparison of Mirasol PRT treatment and culture of platelet products for addressing bacterial contamination concerns. Unpublished manuscript.
24. Foley MJ, Doherty C, Salami A, et al. Re-testing platelet concentrations for bacterial contamination at day 4 prior to extended storage for 7 days. *Transfusion*. 2007;47(Suppl):38A(S96–04CE).
25. Eder AF, Kennedy JM, Dy B, et al. Bacterial screening of apheresis platelets and the residual risk of septic transfusion reactions: the American Red Cross experience (2004–2006). *Transfusion*. 2007;47:1134–1142.
26. Blajchman MA, Beckers EAM, Dickmeiss E, et al. Bacterial detection of platelets: current problems and possible resolutions. *Transfus Med News*. 2006;19:259–272.
27. Leiby DA, Rentas FJ, Nelson KE, et al. Evidence of *Trypanosoma cruzi* infection (Chagas' disease) among patients undergoing cardiac surgery. *Circulation*. 2000;102:2978–2982.
28. Kelly DJ, Richards AL, Temenak J, et al. The past and present threat of rickettsial diseases to military medicine and international public health. *Clin Infect Dis*. 2002;34:S145–S169.
29. Young C, Losikoff P, Chawla A, et al. Transfusion-acquired *Trypanosoma cruzi* infection. *Transfusion*. 2007;47:540–544.
30. Dey A, Singh S. Transfusion transmitted leishmaniasis: a case report and review of literature. *Indian J Med Microbiol*. 2006;24:165–170.
31. Kitchen AD, Chiodini PL. Malaria and blood transfusion. *Vox Sang*. 2006;90:77–84.
32. Cardo LJ, Rentas FJ, Kelchum L, et al. Pathogen inactivation of *Leishmania donovani infantum* in plasma and platelet concentrates using riboflavin and ultraviolet light. *Vox Sang*. 2006;90:85–91.
33. Lippert L, Watson R, Doane S, et al. Inactivation of *P. falciparum* by riboflavin and light. *Vox Sang*. 2002;83:163.
34. Rentas F, Harman R, Gomez C, et al. Inactivation of *Orientia tsutsugamushi* in red blood cells, plasma, and platelets with riboflavin and light, as demonstrated in an animal model. *Transfusion*. 2007;47:240–247.

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References

1. Schroeder ML. Transfusion-associated graft-versus-host disease. *Br J Haematol.* 2002;117:275-287.
2. Lee TH, Paglieroni T, Utter GH, et al. High-level long-term white blood cell microchimerism after transfusion of leukoreduced blood components to patients resuscitated after severe traumatic injury. *Transfusion.* 2005;45:1280-1290.
3. Utter GH, Reed WF, Lee TH, Busch MP. Transfusion-associated microchimerism. *Vox sang.* 2007;93:188-195.
4. Toy P, Popovsky MA, Abraham E, et al. Transfusion-related acute lung injury: definition and review. *Crit Care Med.* 2005;33:721-726.
5. Silliman CC. The two-event model of transfusion-related acute lung injury. *Crit Care Med.* 2006;34(5 Suppl):S124-S131.
6. Heddle NM, Klama LN, Griffith L, et al. A prospective study to identify the risk factors associated with acute reactions to platelet and red cell transfusions. *Transfusion.* 1993;33:794-797.
7. Yazer MH, Podlosky L, Clarke G, et al. The effect of prestorage WBC reduction on the rates of febrile nonhemolytic transfusion reactions to platelet concentrates and RBC. *Transfusion.* 2004;44:10-15.
8. Rogers MA, Blumberg N, Saint SK, et al. Allogeneic blood transfusions explain increased mortality in women after coronary artery bypass graft surgery. *Am Heart J.* 2006;152:1028-1034.
9. Council of Europe Expert Committee in Blood transfusion Study Group on Pathogen Inactivation of Labile Blood Components. Pathogen inactivation of labile blood components. *Transfusion Med.* 2001;11:149-175.
10. Ibojje J, Greiss MA, Urbaniak SJ. Limited efficacy of universal leucodepletion in reducing the incidence of febrile nonhaemolytic reactions in red cell transfusions. *Transfusion Med.* 2002;12:181-185.
11. Vamvakas EC. Is white blood cell reduction equivalent to antibody screening in preventing transmission of cytomegalovirus by transfusion? A review of the literature and meta-analysis. *Transfus Med Rev.* 2005;19:181-199.
12. Akahoshi M, Takanashi M, Masuda M, et al. A case of transfusion-associated graft-versus-host disease not prevented by white cell-reduction filters. *Transfusion.* 1992;32:169-172.
13. Fast LD, DiLeone G, Li J, et al. Functional inactivation of white blood cells by Mirasol treatment. *Transfusion.* 2006;46:642-648.
14. Fast LD, Marschner S, DiLeone G, et al. Inactivation of human WBCs in RBC products using the Mirasol® system for whole blood. In: Abstracts of the American Society of Hematology 49th Annual Meeting and Exposition; December 8-11, 2007; Atlanta, GA, USA.
15. Fast LD, DiLeone G, Cardarelli G, et al. Mirasol PRT treatment of donor white blood cells prevents the development of xenogeneic graft-versus-host disease in Rag2-/-gcl-/- double knockout mice. *Transfusion.* 2006;46:1553-1560.
16. Data on file. Navigant Biotechnologies LLC.
17. Sachs UJ, Hattar K, Weissmann N, et al. Antibody-induced neutrophil activation as a trigger for transfusion-related acute lung injury in an ex vivo rat lung model. *Blood.* 2006;107:1217-1219.
18. Ambuso DR, Thurman G, Marschner S, et al. Mirasol™ pathogen reduction of packed red blood cells (PRBCs) and apheresis platelet (Apl) concentrates decreases neutrophil priming activity generated during storage. In: Abstracts of the American Society of Hematology 49th Annual Meeting and Exposition; December 8-11, 2007; Atlanta, GA, USA.
19. Pelszynski MM, Moroff G, Luban NL, et al. Effect of gamma irradiation of red blood cell units on T-cell inactivation as assessed by limiting dilution analysis: implications for preventing transfusion-associated graft-versus-host disease. *Blood.* 1994;83:1683-1689.
20. Luban NL, Drothler D, Moroff G, Quinones R. Irradiation of platelet components: inhibition of lymphocyte proliferation assessed by limiting-dilution analysis. *Transfusion.* 2000;40:348-352.

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