Table 5. RB and photoproducts in normal, untreated blood

Compound	Concentration (nM) average (range)
RB	23.9 (8.6–79.6)
LC	11.4 (0.0–75.3)
2KF	2.0 (0.9–3.4)
4KF	1.2 (0.5–2.4)
FMF	5.0 (1.7–11.5)

Hustad has reported a value of 6.9 nM (2.7–42.5 nM) for RB concentration in plasma (n = 63) (37). In a later report, Hustad reported the average RB concentration in the plasma of 124 senior citizens as 15.3 nM (5.4–28.4 nM) (25). Capo-chichi reported in 2000 that the median concentration of RB in 10 infants and 10 adolescents to be 20.9 nM (12.7–53.4) and 18.5 nM (8.2–57.8), respectively (25). We now report for n = 30, an average value of 23.9 nM (8.6–79.6 nM) (Fig. 3). The higher mean value observed in this study might be a reflection of the fact that no additional protein precipitation or extraction step was required in this procedure resulting in higher retention of starting levels of RB, or might be due to partial leakage of erythrocytes still present in the buffy coat platelet products.

Photolysis of RB during the MIRASOL process results in the formation of four photoproducts: 2KF, 4KF, FMF and LC. These four photoproducts were found to be present in apheresis platelets that had not undergone any photochemical treatment, although at a much lower concentration (Table 5, Fig. 4). It is also important to note in this context that the concentrations measured in the platelet products will be diluted by a factor of 16- to 20-fold upon infusion of the products into a patient's blood stream, thus lowering the difference between concentrations determined in these products and the levels naturally circulating in blood. The demonstration of the existence of these agents in naturally occurring blood products suggests that the introduction of a RB-based PRT process will not introduce new agents into the blood supply, which are not already present to some extent. The consequence of increased levels of these agents in blood products resulting from this process is being evaluated in separate toxicology studies, which include short-term and long-term exposure to each of the agents described and

characterized in this study (38–45). The presence of these agents in our blood, the ubiquitous nature of RB exposure, its presence in our diets and our ability to metabolize it and manage its inherent photochemistry suggests a low risk profile for this product. These features in combination with the pathogen reduction capacity of this system have the potential to greatly enhance the safety of blood products presently offered in routine clinical practice.

Acknowledgements—The authors acknowledge support in part by the following: Laboratories at Bonfils, Inc. (717 Yosemite Street, Denver, CO 80230) for the collection of apheresis platelets; Gambro BCT (10811 West Collins Avenue, Lakewood, CO 80215); Dr. Donald B. McCormick, Ph.D., Fuller E. Callaway Professor of Biochemistry, Department of Biochemistry, School of Medicine, Emory University, Atlanta, GA 30322-3050. Kristin E. Bjorgo, Navigant Biotechnologies, for help in preparing the manuscript.

REFERENCES

- Curran, J. W. and D. N. Lawrence (1984) AIDS associated with transfusion. *N. Engl. J. Med.* **310**, 69–75.
- Warburg, O. and W. Christian (1932) Über das neue oxydationsferment. *Naturwissenschaften* **20**, 980.
- Speck, W. T., S. Rosenkranz and H. S. Rosenkranz (1976) Further observations on the photooxidation of DNA in the presence of Riboflavin. *Biochim. Biophys. Acta* **435**, 39–44.
- Kuratomi, K. and Y. Kobayashi (1977) Studies on the interactions between DNA and flavins. *Biochim. Biophys. Acta* **476**, 207–217.
- Olsen, J. H., H. Hertz, S. K. Kjaer, A. Bautz, L. Møllekjær and J. D. Boice Jr. (1996) Childhood leukemia following phototherapy for neonatal hyperbilirubinemia (Denmark). *Cancer Causes Control* **7**, 411–414.
- Sisson, T. R., N. Kendall, E. Shaw and L. Kechavarz-Oliai (1972) Phototherapy of jaundice in the newborn infant. II. Effect of various light intensities. *J. Pediatr.* **81**, 35–38.
- Sisson, T. R. (1987) Photodegradation of Riboflavin in neonates. *Fed. Proc.* **46**, 1883–1885.
- Sisson, T. R., B. Slaven and P. B. Hamilton (1976) Effect of broad and narrow spectrum fluorescent light on blood constituents. *Birth Defects Orig. Artic. Ser.* **12**, 122–133.
- Boadi, W. Y., L. Thaire, D. Kerem and S. Yannai (1991) Effects of dietary supplementation with vitamin E, Riboflavin and selenium on central nervous system oxygen toxicity. *Pharmacol. Toxicol.* **68**, 77–82.
- Webster, R. P., M. D. Gawde and R. K. Bhattacharya (1996) Modulation of carcinogen-induced DNA damage and repair enzyme activity by dietary Riboflavin. *Cancer Lett.* **98**, 129–135.
- Bates, C. J. (1987) Human Riboflavin requirements, and metabolic consequences of deficiency in man and animals. *World Rev. Nutr. Diet* **50**, 215–265.
- Corbin, F. III. (2002) Pathogen inactivation of blood components: current status and introduction of an approach using Riboflavin as a photosensitizer. *Int. J. Hematol.* **76**(Suppl. 2), 253–257.
- Goodrich, R. and M. S. Platz (1997) The design and development of selective, photoactivated drugs for sterilization of blood products. *Drugs Future* **22**, 159–171.
- Ruane, P. H., R. Edrich, D. Gampp, S. D. Keil, R. L. Leonard and R. P. Goodrich (2004) Photochemical inactivation of selected viruses and bacteria in platelet concentrates using Riboflavin and light. *Transfusion* **44**, 877–885.
- Treadwell, G. E., W. L. Cairns and D. E. Metzler (1968) Photochemical degradation of flavins V. Chromatographic studies of the products of photolysis of Riboflavin. *J. Chromatogr.* **35**, 376–388.
- Cairns, W. L. and D. E. Metzler (1971) Photochemical degradation of flavins. VI. A new photoproduct and its use in studying photolytic mechanism. *J. Am. Chem. Soc.* **93**, 2772–2777.
- Oka, M. and D. B. McCormick (1985) Urinary lumichrome-level catabolites of Riboflavin are due to microbial and photochemical events and not rat tissue enzymatic cleavage of the ribityl side chain. *J. Nutr.* **115**, 496–499.
- Smith, E. C. and D. E. Metzler (1963) The photochemical degradation of Riboflavin. *J. Am. Chem. Soc.* **85**, 3285–3288.


Figure 4. HPLC photochemical profile generated during treatment of 50 μ M RB in apheresis platelet concentrates (MIRASOL-treated apheresis platelet concentrates) overlaying the photochemical profile generated by direct measurement of a sample of untreated apheresis platelets.

19. Song, P. S. and D. E. Metzler (1967) Photochemical degradation of flavins-IV. Studies of the anaerobic photolysis of Riboflavin. *Photochem. Photobiol.* **6**, 691–709.
20. Halwer, M. (1951) The photochemistry of Riboflavin and related compounds. *Dissertation* **73**, 4870–4874.
21. Speck, A. J., F. van Schaik, J. Schrijver and W. H. Schreurs (1982) Determination of the B2 vitamer flavin-adenine dinucleotide in whole blood by high-performance liquid chromatography with fluorometric detection. *J. Chromatogr.* **228**, 311–316.
22. Ichinose, N., K. Adachi and G. Schwedt (1985) Determination of B2 vitamers in serum of fish using high-performance liquid chromatography with fluorescence detection. *Analyst* **110**, 1505–1508.
23. Zemleni, J., J. R. Galloway and D. B. McCormick (1996) The identification and kinetics of 7 α -hydroxyriboflavin (7-hydroxymethylriboflavin) in blood plasma from humans following oral administration of Riboflavin supplements. *Int. J. Vitam. Nutr. Res.* **66**, 151–157.
24. Zemleni, J. (1995) Determination of Riboflavin and flavocoenzymes in human blood plasma by high-performance liquid chromatography. *Ann. Nutr. Metab.* **39**, 224–226.
25. Hustad, S., M. C. McKinley, H. McNulty, J. Schneede, J. J. Strain, J. M. Scott and P. M. Ueland (2002) Riboflavin, flavin mononucleotide, and flavin adenine dinucleotide in human plasma and erythrocytes at baseline and after low-dose Riboflavin supplementation. *Clin. Chem.* **48**, 1571–1577.
26. Lopez-Anaya, A. and M. Mayersohn (1987) Quantification of Riboflavin, Riboflavin 5'-phosphate and flavin adenine dinucleotide in plasma and urine by high-performance liquid chromatography. *J. Chromatogr.* **423**, 105–113.
27. Capo-chichi, C. D., J. L. Gueant, F. Feillet, F. Namour and M. Vidailhet (2000) Analysis of Riboflavin and Riboflavin cofactor levels in plasma by high-performance liquid chromatography. *J. Chromatogr. B Biomed. Sci. Appl.* **739**, 219–224.
28. Floridi, A., C. A. Palmerini, C. Fini, M. Pupita and F. Fidanza (1985) High performance liquid chromatographic analysis of flavin adenine dinucleotide in whole blood. *Int. J. Vitam. Nutr. Res.* **55**, 187–191.
29. Johansen, K. and P. O. Edlund (1990) Determination of water-soluble vitamins in blood and plasma by coupled-column liquid chromatography. *J. Chromatogr.* **506**, 471–479.
30. Roughead, Z. K. and D. B. McCormick (1990) Flavin composition of human milk. *Am. J. Clin. Nutr.* **52**, 854–857.
31. Roughead, Z. K. and D. B. McCormick (1990) Qualitative and quantitative assessment of flavins in cow's milk. *J. Nutr.* **120**, 382–388.
32. Oka, M. and D. B. McCormick (1985) Urinary lumichrome-level catabolites of Riboflavin are due to microbial and photochemical events and not rat tissue enzymatic cleavage of the ribityl chain. *J. Nutr.* **115**, 496–499.
33. Zemleni, J., G. Link and W. Kubler (1992) The transport of thiamine, Riboflavin and pyridoxal 5'-phosphate by human placenta. *Int. J. Vitam. Nutr. Res.* **62**, 165–172.
34. Zemleni, J., J. R. Galloway and D. B. McCormick (1996) Pharmacokinetics of orally and intravenously administered Riboflavin in healthy humans. *Am. J. Clin. Nutr.* **63**, 54–66.
35. Fall, H. F. and H. G. Petering (1956) Metabolite inhibitors. I. 6,7-dimethyl-9-formylmethylisoalloxazine, 6,7-deimethyl-9-(2'-hydroxyethyl)-isoalloxazine and derivatives. *J. Am. Chem. Soc.* **78**, 377–380.
36. [USP] United States Pharmacopeia (2003) MD Validation of Compendial Methods. In *USP-NF26*. (Edited by USP Compendial Committee), pp. 2439–2442. USP, Rockville, MD.
37. Hustad, S., P. M. Ueland and J. Schneede (1999) Quantification of Riboflavin, flavin mononucleotide, and flavin adenine dinucleotide in human plasma by capillary electrophoresis and laser-induced fluorescence detection. *Clin. Chem.* **45**, 862–868.
38. Piper, J., S. Murphy, R. Schuyler, R. Haynes and A. Wilson (2002) Assessment of the acute toxicity and genotoxicity risks associated with lumichrome, the primary photoproduct of Riboflavin. *Vox Sang.* **83**, 192.
39. Piper, J. T., S. E. Murphy, R. S. Schuyler, R. Haynes and A. Wilson (2001) Evaluation of the acute toxicity and genotoxicity risks associated with the Riboflavin photoproduct lumichrome. *Transfusion* **41**, 90S.
40. Piper, J., E. Hansen, M. Woolum, D. Gampp, L. Goodrich and G. Goodrich (2001) Evaluation of the genotoxicity and acute toxicity risks associated with a Riboflavin based pathogen inactivation process. *Int. J. Toxicol.* **20**, 404.
41. Piergiorgio, P. and A. Calatroni (1982) Hydrolysis of Riboflavin nucleotides in plasma monitored by high performance liquid chromatography. *J. Chromatogr.* **229**, 445–449.
42. Baeckert, P. A., H. L. Greene, I. Fritz, D. G. Oelberg and E. W. Adcock (1988) Vitamin concentrations in very low birth weight infants given vitamins intravenously in a lipid emulsion: measurement of vitamins A, D, and E and Riboflavin. *J. Pediatr.* **113**, 1057–1065.
43. Chastain, J. L. and D. B. McCormick (1987) Clarification and quantitation of primary (tissue) and secondary (microbial) catabolites of Riboflavin that are excreted in mammalian (rat) urine. *J. Nutr.* **117**, 468–475.
44. Chastain, J. L. and D. B. McCormick (1988) Characterization of a new flavin metabolite from human urine. *Biochim. Biophys. Acta* **967**, 131–134.
45. Chastain, J. L. and D. B. McCormick (1987) Flavin catabolites: identification and quantitation in human urine. *Am. J. Clin. Nutr.* **46**, 830–834.

Toxicity Testing of a Novel Riboflavin-Based Technology for Pathogen Reduction and White Blood Cell Inactivation

Heather L. Reddy, Anthony D. Dayan, Joy Cavagnaro, Shayne Gad, Junzhi Li, and Raymond P. Goodrich

The Mirasol PRT System (Gambro BCT, Lakewood, CO) for platelets and plasma uses riboflavin and UV light to reduce pathogens and inactivate white blood cells in donated blood products. An extensive toxicology program, developed in accordance with International Organisation for Standardisation (ISO) 10993 guidelines, was performed for the Mirasol PRT system. Test and control articles for most of the reported studies were treated (test) or untreated (control) blood products. For some studies, pure lumichrome (the major photoproduct of riboflavin) or photolyzed riboflavin solution was used. Systemic toxicity was evaluated with in vivo animal studies in the acute and subchronic

settings. Developmental toxicity was evaluated with an in vivo animal study. Genotoxicity and neoantigenicity were evaluated with in vitro and in vivo tests. Hemocompatibility and cytotoxicity were assessed with standard, in vitro assays. The pharmacokinetics, excretion, and tissue distribution of ^{14}C -riboflavin and its photoproducts was evaluated with an in vivo animal study. The possible presence of leachable or extractable compounds (from the disposable set) was evaluated with novel assays for measuring these compounds in blood. No treatment-related toxicity was observed in any of the studies.

© 2008 Elsevier Inc. All rights reserved.

THE COLLECTION, SEPARATION, and transfusion of red blood cells, platelets, plasma, and fractionated plasma components are essential medical practices, often required for the preservation of life and for the treatment of disease. Although the transfusion of these components is a vital therapy, transfusions are still associated with some risk for transmission of disease to the patient. Bacteria, viruses, and parasites are all potential sources of infection that have been transmitted by allogeneic transfusions.¹ Residual donor white blood cells (WBCs) can cause a series of severe immune responses in a transfusion recipient, resulting in the failure of the transfusion or other complications as a result.²⁻⁴ Prion-transmitted diseases are also a recent and growing concern owing to the potential transfusion-mediated transmission of variant Creutzfeldt-Jakob disease (vCJD).⁵⁻⁸

A number of safety measures have been taken to reduce the risk of disease transmission. They include increased diligence in preparing the donor's arm before donation and diversion of the first portion of blood from the donor to decrease the risks of bacterial contamination of the blood product.⁹ Decreasing the risk of bacteria-related morbidity and mortality in recipients was also the driver for the reduction in storage time for platelets from 7 to 5 days.¹⁰ Additional approaches to prevent disease transmission include the use of questionnaires to screen donors and thereby minimize the risks of collecting blood with viral or other diseases, the use of antibody testing and nucleic acid testing to detect viruses, and the use of culture techniques to detect bacteria.^{1,11} γ -Irradiation has

been used to reduce the possibility of transfusion-associated graft-vs-host disease (GVHD) in susceptible patient populations.¹² Leukoreduction has also been adopted in some parts of the world to reduce the possibility of prion transmission.¹³

Despite the use of the various approaches, as outlined, transfusion-transmitted infection from agents such as HIV, hepatitis B virus, hepatitis C virus, and bacteria continue.¹⁴⁻¹⁶ In part, these events arise from the inability of available tests to detect agents in blood at low levels, which are nonetheless infectious, or during a latent period. Risks related to transfusion-associated GVHD remain, primarily because γ -irradiation of blood components is not universally used. There are some infectious agents for which screening tests are not routinely used, notably Parvovirus B19; cytomegalovirus; and the causative agents of parasitic diseases such as Chagas disease, malaria, babesiosis, and leishmaniasis.¹⁷⁻¹⁹ The recent entry of West Nile virus (WNV) into the blood supply in North America is an example of the emergence of a new

From Navigant Biotechnologies, Inc, Lakewood, CO; Access BIO, Boyce, VA; and Gad Consulting Services, Cary, NC. A.D.D. is an Independent Consultant in London, UK

Current affiliation for J.L.: Thermogenesis, Sacramento, CA.

Address reprint requests to: Raymond Goodrich, PhD, Navigant Biotechnologies, Inc, 1215 Quail Street, Lakewood, CO 80215. Current affiliation for J.L.: Thermogenesis, Sacramento, CA.

E-mail: ray.goodrich@navbio.com

0887-7963/08/\$ - see front matter

© 2008 Elsevier Inc. All rights reserved.

doi: 10.1016/j.tmr.2007.12.003

Table 1. Reported Toxicity Findings After Riboflavin Administration

Study type	Species	Route	Results	Ref
Acute	Mouse	SC	1200 mg/kg: slight reduction in body temperature 2400 mg/kg: neuromuscular effects and mortalities	27
Acute	Mouse	IP	LD ₅₀ = 340 mg/kg LD ₀ = 750 mg/kg	28,29
Acute	Rat	IP	LD ₅₀ = 560 mg/kg	27
Acute	Rat	SC	LD ₀ > 5000 mg/kg	27
Acute	Rat	Oral	LD ₅₀ > 10000 mg/kg	27
Acute	Mouse	IV	LD ₅₀ > 50 mg/kg	30
Acute	Dog	Oral	NOAEL of 2000 mg/kg	27
Subacute 4 d	Mouse	SC	2/6 Mortalities at 1000 mg × kg ⁻¹ × d ⁻¹	28
Subchronic 4-6 wk	Rat	Oral	Increased liver enzymes at 2.0 mg × kg ⁻¹ × d ⁻¹	31
Chronic 22 mo	Rat	Oral	Mild liver cell hypertrophy at 1.5 mg × kg ⁻¹ × d ⁻¹	32
Chronic 140 d	Rat	Oral	50 mg × animal × d ⁻¹ , no adverse effects	27
Chronic 5 mo	Dog	Oral	25 mg × kg ⁻¹ × d ⁻¹ , no adverse effects	27
Developmental	Rat	Oral	During pregnancy and lactation NOAEL of 40 mg × kg ⁻¹ × d ⁻¹	33
Reproduction	Rat	Oral	3 Generations, NOAEL of 50 mg × animal × d ⁻¹	27
Carcinogenicity	Rat	Oral	1.5 mg × kg ⁻¹ × d ⁻¹ for 22 mo, no evidence of carcinogenicity	32
Micronucleus	Mouse	IP	Four doses at 1000 mg/kg, negative	28
Micronucleus	Mouse	Perioral	100 mg/kg, negative	34
Micronucleus (oral mucosa)	Human	Perioral	200 mg × kg ⁻¹ × wk ⁻¹ for 1 y, negative	35
Clastogenesis	CHL	Preincubation	Negative	36
Gene mutation	<i>E coli</i>	Preincubation	Negative	37
Ames test	TA1535, TA1537, TA1538, TA98, TA100 ± S9 activation	Preincubation	10 µmol/plate, negative results	38
Ames test	TA100, TA102, TA98, TA97A; ± S9 activation, cecal cell-free extract	Preincubation	25-100 µg/mL, negative results	39
Pharmacokinetics	Human	IV	11.6 mg/subject in bolus (242 mL) C _{max} = 1209 nmol/L Disposition t _{1/2} = 0.1006 h Renal clearance t _{1/2} = 0.7385 h	40

Abbreviations: SC, subcutaneous; IP, intraperitoneal; LD₅₀, lethal dose at which 50% of the test population survives; LD₀, lethal dose at which none of the test population survives; NOAEL, no observable adverse effect level; CHL, Chinese hamster lunch fibroblast cell line; C_{max}, maximum concentration in the bloodstream; t_{1/2}, time to reduction in concentration by 50%.

agent for which testing was not readily available. Since the appearance of WNV in the United States, several cases of transfusion-related WNV infections have occurred.²⁰ The low levels of WNV present in blood during a subject's subclinical viremia makes its detection difficult even with the most sensitive of the available nucleic acid testing technologies.²¹

Pathogen reduction via a chemical or photochemical process has also been used as a blood safety measure. Methylene blue-treated and solvent/detergent-treated plasmas have been used as transfusion products, such as fresh frozen plasma (FFP), and in the manufacture of fractionated and purified blood-derived products such as factor VIII

concentrate. These treatments are effective against enveloped viruses and yield only slight reductions in protein quality.^{22,23} Recently, a technology using a psoralen compound and UV light has been developed for the reduction of pathogens as well as the inactivation of leukocytes in platelets and in plasma.²⁴⁻²⁶

Similarly, the Mirasol PRT System (Gambro, BCT, Lakewood, CO) was developed to provide reduction of a wide variety of pathogens and inactivation of leukocytes. This device uses UV light and riboflavin to reduce pathogens. Riboflavin was chosen as a photochemical sensitizer because of its well-documented safety profile, as evidenced by


Fig. 1. Diagram of the steps involved with the use of the Mirasol PRT System. Platelets or plasma product in an illumination/storage bag is combined with riboflavin solution to achieve approximately 50 $\mu\text{mol/L}$ and then exposed to 6.2 J/mL UV light. The treated product is then ready for clinical use or blood bank storage. Riboflavin remains in the treated blood product, no removal device is used.

the reports referenced in Table 1. It is an essential nutrient for all living organisms, including humans. The novel nature of this pathogen reduction treatment (PRT) required comprehensive toxicological evaluation to prove that use of this technology to treat blood components would not introduce new hazards. The strategy, nature, and results of these toxicological studies are the subject of this review.

MIRASOL PRT SYSTEM—DESCRIPTION

The Mirasol PRT System consists of an illuminator, the illumination/storage bag, and riboflavin solution (Fig 1). The Illumination/Storage bag is made of biocompatible materials and is adapted from the container used to collect and store platelet products obtained with the Trima automated blood collection device (Gambro, BCT). The sterile riboflavin solution contains riboflavin (500 $\mu\text{mol/L}$) in a 0.9% sodium chloride (NaCl) solution with pH adjusted with hydrochloric acid to be in the range of 4.0 to 5.0. A volume of 35 ± 5 mL of this solution is added to a blood component product to produce a final concentration of approximately 50 $\mu\text{mol/L}$. The illuminator delivers the necessary dose of UV light (6.2 J/mL) to the contents of the illumination bag based on product volume (derived from weight) and measured flux rate. Robustness studies demonstrated that the system is suitable for use with initial platelet

or plasma product volumes ranging from 170 to 365 mL and platelet product cell concentrations ranging from 1100 to $2100 \times 10^3/\text{mL}$.

The Mirasol PRT System reduces the infectivity of pathogens by 3 combined mechanisms. The first is the direct damage of nucleic acids of pathogens with the UV light source used in this system. The second is the damage of pathogen nucleic acids, proteins, and membranes by reactive oxygen species generated when riboflavin absorbs light and interacts with dissolved oxygen in solution. The Mirasol process was designed to minimize the contributions of this mechanism because of its non-specific nature. The third is the damage of pathogen nucleic acid by the interaction of riboflavin with nucleic acids. Exposure of riboflavin to light when it is associated with or in proximity to nucleic acids leads to a Type I photochemical reaction that can induce oxidation of guanine residues and strand breaks. The result of this chemistry is irreversible damage to the nucleic acid.⁴¹

MIRASOL PRT PHOTOCHEMISTRY

Photolysis of riboflavin during treatment results in the formation of 4 major, quantifiable photo-products: 2'-ketoriboflavin, 4'-ketoriboflavin, formylmethylflavin, and lumichrome. Several minor components are present in negligible quantities and

Table 2. Concentration of Riboflavin and Photoproducts in Treated and Untreated Platelets⁴²

	Mirasol PRT platelets* (n = 30)			Untreated platelets (n = 30)
	Pretreatment [†] (μmol/L)	Posttreatment [‡] (μmol/L)	Day 5, treated [§] (μmol/L)	nmol/L
Riboflavin	48.33 ± 1.13 (47.1-50.4) [¶]	38.68 ± 0.68 (38.0-39.9)	39.65 ± 1.22 (37.0-41.3)	23.9 [#] (8.6-79.6)
Lumichrome	0.23 ± 0.08 (0.1-0.3)	2.90 ± 0.30 (2.5-3.0)	3.60 ± 0.30 (3.2-4.1)	11.4 (0.0-75.3)
2-Ketoflavin	ND	1.27 ± 0.17 (1.10-1.49)	0.71 ± 0.12 (0.59-0.90)	2.0 (0.9-3.4)
4-Ketoflavin	ND	0.53 ± 0.03 (0.48-0.58)	0.16 ± 0.04 (0.11-0.22)	1.2 (0.5-2.4)
Formylmethylflavin	ND	1.99 ± 0.24 (1.63-2.32)	1.72 ± 0.22 (1.42-2.07)	5.0 (1.7-11.5)

Abbreviation: ND, not detected.

* All studies for treated products were performed with and without γ-irradiation. No significant difference (Student *t* test, *P* < .05) was observed.

[†] Platelet product with Mirasol riboflavin solution, before UV light exposure.

[‡] Platelet product with Mirasol riboflavin solution, immediately after UV light exposure.

[§] Platelet product with Mirasol riboflavin solution and UV light exposure, stored for 5 days at room temperature with agitation.

^{||} Mean ± 1 SD.

[¶] Range of measured concentrations.

[#] Mean.

represent less than 1% of the total photoproduct peak areas. These minor components have been identified as isomers of the major components by analysis using high-pressure liquid chromatography followed by 2 stages of mass spectroscopy. Each of these 4 major photoproducts is a normal metabolite of riboflavin, and each has been detected in untreated apheresis platelets (Table 2), although at much lower concentrations than observed in Mirasol PRT-treated platelets.⁴² The concentrations measured in Mirasol PRT-treated blood products are diluted by a factor of 16- to 20-fold upon infusion of the products into a patient's blood stream (based on 250-300 mL infused into a subject of 5600 mL blood volume on average), thus reducing the concentrations toward the levels naturally circulating in the bloodstream. The existence of these photoproducts in freshly donated, nonilluminated blood products shows that a riboflavin-based PRT system will not introduce new chemicals into the blood supply.

These photoproducts have been examined in earlier studies of riboflavin photolysis in various solutions. After exposure to light, decreased reactivity of the solutions was demonstrated.⁴³⁻⁴⁵ The results clearly show that conversion of riboflavin from its native form into the photoproducts leads to decreased reactivity of the solution owing to decreases in the extinction coefficients of the resulting solution (photobleaching).

MIRASOL PRT SYSTEM—EFFICACY

Pathogen reduction with the Mirasol PRT System has been tested with a wide variety of enveloped and

nonenveloped viruses as well as gram-negative and gram-positive bacteria.⁴⁶ Prevention of GVHD due to WBC inactivation by the Mirasol PRT System has been studied both in vitro and in vivo.^{47,48} Parasite reduction studies have been conducted in collaboration with the Walter Reed Army Institute of Research and are the subject of separate publications.^{49,50} The results of these studies with bacteria, viruses, parasites, and WBCs demonstrate broad efficacy and capacity for the process in the reduction of pathogenic agents in blood components.

Blood components treated with riboflavin and UV light exhibit satisfactory quality, as assessed by in vitro measurements.^{46,51,52} Initial studies on protein quality after treatment demonstrated adequate retention of coagulation factors, coagulation inhibitors, and immunoglobulins, all at levels required for therapeutic efficacy.^{53,54} Some platelet quality indices were reduced in treated products in comparison with untreated controls but remained within acceptable transfusion limits.^{46,51,52} These studies were repeated after γ-irradiation of the products and showed no significant additional effect (C. Hardwick and R. Goodrich, unpublished data, 2003). Despite the differences in in vitro measurements, Mirasol PRT-treated platelets exhibited normal function in models of adhesion, aggregation, and thrombus formation.⁵⁵

Data from a clinical trial in healthy human subjects illustrated that Mirasol PRT-treated platelets exhibited decreased recovery and survival time in comparison with untreated platelets. Nevertheless, recovery values and cell quality indices for Mirasol PRT-treated platelets remained well

within the acceptable range for clinical utility of platelet products currently licensed and in routine use.^{56,57} Additional studies involving evaluation of the efficacy of treated platelet products in the clinical setting with thrombocytopenic patients are in progress.

THERAPEUTIC USE AND TOXICITY OF RIBOFLAVIN IN HUMANS

Riboflavin was chosen as the photosensitizing agent for this pathogen reduction process because it can be photoactivated to damage pathogens; it has a substantial history of safe clinical use; and it is both a naturally occurring vitamin (B₂) and an essential dietary nutrient. The latter features afforded it the advantage of a well-known and well-characterized safety profile on which to base its evaluation in transfusion medicine. The metabolic role of riboflavin has been well characterized, principally as a source of the essential coenzymes flavin adenine mono- and dinucleotides. The recommended daily intake of riboflavin is about 1.3 mg/d for the average adult⁵⁸⁻⁶⁰ and up to 1.6 mg/d for the lactating woman. Riboflavin is administered in medical practice by the oral, intramuscular, and intravenous (IV) routes to treat deficiency states, and it is widely taken orally as a dietary supplement. Recommended therapeutic doses by any route are up to 30 mg/d for an adult,⁶¹ but there are records of people taking up to 200 mg/d or more orally for 6 months with no adverse effect.^{59,62}

Riboflavin is very widely used as a permitted food coloring in the United States, where it has “generally regarded as safe” status,⁶³ and in Europe, where it has been approved by the Scientific Committee on Food.⁶⁴ The Joint Food Agriculture Organization/World Health Organization (FAO/WHO) Expert Committee on Food Additives has reviewed its use in food and it has been given a value for acceptable daily intake (ADI) of 0 to 0.5 mg · kg⁻¹ · d⁻¹, which is a maximum ADI of 35 mg/d for a subject of average weight (70 kg).⁵⁸

Riboflavin is given to newborn babies undergoing phototherapy for neonatal jaundice.^{65,66} Administration of supplemental riboflavin to these vulnerable patients has become standard practice because the treatment-accelerated clearance of bilirubin via phototherapy could also make them at least temporarily deficient in riboflavin owing to overlap of the light absorption spectra of riboflavin and bilirubin.⁶⁶⁻⁶⁹ To avoid the risk of the potential

riboflavin photochemistry-induced deficit, riboflavin (3 mg/kg perioral) is conventionally given to infants during the phototreatment regimen. Concern was once expressed about the toxic potential of riboflavin with direct phototherapy of newborns, based on theoretical consideration of its possible effect on DNA in the presence of light and on the results of experiments in simplified *in vitro* systems,^{70,71} but there have been no reports of adverse effects in the clinical setting. This includes results from large retrospective analyses, which directly evaluated potential tumorigenic consequences in more than 50 000 infants undergoing the therapy for more than 10 years after phototherapy.⁷⁰ This type of therapy continues in wide clinical use today.⁷²

There are a few claims of uncertain quality in individual subjects of “sensitization” and possible enhancement of epileptiform seizures after several months of high daily doses of riboflavin (eg, >5 mg · kg⁻¹ · d⁻¹ for >6 months), but those reports remain unsubstantiated.^{61,62} In no other reports of high oral and parenteral doses of riboflavin in humans have toxic effects been described,⁷³⁻⁷⁸ probably because of its physiologic nature and the rapid excretion of an excessive dose.⁴⁰

There is an abundance of references on the toxicity of riboflavin administered by various routes. These data support its use in total parenteral nutrition regimens. Table 1 summarizes representative published toxicity data for oral, subcutaneous, perioral, intraperitoneal, and IV exposures to riboflavin.

Although the safety of riboflavin has been extensively studied, there were no reports that directly supported its use in the Mirasol PRT System. Therefore, Navigant Biotechnologies, Lakewood, Colo, conducted a comprehensive preclinical safety evaluation program in support of the Mirasol PRT System, designed to investigate all potential sources of concern. The work was done in compliance with International Organisation for Standardisation (ISO) 10993.⁷⁹

METHODS

Test Article Definition

To establish the safety of Mirasol PRT-treated blood products, *in vivo* animal toxicity studies were performed with treated platelets (and in some instances with treated plasma), as well as with pure lumichrome, and with photolyzed riboflavin in

Table 3. Studies of Systemic Toxicity

	Acute toxicity			Subchronic toxicity
	Mirasol PRT–treated rat platelets	Mirasol PRT–treated dog platelets	Photolyzed riboflavin solution	Mirasol PRT–treated dog plasma
No. of animals/group	10 Male rats 10 Female rats	4 Male dogs 4 Female dogs	10 Male rats	6 Male dogs; 2 in recovery 6 Female dogs; 2 in recovery
Duration of study	14 d; dosing on d 1	14 d; dosing on day 1	9 d; dosing on day 1	13 wk with 4 wk for recovery group; dosing daily, 6 d/wk
Dose route	IV injection via tail vein	IV infusion via catheter	IV injection via tail vein	IV infusion via catheter
Group 1				
Dose material	Untreated rat platelets	Untreated dog platelets	Sterile 0.9% NaCl	Untreated dog plasma
Dose frequency	Three 2.0-mL injections, 2 h apart	Three 40-mL infusions, 2 h apart	Three 2.0-mL injections, 2 h apart	Daily 40-mL injections via catheter
Group 2				
Dose material	Mirasol PRT–treated rat platelets	Mirasol PRT–treated dog platelets	Photolyzed riboflavin solution	Mirasol PRT–treated dog plasma
Dose frequency	Three 2.0-mL injections, 2 h apart	Three 40-mL infusions, 2 h apart	Three 2.0-mL injections, 2 h apart	Daily 40-mL injections via catheter
Group 3				
Dose material	Mirasol PRT–treated, γ -irradiated rat platelets	Mirasol PRT–treated γ -irradiated dog platelets	Photolyzed riboflavin solution	NA
Dose frequency	2.0-mL injections	40-mL infusions	2.0-mL injections	
Group 4				
Dose material	Mirasol PRT–treated, γ -irradiated rat platelets	Mirasol PRT–treated γ -irradiated dog platelets	NA	NA
Dose frequency	Three 2.0-mL injections, 2 h apart	Three 40-mL infusions, 2 h apart		

Abbreviation: NA, not applicable.

the absence of platelets and plasma. The Mirasol PRT System parameters used to prepare test articles were always within the range verified for clinical use of the device.

Wherever possible, Mirasol PRT–treated platelets in autologous plasma were used as the test article so that any toxicity from the individual components and any interaction between them (storage bag, riboflavin, light, blood product) would be detected. The doses used in the *in vivo* studies were designed to maximize exposure while minimizing volume-associated effects and discomfort to the animals. In some cases, pure compounds were used to yield higher doses of the compounds of interest in a study (such as lumichrome and photolyzed riboflavin). In certain animal studies, such as embryo-fetal development in rats and subchronic repeated-dose test in dogs, the volume of animal platelets required to prepare test and control articles for the study would have been very large and would have been required over many weeks. To obtain a consistent test article in as humane a fashion as possible for those studies, species-specific plasma was used rather

than platelets. The photochemistry of riboflavin yields equivalent photoproduct profiles in plasma products and in platelet products (which consist mainly of plasma). The absence of platelets eliminates the possibility of detecting toxic alterations to the platelet surface; however, that issue was addressed in the neoantigenicity and ^{14}C -riboflavin binding studies. Therefore, when necessary and appropriate, plasma was regarded as an appropriate substitute for platelets and was used as the test article. Treatment of plasma with the Mirasol PRT System is identical to the treatment of platelets.

Sample Size Determination, Statistics, and Data Analysis

Sample size determination for studies was based on historical experience or published requirements where appropriate. In all cases, data were maintained in validated databases allowing archiving and analysis of data. All studies, assays, and data analysis were performed at accredited contract testing laboratories in accordance with Good Laboratory Practice (GLP) requirements, except

where noted. Comparisons of data between study groups were analyzed using Student *t* test, Wilcoxon distributions and other suitable analyses as appropriate for the particular study data being analyzed.

Systemic Toxicity Studies

Acute toxicity. Two acute toxicity studies, 1 in rats and 1 in dogs, evaluated the acute toxicity of Mirasol PRT–treated, species-specific platelet products, with or without γ -irradiation. In those studies, both males and females were used, to reveal any sex-dependent effects. In addition, an acute toxicity study of photolyzed riboflavin in saline was performed to test the toxicity of higher doses of riboflavin and its photoproducts than are attainable in a blood product. The dose volume used in all rat studies (2.0 mL) corresponded to a range of doses per body weight (mean dose of 9.0 ± 0.7 mL/kg for males and 11.5 ± 1.3 mL/kg for females in the study of treated rat platelets and 9.4 ± 0.3 mL/kg for the male rats in the study of photolyzed riboflavin). The dose volume used in the acute dog studies (40 mL) yielded mean doses per body weight of 3.9 ± 0.3 mL/kg for males and 5.1 ± 0.4 mL/kg for females. Table 3 provides a summary of the number of animals and of the dose frequency, route, and material for the acute toxicity studies. Each study was performed in accordance with ISO 10993^{80,81} and International Conference on Harmonization guidelines.⁸² These studies evaluated any immediate, systemic toxicity due to administration of Mirasol PRT–treated platelet products by observations of mortality, clinical signs of toxicity, hematology, clinical chemistry, food consumption, general behavior, and weight gain. Complete necropsies were conducted for all animals, and macroscopic observations were made. In the study performed with dogs, electrocardiographic examinations were conducted before randomization, at the conclusion of dosing, and before necropsy. The studies provided observations of effects on major physiologic systems, such as functioning of the nervous, cardiovascular, and respiratory systems, to exclude any pharmacologic effects of test article administration.

Subchronic toxicity. The repeated dose study for subchronic toxicity was done in the dog to support more prolonged and repeated administration of Mirasol PRT–treated platelets to patients, some of whom may require intermittent transfusions over several weeks. Parameters evaluated in the repeated dose toxicity study included morbidity

and mortality, outward signs of toxicity (behavioral and physical examinations), clinical chemistry of blood and urine, hematologic indices, body weights and food consumption, and pathologic and histopathologic evaluation of tissues. Test and control article administration is described in Table 3. The daily dose (40 mL) was equivalent to the human dose; the range of doses relative to body weight was 4.7 ± 0.3 mL/kg in the first week of dosing and 4.1 ± 0.3 mL/kg in the last week of dosing for the males and 5.6 ± 0.5 mL/kg in the first week of dosing and 4.8 ± 0.4 mL/kg in the last week of dosing for the females. On days 0 and 86 of the study, blood samples were removed from the animals to measure riboflavin and lumichrome and thereby assess any accumulation of riboflavin or photoproduct during the course of this experiment. Two animals per sex per group were observed for a 4-week recovery period after the 13-week dosing regimen. At the end of the terminal and recovery periods, complete necropsies were performed. Organ weights were recorded, and select tissues were microscopically examined for the non-recovery group animals.

Reproductive Toxicity Study

Embryo-fetal development in rats exposed to Mirasol PRT–treated products. The developmental toxicity potential of rat plasma treated with the Mirasol PRT System was evaluated in female Sprague-Dawley rats (mean gestation interval = 22 days). Time-mated, pregnant animals in the test and control groups ($n = 25$ for control, $n = 24$ for test group) were given IV PRT-treated or untreated products, respectively, at a dose of 2.0 mL per animal (7.7 ± 0.4 mL/kg on the first day of dosing, 6.0 ± 0.3 mL/kg on the last) on gestation days 6 through 17. The animals were observed over the course of the study for clinical signs, body weight, and food consumption. On gestation day 20, the animals were necropsied, and the ovaries and uterus were investigated for the total number of implantations, total number of corpora lutea, early and late resorptions, and viable and nonviable fetuses. Individual body weights and sex of the fetuses were noted. Gravid uterine weights were recorded, and adjusted body weight changes were calculated. All fetuses were given an external examination and processed for visceral or skeletal examination, and malformations and developmental variations were recorded.

Table 4. Mammalian Erythrocyte Micronucleus Test

Dose (20 mL/kg)	No. of mice/sex dosed	No. of mice/sex used for bone marrow collection after dose administration	
		24 h	48 h
Vehicle control: 0.9% NaCl for injection	10	5	5
Negative control: untreated human platelets	10	5	5
Test article: Mirasol PRT-treated platelets			
Low dose (diluted in saline 1:3)	5	5	0
Middle dose (diluted in saline 1:1)	5	5	0
High dose (undiluted)	10	5	5
Positive control: cyclophosphamide (50 mg/kg)	5	5	0

Genotoxicity Studies

Three different tests were performed to evaluate the genotoxicity of Mirasol PRT-treated human platelet products: an in vitro test for gene mutations in bacteria, an in vitro test for clastogenicity in mammalian cells, and an in vivo bone marrow micronucleus test.

Bacterial reverse mutation assay (Ames test).

Mirasol PRT-treated, day 5 platelets (with and without γ -irradiation) and control, untreated day 5 platelets were tested using *Salmonella typhimurium* tester strains TA98, TA100, TA1535, and TA1537 and *Escherichia coli* tester strain WP2 uvrA in the presence and absence of rat-liver derived S9, prepared by the contract research organization. The assay was performed in 2 phases, using the plate incorporation method and normal saline as the diluent. The first phase established the dose range for the confirmatory mutagenicity assay and provided a preliminary mutagenicity evaluation. In the first phase, the dose levels tested were 0.33, 1.0, 3.3, 10, 33, 100, 333, and 1000 μ L per plate. The second phase, the confirmatory mutagenicity assay, was used to evaluate and confirm the mutagenic potential of the test article with dose levels of 3.3, 10, 33, 100, 333, and 1000 μ L per plate.^{83,84}

The in vitro test for gene mutations in bacteria (Ames test) was also performed with lumichrome, the major product of riboflavin photolysis. The in vitro mutagenicity of lumichrome was evaluated in the dark to test its mutagenic potential in the absence of light activation. The maximum concentration tested was at the limit of solubility of lumichrome.

Chromosomal aberration in Chinese hamster ovary cells. Mirasol PRT-treated platelet products and control, untreated platelet products were tested to assess potential clastogenicity (chromosome breakage) in vitro with the mammalian

chromosome aberration test. The assay was performed using Chinese hamster ovary (CHO) cells.⁸⁵⁻⁸⁷ Based on the findings from the preliminary toxicity assay, the doses chosen for the chromosome aberration assay using the test article (Mirasol PRT-treated platelets) and negative control (untreated platelets) ranged from 0.005 to 0.1 mL/mL for all 3 exposure groups. The test and control articles were soluble in saline at a concentration of 1 mL/mL, the maximum concentration used. Both test and control article were tested in the absence and the presence of an Aroclor-induced S9 activation system.

The CHO cells were exposed to test or control article for 4 and 20 hours in the nonactivated test system and for 4 hours in the S9 activated test system. All cells were harvested 20 hours after treatment initiation. Selection of doses for microscopic analysis was based on toxicity of the test article (the lowest dose with at least 50% reduction in cell growth and 2 lower doses) in all harvests.

Mammalian erythrocyte micronucleus test.

Mirasol PRT-treated platelets and control, untreated platelets were tested to assess clastogenicity in vivo with the mammalian erythrocyte micronucleus test.⁸⁸⁻⁹⁰ An increase in the frequency of micronucleated polychromatic erythrocytes is expected in the bone marrow of animals exposed to genotoxic compounds. The assay was performed in 2 phases. In the pilot toxicity study, neat test article was administered intraperitoneally to 5 male and 5 female mice at a volume of 20 mL/kg (the highest optimal single-dose volume). No mortality was observed during the course of the study. All mice were observed for clinical signs of toxicity after dose administration and daily for 3 days afterward. Body weights were recorded before dose administration, 1 day after, and 3 days after dose administration. Results from the pilot study were

used to assess toxicity of the test article and to set dose levels for the definitive study.

In the definitive micronucleus study, 1 intraperitoneal injection was administered to the male and female mice in each group. The animals were observed for clinical signs of toxicity after dose administration. Table 4 describes the study groups, the number of animals per group, and test or control articles given to each group, as well as the sacrifice times. Immediately after sacrifice, bone marrow was aspirated from the femurs into a syringe. Two slides of bone marrow suspension were prepared and fixed for each mouse. The slides were analyzed for the presence of micronuclei in bone marrow cells, polychromatic erythrocytes, and normochromatic erythrocytes.

Neoantigenicity Studies

The potential for Mirasol PRT-treated platelets to exhibit neoantigenicity (the production of new antigenic determinants due to the treatment) was evaluated with 2 *in vitro* and 1 *in vivo* test. Because the photoproducts associated with the Mirasol PRT System are the same as metabolites, catabolites, and photoproducts of riboflavin normally present in humans under routine conditions, there was no reason to consider that they might themselves affect the immune system, and there were no indications of any such effect in the toxicity tests.

¹⁴C-riboflavin binding evaluation. The potential for neoantigen formation was first assessed *in vitro* by measuring the binding of radiolabeled riboflavin to platelets and plasma proteins after Mirasol PRT treatment and comparing the results with those obtained with untreated controls. This method assesses binding of a radiolabeled compound to proteins and cellular membranes. The presence of covalently bound compounds on cellular membranes or proteins indicates that chemical modification has occurred. This modification could be seen as foreign in the bloodstream and result in an immune response. Immune responses to antibiotics, bound to the surface of red blood cells, have been observed.⁹¹ Quinacrine mustard, which binds to proteins and blood cells, also results in the production of antibodies.⁹²⁻⁹⁴ If Mirasol PRT treatment were to lead to chemical modification of blood constituents *in vitro*, the *in vivo* generation of an immune response might be possible. Prior studies have identified the potential for binding of

riboflavin to albumin.⁹⁵ Reports of additional studies, carried out in nonphysiologic conditions, including the use of artificial media preparations, anaerobic conditions, and nonphysiologic pH values, claimed to demonstrate the formation of riboflavin-protein adducts via tryptophan residues in the protein.⁹⁵ This work also extended further to show the formation of possible antibodies against such complexes.⁹⁷ Other studies with human plasma samples were unable to reproduce these results and instead showed riboflavin interactions to occur primarily with immunoglobulins, which act as natural carriers for the compound *in vivo* and are a potential source of claims for antibody generation.^{98,99} The studies of ¹⁴C-riboflavin binding before and after photolysis were thus undertaken to assess the potential for adduct formation under the specific conditions used for Mirasol treatment of blood products.

To determine the distribution of riboflavin and its photoproducts after Mirasol PRT treatment, a known amount of radiolabeled riboflavin (¹⁴C) was added to unlabeled riboflavin solution, and the resulting solution was used to treat platelet products from 6 different donors with the Mirasol PRT System. Preillumination samples from each platelet product were used as controls, and test articles were postillumination samples. Cell concentrations, protein-associated radioactivity, and protein concentrations were determined for all samples. The method was sensitive enough to determine differences of less than 1% in binding of ¹⁴C-labeled compound, at micromolar concentrations per milligram of protein, if such effects were occurring. Given prior reports of the levels and extent of photochemical modification of proteins under specific conditions of treatment, if such modifications were occurring with Mirasol treatment, this method would be able to detect it.^{96,97}

Evaluation of nonimmunologic alterations. To assess the potential for nonimmunologic alterations to cells that could subsequently induce an immunologic response, another test, the Capture-P assay, (Immucor, Inc, Norcross, GA) was used to evaluate the binding of immunoglobulin G to Mirasol PRT-treated platelets, and to untreated platelets, after exposure to allogeneic and autologous plasma. Platelet products, and the associated plasma, from 2 different donors were used to provide the test and control articles for this study. Platelets were tested with the Capture-P assay after 1 and 5 days of

storage. This study was non-GLP but was well-documented, performed in a respected clinical laboratory, and included appropriate controls.

Neoantigenicity evaluated with the Ouchterlony assay. The formation of new antigens was evaluated in vivo by comparing the immune response between animals given control (untreated) human platelets and animals given Mirasol PRT-treated human platelets. Rabbits were immunized with control article or with 1 of 2 test articles (Mirasol PRT-treated, day 5 human platelets, with or without γ -irradiation). The immunizations occurred on days 0, 21, and 42 of the in-life portion of the study. On day 49, sera from the immunized animals in the test and control groups were collected and tested for antibody production. Antibodies were collected and used in the Ouchterlony assay to incubate with either control or test article. Lines of immunoprecipitation were examined to determine reactions of identity, partial identity, or nonidentity. Partial identity immunoprecipitation lines would indicate neoantigenicity if the test article contained new antigens arising as a result of the treatment. Sensitivity of the assay was evaluated with chicken albumin and serum from rabbits immunized with that compound. The chicken albumin could be detected by weak precipitin lines at 0.15 $\mu\text{g}/40 \mu\text{g}$ plasma protein and at 0.15 $\mu\text{g}/200 \mu\text{g}$ plasma protein, with antibody present at 10, 15, 20, or 25 μg per well in the assay plate.

Cytotoxicity Studies

Agar diffusion test. Mirasol PRT-treated human platelets and control, untreated human platelets were tested to assess cytotoxicity with the agar diffusion assay. The indicator cells for this assay were mouse fibroblast L929 cells, which are classically used for cytotoxicity studies because they demonstrate sensitivity to leachable cytotoxic compounds.¹⁰⁰

Cultures of L929 cells were used in the assay after replacement of the liquid medium with a serum-supplemented medium/agar mixture that is stained with a vital dye, neutral red. The culture was protected from light for the duration of the assay to prevent cell damage elicited by photoactivation of the stain. The test article and control article were applied directly to a filter paper disc (surface area 100 mm²) at a volume of 100 μL and placed on the surface of the agar. All plates were incubated for 48 hours, at 37°C \pm 1°C, in a humidified

atmosphere containing 5% \pm 1% carbon dioxide. The negative control article (negative control plastic) and the positive control article (natural rubber) were prepared and exposed similarly to the test and control articles. The extent of decolorization was evaluated at times 0, 24, and 48 hours.

Minimal essential medium elution study. The minimal essential medium elution study¹⁰⁰ was conducted to evaluate the potential cytotoxicity of Mirasol PRT-treated human platelets. The test was designed for the evaluation of test article extracts with assessment of biologic reactivity of a mammalian cell culture (L929). Aliquots of test article (Mirasol PRT-treated human platelets) and control article (untreated human platelets) were centrifuged, after which the resulting platelet pellets were mixed with the extraction medium (minimal essential medium) and incubated at 37°C \pm 1°C for 24 \pm 2 hours. After extraction, the medium was removed and applied to an L929 cell culture. The L929 cultures were incubated for 48 hours, at 37°C \pm 1°C, and then evaluated for the response of the cell monolayer. Biologic reactivity (cellular degeneration and malformation) was rated on a scale of 0 to 4.

In addition, the primary product of riboflavin photolysis, lumichrome, was tested for in vitro cytotoxicity by direct exposure to the indicator. The lumichrome solution, at the maximum possible lumichrome concentration, was exposed to the indicator cells.

Hemocompatibility Studies

Hemolysis test. Mirasol PRT-treated human platelets and control, untreated human platelets were assessed for hemolytic activity in direct contact with human blood.¹⁰¹⁻¹⁰³ Fresh, whole human blood (type-matched to be compatible with test and control articles) was collected into EDTA-coated Vacutainer tubes. The positive control article solution was obtained by adding 10 mL of sterile water for injection per vial. The test and control article solutions were obtained by adding 5 g of material to 10 mL of sterile 0.9% NaCl. All vials were incubated in a 37°C \pm 2°C water bath for 30 minutes, followed by addition of 0.2 mL of diluted human blood and incubation in a 37°C \pm 2°C water bath for 60 minutes. After incubation, the vials were centrifuged for 5 minutes at approximately 1500g, and the absorbance of each supernatant was determined against a NaCl blank

Table 5. Summary of Samples for Pharmacokinetic Study

Group no.	Urine*	Feces*	Cage rinse	Cage wash and wipe	Blood†	Tissues‡
1 (Intact)	Predose, 0-12 h, 12-24 h, 24-h intervals to 144 h	Predose, 24-h intervals to 144 h	24-h intervals to 120 h	144 h	144 h	144 h
2 (Intact)	NA	NA	NA	NA	1 h, 8 h, 24 h, 48 h, 96 h	1 h, 8 h, 24 h, 48 h, 96 h

* Urine and feces collected on wet ice. Urine assayed directly by liquid scintillation chromatography (LSC). Remaining urine saved for analysis of metabolites.

† Blood collected from all animals in groups 1 and 2 predose and at designated time points. Prep and assay plasma by LSC. Remaining plasma saved for analysis of metabolites.

‡ Tissues collected include the following: group 1—residual carcass, pelt (hair and skin with underlying fat layer), and tail skin; group 2—small intestine, large intestine, spleen, kidneys, liver, lymphatics (lymph nodes, includes splenic, inguinal, and popliteal), and femur bone marrow (right and left).

at 540 to 545 nm. Test and control articles were tested in triplicate.

Platelet function studies. Platelet function and the thrombogenicity of Mirasol PRT-treated platelets were directly evaluated in a series of non-GLP experiments. The hemostatic activity and thrombogenicity of treated platelets were measured in an ex vivo model in which Mirasol PRT-treated human platelets were mixed with whole blood previously depleted of platelets. The details of the experimental system are described by Perez-Pujol et al.⁵⁵ In brief, 8 Trima apheresis platelet collections were performed to collect more than 340 mL at a concentration between 1180 and 2160 × 10⁶ cells/mL. Half of the product volume was treated with the Mirasol PRT System. The other half was used as the paired control. Treated and control platelets were stored in a platelet storage incubator at 22°C ± 2°C for 5 days. Samples were tested using flow cytometric and perfusion methods on days 0, 3, and 5 of storage.⁵⁵

The flow cytometric analysis used commercially available monoclonal antibodies, tagged with fluorescein or phycoerythrin, to detect glycoprotein (GP) IIb-IIIa, GPIb, GPIV, P-selectin, the 53-kd lysosomal membrane protein, and coagulation factor V. The exposure of aminophospholipids on the outer leaflet of the platelet membrane was also quantified using annexin V. The binding of fibrinogen and von Willebrand factor was detected with the corresponding polyclonal antibodies, and nonspecific membrane immunofluorescence was evaluated.

Blood for perfusions was obtained from healthy volunteers and depleted of platelets and leukocytes by filtration. Samples taken from the test and

control articles at different days were incorporated into the thrombocytopenic blood at volumes aimed to increase the platelet count in the perfusate to 150 × 10⁹/L with a final volume of 20 mL. Perfusions were performed at 37°C in annular chambers with enzymatically denuded New Zealand rabbit aorta segments. Blood was recirculated through the chamber for 10 minutes at a shear rate of 800 s⁻¹. At the end of the perfusion, the segments were rinsed with 20 mL of PBS, removed from the chamber, and sliced off and washed with a fixing buffer. Platelets interacting with the subendothelium were evaluated and classified as either adhesion, contact platelets that are attached but not spread and groups of platelets that form aggregates of less than 5 µm in height, or thrombus, platelets that form aggregates of 5 µm or more in height.

Coagulation proteins and complement pathway. Mirasol PRT-treated human platelets and control, untreated human platelets were stored at 22°C for 5 days, after which samples were removed and centrifuged to isolate the plasma fraction. The plasma fraction was analyzed with function-based assays for key protein components of each hemostatic pathway—the contact pathway, the thrombin pathway, the fibrinolytic pathway, the complement pathway, and the coagulation pathway. Factor XI (FXI), factor XII (FXII), high-molecular-weight kininogen (HMWK), prekallikrein (PK), and activated factor XII (FXIIa) were assessed as components of the contact pathway. Prothrombin fragments 1 and 2 (F1+2), thrombin-antithrombin complex, and D-dimer were assessed as components of the thrombin pathway. Plasminogen (PLG), α₂-antiplasmin (A2A), and PLG

Table 6. Summary of the Results of the Systemic and the Reproductive Toxicity Studies

Test	Findings
Systemic toxicity	
Acute toxicity of photolyzed riboflavin in rats	<ul style="list-style-type: none"> • No test-article-related changes in mortality, clinical signs, body weight, food consumption, clinical pathology, organ weights, or macroscopic pathology evaluations were observed.
Acute toxicity of Mirasol PRT-treated platelets in rats*	<ul style="list-style-type: none"> • No test-article-related mortality, clinical observations, body weight changes, food consumption changes, or macroscopic changes were observed.
Acute toxicity of Mirasol PRT-treated platelets in dogs*	<ul style="list-style-type: none"> • No test-article-related mortality, clinical observations, body weight changes, food consumption changes, or macroscopic changes were observed.
Subchronic toxicity of Mirasol PRT-treated plasma in dogs	<ul style="list-style-type: none"> • No toxicologically significant findings were observed in clinical, physical, electrocardiographic, or ophthalmoscopic examinations. • Macroscopic observations did not yield any test-article-related effects in dogs of either sex. • Test-article-related findings for decreased body weight, increased eosinophils, increased chloride, and decrease thymus weights with associated microscopic findings were noted; however, these findings were minor and were not considered to be toxicologically significant.
Reproductive toxicity	
Embryo-fetal development in rats exposed to Mirasol PRT-treated plasma	<ul style="list-style-type: none"> • No test-article-related clinical signs of toxicity were observed during treatment period, no dam had all fetuses resorbed, and no deaths occurred during study. • Treated group exhibited significant increase in gestation body weight on gestation day 18, and body weight changes were significantly increased on gestation day intervals 12-15, 6-18, and 0-20. Increase correlated with significantly increased food consumption in treated-group animals. • No significant differences were observed in fetal weight (g), gravid uterine weight, adjusted weight change, postimplantation loss, nonviable fetuses/, litter size, early and late resorptions, mean no. of corpora lutea/litter, no. of implantation sites, preimplantation loss, no. of viable fetuses/letter, and sex ratio (% males/litter). • No external or test-article-related visceral abnormalities were reported in any of fetuses examined. No significant difference in skeletal malformation and variations were observed in control and treated groups.

* Indicates that study was done with and without γ -irradiation of treated blood products.

activator inhibitor type 1 (PAI-1) were assessed as components of the fibrinolytic pathway. CH50/classical pathway, C3a, C5a, and C1 esterase inhibitor were assessed as components of the complement pathway, and antithrombin III, protein C (PC), and protein S (PS) were assessed as components of the coagulation pathway. Mean values for test and control articles were compared using Student *t* test. In addition, assays for chemotaxis (response to C5a) and assays of myeloperoxidase release from cytochalasin B-treated human polymorphonuclear leukocytes (response to C3a and C5a) were used to further

assess physiologic relevance of any changes after treatment.

Pharmacokinetics of Photolyzed ¹⁴C-Riboflavin in Rats

The in vivo pharmacokinetics and elimination of riboflavin and its photoproducts in rat plasma treated with the Mirasol PRT System were also studied. Sterile-filtered rat plasma was mixed with 500 $\mu\text{mol/L}$ ¹⁴C-riboflavin, treated by the Mirasol PRT System, and administered to male CD rats by IV injection in a single dose. The study consisted of 2 treatment groups of male Crl: CD(SD) IGS BR

Table 7. Summary of the Results for the Mammalian Erythrocyte Micronucleus Study

Test	Findings
Test for clastogenicity in mammalian cells (in vivo; mouse bone marrow micronucleus test)	<ul style="list-style-type: none"> • No mortality occurred at any dose level during course of micronucleus study. • Clinical signs, noted on days after dose administration, included lethargy and piloerection in males and females dosed with untreated platelets and in all animals dosed with test article (neat and diluted). In addition, hunched position was noticed in animals at middle (1:1 dilution) and high (neat) dose. All other mice treated with vehicle or positive control articles appeared normal during course of study. • No appreciable reduction in ratio of polychromatic erythrocytes to total erythrocytes was observed in groups treated with negative control or test article, relative to vehicle control groups. This suggests that negative control and test article did not inhibit erythropoiesis. No significant increase in incidence of micronucleated polychromatic erythrocytes in negative control or test-article-treated groups relative to respective vehicle control groups was observed, regardless of sex, dose, or bone marrow collection time.

rats. A treated group of 3 rats (group 1) and a treated group of 15 rats (group 2) each received the radiolabeled test article at a target radiolabeled riboflavin dose level of 0.48 mg/kg with approximately 555,000 Bq per animal (mean actual dose levels of approximately 0.38 mg/kg and 547,600 Bq per animal). The dose volume for both groups was approximately 4 mL per animal.

All animals were observed twice daily for morbidity, mortality, injury, and the availability of food and water. Urine and feces were collected from group 1 animals predosing and at intervals through 144 hours postdosing. In addition, whole blood was collected, and plasma was prepared, predosing and at the end of the experiment. The tissues collected for postmortem study were residual carcass, pelt, and tail skin, as well as small intestine, large intestine, spleen, kidneys, liver, lymphatics, and right and left femur bone marrow. The concentration of radioactivity in whole blood, plasma, urine, feces, tissues, and cage residues was determined, and analyses of the kinetics of ^{14}C -labeled compound excretion were performed. The plasma and tissue concentrations obtained were for radioactivity associated with riboflavin and its photoproducts, rather than for individual compounds. Table 5 summarizes the samples analyzed for each group. The pharmacokinetic parameters of riboflavin and photoproduct in plasma, small intestine, large intestine, spleen, kidneys, liver, lymphatic, and femur bone marrow tissues were determined from mean concentration-time data in the test animal.

Appropriate validated software was used to calculate these parameters.

Measurement of Leachable and Extractable Compounds in Mirasol PRT-Treated Products

In addition to standard medical device testing for plastic components, several novel methods were applied to evaluate the levels of leachable and extractable compounds in human blood products after treatment with the Mirasol PRT System. This non-GLP analysis was performed to identify and quantify compounds that might migrate directly into the platelet product as a result of treatment with the Mirasol PRT System. Untreated and Mirasol PRT-treated products were analyzed by gas chromatography and mass spectrometry for volatile and semivolatile organic compounds, by inductively coupled plasma-atomic emission spectroscopy (ICP) for metals, and by Fourier transform infrared spectroscopy (FTIR) for polymeric materials.

RESULTS

In all studies of toxicity, γ -irradiated Mirasol PRT-treated products yielded the same results as non- γ -irradiated, Mirasol PRT-treated products.

Systemic Toxicity

No toxicologically significant findings were observed in any of the studies of acute toxicity (Table 6). No toxicologically significant findings were observed in the repeated-dose toxicity study

Table 8. Clotting Factor Levels in Treated and Untreated, Control Platelet Preparations After 5 Days of Storage

	Reference range ¹⁰⁴	Control (n = 6)	Treated* (n = 6)	P
FXI (IU/mL)	0.42-1.44	0.86 ± 0.13	0.59 ± 0.12	<.05
FXII (IU/mL)	0.40-1.52	1.40 ± 2.0	1.07 ± 0.23	<.05
PK (IU/mL)	0.65-1.35 [†]	1.11 ± 0.06	0.73 ± 0.10	<.05
HMWK (IU/mL)	0.65-1.35 [†]	0.78 ± 0.12	0.68 ± 0.29	NS [‡]
Antithrombin III (IU/mL)	0.72-1.45	0.82 ± 0.01	0.84 ± 0.11	NS
PC (IU/mL)	0.58-1.64	1.02 ± 0.19	0.75 ± 0.8	<.05
PS (IU/mL)	0.56-1.68	0.16 ± 0.08 (3 of 6 @ <0.11)	0.13 ± 0.03 (4 of 6 @ <0.11)	NS
PAI-1 (IU/mL)	0-31 [†]	3.6 ± 1.3	17 ± 4	<.05
PLG (IU/mL)	0.68-1.44	1.04 ± 0.10	0.88 ± 0.08	<.05
A2A (IU/mL)	0.72-1.32	0.88 ± 0.05	0.82 ± 0.07	<.05
D-D (ng/mL)	0-255 [†]	5/6 <110, 1 at 139	<110	NA [§]
FXIIa (ng/mL)	0-3.6 [†]	1.5 ± 0.2	4.8 ± 0.9	<.05
F1 + 2 (nmol/L)	0.4-1.8 [†]	0.5 ± 0.0	1.0 ± 0.1	<.05
TAT (ng/mL)	0.0-5.0 [†]	3 of 6 <0.6; 0.7 ± 0.1	1.1 ± 0.4	NA
C3a (ng/mL)	0-940 [†]	5510 ± 552	16594 ± 4939	<.05
C5a (ng/mL)	4.7-9.5 [†]	8 ± 1	16 ± 3	<.05
C1 esterase Inhibitor (mg/mL)	8.0-19.5 [†]	11 ± 2	7 ± 1	<.05
CH50/classic pathway (U/mL)	176-282 [†]	281 ± 56	21 ± 18	<.05

Abbreviations: PC, protein C; PS, protein S; A2A, alpha-2-plasmin; D-D, D-dimer; TAT, thrombin-antithrombin complexes.

* Values reported for plasma from Mirasol PRT-treated platelets are corrected for the dilution due to the addition of riboflavin.

[†] Reference range specific to laboratory or assay kit.

[‡] $P > .05$.

[§] Most data below detection limit of assay.

(Table 6). In the repeated-dose toxicity study, the levels of riboflavin and lumichrome in blood samples from animals receiving Mirasol PRT-treated products were below the limits of quantification, as were the levels in blood samples from animals receiving untreated plasma. These results were consistent with the observed rapid clearance of riboflavin after IV administration, both in the literature⁴⁰ and in the pharmacokinetic study with ¹⁴C-riboflavin in Mirasol PRT-treated products.

Developmental Toxicity and Genotoxicity

No developmental toxicity was observed in the embryo-fetal development study. All fetuses were

examined for malformations and developmental variations (Table 6). No mutagenicity was observed in the Ames test (with and without S9 metabolic activating fraction) for treated or control human platelets, or for lumichrome. The in vitro and in vivo tests for clastogenicity in mammalian cells (chromosomal aberration in cultured CHO cells and micronucleus test in mouse bone marrow cells, respectively) were also performed with Mirasol PRT-treated products. Human platelets treated with the Mirasol PRT System gave negative results in all genotoxicity experiments. Detailed results of the in vivo genotoxicity experiment are summarized in Table 7.

Neoantigenicity and Cytotoxicity

Results of studies using ¹⁴C-labeled riboflavin and exposure of platelets and plasma to UV light did not demonstrate any detectable binding of riboflavin or its photoproducts to platelets or to

Table 9. Clotting Factor Levels in Mirasol PRT-Treated FFP

	Reference range ¹⁰⁴	Mirasol FFP	
		Treated*, assayed postillumination (n = 27)	Treated*, stored for 1 y (n = 27)
Fibrinogen (mg/dL)	145-385	227 ± 60	216 ± 32
FVIII activity (IU/mL)	0.52-1.55	0.6 ± 0.2	0.9 ± 0.2
FII activity (IU/mL)	0.65-1.54	0.8 ± 0.1	0.8 ± 0.1
FV activity (IU/mL)	0.54-1.45	0.7 ± 0.1	0.6 ± 0.1
FVII activity (IU/mL)	0.62-1.65	0.9 ± 0.2	0.6 ± 0.1
FIX activity (IU/mL)	0.45-1.48	0.8 ± 0.2	0.7 ± 0.1
FX activity (IU/mL)	0.68-1.48	0.8 ± 0.1	0.8 ± 0.1
FXI activity (IU/mL)	0.42-1.44	0.8 ± 0.2	0.6 ± 0.1
FXII activity (IU/mL)	0.40-1.52	1.0 ± 0.3	0.8 ± 0.1
Total protein (g/L)	48-64	57 ± 2	53 ± 3
Antithrombin III (IU/mL)	0.72-1.45	0.8 ± 0.1	1.0 ± 0.1
PC (IU/mL)	0.58-1.64	1.1 ± 0.2	0.9 ± 0.2
PS (IU/mL)	0.56-1.68	0.8 ± 0.1	1.0 ± 2
PLG (IU/mL)	0.68-1.44	0.8 ± 0.2	0.9 ± 0.2
A2A (IU/mL)	0.72-1.32	1.0 ± 0.2	0.9 ± 0.1
HMWK (IU/mL)	0.65-1.35	0.8 ± 0.2	0.7 ± 0.2

* Values reported for Mirasol PRT-treated plasma are corrected for the dilution due to the addition of riboflavin.

Table 10. Complement Test Results in Treated and Untreated, Control Platelet Preparations

Sample	Chemotaxis assay	Enzyme release assay
Reference range	Positive control = 60.4 Negative control = 16.8	Positive control = 1.13 Negative control = 0.32
Control (n = 6)	22.8 ± 1.8	0.41 ± 0.03
Treated (n = 6)	24.2 ± 5.0	0.45 ± 0.14
P	NS	NS

NOTE: NS, $P > .05$.

plasma proteins. No evidence of neoantigenicity was observed with the Ouchterlony assay, indicating that no new antigens were formed during treatment with the Mirasol PRT System. Treatment with the Mirasol PRT System did not result in greater immunoglobulin G binding than what was observed in comparison with untreated controls, when assessed with the Capture-P assay. In the tests of lumichrome cytotoxicity, and of the cytotoxicity of Mirasol PRT-treated products, no cytotoxicity was observed.

Hemocompatibility

In tests of hemocompatibility, no hemolysis was observed. In functional assessments, when mixed with thrombocytopenic whole blood, the function of Mirasol PRT-treated platelets was well preserved, in comparison with controls.⁵⁵ Treated platelets displayed no evidence of hyperactivation or hypercoagulability. In evaluation of plasma

removed from platelet products stored for 5 days, statistically significant changes in PRT-treated samples were observed in FXI, FXII, PK, PC, PS, PAI-1, PLG, A2A, FXIIa, and F1+2. The activities of FXII, PK, PC, PAI-1, PLG, and F1+2 were within the reference range. The activities of the remaining proteins (FXI, PS, A2A, HMWK, and FXIIa) were outside the reference range, but the changes were not deemed to be clinically relevant compared with the levels in normal, untreated blood (Table 8). Similar studies on platelets after γ -irradiation yielded similar findings. Because analysis of plasma protein levels in platelet products stored at room temperature for extended periods shows degradation in protein quality due to storage alone, Mirasol PRT-treated FFP was assessed for factor activity. Results are presented in Table 9 for samples held for 8 hours at room temperature after collection and before treatment with Mirasol PRT to simulate worst-case conditions for routine treatment of plasma products. Samples were then treated and stored frozen until analysis of factors in standard coagulation assays. Results for all samples analyzed remained within historical reference ranges based on untreated, historical, and concurrent controls. Results for Mirasol PRT-treated FFP are representative of values obtained for plasma from Mirasol PRT-treated platelets immediately posttreatment.

The plasma from Mirasol PRT-treated human products contained significantly more C3a desarg

Table 11. Results of Pharmacokinetic Study and Test for Leachables and Extractables

Test	Findings
Pharmacokinetics of photolyzed ¹⁴ C-riboflavin in Mirasol PRT-treated rat plasma	Approximately 95% of total radioactivity administered to animals was eliminated within approximately 260-275 h postadministration Based upon observed areas under the curve, tissues exhibiting highest to lowest overall exposure were liver > kidneys > large intestine > small intestine > spleen > femur bone marrow > lymphatics Observed mean $t_{1/2}$ estimates determined from radioactivity excretion rate data in urine (55.3 h) and determined from radioactivity excretion rate data in feces (55.1 h) were in excellent agreement with $t_{1/2}$ estimate (52.1 h) for plasma
Assessment of leachable and extractable compounds using standard medical device biocompatibility testing	No significant differences before and after exposure to treatment conditions
Assessment of leachable and extractable compounds in Mirasol PRT-treated human platelets and untreated human platelets	No polymeric material was detected by FTIR in either treated or untreated platelet extracts All metals detected in treated-platelet extracts were present in untreated-platelet extracts in similar amounts.

Abbreviation: FTIR, Fourier transform infrared spectroscopy.

and C5a desarg (2.7- and 1.8-fold higher, respectively) and significantly lower values for C1 esterase (1.8-fold lower) and CH50 (14-fold lower) than in plasma from untreated controls (Table 8). The desarg forms of C3a and C5a are biologically inactive forms of the anaphylatoxins produced in the common pathway of complement activation. The increase in these inactive forms in treated products and decrease in the active component, CH50 (total hemolytic complement), are consistent. Little to no complement-derived activity was observed in plasma from either Mirasol PRT-treated platelets or untreated controls, and there were no significant differences between treated and control results (Table 10). Similar studies performed on products after γ -irradiation yielded identical results (data not shown).

Pharmacokinetics of Photolyzed ^{14}C -Riboflavin in Rats

After a single IV administration of Mirasol-treated plasma containing photolyzed ^{14}C -riboflavin, the radioactivity was well distributed from the whole blood to tissues selected for assay within the first hour postdose. Most of the excreted urinary radioactivity was recovered by 12 hours postdose, and more than half of all radioactivity was excreted in urine. Blood levels of radioactivity declined rapidly post-dose, as expected from studies of riboflavin metabolism and excretion in humans.⁴⁰ Measurements of the radioactivity associated with the ^{14}C -riboflavin-treated plasma indicated rapid initial apparent distribution (and/or clearance) from the systemic circulation that appeared to be complete within the first 8 to 48 hours postdose. Table 11 contains a summary of the key study results.

Leachables and Extractables

The leachables and extractables analyses detected no polymeric material in either test or control platelet products (Table 11). The Mirasol illumination/storage bag does not contain the plasticizer di(2-ethylhexyl)phthalate (DEHP), and testing verified that this plasticizer was not present in treated and stored products. No toxicologically relevant concentrations of metals were found. These results correlate with those from the biocompatibility testing of the Mirasol illumination/storage bag elements—all elements are biocompatible.

DISCUSSION

The design of the Mirasol PRT device, the strategy of the investigative work, and the specific experiments that were conducted to evaluate toxicity were based first on scientific considerations and then on the need to comply with regulatory requirements. Consideration was also given to minimizing the need for animal experimentation as far as permitted by the information required and the need for reliable scientific experiments. The system design and the experimental program were intended to answer potential scientific and medical concerns and to permit a comprehensive safety assessment of the Mirasol PRT System. The design of the preclinical safety studies was also intended to conform to specific regulatory requirements in different countries, adapted to the nature of this innovative device. The official guidelines used for the study design included ISO 10993⁷⁹ for medical devices and the general toxicity and pharmacokinetic guidelines for drugs of the International Conference on Harmonization of Technical Requirements for Registration of Pharmaceuticals for Human Use.^{82,105} The nonclinical toxicity program was based throughout on studying preparations that had undergone Mirasol PRT treatment, a process shown to be effective in reducing levels of several types of infectious agents.^{46,49,50,106,107}

The testing of treated blood products has examined the potential for toxicity resulting from the entire process (from addition of riboflavin and controlled illumination in the Mirasol illumination bag) and has examined the potential for harm to the blood products themselves. Assessments included effects on the cells and plasma proteins, the toxicity of riboflavin itself and its photoproducts, and the effects on other system components comprising the medical grade plastic connectors and bags used in this process. The results obtained have been compared with extensive literature on the effects of riboflavin and its photoproducts in animals and humans arising from its medical use as a vitamin and more widely as a food colorant or additive. The methods used have been those recommended by current official guidelines interpreted according to the special nature and properties of platelets and plasma and designed to provide as much information as possible from focused in vitro and in vivo studies. The difference in size between human patients and in vivo animal

test systems has restricted application of the standard toxicological strategy of administering a range of doses far exceeding the human dose to reveal possible toxic hazards. In studying Mirasol PRT-treated blood products, the largest volumes or doses have been given consistent with maintaining the welfare of any animals used and to separate effects due to toxicities of the components from toxicities due to volume overload alone.

The results of the *in vivo* and *in vitro* toxicity tests have reproduced previously reported findings (Table 1) for riboflavin and its photoproducts in animals and humans under relevant conditions, despite the special circumstances of the presence of blood cells and proteins and the use of various medical-grade plastic connectors and storage containers. There has been no indication, even after repeated administration of the maximum achievable exposure to Mirasol PRT-treated products (6 times per week for 13 weeks in the subchronic test in the dog), of any untoward effect on the functions of the major physiologic systems or in producing any target organ toxicity. Furthermore, a variety of special tests have shown that no detectable new antigens are present on Mirasol PRT System-treated cells or proteins and that such preparations did not stimulate antibody formation. The wide range of toxicity tests has also excluded any genotoxic risk and any toxic effect on pregnant animals and on embryo-fetal development. Local tolerance at the site of IV infusion has been shown to be good. The treated platelet preparations have not shown any cytotoxicity, nor did they leach unwanted substances from the containers and connectors used in the Mirasol system.

These results indicate the lack of specific toxicity of Mirasol PRT-treated products, and because of the experiments done, they give a strong indication that the same medical acceptability is likely also to apply to both human platelets and to human plasma treated in the same way. The absence of toxicity was considered likely given the extensive prior experience in humans and the laboratory experience of the safety of the chemicals involved (riboflavin, lumichrome, and the other riboflavin photoproducts), but it has been very important to confirm that prediction and to extend it to include the special cells and plasma proteins also involved in the Mirasol System. To this should be added the evidence of retention of acceptable platelet function and survival after

treatment, which has been the subject of additional peer-reviewed publications.^{46,51,55,57}

Thus, there is evidence of the safety of this process in a broad range of laboratory tests from which a qualitative extrapolation of safety in patients can be made. Quantitative prediction of safety, or at least of the margin of safety represented by the difference between the relative exposures of man and animals and the observed “no effect level” in the laboratory, must be more cautious, but there are relevant data on which to base some realistic predictions.

The likely exposure of humans to riboflavin and its photoproducts is 0.077 mg/kg for each unit of product transfused. Likely exposure of recipients in a clinical setting was calculated on the assumption of a mean recipient weight of 70 kg, nominal riboflavin solution concentration of 500 $\mu\text{mol/L}$, mean riboflavin photoconversion of 18%, and nominal riboflavin solution volume of 35 mL. This exposure level can be compared with the lethal dose at which 50% of the population survives reported for IV riboflavin in mice (50–100 mg/kg)³⁰ yielding a safety factor of at least $50/0.077 = 649$. However, this safety margin calculation is somewhat illusory because there is no reliable indication from the clinical literature, or from appropriate animal experiments, of the toxic dose of riboflavin in man, and so, the upper bound of the “no toxic effect level” remains unknown. It is reasonable only to conclude that there is no reason to predict toxicity in man due to parenteral exposure to riboflavin and its photoproducts in Mirasol PRT-treated products until much higher exposures are reached than would ever be feasible from platelet or plasma transfusions. For similar reasons, comparison of exposure from treated platelets with the ADI for riboflavin of $0 \text{ to } 0.54 \text{ mg} \times \text{kg}^{-1} \times \text{d}^{-1}$ is no less misleading because the ADI refers to life-long, daily exposure to the substance in the diet and is not an appropriate comparator for a short-term parenteral exposure even if the latter is repeated up to several times per week, at least for a few months.

In these circumstances, which are not unusual for physiologically essential substances, the best that the toxicologist can offer is the prediction that riboflavin and its photoproducts have low experimental toxicity *in vivo*. Published work suggests that administration of the agent in the extreme may result in renal damage due to crystalluria as the most prominent toxic effect in rodents. This work is

derived from studies in which animals were given very large IV doses, and the harmful effect is very probably a consequence of the limited aqueous solubility of riboflavin. This is far removed from any circumstance associated with administration of Mirasol PRT-treated products in the transfusion setting. There is no realistic indication that chemical toxicity will be a limiting factor in any likely course of treatment with Mirasol PRT-treated products.

Results from all genotoxicity studies with the Mirasol PRT System have been negative. This has included evaluation of all photoproducts and photolyzed platelets and plasma at the maximum tolerable or maximum soluble limits for these assays. The absence of any indication of potential tumorigenic effects in the subchronic and other tests, the lack of genotoxicity, the chemical nature of riboflavin and its photoproducts, and the lack of suspicion arising from the physiologic role and human experience of this vitamin obviated the need for further carcinogenicity testing. Prior carcinogenicity testing of riboflavin in historical studies (Table 1) is consistent with these findings and also supportive of the lack of carcinogenic potential of this compound.

The Mirasol PRT System has been evaluated in an extensive preclinical safety evaluation program. The historical literature available on riboflavin and its photoproducts as well as the results of this work support its safety profile in man. The results obtained from the studies conducted and reported here are consistent with these data and indicate that the risk of toxicity due to the use of this system in the transfusion setting should be low.

ACKNOWLEDGMENTS

We would like to thank Mr Jon White, Dr Ali Faqi, Dr Mary Sherman, Ms Sheryl Loux, Ms Leigh Bergeron, Mr Devaki Sadhu, Dr Florence Burleson, and Dr Ramadevi Gudi for conducting the toxicology research work described in this review article. We also would like to acknowledge the work of Dr George Garratty's and Dr Patsy Giclas' laboratories for their analysis of platelet immunologic properties. We are grateful to Deanna Gamp, Nick Hovenga, Suzann Doane, Denise Gilmour, Shawn Keil, and Christopher Hardwick for their laboratory support work for Navigant's toxicology program.

REFERENCES

1. Busch MP, Kleinman SH, Nemo GJ: Current and emerging infectious risks of blood transfusions. *JAMA* 289:959-962, 2003
2. Vamvakas EC, Blajchman MA: Deleterious clinical effects of transfusion-associated immunomodulation: Fact or fiction? *Blood* 97:1180-1195, 2001
3. Lee JH, Klein HG: From leukocyte reduction to leukocyte transfusion: The immunological effects of transfused leukocytes. *Baillieres Best Pract Res Clin Haematol* 13:585-600, 2000
4. Gould DS, Auchincloss HJ: Direct and indirect recognition: The role of MHC antigens in graft rejection. *Immunol Today* 20:77-82, 1999
5. Llewelyn CA, Hewitt PE, Knight RS, et al: Possible transmission of variant Creutzfeldt-Jakob disease by blood transfusions. *Lancet* 363:417-421, 2004
6. Hewitt PE, Llewelyn CA, Mackenzie J, et al: Creutzfeldt-Jakob disease and blood transfusion: Results of the UK Transfusion Medicine Epidemiological Review study. *Vox Sang* 91:221-230, 2006
7. Wroe SJ, Pal S, Siddique D, et al: Clinical presentation and pre-mortem diagnosis of variant Creutzfeldt-Jakob disease associated with blood transfusion: A case report. *Lancet* 368:2061-2067, 2006
8. Peden AH, Head MW, Ritchie DL, et al: Preclinical vCJD after blood transfusion in a PRNP codon 129 heterozygous patient. *Lancet* 364:527-529, 2004
9. Blajchman MA: Bacterial contamination of cellular blood components: Risks, sources and control. *Vox Sang* 87 (Suppl):98-103, 2004
10. US Food and Drug Administration. Minutes of 36th Meeting of the Blood Product Advisory Committee; 1992, May 28-29; Bethesda, Maryland. Washington (DC), USFDA; 1992
11. Yomtovian R, Lazarus HM, Goodnough LT, et al: A prospective microbiologic surveillance program to detect and prevent the transfusion of bacterially contaminated platelets. *Transfusion* 33:902-909, 1993
12. Schroeder ML: Transfusion-associated graft-versus-host disease. *Br J Haematol* 117:275-287, 2002
13. Saunders C, Herbert P, Rowe G, et al: In vitro evaluation of the PALL Leukotrap Affinity Prion Reduction Filter as a secondary device following primary leucoreduction. *Vox Sang* 89:220-228, 2005
14. Blajchman MA, Goldman M, Baeza F: Improving the bacteriological safety of platelet transfusions. *Transfus Med Rev* 18:11-24, 2004
15. Cohen ND, Munoz A, Reitz BA: Transmission of retroviruses by transfusion of screened blood in patients undergoing cardiac surgery. *N Engl J Med* 320:1172-1176, 1989
16. Nelson KE, Donahue JG, Munoz A, et al: Transmission of retroviruses from seronegative donors by transfusion during cardiac surgery. A multicenter study of HIV-1 and HTLV-I/II infections. *Ann Intern Med* 117:554-559, 1992

17. Grant IH, Gold JW, Wittner M, et al: Transfusion-associated acute Chagas disease acquired in the United States. *Ann Intern Med* 111:849-851, 1989
18. Leiby DA, Lenes BA, Tibbals MA, et al: Prospective evaluation of a patient with *Trypanosoma cruzi* infection transmitted by transfusion. *N Engl J Med* 341:1237-1239, 1999
19. Wu DG, Mason B, Jong J, et al: Parvovirus B19 transmission by a high-purity factor VIII concentrate. *Transfusion* 45:1003-1010, 2005
20. Biggerstaff BJ, Petersen LR: Estimated risk of transmission of the West Nile virus through blood transfusion in the US. *Transfusion* 43:1007-1017, 2002
21. Custer B, Tomasulo PA, Murphy EL, et al: Triggers for switching from minipool testing by nucleic acid technology to individual-donation nucleic acid testing for West Nile virus: Analysis of 2003 data to inform 2004 decision making. *Transfusion* 44:1547-1554, 2004
22. Hellstern P: Solvent/detergent-treated plasma: Composition, efficacy, and safety. *Curr Opin Hematol* 11:346-350, 2004
23. Williamson LM, Cardigan R, Prowse CV: Methylene blue-treated fresh-frozen plasma: What is its contribution to blood safety? *Transfusion* 43:1322-1339, 2003
24. Corash L, Lin L: Novel processes for inactivation of leukocytes to prevent transfusion-associated graft-versus-host disease. *Bone Marrow Transplant* 33:1-7, 2004
25. Lin L, Dikeman R, Molini B, et al: Photochemical treatment of platelet concentrates with amotosalen and long-wavelength ultraviolet light inactivates a broad spectrum of pathogenic bacteria. *Transfusion* 44:1496-1504, 2004
26. Lin L, Hanson CV, Alter HJ, et al: Inactivation of viruses in platelet concentrates by photochemical treatment with amotosalen and long-wavelength ultraviolet light. *Transfusion* 45:580-590, 2005
27. Unna K, Greslin JG: Studies on the toxicity and pharmacology of riboflavin. *J Pharmacol* 76:75-80, 1942
28. Hayashi M, Kishi M, Sofuni T, et al: Micronucleus tests in mice on 39 food additives and eight miscellaneous chemicals. *Food Chem Toxicol* 26:487-500, 1988
29. Kuhn R, Boulanger P: Über die Giftigkeit der Flavine. *Hoppe Seylers Z Physiol Chem* 241:233-238, 1936
30. Studer A, Zbinden G, Uehlinger E: Die pathologie der avitaminosen und hypervitaminosen, in Buchner F, Letterer E, Roulet F (eds): *Handbuch der Allgemeinen Pathologie*, Band II. Berlin, Springer Verlag, 1962, pp 734-987
31. Williams Jr JR, Grantham PH, Yamamoto RS, et al: Effect of dietary riboflavin on azo dye reductase in liver and in bacteria of cecal contents of rats. *Biochem Pharmacol* 19:2523-2525, 1970
32. Scotto JM, Stralin HG, Lageron A, et al: Influence of carbon tetrachloride or riboflavin on liver carcinogenesis with a single dose of aflatoxin b1. *Br J Exp Pathol* 56:133-138, 1975
33. Leclerc J: Influence of thiamine, riboflavin and vitamin B6 content in food on tissue content of these vitamins in the female lactating rat and young rats. *Ann Nutr Aliment* 28:11-20, 1974
34. Lim-Sylianco C: Genotoxic studies of street soils in Quezon City, Philippines, in Sumino K (ed): *Environmental and occupational chemical hazards*, no. 8. Proceedings of the Asia-Pacific Symposium in Environmental and Occupational Toxicology, 1987, Oct 4-10; Singapore. Kobe, Japan, Kobe University School of Medicine
35. Munoz N, Hayashi M, Bang LJ, et al: Effect of riboflavin, retinol, and zinc on micronuclei of buccal mucosa and of esophagus: A randomized double-blind intervention study in China. *J Natl Cancer Inst* 79:687-691, 1987
36. Sugiyama M, Ando A, Nakao K, et al: Influence of vitamin B2 on formation of chromium(V), alkali-labile sites, and lethality of sodium chromate(VI) in Chinese hamster V-79 cells. *Cancer Res* 49:6180-6184, 1989
37. Haveland-Smith RB: Evaluation of the genotoxicity of some natural food colours using bacterial assays. *Mutat Res* 91:285-290, 1981
38. De Flora S: Study of 106 organic and inorganic compounds in the *Salmonella*/microsome test. *Carcinogenesis* 2:283-298, 1981
39. Kale H, Harikumar P, Nair PM, et al: Assessment of the genotoxic potential of riboflavin and lumiflavin. *Mutat Res* 298:9-16
40. Zemleni J, Galloway JR, McCormick DB: Pharmacokinetics of orally and intravenously administered riboflavin in healthy humans. *Am J Clin Nutr* 63:54-66, 1996
41. Kumar V, Lockerbie O, Keil SD, et al: Riboflavin and UV-light based pathogen reduction: Extent and consequence of DNA damage at the molecular level. *J Photochem Photobiol* 80:15-21, 2004
42. Hardwick CC, Herivel TR, Hernandez SC, et al: Separation, identification and quantification of riboflavin and its photoproducts in blood products using high-performance liquid chromatography with fluorescence detection: a method to support pathogen reduction technology. *J Photochem Photobiol* 80:609-615, 2004
43. Holmstrom B: Spectral studies of the photobleaching of riboflavin phosphate. *Arkiv for Kemi* 22:281-301, 1964
44. Dekker RH, Srinivasan BN, Huber JR, et al: Photochemistry of flavins—I. Conventional and laser flash photolysis study of alloxazines. *J Photochem Photobiol* 18:457-466, 1973
45. Koziol J: Studies of flavins in organic solvents—I. Spectral characteristics of riboflavin, riboflavin tetrabutylate and lumichrome. *J Photochem Photobiol* 5:41-54, 1966
46. Ruane PH, Edrich R, Gampp D, et al: Photochemical inactivation of selected viruses and bacteria in platelet concentrates using riboflavin and light. *Transfusion* 44:877-885, 2004
47. Fast LD, DiLeone G, Cardarelli G, et al: Mirasol PRT treatment of donor white blood cells prevents the development of xenogeneic graft-versus-host disease in Rag 2^{-/-}gc^{-/-} double knockout mice. *Transfusion* 46:1553-1560, 2006
48. Fast LD, DiLeone G, Li J, et al: Functional inactivation of white blood cells by Mirasol treatment. *Transfusion* 46:642-648, 2006
49. Cardo LJ, Rentas FJ, Ketchum L, et al: Pathogen inactivation of *Leishmania donovani* infantum in plasma and platelet concentrates using riboflavin and ultraviolet light. *Transfusion* 90:85-91, 2006
50. Rentas FJ, Harman R, Salata J, et al: Inactivation of *Orientia tsutsugamushi* in RBCs, plasma, and platelets using riboflavin and light, as demonstrated in an animal model. *Transfusion* 47:240-247, 2007
51. Li J, de Korte D, Woolum MD, et al: Pathogen reduction of buffy coat platelet concentrates using riboflavin and light: Comparisons with pathogen-reduction technology-treated apheresis platelet products. *Vox Sang* 87:82-90, 2004

52. Li J, Lockerbie O, de Korte D, et al: Evaluation of platelet mitochondria integrity after treatment with Mirasol pathogen reduction technology. *Transfusion* 45:920-926, 2005
53. Kumar V, Motheral T, Luzniak G, et al: Mirasol (TM) plasma pathogen reduction technology: Riboflavin-based process conserves, protein C, protein S and antithrombin activities. *Transfusion* 43(Suppl):80A, 2003 (abstr)
54. Kumar V, Motheral T, Luzniak G, et al: Mirasol (TM) pathogen reduction technology: Riboflavin-based process retains functional activities of immunoglobulins. *Transfusion* 43 (Suppl):80A, 2003 (abstr)
55. Perez-Pujol S, Tonda R, Lozano M, et al: Effects of a new pathogen-reduction technology (Mirasol PRT) on functional aspects of platelet concentrates. *Transfusion* 45:911-919, 2005
56. Goodrich RP, Li J, Pieters H, et al: Correlation of in vitro platelet quality measurements with in vivo platelet viability in human subjects. *Vox Sang* 90:279-285, 2006
57. AuBuchon JP, Herschel L, Roger J, et al: Efficacy of apheresis platelets treated with riboflavin and ultraviolet light for pathogen reduction. *Transfusion* 45:1335-1341, 2005
58. WHO: Twenty-fifth report of the Joint FAO/WHO Expert Committee on Food Additives. Evaluation of certain food additives. *Tech Rep Ser* 669:1-48
59. EVM. Expert Group on Vitamins and Minerals; revised review of riboflavin. 2002. 1-48
60. Federation of American Societies for Experimental Biology: Evaluation of the health aspects of riboflavin and riboflavin-5'-phosphate as food ingredients, SCOGS-114. Life Sciences Research Office, Federation of American Societies for Experimental Biology, Bethesda, Maryland. Prepared for US Food and Drug Administration, Washington, DC, 1979, pp 1-25
61. Martindale W, Sweetman SC (eds): The complete drug reference. 33rd ed. London, UK, Pharmaceutical Press, 2002, p 1386
62. BIBRA Working Group: Riboflavin and its derivatives: toxicity profile. Epscm, Surrey, UK, BIBRA, 1990, pp 1-7
63. Direct Food Substances Affirmed as Generally Recognized as Safe; Riboflavin. 21 CFR Sect. 184.1695; 2001
64. European Commission Health and Consumer Protection Directorate-General, Scientific Committee on Food: Opinion of the Scientific Committee on Food on the tolerable upper intake level of vitamin B2. Brussels, Belgium, European Commission, 2000, pp 1-10
65. Sisson TR, Kendall N, Shaw E, et al: Phototherapy of jaundice in the newborn infant. II. Effect of various light intensities. *J Pediatr* 81:35-38, 1972
66. Sisson TR: Photodegradation of riboflavin in neonates. *Fed Proc* 46:1883-1885, 1987
67. Bates CJ: Human requirements for riboflavin. *Am J Clin Nutr* 46:122-123, 1987
68. Kostenbauder HB, DeLuca PP, Kowarski CR: Photo-binding and photoreactivity of riboflavin in the presence of macromolecules. *J Pharm Sci* 54:1243-1251, 1965
69. Yurdakok M, Erdem G, Tekinalp G: Riboflavin in the treatment of neonatal hyperbilirubinemia. *Turk J Pediatr* 30:159-161, 1988
70. Olsen JH, Hertz H, Kjaer SK, et al: Childhood leukemia following phototherapy for neonatal hyperbilirubinemia (Denmark). *Cancer Causes Control* 7:411-414, 1996
71. Speck WT, Rosenkranz HS: Phototherapy for neonatal hyperbilirubinemia—A potential environmental health hazard to newborn infants: A review. *Environ Mutagen* 1:321-336, 1979
72. Maisels MA, Kring EA: Does intensive phototherapy produce hemolysis in newborns of 35 or more weeks of gestation? *J Perinatol* 26:498-500, 2006
73. Stein G, Sperschneider H, Koppe S: Vitamin levels in chronic renal failure and need for supplementation. *Blood Purif* 3:52-62, 1985
74. Allman MA, Truswell AS, Tiller DJ, et al: Vitamin supplementation of patients receiving haemodialysis. *Med J Aust* 150:130-133, 1989
75. Gentile MG, Porrini M, Manna GM, et al: Water- and fat-soluble vitamin status in chronic renal insufficiency patients. *Contrib Nephrol* 98:89-97, 1992
76. Fouty B, Ferman F, Reves R: Riboflavin to treat nucleoside analogue-induced lactic acidosis. *Lancet* 352: 291-292, 1998
77. Schoenen J, Jacquy J, Lenaerts M: Effectiveness of high-dose riboflavin in migraine prophylaxis. *Neurology* 50:466-470, 1998
78. Luzzati R, Del Bravo P, Di Perri G, et al: Riboflavine and severe lactic acidosis. *Lancet* 353:901-902, 1999
79. International Organisation for Standardisation: International Standard ISO 10993 biological evaluation of medical devices—Part 1: Evaluation and testing. Geneva, Switzerland, International Organisation for Standardisation, 1993
80. International Organisation for Standardisation: International Standard ISO 10993 biological evaluation of medical devices—Part 11: Tests for systemic toxicity. Geneva, Switzerland, International Organisation for Standardisation, 1993
81. International Organisation for Standardisation. International Standard ISO 10993 biological evaluation of medical devices—Part 10: Tests for irritation and sensitization. Geneva, Switzerland, International Organisation for Standardisation, 1995
82. International Conference on Harmonisation of Technical Requirements for Pharmaceuticals (ICH): M3(R1) nonclinical safety studies for the conduct of human clinical trials for pharmaceuticals. Geneva, Switzerland, ICH, 1997
83. Ames B, McCann J, Yamasaki E: Methods for detecting carcinogens and mutagens with the Salmonella/mammalian-microsome mutagenicity test. *Mutat Res* 31:347-364, 1975
84. Green MHL, Muriel WJ: Mutagen testing using trp+ reversion in *Escherichia coli*. *Mutat Res* 38:3-32, 1976
85. Evans HJ: Cytological methods for detecting chemical mutagens, in Hollaender A (ed): Chemical mutagens, principles and methods for their detection, vol. 4. New York (NY), Plenum Press, 1976
86. Swierenga SHH, Heddle JA, Sigal EA, et al: Recommended protocols based on a survey of current practice in genotoxicity testing laboratories, IV. Chromosome aberration and sister-chromatid exchange in Chinese hamster ovary, V79 Chinese lung and human lymphocyte cultures. *Mutat Res* 246: 301-322, 1991
87. Galloway SM, Ardema AMJ, Ishidate MJ, et al: Report from working group on in vitro tests for chromosomal aberrations. *Mutat Res* 312:241-261, 1994
88. Heddle JA: A rapid in vivo test for chromosomal damage. *Mutat Res* 18:187-190, 1973
89. Mavournin KH, Blakey DH, Cimino, et al: The in vivo micronucleus assay in mammalian bone marrow and peripheral

blood. A report of the U.S. Environmental Protection Agency Gene-Tox Program. *Mutat Res* 239:29-80, 1990

90. Hayashi M, Tice RR, Macgregor JT, et al: In vivo rodent erythrocyte micronucleus assay. *Mutat Res* 312:293-304, 1994

91. Arndt PA, Garratty G, Hill J, et al: Two cases of immune haemolytic anaemia, associated with anti-piperacillin, detected by the 'immune complex' method. *Vox Sang* 83:273-278, 2002

92. Creech HJ, O'Connell AP: Immunochemistry of conjugates prepared from serum albumins and acridine nitrogen mustards (ICR mutagens). *Cancer Res* 41:3844-3851, 1981

93. Dise CA, Burch JW, Goodman DB: Direct interaction of mepacrine with erythrocyte and platelet membrane phospholipid. *J Biol Chem* 257:4701-4704, 1982

94. Chen RF: Fluorescence of quinacrine mustard conjugated to proteins. *Arch Biochem Biophys* 172:39-50, 1976

95. Jusko WJ, Levy G: Plasma protein binding of riboflavin and riboflavin 5'-phosphate in man. *J Pharm Sci* 58:58-62, 1969

96. Tapia G, Silva E: Photo-induced binding to the tryptophan residues of bovine and human serum albumin. *Radiat Environ Biol* 30:131-138, 1991

97. Diaz M, Ines Becker M, De Ioannes AE, et al: Development of monoclonal antibodies against a riboflavin-tryptophan photoinduced adduct: Reactivity to eye lens proteins. *J Photochem Photobiol* 63:762-767, 1996

98. Innis WSA, McCormick DB, Merrill AH: Variations in riboflavin binding by human plasma: Identification of immunoglobulins as the major proteins responsible. *Biochem Med* 34:151-156, 1985

99. Merrill AH, Innis-Whitehouse WSA, McCormick DB: Characterization of human riboflavin-binding immunoglobulins, in Edmondson DE, McCormick DB (eds): *Flavins and flavoproteins*. Berlin, Walter de Gruyter & Co, 1987, pp 445-448

100. International Organisation for Standardisation: International Standard ISO 10993 biological evaluation of medical devices—Part 5: Tests for in vitro cytotoxicity. Geneva, Switzerland, International Organisation for Standardisation, 1999

101. International Organisation for Standardisation: International Standard ISO 10993 biological evaluation of medical devices—Part 4: Selection of tests for interactions with blood. Geneva, Switzerland, International Organisation for Standardisation, 2002

102. Material Sciences Toxicology Laboratories: Autian method, ATP-I. Memphis (TN), University of Tennessee Center for the Health Sciences, 1977

103. Hemolysis—Rabbit blood, evaluation of hemodialyzers and dialysis membranes. 77-1294, 213. 1977. DHEW Publication (NIH).

104. Hellstern P, Sachse H, Schwinn H, et al: Manufacture and in vitro characterization of a solvent/detergent-treated human plasma. *Vox Sang* 63:178-185, 1992

105. International Conference on Harmonisation of Technical Requirements for Pharmaceuticals (ICH): S3B Pharmacokinetics: Guidance for Repeat Dose Tissue Distribution Studies. Geneva, Switzerland, ICH, 1995

106. Goodrich RP, Edrich RA, Li J, et al: The Mirasol PRT System for pathogen reduction of platelets and plasma: An overview of current status and future trends. *Transfus Apheresis Sci* 35:5-17, 2006

107. Goodrich RP, Edrich RA, Goodrich L, et al: The antiviral and antibacterial properties of riboflavin and light: Applications to blood safety and transfusion medicine, in Silva E, Edwards AM (eds): *Flavins: Photochemistry and photobiology*. Cambridge (UK), The Royal Society of Chemistry, 2006, pp 83-113

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市場投入およびサービス開始(Placing on the market and putting into service)

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C 98/8

EN

Official Journal of the European Union

23.4.2004

Commission communication in the framework of the implementation of the Council directive 89/336/EEC

(2004/C 98/05)

(Text with EEA relevance)

(Publication of titles and references of harmonized standards under the directive)

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ISO (J)	Reference and title of the standard	Reference documents	Reference of the superseded standard	Date of cessation of presumption of conformity of the superseded standard Note 1
Cenlec	EN 50065-1:2001 Signalling on low-voltage electrical installations in the frequency range 3 kHz to 148.5 kHz — Part 1: General requirements, frequency bands and electromagnetic disturbances		EN 50065-1:1991 +A1:1992 +A2:1995 +A3:1996 Note 2.1	Date expired (1.4.2003)
Cenlec	EN 50065-2:2003 Signalling on low-voltage electrical installations in the frequency range 3 kHz to 148.5 kHz			

	immunity for residential, commercial and light-industrial environments	IEC 61000-6-1:1997 (Modified)	EN 50082-1:1997 Note 2.1	1.7.2004
Cenlec	EN 61000-6-2:1999 Electromagnetic compatibility (EMC) — Part 6-2: Generic standards — Immunity for industrial environments	IEC 61000-6-2:1999	EN 50082-2:1995 Note 2.1	Date expired (1.4.2002)
Cenlec	EN 61000-6-2:2001 Electromagnetic compatibility (EMC) — Part 6-2: Generic standards — Immunity for industrial environments	IEC 61000-6-2:1999 (Modified)	EN 61000-6-2:1999 Note 2.1	1.7.2004
Cenlec	EN 61000-6-3:2001 Electromagnetic compatibility (EMC) — Part 6-3: Generic standards — Emission standard for residential, commercial and light-industrial environments	CISPR/IEC 61000-6-3:1996 (Modified)	EN 50081-1:1992 Note 2.1	1.7.2004
Cenlec	EN 61000-6-4:2001 Electromagnetic compatibility (EMC) — Part 6-4: Generic standards — Emission standard for industrial environments	IEC 61000-6-4:1997 (Modified)	EN 50081-2:1993 Note 2.1	1.7.2004

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EC Declaration of Conformity

Black Monolith Inc.

2001 Monolith Road, Sea of Tranquility, The Moon

We declare, under our solo responsibility, that the product:

Product Super computer

Model HAL 9000

to which this declaration relates complies with the provisions of following European Directives:

73/23/EEC as amended by 93/68/EEC — Directive on the harmonization of the laws of Member States relating to electrical equipment designed for use within certain voltage limits;

89/336/EEC as amended by 92/31/EEC and 93/68/EEC — Directive on the approximation of the laws of the Member States relating to electromagnetic compatibility.

Applied Harmonized Standards

EN 60950-1:2001

Information technology equipment – Safety – Part 1: General requirements

EN 55022:1998 + A1:2000 + A2:2003

Limits and methods of measurement of radio disturbance characteristics of information technology equipment

EN 61000-3-2:2000

Limits for harmonic current emissions (equipment input current up to and including 16 A per phase)

EN 61000-3-3:1995 + A1:2001

Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current ≤ 16 A per phase and not subject to conditional connection

Authorized Representative

Black Monolith Europe GmbH
2010 Discovery Avenue, Neverland, Germany

Signature: *Bowman*

David Bowman

Technical Director, Black Monolith Inc.

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