SUPPLEMENT II TO THE JAPANESE PHARMACOPOEIA SEVENTEENTH EDITION

Official from June 28, 2019

English Version

THE MINISTRY OF HEALTH, LABOUR AND WELFARE

Notice: This *English Version* of the Japanese Pharmacopoeia is published for the convenience of users unfamiliar with the Japanese language. When and if any discrepancy arises between the Japanese original and its English translation, the former is authentic.

The Ministry of Health, Labour and Welfare Ministerial Notification No. 49

Pursuant to Paragraph 1, Article 41 of Act on Securing Quality, Efficacy and Safety of Products Including Pharmaceuticals and Medical Devices (Act No. 145, 1960), this notification stated that a part of the Japanese Pharmacopoeia was revised as follows*.

NEMOTO Takumi

The Minister of Health, Labour and Welfare

June 28, 2019

A part of the Japanese Pharmacopoeia (Ministerial Notification No. 64, 2016) was revised as follows*.

(The text referred to by the term "as follows" are omitted here. All of the revised Japanese Pharmacopoeia in accordance with this notification (hereinafter referred to as "new Pharmacopoeia" in Supplement 2) are made available for public exhibition at the Pharmaceutical Evaluation Division, Pharmaceutical Safety and Environmental Health Bureau, Ministry of Health, Labour and Welfare, at each Regional Bureau of Health and Welfare, and at each Prefectural Office in Japan).

Supplementary Provisions

(Effective Date)

Article 1 This Notification is applied from June 28, 2019.

(Transitional measures)

Article 2 In the case of drugs which are listed in the Japanese Pharmacopoeia (hereinafter referred to as "previous Pharmacopoeia") [limited to those listed in new Pharmacopoeia] and drugs which have been approved as of June 28, 2019 as prescribed under Paragraph 1, Article 14 of Act on Securing Quality, Efficacy and Safety of Products Including Pharmaceuticals and Medical Devices [including drugs the Minister of Health, Labour and Welfare specifies (the Ministry of Health and Welfare Ministerial Notification No. 104, 1994) as of June 27, 2019 as those exempted from marketing approval pursuant to Paragraph 1, Article 14 of the same law (hereinafter referred to as "drugs exempted from approval")], the Standards established in the previous Pharmacopoeia (limited to part of the Standards for the drugs concerned) may be accepted to conform to the Standards established in the new Pharmacopoeia before and on December 31, 2020. In the case of drugs which are listed in the new Pharmacopoeia (excluding those listed in the previous Pharmacopoeia) and drugs which have been approved as of June 28, 2019 as prescribed under the Paragraph 1, Article 14 of the same law (including those exempted from approval), they may be accepted as those being not listed in the new Pharmacopoeia before and on December 31, 2020.

(Partial revision of Products Delivery Rules of National Institute of Infectious Diseases)

Article 3 Omitted.

^{*}The term "as follows" here indicates the content of Supplement II to the Japanese Pharmacopoeia Seventeenth Edition from General Notice to Ultraviolet-visible Reference Spectra (pp. 2877 – 2991).

CONTENTS

Preface 1	General Information
Supplement II to The Japanese Pharmacopoeia,	G1 Physics and Chemistry
Seventeenth Edition 2877–2991	Control of Elemental Impurities in
General Notices	Drug Products
General Rules for Preparations 2879	G3 Biotechnological/Biological Products
General Tests, Processes and Apparatus 2881	Host Cell Protein Assay 2997
2.01 Liquid Chromatography	Total Protein Assay 3000
2.26 Raman Spectroscopy 2882	G4 Microorganisms
2.46 Residual Solvents	Media Fill Test (Process Simulation) 3000
2.51 Conductivity Measurement 2891	Microbiological Environmental Moni-
2.66 Elemental Impurities—Procedures 2893	toring Methods of Processing Areas
6.16 Rheological Measurements for Semi-	for Sterile Pharmaceutical Products 3000
solid Preparations	Parametric Release of Terminally Steri-
6.17 Insoluble Particulate Matter Test for	lized Pharmaceutical Products 3000
Therapeutic Protein Injections 2900	G5 Crude Drugs
9.01 Reference Standards 2901	Purity Tests on Crude Drugs using
9.21 Standard Solutions for Volumetric	Genetic Information 3000
Analysis	On the Scientific Names of Crude Drugs
9.22 Standard Solutions	listed in the JP 3005
9.41 Reagents, Test Solutions	G8 Water
9.42 Solid Supports/Column Packings	Quality Control of Water for Pharma-
for Chromatography2908	ceutical Use
000	G10 Others
Official Monographs	Basic Concepts for Quality Assurance
Crude Drugs and Related Drugs 2973	of Drug Substances and Drug Prod-
1 6 1 D 6 G 4 2007 2007	ucts
Infrared Reference Spectra 2985-2987	Concept on Impurities in Chemically
THE 114 111 D 6 C 4 2000 2001	synthesized Drug Substances and
Ultraviolet-visible Reference Spectra 2989–2991	Drug Products
	Glossary for Quality by Design (QbD),
	Quality Risk Management (QRM)
	and Pharmaceutical Quality System
	(PQS)
	International Harmonization Imple-
	mented in the Japanese Pharmaco-
	poeia Seventeenth Edition 3014
	Index 3019
	Index in Latin Name 3039
	Index in Japanese 3041

PREFACE

The 17th Edition of the Japanese Pharmacopoeia (JP) was promulgated by Ministerial Notification No.64 of the Ministry of Health, Labour and Welfare (MHLW) on March 7, 2016.

In July 2016, the Committee on JP established the basic principles for the preparation of the JP 18th Edition, setting out the roles and characteristics of the JP, the definite measures for the revision, and the date of the revision.

It was agreed that the JP should provide an official standard, being required to assure the quality of medicines in Japan in response to the progress of science and technology and medical demands at the time. It should define the standards for specifications, as well as the methods of testing to assure overall quality of all drugs in principle, and it should have a role in clarifying the criteria for quality assurance of drugs that are recognized to be essential for public health and medical treatment.

The JP has been prepared with the aid of the knowledge and experience of many professionals in the pharmaceutical field. Therefore, the JP should have the characteristics of an official standard, which might be widely used by all parties concerned, and it should play an appropriate role of providing information and understanding about the quality of drugs to the public. Moreover, as a pharmaceutical quality standard, it should contribute promoting and maintaining of advancedness as well as international consistency and harmonization of technical requirements in the international community.

At the Committee, the five basic principles of JP, which we refer to as the "five pillars", were established as follows: 1) Including all drugs which are important from the viewpoint of health care and medical treatment; 2) Making qualitative improvement by introducing the latest science and technology; 3) Further promoting internationalization in response to globalization of drug market; 4) Making prompt partial revision as necessary and facilitating smooth administrative operation; and 5) Ensuring transparency regarding the revision, and disseminating the JP to the public. It was agreed that the Committee on JP should make efforts, on the basis of these principles, to ensure that the JP is used more effectively in the fields of health care and medical treatment by taking appropriate measurements, including getting the understanding and cooperation of other parties concerned.

It was also agreed that JP articles should cover drugs, which are important from the viewpoint of health care and medical treatment, clinical performance or merits and frequency of use, as soon as possible after they reach the market.

The target date for the publication of JP 18th Edition (the Japanese edition) was set as April 2021.

JP drafts are discussed in the following committees that were established in the Pharmaceuticals and Medical Devices Agency: Expert Committee; Sub-expert Committee; Sub-committee on Manufacturing Process-related Matters; Committee on Chemicals; Committee on Antibiotics; Committee on Biologicals; Committee on Crude Drugs; Committee on Pharmaceutical Excipients; Committee on Physico-Chemical Methods; Committee on Drug Formulation; Committee on Physical Methods; Committee on Biological Methods; Committee on Nomenclature for Pharmaceuticals; Committee on International Harmonization; and Committee on Reference Standards. Furthermore, working groups are established under the Expert Committee, Committee on Pharmaceutical Excipients, Committee on Physico-Chemical Methods and Committee on Drug Formulation.

The committees initiated deliberations on several revisions.

Draft revisions covering subjects in General Notices, General Rules for Preparations, General Tests, Monographs, Ultraviolet-visible Reference Spectra and Infrared Reference Spectra for which discussions were finished between April 2017 and November 2018, were prepared for a supplement to the JP 17.

Numbers of discussions in the committees to prepare the supplement drafts were as follows: Expert Committee (11, including a working group); Sub-committee on Manufacturing Process-related Matters (6), Committee on Chemicals (18), Committee on Antibiotics (3); Committee on Biologicals (7); Committee on Crude Drugs (15); Committee on Pharmaceutical Excipients (10); Committee on Physico-Chemical Methods (11, including a working group); Committee on Drug Formulation (19, including working groups); Committee on Physical Methods (7); Committee on Biological Methods (6); Committee on Nomenclature for Pharmaceuticals (5); Committee on International Harmonization (4).

It should be noted that in the preparation of the drafts for the supplement, generous cooperation was

given by the Pharmaceutical Technology Committee of the Kansai Pharmaceutical Industries Association, the Pharmacopeia and CMC Committee of the Pharmaceutical Manufacturers' Association of Tokyo, the Tokyo Crude Drugs Association, the International Pharmaceutical Excipients Council Japan, the Home Medicine Association of Japan, the Japan Kampo Medicines Manufacturers Association, the Japan Flavor and Fragrance Materials Association, the Japan Medical Plants Federation, the Japan Pharmaceutical Manufacturers Association, the Federation of Pharmaceutical Manufacturers' Association of Japan, the Parenteral Drug Association Japan Chapter, the Japan Reagent Association, the Japan Oilseeds Processors Association, the Japan Analytical Instruments Manufacturers' Association, and the Asian Society of Innovative Packaging Technology.

The draft revisions were examined by the Committee on JP in January 2019, followed by the Pharmaceutical Affairs and Food Sanitation Council (PAFSC) in March 2019, and then submitted to the Minister of Health, Labour and Welfare. In the Committee on JP, Mitsuru Hashida took the role of chairman from January 2011 to June 2019.

In consequence of this revision, the JP 17th Edition carries 2008 articles, owing to the addition of 34 articles and the deletion of 3 articles.

The principles of description and the salient points of the revision in this volume are as follows:

- 1. The Supplement II to JP 17th Edition comprises the following items, in order: Notification of MHLW; Contents; Preface; General Notices, General Rules for Preparations; General Tests, Processes and Apparatus; Official Monographs; then followed by Infrared Reference Spectra and Ultraviolet-visible Reference Spectra; General Information; and as an appendix a Cumulative Index containing references to the main volume, the Supplement I and the Supplement II.
- 2. The articles in Official Monographs, Infrared Reference Spectra and Ultraviolet-visible Reference Spectra are respectively placed in alphabetical order in principle.
- 3. The following items in each monograph are put in the order shown below, except that unnecessary items are omitted depending on the nature of the drug:
- (1) English title
- (2) Commonly used name(s)
- (3) Latin title (only for crude drugs)
- (4) Title in Japanese

- (5) Structural formula or empirical formula
- (6) Molecular formula and molecular mass
- (7) Chemical name
- (8) Chemical Abstracts Service (CAS) Registry Number
- (9) Origin
- (10) Limits of the content of the ingredient(s) and/or the unit of potency
- (11) Labeling requirements
- (12) Method of preparation
- (13) Manufacture
- (14) Description
- (15) Identification tests
- (16) Specific physical and/or chemical values
- (17) Purity tests
- (18) Potential adulteration
- (19) Loss on drying or Ignition, or Water
- (20) Residue on ignition, Total ash or Acid-insoluble ash
- (21) Tests being required for pharmaceutical preparations
- (22) Other special tests
- (23) Assay
- (24) Containers and storage
- (25) Shelf life
- (26) Others
- 4. In each monograph, the following physical and chemical values representing the properties and quality of the drug are given in the order indicated below, except that unnecessary items are omitted depending on the nature of drug:
- (1) Alcohol number
- (2) Absorbance
- (3) Congealing point
- (4) Refractive index
- (5) Osmolar ratio
- (6) Optical rotation
- (7) Constituent amino acids
- (8) Viscosity
- (9) pH
- (10) Content ratio of the active ingredients
- (11) Specific gravity
- (12) Boiling point
- (13) Melting point
- (14) Acid value
- (15) Saponification value
- (16) Ester value
- (17) Hydroxyl value
- (18) Iodine value
- 5. Identification tests comprise the following items, which are generally put in the order given

below:

- (1) Coloration reactions
- (2) Precipitation reactions
- (3) Decomposition reactions
- (4) Derivatives
- (5) Infrared and/or ultraviolet-visible absorption spectrometry
- (6) Nuclear magnetic resonance spectrometry
- (7) Chromatography
- (8) Special reactions
- (9) Cations
- (10) Anions
- 6. Purity tests comprise the following items, which are generally put in the order given below, except that unnecessary items are omitted depending on the nature of drug:
- (1) Color
- (2) Odor
- (3) Clarity and/or color of solution
- (4) Acidity or alkalinity
- (5) Acidity
- (6) Alkalinity
- (7) Chloride
- (8) Sulfate
- (9) Sulfite
- (10) Nitrate
- (11) Nitrite
- (12) Carbonate
- (13) Bromide
- (14) Iodide
- (15) Soluble halide
- (16) Thiocyanate
- (17) Selenium
- (18) Cationic salts
- (19) Ammonium
- (20) Heavy metals
- (21) Iron
- (22) Manganese
- (23) Chromium
- (24) Bismuth
- (25) Tin
- (26) Aluminum
- (27) Zinc
- (28) Cadmium
- (29) Mercury
- (30) Copper
- (31) Lead
- (32) Silver
- (33) Alkaline earth metals
- (34) Arsenic
- (35) Free phosphoric acid
- (36) Foreign matters

- (37) Related substances
- (38) Isomer
- (39) Optical isomer
- (40) Polymer

(44)

- (41) Residual solvent
- (42) Other impurities
- (43) Residue on evaporation

Readily carbonizable substances

- 7. The following paragraphs of General Notices were revised:
 - (1) Paragraph 5: "Shelf life" in the monographs on preparations were removed from the standards for conformity, because that the stability of a pharmaceutical preparation differs depending on the formulation, type of the container and packaging and the control of the storage temperature, etc.
 - (2) Paragraph 13: As the cases for drug approval with the quality control using the real time-release testing (RTRT) have been accumulated, the rules of RTRT for drug release were added when RTRT is adopted as an alternative method.
 - (3) Paragraph 46: "Expiration Date" in Official Monograph was changed to "Shelf life" in the Supplement I to JP 17th Edition. In this revision, the rule for expression of expiration date was removed.
- **8.** The following paragraphs were newly added to the General Rules for Preparations:
 - (1) The following item was newly added to the [3] Monographs for Preparations basing on the discussion about needs of this item following release of the guideline for developing liposomal formulations:
 - 3-1-4. Liposome Injections
- **9.** The General Rules for Preparations was revised as follows in general:
 - (1) [3] Monographs for Preparations: "3-1. Injections": the definition of "Liposome Injections" introduced in this revision was newly added, and "6.17 Insoluble Particulate Matter Test for Therapeutic Protein Injections" added to General Tests in this revision was referred and prescribed.
- 10. The following items were newly added to General Tests, Processes and Apparatus:
 - (1) 2.26 Raman Spectroscopy
 - (2) 2.66 Elemental Impurities Procedures

- (3) 6.16 Rheological Measurements for Semisolid Preparations
- (4) 6.17 Insoluble Particulate Matter Test for Therapeutic Protein Injections
- 11. The following items in General Tests, Processes and Apparatus were revised:
 - (1) 2.01 Liquid Chromatography
 - (2) 2.46 Residual Solvents
 - (3) 2.51 Conductivity Measurement
- **12.** The following Reference Standards were newly added:

Bromfenac Sodium RS

L-Carnosine RS

Doripenem RS

Gatifloxacin RS

Hydroxyethylcellulose RS for Identification

Lanoconazole RS

Microcrystalline Cellulose RS for Identification

Residual Solvents Class 2C RS

Sitagliptin Phosphate RS

Sitagliptin Phosphate RS for System Suitability

13. The following Reference Standards were removed from "9.01 (2) The reference standards which are prepared by National Institute of Infectious Diseases." and added to "9.01 (1) The reference standards which are prepared by those who have been registered to prepare them by the Minister of Health, Labour and Welfare, according to the Ministerial ordinance established by the Minister separately":

Ampicillin RS

Azithromycin RS

Cefazolin RS

Cefmetazole RS

Fradiomycin Sulfate RS

Meropenem RS

14. The following substances were newly added to the Official Monographs:

Bromfenac Sodium Hydrate

Bromfenac Sodium Ophthalmic Solution

Cefalotin Sodium for Injection

Cefixime Fine Granules

Clarithromycin for Syrup

Diclofenac Sodium Suppositories

Doripenem Hydrate

Doripenem for Injection

Ethylcellulose

Felodipine

Felodipine Tablets

Gatifloxacin Hydrate

Gatifloxacin Ophthalmic Solution

Gentamicin Sulfate Injection

Gentamicin Sulfate Ointment

Hydroxyethylcellulose

Irinotecan Hydrochloride Hydrate

Lanoconazole

Lanoconazole Cream

Lanoconazole Cutaneous Solution

Lanoconazole Ointment

Minocycline Hydrochloride Granules

Nortriptyline Hydrochloride Tablets

Polaprezinc

Polaprezinc Granules

Ritodrine Hydrochloride Injection

Sitagliptin Phosphate Hydrate

Sitagliptin Phosphate Tablets

Sodium Valproate Extended-release Tablets A

Sodium Valproate Extended-release Tablets B

Telmisartan and Hydrochlorothiazide Tablets

Valsartan and Hydrochlorothiazide Tablets

Verapamil Hydrochloride Injection

Goshuyuto Extract

15. The following monographs were revised:

Amphotericin B Tablets

Beclometasone Dipropionate

Betamethasone Dipropionate

Anhydrous Dibasic Calcium Phosphate

Microcrystalline Cellulose

Chloramphenicol

Chlorpromazine Hydrochloride

Cholesterol

Cloperastine Hydrochloride

Dehydrocholic Acid

Purified Dehydrocholic Acid

Epirubicin Hydrochloride

Estriol

Etizolam

Haloperidol

Hydrocortisone

II 1

Hydrocortisone Acetate

Hydrocortisone and Diphenhydramine Ointment

Hydroxypropylcellulose

Hypromellose

Imipramine Hydrochloride

Imipramine Hydrochloride Tablets

Isomalt Hydrate

Mestranol

Methylcellulose

Methyl prednisolone

Pioglitazone Hydrochloride and Glimepiride

Tablets

Saccharin

Saccharin Sodium Hydrate

Light Anhydrous Silicic Acid

Sodium Fusidate Teicoplanin

Testosterone Enanthate Tipepidine Hibenzate Tablets Triamcinolone Acetonide

Ursodeoxycholic Acid

Amomum Seed

Artemisia Capillaris Flower

Belladonna Root Bitter Tincture

Bofutsushosan Extract

Boiogito Extract

Calumba

Powdered Calumba

Cassia Seed

Chrysanthemum Flower

Eucalyptus Oil Gastrodia Tuber Glycyrrhiza Extract Crude Glycyrrhiza Extract Hochuekkito Extract Japanese Angelica Root

Powdered Japanese Angelica Root

Juzentaihoto Extract

Kakkontokasenkyushin'i Extract

Kamikihito Extract Kamishoyosan Extract Lithospermum Root Lycium Bark Nux Vomica

Oriental Bezoar Otsujito Extract Platycodon Root

Powdered Platycodon Root Platycodon Fluidextract

Polygala Root

Powdered Polygala Root

Safflower

Saposhnikovia Root and Rhizome

Scopolia Rhizome

Sinomenium Stem and Rhizome

Swertia Herb

Powdered Swertia Herb

Swertia and Sodium Bicarbonate Powder

Toad Cake

Tokishakuyakusan Extract Yokukansan Extract

16. The following monographs were deleted:

Adsorbed Habu-venom Toxoid

Freeze-dried Tetanus Antitoxin, Equine

Compound Vitamin B Powder

17. The following articles were newly added to Ul-

traviolet-visible Reference Spectra:

Bromfenac Sodium Hydrate

Doripenem Hydrate

Felodipine

Gatifloxacin Hydrate

Irinotecan Hydrochloride Hydrate

Lanoconazole

Sitagliptin Phosphate Hydrate

18. The following articles were newly added to In-

frared Reference Spectra:

Bromfenac Sodium Hydrate

Doripenem Hydrate

Ethyl cellulose

Felodipine

Gatifloxacin Hydrate

Irinotecan Hydrochloride Hydrate

Lanoconazole Polaprezinc

Sitagliptin Phosphate Hydrate

Those who were engaged in the preparation of the

Supplement II to JP 17 are as follows

ABE Yasuhiro
AITA Youhei
AKAO Kenichi
ARATO Teruyo
ARIMOTO Keiko
ARUGA Naoki
ASAMA Hiroshi

ASHIZAWA Kazuhide

ASO Shinichiro
ASO Yukio
DEMIZU Yosuke
EMURA Makoto
FUCHINO Hiroyuki
FUJII Makiko

FUKAZAWA Hidesuke

FUKUDA Shinji FUKUHARA Kiyoshi FURUICHI Harumi FURUKAWA Hiromitsu

GODA Yukihiro GOTO Takashi GOTO Tamami HAISHIMA Yuji HAKAMATA Hideki

HAKAMATSUKA Takashi

HANADA Kentaro HANAJIRI Ruri HARAZONO Akira
HASHIDA Mitsuru*
HASHII Noritaka
HATANO Rika
HAYASHI Ai
HAYASHI Akira
HAYASHI Yoshinori
HIGANO Taro
HIGUCHI Yasuhiko
HIYAMA Yukio
HORI Masatoshi
HYUGA Masashi
ICHINOSE Koji

IKARASHI Yoshiaki IKEDO Shingo IMAMOTO Tsuyoshi ISHIDA Masato ISHII Akiko ITAI Shigeru ITO Michiho ITO Ryoichi ITO Yuii

IIMORI Jumpei

IZUTANI Yusuke
IZUTSU Kenichi
KAMMOTO Toshihiro
KATAYAMA Hirohito
KATO Kumiko

KATORI Noriko KAWAHARA Nobuo KAWAI Tamotsu KAWAMATA Tomomi KAWANISHI Toru

KAWARASAKI Yoshihiko KAWASAKI Nana** KIJIMA Keiji

KIKUCHI Yutaka
KIKUCHI Yuuichi
KISHIMOTO Yasuhiro
KITAJIMA Akihito
KITTAKA Atsushi
KIUCHI Fumiyuki

KOBAYASHI Takashi

KOCHI Rika KOHAMA Ai KOIDE Tatsuo KOJIMA Takashi KOKUBO Hiroyasu KOMATSU Katsuko KONDA Toshifumi

KONDO Seizo KUBOSAKI Atsutaka KUBOTA Kiyoshi KUMASAKA Kenichi KURIHARA Masaaki KUSUNOKI Hideki MAKIURA Toshinobu MARUYAMA Takuro MASADA Sayaka

MATSUMOTO Kazuhiro MATSUMOTO Makoto MATSUMURA Hajime MIKAMI Eiichi MISAWA Takashi MITSUHASHI Takao

MIYATA Naoki
MIYAZAKI Tamaki
MIZUNO Takeshi
MORI Mitsuo
MORIBE Kunikazu
MORIMOTO Seiki
MORIMOTO Takashi
MORISAKI Takahito
MORIYASU Takako
MUROI Masashi
NAKA Nobuyuki

NAKAGAWA Hidehiko NAKAGAWA Shinsaku NAKAGAWA Yukari NAKAKO Mayumi NAKANO Tatsuya NANAURA Mitsuo NASU Masao

NISHIKAWA Noriaki NOGUCHI Shuji **OBA** Sumiaki **OBARA** Sakaé **OGAWA** Toru OGURA Yasumitsu OGURI Kazuki OHGAMI Yasutaka OKUBO Tsuneo **OKUDA** Akihiro OKUDA Haruhiro **OMURA** Koichi **ONO Makoto** ONODA Hiroshi OSUMI Yuko OTSUKA Masami **OUCHI** Tadashi

SAKAI Eiji SAKAMOTO Tomoaki SANTA Tomofumi SASAKI Hiroshi SASAKI Kunio

SAITO Hideyuki

SASAKI Satoshi SASAKI Yuko SATO Koji SATO Kyoko

SATO Norihisa

SEGAWA Masashi SEKIGAWA Fujio

SEKIGUCHI Michiko

SHIBATA Hiroko SHIBAYAMA Keigo

SHIBAZAKI Keiko

SHIDA Shizuka

SHIMADA Yasuo

SHIMAZAWA Rumiko

SHIMIZU Masaro

SHIMOKAWA Sayuri

SHINOHARA Katsuaki

SHIRATORI Makoto

SHIROTA Osamu

SHODA Takuji

SUDO Hirotaka

SUDO Keiichi

SUGIMOTO Naoki

SUZUKI Mikio

SUZUKI Noriyuki

SUZUKI Ryoji

SUZUKI Shigeo

SUZUKI Tomoyuki

TADA Minoru

TADAKI Shinichi

TAKANO Akihito

TAKAO Masaki

TAKEDA Osami

TAKEDA Tomoko

TAKEUCHI Hirohumi

TAKEUCHI Hisashi

TANABE Toyoshige

TANAKA Masakazu

TANAKA Tomohide

TANAKA Satoshi

TANIMOTO Tsuyoshi

TERABAYASHI Susumu

TERADA Katsuhide

TERAOKA Reiko

TOKUMOTO Hiroko

TOKUNAGA Hiroshi

TOKUOKA Shogo

TOMIOKA Kiyoshi

TOYODA Taichi

TSUDA Shigeki

UCHIYAMA Nahoko

UETAKE Atsuhiro

WATANABE Takumi

YAMADA Rumiko
YAMAGUCHI Tetsuji
YAMAJI Hiroki
YAMAMOTO Keiji
YAMASHITA Chikamasa
YASUHARA Masato
YOMOTA Chikako
YONEDA Sachiyo
YONEMOCHI Etsuo
YOSHIDA Hiroyuki
YOSHIDA Naoya

^{*:} Chairman, the Committee on JP

^{**:} Acting Chairman, the Committee on JP

Supplement II to
The Japanese
Pharmacopoeia
Seventeenth Edition

GENERAL NOTICES

Change the paragraphs 5, 13 and 46 as follows:

- 5. The JP Drugs are to be tested according to the provisions given in the pertinent monographs, General Notices, General Rules for Crude Drugs, General Rules for Preparations, and General Tests for their conformity to the Japanese Pharmacopoeia. However, the headings of "Description" and in addition "Containers and storage" and "Shelf life" in the monographs on preparations are given for information, and should not be taken as indicating standards for conformity. Nevertheless, Containers under "Containers and storage" in the monograph on preparations containing crude drugs as main active ingredients are the standards for conformity.
- 13. When an assurance that a product is of the JP Drug quality is obtained consistently from data derived from the manufacturing process validation studies, and from the records of appropriate manufacturing process control and of the test results of the quality control, the performance of some test items in the monograph at release on a product may be omitted as occasion demands. Moreover, the quality evaluation of final products (drug substances and drug products) based on in-process data including in-process testing results and monitoring data on process parameters can replace specifications and test methods in the monograph or performing the test methods, if appropriate.
- **46.** For the JP Drugs, the contents or potency in terms of units of the active ingredient(s) in the monographs have to be shown on the immediate container or wrapping of them.

GENERAL RULES FOR PREPARATIONS

[3] Monographs for Preparations

Change the following paragraphs:

3. Preparations for Injection

3-1. Injections

(1) Injections are sterile preparations to be administered directly into the body through skin, muscle or blood vessel, usually in form of a solution, a suspension or an emulsion of active substance(s), or of a solid that contains active substance(s) to be dissolved or suspended before use.

Parenteral Infusions, Implants/Pellets, Prolonged-Release Injections and Liposome Injections are included in this category.

(14) Unless otherwise specified, Injections and vehicles attached to preparations meet the requirements of Insoluble Particulate Matter Test for Injections $\langle 6.07 \rangle$ or Insoluble Particulate Matter Test for Therapeutic Protein Injections $\langle 6.17 \rangle$.

Add the following next to 3-1-3. Prolonged Release Injections:

3-1-4. Liposome Injections

- (1) Liposome Injections are injections to be used for intravenous administration, which are intended for improvement of in vivo stability, delivery to a target region and control of release, of active substance(s).
- (2) Liposome Injections are usually prepared by using amphipathic lipid, etc. to make aqueous injections or freeze-dried injections in which closed microvesicles composed of a lipid bilayer membrane are dispersed.
- (3) Liposome Injections have an appropriate function of controlled release.
- (4) Liposome Injections have an appropriate particle size.

GENERAL TESTS, PROCESSES AND APPARATUS

Change the introduction to read:

General Tests, Processes and Apparatus includes common methods for tests, useful test methods for quality recognition of drugs and other articles related to them. Unless otherwise specified, acid-neutralizing capacity determination of gastrointestinal medicines, aerodynamic particle size measurement for inhalations, alcohol number determination, amino acid analysis of proteins, ammonium determination, arsenic determination, atomic absorption spectrophotometry, boiling point determination, chloride determination, conductivity measurement, congealing point determination, determination of bulk and tapped densities, digestion test, disintegration test, dissolution test, distilling range determination, elemental impurities procedures, endpoint determination in titrimetry, flame coloration, fluorometry, foreign insoluble matter test for injections, foreign insoluble matter test for ophthalmic liquids and solutions, gas chromatography, glycosylation analysis of glycoprotein, heavy metal determination, inductively coupled plasma-atomic emission spectrometry and inductively coupled plasma-mass spectrometry, infrared spectrophotometry, insoluble particulate matter test for injections, insoluble particulate matter test for ophthalmic liquids and solutions, insoluble particulate matter test for therapeutic protein injections, iron determination, laser diffraction measurement of particle size, liquid chromatography, loss on drying determination, loss on ignition determination, mass spectrometry, melting point determination, methanol determination, methods for color matching, methods of adhesion testing, microbial assay for antibiotics, mineral oil determination, nitrogen determination, nuclear magnetic resonance spectroscopy, optical rotation determination, osmolarity determination, oxygen flask combustion method, particle size determination, particle size distribution test for preparations, pH determination, powder particle density determination, qualitative test, raman spectroscopy, refractive index determination, release test for preparations for cutaneous application, residual solvents, residue on ignition determination, rheological measurements for semi-solid preparations, specific gravity and density determination, specific surface area determination, sulfate determination, test for bacterial endotoxins, test for glass containers for injections, test for metal particles in ophthalmic ointments, test for microbial limit, test for microbial limit for crude drugs, test for plastic containers, test for pyrogen, test for readily carbonizable substances, test for rubber closure for aqueous infusions, test for sterility, test for total organic carbon, test of extractable volume for injection, thermal analysis, thin-layer chromatography, turbidity measurement, ultravioletvisible spectrophotometry, uniformity of delivered dose for inhalations, uniformity of dosage units (test for content uniformity, mass variation test), viscosity determination, vitamin A assay, water determination, and X-ray powder diffraction are performed as directed in the corresponding articles under the General Tests, Processes and Apparatus. The tests for melting point of fats, congealing point of fatty acids, specific gravity, acid value, saponification value, ester value, hydroxyl value, unsaponifiable matter and iodine value of fats and fatty oils are performed as directed in the corresponding items under Fats and Fatty Oils Test, and sampling, preparation of sample for analysis, micro-scopic examination, purity test, loss on drying, total ash, acid-insoluble ash, extract content, essential oil content of crude drugs and assay of marker compounds for the assay of crude drugs and extracts of Kampo Formulations utilizing nuclear magnetic resonance (NMR) spectroscopy are performed as directed in the corresponding items under the Crude Drugs Test.

The number of each test method is a category number given individually. The number in blackets (< >) appeared in monograph indicates the number corresponding to the general test method.

2.01 Liquid Chromatography

Change the Section 7. Point to consider on changing the operating conditions as follows:

7. Point to consider on changing the operating conditions

Among the operating conditions specified in the individual monograph, inside diameter and length of the column, particle size of the packing material (pore size in the case of monolithic columns), column temperature, composition ratio of the mobile phase, composition of buffer solutions in the mobile phase, pH of the mobile phase, concentration of ion-pair forming agents in the mobile phase, ionic strength of the mobile phase, flow rate of the mobile phase, number and timing of mobile phase composition changes in gradient program, flow rate of mobile phase in gradient program, composition and flow rate of derivatizing reagents, and reaction time and chamber temperature in chemical reaction may be modified within the ranges in which the liquid chromatographic system used conforms to the requirements of system suitability.

Add the following:

2.26 Raman Spectroscopy

Raman spectroscopy is a vibrational spectroscopic technique, which evaluates a sample to be examined qualitatively or quantitatively by analyzing a spectrum obtained by dispersing very weak scattered light, having different wavelengths from irradiation light, generated when the sample is irradiated with the light. Raman scattering is observed when the polarizability of molecules changes with the vibration of chemical bonds of molecules in a sample.

Raman spectroscopy generally uses monochromatic laser light as excitation light. When the laser light is irradiated to the sample to be examined, the molecules in the sample is excited and the light with the same wavelength of the irradiation light, known as Rayleigh scattering, is scattered. The scattered light detected in the shorter wavelength side than the Rayleigh scattering is referred to as anti-Stokes scattering. The scattered light detected in the longer wavelength side than the Rayleigh scattering is referred to as Stokes scattering. Generally Stokes scattering with strong Raman scattering intensity is used for analysis. Raman spectra are usually indicated by Raman shift on the horizontal axis and Raman scattering intensity on the vertical axis.

Raman spectroscopy is capable of measuring samples (solid, semi-solid, liquid, gas, etc.) rapidly and nondestructively without pre-treatment. Application of Raman spectroscopy in the pharmaceutical field includes qualitative or quantitative evaluation of the active pharmaceutical ingredients and additives in drug substances or drug products. Raman spectroscopy can also be used for the evaluation of the physical conditions of substances, such as crystal form and crystallinity. Raman micro-spectroscopy can also be used for the evaluation of the distributions of active pharmaceutical ingredients and additives in the drug products. Furthermore, using an optical fiber probe enables it to measure the spectra of samples at a location remote from the equipment body without sampling, so that it can be used to perform pharmaceutical manufacturing process control online (or in-line).

1. Apparatus

Raman spectrometers are composed of a light source unit, a sample unit, a spectrometry unit, a detector unit, a signal processing unit, a data processing unit and a display-record-output unit. Raman spectrometers are classified into dispersive Raman spectrometers and Fourier transform Raman spectrometers according to their spectroscopic systems.

1.1. Light source

The laser which stably emits monochromatic light as excitation light to samples is used for the light source. The lasers include gas lasers such as a He-Ne laser and solid state lasers, and select a laser with wavelengths and output power according to the purpose. Pay attention to safety standards relating to a laser, when this test is performed.

1.2. Sample unit

The sample unit is composed of an optical system for col-

lecting Raman scattering light generated by irradiation of excitation light and a sample cell. Combination of these take the form of a sample chamber, while there are apparatuses with no sample chamber, such as optical fiber probes and portable Raman spectrometers that can be carried. Representative sample chambers are macroscopic sample chambers and microscopic sample chambers. The components of these optical systems are different, respectively.

1.3. Monochromator and detector

Many dispersive Raman spectrometers use an optical filter to eliminate excitation light and use a single monochromator combined with a multichannel detector, since the configuration is simple and high sensitivity can be obtained. Detectors include multiple elements detectors and single element detectors, and general dispersive Raman spectrometers use a multiple elements detector such as a CCD detector.

Fourier transform (FT) Raman spectrometers obtain spectra by Fourier transformation of interference waveforms using an interferometer. FT-Raman spectrometers are mainly used for near infrared Raman measurement.

2. Methods Used for Measurement

The Raman spectroscopy is applicable to solid samples having a complicated shape in addition to gas/solution samples inside a glass sample cell being transparent in the visible region, using mainly light in the visible region as excitation light. In the view of the size of a measurement region and Raman scattering efficiency, an optimum optical system is selected according to the sample. The excitation wavelength, the measurement mode of the apparatus, etc. are selected and set.

2.1. Macroscopic measurement

Since the macroscopic sample chamber has a high degree of freedom in the scattering configuration, samples can be measured irrespective of solid, liquid, gas, size and shape. It is also applicable to Raman measurement under low temperature, high temperature and high pressure which require the setting of a large sample cell. Usually, in the macroscopic sample chamber three configurations: forward scattering (transmission), 90° scattering and back scattering configurations, can be usable and an appropriate scattering configuration can be selected depending on a sample.

2.2. Microscopic measurement

The microscopic sample chamber is based on an optical microscope and applicable to local analysis. In the optical system of the microscopic sample chamber a microscope objective lens works simultaneously as an excitation light converging lens and a Raman scattering light condensing lens.

Mapping measurement repeats local measurements by moving a sample or laser light position to generate a Raman image showing the two or three dimensional distribution of Raman scattering intensity. Raman images are made by using various spectral information such as a ratio of the intensity of two bands.

2.3. Probe measurement

The optical fiber probe is the collective term of the apparatus of which sample section is separated from a Raman spectrometer body by using an optical fiber and is applica-

ble to *in situ* measurement and on-line (or in-line) measurement.

2.4. Measurement by portable apparatus

The portable Raman spectrometer is possible to carry and perform analysis using Raman spectroscopy outside of laboratory. Main application of this apparatus is judgement on acceptance of pharmaceutical materials. It is used for rather simple measurement.

2.5. Points to note in measurement

Note the following points for solid, liquid and suspended samples.

- (i) Measurement of solid sample: There is a possibility that the filling status, the difference in the particle diameter and the roughness of the surface of the sample could affect the scattering intensity. When measuring a crystalline sample, be careful about the effect of crystal shape. There is also a possibility that the light transparency of the sample affect the spectrum intensity. When a sample is physically and chemically inhomogeneous, it might be recommended to enlarge the spot size of laser irradiation, measure plural samples, measure the plural points of the same sample or crush the sample to homogenize.
- (ii) Measurement of liquid sample: It is possible to subtract the spectrum of the solvent if there is no interaction between the solvent and the sample. When there are insoluble matters in solution, remove the matters using a filter before measurement not to obtain the Raman scattering of the matters. When a sample shows high reactivity by laser irradiation in solution, measure the sample by stirring carefully not to irradiate the same place.
- (iii) Measurement of suspended sample: A suspended sample may settle, so be careful about the positioning of laser irradiation. For samples that are prone to settle, devising measurement such as optimizing the irradiation time and stirring might be helpful. When the Raman scattering of a suspended sample is weak, it is also possible to subtract the spectrum of the solvent likewise the case of measurement of a liquid sample.

3. Factors that affect spectrum

When Raman spectroscopy is applied, note the following items as factors affecting spectra.

3.1. Temperature of sample

Sample heating by laser irradiation can cause a variety of effects, such as physical form change (melting and burning) and polymorph transform. Since the chance of sample heating is increased when the spot size of laser irradiation at the sample is squeezed, be careful not to damage sample when microscopic measurement is carried out. To prevent the sample overheating, a variety of methods can be employed such as suppressing laser output, irradiating a laser without focusing and cooling a sample.

3.2. Sample characteristics

Since Raman signals are very weak, the fluorescence of a sample itself and minute impurities may interfere with Raman scattering light. Fluorescence can be reduced by choosing an excitation light source with a longer wavelength, however it should be noted that it generally decreases the

intensity of the Raman scattering. Photobleaching resulted by laser irradiation before measurement, appropriate irradiation time and accumulation count may mitigate the fluorescence.

When measuring a colored sample, select the wavelength of an excitation laser depending on the absorption characteristics of the sample. When measuring a sample in a container such as, a cell for measurement, a bag or a bottle, take careful note of the spectral characteristics derived from the container in addition to the sample.

4. Control of apparatus performance

Estimate the accuracy of the wave number of Raman shift after adjusting a Raman spectrometer. Measure Raman spectra using an excitation laser utilized for actual measurement and an appropriate standard substance. Polystyrene is an example.

In the cases of 2.1., 2.2. and 2.3., make correction using at least three wave numbers among the below peak wave numbers (cm⁻¹) obtained from the spectrum of polystyrene. The number in parentheses indicates the permissible range.

```
620.9 (± 1.5)
1001.4 (± 1.5)
1031.8 (± 1.5)
1602.3 (± 1.5)
```

3054.3 (\pm 3.0) (Note: 3054.3 cm⁻¹ cannot be measured depending on an excitation wavelength.)

In the case of 2.4., make correction in the same manner using at least three wave numbers among the below peak wave numbers (cm⁻¹).

```
620.9 (± 2.5)
1001.4 (± 2.0)
1031.8 (± 2.0)
1602.3 (± 3.0)
```

Other substance such as cyclohexane can be used as a standard substance, if it is validated.

5. Qualitative and quantitative analysis

5.1. Qualitative analysis

As Raman spectroscopy observe the vibrational energy of a molecule and can obtain a characteristic spectrum depending on the structure of a substance to be analyzed, qualitative analysis based on chemical structural information can be performed.

When the Raman spectra of a sample and the Reference Standard of a substance to be identified are compared and both spectra exhibit similar scattering intensities at the same Raman shifts, the identity of those can be confirmed.

When a sample treatment method for a solid sample is indicated in the monograph in the case of nonconformity of the scattering spectrum with that of the Reference Standard, treat the sample and the Reference Standard under the condition as directed in the monograph, then repeat the measurement.

When the characteristic scattering wave numbers of a substance to be identified are specified in the monograph, the clear appearance of the scattering of a sample at all the specified scattering wave numbers can confirm the identity of the sample with the substance to be identified.

Raman spectroscopy is also applicable to the process control of drug substances or drug products by using a score obtained from a Raman spectrum by a chemometric methodology such as principal component analysis, and characteristic peak wave numbers of the substance to be examined, as indices. Chemometrics usually means mathematical technique and statistic technique for quantization and informatization of chemical data.

5.2. Quantitative analysis

The concentrations of components of a sample to be analyzed can be calculated by calibration curves plotting the relationship between scattering intensity at a specified wave number and concentration. In the case where the composition of components of a sample is complicated, the concentrations of components in the sample can also be calculated by developing a calibration model about a spectrum measured using an existing standard sample by a chemometric methodology and applying the model to the spectrum of the sample to be examined. Chemometric methodologies for obtaining a calibration model include multiple regression analysis method and PLS (Partial least squares) regression analysis method.

The variation of peak intensity at around the reference values of wave numbers using a standard sample, polystyrene etc. used in 4., is preferable to be within \pm 10% compared to that obtained in the last measurement.

2.46 Residual Solvents

Change as follows:

The chapter of residual solvents describes the control, identification and quantification of organic solvents remaining in drug substances, excipients and drug products.

I. Control of residual solvents

1. Introduction

Residual solvents in pharmaceuticals (except for crude drugs and their preparations) are defined here as organic volatile chemicals that are used or produced in the manufacture of drug substances or excipients, or in the preparation of drug products. The solvents are not completely removed by practical manufacturing techniques. Appropriate selection of the solvent for the synthesis of drug substance may enhance the yield, or determine characteristics such as crystal form, purity, and solubility. Therefore, the solvent may sometimes be a critical parameter in the synthetic process. The test method described in this chapter does not address solvents deliberately used as excipients nor does it address solvates. However, the content of solvents in such products should be evaluated and justified.

Since there is no therapeutic benefit from residual solvents, all residual solvents should be reduced to the extent possible to meet product specifications, good manufacturing practices, or other quality-based requirements. Drug products should contain no higher levels of residual solvents than can be supported by safety data. Some solvents that

are known to cause unacceptable toxicities (Class 1, Table 2.46-1) should be avoided in the production of drug substances, excipients, or drug products unless their use can be strongly justified in a risk-benefit assessment. Some solvents associated with less severe toxicity (Class 2, Table 2.46-2) should be limited in order to protect patients from potential adverse effects. Ideally, less toxic solvents (Class 3, Table 2.46-3) should be used where practical.

Testing should be performed for residual solvents when production or purification processes are known to result in the presence of such solvents. It is only necessary to test for solvents that are used or produced in the manufacture or purification of drug substances, excipients, or drug products. Although manufacturers may choose to test the drug product, a cumulative method may be used to calculate the residual solvent levels in the drug product from the levels in the ingredients used to produce the drug product. If the calculation results in a level equal to or below that recommended in this chapter, no testing of the drug product for residual solvents needs to be considered. If, however, the calculated level is above the recommended level, the drug product should be tested to ascertain whether the formulation process has reduced the relevant solvent level to within the acceptable amount. Drug product should also be tested if a solvent is used during its manufacture.

The limit applies to all dosage forms and routes of administration. Higher levels of residual solvents may be acceptable in certain cases such as short term (30 days or less) or topical application. Justification for these levels should be made on a case by case basis.

2. General principles

2.1. Classification of residual solvents by risk assessment

The term "PDE" (Permitted Daily Exposure) is defined in this chapter as a pharmaceutically acceptable daily intake of residual solvents. Residual solvents assessed in this chapter were evaluated for their possible risk to human health and placed into one of three classes as follows:

(i) Class 1 solvents: Solvents to be avoided in manufacture of pharmaceuticals

Known human carcinogens, strongly suspected human carcinogens, and environmental hazards. Class 1 solvents are listed in Table 2.46-1.

(ii) Class 2 solvents: Solvents to be limited in pharmaceuticals

Non-genotoxic animal carcinogens or possible causative agents of other irreversible toxicity such as neurotoxicity or teratogenicity. Solvents suspected of other significant but reversible toxicities. Class 2 solvents are listed in Table 2.46-2.

(iii) Class 3 solvents: Solvents with low toxic potential Solvents with low toxic potential to human; no health-based exposure limit is needed. Class 3 solvents are listed in Table 2.46-3 and have PDEs of 50 mg or more per day.

2.2. Option for describing limits of Class 2 solvents

Two options are available when setting limits for Class 2 solvents.

2.2.1. Option 1

The concentration limits in ppm can be calculated using equation (1) below by assuming a product mass of 10 g administered daily.

Concentration limit (ppm) =
$$\frac{1000 \times PDE}{\text{dose}}$$
 (1)

Here, PDE is given in terms of mg per day and dose is given in g per day.

These limits are considered acceptable for all substances, excipients, or products. Therefore this option may be applied if the daily dose is not known or fixed. If all excipients and drug substances in a formulation meet the limits given in Option 1, then these components may be used in any proportion. No further calculation is necessary provided the daily dose does not exceed 10 g. Products that are administered in doses greater than 10 g per day should be considered under Option 2.

2.2.2. Option 2

It is not considered necessary for each component of the drug product to comply with the limits given in Option 1. The PDE in terms of mg per day as stated in Table 2.46-2 can be used with the known maximum daily dose and equation (1) above to determine the concentration of residual solvent allowed in drug product. Such limits are considered acceptable provided that it has been demonstrated that the residual solvent has been reduced to the practical minimum. The limits should be realistic in relation to analytical precision, manufacturing capability, reasonable variation in the manufacturing process, and the limits should reflect contemporary manufacturing standards.

Option 2 may be applied by adding the amounts of a residual solvent present in each of the components of the drug product. The sum of the amounts of solvent per day should be less than that given by the PDE.

3. Analytical procedures

Residual solvents are typically determined using chromatographic techniques such as gas chromatography. If only Class 3 solvents are present, a nonspecific method such as loss on drying may be used. The analytical method should be validated adequately.

4. Reporting levels of residual solvents

Manufacturers of drug products need certain information about the content of residual solvents in excipients or drug substances. The following statements are given as acceptable examples of the information.

- (i) Only Class 3 solvents are likely to be present. Loss on drying is not more than 0.5%.
- (ii) Only Class 2 solvents are likely to be present. Name the Class 2 solvents that are present. All are not more than the Option 1 limit.
- (iii) Only Class 2 solvents and Class 3 solvents are likely to be present. Residual Class 2 solvents are not more than the Option 1 limit and residual Class 3 solvents are not more than 0.5%.

If Class 1 solvents are likely to be present, they should be identified and quantified. "Likely to be present" refers to the solvents that were used in the final manufacturing step

and to the solvents that were used in earlier manufacturing steps and not always possible to be excluded even in a validated process.

If solvents of Class 2 or Class 3 are present at greater than their Option 1 limits or 0.5%, respectively, they should be identified and quantified.

5. Limits of residual solvents

5.1. Solvents to be avoided in manufacture of pharmaceuticals

Solvents in Class 1 should not be employed in the manufacture of drug substances, excipients, and drug products because of their unacceptable toxicity or their deleterious environmental effect. However, if their use is unavoidable in order to produce a drug product with a significant therapeutic advance, then their levels should be restricted as shown in Table 2.46-1, unless otherwise justified. 1,1,1-Trichloroethane is included in Table 2.46-1 because it is an environmental hazard. The stated limit of 1500 ppm shown in Table 2.46-1 is based on a review of the safety data.

Table 2.46-1 Class 1 solvents in drug products (solvents that should be avoided)

Solvent	Concentration Limit (ppm)	Concern
Benzene	2	Carcinogen
Carbon tetrachloride	4	Toxic and environ-
		mental hazard
1,2-Dichloroethane	5	Toxic
1,1-Dichloroethene	8	Toxic
1,1,1-Trichloroethane	1500	Environmental
		hazard

5.2. Solvents to be limited in pharmaceuticals

Solvents in Table 2.46-2 should be limited in drug products because of their inherent toxicity.

PDEs are given to the nearest 0.1 mg per day, and concentrations are given to the nearest 10 ppm. The stated values do not reflect the necessary analytical precision of determination. Precision should be determined as part of the validation of the method.

5.3. Solvents with low toxic potential

Solvents in Class 3 shown in Table 2.46-3 may be regarded as less toxic and of lower risk to human health. Class 3 includes no solvent known as a human health hazard at levels normally accepted in pharmaceuticals. The amounts of these residual solvents of 50 mg per day or less (corresponding to 5000 ppm or 0.5% under Option 1) would be acceptable without justification. Higher amounts may also be acceptable provided they are realistic in relation to manufacturing capability and good manufacturing practice.

5.4 Solvents for which no adequate toxicological data was found

The following solvents (Table 2.46-4) may also be of interest to manufacturers of drug substances, excipients, or drug products. However, no adequate toxicological data on

Table 2.46-2 Class 2 Solvents which residual amount should be limited in drug products

should be infilted in drug products			
Solvent	PDE (mg/day)	Concentration limit (ppm)	
Acetonitrile	4.1	410	
Chlorobenzene	3.6	360	
Chloroform	0.6	60	
Cumene	0.7	70	
Cyclohexane	38.8	3880	
1,2-Dichloroethene	18.7	1870	
Dichloromethane	6.0	600	
1,2-Dimethoxyethane	1.0	100	
N,N-Dimethylacetamide	10.9	1090	
N,N-Dimethylformamide	8.8	880	
1,4-Dioxane	3.8	380	
2-Ethoxyethanol	1.6	160	
Ethylene glycol	3.1	310	
Formamide	2.2	220	
Hexane	2.9	290	
Methanol	30.0	3000	
2-Methoxyethanol	0.5	50	
Methyl butyl ketone	0.5	50	
Methylcyclohexane	11.8	1180	
Methyl isobutyl ketone	45	4500	
N-Methylpyrrolidone	5.3	530	
Nitromethane	0.5	50	
Pyridine	2.0	200	
Sulfolane	1.6	160	
Tetrahydrofuran	7.2	720	
Tetralin	1.0	100	
Toluene	8.9	890	
1,1,2-Trichloroethene	0.8	80	
Xylene*	21.7	2170	

^{*} Usually 60% *m*-xylene, 14% *p*-xylene, 9% *o*-xylene with 17% ethylbenzene

Table 2.46-3 Class 3 solvents which should be limited by GMP or other quality-based requirements

	Acetic acid	Heptane
	Acetone	Isobutyl acetate
	Anisole	Isopropyl acetate
	1-Butanol	Methyl acetate
	2-Butanol	3-Methyl-1-butanol
	n-Butyl acetate	Methyl ethyl ketone
	tert-Butyl methyl ether	2-Methyl-1-propanol
	Dimethylsulfoxide	Pentane
	Ethanol	1-Pentanol
	Ethyl acetate	1-Propanol
	Diethyl ether	2-Propanol
	Ethyl formate	Propyl acetate
	Formic acid	Triethylamine
-		

which to base a PDE was found. Manufacturers should supply justification for residual levels of these solvents in

Table 2.46-4 Solvents for which no adequate toxicological data was found

1,1-Diethoxypropane 1,1-Dimethoxymethane 2,2-Dimethoxypropane	Methyl isopropyl ketone Methyltetrahydrofuran Petroleum ether
Isooctane Isopropyl ether	Trichloroacetic acid Trifluoroacetic acid

drug products.

II. Identification and quantification of residual solvents

Whenever possible, the substance under test needs to be dissolved to release the residual solvent. Because drug products, as well as active ingredients and excipients are treated, it may be acceptable that in some cases, some of the components of the formulation will not dissolve completely. In those cases, the drug product may first need to be pulverized into a fine powder so that any residual solvent that may be present can be released. This operation should be performed as fast as possible to prevent the loss of volatile solvents during the procedure.

In the operating conditions of gas chromatography and headspace described below, parameters to be set and their description may be different depending on the apparatus. When setting these conditions, it is necessary to change them according to the apparatus used, if it is confirmed that they meet the system suitability.

In addition to the reagents specified to be used for the test, those that meet the purpose of the test can be used.

1. Class 1 and Class 2 residual solvents

The following procedures are useful to identify and quantify residual solvents when the information regarding which solvents are likely to be present in the material is not available. When the information about the presence of specific residual solvents are available, it is not necessary to perform Procedure A and Procedure B, and only Procedure C or other appropriate procedure is needed to quantify the amount of residual solvents.

A flow chart for the identification of residual solvents and the application of limit and quantitative tests is shown in Fig. 2.46-1.

1.1. Water-soluble articles

1.1.1. Procedure A

The test is performed by gas chromatography $\langle 2.02 \rangle$ according to the following conditions.

Class 1 standard stock solution: Pipet 1 mL of Residual Solvents Class 1 RS, dissolve in about 9 mL of dimethylsulfoxide, and add water to make exactly 100 mL. Pipet 1 mL of this solution in a volumetric flask, previously filled with about 50 mL of water and add water to make exactly 100 mL. Pipet 10 mL of this solution in a volumetric flask, previously filled with about 50 mL of water and add water to make exactly 100 mL.

Class 1 standard solution: Pipet 1 mL of Class 1 standard stock solution in an appropriate head space vial containing exactly 5 mL of water, apply the stopper, cap, and mix.




Fig. 2.46-1 Flow chart for the identification of residual solvents and the application of limit and qualification tests.

Class 2 standard stock solution A: Pipet 1 mL of Residual Solvents Class 2A RS, add water to make exactly 100 mL.

Class 2 standard stock solution B: Pipet 1 mL of Residual Solvents Class 2B RS, add water to make exactly 100 mL.

Class 2 standard stock solution C: Pipet 1 mL of Residual Solvents Class 2C RS, add water to make exactly 100 mL.

Class 2 standard solution A: Pipet 1 mL of Class 2 standard stock solution A in an appropriate headspace vial, add exactly 5 mL of water, apply the stopper, cap, and mix.

Class 2 standard solution B: Pipet 5 mL of Class 2 standard stock solution B in an appropriate headspace vial, add exactly 1 mL of water, apply the stopper, cap, and mix.

Class 2 standard solution C: Pipet 1 mL of Class 2 standard stock solution C in an appropriate headspace vial, add exactly 5 mL of water, apply the stopper, cap, and mix.

Test stock solution: Dissolve 0.25 g of the article under test in water, and add water to make exactly 25 mL.

Test solution: Pipet 5 mL of test stock solution in an appropriate headspace vial, add exactly 1 mL of water, and apply the stopper, cap, and mix.

Class 1 system suitability solution: Pipet 1 mL of Class 1 standard stock solution in an appropriate headspace vial, add exactly 5 mL of test stock solution, and apply the stopper, cap, and mix.

Operating conditions—

Detector: A hydrogen flame-ionization detector.

Column: A fused silica column (or a wide-bore column) 0.32 mm (or 0.53 mm) in inside diameter and 30 m in length, coated with 6% cyanopropylphenyl-94% dimethyl silicon polymer for gas chromatography in 1.8 μ m (or 3.0 μ m) thickness.

Column temperature: Maintain at 40°C for 20 minutes after injection, raise the temperature to 240°C at a rate of 10°C per minute, and maintain at 240°C for 20 minutes.

Injection port temperature: 140°C. Detector temperature: 250°C. Carrier gas: Nitrogen or Helium. Flow rate: About 35 cm per second.

Split ratio: 1:5. (Note: The split ratio can be modified in order to optimize sensitivity.)

System suitability-

Test for required detectability: When the procedure is run with Class 1 standard solution and Class 1 system suitability solution under the above operating conditions, the SN ratio of the peak of 1,1,1-trichloroethane in Class 1 standard solution is not less than 5, and the SN ratio of each peak in Class 1 system suitability solution is not less than 3, respectively.

System performance: When the procedure is run with Class 2 standard solution A or the solution for system suitability under the above operating conditions, the resolution between acetonitrile and dichloromethane is not less than 1.0. Pipet 1 mL of a solution of Residual Solvents RS for System Suitability (1 in 100) in an appropriate headspace vial, add exactly 5 mL of water, and apply the stopper, cap, and mix, and use this solution as the solution for system suitability.

System repeatability: When the test is repeated 6 times with Class 1 standard solution under the above operating conditions, the relative standard deviation of each peak area is not more than 15%.

Separately inject (following one of the headspace operating parameter sets described in Table 2.46-5) equal volumes of headspace (about 1.0 mL) of Class 1 standard solution, Class 2 standard solution A, Class 2 standard solution B, Class 2 standard solution C and test solution into the chromatograph, record the chromatograms, and measure the responses for the major peaks. If a peak response of any peak, other than a peak for 1,1,1-trichloroethane, in the test solution is greater than or equal to a corresponding peak in either Class 1 standard solution, Class 2 standard solution A, Class 2 standard solution B or Class 2 standard solution C, or a peak response of 1,1,1-trichloroethane is greater than or equal to 150 times the peak response corresponding to 1,1,1-trichloroethane in Class 1 standard solution, proceed to Procedure B to verify the identity of the peak; otherwise the article meets the requirements of this test.

1.1.2. Procedure B

The test is performed by gas chromatography $\langle 2.02 \rangle$ according to the following conditions.

Class 1 standard stock solution, Class 1 standard solution, Class 1 system suitability solution, Class 2 standard

stock solutions A, B and C, Class 2 standard solutions A, B and C, test stock solution and test solution: Prepare as directed for Procedure A.

Operating conditions—

Detector: A hydrogen flame-ionization detector.

Column: A fused silica column (or a wide-bore column) 0.32 mm (or 0.53 mm) in inside diameter and 30 m in length, coated with polyethylene glycol for gas chromatography in $0.25 \, \mu \text{m}$ thickness.

Column temperature: Maintain at 50°C for 20 minutes after injection, raise the temperature to 165°C at a rate of 6°C per minute, and maintain at 165°C for 20 minutes.

Injection port temperature: 140°C. Detector temperature: 250°C.

Carrier gas: Nitrogen or Helium. Flow rate: About 35 cm per second.

Split ratio: 1:5. (Note: The split ratio can be modified in order to optimize sensitivity.)

System suitability—

Test for required detectability: When the procedure is run with Class 1 standard solution and Class 1 system suitability solution under the above operating conditions, the SN ratio of the peak of benzene in Class 1 standard solution is not less than 5, and the SN ratio of each peak in Class 1 system suitability solution is not less than 3, respectively.

System performance: When the procedure is run with Class 2 standard solution A or the solution for system suitability under the above operating conditions, the resolution between acetonitrile and *cis*-1,2-dichloroethene is not less than 1.0. Pipet 1 mL of a solution of Residual Solvents RS for System Suitability (1 in 100) in an appropriate headspace vial, add exactly 5 mL of water, and apply the stopper, cap, and mix, and use this solution as the solution for system suitability.

System repeatability: When the test is repeated 6 times with Class 1 standard solution under the above operating conditions, the relative standard deviation of each peak area is not more than 15%.

Separately inject (following one of the headspace operating parameter sets described in Table 2.46-5) equal volumes of headspace (about 1.0 mL) of Class 1 standard solution, Class 2 standard solution A, Class 2 standard solution B, Class 2 standard solution C and test solution into the chromatograph, record the chromatograms, and measure the responses for the major peaks. If the peak response(s) in the test solution of the peak(s) identified in Procedure A is/are greater than or equal to a corresponding peak(s) in either Class 1 standard solution, Class 2 standard solution A, Class 2 standard solution B or Class 2 standard solution C, proceed to Procedure C to quantify the peak(s); otherwise the article meets the requirements of this test.

1.1.3. Procedure C

The test is performed by gas chromatography $\langle 2.02 \rangle$ according to the following conditions.

Class 1 standard stock solution, Class 1 standard solution, Class 2 standard stock solution A, Class 2 standard solution C, Class 2 standard solution C and Class 1 system suitability solution:

Prepare as directed for Procedure A.

Standard stock solution (Note: Prepare a separate standard stock solution for each peak identified and verified by Procedures A and B. For Class 1 solvents other than 1,1,1-trichloroethane, prepare the first dilution as directed for the first dilution under Class 1 standard stock solution in Procedure A.): Transfer an accurately measured volume of each individual solvent corresponding to each residual solvent peak identified and verified by Procedures A and B to a suitable container, and dilute quantitatively and stepwise if necessary, with water to obtain a solution having a final concentration of 1/20 of the value stated in Table 2.46-1 or Table 2.46-2.

Standard solution: Pipet 1 mL of standard stock solution in an appropriate headspace vial, add exactly 5 mL of water, and apply the stopper, cap, and mix.

Test stock solution: Weigh accurately about 0.25~g of the article under test, dissolve in water, and add water to make exactly 25~mL.

Test solution: Pipet 5 mL of test stock solution in an appropriate headspace vial, add exactly 1 mL of water, and apply the stopper, cap, and mix.

Spiked test solution (Note: prepare a separate spiked test solution for each peak identified and verified by Procedure A and B): Pipet 5 mL of test stock solution in an appropriate headspace vial, add exactly 1 mL of standard stock solution, and apply the stopper, cap, and mix.

Operating conditions and system suitability fundamentally follow the procedure A. Test for required detectability is unnecessary, and use Standard solution instead of Class 1 standard solution for system repeatability. If the results of the chromatography from Procedure A are found to be inferior to those found with Procedure B, the operating conditions from Procedure B may be substituted.

Separately inject (following one of the headspace operating parameters described in Table 2.46-5) equal volumes of headspace (about 1.0 mL) of the standard solution, test solution, and spiked test solution into the chromatograph, record the chromatograms, and measure the responses for the major peaks. Calculate the amount of each residual solvent found in the article under test by the formula:

Residual solvent (ppm) = 5 $(C/M) \{A_T/(A_S - A_T)\}$

- C: Concentration (µg/mL) of the appropriate Reference Standard in the standard stock solution
- M: Amount (g) of the article under test taken to prepare the test stock solution
- $A_{\rm T}$: Peak responses of each residual solvent obtained from the test solution
- A_S: Peak responses of each residual solvent obtained from the spiked test solution

1.2. Water-insoluble articles

1.2.1. Procedure A

The test is performed by gas chromatography $\langle 2.02 \rangle$ according to the following conditions. Dimethylsulfoxide may be substituted as an alternative solvent to N,N-dimethylformamide.

Class 1 standard stock solution: Pipet 1 mL of Residual Solvents Class 1 RS, dissolve in about 80 mL of N,N-dimethylformamide, and add N,N-dimethylformamide to make exactly 100 mL. Pipet 1 mL of this solution in a volumetric flask, previously filled with about 80 mL of N,N-dimethylformamide and add N,N-dimethylformamide to make exactly 100 mL (this solution is the intermediate dilution prepared from Residual Solvents Class 1 RS and use it for preparation of Class 1 system suitability solution). Pipet 1 mL of this solution, and add N,N-dimethylformamide to make exactly 10 mL.

Class 1 standard solution: Pipet 1 mL of Class 1 standard stock solution in an appropriate head space vial containing exactly 5 mL of water, apply the stopper, cap, and mix.

Class 2 standard stock solution A: Pipet 1 mL of Residual Solvents Class 2A RS, dissolve in about 80 mL of *N*,*N*-dimethylformamide, and add *N*,*N*-dimethylformamide to make exactly 100 mL.

Class 2 standard stock solution B: Pipet 0.5 mL of Residual Solvents Class 2B RS, add *N*,*N*-dimethylformamide to make exactly 10 mL.

Class 2 standard stock solution C: Pipet 1 mL of Residual Solvents Class 2C RS, dissolve in about 80 mL of *N*,*N*-dimethylformamide, and add *N*,*N*-dimethylformamide to make exactly 100 mL.

Class 2 standard solution A: Pipet 1 mL of Class 2 standard stock solution A in an appropriate headspace vial, add exactly 5 mL of water, and apply the stopper, cap, and mix.

Class 2 standard solution B: Pipet 1 mL of Class 2 standard stock solution B in an appropriate headspace vial, add exactly 5 mL of water, and apply the stopper, cap, and mix.

Class 2 standard solution C: Pipet 1 mL of Class 2 standard stock solution C in an appropriate headspace vial, add exactly 5 mL of water, and apply the stopper, cap, and mix.

Test stock solution: Dissolve 0.5 g of the article under test in N,N-dimethylformamide, and add N,N-dimethylformamide to make exactly 10 mL.

Test solution: Pipet 1 mL of test stock solution in an appropriate headspace vial, add exactly 5 mL of water, and apply the stopper, cap, and mix.

Class 1 system suitability solution: Pipet 5 mL of test stock solution and 0.5 mL of the intermediate dilution prepared from Residual Solvents Class 1 RS, and mix. Pipet 1 mL of this solution in an appropriate headspace vial, add exactly 5 mL of water, and apply the stopper, cap, and mix. Operating conditions—

Detector: A hydrogen flame-ionization detector.

Column: A wide-bore column 0.53 mm in inside diameter and 30 m in length, coated with 6% cyanopropylphenyl-94% dimethyl silicon polymer for gas chromatography in $3.0 \, \mu \text{m}$ thickness.

Column temperature: Maintain at 40°C for 20 minutes after injection, raise the temperature to 240°C at a rate of 10°C per minute, and maintain at 240°C for 20 minutes.

Injection port temperature: 140°C.

Detector temperature: 250°C.

Carrier gas: Helium.

Flow rate: About 35 cm per second.

Split ratio: 1:3. (Note: The split ratio can be modified in order to optimize sensitivity.)

System suitability—

Test for required detectability: When the procedure is run with Class 1 standard solution and Class 1 system suitability solution under the above operating conditions, the SN ratio of the peak of 1,1,1-trichloroethane in Class 1 standard solution is not less than 5, and the SN ratio of each peak in Class 1 system suitability solution is not less than 3, respectively.

System performance: When the procedure is run with Class 2 standard solution A or the solution for system suitability under the above operating conditions, the resolution between acetonitrile and dichloromethane is not less than 1.0. Pipet 1 mL of the *N*,*N*-dimethylformamide solution of Residual Solvents RS for System Suitability (1 in 100) in an appropriate headspace vial, add exactly 5 mL of water, and apply the stopper, cap, and mix, and use this solution as the solution for system suitability.

System repeatability: When the test is repeated 6 times with Class 1 standard solution under the above operating conditions, the relative standard deviation of each peak areas is not more than 15%.

Separately inject (use headspace operating parameters in column 3 of Table 2.46-5) equal volumes of headspace (about 1.0 mL) of Class 1 standard solution, Class 2 standard solution A, Class 2 standard solution B, Class 2 standard solution C, and test solution into the chromatograph, record the chromatograms, and measure the responses for the major peaks. If a peak response of any peak, other than a peak for 1,1,1-trichloroethane, in the test solution is greater than or equal to a corresponding peak in either Class 1 standard solution, Class 2 standard solution A, Class 2 standard solution B or Class 2 standard solution C, or a peak response of 1,1,1-trichloroethane is greater than or equal to 150 times the peak response corresponding to 1,1,1-trichloroethane in Class 1 standard solution, proceed to Procedure B to verify the identity of the peak; otherwise, the article meets the requirements of this test.

1.2.2. Procedure B

The test is performed by gas chromatography $\langle 2.02 \rangle$ according to the following conditions.

Class 1 standard stock solution, Class 1 standard solution, Class 1 system suitability solution, Class 2 standard stock solutions A, B and C, Class 2 standard solutions A, B and C, test stock solution, and test solution: Proceed as directed for Procedure A.

Proceed as directed for Procedure B under Water-soluble articles with a split ratio of 1:3. (Note: The split ratio can be modified in order to optimize sensitivity.) The solution for system suitability: Proceed as directed for Procedure A.

Separately inject (use headspace operating parameters in Table 2.46-5) equal volumes of headspace (about 1.0 mL) of Class 1 standard solution, Class 2 standard solution A, Class 2 standard solution B, Class 2 standard solution C and test solution into the chromatograph, record the chromatograms, and measure the responses for the major peaks. If the peak response(s) in test solution of the peak(s) identi-

fied in Procedure A is/are greater than or equal to a corresponding peak(s) in either Class 1 standard solution, Class 2 standard solution A, Class 2 standard solution B or Class 2 standard solution C, proceed to Procedure C to quantify the peak; otherwise, the article meets the requirements of this test.

1.2.3 Procedure C

The test is performed by gas chromatography $\langle 2.02 \rangle$ according to the following conditions.

Class 1 standard stock solution, Class 1 standard solution, Class 1 system suitability solution, Class 2 standard stock solution A, Class 2 standard solution A, Class 2 standard stock solution C and Class 2 standard solution C: Proceed as directed for Procedure A.

Standard stock solution (Note: Prepare a separate standard stock solution for each peak identified and verified by Procedures A and B. For Class 1 solvents other than 1,1,1-trichloroethane, prepare the first dilution as directed for the first dilution under Class 1 standard stock solution in Procedure A.): Transfer an accurately measured volume of each individual solvent corresponding to each residual solvent peak identified and verified by Procedures A and B to a suitable container, and dilute quantitatively and stepwise if necessary, with water to obtain a solution having a final concentration of 1/20 of the value stated in Table 2.46-1 or Table 2.46-2.

Standard solution: Pipet 1 mL of standard stock solution in an appropriate headspace vial, add exactly 5 mL of water, and apply the stopper, cap, and mix.

Test stock solution: Weigh accurately about $0.5~{\rm g}$ of the article under test, and add N,N-dimethylformamide to make exactly $10~{\rm mL}$.

Test solution: Pipet 1 mL of test stock solution in an appropriate headspace vial, add exactly 5 mL of water, and apply the stopper, cap, and mix.

Spiked test solution (Note: prepare a separate spiked test solution for each peak identified and verified by Procedure A and B): Pipet 1 mL of test stock solution in an appropriate headspace vial, add exactly 4 mL of water, and apply the stopper, cap, and mix.

Operating conditions and system suitability fundamentally follow the procedure A. Test for required detectability is unnecessary, and use Standard solution instead of Class 1 standard solution for system repeatability. If the results of the chromatography from Procedure A are found to be inferior to those found with Procedure B, the operating conditions from Procedure B may be substituted.

Separately inject (use headspace operating parameters in Table 2.46-5) equal volumes of headspace (about 1.0 mL) of the standard solution, test solution, and spiked test solution into the chromatograph, record the chromatograms, and measure the responses for the major peaks. Calculate the amount of each residual solvent found in the article under test by the formula:

Residual solvent (ppm) = 10 $(C/M)\{A_T/(A_S - A_T)\}$

C: Concentration (µg/mL) of the appropriate Reference Standard in the standard stock solution

Table 2.46-5 Headspace operating parameters

	Headspace Operating Parameter Sets		
	1	2	3
Equilibration temperature (°C)	80	105	80
Equilibration time (min.)	60	45	45
Transfer-line temperature (°C)	85	110	105
Syringe temperature (°C)	80 – 90	105 – 115	80 – 90

Carrier gas: nitrogen or helium at an appropriate pressure

Pressurization time (s)	≥ 60	≥ 60	≥ 60
Injection volume (mL)*	1	1	1

- * Or follow the instrument manufacture's recommendations, as long as the method criteria are met. Injecting less than this amount is allowed as long as adequate sensitivity is achieved.
 - M: Amount (g) of the article under test taken to prepare the test stock solution
 - A_T : Peak responses of each residual solvent obtained from the test solution
 - A_s: Peak responses of each residual solvent obtained from the spiked test solution

1.3. Headspace operating parameters and other considera-

Examples of headspace operating parameters are shown in Table 2.46-5.

These test methods describe the analytical methods using the headspace gas chromatography. The following Class 2 residual solvents are not readily detected by the headspace injection conditions: 2-ethoxyethanol, ethylene glycol, formamide, 2-methoxyethanol, N-methylpyrrolidone, and sulfolane. Other appropriate validated procedures are to be employed for the quantification of these residual solvents. In the headspace methods, N,N-dimethylformamide and N,N-dimethylacetoamide are often used as solvents. As not only 6 solvents described above but these two solvents are not included in the Residual Solvents Class 2A RS, the Residual Solvents Class 2B RS and/or the Residual Solvents Class 2C RS, appropriate validated procedures are to be employed for these residual solvents as necessary.

2. Class 3 residual solvents

Perform the test according to 1. Otherwise an appropriate validated procedure is to be employed. Prepare appropriately standard solutions, etc. according to the residual solvent under test.

If only Class 3 solvents are present, the level of residual solvents may be determined by Loss on Drying <2.41>. However when the value of the loss on drying is more than 0.5%, or other solvents exist, the individual Class 3 residual solvent or solvents present in the article under test should be identified using the procedures as described above or other appropriate procedure, and quantified as necessary.

3. Reference Standards

- (i) Residual Solvents Class 1 RS (A mixture of benzene, carbon tetrachloride, 1,2-dichloroethane, 1,1-dichloroethene and 1,1,1-trichloroethane)
- (ii) Residual Solvents Class 2A RS (A mixture of acetonitrile, chlorobenzene, cumene, cyclohexane, 1,2-dichloroethene (cis-1,2-dichloroethene, trans-1,2-dichloroethene), dichloromethane, 1,4-dioxane, methanol, methylcyclohexane, tetrahydrofuran, toluene and xylene (m-xylene, p-xylene, o-xylene, ethylbenzene))
- (iii) Residual Solvents Class 2B RS (A mixture of chloroform, 1,2-dimethoxyethane, hexane, methyl butyl ketone, nitromethane, pyridine, tetralin and 1,1,2-trichloroethene)
- (iv) Residual Solvents Class 2C RS (Methyl isobutyl ketone)
- (v) Residual Solvents RS for System Suitability (A mixture of acetonitrile, *cis*-1,2-dichloroethene and dichloromethane)

2.51 Conductivity Measurement

Change the following as follows:

This test is harmonized with the European Pharmacopoeia and the U. S. Pharmacopeia.

The parts of the text that are not harmonized among the targeted texts for the harmonization are marked with symbols ($^{\bullet}$ $_{\bullet}$), and the texts that are uniquely specified by the JP other than the targeted texts for the harmonization are marked with symbols ($^{\diamond}$ $_{\diamond}$).

This chapter provides information on how to apply electrical conductivity measurements (hereafter referred to as "conductivity") of fluid solutions, including pure fluids. This chapter is intended for fluid applications when conductivity is used to measure, monitor, or control chemical dispensing, chemical purity, ionic concentration, and other applications where the ionic character of the fluid needs to be known or controlled.

Applications include, but are not limited to, solutions that may be used in clean-in-place, chromatography detection, ionic solution preparations, end point detection, dosing, fermentation, and buffer production. In some cases, conductivity measurements can be extended to pure organic fluids such as alcohols and glycols where a weak conductivity signal exists, and the signal can be significantly increased if the organics become contaminated with water or salts.

Conductivity is the measurement of the ability of a fluid to conduct electricity via its chemical ions. The ability of any ion to electrically conduct is directly related to its ion mobility. Conductivity is directly proportional to the concentrations of ions in the fluid, according to Equation (1):

$$\kappa = 100 \sum_{i}^{all} C_i \cdot \lambda_i$$
 (1)

 $\kappa = \text{conductivity (S/cm)}$

 C_i = concentration of chemical ion i (mol/L)

 λ_i = specific molar conductance of ion i (S·cm²/mol)

Although the SI unit S/m is the appropriate SI unit for conductivity, historically the unit S/cm has been selected by industry as the accepted unit.

On the basis of Equation (1), conductivity is not ion selective because it responds to all ions. Furthermore, the specific molar conductance of each ion is different. As a result, unless the percentage composition of ions of the solution is limited and known, the precise concentrations of ionic species cannot be determined from conductivity measurements. However, for examples such as a solution of a single salt or acid or base, such as a caustic solution used in cleaning, the precise concentration can be directly determined. Despite the lack of ionic specificity, conductivity is a valuable laboratory and process tool for measurement and control of total ionic content because it is proportional to the sum of the concentrations of all ionic species (anions and cations) for diluted solutions as described in Equation (1). At higher concentrations, conductivity measurements are not perfectly linear with concentration. Conductivity measurements cannot be applied to solids or gases, but they can be applied to the condensate of gases.

Another variable that influences conductivity measurements is the fluid temperature. As the fluid temperature increases, the ion conductance increases, making this physicochemical phenomenon the predominant reason for the temperature-compensation requirement when testing conductive fluids.

The conductivity, κ , is proportional to the conductance, G (S), of a fluid between two electrodes (Equation (2)):

$$\kappa = G \times (d/A) = G \times K \tag{2}$$

 $\kappa = \text{conductivity (S/cm)}$

G = conductance (S)

d = distance between the electrodes (cm)

A =area of the conducting electrodes (cm²)

 $K = \text{cell constant (cm}^{-1})$, which also equals the ratio of d/A

The resistivity p ($\Omega \cdot \text{cm}$) of the fluid is, by definition, the reciprocal of the conductivity (Equation (3)):

$$\rho = 1/\kappa = 1/(G \times K) = R/K \tag{3}$$

 $\rho = \text{resistivity } (\Omega \cdot \text{cm})$

 $\kappa = \text{conductivity (S/cm)}$

G = conductance (S)

 $K = \text{cell constant (cm}^{-1})$

 $R = \text{resistance } (\Omega)$, which is the reciprocal of the conductance, G

1. Apparatus

An electrical conductivity measurement consists of the determination of resistance of the fluid between and around the electrodes of the conductivity sensor. To achieve this measurement, the primary instrumentation is the resistance-measuring circuit and the conductivity sensor, and they are

usually connected by a cable when the sensor and the user interface are separated.

The resistance measurement is made by applying an alternating current (AC, meaning the flow of electric charge periodically reverses direction) voltage (or current) to the electrodes, measuring the current (or voltage), and calculating the resistance according to Ohm's Law. The alternating source is used to prevent the polarization (collection of ions) at the electrodes. Depending on the instrument, the measuring frequency of the measuring system adjusts automatically according to the measuring conditions of the instrument, and there may be multiple resistance-measuring circuits embedded in the measuring system. The resistance-measurement circuit may be embedded in the transmitter or in the sensor.

The conductivity sensor consists of at least two electrical conductors of a fixed size and geometry, separated by an electrical insulator. The electrodes, insulator, and any other wetted materials should be constructed of materials that are unreactive to fluids with which they may come into contact. Also, the sensor construction should withstand the environmental conditions (process or ambient temperature, pressure, cleaning applications) that it would be subjected to.

Most conductivity sensors have temperature devices such as a platinum resistance temperature device (RTD) or negative temperature coefficient (NTC) thermistor embedded inside the sensor, although external temperature measurement is possible. The purpose of the temperature measurement is for temperature compensation of the conductivity measurement.

2. Cell Constant Determination

The purpose of the sensor's cell constant is to normalize the conductance (or resistance) measurement for the geometrical construction of the two electrodes.

The cell constant is determined by immersing the conductivity sensor in a solution of known conductivity. Solutions of known conductivity can be obtained by preparation of specific mixtures according to national authoritative sources or procurement of commercially available certified and traceable standard solutions. These recipes or certified solutions can range from 5 to $200,000\,\mu\text{S/cm}$, depending on the level of accuracy desired. Alternatively the cell constant is determined by comparison to other reference conductivity measuring systems (also available as an accredited calibration service). [Note—Conductivity measurements are not perfectly linear with concentration.]

The measured cell constant of the conductivity sensor must be within 5% of the nominal value indicated by the sensor certificate, unless otherwise prescribed.

Modern conductivity sensors normally do not change their cell constant over their lifetime. If a change of the cell constant is detected during calibration, a cleaning of the sensor is appropriate according to the manufacturer's recommendations. Following that, the calibration procedure should be repeated. Sometimes "memory effects" appear, particularly when changing from high to low concentrations if the sensor is not well flushed.

3. Calibration of Temperature

In addition to verifying the sensor's cell constant, the embedded temperature device (or external temperature device) should be appropriately calibrated for the application to apply the temperature compensation algorithm accurately. The temperature accuracy that is required depends on the criticality of the temperature to the application. An accuracy of \pm 1°C typically suffices.

4. Calibration of Measurement Electronics

The measurement circuit of the system is fundamentally an AC resistance measuring device. Appropriate verification and/or calibration of the measuring circuit is required for measurement systems with signal transfer via analog cable. This is accomplished by disconnecting the measuring circuit from the sensor's electrodes, attaching traceable resistors of known value with the cable of the measurement system to the measuring circuit, and verifying that the measured resistance agrees with the resistor value to an acceptable level. A typical acceptance criterion for the resistance accuracy is <2% of the reading at resistances>100 Ω , and increasing to 5% at lower resistances. However, the application criticality should ultimately determine the desired accuracy.

For conductivity systems that cannot have the resistance-measuring circuit disconnected from the electrodes (e.g., measurement circuit and electrodes in one mutual housing), it may be difficult to directly adjust or verify the circuit accuracy, depending on the sensor design. An alternative method of verifying the measurement system integrity is a system calibration according to the procedures for the cell constant determination for each measuring circuit that is intended to be used.

If verification/calibration of the sensor's cell constant, temperature device, and measuring circuit are done at the same service interval, the measuring circuit should be verified first, the temperature device next, and the cell constant last. Because all of these parameters are typically very stable due to modern electronics and stable sensor construction, frequent calibration (such as daily) is not usually required. Comparison to qualified reference systems is also a proper means of calibration. Calibration is performed at appropriate intervals as defined in the quality management system.

5. Temperature Compensation

Because the conductivity of a fluid is temperature dependent, temperature compensation of the conductivity measurement is necessary unless otherwise prescribed. An appropriate temperature compensation algorithm will ensure that changes in the conductivity measurement can be ascribed to concentration changes and not temperature changes. Conductivity measurements are normally referenced to 25°C. A common form of linear temperature compensation uses Equation (4):

$$\kappa_{25} = \frac{\kappa_T}{[1 + \alpha(T - 25)]} \tag{4}$$

 κ_{25} = conductivity compensated to 25°C

 κ_T = conductivity at T

 α = temperature coefficient of the conductivity

T = measured temperature

A temperature coefficient of 2.1% per 1°C is commonly used for many salt solutions. Most salt-based solutions have linear compensation factors ranging from 1.9% to 2.2% per 1°C. Depending on the fluid samples, other forms of temperature compensation may be appropriate. Non-linear temperature compensation will carry out temperature compensation using preprogrammed data in the instrument. Non-linear temperature compensation data for a variety of solutions is widely available, e.g. for natural waters, and for ultrapure water with traces of ammonia. \(\) In cases of very low conductivity ($< 10 \,\mu\text{S/cm}$), such as purified pharmaceutical waters, two compensations need to be made. One is for the intrinsic conductivity of water, and the other is for the other ionic species in water. These compensations are normally combined and embedded in the microprocessorcontrolled conductivity measurement systems. This is not supplied in all conductivity measurement technologies.

6. Conductivity Measurement of Fluids

For off-line batch measurements, rinse the cleaned sensor with the fluid to be measured. Then immerse the sensor in the fluid to be measured, and record the temperature and the temperature-compensated conductivity as required. Be sure that the position of the sensor in the container does not affect the conductivity measurement, because the container walls can affect the measurement for some electrode designs.

For continuous on-line or on-line measurements, install the cleaned sensor into the pipe, tank, or other containment vessel, and flush, if necessary. Make sure proper installation procedures are applied to prevent bubbles or particles from collecting between the electrodes. Be sure that the position of the sensor in the pipe or tank does not affect the conductivity measurement, because the nearby surfaces can affect the measurement for some electrode designs.

Record the temperature and the temperaturecompensated conductivity as required.

For all batch or continuous measurements, ensure that the wetted components of the sensor are compatible with the fluid and the temperature to be measured.

Add the following:

2.66 Elemental Impurities— Procedures

Procedures of Elemental Impurities are methods to control elemental impurities contained in drug products and their components, etc. This chapter describes two analytical procedures (Procedures 1 and 2) and validation criteria for the evaluation of the levels of elemental impurities. The chapter permits the use of any procedure that meets the validation criteria specified in this chapter. As the chemical composition of the considered substances and the specification limits for the element(s) of interest vary considerably, it

is difficult to describe all suitable sample preparation and measurement methods. By means of validation studies, analysts will confirm that the analytical procedure is suitable for use on specified material. It is not necessary to cross validate against either procedure 1 or 2 provided that requirements for procedure validation are met. As elemental impurities may be ubiquitous they have the potential to be present in trace amounts therefore special precautions may be necessary to avoid sample contamination. (Note: Methods such as atomic absorption spectrometry other than methods described in this chapter, if validated, can also be used without cross validation against analytical procedure 1 or 2.)

1. Sample Preparation

Forms of sample preparation include Neat, Direct aqueous solution, Direct organic solution, and Indirect solution. The selection of the appropriate sample preparation depends on the material under test and is the responsibility of the analyst. When a sample preparation is not indicated in the monograph, an analyst may use any appropriately validated sample preparation procedure, including but not limited to procedures described below. In cases where spiking of a material under test is necessary to provide an acceptable signal intensity, the blank should be spiked with the same Target elements, and where possible, using the same spiking solution. The material or mixture under test must be spiked before any sample preparation steps are performed. Standard solutions may contain multiple Target elements. (Note: If intended for a quantitative test, appropriate material handling procedures should be followed e.g. volatile liquids should be pipetted, viscous liquids should be weighed.)

Neat: Used for liquids or samples measurable without addition of solvent.

Direct aqueous solution: Used when the sample is soluble in an aqueous solvent.

Direct organic solution: Used when the sample is soluble in an organic solvent.

Indirect solution: Generally, an indirect solution is obtained when a material is not directly soluble in aqueous or organic solvents. Total metal extraction is the preferred sample preparation approach to obtain an indirect solution. Digest the sample using the Closed vessel digestion procedure provided below or one similar to it.

Closed vessel digestion: This sample preparation procedure is designed for samples that must be digested in a Concentrated acid using a closed vessel digestion apparatus. Closed vessel digestion minimizes the loss of volatile impurities. The choice of a Concentrated acid depends on the sample matrix. The use of any of the Concentrated acids may be appropriate, but each introduces inherent safety risks. Therefore, appropriate safety precautions should be used at all times. (Note: Weights and volumes provided may be adjusted to meet the requirements of the digestion apparatus used.)

An example procedure that has been shown to have broad applicability is the following. Dehydrate and predigest 0.5 g

of material under test in 5 mL of freshly prepared Concentrated acid. Allow to sit loosely covered for 30 minutes in a fume hood. Add an additional 10 mL of Concentrated acid, and digest, using a closed vessel technique, until digestion or extraction is complete. Repeat, if necessary, by adding an additional 5 mL of Concentrated acid. (Note: Where closed vessel digestion is necessary, follow the manufacturer's recommended procedures to ensure safe use.)

Clear solutions are expected in the validation. In those cases where a clear solution cannot be obtained, appropriate studies should ensure that the recovery is suitable for the intended use.

Reagents: All reagents used for the preparation of sample and standard solutions should be sufficiently pure for the intended purpose.

2. Analytical Procedures 1 and 2

System standardization and suitability evaluation using applicable reference materials should be performed for each analytical sequence.

2.1. Procedure and Detection Technique

Procedure 1 can be used for elemental impurities generally amenable to detection by inductively coupled plasma-atomic (optical) emission spectroscopy (ICP – AES or ICP – OES). Procedure 2 can be used for elemental impurities generally amenable to detection by inductively coupled plasma-mass spectrometry (ICP – MS). Before initial use, the analyst should verify that the procedure is appropriate for the instrument and sample used by meeting the procedure validation requirements below.

2.2. Procedure 1: ICP - OES

Standard solution 1: 1.5J of the Target element(s) in a Matrix matched solution

Standard solution 2: 0.5J of the Target element(s) in a Matrix matched solution

Sample stock solution: Proceed as directed in Sample Preparation above. Allow the sample to cool, if necessary. For mercury determination, add an appropriate stabilizer, if necessary.

Sample solution: Dilute the Sample stock solution with an appropriate solvent to obtain a final concentration of the Target element(s) within the calibrated range.

Blank: Matrix matched solution **Elemental spectrometric system**

Mode: ICP

Detector: Optical detection system **Rinse:** Generally, diluent used

Calibration: Standard solution 1, Standard solution 2, and

Blank

System suitability Sample: Standard solution of the Target element(s) in a Matrix matched solution at a concentration within the calibrated range

Suitability requirements

Short term Instrumental Stability: Compare results obtained from System suitability sample before and after the analysis of the Sample solution.

Suitability criteria: NMT 20% deviation between both samples for each Target element. (Note: If samples are high

in mineral content, rinse the system well in order to minimize carryover and check it by measuring a blank solution before introducing the System Suitability Sample.)

Analysis: Analyze according to manufacture's suggestions for wavelength. Calculate and report results on the basis of the original sample size. [Note: Appropriate measures must be taken to correct for matrix-induced interferences (e.g., wavelength overlaps).]

2.3. Procedure 2: ICP - MS

Standard solution 1: 1.5J of the Target element(s) in a Matrix matched solution

Standard solution 2: 0.5J of the Target element(s) in a Matrix matched solution

Sample stock solution: Proceed as directed in Sample Preparation above. Allow the sample to cool, if necessary. For mercury determination, add an appropriate stabilizer, if necessary.

Sample solution: Dilute the Sample stock solution with an appropriate solvent to obtain a final concentration of the Target element(s) within the calibrated range.

Blank: Matrix matched solution

Elemental spectrometric system

Mode: ICP. [Note: An instrument with a cooled spray chamber is recommended. (A collision cell or reaction cell may also be beneficial.)]

Detector: Mass spectrometer **Rinse:** Generally, diluent used

Calibration: Standard solution 1, Standard solution 2, and

Blanl

System suitability Sample: Standard solution of the Target element(s) in a Matrix matched solution at a concentration within the calibrated range

Suitability requirements

Short term Instrumental Stability: Compare results obtained from system suitability sample before and after the analysis of the Sample solution.

Suitability criteria: NMT 20% deviation between both samples for each Target element. (Note: If samples are high in mineral content, rinse the system well in order to minimize carryover and check it by measuring a blank before introducing the System suitability sample.)

Analysis: Analyze according to the manufacturer's suggestions for program and m/z. Calculate and report results based on the original sample size. [Note: Appropriate measures must be taken to correct for matrix-induced interferences (e.g., argon chloride interference with arsenic determinations).]

3. Requirements for Procedure Validation

All procedures must be validated and shown to be acceptable, in accordance with the validation requirements described below. The level of validation necessary to ensure that a procedure is acceptable depends on whether a limit test or a quantitative determination is used. Any procedure that has been validated and meets the acceptance criteria that follow is considered to be suitable for use. If appropriate, the validation method and criteria may be changed according to the purpose of evaluating the levels of the con-

tent of elemental impurities. They may differ from the requirements to meet the system suitability criteria described in Inductively Coupled Plasma-Atomic Emission Spectrometry and Inductively Coupled Plasma-Mass Spectrometry <2.63>.

3.1. Procedures for Limits Tests

The following section defines the validation parameters for the acceptability of limit tests. Meeting these requirements must be demonstrated experimentally using an appropriate system suitability test and reference materials.

The suitability of the method must be determined by conducting studies with the material or mixture under test spiked with known concentrations of each Target element of interest at the appropriate Target concentration.

3.1.1. Detectability

Standard solution: A preparation of reference materials for the Target element(s) at 1.0*J* in a Matrix matched solution. **Spiked sample solution 1:** Prepare a solution of the sample under test, spiked with appropriate reference materials for the Target element(s) at the Target concentration, solubilized or digested as described in Sample Preparation.

Spiked sample solution 2: Prepare a solution of the sample under test, spiked with appropriate reference materials for the Target element(s) at 80% of the Target concentration, solubilized or digested as described in Sample Preparation. Unspiked sample solution: A sample of material under test, solubilized or digested in the same manner as the spiked Sample solutions

Acceptance criteria

Non-instrumental procedures: Spiked sample solution 1 provides a signal or intensity equivalent to or greater than that of the Standard solution. Spiked sample solution 2 must provide a signal or intensity less than that of Spiked sample solution 1. (Note: The signal from each Spiked sample solution is NLT the Unspiked sample solution determination.)

Instrumental procedures: The average value of the three replicate measurements of Spiked sample solution 1 is within \pm 15% of the average value obtained for the replicate measurements of the Standard solution. The average value of the replicate measurements of Spiked sample solution 2 must provide a signal intensity or value less than that of the Standard solution. (Note: Correct the values obtained for each of the spiked solutions using the Unspiked sample solution.)

3.1.2. Specificity

The procedure must be able to unequivocally assess each Target element in the presence of components that may be expected to be present, including other Target elements, and matrix components.

3.1.3. Precision, only for Instrumental Methods (Repeatability)

Sample solutions: Six independent samples of the material under test, spiked with appropriate reference materials for the Target elements at the Target concentration

Acceptance criteria

Relative standard deviation: NMT 20% for each Target element

3.2. Procedures for Quantitative Tests

The following section defines the validation parameters for the acceptability of procedures for quantitative tests. Meeting these requirements must be demonstrated experimentally, using an appropriate system suitability test and reference materials.

3.2.1. Accuracy

Standard solutions: Prepare solutions containing the Target element(s) at three concentrations ranging from 0.5 to 1.5J, using appropriate reference materials, in a Matrix matched solution and blank.

Test samples: Prepare samples of the material under test spiked with appropriate reference materials for the Target element(s) before any sample preparation steps (digestion or solubilization) at 3 concentrations ranging from 50% to 150% of the Target concentration. Spike concentrations should range from 0.5 to 1.5J, and should contain at least three individual concentrations.

Acceptance criteria

Spike recovery: 70% – 150% for the mean of three replicate preparations at each concentration

3.2.2. Precision

Repeatability

Test samples: Six independent samples of material under test (taken from the same lot) spiked with appropriate reference materials for the Target element(s) at the Target concentration. Or at least 9 determinations (e.g., 3 replicates of 3 concentrations) covering the specified range.

Acceptance criteria

Relative standard deviation: NMT 20% (n = 6) for each Target element

Intermediate precision (ruggedness)

Perform the Repeatability analysis again at least once either on a different day, with a different instrumentation, with a different analyst, or a combination thereof. Combine the results of this analysis with the Repeatability analysis so the total number of samples is at least 12.

Acceptance criteria

Relative standard deviation: NMT 25% (n = 12) for each Target element

3.2.3. Specificity

The procedure must be able to unequivocally assess each Target element in the presence of components that may be expected to be present, including other Target elements, and matrix components.

3.2.4. Range and Linearity

Demonstrated by meeting the Accuracy requirement.

3.2.5. Limit of Quantification

The Limit of Quantification is confirmed when the accuracy acceptance criteria for the corresponding spiked solution is met. The Limit of Quantification is smaller or equal to 50% of Target concentration.

4. Glossary

- (i) Concentrated acid: Concentrated ultra-pure nitric, sulfuric, hydrochloric, or hydrofluoric acids or any other acid or mixture of acids that is demonstrated suitable.
 - (ii) Matrix matched solution: Solutions having the same

solvent composition as the Sample solution. In the case of an aqueous solution, Matrix matched solution would indicate that the same acids, acid concentrations and mercury stabilizer are used in both preparations.

- (iii) Target elements: Elements whose levels in the drug product must be controlled within acceptable limits.
- (iv) Target limit or Target concentration: The acceptance value for the elemental impurity being evaluated. Exceeding the Target limit indicates that a material under test exceeds the acceptable value. Target limits in the final drug product can be approximated by dividing the permitted daily exposures (PDEs) by the maximum daily dose of the drug product. When evaluating the significance of elemental impurity levels, it is possible to set the Target limits to the values obtained by dividing 30% of PDEs by the maximum daily dose of the drug product. Furthermore, when the permitted concentration limit of each element in the individual components of the drug product is set, it can be set as the Target concentration.
- (v) J: The concentration (w/v) of the Target element(s) at the Target limit, appropriately diluted to the working range of the instrument. If a dilution is not necessary, J is equal to the Target concentration. For example, if the target elements are lead and arsenic for an analysis of an oral solid drug product with a daily dose of 10 g/day using inductively coupled plasma mass spectrometry (ICP MS), the target limit for these elements would be $0.5 \,\mu\text{g/g}$ and $1.5 \,\mu\text{g/g}$. However, in both cases, the linear dynamic range of the ICP MS is known to extend from $0.01 \,\text{ng/mL}$ to $0.1 \,\mu\text{g/mL}$ for these elements. Therefore, a dilution factor of at least 1:100 is required to ensure that the analysis occurs in the linear dynamic range of the instrument. J would thus equal 5 $\,\text{ng/mL}$ and 15 $\,\text{ng/mL}$ for lead and arsenic, respectively.
- (vi) Appropriate reference materials: In principle, where Appropriate reference materials are specified in the chapter, certified reference materials (CRM) from a national metrology institute (NMI), or reference materials that are traceable to the CRM of an NMI should be used.
- (vii) Cross validation: Verification whether or not the same result can be obtained from the corresponding analyses for the same sample.

Add the following:

6.16 Rheological Measurements for Semi-solid Preparations

Rheological measurements for semi-solid preparations are methods to measure fluidity and deformation by adding force to semi-solid preparations such as Semi-solid Preparations for Oro-mucosal Application, Ophthalmic Ointments, Ointments, Creams and Gels.

There are two methods, the first (spreadability test) and the second (penetrometry).

These methods are more practical to determine rheologi-

cal properties of semi-solid preparations, though Method II Viscosity measurement by rotational viscometer under Viscosity Determination <2.53> can evaluate precisely rheological properties of semi-solids.

1. Method 1 Spreadability test

Spreadability test is a method for measuring flowability of semi-solid preparations using a spread meter (also known as a parallel plate viscometer).

A spread meter measures a spreading diameter etc. of a sample by observing the characteristic of the concentric spreading over time when the sample is sandwiched between two parallel plates [glass (or plastic) plate and fixed plate] placed horizontally and pushed outward with the own weight of the glass plate.

Generally there is a reciprocal relationship between fluidity as an index of flowability and viscosity, however flowability measured by this method does not necessarily correlate with apparent viscosity (mPa·s) measured by the viscosity determination.

This method targets relatively soft preparations in particular among semi-solid preparations.

1.1. Apparatus

An example of spread meters is shown in Fig. 6.16-1.

There are two parallel plates placed horizontally. The fixed plate installed in the lower part has engraved scales at 1 mm intervals to measure the distance from the center and has a cylindrical hole (0.5 mL) at the center for inserting a sample. The weight-loading glass plate positioned in the upper part is made of transparent glass or acrylate resin etc. and has a mass of 115 g. The glass plate is supported by rods. Pushing the piston up raise the bottom of the hole to push the sample above the fixed plate. In connection with it the glass plate falls (20 mm, commonly) and the sample spreads concentrically on the fixed plate with the own weight of the glass plate. Measuring the extent of the spread provides the flowability of the sample.

1.2. Procedure and measurement conditions

Before measurement remove a glass plate from an apparatus and clean the glass plate, a fixed plate and a specimen hole. Fill a sample in the specimen hole at the center of the fixed plate and flatten the sample by a spatula etc. so that the upper face of the sample is flush with the fixed plate and become flat. Wipe the sample protruded off. Make sure that the apparatus is horizontal by a level and install the glass plate to the support rod. Push the piston up and start time measurement simultaneously. Measure the spread of the sample on the fixed plate by the scale in mm unit over time and record them.

A constant temperature is needed during measurement and it is preferable to perform the measurement in the room controlled at $25 \pm 2^{\circ}$ C.

Bring the temperature of samples equal to that of the measurement environment by allowing them to stand in the measurement environment.

1.3. Analysis

Characteristics of flowability obtained with the measurement using a spread meter are expressed as spread meter

diameter D and spread meter yield value YV in the case of the single-point method and as spread meter slope S and spread meter intercept IC in the case of the multi-point method. They are calculated by the following methods, respectively.

1.3.1. Single-point method

(i) Spread meter diameter *D*: Expressed as the spreading diameter (mm) after the specified time (usually, 60 seconds).

A larger spreading diameter indicates higher flowability.

(ii) Spread meter yield value YV: Calculate by the following equation.

$$YV = (4.8 \times W \times V \times g_n)/(\pi^2 \times D_{\infty}^5)$$

YV: Yield value (Pa)

W: Mass (kg) of glass plate

V: Sample volume (m³)

 g_n : Standard acceleration due to gravity (m/s²)

 D_{∞} : Maximum of spreading diameter (m)

The above equation is expressed in the International System of Units, however actual measurement is performed in the Centimeter-Gram-Second System of Units.

Many semi-solid preparations have no flowability when leave them as they are, however they flow by adding force. The limit value of the force is the yield value and a larger yield value indicates that a larger force is required for the flow of a sample.

1.3.2. Multi-point method

(i) Spread meter slope S: Calculate by the following equation.

$$S = (D_2 - D_1)/\log_{10}(T_2/T_1)$$

 D_1 : Spreading diameter after T_1 seconds (mm)

 D_2 : Spreading diameter after T_2 seconds (mm)

$$T_1$$
, T_2 : Measurement time (seconds) $T_2 > T_1$, $5 \le T_1$ and $T_2 \le 100$, $\Delta T = (T_2 - T_1) > 40$

In general, plot the spreading diameter measured each time on a semi-logarithmic graph to obtain an almost straight line connecting each point. Spread meter slope *S* corresponds to its slope.

A larger S indicates a larger flow of a sample.

(ii) Spread meter intercept IC: Plot logarithms of times (T_1, T_2) as abscissa against spreading diameters (D_1, D_2) as ordinate and draw a line connecting these two points. Obtain the intersection point (T = 1) of this extension line and the ordinate axis as spread meter IC.


A larger IC indicates higher flowability of a sample.

2. Method 2 Penetrometry

Penetrometry is a method for measuring hardness or softness of semi-solid preparations using a penetrometer.

A penetrometer is an apparatus for measuring the distance traveled by a cone penetrant inside a sample and the consistency is expressed as ten times of the value measured in units of 0.1 mm. A smaller value indicates a harder sample.

This method targets relatively hard preparations in particular among semi-solid preparations.






Fig. 6.16-1 Example of spread meter

2.1. Apparatus

2.1.1. Penetrometer

An example of penetrometers is shown in Fig. 6.16-2.

Adjust a penetrometer so that the position of the tip of a prescribed cone touches the surface of the sample filled in a sample container, penetrate the cone into the sample with its own weight for a certain time and calculate the consistency from the depth measured in units of 0.1 mm.

Adjust exactly the cone part or the sample stage of the penetrometer so that the position of the tip of the cone touches the horizontal surface of the sample keeping the reading of the dial gage zero. Adjust the cone in advance so that the cone falls more than 62 mm smoothly when released from the state fixed to the penetrometer and the top of the cone does not hit the bottom of the sample container after the fall. The penetrometer must be equipped with a screw for horizontal control and a level to maintain a cone holder being vertical.

2.1.2. Cone

A standard cone is a conical body made of magnesium or other suitable metal, which has a removable hardened steel needle. Fig. 6.16-3 shows the specifications of the standard cone.

The upper end of a holder is equipped with a stopper and the lower end is devised to connect a cone. The inner structure of the cone may be changed in order to conform to the prescribed mass as far as there is no change in shape and mass distribution of the cone. The outer surface of the cone must be without scratches and smooth enough.

Fig. 6.16-4, 6.16-5 and 6.16-6 show the specifications of an optional cone, a half-scale cone and a quarter-scale cone. The half-scale and quarter-scale cones are the prescribed cones that are scaled down to a half and a quarter of the






Fig. 6.16-2 Example of penetrometer



Total mass of cone: $102.5\pm0.05~g$ Total mass of holder: $47.5\pm0.05~g$

Fig. 6.16-3 Standard cone



Total mass of cone: $102.5\pm0.05~g$ Total mass of holder: $47.5\pm0.05~g$

Fig. 6.16-4 Optional cone




Fig. 6.16-5 Half-scale cone




Fig. 6.16-6 Quarter-scale cone

standard cone or the optional cone.

Select a suitable cone based on the sample volume and the consistency of the preparation to be tested and convert the measured consistency to the value equivalent to that measured using a standard cone.

An optional cone may be used to measure the sample having a consistency of not more than 400. Use a half-scale cone or a quarter-scale cone, when a sample has a consistency of 175 to 385 and a standard cone is not available because of small sample amounts.

2.2. Procedure and measurement conditions

2.2.1. Sample preparation

Stand a prepared empty container and a container with a lid containing a necessary amount of sample in a water bath at 25°C and bring the temperature of the sample at 25 ± 0.5°C. Transfer the sample in the container with a lid to the empty container at a time, if possible. Remove air bubbles mixed in the sample in the sample container by a suitable method and fill the container with an excessive amount of the sample using a spatula again. Be careful not to stir the sample and not to make vacant space inside the sample in this procedure. Flatten the surface of the sample by removal of the excessive sample by moving a spatula, of which surface is inclined about 45 degrees in the direction of the movement, along the upper edge of the container. Hereafter do not touch the surface of the sample before measurement. Perform the test promptly to keep the temperature of the sample evenly at 25 ± 0.5 °C.

The consistency of soft samples is affected by the container diameter. With soft samples having a consistency of not less than 265, use a container with inside diameter of 76.2 mm (a half-scale cone: 38.1 mm, a quarter-scale cone: 19.0 mm). The consistency of relatively hard samples having a consistency of not more than 265 is hardly affected by the container diameter when using a container with inside diameter of not less than 76.2 mm.

The amount of sample required for the test depends on the sizes of the sample container and the cone, and the procedure specified according to the cosistency.

2.2.2. Procedure

2.2.2.1. Standard cone

Stand a sample container gently on the sample stage of the penetrometer adjusted to horizontal position. After placing the position of the cone to the zero point of the dial gage adjust the cone so that the tip touches the surface of the sample at the position prescribed in (i), (ii) or (iii) by moving either the cone part or the sample stage up and down. Push the metal clamping tool immediately to penetrate the cone for 5 ± 0.1 seconds. The holder must move smoothly in the drop unit. Push the rack for measurement down gently until it stops and read the dial gage in first place after the decimal point.

- (i) With a sample having a consistency of more than 400, only one test may be made in one container by placing the center of the sample container within 0.3 mm from the needle tip. Three tests may be made in three containers.
- (ii) With a sample having a consistency of more than 200 and less than 400, perform carefully the centering of the

cone in the sample container. Only one test may be made with this sample. Three tests may be made in three containers.

(iii) With a sample having a consistency of not more than 200, three tests may be made in one container. The measurement points are the midpoints between the center and the edge of the sample container at about 120 degrees intervals on a concentric circle.

2.2.2.2. Half-scale cone and quarter-scale cone

Place a cone at the center of a sample surface and perform a penetration pretest. The pretest may be omitted if an approximate value of consistency is known. Perform the measurement of consistency according to the procedure in 2.2.2.1. and place the position of the tip according to the following (i) or (ii).

- (i) With a sample having a consistency of more than 97 when using a half-scale cone or more than 47 when using a quarter-scale cone, only one test may be made in one container by placing carefully the needle tip at the center of the sample container. Three tests may be made in three containers.
- (ii) With a sample having a consistency of not more than 97 when using a half-scale cone or not more than 47 when using a quarter-scale cone, three tests may be made in one sample container. The measurement points are the midpoints between the center and the edge of the sample container at about 120 degrees intervals on a concentric circle so that the cone does not collide with the edge of the sample container and the measurement position of previous tests.

2.3. Analysis

2.3.1. Conversion of consistency measured using a half-scale or quarter-scale cone

Convert the consistency obtained using a half-scale or quarter-scale cone to the value equivalent to that obtained using a standard cone by the following equation.

2.3.1.1. Conversion of consistency obtained using a half-scale cone

$$P = 2p_{1/2} + 5$$

P: Consistency converted to the value equivalent to that measured using the standard cone

 $p_{1/2}$: Consistency measured using a half-scale cone

2.3.1.2. Conversion of consistency obtained using a quarter -scale cone

$$P = 3.75p_{1/4} + 24$$

P: Consistency converted to the value equivalent to that measured using a standard cone

 $p_{1/4}$: Consistency measured using a quarter-scale cone

Add the following:

6.17 Insoluble Particulate Matter Test for Therapeutic Protein Injections

Insoluble particulate matters in injections consist of mobile undissolved particles other than gas bubbles in preparations. Extraneous substances, substances derived from manufacturing processes, protein aggregates and so on may be included in therapeutic protein injections. In this chapter, Method 1 (Light Obscuration Particle Count Test) under Insoluble Particulate Matter Test for Injections <6.07> is used for the determination of insoluble particulates in therapeutic protein injections. This test is applied to the injections whose active ingredients are peptides, proteins or their derivatives.

Since this test is a sampling test conducted on a part of samples, the test must be performed under a statistically sound sampling plan in order to estimate correctly the number of particles in the population.

1. Apparatus

Use a suitable apparatus based on the principle of light blockage which allows an automatic determination of the size of particles and the number of particles according to size. Calibration, and verification of sample volume accuracy, sample flow rate accuracy and counting accuracy are performed according to Method 1 under Insoluble Particulate Matter Test for Injections <6.07>. When one measurement is performed with volume less than 1 mL, confirm the sample volume accuracy separately by an appropriate method.

2. General precautions

The test is carried out under conditions limiting particulate contamination, preferably in a laminar-flow cabinet. Very carefully wash the glassware and filtration equipment used, except for the membrane filters, with a warm detergent solution and rinse with abundant amounts of water to remove all traces of detergent. Immediately before use, rinse the equipment from top to bottom, outside and then inside, with particle-free water. Take care not to introduce air bubbles into the preparation to be examined, especially when fractions of the preparation are being transferred to the container in which the determination is to be carried out. In order to check that the environment is suitable for the test, that the glassware is properly cleaned and that the number of particles in the particle-free water to be used is within specifications, the following test is carried out using 5 mL of the particle-free water. When one measurement is performed with volume less than 1 mL, 1 mL of particlefree water may be used. Determine the particulate contamination of 5 samples of particle-free water. If the number of particles of $10 \,\mu m$ or greater size exceeds 1 per 1 mL, the precautions taken for the test are not sufficient. In this case, the preparatory steps must be repeated until the environment, filtration equipment glassware and particle-free water are suitable for the test.

3. Method

Treat a protein solution in an appropriate manner because of its tendency to generate air bubbles. In the case of an injection to be dissolved before use, use a specified solvent. When solvent is not specified, dissolve in particle-free water or use other suitable solvent comparable to particlefree water. Mix the appropriate contents of the sample gently and thoroughly by an appropriate procedure such as swirling the container slowly. If the container is sealed, cautiously remove the sealing closure. Clean the outer surfaces of the container opening using a jet of particle-free water and remove the closure, avoiding any contamination at the contents. For elimination of air bubbles, it is recommended to allow a container to stand under ambient pressure or reduced pressure. Other procedures are applicable if confirmed to be appropriate. Sonicating is not appropriate because it may aggregate or denature proteins. If necessary, after degassing, mix the contents gently and thoroughly by swirling slowly the container so as not to introduce air bubbles, and use it for the test. The measurement volume is 1 to 5 mL. The measurement volume can be reduced to 0.2 mL when the validity of the reduction is confirmed in considering the property of the sample and the tare volume of the apparatus. Set the volume necessary for the test in consideration of counting 4 portions.

In the case of injections where the volume necessary for the test can be obtained from one container, individual containers are tested. For injections with insufficient volume, combine the contents of several containers in one clean container to obtain the necessary volume after mixing the contents of containers gently and thoroughly. Where justified, the volume necessary for the test may be prepared by diluting the test solution with *particle-free water* or an appropriate solvent comparable to *particle-free water*. The validity of the dilution procedure and the solvent used for the dilution is confirmed by, for example, demonstrating consistent result regardless of the dilution. The number of test specimens must be adequate to provide a statistically sound assessment.

Take 4 portions and count the number of particles equal to or greater than 10 μ m and 25 μ m. Disregard the result obtained for the first portion, and calculate the mean number of particles for the preparation to be examined.

4. Evaluation

The preparation complies with the test if the average number of particles meets the following requirement. A—Solutions for injection supplied in containers with a nominal content of equal to or more than 100 mL.

The average number of particles of equal to or greater than $10 \, \mu m$ does not exceed 25 per mL and that of particles of equal to or greater than 25 μm does not exceed 3 per mL. B—Solutions for injection supplied in containers with a nominal content of less than 100 mL.

The average number of particles of equal to or greater than $10 \, \mu m$ does not exceed 6000 per container and that of particles of equal to or greater than $25 \, \mu m$ does not exceed

600 per container.

9.01 Reference Standards

Add the following to section (1):

Bromfenac Sodium RS
L-Carnosine RS
Doripenem RS
Gatifloxacin RS
Hydroxyethylcellulose RS for Identification
Lanoconazole RS
Microcrystalline Cellulose RS for Identification
Residual Solvents Class 2C RS
Sitagliptin Phosphate RS
Sitagliptin Phosphate RS for System Suitability

Move the following from section (2) to section (1):

Ampicillin RS
Azithromycin RS
Cefazolin RS
Cefmetazole RS
Fradiomycin Sulfate RS
Meropenem RS

9.21 Standard Solutions for Volumetric Analysis

Add the following:

Cerium (IV) Sulfate, 0.1 mol/L

1000 mL of this solution contains 40.43 g of cerium (IV) sulfate tetrahydrate [Ce $(SO_4)_2.4H_2O: 404.30$].

Preparation—Dissolve 40.43 g of cerium sulfate (IV) tetrahydrate in water to make 1000 mL, and standardize the solution as follows:

Standardization—Weigh accurately about 0.2 g of sodium oxalate (standard reagent), previously dried between 150°C and 200°C for 1 to 1.5 hours, and allow to cool in a desiccator (silica gel), and dissolve in 75 mL of water. Add a mixture of 5 mL of water and 2 mL of sulfuric acid with stirring, add 10 mL of hydrochloric acid, and warm to 70 – 75°C. Titrate <2.50> the solution with 0.1 mol/L cerium (IV) sulfate VS until the solution shows a persistent slightly yellow color, and calculate the molarity factor.

Each mL of 0.1 mol/L cerium sulfate (IV) VS = 6.700 mg of Na₂C₂O₄

Iodine, 0.025 mol/L

1000 mL of this solution contains 6.345 g of iodine (I:126.90).

Preparation—Before use, dilute 0.05 mol/L iodine VS with water to make exactly twice the initial volume.

9.22 Standard Solutions

Add the following:

Standard Chloride Solution Pipet 10 mL of Standard Chloride Stock Solution, add water to make exactly 1000 mL. Prepare before use. Each mL of this solution contains 5 μ g of chloride (Cl).

Standard Chloride Stock Solution Weigh accurately 0.842 g of sodium chloride, previously dried at 130°C for 2 hours, and dissolve in water to make exactly 1000 mL.

Standard Glyoxal Solution Dilute Standard Glyoxal Stock Solution to 10 times with ethanol (99.5). Prepare before use. Each mL of this solution contains $2 \mu g$ of glyoxal $(C_2H_2O_2)$.

Standard Glyoxal Stock Solution Transfer a quantity of 40% glyoxal TS, equivalent to 0.200 g of glyoxal, in a 100-mL volumetric flask, and dilute to 100 mL with ethanol (99.5). Dilute to 100-fold with ethanol (99.5) before use. Each mL of this solution contains $20 \,\mu g$ of glyoxal $(C_2H_2O_2)$.

9.41 Reagents, Test Solutions

Add the following:

4-Amino-6-chlorobenzene-1,3-disulfonamide

C₆H₈ClN₃O₄S₂ White, crystals or crystalline powder.

Identification—Determine the infrared absorption spectrum of 4-amino-6-chlorobenzene-1,3-disulfonamide as directed in the potassium bromide disk method under Infrared Spectrophotometry <2.25>: it exhibits absorption at the wave numbers of about 3380 cm⁻¹, 3250 cm⁻¹, 1638 cm⁻¹, 1597 cm⁻¹, 1544 cm⁻¹ and 1324 cm⁻¹.

Storage—Preserve in tight containers.

Bilirubin for assay $C_{33}H_{36}N_4O_6$ A red-orange powder. Very slightly soluble in dimethyl sulfoxide, and practically insoluble in water and in ethanol (99.5).

Identification Determine the infrared absorption spectrum of bilirubin for assay as directed in the ATR method under Infrared Spectrophotometry <2.25>: it exhibits absorption at the wave numbers of about 3400 cm⁻¹, 2910 cm⁻¹, 1686 cm⁻¹ and 1643 cm⁻¹.

Absorbance $\langle 2.24 \rangle$ $E_{1 \text{ cm}}^{1\%}$ (453 nm): 970 – 1134 (1 mg, dimethyl sulfoxide, 200 mL). Conduct this procedure without exposure to light, using light-resistant vessels.

Purity Related substances—Conduct this procedure without exposure to light, using light-resistant vessels. The following sample solution and standard solution should be prepared before use. Dissolve 5 mg of bilirubin for assay in 50 mL of a warmed mixture of dimethyl sulfoxide and acetic acid (100) (9:1), cool and use this solution as the sample solution. Pipet 1 mL of the sample solution, add a mixture of dimethyl sulfoxide and acetic acid (100) (9:1) to

make exactly 50 mL, and use this solution as the standard solution. Perform the test with exactly $10\,\mu\text{L}$ each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and determine each peak area by the automatic integration method: the total area of the peaks other than bilirubin obtained from the sample solution is not larger than the peak area of bilirubin from the standard solution. Operating conditions

Detector, column, column temperature, mobile phase and flow rate: Proceed as directed in the operating conditions in the Assay under Oriental Bezoar.

Time span of measurement: About 3 times as long as the retention time of bilirubin, beginning after the solvent peak. System suitability

System performance and system repeatability: Proceed as directed in the system suitability in the Assay under Oriental Bezoar.

Deuterated acetone for nuclear magnetic resonance spectroscopy CD₃COCD₃ Prepared for nuclear magnetic resonance spectroscopy.

Diclofenac sodium for assay $C_{14}H_{10}Cl_2NNaO_2$ [Same as the monograph Diclofenac Sodium. When dried, it contains not less than 99.0% of diclofenac sodium $(C_{14}H_{10}Cl_2NNaO_2)$.]

Diisopropyl 1,3-dithiolan-2-ylidenemalonate

 $C_{12}H_{18}O_4S_2$ White crystals.

Identification—Determine the absorption spectrum of a solution of diisopropyl 1,3-dithiolan-2-ylidenemalonate in methanol (1 in 125,000) as directed under Ultraviolet-visible Spectrophotometry <2.24>: it exhibits a maximum between 304 nm and 308 nm.

Melting point <2.60>: 54 - 57°C

Dilute phenolphthalein TS See phenolphthalein TS, dilute.

Euodia fruit [Same as the namesake monograph]

Evodiamine for assay $C_{19}H_{17}N_3O$ White to light yellow, crystals or crystalline powder. Very slightly soluble in methanol and in ethanol (99.5), and practically insoluble in water. Correct the content based on the amount (%) obtained in the Assay.

Identification—Proceed as directed in the Assay: it exhibits a double doublet-like signal equivalent to one proton around δ 2.82 ppm, signals equivalent to four protons which includes a singlet signal around δ 2.91 ppm and a multiplet signal around δ 2.90 ppm – δ 2.98 ppm, a double triplet-like signal equivalent to one proton around δ 3.23 ppm, a double doublet-like signal equivalent to one proton around δ 4.66 ppm, a singlet signal equivalent to one proton around δ 6.16 ppm, a triplet-like signal equivalent to one proton around δ 7.00 ppm, a triplet-like signal equivalent to one proton around δ 7.05 ppm, a doublet-like signal equivalent to one proton around δ 7.08 ppm, a triplet-like signal equivalent to one proton around δ 7.14 ppm, a doublet-like signal equivalent to one proton around δ 7.39 ppm, a

doublet-like signal equivalent to one proton around δ 7.51 ppm, a multiplet signal equivalent to one proton around δ 7.52 ppm and a double doublet-like signal equivalent to one proton around δ 7.83 ppm.

Unity of peak—Dissolve 1 mg of evodiamine for assay in 20 mL of methanol, and use this solution as the sample solution. Perform the test with $10 \,\mu\text{L}$ of the sample solution as directed under Liquid Chromatography $\langle 2.0I \rangle$ according to the following conditions, and compare the absorption spectra of at least 3 points including the top of evodiamine peak and around the two middle peak heights of before and after the top: no difference in form is observed among their spectra.

Operating conditions

Column, column temperature, mobile phase and flow rate: Proceed as directed in the operating conditions in the Assay (1) under Goshuyuto Extract.

Detector: A photodiode array detector (wavelength: 282 nm: spectrum range of measurement: 220 - 400 nm). System suitability

System performance: When the procedure is run with 10 μ L of the sample solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of evodiamine are not less than 5000 and not more than 1.5, respectively.

Assay—Weigh accurately 5 mg of evodiamine for assay and 1 mg of DSS- d_6 for nuclear magnetic resonance spectroscopy using an ultramicrobalance, dissolve them in 1 mL of deuterated dimethylsulfoxide for nuclear magnetic resonance spectroscopy, and use this solution as the sample solution. Transfer the sample solution into an NMR tube 5 mm in outer diameter, and measure ¹H-NMR as directed under Nuclear Magnetic Resonance Spectroscopy $\langle 2.2I \rangle$ and Crude Drugs Test $\langle 5.0I \rangle$ according to the following conditions, using DSS- d_6 for nuclear magnetic resonance spectroscopy as the internal reference compound. Calculate the resonance intensity A (equivalent to 1 hydrogen) of the signal around δ 6.16 ppm assuming the signal of the internal reference compound as δ 0 ppm.

Amount (%) of evodiamine
$$(C_{19}H_{17}N_3O)$$

= $M_S \times I \times P/(M \times N) \times 1.3521$

M: Amount (mg) of evodiamine for assay taken

 $M_{\rm S}$: Amount (mg) of DSS- d_6 for nuclear magnetic resonance spectroscopy taken

I: Signal resonance intensity A based on the signal resonance intensity of DSS- d_6 for nuclear magnetic resonance spectroscopy as 9.000

N: Number of the hydrogen derived from A

P: Purity (%) of DSS- d_6 for nuclear magnetic resonance spectroscopy

Operating conditions

Apparatus: A nuclear magnetic resonance spectrometer having ¹H resonance frequency of not less than 400 MHz.

Target nucleus: 1H.

Digital resolution: 0.25 Hz or lower.

Measuring spectrum range: 20 ppm or upper, including

2903

between -5 ppm and 15 ppm.

Spinning: off. Pulse angle: 90°. ¹³C decoupling: on.

Delay time: Repeating pulse waiting time not less than 60 seconds.

Integrating times: 8 or more times. Dummy scanning: 2 or more times.

Measuring temperature: A constant temperature between 20°C and 30°C.

System suitability

Test for required detectability: When the procedure is run with the sample solution under the above operating conditions, the SN ratio of the signal around δ 6.16 ppm is not less than 100.

System performance: When the procedure is run with the sample solution under the above operating conditions, the signals around δ 6.16 ppm is not overlapped with any signal of obvious foreign substance.

System repeatability: When the test is repeated 6 times with the sample solution under the above operating conditions, the relative standard deviation of the ratio of the resonance intensity A to that of the internal reference compound is not more than 1.0%.

Felodipine for assay [Same as the monograph Felodipine. It contains not less than 99.5% of ferodipine $(C_{18}H_{19}Cl_2NO_4)$, calculated on the dried basis.]

40% glyoxal TS *Content:* 38 – 42%. Assay—Put 1.000 g of 40% glyoxal TS in a glass-stoppered flask, add 20 mL of a solution of hydroxylammonium chloride (7 in 100) and 50 mL of water. Stopper tightly, allow to stand for 30 minutes, titrate <2.50> with 1 mol/L sodium hydroxide VS (indicator: 1.0 mL of methyl red-methylene blue TS). Perform a blank determination in the same manner, and make any necessary correction.

Each mL of 1 mol/L sodium hydroxide VS = 29.02 mg of $C_2H_2O_2$

Iodoethane for assay C_2H_5I Colorless to slightly yellow liquid, turning brown on exposure to air and light. Miscible with ethanol (95). Specific gravity d_{20}^{20} : about 1.95; Boiling point: about 72°C.

Refractive index $\langle 2.45 \rangle$ n_D^{20} : 1.509 – 1.515.

Content: not less than 99.0%. Assay—Proceed as directed in the Assay under isopropyl iodide for assay.

Each mL of 0.1 mol/L silver nitrate VS = 15.60 mg of C_2H_5I

Storage—Preserve in tight, light-resistant containers.

Iron (III) chloride-amidosulfuric acid TS Dissolve 10 g of iron (III) chloride hexahydrate and 16 g of amidosulfuric acid (standard reagent) in water to make 1000 mL.

(Z)-Ligustilide TS for thin-layer chromatography Dis-

solve 1 mg of (*Z*)-ligustilide for thin-layer chromatography in 10 mL of methanol.

Methyl 4-aminobenzoate H₂NC₆H₄COOCH₃ Pale yellow, crystals or crystalline powder.

Melting point <2.60>: 111 - 114°C

3-Methyl-2-benzothiazolonehydrazone hydrochloride monohydrate $C_8H_{10}ClN_3S.H_2O$ A white to light yellow-white crystalline powder. Melting point: about 270°C (with decomposition).

Methyl red-sodium hydroxide TS Dissolve 50 mg of methyl red in a mixture of 1.86 mL of 0.1 mol/L sodium hydroxide VS and 50 mL of ethanol (95), and add water to make 100 mL.

Nortriptyline hydrochloride $C_{19}H_{21}N.HCl$ [Same as the namesake monograph]

Nortriptyline hydrochloride for assay $C_{19}H_{21}N.HCl$ [Same as the monograph Nortriptyline Hydrochloride. When dried, it contains not less than 99.0% of nortriptyline hydrochloride ($C_{19}H_{21}N.HCl$).]

Peucedanum ledebourielloides for purity Powder of the root and rhizome of *Peucedanum ledebourielloides* K. T. Fu (*Umbelliferae*).

Identification—Place 1.0 g of peucedanum·ledebouriel-loides for purity in a glass-stoppered centrifuge tube, add 5 mL of hexane, shake for 10 minutes, centrifuge, and use the supernatant liquid as the sample solution. Perform the test with the sample solution as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 5 μ L of the sample solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of hexane, ethyl acetate and acetic acid (100) (20:10:1) to a distance of about 7 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 365 nm): blue fluorescent spots at Rf values of about 0.35 (agasyllin) and about 0.4 [xanthalin ($C_{24}H_{26}O_7$)] are observed.

Phenolphthalein TS, dilute Dissolve 0.1 g of phenolphthalein in 80 mL of ethanol (95), and add water to make 100 mL.

Platycodin D for thin-layer chromatography $C_{57}H_{92}O_{28}$ A white powder. Freely soluble in methanol.

Identification Determine the infrared absorption spectrum of platycodin D for thin-layer chromatography as directed in the potassium bromide disk method under Infrared Spectrophotometry $\langle 2.25 \rangle$: it exhibits absorption at the wave numbers of about $1734 \, \mathrm{cm}^{-1}$, $1637 \, \mathrm{cm}^{-1}$, $1385 \, \mathrm{cm}^{-1}$, $825 \, \mathrm{cm}^{-1}$ and $783 \, \mathrm{cm}^{-1}$.

Purity Related substances—Dissolve 2 mg of platycodin D for thin-layer chromatography in 2 mL of methanol, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add methanol to make exactly 50 mL, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 5 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chro-

matography. Develop the plate with a mixture of 1-propanol, ethyl acetate, water and acetic acid (100) (5:3:2:1) to a distance of about 10 cm, and air-dry the plate. Spray dilute sulfuric acid on the plate, and heat the plate at 105° C for 10 minutes: the spots other than the principal spot with an Rf value of about 0.5 from the sample solution are not more intense than the spot from the standard solution.

Sodium stearyl fumarate $C_{22}H_{39}NaO_4$ A white crystalline powder.

Identification (1) Determine the infrared absorption spectrum of sodium stearyl fumarate as directed in the potassium bromide disk method under Infrared Spectrophotometry $\langle 2.25 \rangle$: it exhibits absorption at the wave numbers of about 2950 cm⁻¹, 2920 cm⁻¹, 2850 cm⁻¹, 1720 cm⁻¹, 1610 cm⁻¹, 1313 cm⁻¹, 1186 cm⁻¹, 980 cm⁻¹, and 665 cm⁻¹.

(2) Sodium stearyl fumarate responds to Qualitative Tests $\langle 1.09 \rangle$ (1) for sodium salt.

 $\label{eq:Valsartan} \textbf{Valsartan} \quad C_{24}H_{29}N_5O_3 \ \ [\text{Same as the namesake monograph}]$

Change the following as follows:

Cineol for assay $C_{10}H_{18}O$ Clear and colorless liquid, having a characteristic aroma.

Refractive index $\langle 2.45 \rangle$ $n_{\rm D}^{20}$: 1.457 – 1.459 Specific gravity $\langle 2.56 \rangle$ d_{20}^{20} : 0.920 – 0.930

Purity Related substances—Dissolve 0.10 g of cineol for assay in 25 mL of hexane and use this solution as the sample solution. Perform the test with 2 μ L of the sample solution as directed under Gas Chromatography <2.02> according to the following conditions. Determine each peak area by the automatic integration method and calculate the amount of them by the area percentage method: the total amount of the peaks other than cineol is not more than 1.0%.

Operating conditions

Proceed the operating conditions in the Assay under Eucalyptus Oil except test for required detectability and time span of measurement.

Test for required detectability: To 1 mL of the sample solution add hexane to make 100 mL. Adjust so that the peak height of cineol obtained with $2 \mu L$ of this solution is 40 to 60% of the full scale.

Time span of measurement: About 3 times as long as the retention time of cineol, beginning after the solvent peak.

(*E*)-Ferulic acid for assay $C_{10}H_{10}O_4$ (*E*)-Ferulic acid. It meets the requirements of the following 1) (*E*)-ferulic acid for assay 1 or 2) (*E*)-ferulic acid for assay 2 (Purity value by quantitative NMR). The former is used after drying in a desiccator (silica gel) for 24 hours, and latter is used with correction for its amount based on the result obtained in the Assay.

1) (E)-Ferulic acid for assay 1

Absorbance $\langle 2.24 \rangle$ $E_{1 \text{ cm}}^{1\%}$ (320 nm): 878 – 969 (5 mg, methanol, 1000 mL).

Purity Related substances—Conduct this procedure without exposure to light, using light-resistant vessels. Dis-

solve 5 mg of (E)-ferulic acid for assay 1 in 10 mL of a mixture of water and methanol (1:1), and use this solution as the sample solution. Pipet 1 mL of the sample solution, add a mixture of water and methanol (1:1) to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly $10 \,\mu$ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and determine each peak area by the automatic integration method: the total area of the peaks other than (E)-ferulic acid obtained from the sample solution is not larger than the peak area of (E)-ferulic acid from the standard solution.

Operating conditions

Detector, column, column temperature, mobile phase and flow rate: Proceed as directed in the operating conditions in the Assay (1) under Tokishakuyakusan Extract.

Time span of measurement: About 6 times as long as the retention time of (E)-ferulic acid, beginning after the solvent peak.

System suitability

Test for required detectability: Pipet 1 mL of the standard solution, and add a mixture of water and methanol (1:1) to make exactly 20 mL. Confirm that the peak area of (*E*)-ferulic acid obtained with $10\,\mu\text{L}$ of this solution is equivalent to 3.5 to 6.5% of that with $10\,\mu\text{L}$ of the standard solution.

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of (*E*)-ferulic acid are not less than 5000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of (*E*)-ferulic acid is not more than 1.5%.

2) (E)-Ferulic acid for assay 2 (Purity value by quantitative NMR)

Unity of peak—Conduct this procedure without exposure to light, using light-resistant vessels. Dissolve 5 mg of (E)-ferulic acid for assay 2 in 10 mL of a mixture of water and methanol (1:1), and use this solution as the sample solution. Perform the test with $10 \,\mu\text{L}$ of the sample solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and compare the absorption spectra of at least 3 points including the top of (E)-ferulic acid peak and around the two middle peak heights of before and after the top: no difference in form is observed among their spectra.

Operating conditions

Column, column temperature, mobile phase and flow rate: Proceed as directed in the operating conditions in the Assay (1) under Tokishakuyakusan Extract.

Detector: A photodiode array detector (wavelength: 320 nm, measuring range of spectrum: 220 - 400 nm). System suitability

System performance: When the procedure is run with 10 μ L of the sample solution under the above operating conditions, the number of theoretical plates and the symmetry

2905

factor of the peak of (E)-ferulic acid are not less than 5000 and not more than 1.5, respectively.

Assay—Weigh accurately 5 mg of (E)-ferulic acid for assay 2 and 1 mg of 1,4-BTMSB- d_4 for nuclear magnetic resonance spectroscopy using an ultramicrobalance, dissolve in 1 mL of deuterated methanol for nuclear magnetic resonance spectroscopy, and use this solution as the sample solution. Transfer the sample solution into an NMR tube 5 mm in outer diameter, and measure ¹H-NMR as directed under Nuclear Magnetic Resonance Spectroscopy $\langle 2.21 \rangle$ and Crude Drugs Test $\langle 5.01 \rangle$ according to the following conditions, using 1,4-BTMSB- d_4 for nuclear magnetic resonance spectroscopy as the internal reference compound. Calculate the resonance intensity A (equivalent to 1 hydrogen) of the signal around δ 6.06 ppm assuming the signal of the internal reference compound as δ 0 ppm.

Amount (%) of (*E*)-ferulic acid (
$$C_{10}H_{10}O_4$$
)
= $M_S \times I \times P/(M \times N) \times 0.8573$

M: Amount (mg) of (E)-ferulic acid for assay 2 taken

 $M_{\rm S}$: Amount (mg) of 1,4-BTMSB- d_4 for nuclear magnetic resonance spectroscopy taken

- *I*: Signal resonance intensity A based on the signal resonance intensity of 1,4-BTMSB- d_4 for nuclear magnetic resonance spectroscopy as 18.000
- N: Number of hydrogen derived from A
- P: Purity (%) of 1,4-BTMSB-d₄ for nuclear magnetic resonance spectroscopy

Operating conditions

Apparatus: A nuclear magnetic resonance spectrometer having ¹H resonance frequency of not less than 400 MHz.

Target nucleus: ¹H.

Digital resolution: 0.25 Hz or lower.

Measuring spectrum range: 20 ppm or upper, including between -5 ppm and 15 ppm.

Spinning: off. Pulse angle: 90°. ¹³C decoupling: on.

Delay time: Repeating pulse waiting time not less than 60 econds.

Integrating times: 8 or more times. Dummy scanning: 2 or more times.

Measuring temperature: A constant temperature between 20°C and 30°C.

System suitability

Test for required detectability: When the procedure is run with the sample solution under the above operating conditions, the SN ratio of the signal around δ 6.06 ppm is not less than 100.

System performance: When the procedure is run with the sample solution under the above operating conditions, the signal around δ 6.06 ppm is not overlapped with any signal of obvious foreign substances.

System repeatability: When the test is repeated 6 times with the sample solution under the above operating conditions, the relative standard deviation of the ratio of the resonance intensity A to that of the internal reference com-

pound is not more than 1.0%.

10-Hydroxy-2-(E)-decenoic acid for assay $C_{10}H_{18}O_3$ 10-hydroxy-2-(E)-decenoic acid for thin-layer chromatography. It meets the requirement of the following 1) 10-Hydroxy-2-(E)-decenoic acid for assay 1 or 2) 10-Hydroxy-2-(E)-decenoic acid for assay 2 (Purity value by quantitative NMR). The latter is used with correction for its amount based on the result obtained in the Assay.

1) 10-Hydroxy-2-(E)-decenoic acid for assay 1

Purity Related substances—Dissolve 10 mg of 10-hydroxy-2-(E)-decenoic acid for assay 1 in 100 mL of methanol, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add methanol to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 10 μ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions. Determine each peak area by the automatic integration method: the total area of the peaks other than 10-hydroxy-2-(E)-decenoic acid obtained from the sample solution is not larger than the peak area of 10-hydroxy-2-(E)-decenoic acid from the standard solution.

Operating conditions

Detector, column, column temperature, mobile phase and flow rate: Proceed as directed in the operating conditions in the Assay under Royal Jelly.

Time span of measurement: About 4 times as long as the retention time of 10-hydroxy-2-(*E*)-decenoic acid, beginning after the solvent peak.

System suitability

Test for required detectability: Pipet 1 mL of the standard solution, and add methanol to make exactly 20 mL. Confirm that the peak area of 10-hydroxy-2-(E)-decenoic acid obtained with 10 μ L of this solution is equivalent to 3.5 to 6.5% of that with 10 μ L of the standard solution.

System performance: Dissolve 1 mg of propyl parahydroxybenzoate for resolution check in 10 mL of the sample solution. When the procedure is run with $10 \,\mu\text{L}$ of this solution under the above operating conditions, 10-hydroxy-2-(E)-decenoic acid and propyl parahydroxybenzoate are eluted in this order with the resolution between these peaks being not less than 6.

System repeatability: When the test is repeated 6 times with $10 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of 10-hydroxy-2-(*E*)-decenoic acid is not more than 1.5%.

2) 10-Hydroxy-2-(*E*)-decenoic acid for assay 2 (Purity value by quantitative NMR)

Unity of peak—Dissolve 1 mg of 10-hydroxy-2-(E)-decenoic acid for assay 2 in 50 mL of methanol, and use this solution as the sample solution. Perform the test with 10 μ L of the sample solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and compare the absorption spectra of at least 3 points including the top of 10-hydroxy-2-(E)-decenoic acid peak and around the two middle peak heights of before and after the

top: no difference in form is observed among their spectra.

Operating conditions

Column, column temperature, mobile phase, and flow rate: Proceed as directed in the operating conditions in the Assay under Royal Jelly.

Detector: A photodiode array detector (wavelength: 215 nm: spectrum range of measurement: 200 - 400 nm). System suitability

System performance: Dissolve 1 mg each of 10-hydroxy-2-(E)-decenoic acid for assay 2 and propyl parahydroxybenzoate for resolution check in methanol to make 50 mL. When the procedure is run with 10 μ L of this solution under the above operating conditions, 10-hydroxy-2-(E)-decenoic acid and propyl parahydroxybenzoate are eluted in this order with the resolution between these peaks being not less than 6.

Assay—Weigh accurately 5 mg of 10-hydroxy-2-(E)decenoic acid for assay 2 and 1 mg of 1,4-BTMSB-d4 for nuclear magnetic resonance spectroscopy using an ultramicrobalance, dissolve in 1 mL of deuterated methanol for nuclear magnetic resonance spectroscopy, and use this solution as the sample solution. Transfer the sample solution into an NMR tube 5 mm in outer diameter, and measure ¹H-NMR as directed under Nuclear Magnetic Resonance Spectroscopy <2.21> and Crude Drugs Test <5.01> according to the following conditions, using 1,4-BTMSB- d_4 for nuclear magnetic resonance spectroscopy as the internal reference compound. Calculate the resonance intensities, A_1 (equivalent to 1 hydrogen) and A_2 (equivalent to 1 hydrogen), of the signals around δ 5.54 ppm and δ 6.70 ppm assuming the signal of the internal reference compound as $\delta 0$ ppm.

Amount (%) of 10-hydroxy-2-(*E*)-decenoic acid ($C_{10}H_{18}O_3$) = $M_S \times I \times P/(M \times N) \times 0.8223$

M: Amount (mg) of 10-hydroxy-2-(E)-decenoic acid for assay 2 taken

 $M_{\rm S}$: Amount (mg) of 1,4-BTMSB- d_4 for nuclear magnetic resonance spectroscopy taken

I: Sum of the signal resonance intensities, A_1 and A_2 , based on the signal resonance intensity of 1,4-BTMSB- d_4 for nuclear magnetic resonance spectroscopy as 18.000

N: Sum of the numbers of the hydrogen derived from A_1 and A_2

P: Purity (%) of 1,4-BTMSB-d₄ for nuclear magnetic resonance spectroscopy

Operating conditions

Apparatus: A nuclear magnetic resonance spectrometer having ¹H resonance frequency of not less than 400 MHz.

Target nucleus: 1H.

Digital resolution: 0.25 Hz or lower.

Measuring spectrum range: 20 ppm or upper, including between -5 ppm and 15 ppm.

Spinning: off. Pulse angle: 90°. ¹³C decoupling: on. Delay time: Repeating pulse waiting time not less than 60 seconds.

Integrating times: 8 or more times.

Dummy scanning: 2 or more times.

Measuring temperature: A constant temperature between 20°C and 30°C.

System suitability

Test for required detectability: When the procedure is run with the sample solution under the above operating conditions, the SN ratio of each signal around δ 5.54 ppm and δ 6.70 ppm is not less than 100.

System performance: When the procedure is run with the sample solution under the above operating conditions, the signals around δ 5.54 ppm and δ 6.70 ppm are not overlapped with any signal of obvious foreign substances, and the ratio of the resonance intensities, A_1/A_2 , of each signal around δ 5.54 ppm and δ 6.70 ppm is between 0.99 and 1.01.

System repeatability: When the test is repeated 6 times with the sample solution under the above operating conditions, the relative standard deviations of the ratios of the resonance intensities, A_1 and A_2 , to that of the internal reference are not more than 1.0%.

Limonene $C_{10}H_{16}$ Clear and colorless liquid, having a characteristic aroma and a bitter taste.

Refractive index $\langle 2.45 \rangle$ $n_{\rm D}^{20}$: 1.472 – 1.474 Specific gravity $\langle 2.56 \rangle$ d_{20}^{20} : 0.841 – 0.846

Purity Related substances—Dissolve 0.1 g of limonene in 25 mL of hexane and use this solution as the sample solution. Perform the test with $2\,\mu\text{L}$ of the sample solution as directed under Gas Chromatography $\langle 2.02 \rangle$ according to the following conditions. Determine each peak area by the automatic integration method and calculate the amount of them by the area percentage method: the total amount of the peaks other than limonene is not more than 3.0%.

Operating conditions

Proceed the operating conditions in the Assay under Eucalyptus Oil except test for required detectability and time span of measurement.

Test for required detectability: To 1 mL of the sample solution add hexane to make 100 mL. Adjust so that the peak height of limonene obtained with $2 \mu L$ of this solution is 40 to 60% of the full scale.

Time span of measurement: About 3 times as long as the retention time of limonene, beginning after the solvent peak.

Saikosaponin a for thin-layer chromatography A white, crystalline powder or powder. Freely soluble in methanol and in ethanol (99.5), and practically insoluble in water. Melting point: 225 – 232°C (with decomposition).

Absorbance $\langle 2.24 \rangle$ $E_{1\,\text{cm}}^{1\,\%}$ (206 nm): 65 – 73 (15 mg, methanol, 200 mL). Previously dried in a desiccator (in vacuum, silica gel) for 24 hours.

Purity Related substances—Dissolve 1.0 mg of saikosaponin a for thin-layer chromatography in exactly 1 mL of methanol, and perform the test with 10 μ L of this solution as directed in the Identification (2) under Bupleurum Root: any spot other than the principal spot at the Rf value of about 0.4 does not appear.

Saikosaponin d for assay $C_{42}H_{68}O_{13}$ A white, crystalline powder or powder. Freely soluble in methanol and in ethanol (99.5), and practically insoluble in water. Melting point: about 240°C.

Absorbance $\langle 2.24 \rangle$ $E_{1 \text{ cm}}^{1\%}$ (206 nm): 66 – 74 (15 mg, methanol, 200 mL). Previously dried in a desiccator (in vacuum, silica gel) for 24 hours.

Purity Related substances—

- (1) Dissolve 2.0 mg of saikosaponin d for assay in 2 mL of methanol, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add methanol to make exactly 100 mL, and use this solution as the standard solution. Proceed the test with 10 μ L each of the sample solution and standard solution as directed in the Identification (2) under Bupleurum Root: the spot other than the principal spot around Rf value of 0.4 obtained from the sample solution is not larger and not more intense than the spot from the standard solution.
- (2) Dissolve 10 mg of saikosaponin d for assay in 20 mL of methanol, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add methanol to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine each peak area by the automatic integration method: the total area of the peaks other than saikosaponin d obtained from the sample solution is not larger than the peak area of saikosaponin d from the standard solution.

Operating conditions

Detector and column: Proceed as directed in the operating conditions in the Assay under Bupleurum Root.

Column temperature: A constant temperature of about 40°C.

Mobile phase: A mixture of water and acetonitrile (11:9). Flow rate: Adjust so that the retention time of saikosaponin d is about 13 minutes.

Time span of measurement: About 4 times as long as the retention time of saikosaponin d, beginning after the solvent peak.

System suitability

Test for required detectability: Pipet 1 mL of the standard solution, and add methanol to make exactly 20 mL. Confirm that the peak area of saikosaponin d obtained with 20 μ L of this solution is equivalent to 3.5 to 6.5% of that with 20 μ L of the standard solution.

System performance: Dissolve 6 mg each of saikosaponin d for assay and saikosaponin a for assay in methanol to make 100 mL. When the procedure is run with 20μ L of this solution under the above operating conditions, saikosaponin a and saikosaponin d are eluted in this order with the resolution between these peaks being not less than 1.5.

System repeatability: When the test is repeated 6 times with 20 μ L of the standard solution under the above operat-

ing conditions, the relative standard deviation of the peak area of saikosaponin d is not more than 1.0%.

Sinomenine for assay $C_{19}H_{23}NO_4$ Sinomenine for thinlayer chromatography. It meets the requirements of the following 1) sinomenine for assay 1 or 2) sinomenine for assay 2 (Purity value by quantitative NMR). The former is used after drying in a desiccator (silica gel) for 24 hours, and latter is used with correction for its amount based on the result obtained in the Assay.

1) Sinomenine for assay 1

Identification Determine the absorption spectrum of a solution of sinomenine for assay in methanol (1 in 25,000) as directed under Ultraviolet-visible Spectrophotometry <2.24>: it exhibits the maximum between 259 nm and 263 nm.

Purity Related substances—Dissolve 5 mg of sinomenine for assay 1 in 10 mL of a mixture of water and acetonitrile (7:3), and use this solution as the sample solution. Pipet 1 mL of the sample solution, add the mixture of water and acetonitrile (7:3) to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly $10 \,\mu$ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine each peak area by the automatic integration method: the total area of the peaks other than sinomenine obtained from the sample solution is not larger than the peak area of sinomenine from the standard solution.

Operating conditions

Column, column temperature, mobile phase and flow rate: Proceed as directed in the operating conditions in the Assay (1) under Boiogito Extract.

Detector: An ultraviolet absorption photometer (wavelength: 261 nm).

Time span of measurement: About 3 times as long as the retention time of sinomenine, beginning after the solvent peak.

System suitability

Test for required detectability: Pipet 1 mL of the standard solution, and add a mixture of water and acetonitrile (7:3) to make exactly 20 mL. Confirm that the peak area of sinomenine obtained with $10 \,\mu\text{L}$ of this solution is equivalent to 3.5 to 6.5% of that with $10 \,\mu\text{L}$ of the standard solution

System performance: When the procedure is run with $10 \mu L$ of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of sinomenine are not less than 5000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of sinomenine is not more than 1.5%.

2) Sinomenine for assay 2 (Purity value by quantitative NMR)

Unity of peak—Dissolve 5 mg of sinomenine for assay 2 in 10 mL of a mixture of water and acetonitrile (7:3), and

use this solution as the sample solution. Perform the test with $10\,\mu\text{L}$ of the sample solution as directed under Liquid Chromatography <2.01> according to the following conditions, and compare the absorption spectra of at least 3 points including the top of sinomenine peak and around the two middle peak heights of before and after the top: no difference in form is observed among their spectra.

Operating conditions

Column, column temperature, mobile phase and flow rate: Proceed as directed in the operating conditions in the Assay (1) under Boiogito Extract.

Detector: A photodiode array detector (wavelength: 261 nm, spectrum range of requirement: 220 - 400 nm). System suitability

System performance: When the procedure is run with 10 μ L of the sample solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of sinomenine are not less than 5000 and not more than 1.5, respectively.

Assay—Weigh accurately 5 mg of sinomenine for assay 2 and 1 mg of 1,4-BTMSB- d_4 for nuclear magnetic resonance spectroscopy using an ultramicrobalance, dissolve in 1 mL of deuterated acetone for nuclear magnetic resonance spectroscopy, and use this solution as the sample solution. Transfer the sample solution into an NMR tube 5 mm in outer diameter, and measure ¹H-NMR as directed under Nuclear Magnetic Resonance Spectroscopy $\langle 2.21 \rangle$ and Crude Drugs Test $\langle 5.01 \rangle$ according to the following conditions, using 1,4-BTMSB- d_4 for nuclear magnetic resonance spectroscopy as the internal reference compound. Calculate the resonance intensity A (equivalent to 1 hydrogen) of the signal around δ 5.42 ppm assuming the signal of the internal reference compound as δ 0 ppm.

Amount (%) of sinomenine
$$(C_{19}H_{23}NO_4)$$

= $M_S \times I \times P/(M \times N) \times 1.4543$

M: Amount (mg) of sinomenine for assay 2 taken

 $M_{\rm S}$: Amount (mg) of 1,4-BTMSB- d_4 for nuclear magnetic resonance spectroscopy taken

I: Signal resonance intensity A based on the signal resonance intensity of 1,4-BTMSB- d_4 for nuclear magnetic resonance spectroscopy as 18.000

N: Number of hydrogen derived from A

P: Purity (%) of 1,4-BTMSB-d₄ for nuclear magnetic resonance spectroscopy

Operating conditions

Apparatus: A nuclear magnetic resonance spectrometer having ¹H resonance frequency of not less than 400 MHz.

Target nucleus: 1H.

Digital resolution: 0.25 Hz or lower.

Measuring spectrum range: 20 ppm or upper, including between -5 ppm and 15 ppm.

Spinning: off. Pulse angle: 90°. ¹³C decoupling: on.

Delay time: Repeating pulse waiting time not less than 60 seconds.

Integrating times: 8 or more times. Dummy scanning: 2 or more times.

Measuring temperature: A constant temperature between 20° C and 30° C.

System suitability

Test for required detectability: When the procedure is run with the sample solution under the above operating conditions, the SN ratio of the signal around δ 5.42 ppm is not less than 100.

System performance: When the procedure is run with the sample solution under the above operating conditions, the signal around δ 5.42 ppm is not overlapped with any signal of obvious foreign substances.

System repeatability: When the test is repeated 6 times with the sample solution under the above operating conditions, the relative standard deviation of the ratio of the resonance intensity A to that of the internal reference compound is not more than 1.0%.

Starch, soluble [K 8659, Starch, soluble, First class]

9.42 Solid Supports/Column Packings for Chromatography

Add the following:

Amylose tris-(3,5-dimethylphenylcarbamate)-coated silica gel for liquid chromatography Prepared for liquid chromatography.

Octadecyl-strong anion exchange-silanized silica gel for liquid chromatography Prepared for liquid chromatography.

Porous silica gel for gas chromatography Prepared for gas chromatography.

Official Monographs

Amphotericin B Tablets

アムホテリシン **B** 錠

Add the following next to the Uniformity of dosage units:

Disintegration <6.09> Perform the test using the disk: it meets the requirement.

Beclometasone Dipropionate

ベクロメタゾンプロピオン酸エステル

Change the Description and Optical rotation as follows:

Description Beclometasone Dipropionate occurs as a white to pale yellow powder.

It is soluble in methanol, sparingly soluble in ethanol (99.5), and practically insoluble in water.

Melting point: about 208°C (with decomposition).

It shows crystal polymorphism.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁵: +106 - +114° (after drying, 0.1 g, ethanol (99.5), 10 mL, 100 mm).

Betamethasone Dipropionate

ベタメタゾンジプロピオン酸エステル

Change the Description and Optical rotation as follows:

Description Betamethasone Dipropionate occurs as a white to pale yellow-white crystalline powder. It is odorless.

It is freely soluble in acetone and in chloroform, soluble in methanol and in ethanol (99.5), and practically insoluble in water.

It is gradually affected by light.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁵: $+84 - +89^{\circ}$ (after drying, 50 mg, ethanol (99.5), 10 mL, 100 mm).

Add the following:

Bromfenac Sodium Hydrate

ブロムフェナクナトリウム水和物

 $C_{15}H_{11}BrNNaO_3.1\frac{1}{2}H_2O$: 383.17 Sodium 2-[2-amino-3-(4-bromobenzoyl)phenyl]acetate sesquihydrate [120638-55-3]

Bromfenac Sodium Hydrate contains not less than 97.5% and not more than 101.5% of bromfenac sodium ($C_{15}H_{11}BrNNaO_3$: 356.15), calculated on the anhydrous basis.

Description Bromfenac Sodium Hydrate occurs as a yellow to orange crystalline powder.

It is freely soluble in water, soluble in methanol, and slightly soluble in ethanol (99.5).

It dissolves in a solution of sodium hydrogen carbonate (21 in 2500).

Identification (1) Dissolve 10 mg of Bromfenac Sodium Hydrate in 500 mL of a solution of sodium hydrogen carbonate (21 in 2500). Determine the absorption spectrum of this solution as directed under Ultraviolet-visible Spectrophotometry $\langle 2.24 \rangle$, and compare the spectrum with the Reference Spectrum or the spectrum of a solution of Bromfenac Sodium RS prepared in the same manner as the sample solution: both spectra exhibit similar intensities of absorption at the same wavelengths.

- (2) Determine the infrared absorption spectrum of Bromfenac Sodium Hydrate as directed in the potassium bromide disk method under Infrared Spectrophotometry <2.25>, and compare the spectrum with the Reference Spectrum or the spectrum of Bromfenac Sodium RS: both spectra exhibit similar intensities of absorption at the same wave numbers.
- (3) A solution of Bromfenac Sodium Hydrate (1 in 20) responds to Qualitative Tests <1.09> (1) for sodium salt.

pH $\langle 2.54 \rangle$ Dissolve 1.0 g of Bromfenac Sodium Hydrate in 20 mL of water: the pH of the solution is between 8.4 and 10.2.

Purity (1) Heavy metals <1.07>—Proceed with 1.0 g of Bromfenac Sodium Hydrate according to Method 2, and perform the test. Prepare the control solution with 2.0 mL

of Standard Lead Solution (not more than 20 ppm).

(2) Related substances—Dissolve 50 mg of Bromfenac Sodium Hydrate in 100 mL of methanol, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add methanol to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine each peak area by the automatic integration method: the area of the peak other than bromfenac obtained from the sample solution is not larger than 1/10 times the peak area of bromfenac from the standard solution, and the total area of the peaks other than bromfenac from the sample solution is not larger than the peak area of bromfenac from the standard solution.

Operating conditions—

Detector, column and column temperature: Proceed as directed in the operating conditions in the Assay.

Mobile phase: Dissolve 3.85 g of ammonium acetate in 1000 mL of water, and adjust to pH 4.0 with acetic acid (100). To 570 mL of this solution add 430 mL of acetonitrile.

Flow rate: Adjust so that the retention time of bromfenac is about 8 minutes.

Time span of measurement: About 3 times as long as the retention time of bromfenac, beginning after the solvent peak.

System suitability—

Test for required detectability: Pipet 2 mL of the standard solution, and add methanol to make exactly 20 mL. Confirm that the peak area of bromfenac obtained with 20 μ L of this solution is equivalent to 7 to 13% of that with 20 μ L of the standard solution.

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of bromfenac are not less than 5000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with $20 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of bromfenac is not more than 2.0%.

Water $\langle 2.48 \rangle$ 6.9 - 8.5% (0.15 g, volumetric titration, direct titration. Use a solution of imidazole for water determination in methanol for water determination (1 in 80) instead of methanol for water determination).

Assay Weigh accurately about 30 mg each of Bromfenac Sodium Hydrate and Bromfenac Sodium RS (separately determine the water $\langle 2.48 \rangle$ in the same manner as Bromfenac Sodium Hydrate), dissolve each in methanol to make exactly 50 mL. Pipet 5 mL each of these solutions, add the mobile phase to make them exactly 100 mL, and use these solutions as the sample solution and the standard solution, respectively. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following

conditions, and determine the peak areas, A_T and A_S , of bromfenac in each solution.

Amount (mg) of bromfenac sodium ($C_{15}H_{11}BrNNaO_3$) = $M_S \times A_T/A_S$

 M_S : Amount (mg) of Bromfenac Sodium RS taken, calculated on the anhydrous basis

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 266 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 35° C.

Mobile phase: Dissolve 3.85 g of ammonium acetate in 1000 mL of water. To 600 mL of this solution add 250 mL of methanol and 150 mL of tetrahydrofuran.

Flow rate: Adjust so that the retention time of bromfenac is about 9 minutes.

System suitability-

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of bromfenac are not less than 5000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of bromfenac is not more than 1.0%.

Containers and storage Containers—Tight containers. Storage—Light-resistant.

Add the following:

Bromfenac Sodium Ophthalmic Solution

ブロムフェナクナトリウム点眼液

Bromfenac Sodium Ophthalmic Solution is an aqueous ophthalmic preparation.

It contains not less than 90.0% and not more than 110.0% of the labeled amount of Bromfenac Sodium Hydrate ($C_{15}H_{11}BrNNaO_3.1\frac{1}{2}H_2O$: 383.17).

Method of preparation Prepare as directed under Ophthalmic Liquids and Solutions, with Bromfenac Sodium Hydrate.

Description Bromfenac Sodium Ophthalmic Solution occurs as a clear and yellow liquid.

Identification To a volume of Bromfenac Sodium Ophthalmic Solution, equivalent to 1 mg of Bromfenac Sodium Hydrate, add a solution of sodium hydrogen carbonate (21)

in 2500) to make 50 mL, and determine the absorption spectrum of this solution as directed under Ultraviolet-visible Spectrophotometry <2.24>: it exhibits maxima between 266 nm and 270 nm, and between 377 nm and 381 nm.

pH Being specified separately when the drug is granted approval based on the Law.

Purity Related substances—Being specified separately when the drug is granted approval based on the Law.

Foreign insoluble matter <6.11> It meets the requirement.

Insoluble particulate matter <6.08> It meets the requirement.

Sterility <4.06> Perform the test according to the Membrane filtration method: it meets the requirement.

Assay Pipet a volume of Bromfenac Sodium Ophthalmic Solution, equivalent to about 2 mg of bromfenac sodium hydrate ($C_{15}H_{11}BrNNaO_3.1\frac{1}{2}H_2O$), add the mobile phase to make exactly 20 mL, and use this solution as the sample solution. Separately, weigh accurately about 20 mg of Bromfenac Sodium RS (separately determine the water <2.48> in the same manner as Bromfenac Sodium Hydrate), and dissolve in the mobile phase to make exactly 20 mL. Pipet 2 mL of this solution, add the mobile phase to make exactly 20 mL, and use this solution as the standard solution. Perform the test with exactly $10 \,\mu\text{L}$ each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, A_T and A_S , of bromfenac in each solution.

Amount (mg) of bromfenac sodium hydrate $(C_{15}H_{11}BrNNaO_3.1\frac{1}{2}H_2O)$ = $M_S \times A_T/A_S \times 1/10 \times 1.076$

 $M_{\rm S}$: Amount (mg) of Bromfenac Sodium RS taken, calculated on the anhydrous basis

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 266 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 25 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 40°C.

Mobile phase: Dissolve 1.98 g of diammonium hydrogen phosphate in 750 mL of water, adjust to pH 7.3 with phosphoric acid, and add 250 mL of acetonitrile.

Flow rate: Adjust so that the retention time of bromfenac is about 18 minutes.

System suitability—

System performance: When the procedure is run with $10 \mu L$ of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of bromfenac are not less than 13,000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with $10 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of bromfenac is not more than 1.0%.

Containers and storage Containers—Tight containers.

Change the following as follows:

Anhydrous Dibasic Calcium Phosphate

無水リン酸水素カルシウム

CaHPO₄: 136.06 [7757-93-9]

This monograph is harmonized with the European Pharmacopoeia and the U.S. Pharmacopeia.

The corresponding part of the attributes/provisions which are agreed as non-harmonized within the scope of the harmonization is marked with symbols ($^{\blacklozenge}$ $_{\spadesuit}$), and the corresponding parts which are agreed as the JP local requirement other than the scope of the harmonization are marked with symbols ($^{\diamondsuit}$ $_{\diamondsuit}$).

Anhydrous Dibasic Calcium Phosphate contains not less than 97.5% and not more than 102.5% of dibasic calcium phosphate (CaHPO₄).

◆Description Anhydrous Dibasic Calcium Phosphate occurs as white, crystalline powder or granules.

It is practically insoluble in water and in ethanol (99.5). It dissolves in dilute hydrochloric acid and in dilute nitric acid. ♠

Identification (1) Dissolve 0.1 g of Anhydrous Dibasic Calcium Phosphate in 10 mL of 2 mol/L hydrochloric acid TS by warming, add 2.5 mL of ammonia TS dropwise with shaking, and add 5 mL of ammonium oxalate TS: a white precipitate is produced.

(2) Dissolve 0.1 g of Anhydrous Dibasic Calcium Phosphate in 5 mL of dilute nitric acid, and add 2 mL of hexaammonium heptamolybdate TS after warming at 70°C for 1 to 2 minutes: a yellow precipitate is produced.

Purity (1) Acid-insoluble substances—Dissolve 5.0 g of Anhydrous Dibasic Calcium Phosphate in 40 mL of water and 10 mL of hydrochloric acid, and boil gently for 5 minutes. After cooling, collect the insoluble substance using a filter paper for assay. Wash with water until no more turbidity of the washings is produced when silver nitrate TS is added. Ignite to incinerate the residue and the filter paper at 600 ± 50 °C: the mass is not more than 10 mg (not more than 0.2%).

(2) Chloride—To 0.20 g of Anhydrous Dibasic Calcium Phosphate add 20 mL of water and 13 mL of dilute nitric acid, dissolve by warming, if necessary, add water to make

100 mL, and filter, if necessary. Put 50 mL of this solution in a Nessler tube, and use this as the sample solution. Separately, transfer 0.70 mL of 0.01 mol/L hydrochloric acid VS to another Nessler tube, add 6 mL of dilute nitric acid and water to make 50 mL, and use this solution as the control solution. Add 1 mL of silver nitrate TS to the sample solution and control solution, mix well, and allow to stand for 5 minutes protecting from light. Compare the opalescence developed in both solutions against a black background by viewing down-ward or transversely. The opalescence developed in the sample solution is not more than that of the control solution (not more than 0.25%).

- (3) Sulfate—Dissolve 0.50 g of Anhydrous Dibasic Calcium Phosphate in 5 mL of water and 5 mL of dilute hydrochloric acid, add water to make 100 mL, and filter, if necessary. Put 20 mL of this solution in a Nessler tube, add 1 mL of dilute hydrochloric acid, and add water to make 50 mL, and use this as the sample solution. Transfer 1.0 mL of 0.005 mol/L sulfuric acid VS to another Nessler tube, add 1 mL of dilute hydrochloride acid and water to make 50 mL, and use this solution as the control solution. Add 2 mL of barium chloride TS to the sample solution and control solution, mix well, and allow to stand for 10 minutes. Compare the white turbidity produced in both solutions against a black background by viewing downward or transversely. The turbidity produced in the sample solution is not thicker than that of the control solution (not more than 0.48%).
- (4) Carbonate—Mix 1.0 g of Anhydrous Dibasic Calcium Phosphate with 5 mL of freshly boiled and cooled water, and add immediately 2 mL of hydrochloric acid: no effervescence occurs.
- ◇(5) Heavy metals <1.07>—Dissolve 0.65 g of Anhydrous Dibasic Calcium Phosphate in a mixture of 5 mL of water and 5 mL of dilute hydrochloric acid by warming, cool, and add ammonia TS until precipitates begin to form in the solution. Dissolve the precipitates by adding a small amount of dilute hydrochloric acid dropwise, add 10 mL of hydrochloric acid-ammonium acetate buffer solution (pH 3.5) and water to make 50 mL, and perform the test using this solution as the sample solution. Prepare the control solution as follows: to 10 mL of hydrochloric acid-ammonium acetate buffer solution (pH 3.5) add 2.0 mL of Standard Lead Solution and water to make 50 mL (not more than 31 ppm).
- (6) Barium—Heat 0.5 g of Anhydrous Dibasic Calcium Phosphate with 10 mL of water, add 1 mL of hydrochloric acid dropwise with stirring, and filter after cooling, if necessary. Add 2 mL of potassium sulfate TS to this solution, and allow to stand for 10 minutes: no turbidity forms.
- **♦**(7) Arsenic <1.11>—Dissolve 1.0 g of Anhydrous Dibasic Calcium Phosphate in 5 mL of dilute hydrochloric acid, and perform the test with this solution as the test solution (not more than 2 ppm).**♦**

Loss on ignition $\langle 2.43 \rangle$ Not less than 6.6% and not more than 8.7% (1 g, 800 – 825°C, constant mass).

Assay Weigh accurately about 0.4 g of Anhydrous Dibasic Calcium Phosphate, dissolve in 12 mL of dilute hydrochlo-

ric acid by heating on a water bath, if necessary, and add water to make exactly 200 mL. Pipet 20 mL of this solution, add exactly 25 mL of 0.02 mol/L disodium dihydrogen ethylenediamine tetraacetate VS, 50 mL of water and 5 mL of ammonia-ammonium chloride buffer solution (pH 10.7), and titrate <2.50> the excess disodium dihydrogen ethylenediamine tetraacetate with 0.02 mol/L zinc sulfate VS (indicator: 25 mg of eriochrome black T-sodium chloride indicator). Perform a blank determination in the same manner.

Each mL of 0.02 mol/L disodium dihydrogen ethylenediamine tetraacetate VS = 2.721 mg of CaHPO₄

Containers and storage Containers—Well-closed containers.◆

Add the following:

Cefalotin Sodium for Injection

注射用セファロチンナトリウム

Cefalotin Sodium for Injection is a preparation for injection, which is dissolved before use.

It contains not less than 90.0% and not more than 110.0% of the labeled potency of cefalotin $(C_{16}H_{16}N_2O_6S_2: 396.44)$.

Method of preparation Prepare as directed under Injections, with Cefalotin Sodium.

Description Cefalotin Sodium for Injection occurs as white to light yellow-white, crystals or crystalline powder.

Identification Determine the infrared absorption spectrum of Cefalotin Sodium for Injection as directed in the potassium bromide disk method under Infrared Spectrophotometry <2.25>, and compare the spectrum with the Reference Spectrum of Cefalotin Sodium or the spectrum of Cefalotin Sodium RS: both spectra exhibit similar intensities of absorption at the same wave numbers.

- pH < 2.54 Dissolve an amount of Cefalotin Sodium for Injection, equivalent to 0.5 g (potency) of Cefalotin Sodium, in 5 mL of water: the pH of the solution is between 4.5 and 7.0.
- **Purity** (1) Clarity and color of solution—Dissolve 1.0 g of Cefalotin Sodium for Injection in 10 mL of water: the solution is clear. Perform the test with this solution as directed under Ultraviolet-visible Spectrophotometry <2.24>: the absorbance at 450 nm is not more than 0.20.
- (2) Related substances—Dissolve an amount of Cefalotin Sodium for Injection, equivalent to 25 mg (potency), in the mobile phase to make 25 mL, and use this solution as the sample solution. Pipet 1 mL of the sample solution, and add the mobile phase to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly $10 \, \mu$ L each of the sample solution and standard solu-

tion as directed under Liquid Chromatography <2.01> according to the following conditions, and determine each peak area by the automatic integration method: the area of the peak other than cefalotin obtained from the sample solution is not larger than the peak area of cefalotin from the standard solution, and the total area of the peaks other than cefalotin from the sample solution is not larger than 3 times the peak area of cefalotin from the standard solution.

Operating conditions—

Detector, column, column temperature, mobile phase and flow rate: Proceed as directed in the operating conditions in the Assay under Cefalotin Sodium.

Time span of measurement: About 4 times as long as the retention time of cefalotin.

System suitability—

Proceed as directed in the system suitability in the Purity (4) under Cefalotin Sodium.

Water $\langle 2.48 \rangle$ Not more than 1.0% (0.5 g, volumetric titration, back titration).

Bacterial endotoxins <4.01> Less than 0.2 EU/mg (potency).

Uniformity of dosage units < 6.02> It meets the requirement of the Mass variation test.

Foreign insoluble matter <6.06> Perform the test according to Method 2: it meets the requirement.

Insoluble particulate matter < 6.07> It meets the requirement.

Sterility <4.06> Perform the test according to the Membrane filtration method: it meets the requirement.

Assay Weigh accurately the mass of the contents of not less than 10 containers of Cefalotin Sodium for Injection. Weigh accurately an amount of the contents, equivalent to about 25 mg (potency) of Cefalotin Sodium, dissolve in the mobile phase to make exactly 25 mL, and use this solution as the sample solution. Separately, weigh accurately about 25 mg (potency) of Cefalotin Sodium RS, and dissolve in the mobile phase to make exactly 25 mL, and use this solution as the standard solution. Then, proceed as directed in the Assay under Cefalotin Sodium.

Amount [μ g(potency)] of cefalotin ($C_{16}H_{16}N_2O_6S_2$) = $M_S \times A_T/A_S \times 1000$

 M_S : Amount [mg(potency)] of Cefalotin Sodium RS taken

Containers and storage Containers—Hermetic containers.

Add the following:

Cefixime Fine Granules

セフィキシム細粒

Cefixime Fine Granules contain not less than 90.0% and not more than 105.0% of the labeled potency of cefixime ($C_{16}H_{15}N_5O_7S_2$: 453.45).

Method of preparation Prepare as directed under Granules, with Cefixime Hydrate.

Identification To a quantity of powdered Cefixime Fine Granules, equivalent to 2 mg (potency) of Cefixime Hydrate add 150 mL of 0.1 mol/L phosphate buffer solution (pH 7.0), and shake. If necessary, filter or centrifuge. Determine the absorption spectrum of this solution as directed under Ultraviolet-visible Spectrophotometry <2.24>: it exhibits maximum between 286 nm and 290 nm.

Purity Related substances—To a quantity of powdered Cefixime Fine Granules, equivalent to 0.1 g (potency) of Cefixime Hydrate, add 100 mL of 0.1 mol/L phosphate buffer solution (pH 7.0), shake, filter through a membrane filter with a pore size not exceeding 0.45 μ m, and use the filtrate as the sample solution. Perform the test with 10 μ L of the sample solution as directed under Liquid Chromatography <2.01> according to the following conditions. Determine each peak area by the automatic integration method, and calculate their amounts by the area percentage method: the amount of the peak other than cefixime is not more than 1.0%, and the total amount of the peaks other than cefixime is not more than 2.5%.

Operating conditions—

Detector, column, column temperature, mobile phase and flow rate: Proceed as directed in the operating conditions in the Assay under Cefixime Hydrate.

Time span of measurement: Proceed as directed in the operating conditions in the Purity under Cefixime Hydrate. System suitability—

Test for required detectability: To 1 mL of the sample solution add 0.1 mol/L phosphate buffer solution (pH 7.0) to make 100 mL, and use this solution as the solution for system suitability test. Pipet 1 mL of the solution for system suitability test, and add 0.1 mol/L phosphate buffer solution (pH 7.0) to make exactly 10 mL. Confirm that the peak area of cefixime obtained with 10 μ L of this solution is equivalent to 7 to 13% of that with 10 μ L of the solution for system suitability test.

System performance: When the procedure is run with $10 \mu L$ of the solution for system suitability test under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of cefixime are not less than 4000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the solution for system suitability test under the above operating conditions, the relative standard deviation of the peak area of cefixime is not more than 2.0%.

Water <2.48> Not more than 3.0% (1 g, volumetric titration, direct titration. Use a mixture of formamide for water determination and methanol for water determination (2:1) instead of methanol for water determination).

Uniformity of dosage units <6.02> Perform the test according to the following method: Cefixime Fine Granules in single-dose packages meet the requirement of the Content uniformity test.

To the total content of 1 package of Cefixime Fine Granules add 7 V/10 mL of 0.1 mol/L phosphate buffer solution (pH 7.0), shake, add 0.1 mol/L phosphate buffer solution (pH 7.0) to make exactly V mL so that each mL contains about 1 mg (potency) of Cefixime Hydrate. Centrifuge this solution, pipet 10 mL of the supernatant liquid, add 0.1 mol/L phosphate buffer solution (pH 7.0) to make exactly 50 mL, and use this solution as the sample solution. Separately, weigh accurately an amount of Cefixime RS, equivalent to about 20 mg (potency), dissolve in 0.1 mol/L phosphate buffer solution (pH 7.0) to make exactly 100 mL, and use this solution as the standard solution. Then, proceed as directed in the Assay under Cefixime Hydrate.

Amount [mg (potency)] of cefixime ($C_{16}H_{15}N_5O_7S_2$) = $M_S \times A_T/A_S \times V/20$

M_S: Amount [mg (potency)] of Cefixime RS taken

Dissolution <6.10> When the test is performed at 50 revolutions per minute according to the Paddle method, using 900 mL of 2nd fluid for dissolution test as the dissolution medium, the dissolution rates in 30 minutes of Cefixime Fine Granules is not less than 75%.

Start the test with an accurately weighed amount of Cefixime Fine Granules, equivalent to about 0.1 g (potency) of Cefixime Hydrate, withdraw not less than 20 mL of the medium at the specified minute after starting the test, and filter through a membrane filter with a pore size not exceeding 0.45 µm. Discard not less than 10 mL of the first filtrate, and use the subsequent filtrate as the sample solution. Separately, weigh accurately an amount of Cefixime RS, equivalent to about 28 mg (potency), and dissolve in the dissolution medium to make exactly 50 mL. Pipet 4 mL of this solution, add the dissolution medium to make exactly 20 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, A_T and A_S , of cefixime in each solu-

Dissolution rate (%) with respect to the labeled amount of cefixime ($C_{16}H_{15}N_5O_7S_2$)

 $= M_{\rm S}/M_{\rm T} \times A_{\rm T}/A_{\rm S} \times 1/C \times 360$

M_S: Amount [mg (potency)] of Cefixime RS taken
M_T: Amount (g) of Cefixime Fine Granules taken
C: Labeled amount [mg (potency)] of cefixime (C₁₆H₁₅N₅O₇S₂) in 1 g

Operating conditions—

Proceed as directed in the operating conditions in the Assay under Cefixime Hydrate.

System suitability—

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of cefixime are not less than 4000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of cefixime is not more than 2.0%.

Assay Weigh accurately an amount of powdered Cefixime Fine Granules, equivalent to about 0.1 g (potency) of Cefixime Hydrate, add 70 mL of 0.1 mol/L phosphate buffer solution (pH 7.0), shake, and add 0.1 mol/L phosphate buffer solution (pH 7.0) to make exactly 100 mL. Centrifuge this solution, pipet 10 mL of the supernatant liquid, add 0.1 mol/L phosphate buffer solution (pH 7.0) to make exactly 50 mL, and use this solution as the sample solution. Separately, weigh accurately an amount of Cefixime RS, equivalent to about 20 mg (potency), dissolve in 0.1 mol/L phosphate buffer solution (pH 7.0) to make exactly 100 mL, and use this solution as the standard solution. Then, proceed as directed in the Assay under Cefixime Hydrate.

Amount [mg (potency)] of cefixime $(C_{16}H_{15}N_5O_7S_2)$ = $M_S \times A_T/A_S \times 5$

M_S: Amount [mg (potency)] of Cefixime RS taken

Containers and storage Containers—Tight containers.

Change the following as follows:

Microcrystalline Cellulose

結晶セルロース

[9004-34-6, cellulose]

This monograph is harmonized with the European Pharmacopoeia and the U.S. Pharmacopeia.

The corresponding part of the attributes/provisions which are agreed as non-harmonized within the scope of the harmonization is marked with symbols ($^{\blacklozenge}$ $_{\blacklozenge}$), and the corresponding parts which are agreed as the JP local requirement other than the scope of the harmonization are marked with symbols ($^{\diamondsuit}$ $_{\diamondsuit}$).

Microcrystalline Cellulose is purified, partially depolymerized α -cellulose, obtained as a pulp from fibrous plant material, with mineral acids.

The label indicates the ^{\(\)}mean degree of polymerization, loss on drying, ^{\(\)} and bulk density values with a range.

◆Description Microcrystalline Cellulose occurs as a white crystalline powder having fluidity.

It is practically insoluble in water, in ethanol (95) and in diethyl ether.

It swells with sodium hydroxide TS on heating.

- **Identification** (1) Dissolve 20 g of zinc chloride and 6.5 g of potassium iodide in 10.5 mL of water, add 0.5 g of iodine, and shake for 15 minutes. Place about 10 mg of Microcrystalline Cellulose on a watch glass, and disperse in 2 mL of this solution: the substance develops a blue-violet color.
- (2) Determine the infrared absorption spectrum of Microcrystalline Cellulose as directed in the ATR method under Infrared Spectrophotometry <2.25>, and compare the spectrum with the spectrum of Microcrystalline Cellulose RS for Identification: both spectra exhibit similar intensities of absorption at the same wave numbers. If there are absorptions between 800 cm⁻¹ and 825 cm⁻¹, and between 950 cm⁻¹ and 1000 cm⁻¹, disregard the absorptions.
- (3) Transfer about 1.3 g of Microcrystalline Cellulose, accurately weighed, to a 125-mL conical flask, and add exactly 25 mL each of water and 1 mol/L cupriethylenediamine TS. Immediately purge the solution with nitrogen, insert the stopper, and shake on a suitable mechanical shaker to dissolve. Perform the test with a suitable amount of this solution, taken exactly, according to Method 1 under Viscosity Determination $\langle 2.53 \rangle$ using a capillary viscometer having the viscosity constant (K) of approximately 0.03, at 25 ± 0.1 °C, and determine the kinematic viscosity, ν . Separately, perform the test with a mixture of exactly 25 mL each of water and 1 mol/L cupriethylenediamine TS in the same manner as above, using a capillary viscometer having K of approximately 0.01, and determine the kinematic viscosity, ν_0 .

Calculate the relative viscosity, η_{rel} , of Microcrystalline Cellulose by the following formula:

$$\eta_{\rm rel} = \nu/\nu_0$$

Obtain the product, $[\eta]C$, of intrinsic viscosity $[\eta](\text{mL/g})$ and concentration C (g/100 mL) from the value η_{rel} of the table. When calculate the degree of polymerization, P, by the following formula, P is not more than 350 $^{\lozenge}$ and within the labeled range. $_{\lozenge}$

$$P = 95(\eta)C/M_{\rm T}$$

 $M_{\rm T}$: Amount (g) of the Microcrystalline Cellulose taken, calculated on the dried basis

- **pH** <2.54> Shake 5.0 g of Microcrystalline Cellulose with 40 mL of water for 20 minutes, and centrifuge: the pH of the supernatant liquid is between 5.0 and 7.5.
- **Purity** [⋄](1) Heavy metals <1.07>—Proceed with 2.0 g of Microcrystalline Cellulose according to Method 2, and perform the test. Prepare the control solution with 2.0 mL of Standard Lead Solution (not more than 10 ppm)._⋄
- (2) Water-soluble substances—Shake 5.0 g of Microcrystalline Cellulose with 80 mL of water for 10 minutes,

filter with the aid of vacuum through a filter paper for quantitative analysis (5C) into a vacuum flask. Evaporate the filtrate in a tared evaporating dish to dryness without charring, dry at 105°C for 1 hour, cool in a desiccator, and weigh: the difference between the mass of the residue and the mass obtained from a blank determination does not exceed 12.5 mg.


(3) Diethyl ether-soluble substances—Place 10.0 g of Microcrystalline Cellulose in a column having an internal diameter of about 20 mm, and pass 50 mL of peroxide-free diethyl ether through the column. Evaporate the eluate to dryness in a previously dried and tared evaporation dish. Dry the residue at 105°C for 30 minutes, allow to cool in a desiccator, and weigh: the difference between the mass of the residue and the mass obtained from a blank determination does not exceed 5.0 mg.

Conductivity $\langle 2.51 \rangle$ Perform the test as directed in the Conductivity Measurement with the supernatant liquid obtained in the pH as the sample solution, and determine the conductivity at 25 ± 0.1 °C. Determine in the same manner the conductivity of water used for the preparation of the sample solution: the difference between these conductivities is not more than $75 \,\mu\text{S} \cdot \text{cm}^{-1}$.

Loss on drying $\langle 2.41 \rangle$ Not more than 7.0% \Diamond and within a range as specified on the label \Diamond (1 g, 105°C. 3 hours).

Residue on ignition $\langle 2.44 \rangle$ Not more than 0.1% (2 g).

Bulk density (i) Apparatus—Use a volumeter shown in the figure. Put a No.8.6 sieve $(2000 \, \mu\text{m})$ on the top of the volumeter. A funnel is mounted over a baffle box, having four glass baffle plates inside which the sample powder slides as it passes. At the bottom of the baffle box is a funnel that collect the powder, and allows it to pour into a sample receiving cup mounted directly below it.



(ii) Procedure—Weigh accurately the mass of a brass or stainless steel cup, which has a capacity of 25.0 ± 0.05 mL and an inside diameter of 30.0 ± 2.0 mm, and put the cup directly below the funnel of the volumeter. Slowly pour Microcrystalline Cellulose 5.1 cm height from the upper part of the powder funnel through the sieve, at a rate suitable to prevent clogging, until the cup overflows. If the clogging occurs, take out the sieve. Level the excess powder with the aid of a slide glass, weigh the filled cup, and weigh accurately the content of the cup, and then calculate the bulk

Table for Conversion of Relative Viscosity (η_{rel}) into the Product of Limiting Viscosity and Concentration ($[\eta]C$)

				$[\eta]C$								
η_{rel}	0.00	0.01	0.02	0.03	0.04	0.05	0.06	0.07	0.08	0.09		
1.1	0.098	0.106	0.115	0.125	0.134	0.143	0.152	0.161	0.170	0.180		
1.2	0.189	0.198	0.207	0.216	0.225	0.233	0.242	0.250	0.259	0.268		
1.3	0.276	0.285	0.293	0.302	0.310	0.318	0.326	0.334	0.342	0.350		
1.4	0.358	0.367	0.375	0.383	0.391	0.399	0.407	0.414	0.422	0.430		
1.5	0.437	0.445	0.453	0.460	0.468	0.476	0.484	0.491	0.499	0.507		
1.6	0.515	0.522	0.529	0.536	0.544	0.551	0.558	0.566	0.573	0.580		
1.7	0.587	0.595	0.602	0.608	0.615	0.622	0.629	0.636	0.642	0.649		
1.8	0.656	0.663	0.670	0.677	0.683	0.690	0.697	0.704	0.710	0.717		
1.9	0.723	0.730	0.736	0.743	0.749	0.756	0.762	0.769	0.775	0.782		
2.0	0.788	0.795	0.802	0.809	0.815	0.821	0.827	0.833	0.840	0.846		
2.1	0.852	0.858	0.864	0.870	0.876	0.882	0.888	0.894	0.900	0.906		
2.2	0.912	0.918	0.924	0.929	0.935	0.941	0.948	0.953	0.959	0.965		
2.3	0.971	0.976	0.983	0.988	0.994	1.000	1.006	1.011	1.017	1.022		
2.4	1.028	1.033	1.039	1.044	1.050	1.056	1.061	1.067	1.072	1.078		
2.5	1.083	1.089	1.094	1.100	1.105	1.111	1.116	1.121	1.126	1.131		
2.6	1.137	1.142	1.147	1.153	1.158	1.163	1.169	1.174	1.179	1.184		
2.7	1.190	1.195	1.200	1.205	1.210	1.215	1.220	1.225	1.230	1.235		
2.8	1.240	1.245	1.250	1.255	1.260	1.265	1.270	1.275	1.280	1.285		
2.9	1.290	1.295	1.300	1.305	1.310	1.314	1.319	1.324	1.329	1.333		
3.0	1.338	1.343	1.348	1.352	1.357	1.362	1.367	1.371	1.376	1.381		
3.1	1.386	1.390	1.395	1.400	1.405	1.409	1.414	1.418	1.423	1.427		
3.2	1.432	1.436	1.441	1.446	1.450	1.455	1.459	1.464	1.468	1.473		
3.3	1.477	1.482	1.486	1.491	1.496	1.500	1.504	1.508	1.513	1.517		
3.4	1.521	1.525	1.529	1.533	1.537	1.542	1.546	1.550	1.554	1.558		
3.5	1.562	1.566	1.570	1.575	1.579	1.583	1.587	1.591	1.595	1.600		
3.6	1.604	1.608	1.612	1.617	1.621	1.625	1.629	1.633	1.637	1.642		
3.7	1.646	1.650	1.654	1.658	1.662	1.666	1.671	1.675	1.679	1.683		
3.8	1.687	1.691	1.695	1.700	1.704	1.708	1.712	1.715	1.719	1.723		
3.9	1.727	1.731	1.735	1.739	1.742	1.746	1.750	1.754	1.758	1.762		
4.0	1.765	1.769	1.773	1.777	1.781	1.785	1.789	1.792	1.796	1.800		
4.1	1.804	1.808	1.811	1.815	1.819	1.822	1.826	1.830	1.833	1.837		
4.2	1.841	1.845	1.848	1.852	1.856	1.859	1.863	1.867	1.870	1.874		
4.3	1.878	1.882	1.885	1.889	1.893	1.896	1.900	1.904	1.907	1.911		
4.4	1.914	1.918	1.921	1.925	1.929	1.932	1.936	1.939	1.943	1.946		
4.5	1.950	1.954	1.957	1.961	1.964	1.968	1.971	1.975	1.979	1.982		
4.6	1.986	1.989	1.993	1.996	2.000	2.003	2.007	2.010	2.013	2.017		
4.7	2.020	2.023	2.027	2.030	2.033	2.037	2.040	2.043	2.047	2.050		
4.8	2.053	2.057	2.060	2.063	2.067	2.070	2.073	2.077	2.080	2.083		
4.9	2.087	2.090	2.093	2.097	2.100	2.103	2.107	2.110	2.113	2.116		
5.0	2.119	2.122	2.125	2.129	2.132	2.135	2.139	2.142	2.145	2.148		
5.1	2.151	2.154	2.158	2.160	2.164	2.167	2.170	2.173	2.176	2.180		
5.2	2.183	2.186	2.190	2.192	2.195	2.197	2.200	2.203	2.206	2.209		
5.3	2.212	2.215	2.218	2.221	2.224	2.227	2.230	2.233	2.236	2.240		
5.4	2.243	2.246	2.249	2.252	2.255	2.258	2.261	2.264	2.267	2.270		
5.5	2.273	2.276	2.279	2.282	2.285	2.288	2.291	2.294	2.297	2.300		
5.6	2.303	2.306	2.309	2.312	2.315	2.318	2.320	2.324	2.326	2.329		
5.7	2.332	2.335	2.338	2.341	2.344	2.347	2.350	2.353	2.355	2.358		
5.8	2.361	2.364	2.367	2.370	2.373	2.376	2.379	2.382	2.384	2.387		
5.9	2.390	2.393	2.396	2.400	2.403	2.405	2.408	2.411	2.414	2.417		
6.0	2.419	2.422	2.425	2.428	2.431	2.433	2.436	2.439	2.442	2.444		
6.1	2.447	2.450	2.453	2.456	2.458	2.461	2.464	2.467	2.470	2.472		
6.2	2.475	2.478	2.481	2.483	2.486	2.489	2.492	2.494	2.497	2.500		
6.3	2.503	2.505	2.508	2.511	2.513	2.516	2.518	2.521	2.524	2.526		
6.4	2.529	2.532	2.534	2.537	2.540	2.542	2.545	2.547	2.550	2.553		
6.5	2.555	2.558	2.561	2.563	2.566	2.568	2.571	2.574	2.576	2.579		
6.6	2.581	2.584	2.587	2.590	2.592	2.595	2.597	2.600	2.603	2.605		
6.7	2.608	2.610	2.613	2.615	2.618	2.620	2.623	2.625	2.627	2.630		
6.8 6.9	2.633	2.635	2.637	2.640	2.643	2.645	2.648	2.650	2.653	2.655		
6 ()	2.658	2.660	2.663	2.665	2.668	2.670	2.673	2.675	2.678	2.680		

	$[\eta]C$									
η_{rel}	0.00	0.01	0.02	0.03	0.04	0.05	0.06	0.07	0.08	0.09
7.0	2.683	2.685	2.687	2.690	2.693	2.695	2.698	2.700	2.702	2.705
7.1	2.707	2.710	2.712	2.714	2.717	2.719	2.721	2.724	2.726	2.729
7.2	2.731	2.733	2.736	2.738	2.740	2.743	2.745	2.748	2.750	2.752
7.3	2.755	2.757	2.760	2.762	2.764	2.767	2.769	2.771	2.774	2.776
7.4	2.779	2.781	2.783	2.786	2.788	2.790	2.793	2.795	2.798	2.800
7.5	2.802	2.805	2.807	2.809	2.812	2.814	2.816	2.819	2.821	2.823
7.6	2.826	2.828	2.830	2.833	2.835	2.837	2.840	2.842	2.844	2.847
7.7	2.849	2.851	2.854	2.856	2.858	2.860	2.863	2.865	2.868	2.870
7.8	2.873	2.875	2.877	2.879	2.881	2.884	2.887	2.889	2.891	2.893
7.9	2.895	2.898	2.900	2.902	2.905	2.907	2.909	2.911	2.913	2.915
8.0	2.918	2.920	2.922	2.924	2.926	2.928	2.931	2.933	2.935	2.937
8.1	2.939	2.942	2.944	2.946	2.948	2.950	2.952	2.955	2.957	2.959
8.2	2.961	2.963	2.966	2.968	2.970	2.972	2.974	2.976	2.979	2.981
8.3	2.983	2.985	2.987	2.990	2.992	2.994	2.996	2.998	3.000	3.002
8.4	3.004	3.006	3.008	3.010	3.012	3.015	3.017	3.019	3.021	3.023
8.5	3.025	3.027	3.029	3.031	3.033	3.035	3.037	3.040	3.042	3.044
8.6	3.046	3.048	3.050	3.052	3.054	3.056	3.058	3.060	3.062	3.064
8.7	3.067	3.069	3.071	3.073	3.075	3.077	3.079	3.081	3.083	3.085
8.8	3.087	3.089	3.092	3.094	3.096	3.098	3.100	3.102	3.104	3.106
8.9	3.108	3.110	3.112	3.114	3.116	3.118	3.120	3.122	3.124	3.126
9.0	3.128	3.130	3.132	3.134	3.136	3.138	3.140	3.142	3.144	3.146
9.1	3.148	3.150	3.152	3.154	3.156	3.158	3.160	3.162	3.164	3.166
9.2	3.168	3.170	3.172	3.174	3.176	3.178	3.180	3.182	3.184	3.186
9.3	3.188	3.190	3.192	3.194	3.196	3.198	3.200	3.202	3.204	3.206
9.4	3.208	3.210	3.212	3.214	3.215	3.217	3.219	3.221	3.223	3.225
9.5	3.227	3.229	3.231	3.233	3.235	3.237	3.239	3.241	3.242	3.244
9.6	3.246	3.248	3.250	3.252	3.254	3.256	3.258	3.260	3.262	3.264
9.7	3.266	3.268	3.269	3.271	3.273	3.275	3.277	3.279	3.281	3.283
9.8	3.285	3.287	3.289	3.291	3.293	3.295	3.297	3.298	3.300	3.302
9.9	3.304	3.305	3.307	3.309	3.311	3.313	3.316	3.318	3.320	3.321
•	0.0	0.1	0.2	0.3	0.4	0.5	0.6	0.7	0.8	0.9
10	3.32	3.34	3.36	3.37	3.39	3.41	3.43	3.45	3.46	3.48
11	3.50	3.52	3.53	3.55	3.56	3.58	3.60	3.61	3.63	3.64
12	3.66	3.68	3.69	3.71	3.72	3.74	3.76	3.77	3.79	3.80
13	3.80	3.83	3.85	3.86	3.88	3.89	3.90	3.92	3.93	3.95
14	3.96	3.97	3.99	4.00	4.02	4.03	4.04	4.06	4.07	4.09
15	4.10	4.11	4.13	4.14	4.15	4.17	4.18	4.19	4.20	4.22
16	4.23	4.24	4.25	4.27	4.28	4.29	4.30	4.31	4.33	4.34
17	4.35	4.36	4.37	4.38	4.39	4.41	4.42	4.43	4.44	4.45
18	4.46	4.47	4.48	4.49	4.50	4.52	4.53	4.54	4.55	4.56
19	4.57	4.58	4.59	4.60	4.61	4.62	4.63	4.64	4.65	4.66

density by the following expression: the bulk density is within the labeled specification.

Bulk density $(g/cm^3) = A/25$

A: Measured mass (g) of the content of the cup

Microbial limit <4.05> The acceptance criteria of TAMC and TYMC are 10³ CFU/g and 10² CFU/g, respectively. *Escherichia coli*, *Salmonella*, *Pseudomonas aeruginosa* and *Staphylococcus aureus* are not observed.

◆Containers and storage Containers—Tight containers. ◆

Chloramphenicol

クロラムフェニコール

Delete the Purity (2), move up (3) to (2) and change (2) as follows:

Purity

(2) Related substances—Dissolve 0.10 g of Chloramphenicol in 10 mL of methanol, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add methanol to make exactly 100 mL, and use this solution as the standard solution (1). Pipet 10 mL of the standard solution (1), add methanol to make exactly 20 mL, and use this

solution as the standard solution (2). Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 20 μ L each of the sample solution and standard solutions (1) and (2) on a plate of silica gel with fluorescent indicator for thin-layer chromatography, develop the plate with a mixture of ethyl acetate, methanol and acetic acid (100) (10:1:1) to a distance of about 15 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 254 nm): the spots other than the principal spot and the spot of the starting point obtained from the sample solution are not more intense than the spot from the standard solution (1), and the total amount of these spots from the sample solution is not more than 2.0%.

Change the Assay as follows:

Assay Weigh accurately an amount of Chloramphenicol and Chloramphenicol RS, equivalent to about 50 mg (potency), dissolve each in 10 mL of methanol, and add water to make exactly 50 mL. Pipet 20 mL each of these solutions, and add water to make exactly 100 mL. Pipet 10 mL each of these solutions, add water to make exactly 100 mL, and use these solutions as the sample solution and standard solution. Determine the absorbances, $A_{\rm T}$ and $A_{\rm S}$, of the sample solution and standard solution at 278 nm as directed under Ultraviolet-visible Spectrophotometry $\langle 2.24 \rangle$.

Amount [μ g (potency)] of chloramphenicol ($C_{11}H_{12}Cl_2N_2O_5$) = $M_S \times A_T/A_S \times 1000$

 M_S : Amount [mg (potency)] of Chloramphenicol RS taken

Chlorpromazine Hydrochloride

クロルプロマジン塩酸塩

Change the Melting point as follows:

Melting point <2.60> 196 – 200°C

Cholesterol

コレステロール

Change the Description and Optical rotation as follows:

Description Cholesterol occurs as white to pale yellow, crystals or grains. It is odorless, or has a slight odor. It is tasteless.

It is freely soluble in chloroform and in diethyl ether, sparingly soluble in ethanol (99.5) and in acetone, and practically insoluble in water.

It gradually changes to a yellow to light yellow-brown color by light.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁵: $-29 - -36^{\circ}$ (after drying, 0.2 g, acetone, 10 mL, 100 mm).

Add the following:

Clarithromycin for Syrup

シロップ用クラリスロマイシン

Clarithromycin for Syrup is a preparation for syrup, which is suspended before use.

It contains not less than 90.0% and not more than 110.0% of the labeled potency of clarithromycin ($C_{38}H_{60}NO_{13}$: 747.95).

Method of preparation Prepare as directed under Preparations for Syrups, with Clarithromycin.

Identification To an amount of Clarithromycin for Syrup, equivalent to 0.1 g (potency) of Clarithromycin, add 5 mL of acetone, and sonicate. After cooling with ice, centrifuge, take the supernatant liquid, and evaporate the solvent. Dissolve 10 mg of the residue and 2 mg of Clarithromycin RS in separate 2 mL of acetone, and use these solutions as the sample solution and the standard solution, respectively. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of methanol, ethyl acetate and acetic acid (100) (90:10:1) to a distance of about 15 cm, and air-dry the plate. Spray evenly sulfuric acid on the plate, and heat the plate at 105°C for 10 minutes: the principal spot obtained from the sample solution and the spot from the standard solution show a black-purple color and the same Rf value.

Water $\langle 2.48 \rangle$ Not more than 3.0% (0.5 g, volumetric titration, direct titration).

Uniformity of dosage units <6.02> Perform the test according to the following method: Clarithromycin for Syrup in single-dose packages meet the requirement of the Content uniformity test.

To the total content of 1 package of Clarithromycin for Syrup add 3V/5 mL of ethanol (99.5), add exactly V/10 mL of the internal standard solution, sonicate with occasional vigorous shaking, and add ethanol (99.5) to make V mL so that each mL contains about 0.5 mg (potency) of Clarithromycin. Centrifuge this solution, and filter the supernatant liquid through a membrane filter with a pore size not exceeding 0.45 μ m. Discard the first 3 mL of the filtrate, and use the subsequent filtrate as the sample solution. Then, proceed as directed in the Assay.

Amount [mg (potency)] of clarithromycin ($C_{38}H_{69}NO_{13}$) = $M_S \times Q_T/Q_S \times V/100$

 M_S : Amount [mg (potency)] of Clarithromycin RS taken

Internal standard solution—A solution of butyl parahydroxybenzoate in ethanol (99.5) (1 in 12,500).

Dissolution <6.10> When the test is performed at 50 revolutions per minute according to the Paddle method, using 900 mL of disodium hydrogen phosphate-citric acid buffer

solution (pH 5.5) as the dissolution medium, the dissolution rate in 90 minutes of Clarithromycin for Syrup is not less than 75%.

Start the test with an accurately weighed amount of Clarithromycin for Syrup, equivalent to about 50 mg (potency) of Clarithromycin, withdraw not less than 20 mL of the medium at the specified minute after starting the test, and filter through a membrane filter with a pore size not exceeding 0.45 μ m. Discard not less than 10 mL of the first filtrate, pipet 10 mL of the subsequent filtrate, add the mobile phase to make exactly 20 mL, and use this solution as the sample solution. Separately, weigh accurately an amount of Clarithromycin RS, equivalent to about 28 mg (potency), and dissolve in acetonitrile for liquid chromatography to make exactly 100 mL. Pipet 5 mL of this solution, add the mobile phase to make exactly 50 mL, and use this solution as the standard solution. Perform the test with exactly 100 µL each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, $A_{\rm T}$ and $A_{\rm S}$, of clarithromycin in each solution.

Dissolution rate (%) with respect to the labeled amount of clarithromycin ($C_{38}H_{69}NO_{13}$)

 $= M_{\rm S}/M_{\rm T} \times A_{\rm T}/A_{\rm S} \times 1/C \times 180$

 M_S : Amount [mg (potency)] of Clarithromycin RS taken M_T : Amount (g) of Clarithromycin for Syrup taken

C: Labeled amount of [mg (potency)] of clarithromycin (C₃₈H₆₉NO₁₃) in 1 g

Operating conditions—

Proceed as directed in the operating conditions in the Assay.

System suitability—

System performance: When the procedure is run with 100 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of clarithromycin are not less than 3000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with $100 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of clarithromycin is not more than 2.0%.

Assay Weigh accurately an amount of crushed Clarithromycin for Syrup, equivalent to about 50 mg (potency) of Clarithromycin, add 60 mL of ethanol (99.5), add exactly 10 mL of the internal standard solution, sonicate with occasional vigorous shaking, and add ethanol (99.5) to make 100 mL. Centrifuge this solution, and filter the supernatant liquid through a membrane filter with a pore size not exceeding 0.45 μ m. Discard the first 3 mL of the filtrate, and use the subsequent filtrate as the sample solution. Separately, weigh accurately an amount of Clarithromycin RS, equivalent about 50 mg (potency), and dissolve in ethanol (99.5) to make exactly 50 mL. Pipet 10 mL of this solution, add exactly 2 mL of the internal standard solution, add ethanol (99.5) to make 20 mL, and use this solution as the standard solution. Perform the test with 10 μ L each of the sample so-

lution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and calculate the ratios, $Q_{\rm T}$ and $Q_{\rm S}$, of the peak area of clarithromycin to that of the internal standard.

Amount [mg (potency)] of clarithromycin ($C_{38}H_{69}NO_{13}$) = $M_S \times Q_T/Q_S$

 $M_{\rm S}$: Amount [mg (potency)] of Clarithromycin RS taken

Internal standard solution—A solution of butyl parahydroxybenzoate in ethanol (99.5) (1 in 12,500).

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 210 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 50° C.

Mobile phase: A mixture of diluted 0.2 mol/L potassium dihydrogen phosphate TS (1 in 3) and acetonitrile for liquid chromatography (13:7).

Flow rate: Adjust so that the retention time of clarithromycin is about 8 minutes.

System suitability—

System performance: When the procedure is run with $10 \mu L$ of the standard solution under the above operating conditions, clarithromycin and the internal standard are eluted in this order with the resolution between these peaks being not less than 3.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of clarithromycin to that of the internal standard is not more than 2.0%.

Containers and storage Containers—Tight containers. Storage—Light-resistant.

Cloperastine Hydrochloride

クロペラスチン塩酸塩

Change the Melting point and Purity (2) as follows:

Melting point <2.60> 149 – 153°C

Purity

(2) Related substances—Dissolve 40 mg of Cloperastine Hydrochloride in 50 mL of the mobile phase, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add the mobile phase to make exactly 200 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine each peak area by the automatic integration method: the area of

the peaks, having the relative retention times of about 0.8 and about 3.0 to cloperastine obtained from the sample solution are not larger than the peak area of cloperastine from the standard solution, and the area of the peak having the relative retention time about 2.0 from the sample solution is not larger than 5/3 times the peak area of cloperastine from the standard solution. The area of the peak other than cloperastine and the peaks mentioned above from the sample solution are not larger than 3/5 times the peak area of cloperastine from the standard solution, and the total area of these peaks is not larger than 2 times the peak area of cloperastine from the standard solution.

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 222 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 25°C.

Mobile phase: A mixture of methanol, 0.1 mol/L potassium dihydrogen phosphate TS and perchloric acid (500:250:1).

Flow rate: Adjust so that the retention time of cloperastine is about 7 minutes.

Time span of measurement: About 4 times as long as the retention time of cloperastine, beginning after the solvent peak.

System suitability-

Test for required detectability: Pipet 2 mL of the standard solution, add the mobile phase to make exactly 20 mL. Confirm that the peak area of cloperastine obtained with 20 μ L of this solution is equivalent to 7 to 13% of that with 20 μ L of the standard solution.

System performance: Dissolve 30 mg of Cloperastine Hydrochloride and 40 mg of benzophenone in 100 mL of the mobile phase. To 2 mL of this solution add the mobile phase to make 50 mL. When the procedure is run with 20 μ L of this solution under the above operating conditions, cloperastine and benzophenone are eluted in this order with the resolution between these peaks being not less than 6.

System repeatability: When the test is repeated 6 times with $20 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of cloperastine is not more than 2.0%.

Dehydrocholic Acid

デヒドロコール酸

Change the Description and Optical rotation as follows:

Description Dehydrocholic Acid occurs as a white crystalline powder. It is odorless, and has a bitter taste.

It is sparingly soluble in acetone, slightly soluble in

ethanol (95), and practically insoluble in water.

It dissolves in sodium hydroxide TS.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁵: $+20 - +26^{\circ}$ (after drying, 0.2 g, acetone, 10 mL, 100 mm).

Purified Dehydrocholic Acid

精製デヒドロコール酸

Change the Description and Optical rotation as follows:

Description Purified Dehydrocholic Acid occurs as a white crystalline powder. It is odorless, and has a bitter taste.

It is sparingly soluble in acetone, slightly soluble in ethanol (95), and practically insoluble in water.

It dissolves in sodium hydroxide TS.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁵: $+20 - +26^{\circ}$ (after drying, 0.2 g, acetone, 10 mL, 100 mm).

Add the following:

Diclofenac Sodium Suppositories

ジクロフェナクナトリウム坐剤

Diclofenac Sodium Suppositories contain not less than 93.0% and not more than 107.0% of the labeled amount of diclofenac sodium ($C_{14}H_{10}Cl_2NNaO_2$: 318.13).

Method of preparation Prepare as directed under Suppositories for Rectal Application, with Diclofenac Sodium.

Identification To an amount of Diclofenac Sodium Suppositories, equivalent to 25 mg of Diclofenac Sodium, add 200 mL of a mixture of methanol and 0.01 mol/L sodium hydroxide TS (99:1), and dissolve by warming. Cool while shaking, add a mixture of methanol and 0.01 mol/L sodium hydroxide TS (99:1) to make 250 mL, and filter through a pledget of absorbent cotton if necessary. To 10 mL of this solution add a mixture of methanol and 0.01 mol/L sodium hydroxide TS (99:1) to make 100 mL. Determine the absorption spectrum of this solution as directed under Ultraviolet-visible Spectrophotometry <2.24>: it exhibits a maximum between 280 nm and 284 nm.

Uniformity of dosage unit <6.02> Perform the test according to the following method: it meets the requirement of the Content uniformity test.

To 1 suppository of Diclofenac Sodium Suppositories add 5 mL of tetrahydrofuran, and sonicate to dissolve. Add a mixture of methanol and water (3:2) to make exactly 100 mL, shake, and filter through a membrane filter with a pore size not exceeding $0.5 \,\mu\text{m}$. Discard the first 5 mL of the filtrate, pipet $V \, \text{mL}$ of the subsequent filtrate, add a mixture of methanol and water (3:2) to make exactly $V' \, \text{mL}$ so that

each mL contains about 0.125 mg of diclofenac sodium $(C_{14}H_{10}Cl_2NNaO_2)$, and use this solution as the sample solution. Then, proceed as directed in the Assay.

Amount (mg) of diclofenac sodium ($C_{14}H_{10}Cl_2NNaO_2$) = $M_S \times A_T/A_S \times V'/V \times 1/4$

 $M_{\rm S}$: Amount (mg) of diclofenac sodium for assay taken

Melting behavior of suppositories Perform the test according to Method 2 under Melting Point Determination <2.60>: the melting range is between 33°C and 36°C.

Assay Weigh accurately the mass of not less than 20 Diclofenac Sodium Suppositories, cut into small pieces carefully, and mix uniformly. Weigh accurately a portion of the pieces, equivalent to about 25 mg of diclofenac sodium (C₁₄H₁₀Cl₂NNaO₂), add 5 mL of tetrahydrofuran, and sonicate to dissolve. Add a mixture of methanol and water (3:2) to make exactly 100 mL, shake, and filter through a membrane filter with a pore size not exceeding $0.5 \mu m$. Discard the first 5 mL of the filtrate, pipet 10 mL of the subsequent filtrate, add a mixture of methanol and water (3:2) to make exactly 20 mL, and use this solution as the sample solution. Separately, weigh accurately about 50 mg of diclofenac sodium for assay, previously dried, and dissolve in a mixture of methanol and water (3:2) to make exactly 100 mL. Pipet 5 mL of this solution, add a mixture of methanol and water (3:2) to make exactly 20 mL, and use this solution as the standard solution. Perform the test with exactly 20 µL each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, A_T and A_S , of diclofenac in each solution.

Amount (mg) of diclofenac sodium ($C_{14}H_{10}Cl_2NNaO_2$) = $M_S \times A_T/A_S \times 1/2$

M_S: Amount (mg) of diclofenac sodium for assay taken

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 254 nm).

Column: A stainless steel column 4.0 mm in inside diameter and 12.5 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about $25\,^{\circ}\text{C}$.

Mobile phase: Dissolve 13.6 g of sodium acetate trihydrate in water to make 1000 mL. To 200 mL of this solution add 300 mL of methanol.

Flow rate: Adjust so that the retention time of diclofenac is about 3.5 minutes.

System suitability—

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of diclofenac are not less than 2000 and 0.7 to 1.5, respectively.

System repeatability: When the test is repeated 6 times

with $20 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of diclofenac is not more than 1.0%.

Containers and storage Containers—Tight containers. Storage—In a cold place.

Add the following:

Doripenem Hydrate

ドリペネム水和物

 $C_{15}H_{24}N_4O_6S_2.H_2O$: 438.52 (4R,5S,6S)-6-[(1R)-1-Hydroxyethyl]-4-methyl-7-oxo-3-{(3S,5S-5-[(sulfamoylamino)methyl]pyrrolidin-3-ylsulfanyl}-1- azabicyclo[3.2.0]hept-2-ene-2-carboxylic acid monohydrate [364622-82-2]

Doripenem Hydrate contains not less than 970 μg (potency) and not more than 1020 μg (potency) per mg, calculated on the anhydrous basis. The potency of Doripenem Hydrate is expressed as mass (potency) of doripenem ($C_{15}H_{24}N_4O_6S_2$: 420.50).

Description Doripenem Hydrate occurs as a white to pale yellow-brown-white crystalline powder.

It is sparingly soluble in water, slightly soluble in methanol, and practically insoluble in ethanol (99.5).

It is gradually colored to pale yellow-brown-white by light.

Identification (1) Determine the absorption spectrum of a solution of Doripenem Hydrate (1 in 50,000) as directed under Ultraviolet-visible Spectrophotometry $\langle 2.24 \rangle$, and compare the spectrum with the Reference Spectrum or the spectrum of a solution of Doripenem RS prepared in the same manner as the sample solution: both spectra exhibit similar intensities of absorption at the same wavelengths.

(2) Determine the infrared absorption spectrum of Doripenem Hydrate as directed in the potassium bromide disk method under Infrared Spectrophotometry <2.25>, and compare the spectrum with the Reference Spectrum or the spectrum of Doripenem RS: both spectra exhibit similar intensities of absorption at the same wave numbers.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁰: $+33 - +38^{\circ}$ (0.25 g calculated on the anhydrous basis, water, 25 mL, 100 mm).

pH <2.54> Dissolve 0.3 g of Doripenem Hydrate in 30 mL of water: the pH of the solution is between 4.5 and 6.0.

Purity (1) Clarity and color of solution—Dissolve 0.2 g of Doripenem Hydrate in 20 mL of water, and perform the test with this solution as directed under Turbidity Meas-

urement <2.61>: the solution is clear. Perform the test with this solution according to Method 2 under Methods for Color Matching <2.65>: the solution is not more colored than Matching Fluid Y4.

(2) Heavy metals <1.07>—Moisten 1.0 g of Doripenem Hydrate with sulfuric acid, cover loosely, and heat gently to carbonize. Then proceed according to Method 2, and perform the test. Prepare the control solution with 2.0 mL of Standard Lead Solution (not more than 20 ppm).

(3) Related substances (i)—Dissolve 20 mg of Doripenem Hydrate in 10 mL of water, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add water to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly $20 \,\mu\text{L}$ each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions. Determine each peak area by the automatic integration method: the peak areas of related substance A, having the relative retention time of about 2.2 to doripenem, related substance B, having the relative retention time of about 2.5, and related substance C, having the relative retention time of about 3.2, obtained from the sample solution, are not larger than 1/10 times the peak area of doripenem from the standard solution, and the area of the peak other than doripenem, the peaks mentioned above and the peak having the relative retention time of about 2.1, from the sample solution, is not larger than 1/20 times the peak area of doripenem from the standard solution. Furthermore, the total area of the peaks other than doripenem and the peak having the relative retention time of about 2.1 from the sample solution is not larger than 1/2 times the peak area of doripenem from the standard solution.

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 230 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 30°C.

Mobile phase A: Dissolve 2.04 g of potassium dihydrogen phosphate in water to make 1000 mL, and adjust to pH 5.6 – 5.7 with a solution prepared by dissolving 2.61 g of dipotassium hydrogen phosphate in water to make 1000 mL. To 970 mL of this solution add 30 mL of acetonitrile for liquid chromatography.

Mobile phase B: Dissolve 2.04 g of potassium dihydrogen phosphate in water to make 1000 mL, and adjust to pH 5.6 – 5.7 with a solution prepared by dissolving 2.61 g of dipotassium hydrogen phosphate in water to make 1000 mL. To 700 mL of this solution add 300 mL of acetonitrile for liquid chromatography.

Flowing of mobile phase: Control the gradient by mixing the mobile phases A and B as directed in the following table.

Time after injection of sample (min)	Mobile phase A (vol%)	Mobile phase B (vol%)
0 – 15	100	0
15 – 45	$100 \rightarrow 50$	$0 \rightarrow 50$
45 - 50	$50 \rightarrow 0$	$50 \rightarrow 100$
50 – 55	0	100

Flow rate: 1.0 mL per minute.

Time span of measurement: For 55 minutes after injection, beginning after the peak having the relative retention time of about 0.2 to doripenem.

System suitability—

Test for required detectability: Pipet 1.5 mL of the standard solution, and add water to make exactly 50 mL. Confirm that the peak area of doripenem obtained with 20 μ L of this solution is equivalent to 2.1 to 3.9% of that with 20 μ L of the standard solution.

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of doripenem are not less than 5000 and not more than 1.3, respectively.

System repeatability: When the test is repeated 3 times with $20 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of doripenem is not more than 0.95%.

(ii) Dissolve 20 mg of Doripenem Hydrate in 10 mL of water, and use this solution as the sample solution. Pipet 1 mL of the sample solution, and add water to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions. Determine each peak area by the automatic integration method: the peak area of related substance D, having the relative retention time of about 0.5 to doripenem, obtained from the sample solution is not larger than 2/5 times the peak area of doripenem from the standard solution.

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 215 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecyl-strong anion exchange-silanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about $40\,^{\circ}\mathrm{C}.$

Mobile phase: To 9 mL of phosphoric acid add 200 mL of water, add 20 mL of triethylamine, and add water to make 2000 mL. Adjust to pH 5.7 – 5.9 with phosphoric acid. To 950 mL of this solution add 50 mL of acetonitrile for liquid chromatography.

Flow rate: Adjust so that the retention time of doripenem is about 10 minutes.

System suitability—

Test for required detectability: Pipet 2 mL of the stand-

ard solution, and add water to make exactly 20 mL. Confirm that the peak area of doripenem obtained with 20 μ L of this solution is equivalent to 7 to 13% of that with 20 μ L of the standard solution.

System performance: To 1 mL of the sample solution add 1 mL of 0.1 mol/L hydrochloric acid TS, allow to stand at 25 \pm 5°C for 15 minutes, and add water to make 100 mL. When the procedure is run with 20 μ L of this solution under the above operating conditions, related substance D and doripenem are eluted in this order with the resolution between these peaks being not less than 5. The number of theoretical plates and the symmetry factor of the peak of related substance D are not less than 300 and 0.7 to 1.3, respectively, and those of the peak of doripenem are not less than 5000 and 0.7 to 1.3, respectively.

System repeatability: When the test is repeated 6 times with $20 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of doripenem is not more than 2.0%.

(iii) Dissolve 20 mg of Doripenem Hydrate in 10 mL of water, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add water to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions. Determine each peak area by the automatic integration method: the areas of the peaks, having the relative retention time of about 1.8, about 2.2 and about 2.3 to doripenem, obtained from the sample solution are not larger than 1/20, 7/100 and 1/20 times the peak area of doripenem from the standard solution, respectively.

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 310 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 25 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (3 μ m in particle diameter).

Column temperature: A constant temperature of about 30°C.

Mobile phase A: To 11 mL of perchloric acid add water to make 500 mL. To 100 mL of this solution add water to make 1000 mL. To 600 mL of this solution add 100 mL of water, and adjust to pH 1.9 – 2.0 with a solution prepared by adding water to 2.81 g of sodium perchlorate monohydrate to make 1000 mL. To 900 mL of this solution add 100 mL of acetonitrile.

Mobile phase B: To 11 mL of perchloric acid add water to make 500 mL. To 100 mL of this solution add water to make 1000 mL. To 600 mL of this solution add 100 mL of water, and adjust to pH 1.9-2.0 with a solution prepared by adding water to 2.81 g of sodium perchlorate monohydrate to make 1000 mL. To 300 mL of this solution add 200 mL of acetonitrile.

Flowing of mobile phase: Control the gradient by mixing the mobile phases A and B as directed in the following table.

Time after injection of sample (min)	Mobile phase A (vol%)	Mobile phase B (vol%)
0 – 25	100	0
25 - 55	$100 \rightarrow 0$	$0 \rightarrow 100$
55 - 60	0	100

Flow rate: 0.8 mL per minute.

System suitability—

Test for required detectability: Pipet 2.5 mL of the standard solution, and add water to make exactly 50 mL. Confirm that the peak area of doripenem obtained with 20 μ L of this solution is equivalent to 3.5 to 6.5% of that with 20 μ L of the standard solution.

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of doripenem are not less than 15,000 and not more than 1.3, respectively.

System repeatability: When the test is repeated 3 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of doripenem is not more than 0.95%.

Water $\langle 2.48 \rangle$ 4.0 – 5.0% (0.3 g, volumetric titration, back titration).

Residue on ignition $\langle 2.44 \rangle$ Not more than 0.1% (1 g).

Assay Weigh accurately amounts of Doripenem Hydrate and Doripenem RS (separately determine the water $\langle 2.48 \rangle$ in the same manner as Doripenem Hydrate), equivalent to about 25 mg (potency), dissolve each in water to make exactly 200 mL, and use these solutions as the sample solution and the standard solution, respectively. Perform the test with exactly 10 μ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and determine the peak areas, A_T and A_S , of doripenem in each solution.

Amount [μ g (potency)] of doripenem ($C_{15}H_{24}N_4O_6S_2$) = $M_S \times A_T/A_S \times 1000$

 M_S : Amount [mg (potency)] of Doripenem RS taken, calculated on the anhydrous basis

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 300 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter)

Column temperature: A constant temperature of about 25°C.

Mobile phase: Adjust the pH of 90 mL of 0.02 mol/L potassium dihydrogen phosphate TS to pH 5.6 – 5.7 with a solution prepared by dissolving 3.48 g of dipotassium hydrogen phosphate in water to make 1000 mL. To 100 mL of this solution add water to make exactly 1000 mL. To 970

mL of this solution add 30 mL of acetonitrile.

Flow rate: Adjust so that the retention time of doripenem is about 15 minutes.

System suitability—

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of doripenem are not less than 5000 and not more than 1.3, respectively.

System repeatability: When the test is repeated 6 times with $10 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of doripenem is not more than 1.0%.

Containers and storage Containers—Tight containers. Storage—At a temperature between 2°C and 8°C.

Others

Related substance A:

 $(4R,5S,6S)-3-\{(3S,5S)-5-[(\{N-(E)-(S,5S)-5\}-(\{N$

(Dimethylamino)methylene]sulfamoyl}amino)methyl] pyrrolidin-3-ylsulfanyl}-6-[(1*R*)-1-hydroxyethyl]-4-methyl-7-oxo-1-azabicyclo[3.2.0]hept-2-ene-2-carboxylic acid

Related substance B:

$$\label{eq:continuous} \begin{split} &(1S,4S,5S,6R)\text{-}4\text{-}[(1R)\text{-}1\text{-Hydroxyethyl}]\text{-}8\text{-}[(4R,5S,6S)\text{-}6\text{-}[(1R)\text{-}1\text{-hydroxyethyl}]\text{-}4\text{-methyl}\text{-}7\text{-}oxo\text{-}3\text{-}\{(3S,5S)\text{-}5\text{-}[(sulfamoylamino)methyl]pyrrolidin-3\text{-ylsulfanyl}}\text{-}1\text{-}azabicyclo[3.2.0]hept-2-ene-2-carbonyl]\text{-}6\text{-methyl}\text{-}3\text{-}oxo\text{-}7\text{-}\{(3S,5S)\text{-}5\text{-}[(sulfamoylamino)methyl]pyrrolidin-3\text{-ylsulfanyl}}\text{-}2\text{-}oxa\text{-}8\text{-}azabicyclo[3.2.1]octane-1\text{-}carboxylic acid} \end{split}$$

Related substance C: (4R,5S,6S)-3-[(3S,5S)-5- $(\{[N-(1,1-Dimethylethyl) sulfamoyl]amino\}methyl)pyrrolidin-3-ylsulfanyl]-6-[<math>(1R)$ -1-hydroxyethyl]-4-methyl-7-oxo-1-azabicyclo[3.2.0]hept-2-ene-2-carboxylic acid

Related substance D:

(2S,3R,4S)-2-[(1S,2R)-1-Carboxy-2-hydroxypropyl]-3-methyl-4-{(3S,5S)-5-[(sulfamoylamino)methyl]pyrrolidin-3-ylsulfanyl}-3,4-dihydro-2*H*-pyrrole-5-carboxylic acid

Add the following:

Doripenem for Injection

注射用ドリペネム

Doripenem for Injection is a preparation for injection, which is dissolved before use.

It contains not less than 95.0% and not more than 105.0% of the labeled potency of doripenem $(C_{15}H_{24}N_4O_6S_2: 420.50)$.

Method of preparation Prepare as directed under Injections, with Doripenem Hydrate.

Description Doripenem for Injection occurs as a white to pale yellow-brown-white crystalline powder.

Identification Proceed as directed in the Identification (2) under Doripenem Hydrate.

pH $\langle 2.54 \rangle$ Dissolve an amount of Doripenem for Injection, equivalent to 0.3 g (potency) of Doripenem Hydrate, in 30 mL of water: the pH of the solution is between 4.5 and 6.0.

Purity (1) Clarity and color of solution—Dissolve an amount of Doripenem for Injection, equivalent to 0.2 g (potency) of Doripenem Hydrate in 20 mL of water, and proceed as directed in the Purity (1) under Doripenem Hydrate.

(2) Related substances—(i) Dissolve an amount of Doripenem for Injection, equivalent to 20 mg (potency) of Doripenem Hydrate, in 10 mL of water, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add water to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly $20 \,\mu L$ each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine each peak area by the automatic integration method: the area of the peak other than doripenem and the peak, having the relative retention time of about 2.1 to doripenem, related substance A having the relative retention time of about 2.2, related substance B having the relative retention time of about 2.5 and related substance C having the relative retention time of about 3.2, obtained from the sample solution is not larger than 1/10 times the peak area of doripenem from the standard solution, and the total area of the peaks other than

doripenem and the peak mentioned above from the sample solution is not larger than 1/2 times the peak area of doripenem from the standard solution.

Operating conditions—

Proceed as directed in the operating conditions in the Purity (3) (i) under Doripenem Hydrate.

System suitability—

Test for required detectability: Pipet 2.5 mL of the standard solution, and add water to make exactly 50 mL. Confirm that the peak area of doripenem obtained with 20 μ L of this solution is equivalent to 3.5 to 6.5% of that with 20 μ L of the standard solution.

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of doripenem are not less than 5000 and not more than 1.3, respectively.

System repeatability: When the test is repeated 3 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of doripenem is not more than 0.95%.

(ii) Dissolve an amount of Doripenem for Injection, equivalent to 20 mg (potency) of Doripenem Hydrate, in 10 mL of water, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add water to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions. Determine each peak area by the automatic integration method: the peak area of related substance D, having the relative retention time of about 0.5 to doripenem, obtained from the sample solution is not larger than the peak area of doripenem from the standard solution.

Operating conditions—

Proceed as directed in the operating conditions in the Purity (3) (ii) under Doripenem Hydrate.

System suitability—

System performance: To 1 mL of the sample solution add 1 mL of 0.1 mol/L hydrochloric acid TS, allow to stand at 25 \pm 5°C for 15 minutes, and add water to make 100 mL. When the procedure is run with 20 μ L of this solution under the above operating conditions, related substance D and doripenem are eluted in this order with the resolution between these peaks being not less than 5. The number of theoretical plates and the symmetry factor of the peak of related substance D are not less than 300 and 0.7 to 1.3, respectively, and those of the peak of doripenem are not less than 5000 and 0.7 to 1.3, respectively.

System repeatability: When the test is repeated 6 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of doripenem is not more than 2.0%.

Water $\langle 2.48 \rangle$ 4.0 – 5.0% (0.3 g, volumetric titration, back titration)

Bacterial endotoxins <4.01> Less than 0.25 EU/mg (potency).

Uniformity of dosage units <6.02> It meets the requirement of the Mass variation test.

Foreign insoluble matter <6.06> Perform the test according to Method 2: it meets the requirement.

Insoluble particulate matter < 6.07> It meets the requirement.

Sterility < 4.06> Perform the test according to the Membrane filtration method: it meets the requirement.

Assay Weigh accurately the mass of the contents of not less than 10 containers of Doripenem for Injection. Weigh accurately an amount of the contents, equivalent to about 25 mg (potency) of Doripenem Hydrate, dissolve in water to make exactly 200 mL, and use this solution as the sample solution. Separately, weigh accurately about 25 mg (potency) of Doripenem RS (separately determine the water <2.48> in the same manner as Doripenem Hydrate), dissolve in water to make exactly 200 mL, and use this solution as the standard solution. Then, proceed as directed in the Assay under Doripenem Hydrate.

Amount [μ g (potency)] of doripenem (C₁₅H₂₄N₄O₆S₂) = $M_S \times A_T/A_S \times 1000$

 M_S : Amount [mg (potency)] of Doripenem RS taken, calculated on the anhydrous basis

Containers and storage Containers—Hermetic containers. Plastic containers for aqueous injections may be used.

Others

Related substances A, B, C and D: Refer to them described in Doripenem Hydrate.

Epirubicin Hydrochloride

エピルビシン塩酸塩

Delete the Purity (4):

Estriol

エストリオール

Change the Description and Optical rotation as follows:

Description Estriol occurs as a white crystalline powder. It is odorless.

It is sparingly soluble in methanol, slightly soluble in ethanol (99.5), and practically insoluble in water.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁵: +55 - +65° (after drying, 40 mg, ethanol (99.5), 10 mL, 100 mm).

Add the following:

Ethylcellulose

エチルセルロース

[9004-57-3]

This monograph is harmonized with the European Pharmacopoeia and the U. S. Pharmacopeia.

The corresponding part of the attributes/provisions which are agreed as non-harmonized within the scope of the harmonization is marked with symbols ($^{\blacklozenge}$ $_{\blacklozenge}$), and the corresponding parts which are agreed as the JP local requirement other than the scope of the harmonization are marked with symbols ($^{\diamondsuit}$ $_{\diamondsuit}$).

Ethylcellulose is a partly O-ethylated cellulose.

It contains not less than 44.0% and not more than 51.0% of ethoxy group (-OC₂H₅: 45.06), calculated on the dried basis.

It may contain a suitable antioxidant.

The viscosity of Ethylcellulose is shown in millipascal second (mPa·s) on the label.

◆Description Ethylcellulose occurs as a white to yellowish white, amorphous powder or grains.

It is soluble in dichloromethane.

It forms a slightly white-turbid or white-turbid, viscous liquid upon addition of ethanol (95).

To 1 g of Ethylcellulose add 100 mL of hot water, shake to become turbid, cool to room temperature, and add freshly boiled and cooled water to make 100 mL: the solution is neutral.

Identification Spread 2 drops of a solution of Ethylcellulose in dichloromethane (1 in 25) between sodium chloride plates, then remove one of the plates to evaporate the solvent, and determine the infrared absorption spectrum of the plate as directed in the film method under Infrared Spectrophotometry <2.25>, and compare the spectrum with the Reference Spectrum: both spectra exhibit similar intensities of absorption at the same wave numbers.

Viscosity <2.53> Weigh exactly a quantity of Ethylcellulose, equivalent to 5.00 g calculated on the dried basis, add 95 g of a mixture of 80 g of toluene and 20 g of ethanol (95), and shake to dissolve. Perform the test with this solution at 25°C as directed in Method I: not less than 80.0% and not more than 120.0% of the labeled viscosity for a nominal viscosity more than 6 mPa·s, and not less than 75.0% and not more than 140.0% of the labeled viscosity for a nominal viscosity not more than 6 mPa·s.

Purity (1) Acidity or alkalinity—To 0.5 g of Ethylcellulose add 25 mL of freshly boiled and cooled water, shake for 15 minutes, filter through a glass filter (G3), and use the filtrate as the sample solution. To 10 mL of the sample solution add 0.1 mL of dilute phenolphthalein TS and 0.5 mL of 0.01 mol/L sodium hydroxide VS: a light red color

develops. To 10 mL of the sample solution add 0.1 mL of methyl red-sodium hydroxide TS and 0.5 mL of 0.01 mol/L hydrochloric acid VS: a red color develops.

- (2) Chloride—Disperse 0.250 g of Ethylcellulose in 50 mL of water, and boil with occasional shaking. Allow to cool, and filter. Discard the first 10 mL of the filtrate, to 10 mL of the subsequent filtrate add water to make 15 mL, and use this solution as the sample solution. Separately, to 10 mL of Standard Chloride Solution add 5 mL of water, and use this solution as the control solution. To 15 mL each of the sample solution and control solution add 1 mL of 2 mol/L nitric acid TS, transfer to test tubes containing 1 mL of a solution of silver nitrate (17 in 1000), allow to stand for 5 minutes protecting from light, and compare the opalescence developed in the both solutions against a black background by viewing transversely: the opalescence developed in the sample solution is not more intense than that of the control solution (not more than 0.1%).
- ◇(3) Heavy metal <1.07>—Proceed with 1.0 g of Ethylcellulose according to Method 2, and perform the test. Prepare the control solution with 4.0 mL of Standard Lead Solution (not more than 40 ppm).
- (4) Acetaldehyde—Introduce 3.0 g of Ethylcellulose into a 250-mL glass-stoppered conical flask, add 10 mL of water, and stir for 1 hour. Allow to stand for 24 hours, filter, add water to the filtrate to make 100 mL, and use this solution as the sample solution. Separately, dissolve 1.0 g of acetaldehyde for assay in water to make 100 mL. To 5 mL of this solution add water to make 500 mL. To 3 mL of this solution add water to make 100 mL, and use this solution as the control solution. Transfer 5 mL each of the sample solution and control solution to 25-mL volumetric flasks, add 5 mL of a solution of 3-methyl-2-benzothiazolonehydrazone hydrochloride monohydrate (1 in 2000), and heat in a water bath at 60°C for 5 minutes. Add 2 mL of iron (III) chlorideamidosulfuric acid TS, and warm again at 60°C for 5 minutes. After cooling, add water to make 25 mL, and compare the color of these solutions: the sample solution is not more intensely colored than the control solution (not more than 100 ppm).

Loss on drying $\langle 2.41 \rangle$ Not more than 3.0% (1 g, 105°C, 2 hours).

Residue on Ignition $\langle 2.44 \rangle$ Not more than 0.5% (1 g).

Assay Weigh accurately about 30 mg of Ethylcellulose, transfer to a 5-mL pressure-tight serum vial, add exactly 60 mg of adipic acid, 2 mL of the internal standard solution and 1 mL of hydroiodic acid, seal the vial immediately with a septum coated with fluororesin and an aluminum cap or any other sealing system providing a sufficient air-tightness, and weigh accurately the vial. Take care not to mix the contents in the vial before heating. Place the vial in an oven or heat in a suitable heater with continuous stirring, maintaining an internal temperature of about $115 \pm 2^{\circ}$ C for 70 min. Allow to cool, and weigh accurately the vial. If the difference of the mass between before heating and after heating is more than 10 mg, prepare a new sample solution. If the

difference of the mass between before heating and after heating is not more than 10 mg, after phase separation, pierce through the septum of the vial with a cooled syringe, and withdraw a sufficient volume of the upper phase as the sample solution. Separately, place exactly 60 mg of adipic acid, 2 mL of the internal standard solution and 1 mL of hydroiodic acid in another serum vial, and seal immediately. Weigh accurately the vial, inject 25 μ L of iodoethane for assay through the septum in the vial, and weigh again accurately. Shake thoroughly, after phase separation, pierce through the septum of the vial with a cooled syringe, withdraw a sufficient volume of the upper phase, and use this solution as the standard solution. Perform the test with $1 \mu L$ each of the sample solution and standard solution as directed under Gas Chromatography <2.02> according to the following conditions, and calculate the ratios, Q_T and Q_S , of the peak area of iodoethane to that of the internal standard.

> Amount (%) of ethoxy group (C_2H_5O) = $M_S/M_T \times Q_T/Q_S \times 28.89$

 M_S : Amount (mg) of iodoethane for assay taken

 M_T : Amount (mg) of Ethylcellulose taken, calculated on the dried basis

Internal standard solution—A solution of *n*-octane in *o*-xylene (1 in 200).

Operating conditions—

Detector: A hydrogen flame-ionization detector.

Column: A fused silica column 0.53 mm in inside diameter and 30 m in length, coated with dimethylpolysiloxane for gas chromatography in 3 μ m thickness.

Column temperature: Maintain the temperature at 50°C for 3 minutes, raise to 100°C at a rate of 10°C per minute, then to 250°C at a rate of 35°C per minute, and maintain at 250°C for 8 minutes.

Injection port temperature: A constant temperature of about 250°C.

Detector temperature: A constant temperature of about $280^{\circ}\mathrm{C}$.

Carrier gas: Helium.

Flow rate: 4.2 mL per minute (the retention time of the internal standard is about 10 minutes).

Split ratio: 1:40.

System suitability—

System performance: When the procedure is run with 1 μ L of the standard solution under the above operating conditions, iodoethane and the internal standard are eluted in this order with the relative retention time of iodoethane to the internal standard being about 0.6, and with the resolution between these peaks being not less than 5.0.

System repeatability: When the test is repeated 6 times with $1 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of iodoethane to that of the internal standard is not more than 2.0%.

Containers and storage Containers—Well-closed containers.

Etizolam

エチゾラム

Change the Melting point as follows:

Melting point <2.60> 147 - 151°C

Add the following:

Felodipine

フェロジピン

C₁₈H₁₉Cl₂NO₄: 384.25

Ethyl methyl (4RS)-4-(2,3-dichlorophenyl)-2,6-dimethyl-1,4-dihydropyridine-3,5-dicarboxylate [72509-76-3]

Felodipine contains not less than 99.0% and not more than 101.0% of felodipine ($C_{18}H_{19}Cl_2NO_4$), calculated on the dried basis.

Description Felodipine occurs as pale yellow-white to light yellow-white, crystals or crystalline powder.

It is freely soluble in methanol and in ethanol (99.5), and practically insoluble in water.

A solution of Felodipine in methanol (1 in 20) shows no optical rotation.

Identification (1) Determine the absorption spectrum of a solution of Felodipine in methanol (1 in 62,500) as directed under Ultraviolet-visible Spectrophotometry $\langle 2.24 \rangle$, and compare the spectrum with the Reference Spectrum: both spectra exhibit similar intensities of absorption at the same wavelengths.

(2) Determine the infrared absorption spectrum of Felodipine as directed in the potassium bromide disk method under Infrared Spectrophotometry <2.25>, and compare the spectrum with the Reference Spectrum: both spectra exhibit similar intensities of absorption at the same wave numbers.

Purity (1) Heavy metals—Being specified separately when the drug is granted approval based on the Law.

(2) Related substances—Dissolve 25 mg of Felodipine in 50 mL of the mobile phase, and use this solution as the sample solution. Pipet 1 mL of the sample solution, and add the mobile phase to make exactly 100 mL. Pipet 1 mL of this solution, add the mobile phase to make exactly 10 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard

solution as directed under Liquid Chromatography <2.01> according to the following conditions. Determine each peak area by the automatic integration method: the area of the peak other than felodipine, related substance B, having the relative retention time of about 0.7 to felodipine, and related substance C, having the relative retention time of about 1.4, obtained from the sample solution is not larger than the peak area of felodipine from the standard solution. Furthermore, the total area of the peaks of related substances B and C from the sample solution is not larger than 10 times the peak area of felodipine from the standard solution, and the total area of the peaks other than felodipine and related substances mentioned above from the sample solution is not larger than 3 times the peak area of felodipine from the standard solution. For this calculation the peak area less than 1/5 times the peak area of felodipine from the standard solution is excluded.

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 254 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 25°C.

Mobile phase: Dissolve 3.2 g of sodium dihydrogen phosphate dihydrate in 400 mL of water, adjust to pH 3.0 with phosphoric acid, and add 200 mL of methanol and 400 mL of acetonitrile.

Flow rate: Adjust so that the retention time of felodipine is about 12 minutes.

Time span of measurement: About 2 times as long as the retention time of felodipine, beginning after the solvent peak.

System suitability—

Test for required detectability: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the SN ratio of the peak of felodipine is not less than 30.

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of felodipine are not less than 5000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with $20 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of felodipine is not more than 2.0%.

Loss on drying $\langle 2.41 \rangle$ Not more than 0.2% (1 g, 105°C, 3 hours).

Residue on ignition $\langle 2.44 \rangle$ Not more than 0.1% (1 g).

Assay Weigh accurately about 0.16 g of Felodipine, dissolve in 25 mL of t-butyl alcohol and 25 mL of diluted perchloric acid (17 in 200), and titrate $\langle 2.50 \rangle$ with 0.1 mol/L cerium (IV) sulfate VS (indicator: $50 \,\mu$ L of 1,10-

phenanthroline TS) until the color of the solution changes from orange to colorless. Perform a blank determination in the same manner, and make any necessary correction.

> Each mL of 0.1 mol/L cerium (IV) sulfate VS = 19.21 mg of $C_{18}H_{19}Cl_2NO_4$

Containers and storage Containers—Well-closed containers.

Others

Related substance B:

Dimethyl 4-(2,3-dichlorophenyl)-2,6-dimethyl-1,4-dihydropyridine-3,5-dicarboxylate

Related substance C: Diethyl 4-(2,3-dichlorophenyl)-2,6-dimethyl-1,4-dihydropyridine-3,5-dicarboxylate

Add the following:

Felodipine Tablets

フェロジピン錠

Felodipine Tablets contain not less than 95.0% and not more than 105.0% of the labeled amount of felodipine ($C_{18}H_{19}Cl_2NO_4$: 384.25).

Method of preparation Prepare as directed under Tablets, with Felodipine.

Identification To a quantity of powdered Felodipine Tablets, equivalent to 4 mg of Felodipine, add 200 mL of methanol, shake thoroughly, add methanol to make 250 mL, and centrifuge. Determine the absorption spectrum of the supernatant liquid as directed under Ultraviolet-visible Spectrophotometry <2.24>: it exhibits maxima between 235 nm and 239 nm, and between 357 nm and 363 nm.

Uniformity of dosage units <6.02> Perform the Mass variation test, or the Content uniformity test according to the following method: it meets the requirement.

To 1 tablet of Felodipine Tablets add 1 mL of water per 2.5 mg of felodipine ($C_{18}H_{19}Cl_2NO_4$), and shake thoroughly until the tablet is completely disintegrated. Add exactly 1 mL of the internal standard solution per 2.5 mg of felodipine ($C_{18}H_{19}Cl_2NO_4$), and add methanol to make V mL so

that each mL contains about 0.25 mg of felodipine $(C_{18}H_{19}Cl_2NO_4)$. Centrifuge this solution, filter the supernatant liquid, and use this filtrate as the sample solution. Separately, weigh accurately about 25 mg of felodipine for assay (separately determine the loss on drying $\langle 2.4I \rangle$ in the same conditions as Felodipine), add 10 mL of water, and add exactly 10 mL of the internal standard solution. Add methanol to make 100 mL, and use this solution as the standard solution. Then, proceed as directed in the Assay.

Amount (mg) of felodipine
$$(C_{18}H_{19}Cl_2NO_4)$$

= $M_S \times Q_T/Q_S \times V/100$

 M_S : Amount (mg) of felodipine for assay taken, calculated on the dried basis

Internal standard solution—A solution of butyl parahydroxybenzoate in methanol (1 in 3000).

Dissolution <6.10> When the test is performed at 50 revolutions per minute according to the Paddle method, using 900 mL of a solution prepared by dissolving 1 g of polysorbate 80 in water to make 5000 mL, as the dissolution medium, the dissolution rates in 45 minutes of a 2.5-mg tablet and a 5-mg tablet are not less than 80% and not less than 75%, respectively.

Start the test with 1 tablet of Felodipine Tablets, withdraw not less than 20 mL of the medium at the specified minute after starting the test, and filter through a membrane filter with a pore size not exceeding $0.45 \,\mu\text{m}$. Discard not less than 10 mL of the first filtrate, pipet V mL of the subsequent filtrate, add the dissolution medium to make exactly V' mL so that each mL contains about 2.8 μ g of felodipine (C₁₈H₁₉Cl₂NO₄), and use this solution as the sample solution. Separately, weigh accurately about 28 mg of felodipine for assay (separately determine the loss on drying <2.41> in the same conditions as Felodipine), and dissolve in methanol to make exactly 200 mL. Pipet 2 mL of this solution, add the dissolution medium to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 50 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, A_T and A_S , of felodipine in each solution.

Dissolution rate (%) with respect to the labeled amount of felodipine ($C_{18}H_{19}Cl_2NO_4$)

$$= M_{\rm S} \times A_{\rm T}/A_{\rm S} \times V'/V \times 1/C \times 9$$

 M_S : Amount (mg) of felodipine for assay taken, calculated on the dried basis

C: Labeled amount (mg) of felodipine ($C_{18}H_{19}Cl_2NO_4$) in 1 tablet

Operating conditions—

Column, column temperature, mobile phase and flow rate: Proceed as directed in the operating conditions in the Assay.

Detector: An ultraviolet absorption photometer (wavelength: 238 nm).

System suitability-

System performance: When the procedure is run with 50 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of felodipine are not less than 3000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with $50 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of felodipine is not more than 2.0%.

Assay Weigh accurately the mass of not less than 20 Felodipine Tablets, and powder. Weigh accurately a portion of the powder, equivalent to about 10 mg of felodipine (C₁₈H₁₉Cl₂NO₄), add 20 mL of water, add exactly 4 mL of the internal standard solution, add methanol to make 100 mL and shake. Centrifuge this solution, filter the supernatant liquid through a membrane filter with a pore size not exceeding 0.45 μ m, and use the filtrate as the sample solution. Separately, weigh accurately about 10 mg of felodipine for assay (separately determine the loss on drying <2.41> in the same conditions as Felodipine), add 20 mL of water, add exactly 4 mL of the internal standard solution, add methanol to make 100 mL, and use this solution as the standard solution. Perform the test with 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and calculate the ratios, Q_T and Q_S , of the peak area of felidipine to that of the internal standard.

Amount (mg) of felodipine $(C_{18}H_{19}Cl_2NO_4) = M_S \times Q_T/Q_S$

 $M_{\rm S}$: Amount (mg) of felodipine for assay taken, calculated on the dried basis

Internal standard solution—A solution of butyl parahydroxybenzoate in methanol (1 in 6000).

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 264 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 25°C.

Mobile phase: A mixture of methanol, water, a solution of sodium perchlorate monohydrate (281 in 2000) and diluted perchloric acid (17 in 200) (65:25:8:2).

Flow rate: Adjust so that the retention time of felodipine is about 12 minutes.

System suitability—

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the internal standard and felodipine are eluted in this order with the resolution between these peaks being not less than 5.

System repeatability: When the test is repeated 6 times with $20 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of felodipine is not more than 1.0%.

Containers and storage Containers—Tight containers.

Add the following:

Gatifloxacin Hydrate

ガチフロキサシン水和物

 $C_{19}H_{22}FN_3O_4.1\frac{1}{2}H_2O: 402.42$

1-Cyclopropyl-6-fluoro-8-methoxy-7-[(3RS)-3-

methylpiperazin-1-yl]-4-oxo-1,4-dihydroquinoline-3-carboxylic acid sesquihydrate

[180200-66-2]

Gatifloxacin Hydrate contains not less than 98.5% and not more than 101.5% of gatifloxacin ($C_{19}H_{22}FN_3O_4$: 375.39), calculated on the anhydrous basis.

Description Gatifloxacin Hydrate occurs as white to pale yellow, crystals or crystalline powder.

It is slightly soluble in methanol and in ethanol (99.5), and very slightly soluble in water.

It dissolves in sodium hydroxide TS.

It is gradually colored to pale yellow by light.

A solution of Gatifloxacin Hydrate in dilute sodium hydroxide TS (1 in 100) shows no optical rotation.

Identification (1) Determine the absorption spectrum of a solution of Gatifloxacin Hydrate in dilute sodium hydroxide TS (1 in 100,000) as directed under Ultraviolet-visible Spectrophotometry <2.24>, and compare the spectrum with the Reference Spectrum or the spectrum of a solution of Gatifloxacin RS prepared in the same manner as the sample solution: both spectra exhibit similar intensities of absorption at the same wavelengths.

(2) Determine the infrared absorption spectrum of Gatifloxacin Hydrate as directed in the potassium bromide disk method under Infrared Spectrophotometry <2.25>, and compare the spectrum with the Reference Spectrum or the spectrum of Gatifloxacin RS: both spectra exhibit similar intensities of absorption at the same wave numbers.

Purity (1) Clarity and color of solution—Dissolve 1.0 g of Gatifloxacin Hydrate in 10 mL of sodium hydroxide TS: the solution is clear. Perform the test with the solution as directed under Methods for Color Matching <2.65>: the solution is not more colored than diluted Matching Fluid O (1 in 5).

(2) Heavy metals <1.07>—Proceed with 1.0 g of Gatifloxacin Hydrate according to Method 4, and perform the test. Prepare the control solution with 2.0 mL of Standard Lead Solution (not more than 20 ppm).

(3) Related substances—Dissolve 20 mg of Gatifloxacin Hydrate in 50 mL of the dissolving solution, and use this solution as the sample solution. Pipet 1 mL of the sample solution, and add the dissolving solution to make exactly 100 mL. Pipet 2 mL of this solution, add the dissolving solution to make exactly 20 mL, and use this solution as the standard solution. Perform the test with exactly 20 µL each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions. Determine each peak area by the automatic integration method: the peak area of the related substance A, having the relative retention time of about 1.2 to gatifloxacin, obtained from the sample solution is not larger than 2 times the peak area of gatifloxacin from the standard solution, and the area of the peak other than gatifloxacin and the peak mentioned above from the sample solution is not larger than the peak area of gatifloxacin from the standard solution. Furthermore, the total area of the peaks other than gatifloxacin from the sample solution is not larger than 3 times the peak area of gatifloxacin from the standard solu-

Dissolving solution: A mixture of diluted phosphoric acid (1 in 1000) and acetonitrile (4:1).

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 325 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 35° C.

Mobile phase A: A mixture of diluted triethylamine (1 in 100), adjusted to pH 4.3 with phosphoric acid, and acetonitrile (22:3).

Mobile phase B: A mixture of diluted triethylamine (1 in 100), adjusted to pH 4.3 with phosphoric acid, and acetonitrile (1:1).

Flowing of mobile phase: Control the gradient by mixing the mobile phases A and B as directed in the following table.

Time after injection of sample (min)	Mobile phase A (vol%)	Mobile phase B (vol%)
0 – 15	100	0
15 – 30	$100 \rightarrow 0$	$0 \rightarrow 100$
30 - 40	0	100

Flow rate: 1.0 mL per minute (the retention time of gatifloxacin is about 16 minutes).

Time span of measurement: About 2.5 times as long as the retention time of gatifloxacin, beginning after the solvent peak.

System suitability—

Test for required detectability: Pipet 5 mL of the standard solution, and add the dissolving solution to make exactly 10 mL. Confirm that the peak area of gatifloxacin ob-

tained with 20 μ L of this solution is equivalent to 40 to 60% of that with 20 μ L of the standard solution.

System performance: Dissolve 20 mg of methyl 4-aminobenzoate in 50 mL of the dissolving solution. To 5 mL of this solution add 1 mL of the sample solution and the dissolving solution to make 100 mL. When the procedure is run with $20 \,\mu\text{L}$ of this solution under the above operating conditions, gatifloxacin and methyl 4-aminobenzoate are eluted in this order with the resolution between these peaks being not less than 4.

System repeatability: When the test is repeated 6 times with $20 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of gatifloxacin is not more than 3.0%.

Water $\langle 2.48 \rangle$ 6.0 – 9.0% (0.1 g, volumetric titration, direct titration).

Residue on ignition $\langle 2.44 \rangle$ Not more than 0.1% (1 g).

Assay Weigh accurately about 50 mg of Gatifloxacin Hydrate and Gatifloxacin RS (separately determine the water $\langle 2.48 \rangle$ in the same manner as Gatifloxacin Hydrate), and dissolve each in the dissolving solution to make exactly 100 mL. Pipet 2 mL each of these solution, add exactly 2 mL of the internal standard solution to them, add the dissolving solution to make 25 mL, and use these solutions as the sample solution and standard solutions. Perform the test with 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and calculate the ratios, Q_T and Q_S , of the peak area of gatifloxacin to that of the internal standard.

Amount (mg) of gatifloxacin (
$$C_{19}H_{22}FN_3O_4$$
)
= $M_S \times Q_T/Q_S$

 $M_{\rm S}$: Amount (mg) of Gatifloxacin RS taken, calculated on the anhydrous basis

Internal standard solution—A solution of methyl 4-aminobenzoate in the dissolving solution (1 in 4000).

Dissolving solution: A mixture of diluted phosphoric acid (1 in 1000) and acetonitrile (4:1).

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 280 nm).

Column: A stainless steel column 4 mm in inside diameter and 12.5 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 40°C.

Mobile phase: A mixture of diluted triethylamine (1 in 100), adjusted to pH 4.5 with phosphoric acid, and acetonitrile (87:13).

Flow rate: Adjust so that the retention time of gatifloxacin is about 5 minutes.

System suitability—

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, gatifloxacin and the internal standard are eluted in

this order with the resolution between these peaks being not less than 4.

System repeatability: When the test is repeated 6 times with $20 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of gatifloxacin to that of the internal standard is not more than 1.0%.

Containers and storage Containers—Tight containers. Storage—Light-resistant.

Others

Related substance A:

1-Cyclopropyl-6-fluoro-8-methoxy-7-[(2RS)-2-methylpiper-azin-1-yl]-4-oxo-1,4-dihydroquinoline-3-carboxylic acid

Add the following:

Gatifloxacin Ophthalmic Solution

ガチフロキサシン点眼液

Gatifloxacin Ophthalmic Solution is an aqueous ophthalmic preparation.

It contains not less than 95.0% and not more than 107.0% of the labeled amount of gatifloxacin ($C_{19}H_{22}FN_3O_4$: 375.39).

Method of preparation Prepare as directed under Ophthalmic Liquids and Solutions, with Gatifloxacin Hydrate.

Description Gatifloxacin Ophthalmic Solution is a clear, pale yellow liquid.

Identification To a volume of Gatifloxacin Ophthalmic Solution, equivalent to 6 mg of gatifloxacin ($C_{19}H_{22}FN_3O_4$), add diluted sodium hydroxide TS (1 in 10) to make 30 mL. To 1 mL of this solution add diluted sodium hydroxide TS (1 in 10) to make 20 mL, and determine the absorption spectrum of this solution as directed under Ultraviolet-visible Spectrophotometry $\langle 2.24 \rangle$: it exhibits maxima between 238 nm and 242 nm, between 287 nm and 291 nm, and between 336 nm and 340 nm.

Osmotic pressure ratio Being specified separately when the drug is granted approval based on the Law.

pH Being specified separately when the drug is granted approval based on the Law.

Purity Related substance—To a volume of Gatifloxacin Ophthalmic Solution, equivalent to 6 mg of gatifloxacin $(C_{19}H_{22}FN_3O_4)$, add the diluting solution to make 30 mL, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add the diluting solution to make exactly 100 mL. Pipet 2 mL of this solution, add the diluting

solution to make exactly 20 mL, and use this solution as the standard solution. Perform the test with exactly $40~\mu L$ each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions. Determine each peak area by the automatic integration method: the peak area of the related substance A, having the relative retention time of about 1.2 to gatifloxacin, obtained from the sample solution is not larger than 2 times the peak area of gatifloxacin from the standard solution, and the area of the peak other than gatifloxacin and the peak mentioned above from the sample solution is not larger than the peak area of gatifloxacin from the standard solution. Furthermore, the total area of the peaks other than gatifloxacin from the sample solution is not larger than 3 times the peak area of gatifloxacin from the standard solution.

Diluting solution: A mixture of diluted phosphoric acid (1 in 1000) and acetonitrile (4:1).

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 325 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 40°C.

Mobile phase A: A mixture of diluted triethylamine (1 in 100) and acetonitrile (22:3), adjusted to pH 4.3 with phosphoric acid.

Mobile phase B: A mixture of diluted triethylamine (1 in 100) and acetonitrile (1:1), adjusted to pH 4.3 with phosphoric acid.

Flowing of mobile phase: Control the gradient by mixing the mobile phases A and B as directed in the following table.

Time after injection of sample (min)	Mobile phase A (vol%)	Mobile phase B (vol%)
0 - 15	100	0
15 – 30	$100 \rightarrow 0$	$0 \rightarrow 100$
30 - 40	0	100

Flow rate: 0.9 mL per minute (the retention time of gatifloxacin is about 16 minutes).

Time span of measurement: For 40 minutes after injection, beginning after the solvent peak.

System suitability—

Test for required detectability: Pipet 5 mL of the standard solution, and add the diluting solution to make exactly 10 mL. Confirm that the peak area of gatifloxacin obtained with 40 μ L of this solution is equivalent to 40 to 60% of that with 40 μ L of the standard solution.

System performance: Dissolve 20 mg of methyl 4-aminobenzoate in 100 mL of the diluting solution. To 5 mL of this solution and 1 mL of the sample solution add the diluting solution to make 100 mL. When the procedure is

run with $40 \mu L$ of this solution under the above operating conditions, gatifloxacin and methyl 4-aminobenzoate are eluted in this order with the resolution between these peaks being not less than 4.

System repeatability: When the test is repeated 6 times with 40 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of gatifloxacin is not more than 3.0%.

Foreign insoluble matter <6.11> It meets the requirement.

Insoluble particulate matter <6.08> It meets the requirement

Sterility <4.06> Perform the test according to the Membrane filtration method: it meets the requirement.

Assay Pipet a volume of Gatifloxacin Ophthalmic Solution, equivalent to about 6 mg of gatifloxacin (C₁₉H₂₂FN₃O₄), and add the diluting solution to make exactly 30 mL. Pipet 2 mL of this solution, add exactly 3 mL of the internal standard solution, add the diluting solution to make 20 mL, and use this solution as the sample solution. Separately, weigh accurately about 22 mg of Gatifloxacin RS (separately determine the water <2.48> in the same manner as Gatifloxacin Hydrate), and dissolve in the diluting solution to make exactly 100 mL. Pipet 2 mL of this solution, add exactly 3 mL of the internal standard solution, add the diluting solution to make 20 mL, and use this solutions as the standard solution. Perform the test with 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and calculate the ratios, Q_T and $Q_{\rm S}$, of the peak area of gatifloxacin to that of the internal standard.

Amount (mg) of gatifloxacin (
$$C_{19}H_{22}FN_3O_4$$
)
= $M_S \times Q_T/Q_S \times 3/10$

 M_S : Amount (mg) of Gatifloxacin RS taken, calculated on the anhydrous basis

Internal standard solution—A solution of methyl 4-aminobenzoate in the diluting solution (1 in 10,000).

Diluting solution: A mixture of diluted phosphoric acid (1 in 1000) and acetonitrile (4:1).

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 280 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 40°C.

Mobile phase: A mixture of water, acetonitrile and triethylamine (81:18:1), adjusted to pH 4.5 with phosphoric acid.

Flow rate: Adjust so that the retention time of gatifloxacin is about 6 minutes.

System suitability—

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, gatifloxacin and the internal standard are eluted in this order with the resolution between these peaks being not less than 10.

System repeatability: When the test is repeated 6 times with $20 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of gatifloxacin to that of the internal standard is not more than 1.0%.

Containers and storage Containers—Tight containers.

Others

Related substances A: Refer to it described in Gatifloxacin Hydrate.

Add the following:

Gentamicin Sulfate Injection

ゲンタマイシン硫酸塩注射液

Gentamicin Sulfate Injection is an aqueous injection.

It contains not less than 90.0% and not more than 110.0% of the labeled potency of gentamicin C_1 ($C_{21}H_{43}N_5O_7$: 477.60).

Method of preparation Prepare as directed under Injections, with Gentamicin Sulfate.

Description Gentamicin Sulfate Injection is a clear and colorless liquid.

Identification To a volume of Gentamicin Sulfate Injection, equivalent to 40 mg (potency) of Gentamicin Sulfate, add water to make 10 mL, and use this solution as the sample solution. Separately, dissolve an amount of Gentamicin Sulfate RS, equivalent to 20 mg (potency), in 5 mL of water, and use this solution as the standard solution. Perform the test with these solutions as directed under Thinlayer Chromatography $\langle 2.03 \rangle$. Spot 5 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with the lower layer of a mixture of chloroform, ammonia solution (28) and methanol (2:1:1) to a distance of about 15 cm, and air-dry the plate. Spray evenly 0.2% ninhydrin-water saturated 1-butanol TS on the plate, and heat the plate at 100°C for 10 minutes: three principal spots obtained from the sample solution are the same with the corresponding spots from the standard solution in color tone and the Rf value, respectively.

Osmotic pressure ratio Being specified separately when the drug is granted approval based on the Law.

pH <2.54> 4.0 – 6.0

Bacterial endotoxins <4.01> Less than 0.50 EU/mg (po-

tency).

Extractable volume <6.05> It meets the requirement.

Foreign insoluble matter <6.06> Perform the test according to Method 1: it meets the requirement.

Insoluble particulate matter < 6.07> It meets the requirement.

Sterility <4.06> Perform the test according to the Membrane filtration method: it meets the requirement.

Assay Perform the test according to the Cylinder-plate method as directed under Microbial Assay for Antibiotics <4.02> according to the following conditions.

- (i) Test organism, agar media for base and seed layer, agar medium for transferring test organisms, and standard solutions—Proceed as directed in the Assay under Gentamicin Sulfate.
- (ii) Sample solutions—Pipet a volume of Gentamicin Sulfate Injection, equivalent to about 40 mg (potency) of Gentamicin Sulfate, add 0.1 mol/L phosphate buffer solution (pH 8.0) to make exactly 200 mL. Pipet a suitable volume of this solution, add 0.1 mol/L phosphate buffer solution (pH 8.0) to make solutions so that each mL contains 4 μ g (potency) and 1 μ g (potency), and use these solutions as the high concentration sample solution and the low concentration sample solution, respectively.

Containers and storage Containers—Hermetic containers.

Add the following:

Gentamicin Sulfate Ointment

ゲンタマイシン硫酸塩軟膏

Gentamicin Sulfate Ointment contains not less than 90.0% and not more than 110.0% of the labeled potency of gentamicin C_1 ($C_{21}H_{43}N_5O_7$: 477.60).

Method of preparation Prepare as directed under Ointments, with Gentamicin Sulfate.

Identification To an amount of Gentamicin Sulfate Ointment, equivalent to 5 mg (potency) of Gentamicin Sulfate, add 10 mL of diethyl ether, and shake in lukewarm water, if necessary, to dissolve. Add 5 mL of water, shake for 10 minutes, centrifuge, and use the water layer as the sample solution. Separately, dissolve an amount of Gentamicin Sulfate RS, equivalent to 10 mg (potency), in 10 mL of water, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with the lower layer of a mixture of chloroform, ammonia solution (28) and methanol (2:1:1) to a distance of about 15 cm, and air-dry the plate. Spray evenly 0.2% ninhydrin-water saturated 1butanol TS on the plate, and heat the plate at 100°C for 10 minutes: three principal spots obtained from the sample solution are the same with the corresponding spots from the standard solution in color tone and the Rf value, respectively.

Assay Perform the test according to the Cylinder-plate method as directed under Microbial Assay for Antibiotics <4.02> according to the following conditions.

- (i) Test organism, agar media for base and seed layer, agar medium for transferring test organisms, and standard solutions—Proceed as directed in the Assay under Gentamicin Sulfate.
- (ii) Sample solutions—Weigh accurately an amount of Gentamicin Sulfate Ointment, equivalent to about 1 mg (potency) of Gentamicin Sulfate, transfer to a separator, add 50 mL of diethyl ether, and shake until the solution becomes uniform. Add 25 mL of 0.1 mol/L phosphate buffer solution (pH 8.0), shake, and collect the water layer. Repeat the same procedure with 25 mL of 0.1 mol/L phosphate buffer solution (pH 8.0), and combine the water layers. To this solution add 0.1 mol/L phosphate buffer solution (pH 8.0) to make exactly 100 mL. Pipet a suitable volume of this solution, add 0.1 mol/L phosphate buffer solution (pH 8.0) to make solutions so that each mL contains 4 μ g (potency) and 1 μ g (potency), and use these solutions as the high concentration sample solution, respectively.

Containers and storage Containers—Tight containers.

Delete the following Monograph:

Adsorbed Habu-venom Toxoid

沈降はぶトキソイド

Haloperidol

ハロペリドール

Change the Melting point as follows:

Melting point <2.60> 150 - 154°C

Hydrocortisone

ヒドロコルチゾン

Change the Description and Optical rotation as follows:

Description Hydrocortisone occurs as a white crystalline powder.

It is sparingly soluble in methanol and in ethanol (99.5), and very slightly soluble in water.

Melting point: 212 – 220°C (with decomposition). It shows crystal polymorphism.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁵: +160 - +170° (after drying, 0.1 g, ethanol (99.5), 10 mL, 100 mm).

Hydrocortisone Acetate

ヒドロコルチゾン酢酸エステル

Change the Description and Optical rotation as follows:

Description Hydrocortisone Acetate occurs as white, crystals or crystalline powder.

It is freely soluble in dimethylsulfoxide, slightly soluble in methanol and in ethanol (95), and practically insoluble in water

Melting point: about 220°C (with decomposition). It shows crystal polymorphism.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁵: +154 - +164° (after drying, 50 mg, dimethylsulfoxide, 10 mL, 100 mm).

Hydrocortisone and Diphenhydramine Ointment

ヒドロコルチゾン・ジフェンヒドラミン軟膏

Change the Identification (3) as follows:

Identification

(3) To 0.2 g of Hydrocortisone and Diphenhydramine Ointment add 0.5 mL of methanol, warm, and shake. After cooling, separate the methanol layer, and use this layer as the sample solution. Dissolve 10 mg each of hydrocortisone acetate and diphenhydramine in 10 mL each of methanol, and use these solutions as standard solutions (1) and (2). Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 5 μ L each of the sample solution and standard solutions (1) and (2) on a plate of silica gel with fluorescent indicator for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate and diethyl ether (4:1) to a distance of about 5 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 254 nm): two spots obtained from the sample solution show the same Rf value as the corresponding spots from standard solutions (1) and (2).

Add the following:

Hydroxyethylcellulose

ヒドロキシエチルセルロース

[9004-62-0]

This monograph is harmonized with the European Pharmacopoeia and the U. S. Pharmacopeia.

The corresponding part of the attributes/provisions which are agreed as non-harmonized within the scope of the harmonization is marked with symbols ($^{\blacklozenge}$ $_{\spadesuit}$), and the corresponding parts which are agreed as the JP local requirement other than the scope of the harmonization are marked with symbols ($^{\diamondsuit}$ $_{\diamondsuit}$).

Hydroxyethylcellulose is partly O-(2-hydroxyethylated) cellulose.

It contains not less than 30.0% and not more than 70.0% of hydroxyethoxy group (-OC₂H₄OH: 61.06), calculated on the dried basis.

It may contain suitable pH-adjusting agents such as phosphates.

- ◆The viscosity is shown in millipascal second (mPa·s) on the label.
- ◆Description Hydroxyethylcellulose occurs as a white to yellowish white, powder or grains.

It is practically insoluble in ethanol (95).

It forms a viscous liquid upon addition of water.

It is hygroscopic.◆

- **Identification** (1) Determine the infrared absorption spectrum of Hydroxyethylcellulose as directed in the ATR method under Infrared Spectrophotometry <2.25> and compare the spectrum with the spectrum of Hydroxyethylcellulose RS for Identification: both spectra exhibit similar intensities of absorption at the same wave numbers.
- (2) Disperse 1.0 g of Hydroxyethylcellulose, calculated on the dried basis, in 50 mL of freshly boiled and cooled water. After 10 minutes, add freshly boiled and cooled water to make 100 mL, stir to dissolve completely, and use this solution as the sample solution. Boil 10 mL of the sample solution: the solution is clear.
- ♦ Viscosity <2.53> Weigh exactly a quantity of Hydroxyethylcellulose, equivalent to 10.00 g calculated on the dried basis, add 400 mL of water, stir to dissolve, and add water to make exactly 500.0 g. Remove air bubbles, and use this solution as the sample solution. Perform the test with the sample solution at 20 ± 0.1 °C as directed in Method II under Viscosity Determination, using a beaker with an inside diameter of not less than 70 mm and a single cylinder-type rotational viscometer, according to the following operating conditions: not less than 75% and not more than 140% of the labeled viscosity.

Operating conditions—

Apparatus: Brookfield type viscometer LV or RV model. Rotor No., rotation frequency and calculation multiplier: According to the following table, depending on the labeled viscosity.

L		Model	Rotor No.	Rotation frequency /min	Calculation multiplier		
		less	than 200	LV	1	30	2
Not less than 200 and less than 4000				LV	3	30	40
//	4000	//	10,000	LV	4	30	200
//	10,000	//	50,000	RV	6	20	500
//	50,000			RV	7	20	2000

Procedure of apparatus: Read a value after 2 minutes of rotation, and stop the rotation for at least 2 minutes. Repeat this procedure more two times, and average the three observed values.

pH $\langle 2.54 \rangle$ The pH of the sample solution obtained in the Identification (2) is between 5.5 and 8.5.

Purity (1) Chloride—To 1 mL of the sample solution obtained in the Identification (2) add water to make 30 mL, and use this solution as the sample solution. Separately, to 10 mL of Standard Chloride Solution add 5 mL of water, and use this solution as the control solution. To 15 mL each of the sample solution and control solution add 1 mL of diluted nitric acid (1 in 5), transfer to test tubes containing 1 mL of a solution of silver nitrate (17 in 1000), allow to stand for 5 minutes protecting from light, and compare the opalescence developed in the both solutions against a black background by viewing transversely: the opalescence developed in the sample solution is not more intense than that of the control solution (not more than 1.0%).

(2) Nitrate—Prepare the solutions before use. Dissolve 0.50 g of Hydroxyethylcellulose in the diluting solution to make exactly 100 mL, and use this solution as the sample solution. Separately, dissolve 0.8154 g of potassium nitrate in the diluting solution to make 1000 mL, and use this solution as the standard nitrate stock solution. If the viscosity of Hydroxyethylcellulose is not more than 1000 mPa·s, pipet 10 mL, 20 mL and 40 mL of the standard nitrate stock solution, add the diluting solution to each to make exactly 100 mL, and use these solutions as the standard solutions. If the viscosity of Hydroxyethylcellulose is more than 1000 mPa·s, pipet 1 mL, 2 mL and 4 mL of the standard nitrate stock solution, add the diluting solution to each to make exactly 100 mL, and use these solutions as the standard solutions. Perform the test with the sample solution and standard solutions using a nitrate-selective electrode as an indicator electrode, a silver-silver chloride electrode as a reference electrode and diluted ammonium sulfate TS (1 in 30) as reference electrolyte. Calculate the concentration of nitrates in the sample solution using a calibration curve obtained from the potential differences of the standard solutions: not more than 3.0%, calculated on the dried basis, if Hydroxyethylcellulose has a viscosity of not more than 1000 mPa·s, and not more than 0.2%, calculated on the dried basis, if Hydroxyethylcellulose has a viscosity of more than 1000 mPa·s.

Diluting solution: To a mixture of 50 mL of 1 mol/L sulfuric acid TS and 800 mL of water add 135 g of potassium

dihydrogen phosphate, and add water to make 1000 mL. To this solution add water to make exactly 25 times the initial volume.

In order to determine the applicable limit, determine the viscosity using the following procedure.

Introduce a quantity of Hydroxyethylcellulose, equivalent to 2.00 g calculated on the dried basis, into 50 g of water, stir, add water to make 100 g, and stir to dissolve completely. Determine the viscosity using a rotating viscometer at 25°C and at a shear rate of $100 \, \text{s}^{-1}$ for substances with an expected viscosity less than $100 \, \text{mPa·s}$, at a shear rate of $10 \, \text{s}^{-1}$ for substances with an expected viscosity not less than $100 \, \text{mPa·s}$, and not more than $20,000 \, \text{mPa·s}$, and at a shear rate of $1 \, \text{s}^{-1}$ for substances with an expected viscosity more than $20,000 \, \text{mPa·s}$. If it is impossible to obtain a shear rate of exactly $10 \, \text{s}^{-1}$ or $100 \, \text{s}^{-1}$ respectively, use a rate slightly higher and a rate slightly lower and interpolate.

- ◇(3) Heavy metal <1.07>—Proceed with 1.0 g of Hydroxyethylcellulose according to Method 2, and perform the test. Prepare the control solution with 2.0 mL of Standard Lead Solution (not more than 20 ppm).
- (4) Aldehydes—Introduce 1.0 g of Hydroxyethylcellulose into a glass-stoppered test tube, add 10 mL of ethanol (99.5), stopper the tube tightly, and stir for 30 minutes. Centrifuge, and use the supernatant liquid as the sample solution. Use Standard Glyoxal Solution as the control solution. Pipet 2 mL each of the sample solution and control solution, to each add 5 mL of a solution prepared by dissolving 4 g of 3-methyl-2-benzothiazolonehydrazone hydrochloride monohydrate in diluted acetic acid (100) (4 in 5) to make 1000 mL, shake to homogenize, and allow to stand for 2 hours. Compare the color of these solutions: the sample solution is not more intensely colored than the control solution (not more than 20 ppm).

Loss on drying $\langle 2.41 \rangle$ Not more than 10.0% (1 g, 105°C, 3 hours).

Residue on Ignition <2.44> Not more than 4.0% if the viscosity of Hydroxyethylcellulose is not more than 1000 mPa·s, and not more than 1.0% if the viscosity of Hydroxyethylcellulose is more than 1000 mPa·s (1 g). In order to determine the applicable limit, determine the viscosity according to the method in the Purity (2).

Assay Weigh accurately about 30 mg of Hydroxyethylcellulose, transfer to a 5-mL pressure-tight serum vial, add exactly 60 mg of adipic acid, 2 mL of the internal standard solution and 1 mL of hydroiodic acid, seal the vial immediately with a septum coated with fluororesin and an aluminum cap or any other sealing system providing a sufficient air-tightness, and weigh accurately the vial. Take care not to mix the content of the vial before heating. Place the vial in an oven or heat in a suitable heater with continuous stirring, maintaining the internal temperature of about $165 \pm 2^{\circ}$ C for 2.5 hours. Allow to cool and weigh accurately the vial. If the difference of the mass between before heating and after heating is more than 10 mg, prepare a new sample solution. If the difference of the mass between before heating

and after heating is not more than 10 mg, after phase separation, pierce through the septum of the vial with a cooled syringe, and withdraw a sufficient volume of the upper phase as the sample solution. Separately, place exactly 60 mg of adipic acid, 2 mL of the internal standard solution and 1 mL of hydroiodic acid in another serum vial, and seal immediately. Weigh accurately the vial, inject 55 μ L of iodoethane for assay through the septum in the vial, and weigh again accurately. Shake thoroughly, after phase separation, pierce through the septum of the vial with a cooled syringe, and withdraw a sufficient volume of the upper phase as the standard solution. Perform the test with $1 \mu L$ each of the sample solution and standard solution as directed under Gas Chromatography <2.02> according to the following conditions, and calculate the ratios, Q_T and Q_S , of the peak area of iodoethane to that of the internal stand-

Amount (%) of hydroethoxy group (
$$C_2H_5O_2$$
)
= $M_S/M_T \times Q_T/Q_S \times 39.15$

 $M_{\rm S}$: Amount (mg) of iodoethane for assay taken $M_{\rm T}$: Amount (mg) of Hydroxyethylcellulose taken, calculated on the dried basis

Internal standard solution—A solution of *n*-octane in *o*-xy-lene (1 in 200).

Operating conditions—

Detector: A hydrogen flame-ionization detector.

Column: A fused silica column 0.53 mm in inside diameter and 30 m in length, coated with dimethylpolysiloxane for gas chromatography in 3 μ m thickness.

Column temperature: Maintain the temperature at 50°C for 3 minutes, raise to 100°C at a rate of 10°C per minute, then raise to 250°C at a rate of 35°C per minute, and maintain at 250°C for 8 minutes.

Injection port temperature: A constant temperature of about 250°C.

Detector temperature: A constant temperature of about 280°C.

Carrier gas: Helium.

Flow rate: 4.2 mL per minutes (the retention time of the internal standard is about 10 minutes).

Split ratio: 1:40. System suitability—

System performance: When the procedure is run with 1 μ L of the standard solution under the above operating conditions, iodoethane and the internal standard are eluted in this order with the relative retention time of iodoethane to the internal standard being about 0.6 and the resolution between these peaks being not less than 5.0.

System repeatability: When the test is repeated 6 times with $1 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of iodoethane to that of the internal standard is not more than 2.0%.

◆Containers and storage Containers—Tight containers. ◆

Change the following as follows:

Hydroxypropylcellulose

ヒドロキシプロピルセルロース

[9004-64-2]

This monograph is harmonized with the European Pharmacopoeia and the U. S. Pharmacopeia.

The corresponding part of the attributes/provisions which are agreed as non-harmonized within the scope of the harmonization is marked with symbols ($^{\blacklozenge}$ $_{\spadesuit}$), and the corresponding parts which are agreed as the JP local requirement other than the scope of the harmonization are marked with symbols ($^{\diamondsuit}$ $_{\diamondsuit}$).

Hydroxypropylcellulose is partially *O*-(2-hydroxypropylated) cellulose.

It contains not less than 53.4% and not more than 80.5% of hydroxypropoxy group (-OC₃H₆OH: 75.09), calculated on the dried basis.

It may contain silicon dioxide as anti-caking agent.

- ◆The label states the addition in the case where silicon dioxide is added as anti-caking agent. ◆
- **◆Description** Hydroxypropylcellulose occurs as a white to yellowish white powder.

It forms a viscous liquid upon addition of water or ethanol (95).◆

Identification (1) Dissolve 1 g of Hydroxypropylcellulose in 100 mL of water, transfer 1 mL of the solution to a glass plate, and allow the water to evaporate: a thin film is formed.

- (2) Determine the infrared absorption spectrum of Hydroxypropylcellulose as directed in the potassium bromide disk method under Infrared Spectrophotometry <2.25>, and compare the spectrum with the Reference Spectrum: both spectra exhibit similar intensities of absorption at the same wave numbers. If there are an absorption at about 1719 cm⁻¹, disregard the absorption.
- **pH** <2.54> Disperse evenly 1.0 g of Hydroxypropylcellulose in 100 mL of freshly boiled water, and allow to cool the mixture while stirring with a magnetic stirrer: the pH of the solution is between 5.0 and 8.0.

Purity

- ♦(1) Heavy metals<1.07>—Proceed with 1.0 g of Hydroxypropylcellulose according to Method 2 and perform the test. Prepare the control solution with 2.0 mL of Standard Lead Solution (not more than 20 ppm).
- (2) Silicon dioxide—Apply to Hydroxypropylcellulose, if the addition of silicon dioxide is stated on the label and if more than 0.2% residue is found in the Residue on ignition test. Weigh accurately the crucible containing the residue tested in the Residue on ignition of Hydroxypropylcellulose (a (g)). Moisten the residue with water, and add 5 mL of hydrofluoric acid, in small portions. Evaporate it on a

steam bath to dryness and cool. Add 5 mL of hydrofluoric acid and 0.5 mL of sulfuric acid, and evaporate to dryness. Slowly increase the temperature until all the acids have been volatilized, and ignite at 1000 ± 25 °C. Cool the crucible in a desiccator, and weigh (b (g)). Calculate the amount of silicon dioxide by the following equation: not more than 0.6%.

Amount (%) of silicon dioxide (SiO₂)
=
$$(a - b)/M \times 100$$

M: Amount (g) of Hydroxypropylcellulose used for residue on ignition test

Loss on drying $\langle 2.41 \rangle$ Not more than 5.0% (1 g, 105°C, 4 hours).

Residue on ignition $\langle 2.44 \rangle$ Not more than 0.8% (1 g, platinum crucible).

Assay Weigh accurately about 30 mg of Hydroxypropylcellulose, transfer to a reaction vial, add exactly 60 mg of adipic acid, 2 mL of the internal standard solution and 1 mL of hydriodic acid, stopper the vial tightly, and weigh accurately. Place the vial in an oven or heat by a suitable heater with continuous stirring, maintaining the internal temperature of 115 ± 2 °C for 70 minutes. Allow the vial to cool and weigh accurately. If the difference of the mass between before heating and after heating is more than 10 mg, prepare a new test solution. If the difference of the mass between before heating and after heating is not more than 10 mg, after phase separation by allowing the vial to stand, pierce through the septum of the vial with a cooled syringe, and withdraw a sufficient volume of the upper phase as the sample solution. Separately, place exactly 60 mg of adipic acid, 2 mL of internal standard solution and 1 mL of hydriodic acid in an another reaction vial, stopper tightly, and weigh accurately. Inject 25 μ L of isopropyl iodide for assay through the septum, and again weigh accurately. Shake the vial thoroughly, and after phase separation by allowing the vial to stand, pierce through the septum of the vial with a cooled syringe, and withdraw a sufficient volume of the upper phase as the standard solution. Perform the test with 2 µL each of the sample solution and standard solution as directed under Gas Chromatography <2.02> according to the following conditions, and calculate the ratios, Q_T and Q_S , of the peak area of isopropyl iodide to that of the internal standard.

Amount (%) of hydroxypropoxy group ($C_3H_7O_2$) = $M_S/M_T \times Q_T/Q_S \times 1.15 \times 44.17$

 $M_{\rm S}$: Amount (mg) of isopropyl iodide for assay taken $M_{\rm T}$: Amount (mg) of Hydroxypropylcellulose taken, calculated on the dried basis

1.15: Correction factor

Internal standard solution—A solution of methylcyclohexan in *o*-xylene (1 in 50).

Operating conditions—

Detector: A hydrogen flame-ionization detector.

Column: A fused silica tube 0.53 mm in diameter and 30

m in length, coated with methylsilicone polymer for gas chromatography in 3 μ m thickness.

Column Temperature: Maintain the temperature at 40°C for 3 minutes, raise to 100°C at a rate of 10°C per minute, then raise to 250°C at a rate of 50°C per minute, and maintain at 250°C for 3 minutes.

Injection port temperature: A constant temperature of about 180°C.

Detector temperature: A constant temperature of about 280°C .

Carrier gas: Helium.

Flow rate: 52 cm per second (the retention time of the internal standard is about 8 minutes).

Split ratio: 1:50.

System suitability—

System performance: When the procedure is run with 2 μ L of the standard solution under the above operating conditions, isopropyl iodide and the internal standard are eluted in this order with the relative retention time of isopropyl iodide to the internal standard being about 0.8, and with the resolution between these peaks being not less than 2.0.

System repeatability: When the test is repeated 6 times with $2 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of isopropyl iodide to that of the internal standard is not more than 2.0%.

Containers and storage Containers—Well-closed containers.◆

Change the following as follows:

Hypromellose

ヒプロメロース

[9004-65-3]

This monograph is harmonized with the European Pharmacopoeia and the U.S. Pharmacopeia.

Hypromellose is a methyl and hydroxypropyl mixed ether of cellulose.

There are four substitution types of Hypromellose, 1828, 2208, 2906 and 2910. They contain methoxy (-OCH₃: 31.03) and hydroxypropoxy (-OC₃H₆OH: 75.09) groups conforming to the limits for the types of Hypromellose shown in the table below, calculated on the dried basis.

The viscosity is shown in millipascal second (mPa·s) on the label, together with the substitution

type.

Substitution Type -		y Group %)	, ,	Hydroxypropoxy Group (%)		
Type -	Min.	Max.	Min.	Max.		
1828	16.5	20.0	23.0	32.0		
2208	19.0	24.0	4.0	12.0		
2906	27.0	30.0	4.0	7.5		
2910	28.0	30.0	7.0	12.0		

◆Description Hypromellose occurs as a white to yellowish white, powder or grains.

It is practically insoluble in ethanol (99.5).

It swells with water and becomes a clear or slightly turbid, viscous solution. ♠

Identification (1) Disperse evenly 1.0 g of Hypromellose over the surface of 100 mL of water in a beaker, while gently tapping the top of the beaker, if necessary, and allow the beaker to stand: it aggregates on the surface of water.

- (2) Add 1.0 g of Hypromellose to 100 mL of hot water, and stir: it becomes a suspension. Cool the suspension to 10°C, and stir: the resulting liquid is a clear or a slightly cloudy, viscous fluid.
- (3) To 0.1 mL of the viscous fluid obtained in (2) add 9 mL of diluted sulfuric acid (9 in 10), shake, heat in a water bath for exactly 3 minutes, and immediately cool in ice water. Add carefully 0.6 mL of ninhydrin TS, shake, and allow to stand at 25°C: the solution shows a red color first, then changes to purple color within 100 minutes.
- (4) Pour and spread out 2 to 3 mL of the viscous fluid obtained in (2) onto a glass plate, and allow the water to evaporate: a transparent film results.
- (5) Pipet 50 mL of water, add exactly 50 mL of the viscous fluid obtained in (2), and warm to raise the temperature at a rate of 2 to 5°C per minute while stirring: the temperature, when a white turbidity of the solution starts to increase, is not less than 50°C.

Viscosity <2.53> (i) Method I: Apply to Hypromellose having a labeled viscosity of less than 600 mPa·s. Put an exact amount of Hypromellose, equivalent to 4.000 g calculated on the dried basis, in a tared, wide-mouth bottle, add water (between 90°C and 99°C) to make 200 g, stopper the bottle, stir by mechanical means at 350 to 450 revolutions per minute for 10 to 20 minutes to get a homogeneous dispersion. If necessary, take off the sample attached on the walls of the bottle, put them in the dispersed solution, and dissolve by continuing the stirring in a water bath at not exceeding 10°C for 20 to 40 minutes. Add cold water, if necessary, to make 200 g, and use this solution as the sample solution. Centrifuge the solution if necessary to expel any entrapped air bubbles. Perform the test with the sample solution at 20 ± 0.1 °C as directed in Method I under Viscosity Determination: not less than 80% and not more than 120% of the labeled viscosity.

(ii) Method II: Apply to Hypromellose having a labeled

viscosity of not less than 600 mPa·s. Put an exact amount of Hypromellose, equivalent to $10.00\,\mathrm{g}$ calculated on the dried basis, in a tared, wide-mouth bottle, add water (between 90°C and 99°C) to make 500 g, and prepare the sample solution in the same manner as directed in Method I. Perform the test with the sample solution at $20\pm0.1^\circ\mathrm{C}$ as directed in Method II under Viscosity Determination, using a single cylinder-type rotational viscometer, according to the following operating conditions: not less than 75% and not more than 140% of the labeled viscosity.

Operating conditions—

Apparatus: Brookfield type viscometer LV model or an equivalent apparatus.

Rotor No., rotation frequency, and calculation multiplier: According to the following table, depending on the labeled viscosity.

	Labeled (mP	viscos 'a·s)	ity		Rotor No.	Rotation frequency /min	Calculation multiplier
Not less	than 600	and le	ess	than 1400	3	60	20
//	1400	/	′/	3500	3	12	100
//	3500	/	′/	9500	4	60	100
//	9500	/	′/	99,500	4	6	1000
//	99,500				4	3	2000

Procedure of apparatus: Read the value after 2 minutes of rotation, and stop the rotation for at least 2 minutes. Repeat this procedure more two times, and average the three observed values.

pH $\langle 2.54 \rangle$ The pH of the sample solution obtained in the Viscosity, measured after 5 minutes immersing the electrode in the sample solution, is between 5.0 and 8.0.

Purity Heavy metals—Put 1.0 g of Hypromellose in a 100-mL Kjeldahl flask, add a sufficient amount of a mixture of nitric acid and sulfuric acid (5:4) to wet the sample, and heat gently. Repeat this procedure until to use totally 18 mL of the mixture of nitric acid and sulfuric acid (5:4). Then boil gently until the solution changes to black. After cooling, add 2 mL of nitric acid, and heat until the solution changes to black. Repeat this procedure until the solution no longer changes to black, and heat strongly until dense white fumes are evolved. After cooling, add 5 mL of water, boil gently until dense white fumes are evolved, then heat until the volume of the solution becomes to 2 to 3 mL. After cooling, if the solution reveals yellow color by addition of 5 mL of water, add 1 mL of hydrogen peroxide (30), and heat until the volume of the solution becomes to 2 to 3 mL. After cooling, dilute the solution with 2 to 3 mL of water, transfer to a Nessler tube, add water to make 25 mL, and use this solution as the sample solution. Separately, put 2.0 mL of Standard Lead Solution in a 100-mL Kjeldahl flask, add 18 mL of the mixture of nitric acid and sulfuric acid (5:4) and an amount of nitric acid equal to that used for preparation of the sample solution, and heat until white fumes are evolved. After cooling, add 10 mL of water. In the case where hydrogen peroxide (30) is added for the preparation of the sample solution, add the same amount of hydrogen peroxide (30), then proceed in the same manner for preparation of the sample solution, and use so obtained solution as the control solution. Adjust the sample solution and the control solution to pH 3.0 to 4.0 with ammonia solution (28), and add water to make 40 mL, respectively. To these solutions add 1.2 mL of thioacetamide-alkaline glycerin TS, 2 mL of acetate buffer solution (pH 3.5) and water to make 50 mL, separately. After allowing to stand for 5 minutes, observe vertically both tubes on a white background: the color obtained with the sample solution is not more intense than that with the control solution (not more than 20 ppm).

Loss on drying $\langle 2.41 \rangle$ Not more than 5.0% (1 g, 105°C, 1 hour).

Residue on ignition $\langle 2.44 \rangle$ Not more than 1.5% (1 g).

Assay (i) Apparatus—Reaction vial: A 5-mL pressuretight serum vial, having 20 mm in outside diameter and 50 mm in height, the neck 20 mm in outside diameter and 13 mm in inside diameter, equipped with a septum of butylrubber processed the surface with fluoroplastics, which can be fixed tightly to vial with aluminum cap, or equivalent.

Heater: A square-shaped aluminum block, having holes 20 mm in diameter and 32 mm in depth, adopted to the reaction vial. Capable of stirring the content of the reaction vial by means of magnetic stirrer or of reciprocal shaker about 100 times per minute.

(ii) Procedure-Weigh accurately about 65 mg of Hypromellose, transfer to the reaction vial, add 60 to 100 mg of adipic acid, 2.0 mL of the internal standard solution and 2.0 mL of hydroiodic acid, immediately stopper the vial tightly, and weigh accurately. Using a magnetic stirrer equipped in the heating module, or using a shaker, stir for 60 minutes while heating so that the temperature of the vial content is 130 ± 2 °C. In the case when a magnetic stirrer or shaker is not available, heat for 30 minutes with repeated shaking at 5-minute intervals by hand, and continue heating for an additional 30 minutes. Allow the vial to cool, and again weigh accurately. If the mass loss is not more than 26 mg and there is no evidence of a leak of the content, use the upper layer of the mixture as the sample solution. Separately, put 60 to 100 mg of adipic acid, 2.0 mL of the internal standard solution and 2.0 mL of hydroiodic acid in a reaction vial, immediately stopper the vial tightly, and weigh accurately. Add 45 μ L of iodomethane for assay and 15 to 22 μ L of isopropyl iodide for assay through the septum using a micro-syringe with weighing accurately every time, shake thoroughly, and use the upper layer of the content as the standard solution. Perform the test with 1 to 2 μ L each of the sample solution and standard solution as directed under Gas Chromatography <2.02> according to the following conditions, and calculate the ratios, Q_{Ta} and Q_{Tb} , of the peak area of iodomethane and isopropyl iodide to that of the internal standard obtained from the sample solution, and Q_{Sa} and Q_{Sb} , of the peak area of iodomethane and isopropyl iodide to that of the internal standard from the standard solution.

Content (%) of methoxy group (CH₃O) = $M_{Sa}/M \times Q_{Ta}/Q_{Sa} \times 21.86$

Content (%) of hydroxypropoxy group ($C_3H_7O_2$) = $M_{Sb}/M \times Q_{Tb}/Q_{Sb} \times 44.17$

 M_{Sa} : Amount (mg) of iodomethane for assay taken M_{Sb} : Amount (mg) of isopropyl iodide for assay taken M: Amount (mg) of Hypromellose taken, calculated on the dried basis

Internal standard solution—A solution of *n*-octane in *o*-xylene (3 in 100).

Operating conditions—

Detector: A thermal conductivity detector or hydrogen flame-ionization detector.

Column: A fused silica column 0.53 mm in inside diameter and 30 m in length, coated with dimethylpolysiloxane for gas chromatography in 3 μ m thickness. Use a guard column, if necessary.

Column temperature: Maintain the temperature at 50°C for 3 minutes, raise to 100°C at a rate of 10°C per minute, then raise to 250°C at a rate of 35°C per minute, and maintain at 250°C for 8 minutes.

Injection port temperature: A constant temperature of about 250°C.

Detector temperature: A constant temperature of about 280°C.

Carrier gas: Helium.

Flow rate: 4.3 mL per minute (the retention time of the internal standard is about 10 minutes).

Split ratio: 1:40. System suitability—

System performance: When the procedure is run with $1-2\,\mu L$ of the standard solution under the above operating conditions, iodomethane, isopropyl iodide and the internal standard are eluted in this order with the resolution between these peaks being not less than 5.

System repeatability: When the test is repeated 6 times with $1-2 \mu L$ of the standard solution under the above operating conditions, the relative standard deviations of the ratio of the peak area of iodomethane and isopropyl iodide to that of the internal standard are not more than 2.0%, respectively.

Containers and storage Containers—Well-closed containers.◆

Imipramine Hydrochloride

イミプラミン塩酸塩

Change the Melting point as follows:

Melting point <2.60> 172 - 176°C (with decomposition).

Imipramine Hydrochloride Tablets

イミプラミン塩酸塩錠

Delete the Identification (3):

Add the following:

Irinotecan Hydrochloride Hydrate

イリノテカン塩酸塩水和物

C₃₃H₃₈N₄O₆.HCl.3H₂O: 677.18

(4*S*)-4,11-Diethyl-4-hydroxy-3,14-dioxo-3,4,12,14-tetrahydro-1*H*-pyrano[3',4':6,7]indolizino[1,2-*b*]quinolin-9-yl [1,4'-bipiperidine]-1'-carboxylate monohydrochloride trihydrate [136572-09-3]

Irinotecan Hydrochloride Hydrate contains not less than 99.0% and not more than 102.0% of irinotecan hydrochloride ($C_{33}H_{38}N_4O_6$.HCl: 623.14), calculated on the anhydrous basis.

Description Irinotecan Hydrochloride Hydrate occurs as pale yellow to light yellow, crystals or crystalline powder.

It is sparingly soluble in methanol, and slightly soluble in water and in ethanol (99.5).

It is gradually colored to yellow-brown and decomposed by light.

Melting point: about 255°C (with decomposition). It shows crystal polymorphism.

Identification (1) Determine the absorption spectrum of a solution of Irinotecan Hydrochloride Hydrate in methanol (1 in 100,000) as directed under Ultraviolet-visible Spectrophotometry <2.24>, and compare the spectrum with the Reference Spectrum: both spectra exhibit similar intensities of absorption at the same wavelengths.

(2) Determine the infrared absorption spectrum of Irinotecan Hydrochloride Hydrate as directed in the potassium bromide disk method under Infrared Spectrophotometry <2.25>, and compare the spectrum with the Reference Spectrum: both spectra exhibit similar intensities of absorption at the same wave numbers.

(3) To 1 g of Irinotecan Hydrochloride Hydrate add 50 mL of water, dissolve by heating, and cool: the solution responds to Qualitative Tests <1.09> (2) for chloride.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁰: + 64 - + 69° (0.5 g calculated on the anhydrous basis, water, heat, after cooling, 50 mL, 100 mm).

pH <2.54> Dissolve 1 g of Irinotecan Hydrochloride Hy-

drate in 50 mL of water by heating, and cool: the pH of this solution is between 3.5 and 4.5.

Purity (1) Heavy metals <1.07>—Proceed with 2.0 g of Irinotecan Hydrochloride Hydrate according to Method 2, and perform the test. Prepare the control solution with 2.0 mL of Standard Lead Solution (not more than 10 ppm).

(2) Related substances—Dissolve 50 mg of Irinotecan Hydrochloride Hydrate in a suitable amount of a mixture of diluted 0.1 mol/L potassium dihydrogen phosphate TS (1 in 10), methanol and acetonitrile (6:4:3) and 1 mL of 1 mol/L hydrochloric acid TS, and add a mixture of diluted 0.1 mol/L potassium dihydrogen phosphate TS (1 in 10), methanol and acetonitrile (6:4:3) to make 20 mL, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add a mixture of diluted 0.1 mol/L potassium dihydrogen phosphate TS (1 in 10), methanol and acetonitrile (6:4:3) to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions. Determine each peak area by the automatic integration method: the peak areas of related substances A and B, having the relative retention time of about 0.8 to irinotecan, and related substances C and D, having the relative retention time of about 1.6, obtained from the sample solution are not larger than 1/5 times the peak area of irinotecan from the standard solution, and the area of the peak other than irinotecan and the peaks mentioned above from the sample solution is not larger than 1/10 times the peak area of irinotecan from the standard solution. Furthermore, the total area of the peaks other than irinotecan from the sample solution is not larger than 4/5 times the peak area of irinotecan from the standard solu-

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 254 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 25 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 40°C.

Mobile phase: Dissolve 1.22 g of sodium 1-decanesulfonate in a mixture of diluted 0.1 mol/L potassium dihydrogen phosphate TS (1 in 10), methanol and acetonitrile (6:4:3) to make 1000 mL.

Flow rate: Adjust so that the retention time of irinotecan is about 12 minutes.

Time span of measurement: About 3 times as long as the retention time of irinotecan.

System suitability—

Test for required detectability: Pipet 1 mL of the standard solution, and add a mixture of diluted 0.1 mol/L potassium dihydrogen phosphate TS (1 in 10), methanol and acetonitrile (6:4:3) to make exactly 20 mL. Confirm that the peak area of irinotecan obtained with $20 \,\mu\text{L}$ of this solution

is equivalent to 3.5 to 6.5% of that with 20 μ L of the standard solution.

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of irinotecan are not less than 6000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with $20 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of irinotecan is not more than 2.0%.

(3) Enantiomer—Being specified separately when the drug is granted approval based on the Law.

Water $\langle 2.48 \rangle$ 7.5 – 9.5% (0.1 g, volumetric titration, direct titration).

Residue on ignition $\langle 2.44 \rangle$ Not more than 0.1% (1 g).

Assay Weigh accurately about 0.44 g of Irinotecan Hydrochloride Hydrate, dissolve in 120 mL of a mixture of acetic anhydride and acetic acid (100) (7:3), and titrate <2.50> with 0.1 mol/L perchloric acid VS (potentiometric titration). Perform a blank determination in the same manner, and make any necessary correction.

Each mL of 0.1 mol/L perchloric acid VS = 31.16 mg of $C_{33}H_{38}N_4O_6$.HCl

Containers and storage Containers—Tight containers.
Storage—Light-resistant.

Others

Related substance A:

(4*S*)-4,11-Diethyl-4,12-dihydroxy-3,14-dioxo-3,4,12,14-tetrahydro-1*H*-pyrano[3',4':6,7]indolizino[1,2-*b*]quinolin-9-yl [1,4'-bipiperidine]-1'-carboxylate

Related substance B:

(4*S*)-4,11-Diethyl-4-hydroxy-3,14-dioxo-3,4,12,14-tetra-hydro-1*H*-pyrano[3',4':6,7]indolizino[1,2-*b*]quinolin-9-yl 2'-hydroxy-[1,4'-bipiperidine]-1'-carboxylate

Related substance C:

(4*S*)-4,8,11-Triethyl-4-hydroxy-3,14-dioxo-3,4,12,14-tetra-hydro-1*H*-pyrano[3',4':6,7]indolizino[1,2-*b*]quinolin-9-yl [1,4'-bipiperidine]-1'-carboxylate

Related substance D:

(4S)-4,11-Diethyl-4-hydroxy-3,14-dioxo-3,4,12,14-tetra-hydro-1H-pyrano[3',4':6,7]indolizino[1,2-b]quinolin-9-yl 2'-ethoxy-[1,4'-bipiperidine]-1'-carboxylate

Isomalt Hydrate

Isomalt

イソマル水和物

Change as follows:

6-O-α-D-Glucopyranosyl-D-glucitol $C_{12}H_{24}O_{11}$: 344.31 1-O-α-D-Glucopyranosyl-D-mannitol dihydrate $C_{12}H_{24}O_{11}$.2 H_2O : 380.34 6-O-α-D-Glucopyranosyl-D-glucitol—1-O-α-D-glucopyranosyl-D-mannitol dihydrate [64519-82-0]

This monograph is harmonized with the European Pharmacopoeia and the U.S. Pharmacopeia.

The corresponding part of the attributes/provisions which are agreed as non-harmonized within the scope of the harmonization is marked with symbols ($^{\blacklozenge}$ $_{\spadesuit}$), and the corresponding parts which are agreed as the JP local requirement other than the scope of the harmonization are marked with symbols ($^{\diamondsuit}$ $_{\diamondsuit}$).

Isomalt Hydrate is a mixture of 6-O- α -D-glucopyranosyl-D-sorbitol and 1-O- α -D-glucopyranosyl-D-mannitol.

It contains not less than 98.0% and not more than 102.0% as the mixture of 6-O- α -D-glucopyranosyl-D-sorbitol ($C_{12}H_{24}O_{11}$) and 1-O- α -D-glucopyranosyl-D-mannitol ($C_{12}H_{24}O_{11}$), calculated on the anhydrous basis, and the amount of each component is not less

than 3.0%, respectively.

The label states the contents (%) of $6-O-\alpha$ -D-glucopyranosyl-D-sorbitol and $1-O-\alpha$ -D-glucopyranosyl-D-mannitol.

◆Description Isomalt Hydrate occurs as a white, powder or grains.

It is freely soluble in water, and practically insoluble in ethanol (95).

Optical rotation $[\alpha]_D^{20}$: about + 92° (1 g calculated on the anhydrous basis, water, 100 mL, 100 mm).

Identification $^{\diamondsuit}$ (1) To 1 mL of a solution of Isomalt Hydrate (1 in 100) add 1 mL of a solution of catechol (1 in 10) prepared before use, shake thoroughly, add 2 mL of sulfuric acid rapidly, and shake: a reddish purple to red-purple color develops. $^{\diamondsuit}$

(2) Perform the test with $20 \,\mu\text{L}$ each of the sample solution and standard solution obtained in the Assay as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and compare the chromatograms obtained from these solutions: the two principal peaks in the chromatogram obtained from the sample solution are similar in retention time to respective two peaks from the standard solution.

Operating conditions—

Proceed as directed in the operating conditions in the Assay.

System suitability—

Proceed as directed in the system suitability in the Assay.

Purity [⋄](1) Heavy metals <1.07>—Proceed with 2.0 g of Isomalt Hydrate according to Method 1, and perform the test. Prepare the control solution with 2.0 mL of Standard Lead Solution (not more than 10 ppm)._⋄

(2) Nickel—Weigh exactly an amount of Isomalt Hydrate, equivalent to 10.0 g calculated on the anhydrous basis, dissolve in 30 mL of 2 mol/L acetic acid TS, and add water to make exactly 100 mL. Add exactly 2 mL of a solution of ammonium pyrrolidinedithiocarbamate (1 in100) and exactly 10 mL of water-saturated 4-methyl-2-pentanone, and shake for 30 seconds protected from light. Allow the layers to separate, and use the 4-methyl-2-pentanone layer as the sample solution. Separately, take in three vessels three exact portions of Isomalt Hydrate, each equivalent to 10.0 g calculated on the anhydrous basis, dissolve in 30 mL of 2 mol/L acetic acid TS, then add exactly 0.5 mL, 1.0 mL and 1.5 mL respectively of Standard Nickel Solution for Atomic Absorption Spectrophotometry, and add water to make them exactly 100 mL. Then, proceed in the same manner as the sample solution, and use the solutions so obtained as the standard solutions. Separately, prepare 4methyl-2-pentanone layer by proceeding in the same manner as the sample solution but omitting the substance to be examined, and use this solution as the blank solution. Perform the test with the sample solution and standard solution as directed in the standard addition method under Atomic Absorption Spectrophotometry <2.23> according to the following conditions. The blank solution is used to set the zero of the instrument, and to ascertain that the readings return to zero after rinsing the sample introduction system with water between each measurement: the amount of nickel is not more than 1 ppm.

Gas: Combustible gas—Acetylene.

Supporting gas—Air.

Lamp: Nickel hollow-cathode lamp.

Wavelength: 232.0 nm.

(3) Related substances—Weigh exactly 0.20 g of Isomalt Hydrate, dissolve in water to make exactly 10 mL, and use this solution as the sample solution. Separately, weigh exactly 10.0 mg of D-sorbitol and 10.0 mg of D-mannitol, dissolve in water to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine each peak area by the automatic integration method: the peak areas of D-mannitol, having a relative retention time of about 1.6 to 1-O- α -D-glucopyranosyl-D-mannitol, and D-sorbitol, having a relative retention time of about 2.0, obtained from the sample solution are not larger than the area of the corresponding peak from the standard solution (not more than 0.5%), and the area of the peak other of 1-O- α -D-glucopyranosyl-Dmannitol and $6-O-\alpha$ -D-glucopyranosyl-D-sorbitol having a relative retention time of about 1.2 and the peaks mentioned above from the sample solution is not larger than the peak area of D-sorbitol from the standard solution (not more than 0.5%). In addition, the total area of the peaks other than 1-O-α-D-glucopyranosyl-D-mannitol and 6-O-α-D-glucopyranosyl-D-sorbitol from the sample solution is not larger than 4 times the peak area of D-sorbitol from the standard solution (not more than 2.0%). However, the peaks which area is not larger than 1/5 times the peak area of D-sorbitol from the standard solution are disregarded (not more than 0.1%).

Operating conditions—

Detector, column, column temperature, mobile phase, and flow rate: Proceed as directed in the operating conditions in the Assay.

Time span of measurement: About 2.5 times as long as the retention time of 1-O- α -D-glucopyranosyl-D-mannitol. *System suitability*—

System performance: Proceed as directed in the system suitability in the Assay.

 $^{\circ}$ Test for required detectability: Pipet 2 mL of the standard solution, and add water to make exactly 10 mL. Confirm that the peak area of D-sorbitol obtained with 20 μ L of this solution is equivalent to 14 to 26% of that with 20 μ L the standard solution. $_{\circ}$

 $^{\circ}$ System repeatability: When the test is repeated 6 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviations of the peak area of D-mannitol and D-sorbitol are not more than 2.0%, respectively. $_{\circ}$

(4) Reducing sugars—Dissolve 3.3 g of Isomalt Hydrate in 10 mL of water with the aid of gentle heat, cool, and add 20 mL of copper (II) citrate TS. Add a few amount of boil-

ing chips, heat so that the boiling begins after 4 minutes, and maintain boiling for 3 minutes. Cool rapidly, add 100 mL of a solution of acetic acid (100) (3 in 125) and exactly 20 mL of 0.025 mol/L iodine VS. With continuous shaking, add 25 mL of a mixture of water and hydrochloric acid (47:3). When the precipitate has dissolved, titrate <2.50> the excess of iodine with 0.05 mol/L sodium thiosulfate VS, until the blue color due to 1 mL of soluble starch TS added at near of the end point disappears: not less than 12.8 mL of 0.05 mol/L sodium thiosulfate VS is required (not more than 0.3% as glucose).

Conductivity $\langle 2.51 \rangle$ Dissolve 20 g of Isomalt Hydrate in a suitable amount of freshly boiled and cooled water with the aid of gentle heat at 40 - 50°C, cool, add the same water to make exactly 100 mL, and use this solution as the sample solution. Measure the conductivity (25°C) of the sample solution at 25 ± 0.1 °C while gently stirring with a magnetic stirrer: not more than $20 \,\mu\text{S} \cdot \text{cm}^{-1}$.

Water $\langle 2.48 \rangle$ Not more than 7.0% (0.3 g, volumetric titration, direct titration. Use a mixture of methanol for water determination and formamide for water determination (1:1) heated at 50 \pm 5°C instead of methanol for water determination).

Assay Weight accurately about 0.2 g of Isomalt Hydrate, dissolve in water to make exactly 10 mL, and use this solution as the sample solution. Separately, weigh accurately about 0.2 g of Isomalt RS (separately determine the water $\langle 2.48 \rangle$ in the same manner as Isomalt Hydrate), dissolve in water to make exactly 10 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and determine the peak areas, $A_{\rm Ta}$ and $A_{\rm Tb}$, and $A_{\rm Sa}$ and $A_{\rm Sb}$, of 1-O- α -D-glucopyranosyl-D-mannitol and 6-O- α -D-glucopyranosyl-D-sorbitol in each solution.

Amount (g) of 1-O- α -D-glucopyranosyl-D-mannitol ($C_{12}H_{24}O_{11}$)

 $=M_{\rm S} \times K_{\rm a}/100 \times A_{\rm Ta}/A_{\rm Sa}$

Amount (g) of 6-O- α -D-glucopyranosyl-D-sorbitol ($C_{12}H_{24}O_{11}$)

 $= M_{\rm S} \times K_{\rm b}/100 \times A_{\rm Tb}/A_{\rm Sb}$

 $M_{\rm S}$: Amount (g) of Isomalt RS taken, calculated on the anhydrous basis

 K_a : Content (%) of 1-O- α -D-glucopyranosyl-D-mannitol ($C_{12}H_{24}O_{11}$) in Isomalt RS

 K_b : Content (%) of 6-O- α -D-glucopyranosyl-D-sorbitol ($C_{12}H_{24}O_{11}$) in Isomalt RS

Operating conditions—

Detector: A differential refractometer maintained at a constant temperature (40°C for example).

Column: Two stainless steel columns, 4.6 mm in inside diameter and 3 cm in length, and 7.8 mm in inside diameter and 30 cm in length, both packed with strongly acidic ion-

exchange resin (Ca type) for liquid chromatography with sulfonic acid groups attached to a polymer lattice consisting of polystyrene cross-linked with divinylbenzene (degree of cross-linkage: 8%) (9 μ m in particle diameter). These are used as the pre-column and the separation column, respectively.

Column temperature: 80 ± 3 °C.

Official Monographs

Mobile phase: Water.

Flow rate: 0.5 mL per minute (retention time of 1-O- α -D-glucopyranosyl-D-mannitol is about 12 minutes). System suitability—

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, 1-O- α -D-glucopyranosyl-D-mannitol and 6-O- α -D-glucopyranosyl-D-sorbitol are eluted in this order with the resolution between these peaks being not less than 2.0.

 $^{\circ}$ System repeatability: When the test is repeated 6 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of 1-O- α -D-glucopyranosyl-D-mannitol and 6-O- α -D-glucopyranosyl-D-sorbitol is not more than 2.0%, respectively. $_{\odot}$

◆Containers and storage Containers—Well-closed containers.

Add the following:

Lanoconazole

ラノコナゾール

 $C_{14}H_{10}ClN_3S_2$: 319.83 (2*E*)-2-[(4*RS*)-4-(2-Chlorophenyl)-1,3-dithiolan-2-ylidene]-2-(1*H*-imidazol-1-yl)acetonitrile [101530-10-3]

Lanoconazole, when dried, contains not less than 98.0% and not more than 102.0% of lanoconazole $(C_{14}H_{10}ClN_3S_2)$.

Description Lanoconazole occurs as white to pale yellow, crystals or crystalline powder.

It is soluble in acetone, sparingly soluble in methanol and in ethanol (99.5), and practically insoluble in water.

It is gradually colored to yellow by light.

A solution of Lanoconazole in acetone (1 in 25) shows no optical rotation.

Identification (1) To 0.1 g of Lanoconazole add 0.5 g of sodium hydroxide, heat gradually to melt, and carbonize. After cooling, add 10 mL of dilute hydrochloric acid: the gas evolved darkens moistened lead (II) acetate paper.

(2) Perform the test with Lanoconazole as directed

under Flame Coloration Test $\langle 1.04 \rangle$ (2): a green color appears.

- (3) Determine the absorption spectrum of a solution of Lanoconazole in methanol (1 in 100,000) as directed under Ultraviolet-visible Spectrophotometry <2.24>, and compare the spectrum with the Reference Spectrum or the spectrum of a solution of Lanoconazole RS prepared in the same manner as the sample solution: both spectra exhibit similar intensities of absorption at the same wavelengths.
- (4) Determine the infrared absorption spectrum of Lanoconazole, previously dried, as directed in the potassium bromide disk method under Infrared Spectrophotometry <2.25>, and compare the spectrum with the Reference Spectrum or the spectrum of previously dried Lanoconazole RS: both spectra exhibit similar intensities of absorption at the same wave numbers.

Melting point <2.60> 141 - 146°C

Purity (1) Heavy metals <1.07>—Proceed with 2.0 g of Lanoconazole according to Method 4, and perform the test. Prepare the control solution with 2.0 mL of Standard Lead Solution (not more than 10 ppm).

(2) Related substances—Conduct this procedure using light-resistant vessels. Dissolve 0.10 g of Lanoconazole in 100 mL of methanol, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add methanol to make exactly 50 mL, and use this solution as the standard solution. Perform the test with exactly $5\,\mu\text{L}$ each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions. Determine each peak area by the automatic integration method: the total area of the peaks other than lanoconazole obtained from the sample solution is not larger than 1/2 times the peak area of lanoconazole from the standard solution.

Operating conditions—

Detector, column, and column temperature: Proceed as directed in the operating conditions in the Assay.

Mobile phase: Dissolve 0.576 g of sodium 1-nonanesulfonate in 1000 mL of a mixture of methanol, water and acetic acid (100) (55:44:1).

Flow rate: Adjust so that the retention time of lanoconazole is about 7 minutes.

Time span of measurement: About 3 times as long as the retention time of lanoconazole, beginning after the solvent peak.

System suitability—

Test for required detectability: Pipet 2.5 mL of the standard solution, and add methanol to make exactly 50 mL. Confirm that the peak area of lanoconazole obtained with 5 μ L of this solution is equivalent to 3.5 to 6.5% of that with 5 μ L of the standard solution.

System performance: Put 20 mL of the sample solution in a colorless vessel, and expose to ultraviolet light (main wavelength: 365 nm) for 30 minutes. When the procedure is run with 5μ L of this solution under the above operating conditions, the resolution between the peak having the relative retention time of about 0.8 to lanoconazole and the

peak of lanoconazole is not less than 1.5.

System repeatability: When the test is repeated 6 times with $5 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of lanoconazole is not more than 1.0%.

Loss on drying $\langle 2.41 \rangle$ Not more than 0.4% (1 g, 105°C, 2 hours).

Residue on ignition $\langle 2.44 \rangle$ Not more than 0.1% (1 g).

Assay Conduct this procedure using light-resistant vessels. Weigh accurately about 50 mg each of Lanoconazole and Lanoconazole RS, both previously dried, and dissolve each in methanol to make exactly 50 mL. Pipet 5 mL each of these solutions, add exactly 5 mL of the internal standard solution, add methanol to make 50 mL, and use these solutions as the sample solution and the standard solution, respectively. Perform the test with 5 μ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and calculate the ratios, Q_T and Q_S , of the peak area of lanoconazole to that of the internal standard.

Amount (mg) of lanoconazole ($C_{14}H_{10}ClN_3S_2$) = $M_S \times Q_T/Q_S$

 M_S : Amount (mg) of Lanoconazole RS taken

Internal standard solution—A solution of diisopropyl 1,3-dithiolan-2-ylidenemalonate in methanol (1 in 1000).

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 295 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 50° C.

Mobile phase: A mixture of methanol and water (11:9).

Flow rate: Adjust so that the retention time of lanoconazole is about 9 minutes.

System suitability—

System performance: When the procedure is run with 5 μ L of the standard solution under the above operating conditions, lanoconazole and the internal standard are eluted in this order with the resolution between these peaks being not less than 3.

System repeatability: When the test is repeated 6 times with $5 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of lanoconazole to that of the internal standard is not more than 1.0%.

Containers and storage Containers—Well-closed containers

Storage—Light-resistant.

Add the following:

Lanoconazole Cream

ラノコナゾールクリーム

Lanoconazole Cream contains not less than 95.0% and not more than 105.0% of the labeled amount of lanoconazole ($C_{14}H_{10}ClN_3S_2$: 319.83).

Method of preparation Prepare as directed under Creams, with Lanoconazole.

Identification Warm Lanoconazole Cream to soften, if necessary. To a quantity of Lanoconazole Cream, equivalent to 50 mg of Lanoconazole, add 10 mL of diluted hydrochloric acid (1 in 6) saturated with sodium chloride, previously warmed, shake vigorously for 15 minutes to disperse, and centrifuge. Filter the supernatant liquid, wash the residue with 1.5 mL of diluted hydrochloric acid (1 in 6) saturated with sodium chloride, filter, and combine the washing with the filtrate. To the combined filtrate add 2.5 g of sodium hydrogen carbonate to dissolve, and extract with 10 mL of diethyl ether. Wash the diethyl ether layer with three 10-mL portions of water, and dry under reduced pressure. Dissolve the residue in 15 mL of acetone, and use this solution as the sample solution. Separately, dissolve 10 mg of lanoconazole in 10 mL of acetone, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography <2.03>. Spot 10 µL each of the sample solution and standard solution on a plate of silica gel with fluorescent indicator for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, toluene, methanol and ammonia solution (28) (400:400:20:1) to a distance of about 15 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 254 nm): the principal spot obtained from the sample solution and the spot from the standard solution show the same Rf value.

Assay Conduct this procedure using light-resistant vessels. Weigh accurately a quantity of Lanoconazole Cream, equivalent to about 15 mg of lanoconazole (C₁₄H₁₀ClN₃S₂), add 80 mL of methanol, sonicate to disperse, and add exactly 10 mL of the internal standard solution. Add methanol to make 100 mL, filter through a membrane filter with a pore size of $0.45 \,\mu\mathrm{m}$ if necessary, and use this solution as the sample solution. Separately, weigh accurately about 15 mg of Lanoconazole RS, previously dried at 105°C for 2 hours, dissolve in methanol, and add exactly 10 mL of the internal standard solution. Add methanol to make 100 mL, and use this solution as the standard solution. Perform the test with $10 \,\mu\text{L}$ each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and calculate the ratios, Q_T and Q_S , of the peak area of lanoconazole to that of the internal standard.

> Amount (mg) of lanoconazole ($C_{14}H_{10}ClN_3S_2$) = $M_S \times Q_T/Q_S$

2946

 $M_{\rm S}$: Amount (mg) of Lanoconazole RS taken

Internal standard solution—A solution of diisopropyl 1,3-dithiolan-2-ylidenemalonate in methanol (1 in 1000).

Operating conditions—

Proceed as directed in the operating conditions in the Assay under Lanoconazole.

System suitability—

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, lanoconazole and the internal standard are eluted in this order with the resolution between these peaks being not less than 3.

System repeatability: When the test is repeated 6 times with $10 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of lanoconazole to that of the internal standard is not more than 1.0%.

Containers and storage Containers—Tight containers. Storage—Light-resistant.

Add the following:

Lanoconazole Cutaneous Solution

ラノコナゾール外用液

Lanoconazole Cutaneous Solution is a liquid for external use.

It contains not less than 95.0% and not more than 105.0% of the labeled amount of lanoconazole $(C_{14}H_{10}CIN_3S_2: 319.83)$.

Method of preparation Prepare as directed under Liquids and Solutions for Cutaneous Application, with Lanocona-

Identification To a volume of Lanoconazole Cutaneous Solution, equivalent to 50 mg of Lanoconazole, add water enough to produce precipitate, and shake vigorously. Filter this solution, rinse the vessel with a suitable amount of water, and collect the precipitates. Wash the precipitates with 100 mL of water, dissolve in acetone, and dry under reduced pressure. If there are water droplets in the residue, dissolve the residue in 40 mL of acetone, and dry again under reduced pressure. Dissolve the residue in 30 mL of acetone, and use this solution as the sample solution. Separately, dissolve 10 mg of lanoconazole in 10 mL of acetone, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L each of the sample solution and standard solution on a plate of silica gel with fluorescent indicator for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, toluene, methanol and ammonia solution (28) (400:400:20:1) to a distance of about 15 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 254 nm): the principal spot obtained from the sample solution and the spot from the standard solution show the same Rf value.

Assay Conduct this procedure using light-resistant vessels. Pipet a volume of Lanoconazole Cutaneous Solution, equivalent to about 50 mg of lanoconazole (C₁₄H₁₀ClN₃S₂), and add methanol to make exactly 50 mL. Pipet 15 mL of this solution, add exactly 10 mL of the internal standard solution, add methanol to make 100 mL, and use this solution as the sample solution. Separately, weigh accurately about 15 mg of Lanoconazole RS, previously dried at 105°C for 2 hours, dissolve in methanol, and add exactly 10 mL of the internal standard solution. Add methanol to make 100 mL, and use this solution as the standard solution. Perform the test with $10 \,\mu$ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and calculate the ratios, Q_T and Q_S , of the peak area of lanoconazole to that of the internal standard.

> Amount (mg) of lanoconazole ($C_{14}H_{10}ClN_3S_2$) = $M_S \times Q_T/Q_S \times 10/3$

 $M_{\rm S}$: Amount (mg) of Lanoconazole RS taken

Internal standard solution—A solution of disopropyl 1,3-dithiolan-2-ylidenemalonate in methanol (1 in 1000).

Operating conditions—

Proceed as directed in the operating conditions in the Assay under Lanoconazole.

System suitability—

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, lanoconazole and the internal standard are eluted in this order with the resolution between these peaks being not less than 3.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of lanoconazole to that of the internal standard is not more than 1.0%.

Containers and storage Containers—Tight containers. Storage—Light-resistant.

Add the following:

Lanoconazole Ointment

ラノコナゾール軟膏

Lanoconazole Ointment contains not less than 93.0% and not more than 107.0% of the labeled amount of lanoconazole ($C_{14}H_{10}ClN_3S_2$: 319.83).

Method of preparation Prepare as directed under Ointments, with Lanoconazole.

Identification To a quantity of Lanoconazole Ointment, equivalent to 50 mg of Lanoconazole, add 15 mL of hexane, sonicate to disperse, add 10 mL of methanol, and shake for 10 minutes. Centrifuge this solution, discard the hexane

layer, and take the methanol layer. Wash the residue with a small amount of methanol if necessary, and combine the washing with the methanol layer. Dry the combined methanol layer under reduced pressure, dissolve the residue in 40 mL of acetone, and use this solution as the sample solution. Separately, dissolve 10 mg of lanoconazole in 10 mL of acetone, and use this solution as the standard solution. Perform the test with these solutions as directed under Thinlayer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L each of the sample solution and standard solution on a plate of silica gel with fluorescent indicator for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, toluene, methanol and ammonia solution (28) (400:400:20:1) to a distance of about 15 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 254 nm): the principal spot obtained from the sample solution and the spot from the standard solution show the same Rf value.

Assay Conduct this procedure using light-resistant vessels. Weigh accurately a quantity of Lanoconazole Ointment, equivalent to about 15 mg of lanoconazole (C₁₄H₁₀ClN₃S₂), add 20 mL of tetrahydrofuran, sonicate to disperse, and add exactly 10 mL of the internal standard solution. Add methanol to make 100 mL, and use this solution as the sample solution. Separately, weigh accurately about 15 mg of Lanoconazole RS, previously dried at 105°C for 2 hours, dissolve in methanol, and add exactly 10 mL of the internal standard solution. Add methanol to make 100 mL, and use this solution as the standard solution. Perform the test with $10 \,\mu L$ each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and calculate the ratios, Q_T and $Q_{\rm S}$, of the peak area of lanoconazole to that of the internal standard.

> Amount (mg) of lanoconazole $(C_{14}H_{10}ClN_3S_2)$ = $M_S \times Q_T/Q_S$

 M_S : Amount (mg) of Lanoconazole RS taken

Internal standard solution—A solution of disopropyl 1,3-dithiolan-2-ylidenemalonate in methanol (1 in 1000).

Operating conditions—

Proceed as directed in the operating conditions in the Assay under Lanoconazole.

System suitability—

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, lanoconazole and the internal standard are eluted in this order with the resolution between these peaks being not less than 3.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of lanoconazole to that of the internal standard is not more than 1.0%.

Containers and storage Containers—Tight containers. Storage—Light-resistant.

Mestranol

メストラノール

Change the Description and Optical rotation as follows:

Description Mestranol occurs as a white to pale yellowwhite crystalline powder. It is odorless.

It is freely soluble in chloroform, sparingly soluble in ethanol (99.5), and practically insoluble in water.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁵: +1 - +6° (after drying, 0.1 g, ethanol (99.5), 10 mL, 100 mm).

Change the following as follows:

Methylcellulose

メチルセルロース

[9004-67-5]

This monograph is harmonized with the European Pharmacopoeia and the U.S. Pharmacopeia.

The corresponding part of the attributes/provisions which are agreed as non-harmonized within the scope of the harmonization is marked with symbols ($^{\blacklozenge}$ $_{\blacklozenge}$), and the corresponding parts which are agreed as the JP local requirement other than the scope of the harmonization are marked with symbols ($^{\diamondsuit}$ $_{\diamondsuit}$).

Methylcellulose is a methyl ether of cellulose.

It contains not less than 26.0% and not more than 33.0% of methoxy group (-OCH₃: 31.03), calculated on the dried basis.

The viscosity of Methylcellulose is shown in millipascal second (mPa·s).

◆**Description** Methylcellulose occurs as a white to yellowish white, powder or grains.

It is practically insoluble in ethanol (99.5).

It swells, when water is added, and forms a clear or slightly turbid, viscous liquid.◆

Identification (1) Disperse evenly 1.0 g of Methylcellulose over the surface of 100 mL of water in a beaker, while gently tapping the top of the beaker, if necessary, and allow the beaker to stand: it aggregates on the surface of water.

- (2) Add 1.0 g of Methylcellulose to 100 mL of hot water, and stir: it becomes a suspension. Cool the suspension to 5°C, and stir: the resulting liquid is a clear or a slightly cloudy, viscous fluid.
- (3) To 0.1 mL of the viscous fluid obtained in (2) add 9 mL of diluted sulfuric acid (9 in 10), shake, heat in a water bath for exactly 3 minutes, and immediately cool in ice water. Add carefully 0.6 mL of ninhydrin TS, shake, and allow to stand at 25°C: the solution shows a red color, and

it does not change to purple color within 100 minutes.

- (4) Pour and spread out 2 to 3 mL of the viscous fluid obtained in (2) onto a glass plate, and allow the water to evaporate: a transparent film results.
- (5) Pipet 50 mL of water, add exactly 50 mL of the viscous fluid obtained in (2), and warm to raise the temperature at a rate of 2 to 5°C per minute while stirring: the temperature, when a white turbidity of the solution starts to increase, is not less than 50°C.

Viscosity <2.53> (i) Method I: Apply to Methylcellulose having a labeled viscosity of less than 600 mPa·s. Put an exact amount of Methylcellulose, equivalent to 4.000 g, calculated on the dried basis, in a tared, wide-mouth bottle, add water (between 90°C and 99°C) to make 200 g, stopper the bottle, stir by mechanical means at 350 to 450 revolutions per minute for 10 to 20 minutes to get a homogeneous dispersion. If necessary, take off the sample attached on the walls of the bottle, put them in the dispersed solution, and dissolve by continuing the stirring in a water bath at not exceeding 5°C for 20 to 40 minutes. Add cold water, if necessary, to make 200 g, and use this solution as the sample solution. Centrifuge the solution if necessary to expel any entrapped air bubbles. Perform the test with the sample solution at 20 ± 0.1 °C as directed in Method I under Viscosity Determination: not less than 80% and not more than 120% of the labeled viscosity.

(ii) Method II: Apply to Methylcellulose having a labeled viscosity of not less than 600 mPa·s. Put an exact amount of Methylcellulose, equivalent to $10.00\,\mathrm{g}$, calculated on the dried basis, in a tared, wide-mouth bottle, add water (between $90^{\circ}\mathrm{C}$ and $99^{\circ}\mathrm{C}$) to make $500\,\mathrm{g}$, and prepare the sample solution in the same manner as directed in Method I. Perform the test with the sample solution at $20 \pm 0.1^{\circ}\mathrm{C}$ as directed in Method II under Viscosity Determination, using a single cylinder-type rotational viscometer, according to the following operating conditions: not less than 75% and not more than 140% of the labeled viscosity. Operating conditions—

Apparatus: Brookfield type viscometer LV model or an equivalent apparatus.

Rotor No., rotation frequency, and calculation multiplier: According to the following table, depending on the labeled viscosity.

	Labeled (mP	viscosi 'a·s)	ty		Rotor No.	Rotation frequency /min	Calculation multiplier
Not less	than 600	and les	s than	1400	3	60	20
//	1400	//		3500	3	12	100
//	3500	//		9500	4	60	100
//	9500	//	9	9,500	4	6	1000
//	99,500				4	3	2000

Procedure of apparatus: Read the value after 2 minutes of rotation, and stop the rotation for at least 2 minutes. Repeat this procedure more two times, and average the three observed values.

pH <2.54> The pH of the sample solution obtained in the Viscosity, measured after 5 minutes immersing the electrode in the sample solution, is between 5.0 and 8.0.

Purity Heavy metals—Put 1.0 g of Methylcellulose in a 100-mL Kjeldahl flask, add a sufficient amount of a mixture of nitric acid and sulfuric acid (5:4) to wet the sample, and heat gently. Repeat this procedure until to use totally 18 mL of the mixture of nitric acid and sulfuric acid (5:4). Then boil gently until the solution changes to black. After cooling, add 2 mL of nitric acid, and heat until the solution changes to black. Repeat this procedure until the solution no longer changes to black, and heat strongly until dense white fumes are evolved. After cooling, add 5 mL of water, boil gently until dense white fumes are evolved, then heat until the volume of the solution becomes to 2 to 3 mL. After cooling, if the solution reveals yellow color by addition of 5 mL of water, add 1 mL of hydrogen peroxide (30), and heat until the volume of the solution becomes to 2 to 3 mL. After cooling, dilute the solution with 2 to 3 mL of water, transfer to a Nessler tube, add water to make 25 mL, and use this solution as the sample solution. Separately, put 2.0 mL of Standard Lead Solution in a 100-mL Kjeldahl flask, add 18 mL of the mixture of nitric acid and sulfuric acid (5:4) and an amount of nitric acid equal to that used for preparation of the sample solution, and heat until dense white fumes are evolved. After cooling, add 10 mL of water. In the case where hydrogen peroxide (30) is added for the preparation of the sample solution, add the same amount of hydrogen peroxide (30), then proceed in the same manner for preparation of the sample solution, and use so obtained solution as the control solution. Adjust the sample solution and the control solution to pH 3.0 to 4.0 with ammonia solution (28), and add water to make 40 mL, respectively. To these solutions add 1.2 mL of thioacetamide-alkaline glycerin TS, 2 mL of acetate buffer solution (pH 3.5) and water to make 50 mL, separately. After allowing to stand for 5 minutes, observe vertically both tubes on a white background: the color obtained with the sample solution is not more intense than that with the control solution (not more than 20 ppm).

Loss on drying $\langle 2.41 \rangle$ Not more than 5.0% (1 g, 105°C, 1 hour).

Residue on ignition $\langle 2.44 \rangle$ Not more than 1.5% (1 g).

Assay (i) Apparatus—Reaction vial: A 5-mL pressure-tight serum vial, having 20 mm in outside diameter and 50 mm in height, the neck 20 mm in outside diameter and 13 mm in inside diameter, equipped with a septum of butyl-rubber processed the surface with fluoroplastics, which can be fixed tightly to vial with aluminum cap, or equivalent.

Heater: A square-shaped aluminum block, having holes 20 mm in diameter and 32 mm in depth, adopted to the reaction vial. Capable of stirring the content of the reaction vial by means of magnetic stirrer or of reciprocal shaker about 100 times per minute.

(ii) Procedure—Weigh accurately about 65 mg of Methylcellulose, transfer to the reaction vial, add 60 to 100

mg of adipic acid, 2.0 mL of the internal standard solution and 2.0 mL of hydroiodic acid, immediately stopper the vial tightly, and weigh accurately. Using a magnetic stirrer equipped in the heating module, or using a shaker, stir for 60 minutes while heating so that the temperature of the vial content is 130 \pm 2°C. In the case when a magnetic stirrer or shaker is not available, heat for 30 minutes with repeated shaking at 5-minute intervals by hand, and continue heating for an additional 30 minutes. Allow the vial to cool, and again weigh accurately. If the mass loss is not more than 26 mg and there is no evidence of a leak of the content, use the upper layer of the mixture as the sample solution. Separately, put 60 to 100 mg of adipic acid, 2.0 mL of the internal standard solution and 2.0 mL of hydroiodic acid in a reaction vial, immediately stopper the vial tightly, and weigh accurately. Add 45 μ L of iodomethane for assay through the septum using a micro-syringe, weigh accurately, shake, and use the upper layer of the mixture as the standard solution. Perform the test with 1 to $2 \mu L$ each of the sample solution and standard solution as directed under Gas Chromatography <2.02> according to the following conditions, and calculate the ratios, Q_T and Q_S , of the peak area of iodomethane to that of the internal standard.

> Content (%) of methoxy group (CH₃O) = $M_S/M \times Q_T/Q_S \times 21.86$

M_S: Amount (mg) of iodomethane for assay takenM: Amount (mg) of Methylcellulose taken, calculated on the dried basis

Internal standard solution—A solution of *n*-octane in *o*-xylene (3 in 100).

Operating conditions—

Detector: A thermal conductivity detector or hydrogen flame-ionization detector.

Column: A fused silica column 0.53 mm in inside diameter and 30 m in length, coated with dimethylpolysiloxane for gas chromatography in 3 μ m thickness. Use a guard column, if necessary.

Column temperature: Maintain the temperature at 50°C for 3 minutes, raise to 100°C at a rate of 10°C per minute, then raise to 250°C at a rate of 35°C per minute, and maintain at 250°C for 8 minutes.

Injection port temperature: A constant temperature of about $250\,^{\circ}\text{C}$.

Detector temperature: A constant temperature of about $280^{\circ}\mathrm{C}$.

Carrier gas: Helium.

Flow rate: 4.3 mL per minute (the retention time of the internal standard is about 10 minutes).

Split ratio: 1:40. System suitability—

System performance: When the procedure is run with $1-2\,\mu\text{L}$ of the standard solution under the above operating conditions, iodomethane and the internal standard are eluted in this order with the resolution between these peaks being not less than 5

System repeatability: When the test is repeated 6 times

with $1-2 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of iodomethane to that of the internal standard is not more than 2.0%.

◆Containers and storage Containers—Well-closed containers.▲

Methylprednisolone

メチルプレドニゾロン

Change the Description and Optical rotation as follows:

Description Methylprednisolone occurs as a white crystalline powder. It is odorless.

It is sparingly soluble in methanol and in ethanol (99.5), and practically insoluble in water.

Melting point: 232 – 240°C (with decomposition).

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁵: +93 - +103° (after drying, 0.1 g, ethanol (99.5), 10 mL, 100 mm).

Add the following:

Minocycline Hydrochloride Granules

ミノサイクリン塩酸塩顆粒

Minocycline Hydrochloride Granules contain not less than 90.0% and not more than 110.0% of the labeled potency of minocycline ($C_{23}H_{27}N_3O_7$: 457.48).

Method of preparation Prepare as directed under Granules, with Minocycline Hydrochloride.

Identification To a quantity of Minocycline Hydrochloride Granules, equivalent to 10 mg (potency) of Minocycline Hydrochloride, add 625 mL of a solution of hydrochloric acid in methanol (19 in 20,000), shake thoroughly, and filter. Determine the absorption spectrum of the filtrate as directed under Ultraviolet-visible Spectrophotometry <2.24>: it exhibits maxima between 221 nm and 225 nm, between 261 nm and 265 nm, and between 354 nm and 358 nm.

Purity Related substances—Conduct this procedure within 30 minutes after the preparation of the sample solution. To a quantity of Minocycline Hydrochloride Granules, equivalent to 50 mg (potency) of Minocycline Hydrochloride, add 60 mL of the mobile phase, shake vigorously, add the mobile phase to make 100 mL. Centrifuge this solution, and use the supernatant liquid as the sample solution. Perform the test with 20 μ L of the sample solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions. Determine each peak area by the automatic in-

tegration method, and calculate their amounts by the area percentage method: the amount of the peak of epiminocycline, having the relative retention time of about 0.83 to minocycline, is not more than 4.0%.

Operating conditions—

Detector, column, column temperature, mobile phase, and flow rate: Proceed as directed in the operating conditions in the Assay under Minocycline Hydrochloride.

Time span of measurement: About 2.5 times as long as the retention time of minocycline, beginning after the solvent peak.

System suitability—

System performance: Proceed as directed in the system suitability in the Assay under Minocycline Hydrochloride.

Test for required detectability: To 2 mL of the standard solution obtained in the Assay add the mobile phase to make 100 mL, and use this solution as the solution for system suitability test. Pipet 5 mL of the solution for system suitability test, and add the mobile phase to make exactly 100 mL. Confirm that the peak area of minocycline obtained with 20 μ L of this solution is equivalent to 3.5 to 6.5% of that with 20 μ L of the solution for system suitability test.

System repeatability: When the test is repeated 6 times with $20 \,\mu\text{L}$ of the solution for system suitability test under the above operating conditions, the relative standard deviation of the peak area of minocycline is not more than 2.0%.

Water <2.48> Not more than 2.0% (4 g of powdered Minocycline Hydrochloride Granules, volumetric titration, back titration).

Uniformity of dosage units <6.02> Perform the test according to the following method: Minocycline Hydrochloride Granules in single-dose packages meet the requirement of the Content uniformity test.

To the total content of 1 package of Minocycline Hydrochloride Granules add water to disintegrate, shake thoroughly, add water to make exactly 100 mL, and filter through a membrane filter with a pore size not exceeding $0.45 \,\mu\text{m}$. Discard the first 10 mL of the filtrate, pipet $V \,\text{mL}$ of the subsequent filtrate, add water to make exactly V' mL so that each mL contains about 20 µg (potency) of minocycline (C₂₃H₂₇N₃O₇), and use this solution as the sample solution. Separately, weigh accurately an amount of Minocycline Hydrochloride RS, equivalent to about 20 mg (potency), dissolve in water to make exactly 100 mL. Pipet 5 mL of this solution, add water to make exactly 50 mL, and use this solution as the standard solution. Determine the absorbances, A_T and A_S , of the sample solution and standard solution at 348 nm as directed under Ultravioletvisible Spectrophotometry $\langle 2.24 \rangle$.

Amount [mg (potency)] of minocycline $(C_{23}H_{27}N_3O_7)$ = $M_S \times A_T/A_S \times V'/V \times 1/10$

 M_S : Amount [mg (potency)] of Minocycline Hydrochloride RS taken

Dissolution <6.10> When the test is performed at 50 revo-

lutions per minute according to the Paddle method, using 900 mL of water as the dissolution medium, the dissolution rate in 15 minutes of Minocycline Hydrochloride Granules is not less than 85%.

Start the test with an accurately weighed amount of Minocycline Hydrochloride Granules, equivalent to about 20 mg (potency) of Minocycline Hydrochloride, withdraw not less than 20 mL of the medium at the specified minute after starting the test, and filter through a membrane filter with a pore size not exceeding $0.45~\mu m$. Discard not less than 10~mL of the first filtrate, and use the subsequent filtrate as the sample solution. Separately, weigh accurately an amount of Minocycline Hydrochloride RS, equivalent to about 22 mg (potency), and dissolve in water to make exactly 100~mL. Pipet 5~mL of this solution, add water to make exactly 50~mL, and use this solution as the standard solution. Determine the absorbances, A_T and A_S , of the sample solution and standard solution at 348~mm as directed under Ultraviolet-visible Spectrophotometry <2.24>.

Dissolution rate (%) with respect to the labeled amount of minocycline ($C_{23}H_{27}N_3O_7$)

$$= M_{\rm S}/M_{\rm T} \times A_{\rm T}/A_{\rm S} \times 1/C \times 90$$

 M_S : Amount [mg (potency)] of Minocycline Hydrochloride RS taken

 M_T : Amount (g) of Minocycline Hydrochloride Granules taken

C: Labeled amount [mg (potency)] of minocycline $(C_{23}H_{27}N_3O_7)$ in 1 g

Assay Weigh accurately an amount of powdered Minocycline Hydrochloride Granules, equivalent to about 50 mg (potency) of Minocycline Hydrochloride, add the mobile phase, shake vigorously, and add the mobile phase to make exactly 100 mL. Centrifuge this solution, and use the supernatant liquid as the sample solution. Separately, weigh accurately an amount of Minocycline Hydrochloride RS, equivalent to about 25 mg (potency), dissolve in the mobile phase to make exactly 50 mL, and use this solution as the standard solution. Then, proceed as directed in the Assay under Minocycline Hydrochloride.

Amount [mg (potency)] of minocycline ($C_{23}H_{27}N_3O_7$) = $M_S \times A_T/A_S \times 2$

 M_S : Amount [mg (potency)] of Minocycline Hydrochloride RS taken

Containers and storage Containers—Tight containers. Storage—Light-resistant.

Add the following:

Nortriptyline Hydrochloride Tablets

ノルトリプチリン塩酸塩錠

Nortriptyline Hydrochloride Tablets contain not less than 95.0% and not more than 105.0% of the labeled amount of nortriptyline ($C_{19}H_{21}N$: 263.38).

Method of preparation Prepare as directed under Tablets, with Nortriptyline Hydrochloride.

Identification (1) Determine the absorption spectrum of the sample solution obtained in the Assay as directed under Ultraviolet-visible Spectrophotometry <2.24>, using diluted 0.1 mol/L hydrochloric acid TS (1 in 50) as the blank: it exhibits a maximum between 237 nm and 241 nm.

(2) To a quantity of powdered Nortriptyline Hydrochloride Tablets, equivalent to 10 mg of nortriptyline ($C_{19}H_{21}N$), add 10 mL of ethanol (99.5), shake thoroughly, centrifuge, and use the supernatant liquid as the sample solution. Separately, dissolve 11 mg of nortriptyline hydrochloride in 10 mL of ethanol (99.5), and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L each of the sample solution and standard solution on a plate of silica gel with fluorescent indicator for thin-layer chromatography. Develop the plate with a mixture of 1-butanol, water and acetic acid (100) (4:1:1) to a distance of about 10 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 254 nm): the principal spot obtained from the sample solution and the spot from the standard solution show the same Rf value.

Uniformity of dosage units <6.02> Perform the test according to the following method: it meets the requirement of the Content uniformity test.

To 1 tablet of Nortriptyline Hydrochloride Tablets add a suitable volume of 0.1 mol/L hydrochloric acid TS, disperse the fine particles by sonicating, add a suitable volume of 0.1 mol/L hydrochloric acid TS, sonicate, extract for 15 minutes while occasional shaking. Shake for 15 minutes, and add 0.1 mol/L hydrochloric acid TS to make exactly $V\,\mathrm{mL}$ so that each mL contains about 0.5 mg of nortriptyline (C₁₉H₂₁N). Centrifuge this solution, pipet 2 mL of the supernatant liquid, add water to make exactly 100 mL, and use this solution as the sample solution. Then, proceed as directed in the Assay.

Amount (mg) of nortriptyline ($C_{19}H_{21}N$) = $M_S \times A_T/A_S \times V/50 \times 0.878$

 M_S : Amount (mg) of nortriptyline hydrochloride for assay taken

Dissolution <6.10> When the test is performed at 50 revolutions per minute according to the Paddle method, using 900 mL of water as the dissolution medium, the dissolution

rates in 30 minutes of a 10-mg tablet and a 25-mg tablet are not less than 70% and not less than 80%, respectively.

Start the test with 1 tablet of Nortriptyline Hydrochloride Tablets, withdraw not less than 20 mL of the medium at the specified minute after starting the test, and filter through a membrane filter with a pore size not exceeding 0.45 μ m. Discard not less than 10 mL of the first filtrate, pipet V mL of the subsequent filtrate, add water to make exactly V' mL so that each mL contains about 11 µg of nortriptyline (C₁₉H₂₁N), and use this solution as the sample solution. Separately, weigh accurately about 25 mg of nortriptyline hydrochloride for assay, previously dried at 105°C for 2 hours, and dissolve in water to make exactly 100 mL. Pipet 5 mL of this solution, add water to make exactly 100 mL, and use this solution as the standard solution. Determine the absorbances, $A_{\rm T}$ and $A_{\rm S}$, of the sample solution and standard solution at 239 nm as directed under Ultravioletvisible Spectrophotometry $\langle 2.24 \rangle$.

Dissolution rate (%) with respect to the labeled amount of nortriptyline ($C_{19}H_{21}N$)

 $= M_{\rm S} \times A_{\rm T}/A_{\rm S} \times V'/V \times 1/C \times 45 \times 0.878$

 M_S : Amount (mg) of nortriptyline hydrochloride for assay taken

C: Labeled amount (mg) of nortriptyline (C₁₉H₂₁N) in 1 tablet

Assay Weigh accurately the mass of not less than 20 tablets of Nortriptyline Hydrochloride Tablets, and powder. Weigh accurately a portion of the powder, equivalent to about 50 mg of nortriptyline (C₁₉H₂₁N), add 50 mL of 0.1 mol/L hydrochloric acid TS, sonicate, and extract for 15 minutes while occasional shaking. Shake for 15 minutes, and add 0.1 mol/L hydrochloric acid TS to make exactly 100 mL. Centrifuge this solution, pipet 2 mL of the supernatant liquid, add water to make exactly 100 mL, and use this solution as the sample solution. Separately, weigh accurately about 28 mg of nortriptyline hydrochloride for assay, previously dried at 105°C for 2 hours, and dissolve in 0.1 mol/L hydrochloric acid TS to make exactly 50 mL. Pipet 2 mL of this solution, add water to make exactly 100 mL, and use this solution as the standard solution. Determine the absorbances, A_T and A_S , of the sample solution and standard solution at 239 nm as directed under Ultraviolet-visible Spectrophotometry <2.24>, using diluted 0.1 mol/L hydrochloric acid TS (1 in 50) as the blank.

> Amount (mg) of nortriptyline ($C_{19}H_{21}N$) = $M_S \times A_T/A_S \times 2 \times 0.878$

 M_S : Amount (mg) of nortriptyline hydrochloride for assay taken

Containers and storage Containers—Tight containers.

Pioglitazone Hydrochloride and Glimepiride Tablets

ピオグリタゾン塩酸塩・グリメピリド錠

Change the Identification (2) as follows:

Identification

(2) Wash the membrane filter obtained in (1) with 100 mL of 0.1 mol/L hydrochloric acid TS, and extract with methanol so that each mL contains about $10 \mu g$ of glimepiride ($C_{24}H_{34}N_4O_5S$). Determine the absorption spectrum of this solution as directed under Ultraviolet-visible Spectrophotometry $\langle 2.24 \rangle$: it exhibits a maximum between 227 nm and 231 nm.

Add the following:

Polaprezinc

ポラプレジンク

 $(C_9H_{12}N_4O_3Zn)_n$ catena-Poly{zinc- μ -[β -alanyl-L-histidinato(2-)-N, N^N ,O: N^T]} [107667-60-7]

Polaprezinc contains not less than 98.0% and not more than 102.0% of polaprezinc ($C_9H_{12}N_4O_3Zn$: 289.60), and contains not less than 21.5% and not more than 23.0% of zinc (Zn: 65.38), calculated on the anhydrous basis.

Description Polaprezinc occurs as a white to pale yellowwhite crystalline powder.

It is practically insoluble in water, in methanol and in ethanol (99.5).

It dissolves in dilute hydrochloric acid.

- **Identification** (1) To 2 mL of a solution of Polaprezinc in 0.2 mol/L hydrochloric acid TS (1 in 1000) add 0.5 mL of a solution of sulfanilic acid in 1 mol/L hydrochloric acid TS (1 in 200), 0.5 mL of a solution of sodium nitrite (1 in 20) and 3 mL of sodium carbonate TS: a red color is produced.
- (2) A solution of Polaprezinc in 0.2 mol/L hydrochloric acid TS (1 in 1000) responds to Qualitative Tests <1.09> for zinc salt.
- (3) Determine the infrared absorption spectrum of Polaprezinc as directed in the potassium bromide disk method under Infrared Spectrophotometry <2.25>, and compare the spectrum with the Reference Spectrum: both spectra exhibit similar intensities of absorption at the same wave numbers.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁰: $+8 - +9^{\circ}$ (1 g calculated on the anhydrous basis, 3 mol/L hydrochloric acid TS, 50 mL, 100 mm).

Purity (1) Lead—Weigh accurately about 0.5 g of Polaprezinc, dissolve in 3 mL of dilute nitric acid, add water to make exactly 10 mL, and use this solution as the sample solution. Separately, pipet 0.5 mL, 1.0 mL, 1.5 mL and 2.0 mL of Standard Lead Solution, to each solution add 3 mL of dilute nitric acid and water to make exactly 10 mL, and use these solutions as the standard solutions. Perform the test with the sample solution and standard solutions as directed under Atomic Absorption Spectrophotometry <2.23> according to the following conditions, and calculate the amount of lead in the sample solution using a calibration curve obtained from the absorbances of the standard solutions: not more than 10 ppm.

Gas: Combustible gas—Acetylene.

Supporting gas—Air.

Lamp: Lead hollow-cathode lamp.

Wavelength: 283.3 nm.

(2) Related substances—Dissolve 50 mg of Polaprezinc in 10 mL of 0.1 mol/L hydrochloric acid TS, add the mobile phase to make 100 mL, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add the mobile phase to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 10 μL each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions. Determine each peak area by the automatic integration method: the peak area of L-histidine, having the relative retention time of about 0.38 to Lcarnosine, obtained from the sample solution is not larger than 1/5 times of the peak area of L-carnosine from the standard solution, the area of the peak other than Lcarnosine and the peak mentioned above from the sample solution is not larger than 1/10 times of the peak area of Lcarnosine from the standard solution. Furthermore, the total area of the peaks other than L-carnosine from the sample solution is not larger than the peak area of L-carnosine from the standard solution.

Operating conditions—

Detector, column, column temperature, mobile phase, and flow rate: Proceed as directed in the operating conditions in the Assay.

Time span of measurement: About 4 times as long as the retention time of L-carnosine, beginning after the solvent peak.

System suitability—

Test for required detectability: Pipet 2 mL of the standard solution, add the mobile phase to make exactly 20 mL. Confirm that the peak area of L-carnosine obtained with 10 μ L of this solution is equivalent to 7 to 13% of that with 10 μ L of the standard solution.

System performance: Dissolve 50 mg each of Polaprezinc and L-histidine in 10 mL of 0.1 mol/L hydrochloric acid TS, and add the mobile phase to make 100 mL. When the procedure is run with $10 \,\mu\text{L}$ of this solution under the above

operating conditions, L-histidine and L-carnosine are eluted in this order with the resolution between these peaks being not less than 12.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of L-carnosine is not more than 2.0%.

Water <2.48> Not more than 5.0% (0.2 g, volumetric titration, direct titration, stir for 30 minutes).

Assay (1) Polaprezinc—Weigh accurately about 25 mg of Polaprezinc, dissolve in 5 mL of 0.1 mol/L hydrochloric acid TS, add the mobile phase to make exactly 100 mL, and use this solution as the sample solution. Separately, weigh accurately about 20 mg of L-Carnosine RS, previously dried at $105\,^{\circ}$ C for 3 hours, dissolve in 5 mL of 0.1 mol/L hydrochloric acid TS, add the mobile phase to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly $10\,\mu$ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions. Determine the peak areas, A_T and A_S , of L-carnosine in each solution.

Amount (mg) of polaprezinc (
$$C_9H_{12}N_4O_3Zn$$
)
= $M_S \times A_T/A_S \times 1.292$

M_S: Amount (mg) of L-Carnosine RS taken

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 210 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 45°C.

Mobile phase: Dissolve 1.4 g of potassium dihydrogen phosphate in 1000 mL of water, adjust to pH 3.5 with diluted phosphoric acid (1 in 100). Dissolve 2 g of sodium 1-octane sulfonate in 900 mL of this solution, and add 100 mL of acetonitrile for liquid chromatography.

Flow rate: Adjust so that the retention time of L-carnosine is about 15 minutes.

System suitability—

System performance: Dissolve 5 mg of L-histidine in 20 mL of the standard solution. When the procedure is run with $10\,\mu\text{L}$ of this solution under the above operating conditions, L-histidine and L-carnosine are eluted in this order with the resolution between these peaks being not less than 12.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of L-carnosine is not more than 1.0%.

(2) Zinc—Weigh accurately about 0.2 g of Polaprezinc, dissolve in 3 mL of dilute hydrochloric acid, and add water to make exactly 100 mL. Pipet 25 mL of this solution, add 10 mL of ammonia-ammonium chloride buffer solution (pH

10.7), and titrate <2.50> with 0.01 mol/L disodium dihydrogen ethylenediamine tetraacetate VS (indicator: 40 mg of eriochrome black T-sodium chloride indicator).

Each mL of 0.01 mol/L disodium dihydrogen ethylenediamine tetraacetate VS = 0.6538 mg Zn

Containers and storage Containers—Tight containers.

Add the following:

Polaprezinc Granules

ポラプレジンク顆粒

Polaprezinc Granules contain not less than 95.0% and not more than 105.0% of the labeled amount of polaprezinc $[(C_9H_{12}N_4O_3Zn)_n]$.

Method of preparation Prepare as directed under Granules, with Polaprezinc.

Identification (1) To a quantity of Polaprezinc Granules, equivalent to 20 mg of Polaprezinc, add 20 mL of 0.2 mol/L hydrochloric acid TS, shake for 10 minutes, centrifuge, and use the supernatant liquid as the sample solution. To 2 mL of the sample solution add 0.5 mL of a solution of sulfanilic acid in 1 mol/L hydrochloric acid TS (1 in 200), 0.5 mL of a solution of sodium nitrite (1 in 20) and 3 mL of sodium carbonate TS: a red color develops.

(2) The sample solution obtained in (1) responds to Qualitative Tests $\langle 1.09 \rangle$ for zinc salt.

Uniformity of dosage units <6.02> Perform the test according to the following method: Polaprezinc Granules in single-dose packages meet the requirement of the Content uniformity test.

To the total content of 1 package of Polaprezinc Granules add exactly $V \, \text{mL}$ of $0.2 \, \text{mol/L}$ hydrochloric acid TS so that each mL contains about 5 mg of polaprezinc $[(C_9H_{12}N_4O_3Zn)_n]$, shake vigorously for 10 minutes, and centrifuge. Pipet 5 mL of the supernatant liquid, add exactly 5 mL of the internal standard solution, add the mobile phase to make 50 mL, and use this solution as the sample solution. Then, proceed as directed in the Assay.

Amount (mg) of polaprezinc
$$[(C_9H_{12}N_4O_3Zn)_n]$$

= $M_S \times Q_T/Q_S \times V/5 \times 1.292$

M_S: Amount (mg) of L-Carnosine RS taken

Dissolution <6.10> When the test is performed at 50 revolutions per minute according to the Paddle method, using 900 mL of 0.05 mol/L acetic acid-sodium acetate buffer solution (pH 4.0) as the dissolution medium, the dissolution rate in 15 minutes of Polaprezinc Granules is not less than 80%.

Start the test with an accurately weighed amount of Polaprezinc Granules, equivalent to about 75 mg of polaprezinc $[(C_9H_{12}N_4O_3Zn)_n]$, withdraw not less than 20

mL of the medium at the specified minute after starting the test, and filter through a membrane filter with a pore size not exceeding 0.45 μ m. Discard not less than 10 mL of the first filtrate, pipet 1 mL of the subsequent filtrate, add diluted nitric acid (77 in 10,000) to make exactly 25 mL, and use this solution as the sample solution. Separately, pipet suitable volumes of Standard Zinc Stock Solution, to each solution add diluted nitric acid (77 in 10,000) so that each mL contains 0.4 to 0.8 μ g of zinc (Zn: 65.38), and use these solutions as the standard solutions. Perform the test with the sample solution and standard solutions as directed under Atomic Absorption Spectrophotometry <2.23> according to the following conditions, and calculate the amount of zinc in the sample solution using a calibration curve obtained from the absorbances of the standard solutions.

Dissolution rate (%) with respect to the labeled amount of polaprezinc $[(C_9H_{12}N_4O_3Zn)_n]$

= Content (μ g/mL) of zinc in the sample solution/ M_T × 1/C × 2250 × 4.429

 $M_{\rm T}$: Amount (g) of Polaprezinc Granules taken

C: Labeled amount (mg) of polaprezinc $[(C_9H_{12}N_4O_3Zn)_n]$ in 1 g

Gas: Combustible gas—Acetylene.

Supporting gas—Air.

Lamp: Zinc hollow-cathode lamp.

Wavelength: 213.9 nm.

Assay Weigh accurately an amount of Polaprezinc Granules, equivalent to about 0.1 g of polaprezinc [$(C_9H_{12}N_4O_3Zn)_n$], add exactly 20 mL of 0.2 mol/L hydrochloric acid TS, shake vigorously for 10 minutes, and centrifuge. Pipet 5 mL of the supernatant liquid, add exactly 5 mL of the internal standard solution, add the mobile phase to make 50 mL, and use this solution as the sample solution. Separately, weigh accurately about 20 mg of L-Carnosine RS, previously dried at 105°C for 3 hours, dissolve in 5 mL of 0.2 mol/L hydrochloric acid TS, add exactly 5 mL of the internal standard solution, add the mobile phase to make 50 mL, and use this solution as the standard solution. Perform the test with $5 \mu L$ each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and calculate the ratios, Q_T and Q_S , of the peak area of L-carnosine to that of the internal standard.

Amount (mg) of polaprezinc $[(C_9H_{12}N_4O_3Zn)_n]$ = $M_S \times Q_T/Q_S \times 4 \times 1.292$

M_S: Amount (mg) of L-Carnosine RS taken

Internal standard solution—Dissolve 0.25 g of 4-aminoacetophenone in 5 mL of acetonitrile, and add the mobile phase to make 100 mL.

Operating conditions—

Proceed as directed in the operating conditions in the Assay (1) under Polaprezinc.

System suitability—

System performance: When the procedure is run with 5 μ L of the standard solution under the above operating conditions, 4-aminoacetophenone and L-carnosine are eluted in this order with the resolution between these peaks being not less than 6.

System repeatability: When the test is repeated 6 times with $5 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of L-carnosine to that of the internal standard is not more than 1.0%.

Containers and storage Containers—Tight containers.

Add the following:

Ritodrine Hydrochloride Injection

リトドリン塩酸塩注射液

Ritodrine Hydrochloride Injection is an aqueous injection.

It contains not less than 95.0% and not more than 105.0% of the labeled amount of ritodrine hydrochloride ($C_{17}H_{21}NO_3$.HCl: 323.81).

Method of preparation Prepare as directed under Injections, with Ritodrine Hydrochloride.

Manufacture Ritodrine Hydrochloride Injection is produced by the formulation and the manufacturing method to ensure that the amounts of related substances do not exceed the limit values of related substances under Ritodrine Hydrochloride.

Description Ritodrine Hydrochloride Injection is a clear and colorless liquid.

Identification To a volume of Ritodrine Hydrochloride Injection, equivalent to 50 mg of Ritodrine Hydrochloride, add 0.01 mol/L hydrochloric acid TS to make 100 mL. To 10 mL of this solution add 0.01 mol/L hydrochloric acid TS to make 100 mL. Determine the absorption spectrum of this solution as directed under Ultraviolet-visible Spectrophotometry <2.24>: it exhibits a maximum between 272 nm and 276 nm.

pH Being specified separately when the drug is granted approval based on the Law.

Bacterial endotoxins <4.01> Less than 25 EU/mg.

Extractable volume <6.05> It meets the requirement.

Foreign insoluble matter <6.06> Perform the test according to Method 1: it meets the requirement.

Insoluble particulate matter <6.07> It meets the requirement.

Sterility <4.06> Perform the test according to the Membrane filtration method: it meets the requirement.

Assay Pipet a volume of Ritodrine Hydrochloride Injection, equivalent to about 20 mg of ritodrine hydrochloride ($C_{17}H_{21}NO_3$.HCl), and add a mixture of 0.02 mol/L sodium dihydrogen phosphate dihydrate solution and methanol (7:3) to make exactly 250 mL, and use this solution as the sample solution. Separately, weigh accurately about 20 mg of Ritodrine Hydrochloride RS, previously dried at 105°C for 2 hours, dissolve in a mixture of 0.02 mol/L sodium dihydrogen phosphate dihydrate solution and methanol (7:3) to make exactly 250 mL, and use this solution as the standard solution. Perform the test with exactly 10 μ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions. Determine the peak areas, A_T and A_S , of ritodrine in each solution.

Amount (mg) of ritodrine hydrochloride ($C_{17}H_{21}NO_3.HCl$) = $M_S \times A_T/A_S$

M_S: Amount (mg) of Ritodrine Hydrochloride RS taken

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 220 nm).

Column: A stainless steel column 6 mm in inside diameter and 15 cm in length, packed with octylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about $25\,^{\circ}\mathrm{C}$.

Mobile phase: Dissolve 6.6 g of diammonium hydrogen phosphate and 1.1 g of sodium 1-heptansulfonate in 840 mL of water, add 160 mL of acetonitrile for liquid chromatography, and adjust to pH 3.0 with phosphoric acid.

Flow rate: Adjust so that the retention time of ritodrine is about 19 minutes.

System Suitability-

System performance: Dissolve 10 mg of ritodrine hydrochloride in 50 mL of dilute sulfuric acid. Heat a portion of this solution in a water bath for about 30 minutes, and allow to cool. Measure a portion of this solution, and add the same volume of 2 mol/L sodium hydroxide TS. Dissolve 2 mg of ritodrine hydrochloride in 10 mL of this solution, and add a mixture of 0.02 mol/L sodium dihydrogen phosphate dihydrate solution and methanol (7:3) to make 25 mL. When the procedure is run with 10 μ L of this solution under the above operating conditions, ritodrine and ritodrine threo-isomer are eluted in this order with the resolution between these peaks being not less than 3.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of ritodrine is not more than 1.0%.

Containers and storage Containers—Hermetic containers. Storage—At a temperature between 2°C and 8°C.

Saccharin

サッカリン

Delete the description of harmonization in the beginning and symbols (\blacklozenge \blacklozenge) in the Description, Melting point, Purity (2) and (4), and Containers and storage.

Saccharin Sodium Hydrate

サッカリンナトリウム水和物

Delete the description of harmonization in the beginning and symbols ($^{\blacklozenge}$ $_{\blacklozenge}$) in the Description, Purity (3) and (5), and Containers and storage.

Light Anhydrous Silicic Acid

軽質無水ケイ酸

Delete the Volume test:

Add the following:

Sitagliptin Phosphate Hydrate

シタグリプチンリン酸塩水和物

 $C_{16}H_{15}F_6N_5O.H_3PO_4.H_2O:$ 523.32 (3*R*)-3-Amino-1-[3-(trifluoromethyl)-5,6-dihydro[1,2,4]triazolo[4,3-*a*]pyrazin-7(8*H*)-yl]-4-(2,4,5-trifluorophenyl)butan-1-one monophosphate monohydrate [654671-77-9]

Sitagliptin Phosphate Hydrate contains not less than 98.0% and not more than 102.0% of sitagliptin phosphate ($C_{16}H_{15}F_6N_5O.H_3PO_4$: 505.31), calculated on the anhydrous basis.

Description Sitagliptin Phosphate Hydrate occurs as a white powder.

It is soluble in water, sparingly soluble in methanol, very slightly soluble in acetonitrile and in ethanol (99.5).

Identification (1) Determine the absorption spectrum of a solution of Sitagliptin Phosphate Hydrate (1 in 10,000) as directed under Ultraviolet-visible Spectrophotometry

- <2.24>, and compare the spectrum with the Reference Spectrum or the spectrum of a solution of Sitagliptin Phosphate RS prepared in the same manner as the sample solution: both spectra exhibit similar intensities of absorption at the same wavelengths.
- (2) Determine the infrared absorption spectrum of Sitagliptin Phosphate Hydrate as directed in the paste method under Infrared Spectrophotometry <2.25>, and compare the spectrum with the Reference Spectrum or the spectrum of Sitagliptin Phosphate RS: both spectra exhibit similar intensities of absorption at the same wave numbers. Alternatively, perform the test by the potassium bromide disk method or the ATR method, and compare the spectrum with the spectrum of Sitagliptin Phosphate RS: both spectra exhibit similar intensities of absorption at the same wave numbers.
- (3) A solution of Sitagliptin Phosphate Hydrate (1 in 25) responds to Qualitative Tests <1.09> (1) for phosphate.
- **Purity** (1) Heavy metals—Being specified separately when the drug is granted approval based on the Law.
- (2) Related substances—Use the sample solution obtained in the Assay as the sample solution. Pipet 1 mL of the sample solution, add a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1) to make exactly 1000 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions. Determine each peak area by the automatic integration method: the area of the peak other than sitagliptin obtained from the sample solution is not larger than the peak area of sitagliptin from the standard solution, and the total area of the peaks other than sitagliptin from the sample solution is not larger than 5 times the peak area of sitagliptin from the standard solution. For this calculation the peak area not larger than 1/2 times the peak area of sitagliptin from the standard solution is excluded.

Operating conditions—

Detector, column, column temperature, mobile phase and flow rate: Proceed as directed in the operating conditions in the Assay.

Time span of measurement: About 5.5 times as long as the retention time of sitagliptin, beginning after the solvent peak.

System suitability-

System performance: Proceed as directed in the system suitability in the Assay.

Test for required detectability: Pipet 5 mL of the standard solution, and add a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1) to make exactly 10 mL. When the procedure is run with 20 μ L of this solution under the above operating conditions, the SN ratio of the peak of sitagliptin is not less than 10.

System repeatability: When the test is repeated 6 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak

area of sitagliptin is not more than 2.0%.

(3) Enantiomer—Dissolve 80 mg of Sitagliptin Phosphate Hydrate in a mixture of methanol and water (9:1) to make 10 mL, and use this solution as the sample solution. Perform the test with $10 \,\mu\text{L}$ of the sample solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions. Determine the total peak area, A_{T} , of sitagliptin and related substance A (enantiomer), having the relative retention time of about 0.9 to sitagliptin and the peak area, A_{S} , of related substance A (enantiomer), and calculate the amount of the enantiomer by the following equation: not more than 0.5%.

Amount (%) of enantiomer = $A_S/A_T \times 100$

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 268 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 25 cm in length, packed with amylose tris-(3,5-dimethylphenylcarbamate)-coated silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 35° C.

Mobile phase: A mixture of ethanol (99.5), heptane, water and diethylamine (600:400:1:1).

Flow rate: 0.8 mL per minute.

System suitability—

Test for required detectability: Pipet 1 mL of the sample solution, and dissolve in a mixture of methanol and water (9:1) to make exactly 100 mL. Pipet 1 mL of this solution, add a mixture of methanol and water (9:1) to make exactly 10 mL. When the procedure is run with $10 \,\mu\text{L}$ of this solution under the above operating conditions, the SN ratio of the peak of sitagliptin is not less than 10.

System performance: Dissolve 8 mg of Sitagliptin Phosphate RS for System Suitability in 1 mL of a mixture of methanol and water (9:1). When the procedure is run with 10μ L of this solution under the above operating conditions, the resolution between the peaks of related substance A (enantiomer) and sitagliptin is not less than 1.5.

Water $\langle 2.48 \rangle$ 3.3 - 3.7% (0.3 g, volumetric titration, direct titration).

Residue on ignition <2.44> Not more than 0.2% (1 g, platinum crucible).

Assay Weigh accurately about 20 mg each of Sitagliptin Phosphate Hydrate and Sitagliptin Phosphate RS (separately determine the water $\langle 2.48 \rangle$ in the same manner as Sitagliptin Phosphate Hydrate), dissolve each in a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1) to make exactly 200 mL, and use these solutions as the sample solution and the standard solution, respectively. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and determine the peak areas, A_T and A_S , of sitagliptin in each solution.

Amount (mg) of sitagliptin phosphate $(C_{16}H_{15}F_6N_5O.H_3PO_4)$ $= M_S \times A_T/A_S$

 M_S : Amount (mg) of Sitagliptin Phosphate RS taken, calculated on the anhydrous basis

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 205 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with cyanopropylsilanized silica gel for liquid chromatography (5 μ m in particle diameter)

Column temperature: A constant temperature of about 30°C.

Mobile phase: Dissolve 1.36 g of potassium dihydrogen phosphate in 900 mL of water, adjust to pH 2.0 with phosphoric acid, and add water to make 1000 mL. To 850 mL of this solution add 150 mL of acetonitrile for liquid chromatography.

Flow rate: 1.0 mL per minute.

System suitability—

System performance: Place 10 mg of Sitagliptin Phosphate RS and 1 mg of sodium stearyl fumarate in a vial, and add 1 mL of water. Stopper the vial tightly, and heat at 80°C for 20 to 48 hours. Take out the contents of the vial, wash the vial 3 times with a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1), combine the washings and the content, and add a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1) to make 100 mL. Stir this solution for 1 hour, and centrifuge for 10 minutes or until the solution becomes clear. Use the supernatant liquid as the solution for system suitability test. When the procedure is run with 20 μ L of the solution for system suitability test under the above operating conditions, the resolution between sitagliptin and the peak having the relative retention time of about 1.2 to sitagliptin is not less than 1.5.

System repeatability: When the test is repeated 6 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of sitagliptin is not more than 1.0%.

Containers and storage Containers—Tight containers.

Others

butan-1-one

Related substance A (enantiomer): (3*S*)-3-Amino-1-[3-(trifluoromethyl)-5,6-dihydro[1,2,4] triazolo[4,3-*a*]pyrazin-7(8*H*)-yl]-4-(2,4,5-trifuluorophenyl)

Add the following:

Sitagliptin Phosphate Tablets

シタグリプチンリン酸塩錠

Sitagliptin Phosphate Tablets contain not less than 95.0% and not more than 105.0% of the labeled amount of sitagliptin ($C_{16}H_{15}F_6N_5O$: 407.31).

Method of preparation Prepare as directed under Tablets, with Sitagliptin Phosphate Hydrate.

Manufacture The management strategy of Sitagliptin Phosphate Tablets is based on systematic development methods, which put emphasis on prior setting targets, understanding of products and processes, and process control, and which is based on quality risk management and proven science. In addition when it can be scientifically possible to explain that a disintegration test ensure quality with distinguishability equal or better than a dissolution test, the following disintegration is alternative for the estimation of dissolution.

Disintegration <6.09> Perform the test for 5 minutes: it meets the requirement.

Identification (1) To 1 tablet of Sitagliptin Phosphate Tablets add water so that each mL contains about 0.2 mg of sitagliptin ($C_{16}H_{15}F_6N_5O$), and shake thoroughly to disintegrate. Centrifuge this solution. Determine the absorption spectrum of the supernatant liquid as directed under Ultraviolet-visible Spectrophotometry $\langle 2.24 \rangle$: it exhibits a maximum between 265 nm and 269 nm.

(2) Perform the test with $20 \,\mu\text{L}$ of the sample solution and standard solution obtained in the Assay as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the operating conditions in the Assay: the retention times of the principal peaks in the chromatograms obtained from the sample solution and the standard solution are the same.

Purity Related substances—Use the sample solution obtained in the Assay as the sample solution. Separately, pipet 1 mL of the standard solution obtained in the Assay, add a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1) to make exactly 500 mL, and use this solution as the standard solution. Perform the test with exactly 20 µL each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions. Determine the peak area, A_T , of related substance obtained from the sample solution and the peak area, A_S , of sitagliptin from the standard solution, and calculate the amount of related substances by the following equation: the total amount of related substances is not more than 0.2%. For this calculation the amount of related substance not more than 0.1% is excluded.

Amount (%) of related substance = $M_S \times A_T/A_S \times V'/V \times 1/C \times 1/50 \times 0.806$

 M_S : Amount (mg) of Sitagliptin Phosphate RS taken, cal-

culated on the anhydrous basis

V'/V: Dilution factor for the sample solution in the Assay

C: Labeled amount (mg) of sitagliptin (C₁₆H₁₅F₆N₅O) in 1 tablet

Operating conditions—

Detector, column, column temperature, mobile phase and flow rate: Proceed as directed in the operating conditions in the Assay.

Time span of measurement: About 5.5 times as long as the retention time of sitagliptin, beginning after the solvent peak.

System suitability—

System performance and system repeatability: Proceed as directed in the system suitability in the Assay.

Test for required detectability: To 5 mL of the standard solution add a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1) to make exactly 10 mL. When the procedure is run with 20 μ L of this solution under the above operating conditions, the SN ratio of the peak of sitagliptin is not less than 10.

Uniformity of dosage units < 6.02> Perform the Mass variation test, or the Content uniformity test according to the following method: it meets the requirement.

To 1 tablet of Sitagliptin Phosphate Tablets add a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1) to make exactly 25 mL, and stir thoroughly. Pipet V mL of this solution, and add a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1) to make exactly V' mL so that each mL contains about 80 μ g of sitagliptin (C₁₆H₁₅F₆N₅O). Centrifuge this solution, and use the supernatant liquid as the sample solution. Then, proceed as directed in the Assay.

Amount (mg) of sitagliptin (
$$C_{16}H_{15}F_6N_5O$$
)
= $M_S \times A_T/A_S \times V'/V \times 1/10 \times 0.806$

 $M_{\rm S}$: Amount (mg) of Sitagliptin Phosphate RS taken, calculated on the anhydrous basis

Dissolution <6.10> When the test is performed at 100 revolutions per minute according to the Basket method, using 900 mL of water as the dissolution medium, the dissolution rate in 15 minutes of Sitagliptin Phosphate Tablets is not less than 85%.

Start the test with 1 tablet of Sitagliptin Phosphate Tablets, withdraw not less than 4 mL of the medium at the specified minute after starting the test, and filter through a membrane filter with a pore size of $0.45 \,\mu\text{m}$. Discard not less than 2 mL of the first filtrate, pipet V mL of the subsequent filtrate, add water to make exactly V' mL so that each mL contains about $14 \,\mu\text{g}$ of sitagliptin ($C_{16}H_{15}F_6N_5O$), and use this solution as the sample solution. Separately, weigh accurately about 29 mg of Sitagliptin Phosphate RS (separately determine the water <2.48> in the same manner as Sitagliptin Phosphate Hydrate), and dissolve in a solution of sodium chloride (37 in 25,000) to make exactly 100

mL. Pipet 6 mL of this solution, and add a solution of sodium chloride (37 in 25,000) to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, A_T and A_S , of sitagliptin in each solution.

Dissolution rate (%) with respect to the labeled amount of sitagliptin ($C_{16}H_{15}F_6N_5O$)

$$= M_{\rm S} \times A_{\rm T}/A_{\rm S} \times V'/V \times 1/C \times 54 \times 0.806$$

M_S: Amount (mg) of Sitagliptin Phosphate RS taken, calculated on the anhydrous basis

C: Labeled amount (mg) of sitagliptin ($C_{16}H_{15}F_6N_5O$) in 1 tablet

Operating conditions—

Column, column temperature and flow rate: Proceed as directed in the operating conditions in the Assay.

Detector: An ultraviolet absorption photometer (wavelength: 267 nm).

Mobile phase: Dissolve 1.36 g of potassium dihydrogen phosphate in 900 mL of water, adjust to pH 2.0 with phosphoric acid, and add water to make 1000 mL. To 750 mL of this solution add 250 mL of acetonitrile for liquid chromatography.

System suitability—

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of sitagliptin are not less than 5000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of sitagliptin is not more than 1.0%.

Assay To 10 Sitagliptin Phosphate Tablets add a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1) to make exactly 250 mL, and stir thoroughly. Pipet VmL of this solution, and add a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1) to make exactly V' mL so that each mL contains about 80 μ g of sitagliptin (C₁₆H₁₅F₆N₅O). Centrifuge this solution, and use the supernatant liquid as the sample solution. Separately, weigh accurately about 26 mg of Sitagliptin Phosphate RS (separately determine the water <2.48> in the same manner as Sitagliptin Phosphate Hydrate), dissolve in a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1) to make exactly 250 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, $A_{\rm T}$ and $A_{\rm S}$, of sitagliptin in each solution.

Amount (mg) of sitagliptin ($C_{16}H_{15}F_6N_5O$) in 1 tablet of Sitagliptin Phosphate Tablets

 $= M_{\rm S} \times A_{\rm T}/A_{\rm S} \times V'/V \times 1/10 \times 0.806$

M_S: Amount (mg) of Sitagliptin Phosphate RS taken, calculated on the anhydrous basis

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 205 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with cyanopropylsilanized silica gel for liquid chromatography (5 μ m in particle diameter)

Column temperature: A constant temperature of about 30°C.

Mobile phase: Dissolve 1.36 g of potassium dihydrogen phosphate in 900 mL of water, adjust to pH 2.0 with phosphoric acid, and add water to make 1000 mL. To 850 mL of this solution add 150 mL of acetonitrile for liquid chromatography.

Flow rate: 1.0 mL per minute.

System suitability—

System performance: Proceed as directed in the system suitability in the Assay under Sitagliptin Phosphate Hydrate. The following method can be applied when sodium stearyl fumarate is contained in the additive of the tablet.

Crush 1 tablet of Sitagliptin Phosphate Tablets, transfer to a vial, and add 1 mL of water. Stopper the vial tightly, and heat at 80°C for 20 to 48 hours. Take out the contents of the vial, wash the vial 3 times with a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1), combine the washings and the content, and add a mixture of diluted phosphoric acid (1 in 1000) and acetonitrile for liquid chromatography (19:1) to make 100 mL. Stir this solution for 1 hour, and centrifuge for 10 minutes or until the solution becomes clear. When the procedure is run with 20 μ L of the supernatant liquid under the above operating conditions, the resolution between sitagliptin and the peak having the relative retention time of about 1.2 to sitagliptin is not less than 1.5.

System repeatability: When the test is repeated 6 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of sitagliptin is not more than 1.0%.

Containers and storage Containers—Tight containers.

Sodium Fusidate

フシジン酸ナトリウム

Change the Structural formula as follows:

$$H_3C$$
 CH_3 CO_2Na CO_2N

Change the purity as follows:

Purity (1) Heavy metals <1.07>—Proceed with 2.0 g of Sodium Fusidate according to Method 4, and perform the test. Prepare the control solution with 2.0 mL of Standard Lead Solution (not more than 10 ppm).

(2) Related substances—Dissolve 25 mg of Sodium Fusidate in a mixture of acetonitrile for liquid chromatography, diluted phosphoric acid (3 in 1000) and methanol (5:4:1) to make 10 mL, and use this solution as the sample solution. Pipet 1 mL of the sample solution, add water to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 20 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions. Determine each peak area by the automatic integration method: the peak area of related substance A, having the relative retention time of about 0.4 to fusidic acid, obtained from the sample solution is not larger than 3/10 times the peak area of fusidic acid from the standard solution, the peak area of related substance B, having the relative retention time of about 0.5, from the sample solution is not larger than 2/5 times the peak area of fusidic acid from the standard solution, the peak areas of related substance C having the relative retention time of about 0.6, related substance D having the relative retention time of about 0.63, an unknown substance having the relative retention time of about 0.65, related substance E having the relative retention time of about 0.7, related substance G having the relative retention time of about 0.96 and related substance H having the relative retention time of about 1.18, from the sample solution are not larger than 1/5 times the peak area of fusidic acid from the standard solution, the peak area of related substance F, having the relative retention time of about 0.82, from the sample solution is not larger than 7/10 times the peak area of fusidic acid from the standard solution, the peak area of related substance I, having the relative retention time of about 1.23, from the sample solution is not larger than 1/2 times the peak area of fusidic acid from the standard solution, the peak area of related substance J, having the relative retention time of about 1.4, from the sample solution is not larger than the peak area of

fusidic acid from the standard solution, the area of the peak other than fusidic acid and the peaks mentioned above from the sample solution is not larger than 1/10 times the peak area of fusidic acid from the standard solution. Furthermore, the total area of the peaks other than fusidic acid from the sample solution is not larger than 2 times the peak area of fusidic acid from the standard solution. For the areas of the peaks, related substances C, D, E, G and H, multiply their relative response factors, 0.7, 0.7, 0.3, 0.6 and 0.6, respectively.

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 235 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (3.5 μ m in particle diameter).

Column temperature: A constant temperature of about 30°C .

Mobile phase A: A mixture of diluted phosphoric acid (3 in 1000), acetonitrile for liquid chromatography and methanol (2:2:1).

Mobile phase B: A mixture of acetonitrile for liquid chromatography, methanol and diluted phosphoric acid (3 in 1000) (7:2:1).

Flowing of mobile phase: Control the gradient by mixing the mobile phases A and B as directed in the following table.

Time after injection of sample (min)	Mobile phase A (vol%)	Mobile phase B (vol%)
0 - 3	100	0
3 – 28	$100 \rightarrow 0$	$0 \rightarrow 100$
28 - 33	0	100

Flow rate: 1.0 mL per minute.

Time span of measurement: For 33 minutes after injection, beginning after the solvent peak.

System suitability—

Test for required detectability: Pipet 1 mL of the standard solution, and add a mixture of acetonitrile for liquid chromatography, diluted phosphoric acid (3 in 1000) and methanol (5:4:1) to make exactly 20 mL. Confirm that the peak area of fusidic acid obtained with 20 μ L of this solution is equivalent to 3.5 to 6.5% of that with 20 μ L of the standard solution.

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of fusidic acid are not less than 43,000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with $20 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of fusidic acid is not more than 2.0%.

Add the following next to the Containers and storage:

Others

Related substance A:

(24*RS*,17*Z*)-*ent*-16α-Acetoxy-3 β ,11 β ,24,25-tetrahydroxy-4 β ,8 β ,14 α -trimethyl-18-nor-5 β ,10 α -cholest-17(20)-en-21-oic acid

and epimer at C*

Related substance B:

(17Z)-ent-3 β ,11 β -Dihydroxy-17-[(6SR)-

6-hydroxy-7,7-dimethyl-2-oxooxepan-3-ylidene]- 4β , 8β , 14α -trimethyl-18-nor- 5β , 10α -androstan- 16α -yl acetate

Related substance C:

(17Z)-ent- 3β , 11β -Dihydroxy-17-[(6S)-6-(2-hydroxypropan-2-yl)-2-oxodihydro-2H-pyran-3(4H)-ylidene]- 4β , 8β , 14α -trimethyl-18-nor- 5β , 10α -androstan- 16α -yl acetate

Related substance D:

(17*Z*)-ent-3 β ,11 β -Dihydroxy-17-[(6*R*)-6-(2-hydroxypropan-2-yl)-2-oxodihydro-2*H*-pyran-3(4*H*)-ylidene]-4 β ,8 β ,14 α -trimethyl-18-nor-5 β ,10 α -androstan-16 α -yl acetate

Related substance E:

(17*Z*,24*EZ*)-ent-16 α -Acetoxy-3 β ,11 β -dihydroxy-4 β ,8 β ,14 α -trimethyl-26-oxo-18-nor-5 β ,10 α -cholesta-17(20),24-dien-21-oic acid

Related substance F:

(17*Z*)-*ent*-16 α -Acetoxy-11 β -hydroxy-4 β ,8 β ,14 α -trimethyl-3-oxo-18-nor-5 β ,10 α -cholesta-17(20),24-dien-21-oic acid

Related substance G:

(17*Z*)-*ent*-3 β ,11 β ,16 β -Trihydroxy-4 β ,8 β ,14 α -trimethyl-18-nor-5 β ,10 α -cholesta-17(20),24-dien-21-oic acid

$$H_3C$$
 CH_3 H CO_2H H CH_3 CH_3

Related substance H:

(17*Z*)-ent-3 β ,11 β -Dihydroxy-4 β ,8 β ,14 α -trimethyl-18-nor-5 β ,10 α -cholesta-17(20),24-dieno-21,16 α -lactone

Related substance I:

(17*Z*)-*ent*-16 α -Acetoxy-3 β -hydroxy-4 β ,8 β ,14 α -trimethyl-18-nor-5 β ,10 α -cholesta-9(11),17(20),24-trien-21-oic acid

Related substance J:

(17*Z*)-ent-16 α -Acetoxy-3 β -hydroxy-4 β ,8 β ,14 α -trimethyl-18-nor-5 β ,10 α -cholesta-17(20),24-dien-21-oic acid

Add the following:

Sodium Valproate Extended-release Tablets A

バルプロ酸ナトリウム徐放錠 A

Sodium Valproate Extended-release Tablets A contain not less than 95.0% and not more than 105.0% of the labeled amount of sodium valproate $(C_8H_{15}NaO_2: 166.19)$.

Method of preparation Prepare as directed under Tablets, with Sodium Valproate.

Identification To a quantity of powdered Sodium Valproate Extended-release Tablets A, equivalent to 0.2 g of Sodium Valproate, add 20 mL of water, shake thoroughly, and centrifuge. To 2 mL of the supernatant liquid add 1 mL of a solution of cobalt (II) nitrate hexahydrate (1 in 20), and heat on a water bath: a purple precipitate is formed.

Uniformity of dosage units <6.02> Perform the test according to the following method: it meets the requirement of the Content uniformity test.

Crush 1 tablet of Sodium Valproate Extended-release Tablets A, add exactly V/40 mL of the internal standard solution, add 4V/5 mL of a mixture of methanol and water (3:2), shake vigorously, add a mixture of methanol and water (3:2) to make V mL so that each mL contains about 1 mg of sodium valproate ($C_8H_{15}NaO_2$), and filter through a membrane filter with a pore size not exceeding $0.45 \, \mu \text{m}$. Discard the first 5 mL of the filtrate, and use the subsequent filtrate as the sample solution. Separately, weigh accurately about 0.1 g of sodium valproate for assay, previously dried at $105\,^{\circ}\text{C}$ for 3 hours, add exactly $2.5 \, \text{mL}$ of the internal standard solution, add a mixture of methanol and water

(3:2) to make 100 mL, and use this solution as the standard solution. Then, proceed as directed in the Assay.

Amount (mg) of sodium valproate (
$$C_8H_{15}NaO_2$$
)
= $M_S \times Q_T/Q_S \times V/100$

 $M_{\rm S}$: Amount (mg) of sodium valproate for assay taken

Internal standard solution—A solution of ethyl parahydroxybenzoate in a mixture of methanol and water (3:2) (1 in 5000).

Dissolution <6.10> When the test is performed at 50 revolutions per minute according to the Paddle method, using 900 mL of water as the dissolution medium, the dissolution rates of a 100-mg tablet in 4 hours, in 6 hours and in 12 hours are 15 to 45%, 40 to 70%, and not less than 75%, respectively, and those of a 200-mg tablet in 4 hours, in 6 hours and in 12 hours are 15 to 45%, 35 to 65%, and not less than 75%, respectively.

Start the test with 1 tablet of Sodium Valproate Extended-release Tablets A, withdraw exactly 20 mL of the medium at the specified minutes after starting the test and supply exactly 20 mL of water warmed to 37 ± 0.5 °C immediately after withdrawing of the medium every time. Filter the media through a membrane filter with a pore size not exceeding 0.45 µm. Discard not less than 10 mL of the first filtrate, pipet V mL of the subsequent filtrate, add water to make exactly V' mL so that each mL contains about 0.11 mg of sodium valproate (C₈H₁₅NaO₂), and use these solutions as the sample solutions. Separately, weigh accurately about 56 mg of sodium valproate for assay, previously dried at 105°C for 3 hours, and dissolve in water to make exactly 50 mL. Pipet 5 mL of this solution, add water to make exactly 50 mL, and use this solution as the standard solution. Perform the test with exactly 50 μ L each of the sample solutions and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, $A_{T(n)}$ and A_{S} , of valproic acid in each solution.

Dissolution rate (%) with respect to the labeled amount of sodium valproate ($C_8H_{15}NaO_2$) on the *n*th medium withdrawing (n = 1, 2, 3)

$$= M_{\rm S} \times \left\{ \frac{A_{\rm T(n)}}{A_{\rm S}} + \sum_{i=1}^{n-1} \left(\frac{A_{\rm T(i)}}{A_{\rm S}} \times \frac{1}{45} \right) \right\} \times \frac{V'}{V} \times \frac{1}{C} \times 180$$

M_S: Amount (mg) of sodium valproate for assay taken
 C: Labeled amount (mg) of sodium valproate (C₈H₁₅NaO₂) in 1 tablet

Operating conditions—

Proceed as directed in the operating conditions in the Assay.

System suitability—

System performance: When the procedure is run with 50 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of valproic acid are not less than 3000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times

with 50 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of valproic acid is not more than 1.5%.

Assay Weigh accurately the mass of not less than 20 Sodium Valproate Extended-release Tablets A, and powder. Weigh accurately a portion of the powder, equivalent to about 0.1 g of sodium valproate (C₈H₁₅NaO₂), add about 80 mL of the mobile phase, shake thoroughly, add the mobile phase to make exactly 100 mL, and centrifuge. Pipet 20 mL of the supernatant liquid, add exactly 5 mL of the internal standard solution, and use this solution as the sample solution. Separately, weigh accurately about 0.1 g of sodium valproate for assay, previously dried at 105°C for 3 hours, dissolve in the mobile phase to make exactly 100 mL. Pipet 20 mL of this solution, add exactly 5 mL of the internal standard solution, and use this solution as the standard solution. Perform the test with 10 µL each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and calculate the ratios, Q_T and Q_S , of the peak area of valproic acid to that of the internal standard.

Amount (mg) of sodium valproate (
$$C_8H_{15}NaO_2$$
)
= $M_S \times Q_T/Q_S$

 M_S : Amount (mg) of sodium valproate for assay taken

Internal standard solution—A solution of ethyl parahydroxybenzoate in the mobile phase (1 in 50,000).

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 210 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 25°C.

Mobile phase: A mixture of 0.05 mol/L sodium dihydrogen phosphate TS (pH 3.0) and acetonitrile for liquid chromatography (1:1).

Flow rate: Adjust so that the retention time of valproic acid is about 6 minutes.

System suitability-

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, the internal standard and valproic acid are eluted in this order with the resolution between these peaks being not less than 7.

System repeatability: When the test is repeated 6 times with $10 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of valproic acid to that of the internal standard is not more than 1.0%.

Containers and storage Containers—Tight containers.

Add the following:

Sodium Valproate Extended-release Tablets B

バルプロ酸ナトリウム徐放錠 B

Sodium Valproate Extended-release Tablets B contain not less than 95.0% and not more than 105.0% of the labeled amount of sodium valproate $(C_8H_{15}NaO_2: 166.19)$.

Method of preparation Prepare as directed under Tablets, with Sodium Valproate.

Identification To a quantity of the powdered Sodium Valproate Extended-release Tablets B, equivalent to 1.0 g of Sodium Valproate, add 10 mL of water, heat on a water bath for 30 minutes, and filter. To 2.5 mL of the filtrate add 2.5 mL of water and 1 mL of a solution of cobalt (II) nitrate hexahydrate (1 in 20), and heat on a water bath: a purple precipitate is formed.

Uniformity of dosage units <6.02> Perform the Mass variation test, or the Content uniformity test according to the following method: it meets the requirement.

To 1 tablet of Sodium Valproate Extended-release Tablets B add 150 mL of the mobile phase, allow to stand for not less than 16 hours, shake until the film is disintegrated, and add the mobile phase to make exactly 200 mL. Pipet V mL of this solution, and add the mobile phase to make exactly V' mL so that each mL contains about 1 mg of sodium valproate ($C_8H_{15}NaO_2$). Centrifuge this solution, and filter the supernatant liquid through a membrane filter with a pore size not exceeding 0.45 μ m. Discard the first 5 mL of the filtrate, pipet 20 mL of the subsequent filtrate, add exactly 5 mL of the internal standard solution, and use this solution as the sample solution. Then, proceed as directed in the Assay.

Amount (mg) of sodium valproate (
$$C_8H_{15}NaO_2$$
)
= $M_S \times Q_T/Q_S \times V'/V \times 2$

M_S: Amount (mg) of sodium valproate for assay taken

Internal standard solution—A solution of methyl parahydroxybenzoate in the mobile phase (1 in 50,000).

Dissolution <6.10> When the test is performed at 50 revolutions per minute according to the Paddle method, using 900 mL of water as the dissolution medium, the dissolution rates of a 200-mg tablet in 8 hours, in 11 hours and in 20 hours are 15 to 45%, 35 to 65%, and not less than 70%, respectively, and those of a 400-mg tablet in 9 hours, in 12 hours and in 21 hours are 15 to 45%, 35 to 65%, and not less than 70%, respectively.

Start the test with 1 tablet of Sodium Valproate Extended-release Tablets B, withdraw exactly 20 mL of the medium at the specified minutes after starting the test and supply exactly 20 mL of water warmed to 37 ± 0.5 °C immediately after withdrawing of the medium every time.

Filter the media through a membrane filter with a pore size not exceeding $0.45~\mu m$. Discard not less than 2 mL of the first filtrate, pipet V mL of the subsequent filtrate, add water to make exactly V' mL so that each mL contains about 0.22 mg of sodium valproate ($C_8H_{15}NaO_2$), and use these solutions as the sample solutions. Separately, weigh accurately about 55 mg of sodium valproate for assay, previously dried at 105 °C for 3 hours, dissolve in water to make exactly 250 mL, and use this solution as the standard solution. Perform the test with exactly $20~\mu$ L each of the sample solutions and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, $A_{T(n)}$ and A_S , of valproic acid in each solution.

Dissolution rate (%) with respect to the labeled amount of sodium valproate ($C_8H_{15}NaO_2$) on the *n*th medium withdrawing (n = 1, 2, 3)

$$= M_{\rm S} \times \left\{ \frac{A_{\rm T(n)}}{A_{\rm S}} + \sum_{i=1}^{n-1} \left(\frac{A_{\rm T(i)}}{A_{\rm S}} \times \frac{1}{45} \right) \right\} \times \frac{V'}{V} \times \frac{1}{C} \times 360$$

M_S: Amount (mg) of sodium valproate for assay taken
C: Labeled amount (mg) of sodium valproate
(C₈H₁₅NaO₂) in 1 tablet

Operating conditions—

Proceed as directed in the operating conditions in the Assay.

System suitability—

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of valproic acid are not less than 3000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of valproic acid is not more than 1.0%.

Assay To 20 Sodium Valproate Extended-release Tablets B add 150 mL of the mobile phase, allow to stand for not less than 16 hours, shake until the film is disintegrated, and add the mobile phase to make exactly 200 mL. Pipet 5 mL of this solution, add the mobile phase to make exactly $V \, \text{mL}$ so that each mL contains about 1 mg of sodium valproate (C₈H₁₅NaO₂). Centrifuge this solution, and filter the supernatant liquid through a membrane filter with a pore size not exceeding $0.45 \mu m$. Discard the first 5 mL of the filtrate, pipet 20 mL of the subsequent filtrate, add exactly 5 mL of the internal standard solution, and use this solution as the sample solution. Separately, weigh accurately about 0.1 g of sodium valproate for assay, previously dried at 105°C for 3 hours, dissolve in the mobile phase to make exactly 100 mL. Pipet 20 mL of this solution, add exactly 5 mL of the internal standard solution, and use this solution as the standard solution. Perform the test with 10 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and calculate the ratios, Q_T and Q_S , of the peak area of valproic acid to that of the internal standard.

Amount (mg) of sodium valproate ($C_8H_{15}NaO_2$) in 1 tablet = $M_S \times Q_T/Q_S \times V/50$

 M_S : Amount (mg) of sodium valproate for assay taken

Internal standard solution—A solution of methyl parahydroxybenzoate in the mobile phase (1 in 50,000). Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 210 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter)

Column temperature: A constant temperature of about $25\,^{\circ}\mathrm{C}$.

Mobile phase: A mixture of 0.05 mol/L sodium dihydrogen phosphate TS (pH 3.0) and acetonitrile for liquid chromatography (1:1).

Flow rate: Adjust so that the retention time of valproic acid is about 6 minutes.

System suitability—

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, the internal standard and valproic acid are eluted in this order with the resolution between these peaks being not less than 7.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of valproic acid to that of the internal standard is not more than 1.0%.

Containers and storage Containers—Tight containers.

Teicoplanin

テイコプラニン

Change the Purity (4) as follows:

Purity

(4) Arsenic <1.11>—Prepare the test solution with 1.0 g of Teicoplanin according to Method 3, and perform the test (not more than 2 ppm).

Add the following:

Telmisartan and Hydrochlorothiazide Tablets

テルミサルタン・ヒドロクロロチアジド錠

Telmisartan and Hydrochlorothiazide Tablets contain not less than 95.0% and not more than 105.0% of the labeled amount of telmisartan ($C_{33}H_{30}N_4O_2$: 514.62) and hydrochlorothiazide ($C_7H_8ClN_3O_4S_2$:

297.74).

Method of preparation Prepare as directed under Tablets, with Telmisartan and Hydrochlorothiazide.

Identification (1) Perform the test with $5 \mu L$ each of the sample solution and standard solution obtained in the Assay (1) as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions: the retention times of the peaks of telmisartan in the chromatograms obtained from the sample solution and standard solution are the same, and absorption spectra of these peaks exhibit similar intensities of absorption at the same wavelengths.

Operating conditions—

Column, column temperature, mobile phase, and flow rate: Proceed as directed in the operating conditions in the Assay (1).

Detector: A photodiode array detector (wavelength: 270 nm, spectrum range of measurement: 210 – 400 nm).

System suitability—

System performance: Proceed as directed in the system suitability in the Assay (1).

(2) Perform the test with $5 \mu L$ each of the sample solution and standard solution obtained in the Assay (2) as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions: the retention times of the peaks of hydrochlorothiazide in the chromatograms obtained from the sample solution and standard solution are the same, and absorption spectra of these peaks exhibit similar intensities of absorption at the same wavelengths.

Operating conditions—

Column, column temperature, mobile phase, and flow rate: Proceed as directed in the operating conditions in the Assay (1).

Detector: A photodiode array detector (wavelength: 270 nm, spectrum range of measurement: 210 - 400 nm). System suitability—

System performance: Proceed as directed in the system suitability in the Assay (2).

Purity Related substances—To a quantity of powdered Telmisartan and Hydrochlorothiazide Tablets, equivalent to 12.5 mg of Hydrochlorothiazide, add 40 mL of the dissolving solution, disperse by sonicating, add the dissolving solution to make exactly 50 mL. Centrifuge this solution, and use the supernatant liquid as the sample solution. Pipet 1 mL of the sample solution, add the dissolving solution to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly $20 \,\mu\text{L}$ each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine each peak area by the automatic integration method: the area of the peak having the relative retention time of about 0.9 to hydrochlorothiazide, obtained from the sample solution is not larger than the peak area of hydrochlorothiazide from the standard solution.

Dissolving solution: Dissolve 2 g of ammonium dihydrogen phosphate in 1000 mL of water, and adjust to pH 1.8 with phosphoric acid. To 1000 mL of this solution add 1000

mL of acetonitrile.

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 270 nm).

Column: A stainless steel column 4.0 mm in inside diameter and 15 cm in length, packed with octylsilanized silica gel for liquid chromatography (3 μ m in particle diameter).

Column temperature: A constant temperature of about 40°C.

Mobile phase A: Dissolve 2 g of ammonium dihydrogen phosphate in 1000 mL of water, and adjust to pH 3.5 with phosphoric acid.

Mobile phase B: Acetonitrile.

Flowing of mobile phase: Control the gradient by mixing the mobile phases A and B as directed in the following table.

Time after injection of sample (min)	Mobile phase A (vol%)	Mobile phase B (vol%)
0 - 8	$90 \rightarrow 50$	10 → 50
8 – 12	50	50
12 – 18	$50 \rightarrow 20$	$50 \rightarrow 80$
18 – 20	20	80

Flow rate: 1.0 mL per minute.

System suitability—

Test for required detectability: Pipet 5 mL of the standard solution and add the dissolving solution to make exactly 50 mL. Confirm that the peak area of hydrochlorothiazide obtained with 20 μ L of this solution is equivalent to 7 to 13% of that with 20 μ L of the standard solution.

System performance: When the procedure is run with 20 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of hydrochlorothiazide are not less than 6000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with 20 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of hydrochlorothiazide is not more than 2.0%.

Uniformity of dosage units <6.02> Perform the test according to the following method: it meets the requirement of the Content uniformity test.

(1) Telmisartan—To 1 tablet of Telmisartan and Hydrochlorothiazide Tablets add 4V/5 mL of the dissolving solution, disintegrate by sonicating, add the dissolving solution to make exactly VmL so that each mL contains about 1.6 mg of telmisartan (C₃₃H₃₀N₄O₂). Centrifuge this solution, pipet 5 mL of the supernatant liquid, add the buffer solution to make exactly 25 mL, and use this solution as the sample solution. Proceed as directed in the Assay (1).

Amount (mg) of telmisartan (
$$C_{33}H_{30}N_4O_2$$
)
= $M_S \times A_T/A_S \times V/50$

 $M_{\rm S}$: Amount (mg) of telmisartan for assay taken

Dissolving solution: Dissolve 2 g of ammonium dihydro-

gen phosphate in 1000 mL of water, and adjust to pH 1.8 with phosphoric acid. To 1000 mL of this solution add 1000 mL of acetonitrile.

Buffer solution: Dissolve 2 g of ammonium dihydrogen phosphate in 1000 mL of water, and adjust to pH 1.8 with phosphoric acid.

(2) Hydrochlorothiazide—To 1 tablet of Telmisartan and Hydrochlorothiazide Tablets add 4V/5 mL of the dissolving solution, disintegrate by sonicating, add the dissolving solution to make exactly V mL so that each mL contains about 0.25 mg of hydrochlorothiazide (C₇H₈ClN₃O₄S₂). Centrifuge this solution, pipet 5 mL of the supernatant liquid, add the buffer solution to make exactly 25 mL, and use this solution as the sample solution. Proceed as directed in the Assay (2).

Amount (mg) of hydrochlorothiazide ($C_7H_8ClN_3O_4S_2$) $= M_{\rm S} \times A_{\rm T}/A_{\rm S} \times V/50$

M_S: Amount (mg) of Hydrochlorothiazide RS taken

Dissolving solution: Dissolve 2 g of ammonium dihydrogen phosphate in 1000 mL of water, and adjust to pH 1.8 with phosphoric acid. To 1000 mL of this solution add 1000 mL of acetonitrile.

Buffer solution: Dissolve 2 g of ammonium dihydrogen phosphate in 1000 mL of water, and adjust to pH 1.8 with phosphoric acid.

Dissolution <6.10> (1) Telmisartan—When the test is performed at 50 revolutions per minute according to the Paddle method, using 900 mL of 2nd fluid for dissolution test as the dissolution medium, the dissolution rates in 45 minutes of a telmisartan 40-mg and hydrochlorothiazide 12.5-mg tablet and a telmisartan 80-mg and hydrochlorothiazide 12.5-mg tablet are not less than 85% and not less than 80%, respectively.

Start the test with 1 tablet of Telmisartan and Hydrochlorothiazide Tablets, withdraw not less than 20 mL of the medium at the specified minute after starting the test, and filter through a membrane filter with a pore size not exceeding $0.45 \,\mu\text{m}$. Discard not less than 15 mL of the first filtrate, pipet VmL of the subsequent filtrate, add the dissolution medium to make exactly V'mL so that each mL contains about 44 μ g of telmisartan (C₃₃H₃₀N₄O₂), and use this solution as the sample solution. Separately, weigh accurately about 44 mg of telmisartan for assay, previously dried at 105°C for 4 hours, dissolve in 10 mL of a solution of meglumine in methanol (1 in 250), and add methanol to make exactly 50 mL. Pipet 5 mL of this solution, add water to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 25 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, A_T and A_S , of telmisartan in each solution.

Dissolution rate (%) with respect to the labeled amount of telmisartan (C₃₃H₃₀N₄O₂)

$$= M_{\rm S} \times A_{\rm T}/A_{\rm S} \times V'/V \times 1/C \times 90$$

 M_S : Amount (mg) of telmisartan for assay taken C: Labeled amount (mg) of telmisartan ($C_{33}H_{30}N_4O_2$) in 1 tablet

Operating conditions—

Proceed as directed in the operating conditions in the Assay (1).

System suitability—

System performance: When the procedure is run with 25 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of telmisartan are not less than 25,000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with 25 μ L of the standard solution under the above operating conditions, the relative standard deviation of the peak area of telmisartan is not more than 2.0%.

(2) Hydrochlorothiazide—When the test is performed at 75 revolutions per minute according to the Paddle method, using 900 mL of 2nd fluid for dissolution test as the dissolution medium, the dissolution rate in 45 minutes of Telmisartan and Hydrochlorothiazide Tablets is not less than 80%.

Start the test with 1 tablet of Telmisartan and Hydrochlorothiazide Tablets, withdraw not less than 20 mL of the medium at the specified minute after starting the test, and filter through a membrane filter with a pore size not exceeding $0.45 \,\mu\text{m}$. Discard not less than 15 mL of the first filtrate, pipet VmL of the subsequent filtrate, add the dissolution medium to make exactly V' mL so that each mL contains about 14 µg of hydrochlorothiazide (C₇H₈ClN₃O₄S₂), and use this solution as the sample solution. Separately, weigh accurately about 14 mg of Hydrochlorothiazide RS, previously dried at 105°C for 2 hours, and dissolve in methanol to make exactly 50 mL. Pipet 5 mL of this solution, add water to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly 25 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, $A_{\rm T}$ and $A_{\rm S}$, of hydrochlorothiazide in each solution.

Dissolution rate (%) with respect to the labeled amount of hydrochlorothiazide ($C_7H_8ClN_3O_4S_2$)

$$= M_{\rm S} \times A_{\rm T}/A_{\rm S} \times V'/V \times 1/C \times 90$$

M_S: Amount (mg) of Hydrochlorothiazide RS taken C: Labeled amount (mg) of hydrochlorothiazide (C₇H₈ClN₃O₄S₂) in 1 tablet

Operating conditions—

Proceed as directed in the operating conditions in the Assay (1).

System suitability—

System performance: When the procedure is run with 25 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of hydrochlorothiazide are not less than 1000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with 25 μ L of the standard solution under the above operat-

ing conditions, the relative standard deviation of the peak area of hydrochlorothiazide is not more than 2.0%.

Assay (1) Telmisartan—Weigh accurately the mass of not less than 20 Telmisartan and Hydrochlorothiazide Tablets, and powder. Weigh accurately a portion of the powder, equivalent to about 80 mg of telmisartan (C₃₃H₃₀N₄O₂), add 40 mL of the dissolving solution, sonicate, and add the dissolving solution to make exactly 50 mL. Centrifuge this solution, pipet 5 mL of the supernatant liquid, add the buffer solution to make exactly 25 mL, and use this solution as the sample solution. Separately, weigh accurately about 80 mg of telmisartan for assay, previously dried at 105°C for 4 hours, and dissolve in the dissolving solution to make exactly 50 mL. Pipet 5 mL of this solution, add the buffer solution to make exactly 25 mL, and use this solution as the standard solution. Perform the test with exactly 5 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, $A_{\rm T}$ and $A_{\rm S}$, of telmisartan in each solution.

Amount (mg) of telmisartan
$$(C_{33}H_{30}N_4O_2)$$

= $M_S \times A_T/A_S$

 $M_{\rm S}$: Amount (mg) of telmisartan for assay taken

Dissolving solution: Dissolve 2 g of ammonium dihydrogen phosphate in 1000 mL of water, and adjust to pH 1.8 with phosphoric acid. To 1000 mL of this solution add 1000 mL of acetonitrile.

Buffer solution: Dissolve 2 g of ammonium dihydrogen phosphate in 1000 mL of water, and adjust to pH 1.8 with phosphoric acid.

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 270 nm).

Column: A stainless steel column 3.0 mm in inside diameter and 7.5 cm in length, packed with octylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 40°C.

Mobile phase A: Dissolve 2 g of ammonium dihydrogen phosphate in 1000 mL of water, and adjust to pH 3.5 with phosphoric acid.

Mobile phase B: Acetonitrile.

Flowing of mobile phase: Control the gradient by mixing the mobile phases A and B as directed in the following table.

Time after injection of sample (min)	Mobile phase A (vol%)	Mobile phase B (vol%)
0 - 2	90	10
2 – 7	$90 \rightarrow 20$	$10 \rightarrow 80$
7 - 8	20	80

Flow rate: 0.8 mL per minute.

System suitability-

System performance: When the procedure is run with 5

 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of telmisartan are not less than 15,000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with $5 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of telmisartan is not more than 1.0%.

(2) Hydrochlorothiazide—Weigh accurately the mass of not less than 20 Telmisartan and Hydrochlorothiazide Tablets, and powder. Weigh accurately a portion of the powder, equivalent to about 12.5 mg of hydrochlorothiazide (C₇H₈ClN₃O₄S₂), add 40 mL of the dissolving solution, sonicate, and add the dissolving solution to make exactly 50 mL. Centrifuge this solution, pipet 5 mL of the supernatant liquid, add the buffer solution to make exactly 25 mL, and use this solution as the sample solution. Separately, weigh accurately about 12.5 mg of Hydrochlorothiazide RS, previously dried at 105°C for 2 hours, and dissolve in the dissolving solution to make exactly 50 mL. Pipet 5 mL of this solution, add the buffer solution to make exactly 25 mL, and use this solution as the standard solution. Perform the test with exactly 5 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, A_T and A_S , of hydrochlorothiazide in each solution.

Amount (mg) of hydrochlorothiazide ($C_7H_8ClN_3O_4S_2$) = $M_S \times A_T/A_S$

M_S: Amount (mg) of Hydrochlorothiazide RS taken

Dissolving solution: Dissolve 2 g of ammonium dihydrogen phosphate in 1000 mL of water, and adjust to pH 1.8 with phosphoric acid. To 1000 mL of this solution add 1000 mL of acetonitrile.

Buffer solution: Dissolve 2 g of ammonium dihydrogen phosphate in 1000 mL of water, and adjust to pH 1.8 with phosphoric acid.

Operating conditions—

Proceed as directed in the operating conditions in (1). System suitability—

System performance: When the procedure is run with 5 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of hydrochlorothiazide are not less than 1500 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with $5 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of hydrochlorothiazide is not more than 1.0%.

Containers and storage Containers—Tight containers.

Testosterone Enanthate

テストステロンエナント酸エステル

Change the Description and Optical rotation as follows:

Description Testosterone Enanthate occurs as white to pale yellow, crystals or crystalline powder, or a pale yellow-brown viscous liquid. It is odorless or has a slight, characteristic odor.

It is freely soluble in ethanol (99.5), and practically insoluble in water.

Melting point: about 36°C

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁵: +76 - +86° (after drying, 0.1 g, ethanol (99.5), 10 mL, 100 mm).

Delete the following Monograph:

Freeze-dried Tetanus Antitoxin, Equine

乾燥破傷風ウマ抗毒素

Tipepidine Hibenzate Tablets

チペピジンヒベンズ酸塩錠

Change the Identification (2) as follows:

Identification

(2) To a quantity of powdered Tipepidine Hibenzate Tablets, equivalent to 11 mg of Tipepidine Hibenzate, add 30 mL of ethanol (99.5), and warm for 10 minutes with occasional shaking. After cooling, add ethanol (99.5) to make 50 mL, and filter. To 1 mL of the filtrate add ethanol (99.5) to make 20 mL, and determine the absorption spectrum of this solution as directed under Ultraviolet-visible Spectrophotometry <2.24>: it exhibits a maximum between 280 nm and 286 nm.

Triamcinolone Acetonide

トリアムシノロンアセトニド

Change the Description and Optical rotation as follows:

Description Triamcinolone Acetonide occurs as a white crystalline powder.

It is sparingly soluble in ethanol (99.5) and in acetone, slightly soluble in methanol, and practically insoluble in water.

Melting point: about 290°C (with decomposition).

It shows crystal polymorphism.

Optical rotation $\langle 2.49 \rangle$ [α]_D²⁵: +110 - +120° (after drying, 0.1 g, ethanol (99.5), 10 mL, 100 mm).

Ursodeoxycholic Acid

ウルソデオキシコール酸

Change the Melting point as follows:

Melting point <2.60> 201 – 205°C

Add the following:

Valsartan and Hydrochlorothiazide Tablets

バルサルタン·ヒドロクロロチアジド錠

Valsartan and Hydrochlorothiazide Tablets contain not less than 95.0% and not more than 105.0% of the labeled amount of valsartan ($C_{24}H_{29}N_5O_3$: 435.52) and hydrochlorothiazide ($C_7H_8ClN_3O_4S_2$: 297.74).

Method of preparation Prepare as directed under Tablets, with Valsartan and Hydrochlorothiazide.

Identification (1) To a quantity of powdered Valsartan and Hydrochlorothiazide Tablets, equivalent to 80 mg of Valsartan, add 5 mL of acetone, shake, centrifuge, and use the supernatant liquid as the sample solution. Separately, dissolve 16 mg of valsartan in 1 mL of acetone, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 5 μ L each of the sample solution and standard solution on a plate of silica gel with fluorescent indicator for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, hexane and acetic acid (100) (15:5:2) to a distance of about 10 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 254 nm): one of the two spots obtained from the sample solution and the spot from the standard solution show the same Rf value.

(2) To a quantity of powdered Valsartan and Hydrochlorothiazide Tablets, equivalent to 6.25 mg of Hydrochlorothiazide, add 5 mL of acetone, shake, centrifuge, and use the supernatant liquid as the sample solution. Separately, dissolve 12.5 mg of hydrochlorothiazide in 10 mL of acetone, and use this solution as the standard solution. Perform the test with these solutions as directed under Thinlayer Chromatography $\langle 2.03 \rangle$. Spot 5 μ L each of the sample solution and standard solution on a plate of silica gel with fluorescent indicator for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, hexane and acetic acid (100) (15:5:2) to a distance of about 10 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 254 nm): one of the two spots obtained from

the sample solution and the spot from the standard solution show the same Rf value.

Uniformity of dosage units <6.02> (1) Valsartan—Perform the Mass variation test, or the Content uniformity test according to the following method: it meets the requirement.

To 1 tablet of Valsartan and Hydrochlorothiazide Tablets add 10 mL of water, and shake until the tablet is disintegrated. Add 10 mL of acetonitrile, shake thoroughly, and add a mixture of water and acetonitrile (1:1) to make exactly 50 mL. Centrifuge this solution, pipet V mL of the supernatant liquid, add a mixture of water and acetonitrile (1:1) to make exactly V' mL so that each mL contains about 0.4 mg of valsartan ($C_{24}H_{29}N_5O_3$), and use this solution as the sample solution. Proceed as directed in the Assay (1).

Amount (mg) of valsartan (
$$C_{24}H_{29}N_5O_3$$
)
= $M_S \times A_T/A_S \times V'/V \times 1/2$

 M_S : Amount (mg) of Valsartan RS taken, calculated on the anhydrous and residual solvent-free basis

(2) Hydrochlorothiazide—Perform the test according to the following method: it meets the requirement of the Content uniformity test.

To 1 tablet of Valsartan and Hydrochlorothiazide Tablets add 10 mL of water, and shake until the tablet is disintegrated. Add 10 mL of acetonitrile, shake thoroughly, and add a mixture of water and acetonitrile (1:1) to make exactly 50 mL. Centrifuge this solution, pipet V mL of the supernatant liquid, add a mixture of water and acetonitrile (1:1) to make exactly V' mL so that each mL contains about 31 μ g of hydrochlorothiazide ($C_7H_8ClN_3O_4S_2$), and use this solution as the sample solution. Proceed as directed in the Assay (2).

Amount (mg) of hydrochlorothiazide ($C_7H_8ClN_3O_4S_2$) = $M_S \times A_T/A_S \times V'/V \times 1/8$

 $M_{\rm S}$: Amount (mg) of Hydrochlorothiazide RS taken

Dissolution <6.10> (1) Valsartan—When the test is performed at 50 revolutions per minute according to the Paddle method, using 900 mL of water as the dissolution medium, the dissolution rate in 30 minutes of Valsartan and Hydrochlorothiazide Tablets is not less than 75%.

Start the test with 1 tablet of Valsartan and Hydrochlorothiazide Tablets, withdraw not less than 20 mL of the medium at the specified minute after starting the test, and filter through a membrane filter with a pore size not exceeding 0.45 μ m. Discard not less than 5 mL of the first filtrate, pipet V mL of the subsequent filtrate, and add water to make exactly V' mL so that each mL contains about 89 μ g of valsartan ($C_{24}H_{29}N_5O_3$). Pipet 5 mL of this solution, add methanol to make exactly 10 mL, and use this solution as the sample solution. Separately, weigh accurately about 45 mg of Valsartan RS (separately determine the water <2.48> and the residual solvent in the same manner as Valsartan), and dissolve in methanol to make exactly 50 mL. Pipet 10 mL of this solution, add exactly 100 mL of water, then add

methanol to make exactly 200 mL, and use this solution as the standard solution. Perform the test with exactly $10 \,\mu\text{L}$ each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and determine the peak areas, A_{T} and A_{S} , of valsartan in each solution.

Dissolution rate (%) with respect to the labeled amount of valsartan ($C_{24}H_{29}N_3O_3$)

$$= M_{\rm S} \times A_{\rm T}/A_{\rm S} \times V'/V \times 1/C \times 180$$

 M_S : Amount (mg) of Valsartan RS taken, calculated on the anhydrous and residual solvent-free basis

C: Labeled amount (mg) of valsartan ($C_{24}H_{29}N_5O_3$) in 1 tablet

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 225 nm).

Column: A stainless steel column 3.0 mm in inside diameter and 12.5 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 25°C.

Mobile phase: Dissolve 14.68 g of disodium hydrogen phosphate dodecahydrate and 3.81 g of potassium dihydrogen phosphate in 1000 mL of water. To 4 volumes of this solution add 1 volume of acetonitrile.

Flow rate: Adjust so that the retention time of valsartan is about 6 minutes.

System suitability—

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of valsartan are not less than 500, and not less than 0.7 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with $10 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of valsartan is not more than 1.0%.

(2) Hydrochlorothiazide—When the test is performed at 50 revolutions per minute according to the Paddle method, using 900 mL of water as the dissolution medium, the dissolution rate in 15 minutes of Valsartan and Hydrochlorothiazide Tablets is not less than 85%.

Start the test with 1 tablet of Valsartan and Hydrochlorothiazide Tablets, withdraw not less than 20 mL of the medium at the specified minute after starting the test, and filter through a membrane filter with a pore size not exceeding 0.45 μ m. Discard the first 5 mL of the filtrate, pipet V mL of the subsequent filtrate, add water to make exactly V' mL so that each mL contains about 6.9 μ g of hydrochlorothiazide ($C_7H_8CIN_3O_4S_2$), and use this solution as the sample solution. Separately, weigh accurately about 14 mg of Hydrochlorothiazide RS, previously dried at 105 °C for 2 hours, and dissolve in methanol to make exactly 100 mL. Pipet 5 mL of this solution, add water to make exactly 100 mL, and use this solution as the standard solution. Perform

the test with exactly $10 \,\mu\text{L}$ each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and determine the peak areas, A_{T} and A_{S} , of hydrochlorothiazide in each solution.

Dissolution rate (%) with respect to the labeled amount of hydrochlorothiazide ($C_7H_8CIN_3O_4S_2$)

$$= M_{\rm S} \times A_{\rm T}/A_{\rm S} \times V'/V \times 1/C \times 45$$

M_S: Amount (mg) of Hydrochlorothiazide RS taken
 C: Labeled amount (mg) of hydrochlorothiazide (C₇H₈ClN₃O₄S₂) in 1 tablet

Operating conditions—

Proceed as directed in the operating conditions in (1). System suitability—

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of hydrochlorothiazide are not less than 3000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of hydrochlorothiazide is not more than 1.0%.

Assay (1) Valsartan—Weigh accurately the mass of not less than 20 tablets of Valsartan and Hydrochlorothiazide Tablets, and powder. Weigh accurately a portion of the powder, equivalent to about 80 mg of valsartan $(C_{24}H_{29}N_5O_3)$, add 10 mL of water, and shake. Add 10 mL of acetonitrile, shake thoroughly, add a mixture of water and acetonitrile (1:1) to make exactly 50 mL, and centrifuge. Pipet 5 mL of the supernatant liquid, add a mixture of water and acetonitrile (1:1) to make exactly 20 mL, and use this solution as the sample solution. Separately, weigh accurately about 40 mg of Valsartan RS (separately determine the water <2.48> and the residual solvent in the same manner as Valsartan), dissolve in a mixture of water and acetonitrile (1:1) to make exactly 25 mL, and use this solution as the valsartan standard stock solution. Pipet 5 mL of the valsartan standard stock solution, add a mixture of water and acetonitrile (1:1) to make exactly 20 mL, and use this solution as the standard solution. Perform the test with exactly $10 \,\mu L$ each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, $A_{\rm T}$ and $A_{\rm S}$, of valsartan in each solution.

Amount (mg) of valsartan (
$$C_{24}H_{29}N_5O_3$$
)
= $M_S \times A_T/A_S \times 2$

 M_S : Amount (mg) of Valsartan RS taken, calculated on the anhydrous and residual solvent-free basis

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 271 nm).

Column: A stainless steel column 3.0 mm in inside diameter and 12.5 cm in length, packed with octadecylsilanized

silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 25°C.

Mobile phase A: A mixture of water, acetonitrile and trifluoroacetic acid (900:100:1).

Mobile phase B: A mixture of acetonitrile, water and trifluoroacetic acid (900:100:1).

Flowing of mobile phase: Control the gradient by mixing the mobile phases A and B as directed in the following table.

Time after injection of sample (min)	Mobile phase A (vol%)	Mobile phase B (vol%)
0 - 25	$90 \rightarrow 10$	10 → 90

Flow rate: Adjust so that the retention time of valsartan is about 16 minutes.

System suitability—

System performance: Dissolve 1 mg of 4-amino-6-chlorobenzene-1,3-disulfonamide in a mixture of water and acetonitrile (1:1) to make 200 mL. To 1 mL of this solution, 5 mL of the valsartan standard stock solution and 5 mL of the hydrochlorothiazide standard stock solution in (2) add a mixture of water and acetonitrile (1:1) to make 20 mL. When the procedure is run with $10\,\mu\text{L}$ of this solution under the above operating conditions, 4-amino-6-chlorobenzene-1,3-disulfonamide, hydrochlorothiazide and valsartan are eluted in this order with the resolution between the peaks of 4-amino-6-chlorobenzene-1,3-disulfonamide and hydrochlorothiazide being not less than 1.5.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of valsartan is not more than 1.0%.

(2) Hydrochlorothiazide—Weigh accurately the mass of not less than 20 tablets of Valsartan and Hydrochlorothiazide Tablets, and powder. Weigh accurately a portion of the powder, equivalent to about 6.25 mg of hydrochlorothiazide (C₇H₈ClN₃O₄S₂), add 10 mL of water, and shake. Add 10 mL of acetonitrile, shake thoroughly, add a mixture of water and acetonitrile (1:1) to make exactly 50 mL, and centrifuge. Pipet 5 mL of the supernatant liquid, add a mixture of water and acetonitrile (1:1) to make exactly 20 mL, and use this solution as the sample solution. Separately, weigh accurately about 12.5 mg of Hydrochlorothiazide RS, previously dried at 105°C for 2 hours, and dissolve in a mixture of water and acetonitrile (1:1) to make exactly 50 mL, and use this solution as the hydrochlorothiazide standard stock solution. Pipet 2.5 mL of the hydrochlorothiazide standard stock solution, add a mixture of water and acetonitrile (1:1) to make exactly 20 mL, and use this solution as the standard solution. Perform the test with exactly $10 \,\mu L$ each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, $A_{\rm T}$ and $A_{\rm S}$, of hydrochlorothiazide in each solution.

Amount (mg) of hydrochlorothiazide ($C_7H_8ClN_3O_4S_2$) = $M_S \times A_T/A_S \times 1/2$

 $M_{\rm S}$: Amount (mg) of Hydrochlorothiazide RS taken

Operating conditions—

Proceed as directed in the operating conditions in (1). System suitability—

System performance: Dissolve 1 mg of 4-amino-6-chlorobenzene-1,3-disulfonamide in a mixture of water and acetonitrile (1:1) to make 200 mL. To 1 mL of this solution, 5 mL of the valsartan standard stock solution in (1) and 5 mL of the hydrochlorothiazide standard stock solution add a mixture of water and acetonitrile (1:1) to make 20 mL. When the procedure is run with $10 \,\mu$ L of this solution under the above operating conditions, 4-amino-6-chlorobenzene-1,3-disulfonamide, hydrochlorothiazide and valsartan are eluted in this order with the resolution between the peaks of 4-amino-6-chlorobenzene-1,3-disulfonamide and hydrochlorothiazide being not less than 1.5.

System repeatability: When the test is repeated 6 times with $10 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of hydrochlorothiazide is not more than 1.0%.

Add the following:

Verapamil Hydrochloride Injection

ベラパミル塩酸塩注射液

Verapamil Hydrochloride Injection is an aqueous injection.

It contains not less than 93.0% and not more than 107.0% of the labeled amount of verapamil hydrochloride ($C_{27}H_{38}N_2O_4$.HCl: 491.06).

Method of preparation Prepare as directed under Injections, with Verapamil Hydrochloride.

Description Verapamil Hydrochloride Injection is a clear, colorless liquid.

Identification To 1 mL of the sample solution obtained in the Assay, add 0.02 mol/L hydrochloric acid TS to make 50 mL, and determine the absorption spectrum of this solution as directed under Ultraviolet-visible Spectrophotometry <2.24>: it exhibits maxima between 227 nm and 231 nm, and between 276 nm and 280 nm.

pH Being specified separately when the drug is granted approval based on the Law.

Bacterial endotoxins <4.01> Less than 12 EU/mg.

Extractable volume <6.05> It meets the requirement.

Foreign insoluble matter <6.06> Perform the test according to Method 1: it meets the requirement.

Insoluble particulate matter <6.07> It meets the requirement.

Sterility < 4.06> Perform the test according to the Membrane filtration method: it meets the requirement.

Assay Pipet a volume of Verapamil Hydrochloride Injection, equivalent to about 10 mg of verapamil hydrochloride ($C_{27}H_{38}N_2O_4$.HCl), add 0.02 mol/L hydrochloric acid TS to make exactly 10 mL, and use this solution as the sample solution. Separately, weigh accurately about 50 mg of verapamil hydrochloride for assay, previously dried at 105°C for 2 hours, dissolve in 0.02 mol/L hydrochloric acid TS to make exactly 50 mL, and use this solution as the standard solution. Perform the test with exactly 10 μ L each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and determine the peak areas, A_T and A_S , of verapamil in each solution.

```
Amount (mg) of verapamil hydrochloride

(C_{27}H_{38}N_2O_4.HCl)

= M_S \times A_T/A_S \times 1/5
```

M_S: Amount (mg) of verapamil hydrochloride for assay

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 279 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 40°C.

Mobile phase: A mixture of methanol, water and perchloric acid (550:450:1).

Flow rate: Adjust so that the retention time of verapamil is about 5 minutes.

System Suitability-

System performance: When the procedure is run with $10 \mu L$ of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of verapamil are not less than 2000 and not more than 2.0, respectively.

System repeatability: When the test is repeated 6 times with $10 \,\mu\text{L}$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of verapamil is not more than 1.0%.

Containers and storage Containers—Hermetic containers, and colored containers may be used.

Storage—Light-resistant.

Delete the following Monograph:

Compound Vitamin B Powder

複方ビタミンB散

Crude Drugs and Related Drugs

Amomum Seed

シュクシャ

Change the Origin/limits of content as follows:

Amomum Seed is the seed mass of Amomum villosum Loureiro var. xanthioides T. L. Wu & S. J. Chen, Amomum villosum Loureiro var. villosum or Amomum longiligulare T. L. Wu (Zingiberaceae).

Artemisia Capillaris Flower

インチンコウ

Change the Description as follows:

Description Capitulum, of ovoid to spherical, about 1.5 – 2 mm in length, about 2 mm in diameter, with linear leaves and pedicels. Outer surface of capitulum, light green to light yellow-brown in color; outer surface of leaf, green to green-brown; outer surface of pedicel, green-brown to dark brown. Under a magnifying glass, at the capitulum, involucral scale in 3 – 4 succubous rows; outer scale, of ovate with obtuse; inner scale, of elliptical, 1.5 mm in length, longer than outer one, with keel midrib and thin membranous margin. Floret, tubular; marginal flower, of female; disk flower, of hermaphrodite. Achene, of obovoid, 0.8 mm in length. Light in texture.

Odor, characteristic, slight; taste, slightly acrid, which gives slightly numbing sensation to the tongue.

Belladonna Root

ベラドンナコン

Change the Description as follows:

Description Cylindrical root, usually 10 - 30 cm in length, 0.5 - 4 cm in diameter; often cut crosswise or lengthwise; externally grayish brown to grayish yellow-brown, with longitudinal wrinkles; periderm often removed; fractured surface is light yellow to light yellow-brown in color and is powdery.

Almost odorless.

Bitter Tincture

苦味チンキ

Change the Identification (2) as follows:

Identification

(2) Use Bitter Tincture as the sample solution. Separately, to 0.5 g of pulverized Bitter Orange Peel add 10 mL of methanol, shake for 5 minutes, filter, and use the filtrate as the standard solution (1). Proceed with 0.5 g each of pulverized Swertia Herb and Japanese Zanthoxylum Peel in the same manner, and use the solutions so obtained as the standard solution (2) and the standard solution (3). Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L of the sample solution and 5 μ L each of standard solutions (1), (2) and (3) on a plate of silica gel with fluorescent indicator for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, ethanol (99.5) and water (8:2:1) to a distance of about 7 cm, and air-dry the plate. Examine the plate under ultraviolet light (main wavelength: 254 nm): one of the several spots obtained from the sample solution has the same Rf value with the clear spot at an Rf value of about 0.7 from the standard solution (3). Spray evenly vanillinsulfuric acid-ethanol TS for spraying on the plate, and heat the plate at 105°C for 5 minutes: two of the several spots from the sample solution show the same color tone and Rfwith the clear spot at an Rf value of about 0.4 from the standard solution (1), and the clear spot at an Rf value of about 0.35 from the standard solution (2).

Bofutsushosan Extract

防風通聖散エキス

Change the Identification (1) as follows:

Identification (1) To 2.0 g of the dry extract (or 6.0 g of the viscous extract) add 10 mL of water, shake, then add 10 mL of diethyl ether, shake, and centrifuge. Separate the diethyl ether layer, add 10 mL of sodium hydroxide TS, shake, centrifuge, separate the diethyl ether layer, and use this solution as the sample solution. Separately, use (*Z*)-ligustilide TS for thin-layer chromatography as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 20 μ L of the sample solution and 10 μ L of the standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of butyl acetate and hexane (2:1) to a distance of about 7 cm, and air-dry the plate. Examine

under ultraviolet light (main wavelength: 365 nm): one of the several spots obtained from the sample solution has the same color tone and Rf value with the blue-white fluorescent spot from the standard solution (Japanese Angelica Root; Cnidium Rhizome).

Boiogito Extract

防已黄耆湯エキス

Change the Assay (1) as follows:

Assay (1) Sinomenine—Weigh accurately about 0.5 g of the dry extract (or an amount of the viscous extract, equivalent to about 0.5 g of the dried substance), add 20 mL of diethyl ether, shake, then add 5.0 mL of 0.1 mol/L hydrochloric acid TS, and shake for 10 minutes, centrifuge, and remove the upper layer. Add 20 mL of diethyl ether, proceed in the same manner as described above, and remove the upper layer. To the aqueous layer add 5.0 mL of diluted sodium hydroxide TS (1 in 10) and 10 mL of methanol, shake for 15 minutes, centrifuge, and take the supernatant liquid. To the residue add 20 mL of diluted methanol (1 in 2), shake for 15 minutes, centrifuge, and take the supernatant liquid. Combine all the supernatant liquids, add diluted methanol (1 in 2) to make exactly 50 mL, and use this solution as the sample solution. Separately, weigh accurately about 5 mg of sinomenine for assay, dissolve in diluted methanol (1 in 2) to make exactly 100 mL, and use this solution as the standard solution. Perform the test with exactly $10 \,\mu\text{L}$ each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, $A_{\rm T}$ and $A_{\rm S}$, of sinomenine in each solution.

Amount (mg) of sinomenine = $M_S \times A_T/A_S \times 1/2$

M_S: Amount (mg) of sinomenine for assay taken

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 254 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 30°C.

Mobile phase: To 3 g of sodium lauryl sulfate add 350 mL of acetonitrile, shake, then add 650 mL of water and 1 mL of phosphoric acid to dissolve lauryl sulfate.

Flow rate: 1.0 mL per minute (the retention time of sinomenine is about 18 minutes).

System suitability—

System performance: When the procedure is run with 10 μ L each of the sample solution, the standard solution of sinomenine and the standard solution of glycyrrhizic acid obtained in Assay (2) under the above operating conditions, peaks of sinomenine and glycyrrhizic acid are observed in

the sample solution, glycyrrhizic acid and sinomenine are eluted in this order with the resolution between these peaks being not less than 4.5. Furthermore, except for the peak of glycyrrhizic acid, distinct peaks are observed before and after the peak of sinomenine, and the resolutions between sinomenine and these peaks are respectively not less than 1.5.

System repeatability: When the test is repeated 6 times with $10 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of sinomenine is not more than 1.5%.

Calumba

コロンボ

Change the Purity (1) as follows:

Purity (1) Heavy metals <1.07>—Proceed with 2.0 g of pulverized Calumba according to Method 3, and perform the test. Prepare the control solution with 3.0 mL of Standard Lead Solution (not more than 15 ppm).

Powdered Calumba

コロンボ末

Change the Purity (1) as follows:

Purity (1) Heavy metals <1.07>—Proceed with 2.0 g of Powdered Calumba according to Method 3, and perform the test. Prepare the control solution with 3.0 mL of Standard Lead Solution (not more than 15 ppm).

Cassia Seed

ケツメイシ

Change the Identification as follows:

Identification To 1.0 g of pulverized Cassia Seed add 10 mL of diluted methanol (4 in 5), heat on a water bath for 5 minutes, and filter. Evaporate the solvent of the filtrate, dissolve the residue in 5 mL of water, add 2 mL of ethyl acetate, and shake for 10 minutes. Centrifuge this solution, and use the upper layer as the sample solution. Perform the test with the sample solution as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot $10 \,\mu$ L of the sample solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of diethyl ether, cyclohexane and formic acid (5:5:1) to a distance of about 7 cm, and airdry the plate. Spray evenly potassium hydroxide-ethanol TS on the plate: an orange to yellow-brown spot appears at an Rf value of about 0.35.

Chrysanthemum Flower

キクカ

Change the Description as follows:

Description 1) Chrysanthemum morifolium origin—Capitulum, 15 – 40 mm in diameter; involucre, consisting of 3 to 4 rows of involucral scales, often with peduncle; the outer involucral scale, linear to lanceolate; inner involucral scale, narrow ovate to ovate; ligulate flowers, numerous, white to yellow in color; tubular flowers, small in number, light yellow-brown, occasionally degenerate; outer surface of involucre, green-brown to brown; light in texture and easy to break.

Odor, characteristic; taste, slightly bitter.

2) Chrysanthemum indicum origin—Capitulum, 3 – 10 mm in diameter; involucre, consisting of 3 to 5 rows of involucral scales, often with peduncle; the outer involucral scale, linear to lanceolate; inner involucral scale, narrow ovate to ovate; ligulate flower, in a single circle, yellow to light yellow-brown in color; tubular flowers, numerous, light yellow-brown; outer surface of involucre, yellow-brown to brown; light in texture and easy to break.

Odor, characteristic; taste, slightly bitter.

Eucalyptus Oil

ユーカリ油

Change the Assay as follows:

Assay Weigh accurately about 0.1 g each of Eucalyptus Oil and cineol for assay, and dissolve each in hexane to make exactly 25 mL. Pipet 5 mL each of these solutions, add exactly 5 mL of the internal standard solution to each, then add hexane to make 100 mL, and use these solutions as the sample solution and the standard solution, respectively. Perform the test with $2 \mu L$ each of the sample solution and standard solution as directed under Gas Chromatography $\langle 2.02 \rangle$ according to the following conditions. Calculate the ratios, Q_T and Q_S , of the peak area of cineol to that of the internal standard.

Amount (mg) of cineol ($C_{10}H_{18}O$) = $M_S \times Q_T/Q_S$

 $M_{\rm S}$: Amount (mg) of cineol for assay taken

Internal standard solution—A solution of anisol in hexane (1 in 250).

Operating conditions—

Detector: A hydrogen flame-ionization detector.

Column: A glass column 3 mm in inside diameter and 5 m in length, packed with silanized porous silica gel for gas chromatography coated in 5% with alkylene glycol phthalate ester for gas chromatography (150 to 180 μ m in particle diameter).

Column temperature: A constant temperature of about

120°C.

Carrier gas: Nitrogen.

Flow rate: Adjust so that the retention time of cineol is about 11 minutes.

System suitability—

System performance: Dissolve 0.1 g each of cineol for assay and limonene in 25 mL of hexane. To 1 mL of this solution add hexane to make 20 mL. When the procedure is run with $2 \mu L$ of this solution under the above operating conditions, limonene and cineol are eluted in this order with the resolution between these peaks being not less than 1.5.

System repeatability: When the test is repeated 6 times with $2 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the ratio of the peak area of cineol to that of the internal standard is not more than 1.0%.

Gastrodia Tuber

テンマ

Change the origin/limits of content as follows:

Gastrodia Tuber is the tuber of *Gastrodia elata* Blume (*Orchidaceae*), after being passed through hot water or steamed.

Glycyrrhiza Extract

カンゾウエキス

Change the Method of preparation as follows:

Method of preparation

- 1) To 1 kg of fine cuttings of Glycyrrhiza or the root and stolon of *Glycyrrhiza glabra* Linné (*Leguminosae*) which meets the requirement of Glycyrrhiza add 5 L of Water, Purified Water or Purified Water in Containers, and macerate for 2 days. Filter the macerated solution through a cloth filter. Add 3 L of Water, Purified Water or Purified Water in Containers to the residue, macerate again for 12 hours, and filter through a cloth filter. Evaporate the combined filtrates until the whole volume becomes 3 L. After cooling, add 1 L of Ethanol, and allow to stand in a cold place for 2 days. Filter, and evaporate the filtrate to a viscous extract.
- 2) Take Glycyrrhiza or the root and stolon of Glycyrrhiza glabra Linné (Leguminosae) which meets the requirement of Glycyrrhiza, pulverized to suitable sizes, and prepare the viscous extract as directed under Extracts using Water, Purified Water or Purified Water in Containers as the solvent. Immediately before making a millet jelly-like consistency for the viscous extract, add Ethanol, Anhydrous Ethanol or ethanol (99.5) to the extract, allow it to stand in a cold place, filter, and concentrate the filtrate to prepare.

Crude Glycyrrhiza Extract

カンゾウ粗エキス

Change the Method of preparation as follows:

Method of preparation Take Glycyrrhiza or the root and stolon of *Glycyrrhiza glabra* Linné (*Leguminosae*) which meets the requirement of Glycyrrhiza, pulverized to suitable sizes, and prepare the dry extracts as directed under Extracts using Water, Purified Water or Purified Water in Containers as the solvent.

Add the following:

Goshuyuto Extract

呉茱萸湯エキス

Goshuyuto Extract contains not less than 0.3 mg (for preparation prescribed 3 g of Euodia Fruit) or not less than 0.4 mg (for preparation prescribed 4 g of Euodia Fruit) of evodiamine, and not less than 0.5 mg and not more than 2.0 mg (for preparation prescribed 1 g of Ginger) or not less than 0.7 mg and not more than 2.8 mg (for preparation prescribed 1.5 g of Ginger) of [6]-gingerol, per extract prepared with the amount specified in the Method of preparation.

Method of preparation

	1)	2)	3)
Euodia Fruit	3 g	4 g	3 g
Ginger	1 g	1.5 g	1.5 g
Ginseng	2 g	3 g	2 g
Jujube	4 g	3 g	4 g

Prepare a dry extract or viscous extract as directed under Extracts, according to the prescription 1) to 3), using the crude drugs shown above.

Description Goshuyuto Extract occurs as a light brown to light red-yellow powder, or a black-brown viscous extract. It has a slight odor and a hot and bitter taste.

Identification (1) To 1.0 g of the dry extract (or 3.0 g of the viscous extract) add 10 mL of sodium hydroxide TS, shake, add 5 mL of 1-butanol, shake, centrifuge, and use the supernatant liquid as the sample solution. Separately, to 1 g of pulverized euodia fruit add 10 mL of methanol, shake, centrifuge, and use the supernatant liquid as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 1 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of acetone, 2-propanol, water and formic acid (7:7:1:1) to a distance of about 7 cm, and airdry the plate. Examine under ultraviolet light (main

wavelength: 365 nm): one of the several spots obtained from the sample solution has the same color tone and Rf value with the blue-white fluorescent spot (Rf value: about 0.5) from the standard solution (Euodia Fruit).

- (2) To 1.0 g of the dry extract (or 3.0 g of the viscous extract) add 10 mL of water, shake, add 25 mL of diethyl ether, and shake. Separate the diethyl ether layer, evaporate the solvent under reduced pressure, add 2 mL of diethyl ether to the residue, and use this solution as the sample solution. Separately, dissolve 1 mg of [6]-gingerol for thinlayer chromatography in 1 mL of methanol, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L of the sample solution and 5 μ L of the standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate and hexane (1:1) to a distance of about 7 cm, and air-dry the plate. Spray evenly 4-dimethylaminobenzaldehyde TS for spraying on the plate, heat the plate at 105°C for 5 minutes, allow to cool, and spray water: one of the several spots obtained from the sample solution has the same color tone and Rf value with the blue-green to grayish green spot from the standard solution (Ginger).
- (3) To 1.0 g of the dry extract (or 3.0 g of the viscous extract) add 10 mL of sodium hydroxide TS, shake, add 5 mL of 1-butanol, shake, centrifuge, and use the supernatant liquid as the sample solution. Separately, dissolve 1 mg of ginsenoside Rb₁ for thin-layer chromatography in 1 mL of methanol, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 5 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, 1-propanol, water and acetic acid (100) (7:5:4:1) to a distance of about 7 cm, and air-dry the plate. Spray evenly vanillin-sulfuric acid-ethanol TS for spraying on the plate, heat the plate at 105°C for 5 minutes, and allow to cool: one of the several spots obtained from the sample solution has the same color tone and Rf value with the blue-purple spot from the standard solution (Ginseng).
- **Purity** (1) Heavy metals <1.07>—Prepare the test solution with 1.0 g of the dry extract (or an amount of the viscous extract, equivalent to 1.0 g of the dried substance) as directed under Extracts (4), and perform the test (not more than 30 ppm).
- (2) Arsenic <1.11>—Prepare the test solution with 0.67 g of the dry extract (or an amount of the viscous extract, equivalent to 0.67 g of the dried substance) according to Method 3, and perform the test (not more than 3 ppm).

Loss on drying $\langle 2.4I \rangle$ The dry extract: Not more than 11.0% (1 g, 105°C, 5 hours).

The viscous extract: Not more than 66.7% (1 g, 105°C, 5 hours).

Total ash $\langle 5.01 \rangle$ Not more than 10.0%, calculated on the dried basis.

Assay (1) Evodiamine—Weigh accurately about 0.5 g of the dry extract (or an amount of the viscous extract, equivalent to about 0.5 g of the dried substance), add exactly 50 mL of diluted methanol (7 in 10), shake for 30 minutes, filter, and use the filtrate as the sample solution. Separately, weigh accurately about 10 mg of evodiamine for assay, and dissolve in methanol to make exactly 200 mL, and use this solution as the standard solution. Perform the test with exactly $10 \,\mu\text{L}$ each of the sample solution and standard solution as directed under Liquid Chromatography $\langle 2.01 \rangle$ according to the following conditions, and determine the peak areas, A_T and A_S , of evodiamine in each solution.

Amount (mg) of evodiamine = $M_S \times A_T/A_S \times 1/4$

 M_S : Amount (mg) of evodiamine for assay taken, calculated on the basis of the content obtained by qNMR

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 282 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 40°C.

Mobile phase: A mixture of water, acetonitrile and phosphoric acid (620:380:1).

Flow rate: 1.0 mL per minute (the retention time of evodiamine is about 18 minutes).

System suitability—

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of evodiamine are not less than 5000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with $10 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of evodiamine is not more than 1.5%.

(2) [6]-Gingerol—Weigh accurately about 0.5 g of the dry extract (or an amount of the viscous extract, equivalent to about 0.5 g of the dried substance), add exactly 50 mL of diluted methanol (7 in 10), shake for 30 minutes, filter, and use the filtrate as the sample solution. Separately, weigh accurately about 10 mg of [6]-gingerol for assay, dissolve in methanol to make exactly 100 mL. Pipet 5 mL of this solution, add methanol to make exactly 50 mL, and use this solution as the standard solution. Perform the test with exactly $10 \,\mu$ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, A_T and A_S , of [6]-gingerol in each solution.

Amount (mg) of [6]-gingerol = $M_S \times A_T/A_S \times 1/20$

 M_S : Amount (mg) of [6]-gingerol for assay taken

Operating conditions—

Detector, column, column temperature and mobile phase:

Proceed as directed in the operating conditions in (1).

Flow rate: 1.0 mL per minute (the retention time of [6]-gingerol is about 14 minutes).

System suitability—

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of [6]-gingerol are not less than 5000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with $10 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of [6]-gingerol is not more than 1.5%.

Containers and storage Containers—Tight containers.

Hochuekkito Extract

補中益気湯エキス

Change the Identification (5) as follows:

Identification

(5) To 3.0 g of the dry extract (or 9.0 g of the viscous extract) add 30 mL of water, shake, then add 50 mL of diethyl ether, shake, and separate the diethyl ether layer. Evaporate the layer under reduced pressure, add 1 mL of diethyl ether to the residue, and use this solution as the sample solution. Separately, use (Z)-ligustilide TS for thinlayer chromatography as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate and hexane (1:1) to a distance of about 7 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 365 nm): one of the several spots obtained from the sample solution has the same color tone and Rf value with the blue-white fluorescent spot from the standard solution (Japanese Angelica Root).

Japanese Angelica Root

トウキ

Add the following next to the Description as follows:

Identification To 1.0 g of pulverized Japanese Angelica Root add 5 mL of methanol, shake for 10 minutes, centrifuge, and use the supernatant liquid as the sample solution. Separately, use (Z)-ligustilide TS for thin-layer chromatography as the standard solution (1). Dissolve 1 mg of scopoletin for thin-layer chromatography in 10 mL of methanol, and use this solution as the standard solution (2). Perform the test with these solutions as directed under Thin-layer Chromatography <2.03>. Spot $10~\mu$ L of the sam-

ple solution and $5 \mu L$ each of the standard solutions (1) and (2) on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of hexane, acetone and acetic acid (100) (30:25:1) to a distance of about 7 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 365 nm): two of the several spots obtained from the sample solution have the same color tones and Rf values with the corresponding blue-white fluorescent spots from the standard solutions (1) and (2).

Powdered Japanese Angelica Root

トウキ末

Add the following next to the Description as follows:

Identification To 1.0 g of Powdered Japanese Angelica Root add 5 mL of methanol, shake for 10 minutes, centrifuge, and use the supernatant liquid as the sample solution. Separately, use (Z)-ligustilide TS for thin-layer chromatography as the standard solution (1). Dissolve 1 mg of scopoletin for thin-layer chromatography in 10 mL of methanol, and use this solution as the standard solution (2). Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L of the sample solution and 5 μ L each of the standard solutions (1) and (2) on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of hexane, acetone and acetic acid (100) (30:25:1) to a distance of about 7 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 365 nm): two of the several spots obtained from the sample solution have the same color tones and Rf values with the corresponding blue-white fluorescent spots from the standard solutions (1) and (2).

Juzentaihoto Extract

十全大補湯エキス

Change the Identification (5) as follows:

Identification

(5) Shake 1.0 g of the dry extract (or 3.0 g of the viscous extract) with 15 mL of water and 5 mL of 0.1 mol/L hydrochloric acid TS, add 25 mL of diethyl ether, and shake. Separate the diethyl ether layer, evaporate the layer under reduced pressure, then add 2 mL of diethyl ether to the residue, and use this solution as the sample solution. Separately, use (Z)-ligustilide TS for thin-layer chromatography as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography <2.03>. Spot $10~\mu$ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate and hexane (1:1) to a distance of about 7 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength:

365 nm): one of the several spots obtained from the sample solution has the same color tone and *Rf* value with the bluewhite fluorescent spot from the standard solution (Cnidium Rhizome; Japanese Angelica Root).

Kakkontokasenkyushin'i Extract

葛根湯加川芎辛夷エキス

Change the Identification (7) as follows:

Identification

(7) Shake 1.0 g of the dry extract (or 3.0 g of the viscous extract) with 15 mL of water and 5 mL of 0.1 mol/L hydrochloric acid TS, and then shake with 25 mL of diethyl ether. Separate the diethyl ether layer, evaporate the solvent under reduced pressure, then dissolve the residue in 2 mL of diethyl ether, and use this solution as the sample solution. Separately, use (Z)-ligustilide TS for thin-layer chromatography as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate and hexane (1:1) to a distance of about 7 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 365 nm): one of the several spots obtained from the sample solution has the same color tone and Rf value with the blue-white fluorescent spot from the standard solution (Cnidium Rhizome).

Kamikihito Extract

加味帰脾湯エキス

Change the Identification (5) as follows:

Identification

(5) To 3.0 g of the dry extract (or 9.0 g of the viscous ex tract) add 15 mL of water, shake, then add 25 mL of diethyl ether, and shake. Separate the diethyl ether layer, evaporate the solvent under reduced pressure, then add 2 mL of diethyl ether to the residue, and use the solution as the sample solution. Separately, use (Z)-ligustilide TS for thin-layer chromatography as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L of the sample solution and 2 μ L of the standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of butyl acetate and hexane (2:1) to a distance of about 7 cm, and air-dry the plate. Examine the plate under ultraviolet light (main wavelength: 365 nm): one of the several spots obtained from the sample solution has the same color tone and Rf value with the blue-white fluorescent spot from the standard solution (Japanese Angelica Root).

Kamishoyosan Extract

加味逍遙散エキス

Change the Identification (1) as follows:

Identification (1) To 2.0 g of the dry extract (or 6.0 g of the viscous extract) add 10 mL of water, shake, then add 5 mL of diethyl ether, shake, centrifuge, and use the diethyl ether layer as the sample solution. Separately, use (Z)ligustilide TS for thin-layer chromatography as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate and hexane (1:1) to a distance of about 7 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 365 nm): one of the several spots obtained from the sample solution has the same color tone and Rf value with the blue-white fluorescent spot from the standard solution (Japanese Angelica Root).

Lithospermum Root

シコン

Change the Purity (1) as follows:

Purity (1) Heavy metals <1.07>—Proceed with 2.0 g of pulverized Lithospermum Root according to Method 3, and perform the test. Prepare the control solution with 3.0 mL of Standard Lead Solution (not more than 15 ppm).

Lycium Bark

ジコッピ

Change the Purity (1) as follows:

Purity (1) Heavy metals <1.07>—Proceed with 2.0 g of pulverized Lycium Bark according to Method 3, and perform the test. Prepare the control solution with 3.0 mL of Standard Lead Solution (not more than 15 ppm).

Nux Vomica

ホミカ

Change the Description as follows:

Description Disk, often slightly bent, 1-3 cm in diameter, 0.3-0.5 cm in thickness; externally light grayish yellow-green to light grayish brown, covered densely with lustrous appressed hairs radiating from the center to the circumference; on both sides, the margin and the central part

bulged a little; the dot-like micropyle situated at one point on the margin, and from which usually a raised line runs to the center on one side; extremely hard in texture; when cracked upon soaking in water, the seed coat thin, the interior consisting of two horny, light grayish yellow endosperms, and leaving a central narrow cavity at the center; a white embryo, about 0.7 cm in length, situated at one end between the inner surfaces of the endosperms.

Odorless.

Oriental Bezoar

ゴオウ

Change the origin/limits of content:

Oriental Bezoar is a stone formed in the gall sac of *Bos taurus* Linné var. *domesticus* Gmelin (*Bovidae*). It contains not less than 10.0% of bilirubin.

Delete the Identification (2):

Add the following next to the Total ash:

Assay Conduct this procedure without exposure to light using light-resistant vessels. The following sample solution and standard solution should be prepared before use. Weigh accurately about 10 mg of pulverized Oriental Bezoar, add 10 mL of a mixture of dimethyl sulfoxide and acetic acid (100) (9:1), warm at 60°C for 20 minutes, and add a mixture of dimethyl sulfoxide and acetic acid (100) (9:1) to make exactly 50 mL. Pipet 5 mL of this solution, add the mobile phase to make exactly 50 mL, filter through a membrane filter, and use the filtrate as the sample solution. Separately, weigh accurately about 10 mg of bilirubin for assay, add about 350 mg of L-ascorbic acid, and dissolve in a mixture of dimethyl sulfoxide and acetic acid (100) (9:1) to make exactly 100 mL. Pipet 5 mL of this solution, add the mobile phase to make exactly 50 mL, and use this solution as the standard solution. Perform the test with exactly 10 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, A_T and A_S , of bilirubin in each solution.

Amount (mg) of bilirubin = $M_S \times A_T/A_S \times 1/2$

M_S: Amount (mg) of bilirubin for assay taken

Operating conditions—

Detector: A visible absorption photometer (wavelength: 450 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 30° C.

Mobile phase: A mixture of acetonitrile and diluted acetic

acid (100) (1 in 100) (19:1).

Flow rate: Adjust so that the retention time of bilirubin is about 10 minutes.

System suitability-

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of bilirubin are not less than 5000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with $10 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of bilirubin is not more than 1.5%.

Delete the Content of the active principle:

Otsujito Extract

乙字湯エキス

Change the Identification (1) as follows:

Identification (1) Shake 2.0 g of the dry extract (or 6.0 g of the viscous extract) with 10 mL of water, add 10 mL of diethyl ether, shake, and centrifuge. Separate the diethyl ether layer, add 10 mL of sodium hydroxide TS, shake, centrifuge, separate the diethyl ether layer, and use this layer as the sample solution. Separately, use (Z)-ligustilide TS for thin-layer chromatography as the standard solution. Perform the test with these solutions as directed under Thinlayer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of butyl acetate and hexane (2:1) to a distance of about 7 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 365 nm): one of the several spots obtained from the sample solution has the same color tone and Rf value with the blue-white fluorescent spot from the standard solution (Japanese Angelica Root).

Platycodon Root

キキョウ

Delete the Identification (1), move up (2) to (1) and add the following next to (1):

Identification

(2) To 2.0 g of pulverized Platycodon Root add 20 mL of sodium carbonate TS, and shake. Add 5 mL of 1-butanol, shake for 10 minutes, centrifuge, and use the 1-butanol layer as the sample solution. Separately, dissolve 1 mg of platycodin D for thin-layer chromatography in 1 mL of methanol, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 5 μ L each of the sample solution and standard solution on a plate of silica

gel for thin-layer chromatography. Develop the plate with a mixture of 1-propanol, ethyl acetate, water and acetic acid (100) (5:3:2:1) to a distance of about 7 cm, and air-dry the plate. Spray evenly dilute sulfuric acid on the plate, and heat the plate at 105° C for 5 minutes: one of the several spots obtained from the sample solution and a spot from the standard solution show the same color tone and the same Rf value.

Powdered Platycodon Root

キキョウ末

Delete the Identification (1), move up (2) to (1) and add the following next to (1):

Identification

(2) To 2.0 g of Powdered Platycodon Root add 20 mL of sodium carbonate TS, and shake. Add 5 mL of 1butanol, shake for 10 minutes, centrifuge, and use the 1butanol layer as the sample solution. Separately, dissolve 1 mg of platycodin D for thin-layer chromatography in 1 mL of methanol, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 5 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of 1-propanol, ethyl acetate, water and acetic acid (100) (5:3:2:1) to a distance of about 7 cm, and air-dry the plate. Spray evenly dilute sulfuric acid on the plate, and heat the plate at 105°C for 5 minutes: one of the several spots obtained from the sample solution and a spot from the standard solution show the same color tone and the same Rf value.

Platycodon Fluidextract

キキョウ流エキス

Change the Method of preparation and Identification as follows:

Method of preparation 1) Take coarse powder of Platycodon Root, and prepare the fluidextract as directed under Fluidextracts using 25 vol% ethanol. An appropriate quantity of Ethanol and Purified Water or Purified Water in Containers may be used in place of 25 vol% ethanol.

2) Take Platycodon Root pulverized to suitable sizes, and prepare the fluidextract as directed under Fluidextracts using 25 vol% ethanol or diluted ethanol (1 in 4) as the solvent.

Identification To 2 mL of Platycodon Fluidextract add 20 mL of water and 5 mL of 1-butanol, mix, shake for 10 minutes, centrifuge, and use the 1-butanol layer as the sample solution. Separately, dissolve 1 mg of platycodin D for thin-layer chromatography in 1 mL of methanol, and use

this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 5 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of 1-propanol, ethyl acetate, water and acetic acid (100) (5:3:2:1) to a distance of about 7 cm, and air-dry the plate. Spray evenly dilute sulfuric acid on the plate, and heat the plate at $105\,^{\circ}$ C for 5 minutes: one of the several spots obtained from the sample solution has the same color tone and Rf value with the spot from the standard solution.

Add the following next to the Identification:

Alcohol number $\langle 1.01 \rangle$ Apply to Platycodon Fluidextract prepared by the Method of preparation 2). 2.0 - 3.0 (Method 1).

Change the Content of the active principle as follows:

Content of the active principle Transfer exactly 5 mL of Platycodon Fluidextract to a tared beaker or porcelain dish, evaporate to dryness on a water bath, and dry at 105°C for 5 hours: the mass of the residue is not less than 0.50 g.

Polygala Root

オンジ

Change the Identification (2) as follows:

Identification

(2) To 1.0 g of pulverized Polygala Root add 10 mL of a solution of sodium hydroxide (1 in 10), and heat under a reflux condenser on a water bath for 20 minutes. After cooling, add 10 mL of dilute hydrochloric acid, and shake. After cooling, add 10 mL of 1-butanol, shake for 10 minutes, centrifuge, and use the upper layer as the sample solution. Perform the test with the sample solution as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 5μ L of the sample solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, methanol, water and acetic acid (100) (20:4:2:1) to a distance of about 7 cm, and air-dry the plate. Spray evenly dilute sulfuric acid on the plate, and heat the plate at 105° C for 10 minutes: a red-brown to light brown spot appears at an Rf value of about 0.35.

Powdered Polygala Root

オンジ末

Change the Identification (2) as follows:

Identification

(2) To 1.0 g of Powdered Polygala Root add 10 mL of a solution of sodium hydroxide (1 in 10), and heat under a

reflux condenser on a water bath for 20 minutes. After cooling, add 10 mL of dilute hydrochloric acid, and shake. After cooling, add 10 mL of 1-butanol, shake for 10 minutes, centrifuge, and use the upper layer as the sample solution. Perform the test with the sample solution as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 5 μ L of the sample solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, methanol, water and acetic acid (100) (20:4:2:1) to a distance of about 7 cm, and air-dry the plate. Spray evenly dilute sulfuric acid on the plate, and heat the plate at 105°C for 10 minutes: a red-brown to light brown spot appears at an Rf value of about 0.35.

Safflower

コウカ

Change the Identification as follows:

Identification To 1.0 g of pulverized Safflower add 10 mL of a mixture of acetone and water (4:1), shake for 10 minutes, filter, and use the filtrate as the sample solution. Perform the test with the sample solution as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 5μ L of the sample solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, water, formic acid and methanol (35:15:10:2) to a distance of about 7 cm, and air-dry the plate: a red spot appears at an Rf value of about 0.5.

Saposhnikovia Root and Rhizome

ボウフウ

Change the Purity (1) as follows:

Purity (1) Heavy metals <1.07>—Proceed with 2.0 g of pulverized Saposhnikovia Root and Rhizome according to Method 3, and perform the test. Prepare the control solution with 3.0 mL of Standard Lead Solution (not more than 15 ppm).

Add the following next to the Purity (3):

Purity

(4) Peucedanum ledebourielloides—Place 1.0 g of pulverized Saposhnikovia Root and Rhizome in a glass-stoppered centrifuge tube, add 5 mL of hexane, shake for 10 minutes, centrifuge, and use the supernatant liquid as the sample solution. Separately, place 1.0 g of peucedanum ledebourielloides for purity in a glass-stoppered centrifuge tube, add 5 mL of hexane, shake for 10 minutes, centrifuge, and use the supernatant liquid as the standard solution. Perform the test with these solutions as directed under Thinlayer Chromatography $\langle 2.03 \rangle$. Spot 5 μ L each of the sample solution and standard solution on a plate of silica gel for

thin-layer chromatography. Develop the plate with a mixture of hexane, ethyl acetate and acetic acid (100) (20:10:1) to a distance of about 7 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 365 nm): the sample solution shows no spot corresponding to the blue fluorescent spot at an Rf value of about 0.4 obtained from the standard solution.

Scopolia Rhizome

ロートコン

Change the Description as follows:

Description Chiefly irregularly branched, slightly curved rhizome, about 15 cm in length, about 3 cm in diameter, occasionally longitudinally cut; externally grayish brown, with wrinkles; constrictions make the rhizome appear nodular; rarely, stem base at one end; stem scars at upper side of each node; roots or root scars on both sides and lower surface of rhizome; fractured surface granular, grayish white to light brown in color, with lighter colored cortex.

Odor characteristic.

Under a microscope <5.01>, xylem reveals groups of vessels arranged stepwise, and accompanied with xylem sieve tubes in medullary rays; parenchyma cells contain starch grains, and sometimes sand crystals of calcium oxalate.

Sinomenium Stem and Rhizome

ボウイ

Change the Identification as follows:

Identification To 1.0 g of pulverized Sinomenium Stem and Rhizome add 5 mL of methanol, shake for 10 minutes, filter, and use the filtrate as the sample solution. Separately, dissolve 1 mg of sinomenine for thin-layer chromatography in 1 mL of methanol, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography <2.03>. Spot 10 μL of the sample solution and 5 μ L of the standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl formate, 1-propanol, water and acetic acid (100) (3:3:2:2) to a distance of about 7 cm, and air-dry the plate. Spray evenly Dragendorff's TS for spraying on the plate, air-dry the plate, and spray evenly sodium nitrite TS on the plate: one of the several spots obtained from the sample solution and the spot from the standard solution show the same color tone and the same Rf value. Further, a spot with the same color tone appears immediately below the spot.

Swertia Herb

センブリ

Change the Identification as follows:

Identification To 0.5 g of pulverized Swertia Herb add 10 mL of methanol, shake for 5 minutes, filter, and use the filtrate as the sample solution. Separately, dissolve 1 mg of Swertiamarin RS or swertiamarin for thin-layer chromatography in 1 mL of methanol, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography <2.03>. Spot $5 \mu L$ each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, ethanol (99.5) and water (8:2:1) to a distance of about 7 cm, and air-dry the plate. Spray evenly vanillin-sulfuric acid-ethanol TS for spraying on the plate, and heat the plate at 105°C for 5 minutes: one of the several spots obtained from the sample solution and the spot from the standard solution show the same color tone and Rf value.

Powdered Swertia Herb

センブリ末

Change the Identification as follows:

Identification To 0.5 g of Powdered Swertia Herb add 10 mL of methanol, shake for 5 minutes, filter, and use the filtrate as the sample solution. Separately, dissolve 1 mg of Swertiamarin RS or swertiamarin for thin-layer chromatography in 1 mL of methanol, and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography <2.03>. Spot $5 \mu L$ each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography. Develop the plate with a mixture of ethyl acetate, ethanol (99.5) and water (8:2:1) to a distance of about 7 cm, and air-dry the plate. Spray evenly vanillin-sulfuric acid-ethanol TS for spraying on the plate, and heat the plate at 105°C for 5 minutes: one of the several spots obtained from the sample solution and the spot from the standard solution show the same color tone and Rf value.

Swertia and Sodium Bicarbonate Powder

センブリ·重曹散

Change the Identification (1) as follows:

Identification (1) To 10 g of Swertia and Sodium Bicarbonate Powder add 10 mL of ethanol (95), shake for 15 minutes, filter, and use the filtrate as the sample solution.

Separately, dissolve 2 mg of Swertiamarin RS or swertiamarin for thin-layer chromatography in 1 mL of ethanol (95), and use this solution as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot $10 \,\mu\text{L}$ of the sample solution and 5 μL of the standard solution on a plate of silica gel for thin-layer chromatography. Proceed as directed in the Identification under Powdered Swertia Herb.

Toad Cake

センソ

Change the Description as follows:

Description A round disk with slightly dented bottom and protuberant surface, about 8 cm in diameter, about 1.5 cm in thickness, the mass of one disk being about 80 to 90 g; or a round disk with almost flattened surfaces on both sides, about 3 cm in diameter, and about 0.5 cm in thickness, the mass of one disk being about 8 g; externally red-brown to black-brown, somewhat lustrous, approximately uniform and horny, hard in texture, and difficult to break; fractured surface nearly flat, and edges of broken pieces red-brown and translucent.

Odorless.

Tokishakuyakusan Extract

当帰芍薬散エキス

Change the Identification (1) as follows:

Identification (1) Shake 1.0 g of the dry extract (or 3.0 g of the viscous extract) with 15 mL of water and 5 mL of 0.1 mol/L hydrochloric acid TS, then add 25 mL of diethyl ether, and shake. Separate the diethyl ether layer, evaporate the layer under reduced pressure, add 2 mL of diethyl ether to the residue, and use this solution as the sample solution. Separately, use (Z)-ligustilide TS for thin-layer chromatography as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot 10 μ L each of the sample solution and standard solution on a plate of silica gel for thin-layer chromatography, develop the plate with a mixture of ethyl acetate and hexane (1:1) to a distance of about 7 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 365 nm): one of the several spots obtained from the sample solution has the same color tone and Rf value with the bluewhite fluorescent spot from the standard solution (Japanese Angelica Root; Cnidium Rhizome).

Change the Assay (1) as follows:

Assay (1) (E)-Ferulic acid—Conduct this procedure without exposure to light, using light-resistant vessels. Weigh accurately about 0.5 g of the dry extract (or an

amount of the viscous extract, equivalent to about 0.5 g of the dried substance), add exactly 50 mL of diluted methanol (1 in 2), shake for 15 minutes, filter, and use the filtrate as the sample solution. Separately, weigh accurately about 10 mg of (*E*)-ferulic acid for assay, and dissolve in diluted methanol (1 in 2) to make exactly 100 mL. Pipet 2 mL of this solution, add diluted methanol (1 in 2) to make exactly 50 mL, and use this solution as the standard solution. Perform the test with exactly 10 μ L each of the sample solution and standard solution as directed under Liquid Chromatography <2.01> according to the following conditions, and determine the peak areas, A_T and A_S , of (*E*)-ferulic acid in each solution.

Amount (mg) of (E)-ferulic acid = $M_S \times A_T/A_S \times 1/50$

 M_S : Amount (mg) of (E)-ferulic acid for assay taken

Operating conditions—

Detector: An ultraviolet absorption photometer (wavelength: 320 nm).

Column: A stainless steel column 4.6 mm in inside diameter and 15 cm in length, packed with octadecylsilanized silica gel for liquid chromatography (5 μ m in particle diameter).

Column temperature: A constant temperature of about 40°C.

Mobile phase: Dissolve 7.8 g of sodium dihydrogen phosphate in 1000 mL of water, and add 2 mL of phosphoric acid. To 850 mL of this solution add 150 mL of acetonitrile.

Flow rate: 1.0 mL per minute [the retention time of (*E*)-ferulic acid is about 10 minutes].

System suitability—

System performance: When the procedure is run with 10 μ L of the standard solution under the above operating conditions, the number of theoretical plates and the symmetry factor of the peak of (*E*)-ferulic acid are not less than 5000 and not more than 1.5, respectively.

System repeatability: When the test is repeated 6 times with $10 \mu L$ of the standard solution under the above operating conditions, the relative standard deviation of the peak area of (*E*)-ferulic acid is not more than 1.5%.

Yokukansan Extract

抑肝散エキス

Change the Identification (1) as follows:

Identification (1) To 2.0 g of the dry extract (or 6.0 g of the viscous extract) add 10 mL of water, shake, then add 10 mL of diethyl ether, shake, and centrifuge. Separate the diethyl ether layer, add 10 mL of sodium hydroxide TS, shake, centrifuge, separate the diethyl ether layer, and use this layer as the sample solution. Separately, use (Z)-ligustilide TS for thin-layer chromatography as the standard solution. Perform the test with these solutions as directed under Thin-layer Chromatography $\langle 2.03 \rangle$. Spot $10 \,\mu$ L each of the sample solution and standard solution on a plate of


2984 Crude Drugs and Related Drugs

silica gel for thin-layer chromatography. Develop the plate with a mixture of butyl acetate and hexane (2:1) to a distance of about 7 cm, and air-dry the plate. Examine under ultraviolet light (main wavelength: 365 nm): one of the several spots obtained from the sample solution has the same color tone and Rf value with the blue-white fluorescent spot from the standard solution (Japanese Angelica Root; Cnidium Rhizome).


Infrared Reference Spectra

Add the following spectra:


Bromfenac Sodium Hydrate



Doripenem Hydrate




Ethylcellulose




The JP Drugs are to be tested according to the provisions given in the pertinent monographs, General Notices, General Rules for Crude Drugs, General Rules for Preparations, and General Tests for their conformity to the Japanese Pharmacopoeia. (See the General Notices 5.)


Felodipine



Gatifloxacin Hydrate




Irinotecan Hydrochloride Hydrate





The JP Drugs are to be tested according to the provisions given in the pertinent monographs, General Notices, General Rules for Crude Drugs, General Rules for Preparations, and General Tests for their conformity to the Japanese Pharmacopoeia. (See the General Notices 5.)

Lanoconazole



Polaprezinc






The JP Drugs are to be tested according to the provisions given in the pertinent monographs, General Notices, General Rules for Crude Drugs, General Rules for Preparations, and General Tests for their conformity to the Japanese Pharmacopoeia. (See the General Notices 5.)


Ultraviolet-visible Reference Spectra

Add the following spectra:


Bromfenac Sodium Hydrate



Doripenem Hydrate




Felodipine





The JP Drugs are to be tested according to the provisions given in the pertinent monographs, General Notices, General Rules for Crude Drugs, General Rules for Preparations, and General Tests for their conformity to the Japanese Pharmacopoeia. (See the General Notices 5.)

Gatifloxacin Hydrate




Irinotecan Hydrochloride Hydrate





The JP Drugs are to be tested according to the provisions given in the pertinent monographs, General Notices, General Rules for Crude Drugs, General Rules for Preparations, and General Tests for their conformity to the Japanese Pharmacopoeia. (See the General Notices 5.)

Sitagliptin Phosphate Hydrate



General Information

GENERAL INFORMATION

G1 Physics and Chemistry

Add the following:

Control of Elemental Impurities in Drug Products

1. Introduction

Elemental impurities in drug products may arise from several sources; they may be residues intentionally added such as catalysts in the synthetic process of drug substances, drug substances being components of the drug product, impurities from natural products contained in additives, etc., and contaminants from manufacturing equipment and container/closure systems. The amounts of these impurities in drug products should be controlled within acceptable limits, except when they are stipulated in monographs.

The permitted daily exposures (PDEs) of elemental impurities are established to protect the health of all patients based on the evaluation of the toxic data of elemental impurities, and more strict limits are not needed if elemental impurities in drug products do not exceed the PDEs. In some cases, lower level of elemental impurities may be warranted when it is known that elemental impurities have been shown to have an impact on the quality attributes of the drug product (e.g., element catalyzed degradation of drug substances).

Elemental impurities in drug products are assessed and controlled based on a risk management approach.

2. Scope

The acceptable limit of elemental impurities apply to drug products, and also apply to drug products containing purified proteins and polypeptides (including proteins and polypeptides produced from genetic recombinant or non-recombinant origins), their derivatives, and drug products which they are components (e.g., conjugates) are within the scope of this guideline, as are drug products containing synthetic polypeptides, polynucleotides, and oligosaccharides.

It does not apply to crude drugs, radiopharmaceuticals, vaccines, cell metabolites, DNA products, allergenic extracts, cells, whole blood, cellar blood components, blood derivatives including plasma and plasma preparations, dialysate solutions not intended for systematic circulation, and drug products based on genes (gene therapy), cells (cell therapy) and tissues (tissue engineering). Also, it does not apply to elements that are intentionally included in the drug product for therapeutic benefit.

3. The PDEs for Elemental Impurities for Oral, Parenteral and Inhalation Routes of Administration, and Element Classification

The PDEs of elemental impurities established for preparations for oral, parenteral and inhalation routes of administration are shown in Table 1. If the PDEs for the other administration route are necessary, generally consider the oral PDE as a starting point in the establishment, and assess if the elemental impurity is expected to have local effects when administered by the intended route of administration.

Parenteral drug products with maximum daily volumes up to 2 L may use the maximum daily volume to calculate permissible concentrations from PDEs. For products whose daily volumes or general clinical practice may exceed 2 L (e.g., saline, dextrose, total parenteral nutrition, solutions for irrigation), a 2-L volume is used to calculate permissible concentrations from PDEs.

As shown in Table 1, elemental impurities are divided into three classes based on their toxicity (PDE) and likelihood of occurrence in the drug product. The likelihood of occurrence is judged from several factors, such as probabil-

Table 1 PDEs for Elemental Impurities

Element	Class	Oral PDE (µg/day)	Parenteral PDE (µg/day)	Inhalation PDE (µg/day)
Cd	1	5	2	2
Pb	1	5	5	5
As	1	15	15	2
Hg	1	30	3	1
Co	2A	50	5	3
V	2A 2A	100	10	1
v Ni	2A 2A	200	20	5
TI	2B	8	8	8
Au	2B	100	100	1
Pd	2B	100	10	1
Ir	2B	100	10	1
Os	2B	100	10	1
Rh	2B	100	10	1
Ru	2B	100	10	1
Se	2B	150	80	130
Ag	2B	150	10	7
Pt	2B	100	10	1
Li	3	550	250	25
Sb	3	1200	90	20
Ba	3	1400	700	300
Mo	3	3000	1500	10
Cu	3	3000	300	30
Sn	3	6000	600	60
Cr	3	11000	1100	3

ity of use in pharmaceutical processes, impurities in materials used in pharmaceutical processes, the observed natural abundance and environmental distribution of the element. Class 1: The elements, As, Cd, Hg, and Pb, are classified as this category and are human toxicant elements. As these elements are limited in the manufacture of pharmaceuticals, they are rarely used. Their presence in drug products usually comes from used materials such as mined excipients. These four elements require evaluation during the risk assessment, across all sources and routes of administration having possibility of contamination. Testing should only be applied when the risk assessment identifies it necessary to ensure that the PDE will be met, however it is not necessary for all components to determine for Class 1 elemental impurities. Class 2: Elemental impurities classified as Class 2 have lower toxicity than the elements in Class 1, and are routedependent human toxicants. These elements are further divided in 2A and 2B based on their relative likelihood of occurrence in the drug products. The class 2A elements are Co, Ni and V, which are known to exist naturally. These elements have relatively high probability of occurrence in drug products, and thus require evaluation during the risk assessment, across all sources and routes of administration having possibility of contamination. Because the Class 2B elements have the low probability of their existence in natural, they may be excluded from the risk assessment unless they are intentionally added during the manufacture of drug substances, excipients or other components of the drug product. The elemental impurities in Class 2B include Ag, Au, Ir, Os, Pd, Pt, Rh, Ru, Se and Tl.

Class 3: The elements in this class have relatively low toxicities by the oral route of administration, and their oral PDEs are more than $500 \,\mu\text{g}/\text{day}$. For oral routes of administration, unless these elements are intentionally added, they do not need to be considered during the risk assessment. For parenteral and inhalation products, the potential for inclusion of these elemental impurities should be evaluated even in the case where they are not intentionally added, unless the route specific PDE is above $500 \,\mu\text{g}/\text{day}$. The elements in this class include Ba, Cr, Cu, Li, Mo, Sb and Sn.

4. Risk Assessment and Control of Elemental Impurities

The technique of quality risk management should be considered in controls for elemental impurities in drug products, and the risk assessment should be based on scientific knowledge and principles. The risk assessment would be focused on assessing the levels of elemental impurities in a drug product in relation to the PDEs. Useful information for this risk assessment includes measured data of drug products and components, measured data and the risk assessment result supplied by drug substance and/or excipient manufacturers, and/or data available in published literature, but is not limited to them.

The risk assessment should be performed depending on the level of risk, and do not always require a formal risk management process. The use of informal risk management processes may also be considered acceptable.

4.1. General Principles

The risk assessment process consists of the following three steps.

- 1) Identify known and potential sources of elemental impurities that may find their way into the drug product.
- 2) Evaluate the presence of a particular elemental impurity in the drug product by determining the observed or predicted level of the impurity and comparing with the established PDE.
- 3) Summarize the risk assessment, and identify if controls built into the process are sufficient. Identify additional controls to be considered to limit elemental impurities in the drug product.

In many cases, the steps are considered simultaneously. The risk assessment may be iterated to develop a final approach to ensure the elemental impurities do not exceed the PDE certainly.

4.2. Source of Elemental Impurities

In considering the production of a drug product, there are broad categories of potential sources of elemental impurities.

- Residual impurities resulting from elements intentionally added (e.g., metal catalysts) in the formation of the drug substance, excipients or other components. The risk assessment of the drug substance should be studied about the potential for inclusion of elemental impurities in the drug product.
- Elemental impurities that are not intentionally added and are potentially present in the drug substance, water or excipients used in the preparation of the drug product.
- Elemental impurities that are potentially introduced into the drug substance and/or drug product from manufacturing equipment.
- Elemental impurities that have the potential to be leached into the drug substance and drug product from container closure systems.

During the risk assessment, the potential contributions from each of these sources should be considered to determine the overall contribution of elemental impurities to the drug product.

4.3. Identification of Potential Elemental Impurities

Potential elemental impurities derived from intentionally added catalysts and inorganic reagents: If any element is intentionally added, it should be considered in the risk assessment.

Potential elemental impurities that may be present in drug substances and excipients: While not intentionally added, some elemental impurities may be present in some drug substances and excipients. The possibility for inclusion of these elements in the drug product should be reflected in the risk assessment.

Potential elemental impurities derived from manufacturing equipment: The contribution of elemental impurities from this source may be limited and the subset of elemental impurities that should be considered in the risk assessment will depend on the manufacturing equipment used in the production of the drug product. The specific elemental

impurities of concern should be assessed based on the knowledge of the composition of the components of the manufacturing equipment that come in contact with components of the drug product. The risk assessment of this source of elemental impurities is one that can potentially be utilized for many drug products using similar process trains or processes.

In general, the processes used to prepare a given drug substance are considerably more aggressive than processes used in preparing the drug product when assessed relative to the potential to leach or remove elemental impurities from manufacturing equipment. Contributions of elemental impurities from drug product processing equipment would be expected to be lower than contributions observed for the drug substance. However, when this is not the case based on process knowledge or understanding, the potential for incorporation of elemental impurities from the drug product manufacturing equipment in the risk assessment (e.g., hot melt extrusion) should be considered.

Elemental impurities leached from container closure systems: The identification of potential elemental impurities that may be introduced from container closure systems should be based on a scientific understanding of likely interactions between a particular drug product type and its packaging. When a review of the materials of construction demonstrates that the container closure system does not contain elemental impurities, no additional risk assessment needs to be performed. It is recognized that the probability of elemental leaching into solid dosage forms is minimal and does not require further consideration in the risk assessment. For liquid and semi-solid dosage forms there is a higher probability that elemental impurities could leach from the container closure system during the shelf-life of the drug product. Studies to understand potential leachables from the container closure system (after washing, sterilization, irradiation, etc.) should be performed.

Factors that should be considered (for liquid and semisolid dosage forms) are shown as follows, but are not limited.

 Hydrophilicity/hydrophobicity, Ionic content, pH, Temperature (cold chain vs room temperature and processing conditions), Contact surface area, Container/material composition, Terminal sterilization, Packaging process, Material sterilization, Duration of storage

Table 2 provides recommendations for inclusion of elemental impurities in the risk assessment. This can be applied to all sources of elemental impurities in the drug product.

4.4. Evaluation

As the potential elemental impurity identification process is concluded, there are following two possible outcomes.

- 1) The risk assessment process does not identify any potential elemental impurities.
- 2) The risk assessment process identifies one or more potential elemental impurities. For any elemental impurities identified in the process, the risk assessment should consider if there are multiple sources of the identified elemental im-

Table 2 Elements to be Considered in the Risk Assessment

T	CI.	If intention-	If not intentionally add			
Element	Class	ally added (all routes)	Oral	Parenteral	Inhalation	
Cd	1	0	\circ	0	$\overline{}$	
Pb	1	\bigcirc	\bigcirc	\bigcirc	\bigcirc	
As	1	\bigcirc	0000	0 0 0	\bigcirc	
Hg	1	\bigcirc	\bigcirc	\bigcirc	0 0 0 0	
Co	2A	\bigcirc	\bigcirc	\bigcirc	\bigcirc	
V	2A	\bigcirc	\bigcirc	\bigcirc	\circ	
Ni	2A	\bigcirc		\bigcirc	Ō	
TI	2B	\bigcirc	×	×	×	
Au	2B	\bigcirc	×	×	×	
Pd	2B	\bigcirc	×	×	×	
Ir	2B	\bigcirc	×	×	×	
Os	2B	\bigcirc	×	×	×	
Rh	2B	\bigcirc	×	×	×	
Ru	2B	\bigcirc	×	×	×	
Se	2B	\bigcirc	×	×	×	
Ag	2B	\bigcirc	×	×	×	
Pt	2B	\bigcirc	×	×	×	
Li	3	\bigcirc	×	\bigcirc	\bigcirc	
Sb	3	\bigcirc	×	\bigcirc	\bigcirc	
Ba	3	\bigcirc	×	×	\bigcirc	
Mo	3	\bigcirc	×	×	\bigcirc	
Cu	3	000000000000000000000000000000000000000	×	\bigcirc	\bigcirc	
Sn	3	\bigcirc	×	×	0 0 0 0 0	
Cr	3	0	×	×	0	

: necessary : unnecessary

purity or impurities.

During the risk assessment, a number of factors that can influence the level of the potential elemental impurity in the drug product should be considered.

4.5. Summary of Risk Assessment Process

The risk assessment is summarized by reviewing relevant product or component specific data combined with information and knowledge gained across products or processes to identify the significant probable elemental impurities that may be observed in the drug product.

The significance of the observed or predicted level of the elemental impurity should be considered in relation to the PDE of the elemental impurity. As a measure of the significance of the observed elemental impurity level, a control threshold is defined as a level that is 30% of the established PDE in the drug product. The control threshold may be used to determine if additional controls may be required.

If the total elemental impurity level from all sources in the drug product is expected to be consistently less than 30% of the PDE, then additional controls are not required, provided adequate controls on elemental impurities are demonstrated by the appropriate assessment of the data.

If the risk assessment fails to demonstrate that an elemental impurity level is consistently less than the control threshold, controls should be established to ensure that the

elemental impurity level does not exceed the PDE in the drug product.

The variability of the level of an elemental impurity should be factored into the application of the control threshold to drug products. Sources of variability may include the following.

- · Variability of the analytical method
- Variability of the elemental impurity level in the specific sources
- Variability of the elemental impurity level in the drug product

For some components that have inherent variability (e.g., mined excipients), additional data may be needed to apply the control threshold.

5. Converting between PDEs and Concentration Limits

The PDEs reported in μ g per day (μ g/day) give the maximum permitted quantity of each element that may be contained in the maximum daily dose of a drug product. Because the PDE reflects total exposure from the drug product, it is useful to convert the PDE into concentrations as a tool in evaluating elemental impurities in drug products or their components. Any of the following options may be selectable as long as the resulting permitted concentrations assure that the drug product does not exceed the PDEs. In the choice of a specific option the daily dose of the drug product needs to be determined or assumed.

Option 1: Common permitted concentration limits of elements across drug product components for drug products with daily doses of not more than 10 g: This option is not intended to imply that all elements are present at the same concentration, but rather provides a simplified approach to the calculations. The option assumes the daily dose of the drug product is 10 g or less, and that elemental impurities identified in the risk assessment (the target elements) are present in all components of the drug product. Using Equation (1) below and a daily dose of 10 g of drug product, this option calculates a common permissible target elemental concentration for each component in the drug product.

Concentration
$$(\mu g/g)$$

= $\frac{PDE(\mu g/day)}{\text{daily dose of drug product } (g/day)}$ (1)

This approach, for each target element, allows determination of a fixed common maximum concentration in μ g per g in each component. The permitted concentrations are provided in Table 3.

If all the components in a drug product do not exceed the Option 1 permitted concentrations for all target elements identified in the risk assessment, then all these components may be used in any proportion in the drug product. If the permitted concentrations in Table 3 are not applied, Options 2a, 2b, or 3 should be followed.

Option 2a: Common permitted concentration limits of elements across drug product components for a drug product with a specified daily dose: This option is similar to Option 1, except that the drug daily dose is not assumed to be 10 g. The common permitted concentration of each element is de-

Table 3 Permitted Concentrations of Elemental Impurities for Option 1

E 1 .	CI	Oral	Parenteral	Inhalation
Element	Class		Concentration	
		$(\mu g/g)$	$(\mu g/g)$	$(\mu g/g)$
Cd	1	0.5	0.2	0.2
Pb	1	0.5	0.5	0.5
As	1	1.5	1.5	0.2
Hg	1	3	0.3	0.1
Co	2A	5	0.5	0.3
V	2A	10	1	0.1
Ni	2A	20	2	0.5
TI	2B	0.8	0.8	0.8
Au	2B	10	10	0.1
Pd	2B	10	1	0.1
Ir	2B	10	1	0.1
Os	2B	10	1	0.1
Rh	2B	10	1	0.1
Ru	2B	10	1	0.1
Se	2B	15	8	13
Ag	2B	15	1	0.7
Pt	2B	10	1	0.1
Li	3	55	25	2.5
Sb	3	120	9	2
Ba	3	140	70	30
Mo	3	300	150	1
Cu	3	300	30	3
Sn	3	600	60	6
Cr	3	1100	110	0.3

termined using Equation (1) and the actual maximum daily dose. This approach, for each target element, allows determination of a fixed common maximum concentration in μg per g in each component based on the actual daily dose provided. If all components in a drug product do not exceed the Option 2a permitted concentrations for all target elements identified in the risk assessment, then all these components may be used in any proportion in the drug product. Option 2b: Permitted concentration limits of elements in individual components of a drug product with a specified daily dose: Permitted concentrations based on the distribution of elements in the components (e.g., higher concentrations in components with the presence of an element in question) may be set. For each element identified as potentially present in the components of the drug product, the maximum expected mass of the elemental impurity in the final drug product can be calculated by multiplying the mass of each component material times the permitted concentration pre-established in each material and summing over all components in the drug product, as described in Equation (2). The total mass of the elemental impurity in the drug product should comply with the PDEs unless justified according to other relevant sections of this general information. If the risk assessment has determined that a specific element is not a potential impurity in a specific component,

there is no need to establish a quantitative result for that element in that component. This approach allows that the maximum permitted concentration of an element in certain components of the drug product may be higher than the Option 1 or Option 2a limit, but this should then be compensated by lower allowable concentrations in the other components of the drug product. Equation (2) may be used to demonstrate that component-specific limits for each element in each component of a drug product assure that the PDE will be met.

$$PDE (\mu g/day) \ge \sum_{k=1}^{N} C_k \cdot M_k$$
 (2)

k =an index for each of N components in the drug product

 C_k = permitted concentration of the elemental impurity in component $k (\mu g/g)$

 M_k = mass of component k in the maximum daily dose of the drug product (g)

Option 3: Finished Product Analysis: The concentration of each element may be measured in the final drug product. Equation (1) may be used with the maximum total daily dose of the drug product to calculate a maximum permitted concentration of the elemental impurity.

6. Speciation and Other Considerations

Speciation is defined as the distribution of elements among chemical species based on the difference of molecular structure including ionic element, molecules, or complexes, reflecting isotopic composition, electronic or oxidation state. When the toxicities of different species of the same element are known to be different, the PDE has been established using the toxicity information on the species expected to be in the drug product.

When elemental impurity measurements are used in the risk assessment, total elemental impurity levels in drug products may be used to assess compliance with the PDEs. The identification of speciation is not particularly expected, however such information could be used to justify lower or higher levels when the identified species is more or less toxic, respectively, than the species used for the calculation of the PDEs.

When total elemental impurity levels in components are used in the risk assessment, providing information on release of an elemental impurity from the component in which it is found is not expected. However, such information could be used to justify levels higher than those based on the total elemental impurity content of the drug product.

7. Analytical Procedures

The determination of elemental impurities should be conducted using appropriate procedures suitable for their intended purposes. Unless otherwise justified, the test should be specific for each elemental impurity identified for control during the risk assessment. Elemental Impurities-Procedures <2.66> or suitable alternative procedures for determining levels of elemental impurities should be used.

8. Lifecycle Management

If changes to the drug product or components have the potential to change the elemental impurity content of the drug product, the risk assessment, including established controls for elemental impurities, should be re-evaluated. Such changes could include changes in synthetic routes, excipient suppliers, raw materials, processes, equipment, container closure systems or facilities.

G3 Biotechnological/Biological Products

Add the following:

Host Cell Protein Assay

Host cell protein (HCP) is a general term for proteins derived from host cells used for the production of pharmaceutical products. This general information describes HCP assays for therapeutic proteins produced by recombinant DNA technology (recombinant therapeutic proteins).

Residual HCP in recombinant therapeutic proteins has potential to elicit immune responses against itself and may also act as adjuvants to induce anti-drug antibodies. Therefore, in order to ensure the efficacy and safety of recombinant therapeutic proteins, it is necessary to establish a purification process to reduce HCP to a level that does not affect safety. In addition, residual levels of HCP must be appropriately controlled by verifying that in-process tests can consistently eliminate HCP or by establishing purity testing of drug substance.

1. Selection of Test Methods for HCP assay

HCP assay is usually performed with a sandwich immunoassay using antibodies against HCP (anti-HCP antibodies) and detection systems including enzyme-linked immunosorbent assay (ELISA), electrochemiluminescent immunoassay (ECLIA), and time-resolved fluorescent immunoassay (TRFIA). This general information addresses the sandwich immunoassay, but does not discourage other assays.

Residual HCP from the manufacturing process of recombinant therapeutic proteins consists a large number of proteins and may have profiles that vary from one host cell to another or depending on manufacturing conditions. Based on differences in its intended use or differences in the concept of preparation of anti-HCP antibodies used for testing, HCP assay is classified into generic assay, product-specific assay, and platform assay. Generic assay is intended for wide use in pharmaceutical products manufactured using similar host cells (e.g., CHO-K1 or CHO-DG44 cells derived from CHO [Chinese hamster ovary] cells) and is a test method established using anti-HCP antibodies that are prepared with proteins from all conponents of the host cells (cell lysate or culture supernatant) as immunogens. Commercially available reagents or kits for an HCP assay com-

monly referred to as generic assay and need to be validated before use. Product-specific assay is intended to control HCP in a specific product and developed in consideration of characteristics of the manufacturing process of the product. Platform assay is developed for the application to recombinant therapeutic proteins (e.g., monoclonal antibodies with adequate experience) produced by a manufacturing platform.

Generic assay is intended to comprehensively obtain antibodies against a wide range of HCP by using proteins from all conponents of the host cells as antigens. However, it should be noted that it is difficult to obtain antibodies against all HCPs because of differences in proportions or immunogenicity of individual proteins used as antigens and that residual HCP from actual manufacturing processes may be inadequately covered because of potential different profiles of residual HCP from different manufacturing processes. In contrast, product-specific assay is expected to prepare anti-HCP antibodies that can detect residual HCP from actual manufacturing processes compared with generic assay because it uses potentially residual HCP as antigens. However, it should be noted that profiles of residual HCP may vary with modifications to manufacturing processes. Platform assay has both aspects, namely, the assay has an advantage of being able to be applied to various products prepared by a manufacturing platform; however, the assay may involve issues as generic or product-specific assay, depending on the method of preparation of antigens used for preparing anti-HCP antibodies.

In some products, it is possible that specific HCPs binding to desired products are present or that HCPs markedly increase in production amount with expression of desired products. If the residue of these HCPs is found, the need to establish other test methods for the HCPs is considered.

In light of these characteristics of the test methods for HCP assays and in consideration of properties of host cells, characteristics of manufacturing processes, knowledge about immunogenicity of HCPs, stages of development of products, etc., an appropriate test method is selected.

2. Preparation and Characterization of Reagents

2.1. HCP Antigens/HCP Reference Materials

For antigens to produce antibodies that specifically detect HCPs in products, it is necessary to prepare HCP not containing desired products. Usually, null cells are used for preparation of HCP antigen while keeping in mind that HCP appropriate to the purpose of the HCP assay are comprehensively contained. In addition, HCP are used not only as antigens but also as reference materials for HCP assay and may also be used as ligands for purification of anti-HCP antibodies by affinity chromatography.

2.1.1. Preparation of HCP Antigens/HCP Reference Materials

The method of preparation of HCP antigens/HCP reference materials varies widely depending on the type of test methods. The method of preparation of HCP used as antigens for preparation of anti-HCP antibodies or as HCP reference materials in the test methods is shown below along

with points to note.

HCP used for generic assays are prepared from culture supernatant or lysed or disrupted null cells using minimal procedures such as concentration and dialysis and keeping preservation of component proteins in mind. It should be noted that these HCPs show profiles that are different from those of residual HCP in products because of preparation under conditions different from culture processes at commercial scale.

HCP used for product-specific assays are prepared from null cells using manufacturing processes of products. Usually, the application of the purification process is minimized to obtain a wide range of HCP. However, if a suitable antibody to detect the HCP is not obtained, it may be necessary to prepare appropriate HCP antigens by exploring conditions of preparation of HCP or excluding specific HCP.

HCP used for the platform assays are prepared from null cells using a manufacturing platform that can be used in multiple products. Usually, as with the other test methods, the application of the purification process is minimized to obtain a wide range and an enough amount of HCP. In addition, a mixture of HCPs prepared under multiple conditions can be used to address differences in HCP spectra due to slight differences in manufacturing conditions.

The use of mock cells as null cells used for preparation of HCP for product-specific and platform assays has advantages such as the presence of proteins expressed as selection markers in antigens and ability to culture the cells under similar culture conditions in actual manufacturing processes. In contrast, it should be noted that even a same cell line does not show consistent properties (such as cell growth rate) among different clones and that differences such as the presence or absence of production of desired products may cause different profiles of HCP.

2.1.2. Characterization of HCP Antigens/HCP Reference Materials

Prepared HCPs are analyzed for the following items.

1) Protein Concentrations

Protein concentrations are determined by a suitable measurement method keeping in mind that host cell nucleic acids or culture medium components may be contained depending on the method of preparation of HCP. For information on detailed measurement methods and points to consider, "Protein Assay" in General Information would be helpful.

2) HCP Profiles

Usually, one-dimensional electrophoresis (SDS-PAGE) or two-dimensional electrophoresis is used to confirm that prepared HCP include HCP species that are considered to remain in manufacturing processes or drug substances. Identification of HCP species by mass spectrometry is also a helpful approach.

2.2. Anti-HCP Antibodies

2.2.1. Preparation of Anti-HCP Antibodies

Since HCPs represent a heterogeneous variety of different proteins, polyclonal antibodies are obtained as anti-HCP antibodies used for the assay to comprehensively detect HCPs. The rabbit, goat, and sheep are commonly used animal species for immunization. For immunization, it is

useful to enhance immune response with adjuvants. Due to different degrees of immunogenicity of individual proteins comprising HCPs, the timing of induction of antibodies or the amount of antibodies produced is not constant, regardless of the amount of proteins as antigens. In addition, inter-individual variability in animals used for immunization makes the profile of induced antibodies inconsistent. Several rounds of immunization are usually required and, after determining the reactivity of induced antibodies with HCPs by Western blotting using antisera in each period of immunization, whole serum is collected. Mixing of anti-HCP antibodies derived from multiple individuals is intended to obtain adequate amounts of antibodies and is also expected to contribute to the elimination of imbalance of HCP profiles.

Anti-HCP antibodies are purified from the obtained antiserum by Protein A or Protein G chromatography. In either case, aggregates may be formed from some antibodies due to use of acid conditions for antibody elution from the columns. It is useful to remove antibody aggregates by a suitable method because they may cause interference with measurement.

Anti-HCP antibodies can also be purified by affinity chromatography using HCP as ligands. This purification is expected to eliminate non-specific reactions because of concentration of antibodies specific to HCPs, however, it should be noted that anti-HCP antibodies may get less diverse due to less adsorption of low-affinity antibodies or less elution of very high-affinity antibodies.

2.2.2. Suitability of Anti-HCP Antibodies

Anti-HCP antibodies have to be able to comprehensively recognize HCPs with wide ranges of electric charges and molecular masses that potentially remain in manufacturing processes or drug substances. However, because differences in immunogenicity of each HCP species may make antibodies against some HCPs less likely to be induced, obtained anti-HCP antibodies have to be qualified, usually as measured by antigen coverage. A specific example of assessment methods is shown below. After separation of HCPs by two-dimensional electrophoresis, total protein on the gel is stained. After performing two-dimensional electrophoresis in the same manner, Western blotting using anti-HCP antibodies is performed. Spot patterns obtained from each staining are compared and the proportion of spots detected by Western blotting vs spots obtained in total protein staining is determined as antigen coverage.

2.3. Storage of Reagents

HCP reference materials and anti-HCP antibodies are stored with attention to stability. The stability of these reagents can be confirmed by continuously monitoring parameters of dose-response curves of reference materials.

3. Validation of HCP Assay

When a sandwich immunoassay is used for HCP assay, "Validation of Analytical Procedures" in General Information would be helpful for information on basic requirements for validation. However, unlike a conventional sandwich immunoassay to determine a single antigen, HCP assay is

an approach to use antibodies prepared with a mixture of various HCP species as antigens to determine the HCP species simultaneously. Therefore, changes in concentrations as a function of dilution ratios of samples (dilution linearity) may not be observed in highly purified samples even within the quantitation limits in which linearity has been obtained for HCP reference materials. This phenomenon is likely attributed to insufficient amounts of antibodies to some HCPs that are found in increased proportions in samples for measurement due to the difference of removal rates of individual HCPs in the purification process and may lead to underestimation of HCP concentrations.

Therefore, HCP assay should be validated for accuracy, precision, specificity, standard curve, quantitative range, and dilution linearity.

(1) Accuracy and Precision

Accuracy and precision are expressed as coefficients of variation of recovery rates and quantitative values of HCP reference materials, respectively, by performing spike and recovery tests of HCP reference materials in purification process pools or drug substances to be measured.

(2) Specificity

Because HCP assay involves measurement of trace amounts of HCPs remaining in samples containing large quantities of desired products, it should be confirmed that there is no interference of desired substances or ingredients in sample solutions.

(3) Standard Curve and Quantitative Range

A standard curve is generated using serially diluted HCP reference materials, a regression expression is obtained, and validity is expressed in terms of determination coefficient, etc. Quantitative values of reference materials at each concentration level are determined from the regression expression and a range of concentrations with acceptable levels of accuracy and precision is defined as quantitative range with the lowest concentration within the range being defined as the minimum limit of quantitation.

(4) Dilution Linearity

Purification process pools or drug substances to be measured are examined for the range of dilution ratios of samples in which quantitative values of samples diluted within the quantitative range of a standard curve are linear.

4. Establishment of HCP Assay

HCP assay is used for confirmation of the status of removal of HCP in manufacturing processes or as a purity test of drug substances. For information on basic concepts of procedures or data analysis in HCP assay, "Enzymelinked Immunosorbent Assay" in General Information would be helpful.

Results from HCP assay are usually presented as contents in desired products. Contents of HCP per total protein or desired products can be determined by separately determining concentrations of the proteins in samples. For example, when total protein concentration is 2 mg/mL and an HCP concentration is 20 ng/mL, the content of HCP should be indicated as 10 ng/mg.

5. Others

5.1. Considerations for Modifications to Manufacturing Processes

Since any modifications made to manufacturing processes for recombinant therapeutic proteins may affect profiles of residual HCP, it is necessary to confirm that appropriate measurement of HCP is performed after modifications to manufacturing processes. If HCP profiles are altered and it is considered inappropriate to apply a test method for HCP assay before modifications to manufacturing processes, a test method for HCP assay has to be established again. Changes in HCP profiles can be assessed by an analytical procedure such as two-dimensional electrophoresis or mass spectrometry.

5.2. Considerations for Modifications to Reagents for Tests

It is desirable to secure sufficient quantities of HCP antigens/HCP reference materials and anti-HCP antibodies, all important reagents, whenever possible, in consideration of the life cycle of products. When HCP antigens/HCP reference materials or anti-HCP antibodies are newly prepared, it should be confirmed that their characteristics are comparable before and after renewal by analytical procedures such as two-dimensional electrophoresis, Western blotting, and mass spectrometry. In addition, test methods should be validated again for necessary items and newly prepared reagents should be used after confirming their consistency with the reagents and test methods before renewal

When generic assay is employed, qualified commercially available kits can also be used. However, availability of information on lot renewal of reagents, etc. is necessary to ensure consistency and quality of tests using commercially available kits, and therefore characterization of important reagents and method validation should be performed as needed.

6. Terms

Mock cell: A cell line established from the host cell line by transferring expression vectors that do not contain genes encoding desired products.

Null cell: A host cell that does not express desired products. It includes parent cells or mock cells.

Antigen coverage: The rate of detection of proteins comprising HCPs by anti-HCP antibodies. For example, an antigen coverage can be calculated from the number of spots obtained in total protein staining and the number of spots obtained in Western blotting using anti-HCP antibodies after separation of HCPs by two-dimensional electrophoresis.

Total Protein Assay

Delete the description of harmonization in the beginning and symbols ($^{\spadesuit}$ $_{\spadesuit}$).

G4 Microorganisms

Delete the following items:

Media Fill Test (Process Simulation)

Microbiological Environmental Monitoring Methods of Processing Areas for Sterile Pharmaceutical Products

Parametric Release of Terminally Sterilized Pharmaceutical Products

G5 Crude Drugs

Purity Tests on Crude Drugs using Genetic Information

Change the following as follows:

The first step in the quality assurance of natural products is the use of raw materials from the right part of the right origin. Therefore, it is clearly stated in Article 4 of the General Rules For Crude Drugs that the source of a crude drug is an approval or rejection criterion. There are various methods for differentiating the sources of crude drugs, such as morphological methods, organoleptic tests, and chemical methods, and appropriate methods for each are described in the individual monographs. Morphological methods, organoleptic tests, and chemical methods are discrimination methods for species that are based on the phenotypic characteristics of the crude drugs. On the other hand, together with recent progress in molecular biology techniques and the accumulation of genetic information on plants, differentiating methods of crude drugs based on genotypes is being established. Unlike morphological and other methods that are based on phenotypic characteristics, the genotypic methods are not affected by environmental factors. Also, the methods have several advantages, such as specialized expertise and skill for classification are not needed, and objective results are easily obtained.

The evolution of living organisms is accomplished by genetic mutation, and differences among the nucleotide sequences of genes of closely related species reflect the strain relationships between the species. Based on this theory, methods that classify species phylogenetically using the nucleotide sequence of rDNA that codes for ribosomal RNA (rRNA) on the nuclear genome have recently been adopted for the classification of microorganisms. In the same way, the sequence of this rDNA is also most often

used in the classification of higher plants based on the genotype. In particular, it is very easy to classify closely related species using the internal transcribed spacer (ITS) region of the rDNA, since nucleotide substitution is more often undertaken by comparison with the coded gene region. Furthermore, since the genes on the nuclear genome originate from the parents' genomes, there is an advantage that interspecies hybrids can be detected. Higher plants also have mitochondrial genes and chloroplastic genes. Although the genes on these genomes are also often used for classification, interspecies hybrids cannot be confirmed because the genes are normally uniparental inheritance.

The three methods presented here are, 1) the purity test of Atractylodes Rhizome for Atractylodes Lancea Rhizome, 2) the purity test of Saposhnikovia Root and Rhizome for *Peucedanum ledebourielloides*, which are developed based on the difference of the gene sequence of the ITS region of rDNA recently reported¹⁻⁴), and the inter-laboratory validation study have been completed.

The plant sources for Atractylodes Lancea Rhizome stipulated in the individual monographs are Atractylodes lancea De Candolle and A. chinensis Koidzumi (Compositae), while those for Atractylodes Rhizome are A. japonica Koidzumi ex Kitamura and A. macrocephala Koidzumi (Compositae). The approval or rejection of the both sources is, in principle, determined by the description of the crude drug, including microscopy, together with thin-layer chromatography in identification tests. In the above scientific paper, it was shown that these 4 plant species can be clearly classified by comparing the nucleotide sequences of the ITS region mentioned above, and that the species can be easily classified without performing sequence analysis by performing PCR using a species-specific primer pair or by using a restriction enzyme which recognizes species-specific sequence.

Likewise, the plant source of Saposhnikovia Root and Rhizome is stipulated as *Saposhnikovia divaricata* Schischkin (*Umbelliferae*), and the approval or rejection of the source is determined by the description of the crude drug and thin-layer chromatography in identification tests. According to the report⁴, crude drugs treated as Saposhnikovia Root and Rhizome in Shaanxi and Shanxi Provinces are frequently derived from *Peucedanum ledebourielloides*, and it was shown that the differentiation of the both is possible by using the nucleotide sequence in the ITS region of the rDNA.

In purity tests on crude drugs using genetic information, the simplicity of the test is given maximum consideration. We established methods that observe PCR amplification bands using species-specific primer pair (Mutant Allele Specific Amplification: Method 1) and methods that observe DNA fragments produced by restriction enzyme treatment of the PCR products, which are prepared using a primer pair common to each plant source (PCR—Restriction Fragment Length Polymorphism: Method 2), without nucleotide sequence analyses. In these methods based on PCR, an extremely small amount of template DNA is amplified to billions to hundreds of billions times. Therefore, when

using them as identification tests for powdered crude drugs, the target DNA fragment can be observed even if the vast majority of the crude drug for analysis is not appropriate plant species and there is only a minute amount of powder from a crude drug derived from a suitable plant. Consequently, in identification tests, either a cut or a whole crude drug must be used, as long as one is careful to avoid contamination by powder originating from other crude drugs. On the other hand, when used as a purity test, the form of the crude drug is irrelevant as long as the gene amplification is performed properly and the target gene is not polymorphic, so if DNA fragments of an inappropriate plant to be examined are confirmed in the purity test, regardless of the form of the crude drug, it becomes clear there is contamination by an inappropriate crude drug to be examined.

The methods shown here are general information and at the present stage results obtained using the methods do not affect the approval or rejection of the crude drug in each monograph. Furthermore, by performing the sequence analysis outlined in the previous paper for a crude drug sample derived from a single individual, it goes without saying that more accurate decision concerning the source species can be made.

1. DNA Amplification Equipment

DNA amplification equipment is used to amplify the DNA which is extracted from a crude drug and then purified. Since there are slight differences in the methods of temperature control, and so on depending on the equipment used, there may be differences in the intensity, etc. of the PCR amplification bands even if PCR is carried out under the stipulated conditions. Therefore, when judging results based solely on the presence or absence of PCR amplification bands as in Methods 1, the use of equipment described in the JAS analytical test handbook: genetically modified food quality, labeling analysis manual for individual products⁵⁾ is recommended. When other equipment is used, confirm that only proper amplification bands are obtained by performing PCR using DNA obtained from samples confirmed beforehand to be the source species. If proper amplification bands are not obtained, the PCR temperature conditions should be slightly adjusted. This equipment can be used for the restriction enzyme treatment in Method 2.

2. General precautions

Crude drugs are different from fresh plants in that they are dried products and a certain amount of time has passed since they were harvested. Therefore, in many cases the DNA has undergone fragmentation. Furthermore, various substances that can block or interfere with the PCR reaction may be present in the plant. For these reasons, the extraction and purification of template DNA is the process that should receive the greatest amount of attention. In the case of Atractylodes crude drugs, the periderm should be removed using a clean scalpel or other clean instrument before pulverizing the sample because very often there are inhibitory substances present in the periderm.

The PCR used for this test is the technique that amplify the target DNA more than hundreds of millions times, and a trace of contamination leads an incorrect result. Therefore, careful attention is required to prevent contamination. For treatment to prevent contamination, refer to the prevention of contamination section⁶⁾ in the above manual.

3. Purity test of Atractylodes Rhizome for Atractylodes Lancea Rhizome

3.1. Method 1 (Mutant Allele Specific Amplification Method)

Generally, this method is referred to as Mutant Allele Specific Amplification (MASA) or Amplification Refractory Mutation System (ARMS), and it provides nucleotide sequence information of sample-derived template DNA, based upon the presence or absence of DNA amplification in PCR using a species specific primer pair.

3.1.1. Procedure

The following is an example procedure.

3.1.1.1. Preparation of template DNA

There are various methods with which to extract and purify DNA from the samples. It is recommended that commercially available DNA extraction kits be used when considering their advantages of not using any noxious reagents and not requiring any complicated purification procedures. In this case, attention should be paid to the final amount (concentration) of DNA obtained, and the initial amount of initial sample and the volume of liquid to elute the DNA need to be controlled. When extraction and purification are performed using silica gel membrane type kits stipulated in notifications⁷⁾ related to inspection methods of the foods produced by recombinant DNA techniques, it is appropriate to use 200 mg of sample, 1 mL of AP1 buffer solution, 2 μ L of RNase A, and 325 μ L of AP2 buffer solution. Also, the most important things are that the supernatant loaded on the first column is clear and that there is no need to load 1 mL unreasonably. Furthermore, 50 μL is an appropriate volume used in the final elution of the DNA, and normally the initial eluate is used as the DNA sample stock solution.

3.1.1.2. Confirmation of purity of DNA in DNA sample stock solution and assay of DNA

The purity of the DNA in the stock solution can be confirmed by the $OD_{260 \text{ nm}}/OD_{280 \text{ nm}}$ ratio using a spectrophotometer. A ratio of 1.5 indicates that the DNA has been adequately purified. The amount of DNA is calculated using $1 OD_{260 \text{ nm}} = 50 \,\mu\text{g/mL}$. The measurement mentioned above is performed using appropriately diluted DNA sample stock solution. Based on the results obtained, dilute with water to the concentration needed for the subsequent PCR reactions, dispense the solution into micro tubes as the sample DNA solution, and if necessary store frozen at not over -20°C . The dispensed DNA sample is used immediately after thawing and any remaining solution should be discarded and not refrozen. If the concentration of the DNA sample stock solution does not reach the concentration stipulated in PCR, it is used as a DNA sample solution.

3.1.1.3. PCR

When a commercially available PCR enzyme mentioned in the above notification⁸⁾ is used, it is appropriate that 25 μ L of reaction mixture consisting of 2.5 μ L of the PCR

buffer solution containing magnesium, dNTP (0.2 mmol/L), 5' and 3'primer (0.4 µmol/L), Taq DNA polymerase (1.25 units), and 5 μ L of 10 ng/ μ L sample DNA solution (50 ng of DNA) is prepared on ice. Among them, the PCR buffer solution and dNTP are provided as adjuncts to the enzyme. When conducting purity tests on Atractylodes Lancea Rhizome in Atractylodes Rhizome, the primer sets used are C and D (C is positive with A. lancea, D is positive with A. chinensis) as described in the paper1) mentioned above, however, when primer sets A and B are used, it is possible to confirm the source species of each of the respective specimens. In order to confirm that the DNA has been extracted correctly, the reaction solution containing the positive control primer pair (Pf and Pr) as shown below should be prepared. In addition, the negative control solutions which are not containing DNA sample or either of the primer pair should be prepared and simultaneously conduct PCR.

Pf: 5'-CAT TGT CGA AGC CTG CAC AGC A-3' Pr: 5'-CGA TGC GTG AGC CGA GAT ATC C-3'

The PCR reaction is performed under the following conditions: starting the reaction at 95°C for 10 minutes, 30 cycles of 0.5 minutes at 95°C and 0.75 minutes at 68°C (69°C only when using the primer set C), terminate reaction at 72°C for 7 minutes, and store at 4°C. The resulting reaction mixture is used for the following process as PCR amplification reaction solution.

3.1.1.4. Agarose gel electrophoresis and detection of PCR products

After completion of the PCR reaction, mix $5\,\mu\text{L}$ of the PCR amplification reaction solution with an appropriate volume of gel loading buffer solution, add the mixture to the wells of 2 w/v% agarose gel, and then perform electrophoresis using 1-fold TAE buffer solution (refer to General Information, Rapid Identification of Microorganisms Based on Molecular Biological Method). Carry out the electrophoresis together with an appropriate DNA molecular marker. Electrophoresis is terminated when the bromophenol blue dye in the gel loading buffer has advanced to a point corresponding to 1/2 to 2/3 the length of the gel.

Stain the gel after electrophoresis when not using gel stained in advance with ethidium bromide. Place the gel that has undergone electrophoresis and staining in a gel image analyzer, irradiate with ultraviolet light (312 nm), and detect its electrophoresis pattern. Compare this to the DNA molecular marker and determine the absence or presence of the target amplification band.

3.1.2. Judgment

Confirm at first that a 305 bp band is found with the reaction solution to which the positive control primer pair has been added, and confirm there are no bands in a solution with no primer sets and a solution with no sample DNA solution. Next, if a 226 bp band is confirmed when the primer set C is added or if a 200 bp band is confirmed when the primer set D is added, the sample is judged to be Atractylodes Lancea Rhizome (in the case of cut crude drug, contamination of Atractylodes Lancea Rhizome is observed)

and it is rejected. The sample is judged not to be Atractylodes Lancea Rhizome (in the case of cut crude drug, there is no contamination of Atractylodes Lancea Rhizome) and the purity test is acceptable if a 305 bp band is confirmed with the positive control primer pair, bands are not observed in reaction solution without primer and reaction solution without DNA sample solution, and a 226 bp band is not observed with the primer set C and a 200 bp band is not observed with the primer set D. If a band is not observed with the positive control primer pair, it is to be concluded that the DNA extraction failed and the procedure should be started over again from the DNA extraction step. If bands are confirmed in reaction solutions without primer sets or without DNA sample solution, it should be assumed that there was an error in the PCR procedure and therefore the procedure should be repeated again from the step 3.1.1.3. PCR.

3.2. Method 2 (PCR—Restriction Fragment Length Polymorphism)

Generally, this method is referred to as PCR—Restriction Fragment Length Polymorphism (RFLP), and it provides nucleotide sequence information of sample-derived template DNA, based upon the DNA fragment pattern produced by restriction enzyme treatment of the PCR products, which are amplified by using a primer pair common to the DNA sequence of the objective plant.

The test is performed with 25 samples randomly taken from a lot, and each sample is designated with a number from 1 to 25. Differentiation of the sources is performed by individual PCR—RFLP measurement of the samples, and decision of the acceptability of the purity is dependent on how many nonconforming samples are present in the first 20 samples, taken in numerical order, for which judgment is possible as described below.

3.2.1. Procedure

The following is an example procedure.

3.2.1.1. Preparation of template DNA

There are various methods with which to extract and purify DNA from the samples. It is recommended that commercially available DNA extraction kits be used, when considering their advantages of not using noxious reagents and not requiring complicated purification procedures. Recently, PCR reagents that inhibit the effect of PCR enzymeinhibiting substances present in samples have become commercially available, and by using these reagents, it is possible to prepare the template DNA from the sample simply by incubating the sample with the DNA extraction reagent. Here, a recommended DNA preparing procedure using such PCR reagents is described for the convenience of experimenters.

Cut 20 mg of the sample into small pieces with a clean knife, add $400 \,\mu\text{L}$ of the DNA extraction reagent, and incubate at $55\,^{\circ}\text{C}$ overnight (16 – 18 hours). Then heat at $95\,^{\circ}\text{C}$ for 5 minutes to inactivate the enzyme in the reagent. Centrifuge to precipitate the sample, and use $50\,\mu\text{L}$ of the supernatant liquid as the template DNA solution. The DNA solution prepared in this method can not be used for con-

centration measurement based on $OD_{260\,nm}$, because it contains many foreign substances affecting $OD_{260\,nm}$ value from the sample.

The composition of the DNA extraction reagent is as follows:

2-Amino-2-hydroxymethyl-1,3propanediol-hydrochloric acid (pH 8.0) 20 mmol/L Ethylenediamine tetraacetate 5 mmol/L Sodium chloride 400 mmol/L Sodium dodecyl sulfate 0.3% Proteinase K 200 μ g/mL

3.2.1.2. PCR

In the method using the PCR enzyme and PCR reagent as described³⁾, the reaction mixture is prepared on an ice bath in a total volume of 20 μ L of a solution containing 10.0 μ L of 2-fold concentrated PCR reagent, 5'- and 3'-primers (0.5 μ mol/L), Taq DNA polymerase (0.5 units) and 0.5 μ L of template DNA solution.

The PCR reaction is performed under the following conditions: 95°C for 10 minutes, 40 cycles of 95°C for 0.5 minute, 65°C for 0.25 minute, and 72°C for 0.25 minute and 72°C for 7 minutes. Store the solution at 4°C, and use this solution as the PCR amplified solution. A negative control (containing water instead of the template DNA solution) must be included in the procedure.

The sequence of each primer is as follows:

5'-primer: 5'-GGC ACA ACA CGT GCC AAG GAA
AA-3'

3'-primer: 5'-CGA TGC GTG AGC CGA GAT ATC C-3'

3.2.1.3. Restriction enzyme treatment

The treatment is performed on individual reaction solutions using two enzymes, Fau I and Msp I. In the case of Fau I, to an appropriate amount of the reaction solution, composed of a reaction buffer containing 1.0 unit of enzyme, add $3.0 \,\mu\text{L}$ of PCR products while cooling in an ice bath to make $15.0 \,\mu\text{L}$. In the case of Msp I, to an appropriate amount of the reaction solution, composed of a reaction buffer containing 20.0 units of enzyme, add $3.0 \,\mu\text{L}$ of PCR products while cooling in an ice bath to make $15.0 \,\mu\text{L}$. Incubate these solutions at the temperature recommended by the manufacturer for 2 hours, and then inactivate the enzyme by heating at 72°C for 10 minutes. The negative control of the PCR reaction is also treated in the same manner.

3.2.1.4. Agarose gel electrophoresis and detection of DNA fragments

After the restriction enzyme treatment, mix the total amount of the reaction solution and an appropriate amount of the gel loading buffer solution, place it in a 4 w/v% agarose gel well, and carry out electrophoresis with 1-fold concentrated TAE buffer solution (see "Rapid Identification of Microorganisms Based on Molecular Biological Methods" under General Information). Carry out the electrophoresis together with appropriate DNA molecular marker. Stop the electrophoresis when the bromophenol blue included in the loading buffer solution has moved about 2 cm from the well. The 4 w/v% agarose gel is sticky, difficult

to prepare and hard to handle, so that it is better to use a commercially available precast gel.

After the electrophoresis, stain the gel, if it is not already stained, with ethidium bromide, and observe the gel on an illuminating device under ultraviolet light (312 nm) to confirm the electrophoretic pattern.

3.2.2. Judgment

3.2.2.1. Judgment of each sample

Confirm that no band is obtained with the negative control of the PCR, other than the primer dimer (about 40 bp) band. A sample treated with *Fau* I, showing bands of about 80 bp and 60 bp, or that treated with *Msp* I, showing bands of about 90 bp and 50 bp, is judged as Atractylodes Lancea Rhizome. A sample not showing any band other than a band at about 140 bp and the primer dimer band is judged as Atractylodes Rhizome. If a sample does not show any band other than the primer dimer band, it is considered that PCR products were not obtained, and judgment is impossible for the sample.

3.2.2.2. Judgment of the purity

Judgment of the purity is based on the result of the judgment of each sample. If there is no sample that is judged as Atractylodes Lancea Rhizome among 20 samples taken in order of the numbering, excluding any sample for which judgment is impossible, the lot is acceptable for purity. When there is one sample that is judged as Atractylodes Lancea Rhizome among the 20 samples, perform the same test with 25 newly taken samples from the lot, and if there is no sample that is judged as Atractylodes Lancea Rhizome, the lot is acceptable for purity. When there is a sample that is judged as Atractylodes Lancea Rhizome in the second test, or there is more than one sample that is judged as Atractylodes Lancea Rhizome in the first test, the lot is not acceptable for purity.

4. Purity test of Saposhnikovia Root and Rhizome for Peucedanum ledebourielloides

4.1. Method 1

Similarly as 3.1., this method provides nucleotide sequence information of sample-derived template DNA, based upon the presence or absence of DNA amplification band in PCR using a species specific primer pair.

4.1.1. Procedure

The following is an example procedure.

4.1.1.1. Preparation of template DNA

For Atractylodes crude drugs, a preparation procedure using a silica gel membrane type kit is adopted, however for the test of Saposhnikovia Root and Rhizome, and *Peucedanum ledebourielloides*, the simple preparation procedure shown below is adopted for the convenience of experimenters, because it was confirmed that the PCR product is stably obtained in using a DNA sample solution prepared by the simple preparation procedure shown in 3.2.1.1. as a template.

Cut 10 mg of the sample into small pieces with a clean knife, add 400 μ L of the DNA extraction reagent, and incubate at 55°C overnight (16 - 18 hours). Then heat at 95°C for 5 minutes to inactivate the enzyme in the reagent. Cen-

trifuge to precipitate the sample, and use $50 \,\mu L$ of the supernatant liquid as the template DNA solution. The DNA solution prepared in this method can not be used for concentration measurement based on OD_{260nm} , because it contains many foreign substances affecting OD_{260nm} value from the sample.

The composition of the DNA extraction reagent is as follows:

2-Amino-2-hydroxymethyl-1,3-

propanediol-hydrochloric acid (pH 8.0) 20 mmol/L Ethylenediamine tetraacetate 5 mmol/L Sodium chloride 400 mmol/L Sodium dodecyl sulfate 0.3% Proteinase K 200 μ g/mL

4.1.1.2. PCR

In the method using the PCR enzyme and PCR reagent as described³⁾, the reaction mixture is prepared on an ice bath in a total volume of 20 μ L of a solution containing 10.0 μ L of 2-fold concentrated PCR reagent, 5'- and 3'-primers (0.5 μ mol/L), Taq DNA polymerase (0.5 units) and 0.5 μ L of template DNA solution.

When the purity test of Saposhnikovia Root and Rhizome for *Peucedanum ledebourielloides* is performed, the reaction solution containing the positive control primer pair as shown below should be prepared besides the reaction solution containing a species specific primer pair in order to confirm that the DNA has been extracted correctly. In addition, the negative control solutions which are not containing DNA sample should be prepared and simultaneously conduct PCR.

The PCR reaction is performed under the following conditions: 95°C for 10 minutes, 45 cycles of 95°C for 0.5 minute, 62°C for 0.5 minute, and 72°C for 0.75 minute and 72°C for 7 minutes. Store the solution at 4°C, and use this solution as the PCR amplified solution. The sequence of each primer is as follows. The positive control 3'-primer for PCR and the species specific 3'-primer for PCR have the same sequence.

5'-primer for positive control PCR: 5'-GCG TGG GTG TCA CGC ATC G-3'

3'-primer for positive control PCR: 5'-GTA GTC CCG CCT GAC CTG-3'

5'-primer for species specific PCR: 5'-CTG AGA AGT TGT GCC CGG-3'

3'-primer for species specific PCR: 5'-GTA GTC CCG CCT GAC CTG-3'

4.1.1.3. Agarose gel electrophoresis and detection of PCR products

After completion of PCR reaction, mix 5 μ L of the PCR amplification reaction solution with an appropriate volume of gel loading buffer solution, add the mixture to the wells of 2 w/v% agarose gel, and then perform electrophoresis using 1-fold TAE buffer solution (refer to General Information, Rapid Identification of Microorganisms Based on Molecular Biological Method). Carry out the electrophoresis together with an appropriate DNA molecular marker. Electrophoresis is terminated when the bromophenol blue

dye in the gel loading buffer has moved about 2 cm from the well.

Stain the gel after electrophoresis when not using gel stained in advance with ethidium bromide. Place the gel that has undergone electrophoresis and staining in a gel image analyzer, irradiate with ultraviolet light (312 nm), and confirm its electrophoresis pattern. Compare this to the DNA molecular marker and determine the absence or presence of the target amplification band.

4.2. Judgment

Confirm at first that a 250 bp band is found with the reaction solution to which the positive control primer pair has been added, and confirm there are no bands other than the primer dimer (about 40 bp) in a solution with no sample DNA solution. Next, if a 200 bp band is confirmed when the species specific primer pair is added, the sample is judged to be contaminated with Peucedanum ledebourielloides and it is rejected. The sample is judged not to be contaminated with Peucedanum ledebourielloides and the purity test is acceptable if a 250 bp band is confirmed with the positive control primer pair, bands are not observed in reaction solution without DNA sample solution, and a 200 bp band is not observed with the species specific primer pair. If a band is not observed with the positive control primer pair, it is to be concluded that the DNA extraction failed and the procedure should be started over again from the DNA extraction step. If bands are confirmed in reaction solution without DNA sample solution, it should be assumed that there was an error in the PCR procedure and therefore the procedure should be repeated again from the step 4.1.1.2. PCR.

5. Reference

- 1) Y. Guo, et al., J. Nat. Med. 60, 149-156 (2006).
- 2) K. Kondo, et al., J. Jpn. Bot. 84, 356-359 (2009).
- T. Maruyama, et al., Shoyakugaku Zasshi 64, 96-101 (2010).
- 4) T. Maruyama, et al., J. Nat. Med. in press.
- 5) JAS analytical test handbook: genetically modified food quality, labeling analysis manual for individual products, ver. 3, I Fundamental procedure 4.4.1 PCR, (September 24, 2012). Incorporated Administrative Agency Food and Agricultural Materials Inspection Center
- 6) JAS analytical test handbook: genetically modified food quality, labeling analysis manual for individual products, ver. 3, IV Contamination prevention, (September 24, 2012). Food and Agricultural Materials Inspection Center Independent Administrative Organization
- Notification No. 110, Director of Food Health Department, March 2001; Partial Amendment: Notification No. 0629002, 2.2.1.2, Director of Food Safety Department, June 2006.
- 8) Notification No. 0629002, 2.1.3.1.1, Director of Food Safety Department, June 2006.

On the Scientific Names of Crude Drugs listed in the JP

Change the following as follows:

Scientific Names used in the JP and Those being used Taxonomically

Crude Drug	Scientific names used in the JP = Scientific names being used taxonomically (Combined notation, Standard form for author or authors)	Family
	Scientific names that are different from those written in JP but identical to them taxonomically or being regarded as identical, and typical subclassified groups belonging to their species. The names marked with "*" are those being written together in JP.	
Amomum Seed シュクシャ	Amomum villosum Loureiro var. xanthioides T. L. Wu & S. J. Chen = Amomum villosum Lour. var. xanthioides (Wall. ex Baker) T. L. Wu & S. J. Chen	
	Amomum xanthioides Wallich = Amomum xanthioides Wall. ex Baker	Zingiberaceae
	Amomum villosum Lour. var. nanum H. T. Tsai & S. W. Zhao	
	Amomum villosum Loureiro var. villosum = Amomum villosum Lour. var. villosum	
	Amomum villosum Lour.	
	Amomum longiligulare T. L. Wu	

G8 Water

Quality Control of Water for Pharmaceutical Use

Change the 4.4.2. Media Growth Promotion Test as follows:

4.4.2. Media Growth Promotion Test

In the media growth promotion test with the R2A Agar Medium, use test strains listed below or other test strains considered equivalent to these test strains.

Methylobacterium extorquens: NBRC 15911
Pseudomonas fluorescens: NBRC 15842, ATCC 17386,

Prior to the test, inoculate these test strains into sterile purified water and starve them at 20 – 25°C for 3 days.

Dilute the fluid containing the test strain starved with sterile purified water to prepare microbial suspensions. When inoculating the R2A Agar Medium with the microorganisms ($5 \times 10^1 - 2 \times 10^2$ CFU) and incubating at 20 - 25°C for 4 - 7 days, sufficient proliferation of the inoculated strains must be observed.

In the media growth promotion test with the Standard Agar Medium, use test strains listed below or other test strains considered equivalent to these test strains.

Staphylococcus aureus: ATCC 6538, NCIMB 9518, CIP 4.83 or NBRC 13276

Pseudomonas aeruginosa: ATCC 9027, NCIMB 8626, CIP 82.118 or NBRC 13275

Escherichia coli: ATCC 8739, NCIMB 8545, CIP 53.126 or NBRC 3972

Prepare microbial suspensions containing the test strains according to the procedure prescribed in the Microbiological Examination of Non-sterile Products <4.05>. When inoculating the Standard Agar Medium with a small number (not more than 100 CFU) of the micro-organisms and incubating at 30 – 35°C for 48 hours, sufficient proliferation of the inoculated strains must be observed.

G10 Others

Change as follows:

Basic Concepts for Quality Assurance of Drug Substances and Drug Products

Introduction

Quality of drug substances and products are generally assured through manufacturing and testing under appropriate Good Manufacturing Practice (GMP) conditions reflecting

knowledge obtained from designing and developmental stages and manufacturing stage on management of raw materials and other materials, control of manufacturing process, specifications, etc. As shown in the General Notice 5, JP listed drugs are to be tested according to the provisions given in the pertinent monographs, General Notices, General Rules for Crude Drugs, General Rules for Preparations, and General Tests for their conformity to the *Japanese Pharmacopoeia*. In addition to these, compliance with GMP, management of raw materials and other materials, and management of manufacturing process are fundamental factors required to assure the quality of JP listed products in actual production.

The present chapter summarizes general concepts concerning measures for quality assurance of drug substances and products mainly aimed at chemicals, including chemically synthesized antibiotics and semisynthetic antibiotics, synthetic peptides, oligonucleotides, and biotechnological/biological products, and shows the principle idea of quality assurance in the process listing a drug as an individual monograph in the JP. Although radiopharmaceuticals, crude drugs, herbal products, and crude products of animal or plant origin are excluded from the subjects of the concepts, these concepts are useful for the management of any type of drugs.

Basic Concept

In recent years, the mainstream concept for quality control of drugs has been implemented according to a control strategy that their quality is assured by control of manufacturing process, including management of raw material and other materials, and quality testing of final products (drug substances or drug products) that are conducted mutually complementary. The control strategy is implemented based on Quality Risk Management (QRM). The first and most important step is identifying Critical Quality Attributes (CQAs) which are the attributes or properties required to ensure the desired product quality, and it is necessary to specify physical, chemical, biological, microbiological characteristics or properties of the product which should be within the appropriate limits, ranges and distributions. The next step is to guarantee that the CQA falls within the defined range, limit, and distribution by using specification tests, in-process tests and various measures, for that the quality of the drug will eventually be realized.

The Specification is one of the elements of control strategy and not all the CQA need to be included in the specifications. CQA is (1) included in specifications and confirmed by testing final products (including periodical or skip testing, described later), (2) included in specifications and confirmed by process controls (e.g., real time release testing, described later), or (3) not included in specifications but can be ensured by controlling starting materials, raw materials and manufacturing process. As an example of (3), effective control over robust manufacturing processes can assure that certain impurities are controlled at an acceptable risk level or are efficiently removed below an acceptable level, and sometimes the purity testing for the final product may not

be required and omitted from specifications. However, in the case of a drug listed in the JP monograph, regarding the manufacturing process control related to CQA, if necessary, the control method and control value are indicated in the Manufacture in individual monograph.

What kind of control strategy should be applied to a certain CQA is individually determined by QRM according to the understanding and risk of the manufacturing process.

1. Management of manufacturing process

1.1. Considerations of manufacturing process

Adequate design of manufacturing processes and knowledge of their capacity are important to establish manufacturing processes yielding drug substances or drug products that meet specifications and fulfill CQA, and to perform consistent manufacturing control, quality control, etc. appropriately.

From this standpoint, the limits for control of manufacturing processes should be based on information obtained from the entire process spanning the period from the early development through commercial scale production. The appropriateness of the limits also need to be confirmed by evaluation, verification, review, and other examinations of manufacturing processes based on QRM.

In-process tests are tests that may be performed during the manufacture of either the drug substance or drug product, rather than specification tests for the final product. In-process tests are performed for quality verification during manufacturing processes that are likely to influence drug substance or drug product quality, or for confirmation of proper functioning of the manufacturing process. Inprocess tests may also be used for the evaluation of CQA.

Usually an in-process test is properly designed according to the risk on quality, however, the use of internal action limits by the manufacturer to assess the consistency of the process at less critical steps is also important. Provisional action limits should be set for the manufacturing process based on data obtained during development of the drug and during evaluation and verification of the manufacturing process, and should be further refined based on additional manufacturing experience and data accumulated after product approval for marketing.

1.2. Considerations of raw materials and other materials (starting materials, excipients, packaging materials, etc.)

The raw materials and other materials used in the production of drug substances (or drug products) should meet quality standards, appropriate for their intended use, and appropriate setting of specifications and test methods assuring CQA are required. Especially, biological raw/source materials may require careful evaluation to establish the presence or the absence of deleterious endogenous or adventitious agents. Procedures that make use of affinity chromatography (for example, employing monoclonal antibodies), should be accompanied by appropriate risk management to ensure that such process-related impurities or potential contaminants arising from their production and use do not compromise the quality and safety of the drug substance or drug product.

The quality of the excipients used in the drug product formulation (and in some cases, in the production of drug substance), as well as the primary packaging materials, should be controlled with specifications established based on the characteristics of the drug. If specifications and test procedures for a material are described by the JP, as a rule, at least the JP criteria should be satisfied. Concerning excipients and other materials not listed in the JP, appropriate specifications and test procedures should be established individually.

2. Quality tests of products (specifications)

A specification is defined as a list of tests, references to analytical procedures, and appropriate acceptance criteria which are numerical limits, ranges, or other criteria for the tests described. Specifications and test methods of the JP monograph are defined sets of quality characteristics needed for determination of whether the use of a drug substance or a drug product is appropriate for the intended purpose. "Conformance to the specifications of the JP monograph" means that the JP-listed drug substances and drug products, when tested according to the procedures described in general tests and drug monographs, will meet the all acceptance criteria except criteria of "Description", "Containers and storage (for drug products)" and "Shelf life" in the JP monographs.

However, as described in "Basic Concept" specifications of monographs and test procedures for drug substance and drug product are one part of a total control strategy for assurance of the quality and consistency of the substances/products. Other parts of control strategies include thorough characterization of the drug in developmental stage (specifications and test procedures are established based on the characterization), and management of manufacturing process and products' quality, such as evaluation, verification and review of manufacturing process, and management of raw materials, other materials and manufacturing process, that is to say, compliance with the GMP.

3. Periodic or Skip Testing

Periodic or skip testing is the performance of specified tests at release on preselected batches and/or at predetermined intervals, rather than on a batch-to-batch basis with the understanding that those batches not being tested still must meet all acceptance criteria established for that product. This concept may be applicable to, for example, residual solvents and microbiological testing for solid oral dosage forms. It is recognized that only limited data may be available at the time of submission of an application for marketing approval. Implementation of this concept should therefore generally be considered post-approval. When tested, any failure to meet acceptance criteria established for the periodic test should be handled by proper notification of the appropriate regulatory authorities. If these data demonstrate a need of testing for all lots, then batch-bybatch release testing should be reinstated.

4. Real-time release testing (RTRT) and parametric release RTRT is a type of tests to evaluate the quality of in-

process or final products based on process data (including results of in-process testing and data on process parameters) and to assure that the quality is acceptable. RTRT is a kind of specifications and consists of a valid combination of materials attribute (intermediate products) pre-evaluated and process control. RTRT is used for judgement of products release instead of the release testing of final products when the application containing RTRT is approved by a regulatory authority.

The usage of RTRT does not mean unnecessity of setting tests of a final product directly. Even if the decision of release is made by RTRT, the tests for final products need to be set as specifications. It is because final product testing may be requested for some reasons such as failure of data acquisition due to troubles of equipments used for RTRT and evaluation of stability of final products. The final products, of course, need to meet their specifications, when tested.

Likewise, in the case that the drugs that was approved for marketing with the RTRT is listed in the JP monograph, the RTRT can be continued to use for release judgement. However specification and test procedure that assure the quality as same as the RTRT for final products should be set in the monograph. Even for drugs whose specifications are listed in the monographs, when a new application (or application for partial change) containing RTRT is approved by the regulatory authority, the products release can be judged based on the results of the RTRT instead of the tests prescribed in the monograph. In addition, it is necessary to comply with the specification in the case of conducting the compendial tests. In either case, it is unnecessary to set specification for RTRT in "Manufacture" of the monograph since the control criteria for the target CQA is already shown for the RTRT.

If RTRT results fail or trend toward failure, RTRT should not easily be substituted by final product testing. In this case, it is important to investigate the cause properly and need to take corrective action. Also, if RTRT results fail, the products cannot be released unless they were caused by analysis failure such as equipment failure. If RTRT results are trending toward failure, the products release should be made carefully based on the results of the investigation.

Parametric release can be considered a type of real time release. One example of parametric release is to determine the suitability for release of terminally sterilized drug products based on the data on sterilizing process instead of the results of sterility testing. In this case, the release of each batch is based on satisfactory results from monitoring specific parameters, e.g., temperature, pressure, and time during the terminal sterilization phase(s) of drug product manufacturing. Parametric release based on above parameters is more reliable in predicting sterility assurance than is determination of suitability for release based on sterility testing using limited number of final products. Besides, even if parametric release is applied, the final product testing need to be set because the testing is necessary in stability testing and post-marketing surveillance. If in-process

data used for parametric release are not acceptable, the products cannot be released. The parametric release differs from RTRT in the case, for example where the data of monitoring specific parameters in terminally sterilized process is failed to obtain by a certain reason such as analysis failure by equipment failure and so on. The incomplete data acquisition means no assurance on sterilization process, it is impossible to substitute parametric release by sterility testing of final products in principle.

Add the following:

Concept on Impurities in Chemically synthesized Drug Substances and Drug Products

1. Classification of impurities found in chemically synthesized pharmaceuticals and the guidance to comply with for their control

Impurities found in chemically synthesized pharmaceuticals are roughly classified into organic impurities, inorganic impurities and residual solvents. Those impurities in the new drug substances and the products are controlled by the following guidelines agreed upon at the International Council for Harmonisation of Technical Requirements for Pharmaceuticals for Human Use (hereinafter referred to as "ICH"). More specifically, "Impurities in New Drugs Substances (PAB/PCD Notification No. 877 dated September 25, 1995)" (hereinafter referred to as "ICH Q3A Guideline")1) on specifications for organic impurities in drug substances applies to applications for marketing approval after April 1, 1997, while "Impurities in New Drug Products (PAB/PCD Notification No. 539 dated June 23, 1997" (hereinafter referred to as "ICH Q3B Guideline")2) on specifications for organic impurities in drug products applies to applications for marketing approval after April 1, 1999. Meanwhile, specifications for inorganic impurities were specified by Japanese pharmacopoeial standards and known safety data. Now "Guidelines for Elemental Impurities (PFSB/ELD Notification No. 4 dated September 30, 2015)" apply to applications for marketing approval after April 1, 2017. In regard to residual solvents, "Impurities: Guidelines for Residual Solvents (PAB/ELD Notification No. 307 dated March 30, 1998)" (hereinafter referred to as "ICH Q3C Guideline") applies to applications for marketing approval after April 1, 2000. Especially in regard to DNA-reactive impurities, "Assessment and control of DNA reactive (mutagenic) impurities in pharmaceuticals to limit potential carcinogenic risk (PSEHB/ELD Notification No. 3 dated November 10, 2015)" applies to applications for marketing approval after January 15, 2016. Although ICH Q3A guideline does not cover optical enantiomers, a type of organic impurities, "Specifications: Test Procedures and Acceptance Criteria for New Drug Substances and New Drug Products: Chemical Substances (PMSB/ELD Notification No. 568 dated May 1, 2001)" (hereinafter referred to as "ICH Q6A Guideline"), which was published subsequently, provides that enantiomers are impurities that should be controlled and, if measurable, should be controlled in accordance with the principle of the ICH Q3A guideline.

Control of impurities in accordance with the guidelines mentioned above is expected also for pharmaceuticals other than new drug substances and new drug products. Their applications for marketing (or applications for partial changes) are subject to those guidelines when necessary. The General Notices of the JP 17th Edition states that residual solvents of all JP-listed drugs, in principle, have to be controlled in accordance with specification "Residual Solvents" in General Tests unless otherwise specified in the individual monograph. In regard to elemental impurities, it has been decided in the basic principles for the preparation of the JP 18th Edition to create a roadmap for their incorporation into the JP for listing and to address its implementation.

2. The concept of ICH Q3A and Q3B guidelines for the control of organic impurities

ICH Q3A and Q3B guidelines require setting acceptance criteria for organic impurities based on the information gained from development stages for new drugs. Concerning impurities in drug substances, ICH Q3A guideline refers to the items to be examined from chemical and safety perspectives. ICH Q3B guideline complements Q3A guideline, and have the same basic concept as Q3A. Chemical aspects to be examined include classification and identification of impurities, their reporting method, specification settings and analytical methods. Safety aspects include specific guidelines for qualifying the safety of impurities that were not present, or were present at substantially lower levels, in batches of a drug substance used in safety and clinical studies.

Qualification of the safety is the process of acquiring and evaluating data that establishes the biological safety of an individual impurity or a given impurity profile at the level(s) specified. The applicant should describe a rationale for establishing impurity acceptance criteria that includes safety considerations in attachments when applicated for approval. The level of any impurities present in a new drug substance that has been adequately tested in safety and/or clinical studies would be considered qualified.

Identified impurities, unidentified impurities and total impurities are specified based on the data obtained according to the guidelines. The threshold of unspecified impurities in a drug substance is determined depending on the daily intake of the drug substance. When the maximum daily dosage is not more than 2 g, it is set at 0.10%. The establishment of individual specifications is required for impurities at a level greater than 0.10%.

In regard to drug products, the ICH Q3B guideline cover the degradation products of drug substances or reaction products between the drug substance and additive/primary packaging. Therefore, even if organic impurities other than degradation products (e.g., by-products and synthetic intermediates) in the drug substance are found as impurities in the drug product, they need not be monitored or specified since they have already been controlled as the drug sub-

stance specifications. However, degradation products elevated in the drug product need to be monitored and specified.

3. Principles for controlling organic impurities in the articles listed in the JP

Conventionally in the JP, specified impurities, unspecified impurities and total impurities are specified in accordance with ICH Q3A and Q3B guidelines for pharmaceutical products, whose impurities have been controlled by those guidelines, in the process listing in the JP. (However, this shall not apply to the long-term listed pharmaceutical products which had existed in the JP before these guidelines were applicable. However, when a new application is filed for those JP-listed pharmaceutical products, control of impurities in accordance with ICH Q3A and Q3B guidelines may be required, if necessary.) In order to specify the impurities, analysis data during development submitted from the drafting company and impurity analysis data from commercial production batches after consistent manufacturing is achieved should be assessed. Safety evaluation is not required again for the process listing in the JP since it has been performed at the time of approval.

ICH Q3A and Q3B guidelines cover impurities in the drug substances manufactured by chemical syntheses and the drug products manufactured with those drug substances. Similarly, the following types of products are not covered in the JP: biological/biotechnological products, peptides, oligonucleotides, radiopharmaceuticals, fermentation products and semi-synthetic products derived therefrom, herbal products and crude products of animal or plant origin.

When organic impurities assessed in accordance with the principles of ICH Q3A and Q3B guidelines are listed as JP tests of purity, the operational rationality of the JP is considered and its own modification is added. (i) Except in exceptional circumstances, impurity reference standards are not established. In order to identify an impurity using liquid chromatography, the relative retention time of the impurity to the drug substance is used for identification. (ii) When only unidentified impurities in highly pure pharmaceutical products (not more than 0.1%) are specified, it is generally exempted to set acceptance criteria for total impurities. (iii) When acceptance criteria set based only on actual measured values result in many impurities with slightly different acceptance criteria, consideration can be given so that the purity test consists of a small number of representative acceptance criteria, if possible. (iv) Chemical structural information and the chemical name of the impurities are not disclosed. Those measures enable impurity control without impurity reference standards, and can simplify system suitability test for highly pure pharmaceutical products.

Meanwhile, the method to identify impurities by use of relative retention time is column-dependent and analysis becomes difficult when appropriate columns are not available. Therefore, the JP 17th Edition also allows the use of the analysis method with impurity reference standards when designing purity tests for a drug substance. In addition, the

JP adopted a policy to disclose chemical names and structure formulas as the information on impurities including, in principle, optical enantiomers.

The JP-specific consideration may be given to purity tests for organic impurities in drug products in the process listing in the JP. Also in the JP, impurities derived from the products of the reaction between the drug substance and additive/primary packaging are specified as impurities in the drug product. Those impurities are formulation-dependent and may not be formed in different formulations. Since the JP is an official compendium that allows a wide variety of formulations, when it is not appropriate to specify impurities uniformly in the individual monograph, they are subject to the specifications at the time of approval, along with the statement "Being specified separately when the drug is granted approval based on the Law."

When the specifications for impurities are reviewed for a new entry of a pharmaceutical product in an individual monograph of the JP, acceptance criteria for impurities may be included in the review according to the following concepts. ICH Q6A guideline point out: Data available upon the marketing application are limited and it has to be taken into consideration that the limited data may influence the design of acceptance criteria. Regarding impurities, since impurity profiles gained during the manufacturing stages may sometimes be different from that gained from development stage, it is stated that changes in impurity profiles at the manufacturing stage should be considered as appropriate. According to this concept, for impurities which should be specified in the process listing in the JP, not only information from development stage but also information about impurity profiles if there are changes at the manufacturing stage, and information at the stage after the product manufacturing becomes stable (hereinafter referred to as the "stable production stage") should be taken into consideration.

However, it is undesirable to remove impurities that are present at substantially lower levels, or become undetectable at the stable production stage indiscriminately from the list of candidate compounds to be specified. JP-listed drugs are accepted as drugs by conformance to the specifications in the individual monograph. However, generic drugs, whose manufacturing methods are not necessarily the same as that of the drug substance used for JP monograph, may have different impurity profiles and contain such impurities. Providing information in the process listing in the JP based on the detection results during development stage may result in encompassing impurities found in drug substances and drug products distributed as JP drugs.

Therefore, before the removal of impurities that are present at substantially lower levels or become undetectable at the stable production stage from the JP specification list, the need to establish specifications should be fully examined based on ICH Q3A and Q3B guidelines with respect to safety

For a drug substance that was approved by the method to identify its impurities with impurity reference materials, it is desirable also in the individual JP monograph, in principle, to establish specifications and test methods appropriately so that the specified impurity becomes identifiable. In regard to impurity control during the manufacturing process, impurities can be controlled by establishing an appropriate control strategy including release testing, in-process tests and process parameters control.

References

- 1) ICH:Harmonised Tripartite Guideline, Impurities in New Drug Substances.
- ICH:Harmonised Tripartite Guideline, Impurities in New Drug Products.

Add the following:

Glossary for Quality by Design (QbD), Quality Risk Management (QRM) and Pharmaceutical Quality System (PQS)

1. Introduction

The purpose of this glossary is to define terms, used for developing the new concept of quality assurance in ICH Q8-11 guidelines so-called Q quartet, and to explain the concept. The terms shown here are determined as the result of discussion for long time in ICH, and are most important to understand the concept of systematic quality assurance based on science and quality risk management, as shown by the guidelines. The usage may not necessarily accord with general usage, however it is necessary to keep in mind that the following definition is used in the regulatory application of pharmaceuticals. The terms used in ICH Q8 to Q11 are shown below in their order. For terms explained in more than one guideline, the name of duplicated guideline is described in parentheses at the end of the corresponding sentence.

2. Glossary

[ICH Q8 Guideline]

Control Strategy: A planned set of controls, derived from current product and process understanding, that ensures process performance and product quality. The controls can include parameters and attributes related to drug substance and drug product materials and components, facility and equipment operating conditions, in-process controls, finished product specifications, and the associated methods and frequency of monitoring and control (ICH Q10, Q11). A control strategy is expected irrespective of development approaches. Under the development approach using Quality by Design, testing, monitoring or controlling can be shifted earlier into the process.

Quality by Design (QbD): A systematic approach to development that begins with predefined objectives and emphasizes product and process understanding and process control, based on sound science and quality risk management.

Continuous Process Verification: An alternative approach to process validation in which manufacturing process performance is continuously monitored and evaluated. Process validation protocol can use Continuous Process Verification (CPV) to the process validation protocol for the initial and ongoing commercial production (ICH Q11). Generally, for initial process validation, CPV is more appropriate when QbD approach has been applied. However, it can also be used when extensive process knowledge has been gained through commercial manufacturing experience.

Process Robustness: Ability of a process to tolerate variability of materials and changes of the process and equipment without negative impact on quality (ICH Q11).

Critical Process Parameter (CPP): A process parameter whose variability has an impact on a critical quality attribute and therefore should be monitored or controlled to ensure the process produces the desired quality.

Critical Quality Attribute (CQA): A physical, chemical, biological or microbiological property or characteristic that should be within an appropriate limit, range, or distribution to ensure the desired product quality (ICH Q11). For example, CQAs of solid oral dosage forms are typically those aspects affecting product purity, strength, drug release and stability as described in ICH Q8, however it is usual to include product purity and strength itself in CQAs.

Formal Experimental Design: A structured, organized method for determining the relationship between factors affecting a process and the output of that process. Also known as "Design of Experiments" (DoE). The factors to be studied in a DoE could come from the risk assessment exercise or prior knowledge.

Design Space (DS): The multidimensional combination and interaction of input variables (e.g., material attributes) and process parameters that have been demonstrated to provide assurance of quality. Working within the design space is not considered as a change. Movement out of the design space is considered to be a change and would normally initiate a regulatory post approval change process. Design space is proposed by the applicant and is subject to regulatory assessment and approval (ICH Q10, Q11). Design space can be updated over the lifecycle as additional knowledge is gained. Since Proven Acceptable Range (PAR) from only univariate experimentation may lack an understanding of interactions between process parameters and/or material attributes, it should be noted that a combination of PAR does not constitute a design space.

Quality: The degree to which a set of inherent properties of a product, system or process fulfills requirements (ICH Q6A, Q8, Q10). The suitability of either a drug substance or a drug product for its intended use. This term includes such attributes as the identity, strength, and purity (ICH Q6A, Q8, Q9, Q10).

Process Analytical Technology (PAT): A system for designing, analyzing, and controlling manufacturing through timely measurements (i.e., during processing) of critical quality and performance attributes of raw and in-process materials and processes with the goal of ensuring final product quality.

Quality Target Product Profile (QTPP): A prospective summary of the quality characteristics of a drug product that ideally will be achieved to ensure the desired quality, taking into account safety and efficacy of the drug product. Quality Target Product Profile describes the design criteria for the product, and should therefore form the basis for development of the product (ICH Q8).

Lifecycle: All phases in the life of a product from the initial development through marketing until the product's discontinuation (ICH Q11).

Real Time Release Testing (RTRT): The ability to evaluate and ensure the quality of in-process and/or final product based on process data, which typically include a valid combination of measured material attributes and process controls (ICHQ11). Parametric release is one type of Real Time Release Testing. It is based on process data rather than testing of material and/or a sample for a specific attribute. For details, refer to "Basic Concepts for Quality Assurance of Drug Substances and Drug Products" in General Informa-

Proven Acceptable Range (PAR): A characterized range of a process parameter for which operation within this range, while keeping other parameters constant, will result in producing a material meeting relevant quality criteria.

[ICH Q9 Guideline]

Decision Maker(s): Person(s) with the competence and authority to make appropriate and timely quality risk management decisions.

Harm: Damage to health, including the damage that can occur from loss of product quality or availability.

Trend: A statistical term referring to the direction or rate of change of a variable(s).

Detectability: The ability to discover or determine the existence, presence, or fact of a hazard.

Severity: A measure of the possible consequences of a haz-

Product Lifecycle: All phases in the life of the product from the initial development through marketing until the product's discontinuation.

Hazard: The potential source of harm (ISO/IEC Guide 51).

Quality System: The sum of all aspects of a system that implements quality policy and ensures that quality objectives are met.

Requirements: The explicit or implicit needs or expectations of the patients or their surrogates (e.g., health care professionals, regulators and legislators). In this document (ICH Q9), "requirements" refers not only to statutory, legisla3012

tive, or regulatory requirements, but also to such needs and expectations.

Stakeholder: Any individual, group or organization that can affect, be affected by, or perceive itself to be affected by a risk. Decision makers might also be stakeholders. For the purposes of this guideline, the primary stakeholders are the patient, healthcare professional, regulatory authority, and industry.

Risk: The combination of the probability of occurrence of harm and the severity of that harm (ISO/IEC Guide 51).

Risk Assessment: A systematic process of organizing information to support a risk decision to be made within a risk management process. It consists of the identification of hazards and the analysis and evaluation of risks associated with exposure to those hazards.

Risk Communication: The sharing of information about risk and risk management between the decision maker and other stakeholders.

Risk Control: Actions implementing risk management decisions (ISO Guide 73).

Risk Acceptance: The decision to accept risk (ISO Guide

Risk Reduction: Actions taken to lessen the probability of occurrence of harm and the severity of that harm.

Risk Identification: The systematic use of information to identify potential sources of harm (hazards) referring to the risk question or problem description.

Risk Evaluation: The comparison of the estimated risk to given risk criteria using a quantitative or qualitative scale to determine the significance of the risk.

Risk Analysis: The estimation of the risk associated with the identified hazards.

Risk Management: The systematic application of quality management policies, procedures, and practices to the tasks of assessing, controlling, communicating and reviewing risk.

Risk Review: Review or monitoring of output/results of the risk management process considering (if appropriate) new knowledge and experience about the risk.

[ICH Q10 Guideline]

Innovation: The introduction of new technologies or methodologies.

Pharmaceutical Quality System (PQS): Management system to direct and control a pharmaceutical company with regard to quality (ICH Q10 based upon ISO 9000:2005) .

Outsourced Activities: Activities conducted by a contract acceptor under a written agreement with a contract giver.

State of Control: A condition in which the set of controls consistently provides assurance of continued process performance and product quality.

Performance Indicators: Measurable values used to quantify quality objectives to reflect the performance of an organization, process or system, also known as "performance metrics" in some regions.

Continual Improvement: Recurring activity to increase the ability to fulfil requirements (ISO 9000:2005).

Senior Management: Person(s) who direct and control a company or site at the highest levels with the authority and responsibility to mobilize resources within the company or site (ICH Q10 based in part on ISO 9000:2005).

Capability of a Process: Ability of a process to realize a product that will fulfil the requirements of that product. The concept of process capability can also be defined in statistical terms (ISO 9000:2005).

Product Realization: Achievement of a product with the quality attributes appropriate to meet the needs of patients, health care professionals, and regulatory authorities (including compliance with marketing authorization) and internal customers requirements.

Corrective Action: Action to eliminate the cause of a detected non-conformity or other undesirable situation. NOTE: Corrective action is taken to prevent recurrence whereas preventive action is taken to prevent occurrence (ISO 9000:2005).

Enabler: A tool or process which provides the means to achieve an objective.

Knowledge Management: Systematic approach to acquiring, analyzing, storing, and disseminating information related to products, manufacturing processes and components.

Quality Planning: Part of quality management focused on setting quality objectives and specifying necessary operational processes and related resources to fulfil the quality objectives (ISO 9000:2005).

Quality Policy: Overall intentions and direction of an organization related to quality as formally expressed by senior management (ISO 9000:2005).

Quality Manual: Document specifying the quality management system of an organization (ISO 9000:2005).

Quality Objectives: A means to translate the quality policy and strategies into measurable activities.

Quality Risk Management (QRM): A systematic process for the assessment, control, communication and review of risks to the quality of the drug (medicinal) product across the product lifecycle (ICH Q9, Q10). For details, refer to "Basic Concept of Quality Risk Management" in General Information.

Feedback/Feedforward: The modification or control of a process or system by its results or effects. Feedback/feedforward can be applied technically in process control strategies and conceptually in quality management. Feedback is to

reflect results to a previous process (for example: a control of the supply of materials in a previous process), and feed-forward is to reflect results to a subsequent process (for example: a control of time for drying in a subsequent process).

Change Management: A systematic approach to proposing, evaluating, approving, implementing and reviewing changes.

Preventive Action: Action to eliminate the cause of a potential non-conformity or other undesirable potential situation. NOTE: Preventive action is taken to prevent occurrence whereas corrective action is taken to prevent recurrence (ISO 9000:2005).

[ICH Q11 Guideline]

Chemical Transformation Step: For Chemical Entities, a step involved in the synthesis of the chemical structure of the drug substance from precursor molecular fragments. Typically it involves C—X or C—C bond formation or breaking.

Contaminants: Any adventitiously introduced materials (e.g., chemical, biochemical, or microbial species) not intended to be part of the manufacturing process of the drug substance or drug product (ICH Q6B).

3. References

- ICH: Guideline for Q8(R2), Pharmaceutical Development
- 2) ICH: Guideline for Q9, Quality Risk Management.
- ICH: Guideline for Q10, Pharmaceutical Quality Systems.
- 4) ICH: Guideline for Q11, Development and Manufacture of Drug Substance (Chemical Entities and Biotechnological/Biological Entities).
- 5) ICH: Quality Implementation Working Group, Points to Consider (R2), ICH-Endorsed Guide for ICH Q8/Q9/Q10 Implementation
- 6) ICH: Quality Implementation Working Group on Q8, Q9 and Q10 Questions & Answers (R4)

International Harmonization Implemented in the Japanese Pharmacopoeia Seventeenth Edition

Add the following:

Sep. 2017

Harmonized items	JP 17 (Supplement II)	Remarks
Conductivity	2.51 Conductivity Measurement	
Introduction	(Introduction)	
Apparatus	1. Apparatus	
Cell constant determination	2. Cell constant determination	
Calibration of temperature	3. Calibration of temperature	
Calibration of measurement electronics	4. Calibration of measurement electronics	
Temperature compensation	5. Temperature compensation	Deletion of ISO and examples of individual monograph.
Conductivity measurement of fluids	Conductivity measurement of fluids	

JP's particular description: 5. Temperature compensation; Additional explanation for non-linear temperature compensation data.

Feb. 2002/Oct.2009 (Rev. 1)/May 2016 (Rev. 2)

Harmonized items	JP 17 (Supplement II)	Remarks
Ethylcellulose	Ethylcellulose	
Definition	origin, limits of content	
	Description	Non-harmonized item
Identification IR	Identification	
Acidity/alkalinity	Purity (1) Acidity or alkalinity	
Viscosity	Viscosity	
Acetaldehyde	Purity (4) Acetaldehyde	
Chlorides	Purity (2) Chloride	
Loss on drying	Loss on drying	
Sulfated ash	Residue on ignition	
Assay	Assay	
Labelling	labelling	Not specified for the name and concen-
		tration of antioxidant in labelling
	Containers and storage	Non-harmonized item

JP's particular description: Purity Heavy metals.

July 2016

Harmonized items	JP 17 (Supplement II)	Remarks
Hydroxyethylcellulose	Hydroxyethylcellulose	
Definition	origin, limits of content	
	labeling of viscosity	Non-harmonized item
	Description	Non-harmonized item
Identification A (IR)	Identification (1)	
Identification B	Identification (2)	
	Viscosity	Non-harmonized item
pН	pН	
Chlorides	Purity (1) Chloride	
Nitrates	Purity (2) Nitrate	
Aldehydes	Purity (4) Aldehydes	
Loss on drying	Loss on drying	
Sulfated ash	Residue on ignition	
Assay	Assay	
	Containers and storage	Non-harmonized item

JP's particular description: Purity (3) Heavy metals.

Change the following as follows:

Nov. 2016 (Corr. 1)

Harmonized items	JP 17 (Supplement II)	Remarks
Isomalt	Isomalt Hydrate	
Definition	origin, limits of content	
	Description	Non-harmonized item
Identification	Identification (2)	
Conductivity	Conductivity	
Reducing sugars	Purity (4) Reducing sugars	
Related substances	Purity (3) Related substances	
Nickel	Purity (2) Nickel	
Water	Water	
Assay	Assay	
Labelling	labelling	
-	Containers and storage	Non-harmonized item

JP's particular description: Identification (1) Color reaction; Purity (1) Heavy metals; Purity (3) Related substances Test for required detectability, System repeatability; Assay System repeatability.

Sep. 2017 (Rev. 2)/Oct. 2018 (Rev. 3)

Harmonized items	JP 17 (Supplement II)	Remarks
Microcrystalline Cellulose	Microcrystalline Cellulose	
Definition	origin	
	Description	Non-harmonized item
Identification A (IR)	Identification (2)	
Identification B (wet chemistry)	Identification (1)	
Identification C (degree of polymeri-	Identification (3)	
zation)		
Conductivity	Conductivity	
pН	pH	
Loss on drying	Loss on drying	
Residue on ignition	Residue on ignition	
Bulk density	Bulk density	
Water-soluble substances	Purity (2) Water-soluble substances	
Ether-soluble substances	Purity (3) Diethyl ether-soluble sub-	
	stances	
Microbial limits	Microbial limits	
	Containers and storage	Non-harmonized item

JP's particular description: Labelling (Mean degree of polymerization and loss on drying value); Identification (3) In the range of label; Purity (1) Heavy metals; Loss on drying In the range of label.

Sep. 2017 (Rev. 2)

Harmonized items	JP 17 (Supplement II)	Remarks
Hypromellose	Hypromellose	
Definition	origin, limits of content	
	Description	Non-harmonized item
Labeling	labeling of viscosity	
Identification (1)	Identification (1)	
Identification (2)	Identification (2)	
Identification (3)	Identification (3)	
Identification (4)	Identification (4)	
Identification (5)	Identification (5)	
Viscosity	Viscosity	
Method 1	Method I	
Method 2	Method II	
рH	рН	
Loss on drying	Loss on drying	
Residue on ignition	Residue on ignition	
Assay	Assay	
•	Containers and storage	Non-harmonized item

JP's particular description: Purity Heavy metals.

Sep. 2017 (Rev. 3)

Harmonized items	JP 17 (Supplement II)	Remarks
Methylcellulose	Methylcellulose	
Definition	origin, limits of content	
	Description	Non-harmonized item
Labeling	labeling of viscosity	
Identification (1)	Identification (1)	
Identification (2)	Identification (2)	
Identification (3)	Identification (3)	
Identification (4)	Identification (4)	
Identification (5)	Identification (5)	
Viscosity	Viscosity	
Method 1	Method I	
Method 2	Method II	
pН	pН	
Loss on drying	Loss on drying	
Residue on ignition	Residue on ignition	
Assay	Assay	
	Containers and storage	Non-harmonized item

JP's particular description: Purity Heavy metals.

Sep. 2017 (Rev. 2)

Harmonized items	JP 17 (Supplement II)	Remarks
Anhydrous Dibasic Calcium Phosphate	Anhydrous Dibasic Calcium Phosphate	
Definition	limits of content	
	Description	Non-harmonized item
Identification (1)	Identification (1)	
Identification (2)	Identification (2)	
Acid-insoluble substances	Purity (1) Acid-insoluble substances	
Chloride	Purity (2) Chloride	
Sulfate	Purity (3) Sulfate	
Carbonate	Purity (4) Carbonate	
Barium	Purity (6) Barium	
	Purity (7) Arsenic	Non-harmonized item
Loss on ignition	Loss on ignition	
Assay	Assay	
	Containers and storage	Non-harmonized item

JP's particular description: Purity (5) Heavy metals.

Delete the following:

Saccharin Saccharin Sodium Hydrate Total Protein Assay

INDEX

A	Ajmaline, 381 Tablets, 382 Akebia Stem, 1793	Injection, 410 Amiodarone Hydrochloride, 411 Tablets, 412
Absorptive Cream, 755	Alacepril, 382	Amitriptyline Hydrochloride, 414
Acacia, 1791	Tablets, 383	Tablets, 414
Powdered, 1791	L-Alanine, 384	Amlexanox, 415
Acebutolol Hydrochloride, 359	Albumin Tannate, 386	Tablets, 416
Aceglutamide Aluminum, 359, 2671	Alcohol, 895	Amlodipine Besilate, 417
Acemetacin, 361	Dehydrated, 896	Orally Disintegrating Tablets, 418
Capsules, 361	for Disinfection, 897	Tablets, 420
Tablets, 362	Aldioxa, 386	Ammonia Water, 420
Acetaminophen, 363	Granules, 387	Amobarbital, 421
Acetazolamide, 364	Tablets, 387	Amomum Seed, 1797, 2973
Acetic Acid, 365	Alendronate	Powdered, 1797
Glacial, 365	Sodium Hydrate, 388	Amosulalol Hydrochloride, 421
Acetohexamide, 366	Sodium Injection, 389	Tablets, 422
Acetylcholine Chloride for Injection,	Sodium Tablets, 390	Amoxapine, 424
367	Alimemazine Tartrate, 391	Amoxicillin
Acetylcysteine, 368	Alisma Tuber, 1793	Capsules, 425
Acetylsalicylic Acid, 449	Powdered, 1794, 2747	Hydrate, 424, 2671
Tablets, 450	Allopurinol, 392	Amphotericin B, 426
Achyranthes Root, 1792	Tablets, 392	for Injection, 427
Aciclovir, 369	Alminoprofen, 393	Syrup, 428
for Injection, 371	Tablets, 394	Tablets, 428, 2909
for Syrup, 373	Aloe, 1794	Ampicillin
Granules, 370	Powdered, 1795	Anhydrous, 429
Injection, 371	Alpinia Officinarum Rhizome, 1796	Ethoxycarbonyloxyethyl Hydrochlo-
Ointment, 372	Alprazolam, 395	ride, 468
Ophthalmic Ointment, 372	Alprenolol Hydrochloride, 396	Hydrate, 430, 2671
Syrup, 372	Alprostadil, 396	Sodium, 431
Tablets, 374	Alfadex, 400	Sodium and Sulbactam Sodium for
Aclarubicin Hydrochloride, 375	Injection, 398	Injection, 433
Acrinol	Alum, 403	Sodium for Injection, 432
and Zinc Oxide Oil, 377	Solution, 401	Ampicillinphthalidyl Hydrochloride,
and Zinc Oxide Oil, Compound,	Powder, Salicylated, 1542	1633
377	Aluminum	Ampiroxicam, 434
and Zinc Oxide Ointment, 378	Acetylsalicylate, 450	Capsules, 435
Hydrate, 376	Monostearate, 402	Amyl Nitrite, 436
Actinomycin D, 378	Potassium Sulfate Hydrate, 403	Anemarrhena Rhizome, 1797
Adrenaline, 379	Silicate Hydrate with Silicon Dioxide,	Anesthamine, 901
Injection, 379	1796	Anesthetic Ether, 898
Solution, 380	Silicate, Natural, 404	Angelica Dahurica Root, 1798
Adsorbed	Silicate, Synthetic, 405	Anhydrous
Diphtheria-Purified Pertussis-Teta-	Sucrose Sulfate Ester, 1609	Aminobenzylpenicillin, 429
nus Combined Vaccine, 816	Amantadine Hydrochloride, 405	Ampicillin, 429
Diphtheria-Tetanus Combined Tox-	Ambenonium Chloride, 406	Caffeine, 543
oid, 816	Amidotrizoic Acid, 407	Citric Acid, 716
Diphtheria Toxoid for Adult Use,	Amikacin Sulfate, 408	Dibasic Calcium Phosphate, 557,
816	for Injection, 409	2911
Habu-venom Toxoid, 1004, 2934	Injection, 409	Ethanol, 896, 2693
Hepatitis B Vaccine, 1018	Aminoacetic Acid, 995	Lactose, 1132, 2714
Purified Pertussis Vaccine, 1384	Aminobenzylpenicillin, 430	Light, Silicic Acid, 1552, 2955
Tetanus Toxoid, 1670	Anhydrous, 429	Sodium Sulfate, 1994
Afloqualone, 380	Sodium, 431	Antipyrine, 437
Agar, 1792	Aminophylline	Apricot Kernel, 1798
Powdered, 1793	Hydrate, 410	Water, 1799

Cacao Butter, 1821 Cadralazine, 541

Tablets, 542

Caffeine

Bacampicillin Hydrochloride, 468

Bacitracin, 469, 2674 Baclofen, 470 Tablets, 470

Aprindine Hydrochloride, 437 Capsules, 438 Aralia Rhizome, 1800 Arbekacin Sulfate, 439 Injection, 440 Areca, 1800 Argatroban Hydrate, 441 L-Arginine, 442 Hydrochloride, 443 Hydrochloride Injection, 443 Aromatic Castor Oil, 1826 Arotinolol Hydrochloride, 444 Arsenic Trioxide, 445 Arsenical Paste, 445 Arsenous Acid, 445	Bakumondoto Extract, 1805, 2747 Bamethan Sulfate, 471 Barbital, 472 Barium Sulfate, 473 Bear Bile, 1807 Bearberry Leaf, 1807 Beclometasone Dipropionate, 473, 2909 Beef Tallow, 1808 Beeswax White, 1808 Yellow, 1809 Bekanamycin Sulfate, 474 Belladonna Extract, 1810	Biotin, 510 Biperiden Hydrochloride, 510 Biphasic Isophane Insulin Human (Genetical Recombination) Injectable Aqueous Suspension, 2706 Bisacodyl, 511 Suppositories, 512 Bismuth Subgallate, 512 Subnitrate, 513 Bisoprolol Fumarate, 514 Tablets, 515 Bitter Cardamon, 1812 Orange Peel, 1812
Artemisia Capillaris Flower, 1800, 2747, 2973 Leaf, 1801 Ascorbic Acid, 445 and Calcium Pantothenate Tablets,	Root, 1809, 2973 Total Alkaloids, 1810 Benidipine Hydrochloride, 475 Tablets, 476 Benincasa Seed, 1811	Tincture, 1813, 2973 Bleomycin Hydrochloride, 516 Sulfate, 518 Bofutsushosan Extract, 1813, 2748,
447 Injection, 446 Powder, 446 Asiasarum Root, 1801	Benoxinate Hydrochloride, 1345 Benserazide Hydrochloride, 478 Bentonite, 478 Benzalkonium Chloride, 479	2973 Boiogito Extract, 1817, 2750, 2974 Boric Acid, 520 Bromazepam, 520 Bromfenac Sodium
L-Aspartic Acid, 449 Aspirin, 449 Aluminum, 450 Tablets, 450	Solution, 480 Solution 50, Concentrated, 480 Benzbromarone, 481 Benzethonium Chloride, 481 Solution, 482	Hydrate, 2909 Ophthalmic Solution, 2910 Bromhexine Hydrochloride, 521 Bromocriptine Mesilate, 522
Aspoxicillin Hydrate, 451 Astragalus Root, 1803 Atenolol, 452 Atorvastatin Calcium Hydrate, 453	Benzocaine, 901 Benzoic Acid, 482 Benzoin, 1812 Benzyl Alcohol, 483	Bromovalerylurea, 523 Brotizolam, 523 Tablets, 524 Brown Rice, 1820 Bucillamine, 525
Tablets, 454 Atractylodes Lancea Rhizome, 1803 Lancea Rhizome, Powdered, 1804 Rhizome, 1804	Benzoate, 484 Benzylpenicillin Benzathine Hydrate, 485 Potassium, 486, 2674 Potassium for Injection, 487	Tablets, 526 Bucumolol Hydrochloride, 527 Bufetolol Hydrochloride, 528 Buformin Hydrochloride, 529 Delayed-release Tablets, 529
Rhizome, Powdered, 1805 Atropine Sulfate Hydrate, 456 Injection, 456 Auranofin, 457	Bepotastine Besilate, 488 Tablets, 489 Beraprost Sodium, 490 Tablets, 492 Berberine	Tablets, 531 Bumetanide, 531 Bunazosin Hydrochloride, 532 Bupivacaine Hydrochloride Hydrate, 533
Tablets, 458 Azathioprine, 459 Tablets, 460 Azelastine Hydrochloride, 461	Chloride Hydrate, 493 Tannate, 494 Betahistine Mesilate, 495 Tablets, 496	Bupleurum Root, 1820 Bupranolol Hydrochloride, 534 Buprenorphine Hydrochloride, 535 Burdock Fruit, 1821
Granules, 462 Azelnidipine, 463 Tablets, 464 Azithromycin Hydrate, 465 Azosemide, 2672	Betamethasone, 497 Dipropionate, 499, 2909 Sodium Phosphate, 500 Tablets, 498 Valerate, 502	Burnt Alum, 403 Busulfan, 536 Butenafine Hydrochloride, 536 Cream, 537 Solution, 538
Tablets, 2672 Aztreonam, 466 for Injection, 467 B	Valerate and Gentamicin Sulfate Cream, 502 Valerate and Gentamicin Sulfate Ointment, 504 Betamipron, 505	Spray, 538 Butropium Bromide, 539 Butyl Parahydroxybenzoate, 540
	Betaxolol Hydrochloride, 506	

Bethanechol Chloride, 507

Extended-release Tablets, 508

Bezafibrate, 507

Bifonazole, 509

and Sodium Benzoate, 544 Lansoprazole Delayed-release, 1140 for Syrup, 606 Cefazolin Sodium, 607 Anhydrous, 543 Methotrexate, 1225 Hydrate, 543 Nifedipine Extended-release, 1310 for Injection, 608 Calciferol, 882 Nizatidine, 1320 Hydrate, 609 Calcitonin Salmon, 545 Pilsicainide Hydrochloride, 1402 Cefbuperazone Sodium, 610 Calcium Pullulan, 574 Cefcapene Pivoxil Hydrochloride Chloride Hydrate, 550 Ribavirin, 1512 Fine Granules, 613 Chloride Injection, 551 Rifampicin, 1519 Hydrate, 612 Roxatidine Acetate Hydrochloride Folinate, 551 Tablets, 614 Gluconate Hydrate, 552 Extended-release, 1531 Cefdinir, 615 Hydroxide, 553 Sodium Iodide (123I), 1579 Capsules, 617 Sodium Iodide (131I), 1580 Lactate Hydrate, 553 Fine Granules, 617 Sulpiride, 1621 Leucovorin, 551 Cefditoren Pivoxil, 618 Tacrolimus, 1633 Levofolinate Hydrate, 2675 Fine Granules, 619 Oxide, 554 Teprenone, 1661 Tablets, 620 Pantothenate, 554 Cefepime Dihydrochloride Tranexamic Acid, 1709 Paraaminosalicylate Granules, 556 Tranilast, 1712 for Injection, 622 Paraaminosalicvlate Hydrate, 555 Trientine Hydrochloride, 1727 Hydrate, 621 Polystyrene Sulfonate, 559 Ubenimex, 1742 Cefixime Sodium Edetate Hydrate, 561 Captopril, 575 Capsules, 624 Stearate, 562 Carbamazepine, 575 Fine Granules, 2913 Calumba, 1821, 2974 Carbazochrome Sodium Sulfonate Hy-Hydrate, 623, 2676 Powdered, 1822, 2974 drate, 576 Cefmenoxime Hydrochloride, 625 Camellia Oil, 1822 Carbetapentane Citrate, 1378 Cefmetazole Sodium, 627 Camostat Mesilate, 562 Carbetapentene Citrate, 1378 for Injection, 628 d-Camphor, 563 Carbidopa Hydrate, 577 Cefminox Sodium Hydrate, 629 dl-Camphor, 564 Cefodizime Sodium, 629 L-Carbocisteine, 578 Candesartan Cilexetil, 565 Tablets, 579 Cefoperazone Sodium, 631 and Sulbactam Sodium for Injection, and Amlodipine Besylate Tablets, Carbolic Acid, 1390 Carbon Dioxide, 580 632 and Hydrochlorothiazide Tablets, Carboplatin, 580 for Injection, 2677 570 Injection, 582 Cefotaxime Sodium, 634 Tablets, 566 Carboxymethylcellulose, 583 Cefotetan, 635 Capsicum, 1822 Calcium, 583 Cefotiam and Salicylic Acid Spirit, 1824 Sodium, 584 Hexetil Hydrochloride, 637 Powdered, 1823 Hydrochloride, 639 Cardamon, 1825 Tincture, 1824 Carmellose, 583 Hydrochloride for Injection, 640 Capsules, 574 Calcium, 583 Cefozopran Hydrochloride, 640 Capsules Sodium, 584 for Injection, 641 Acemetacin, 361 Carmofur, 586 Cefpiramide Sodium, 642 Amoxicillin, 425 Carnauba Wax, 1825 Cefpirome Sulfate, 643 Ampiroxicam, 435 Carteolol Hydrochloride, 587 Cefpodoxime Proxetil, 644 Aprindine Hydrochloride, 438 Carumonam Sodium, 587 for Syrup, 646 Tablets, 647 Cefaclor, 593 Carvedilol, 589 Tablets, 590 Cefadroxil, 598 Cefroxadine Cefalexin, 601 Cassia Seed, 1825, 2974 for Syrup, 649 Castor Oil, 1826 Cefdinir, 617 Hydrate, 648 Cefixime, 624 Aromatic, 1826 Cefsulodin Sodium, 650 Cinoxacin, 710 Catalpa Fruit, 1826 Ceftazidime Clindamycin Hydrochloride, 722 Cefaclor, 592 for Injection, 653 Clofibrate, 729 Capsules, 593 Hydrate, 652 Clorazepate Dipotassium, 740 Combination Granules, Cefteram Pivoxil, 654 Diltiazem Hydrochloride Extended-Fine Granules, 596 Fine Granules, 655 release, 807 Cefadroxil, 597 Tablets, 656 Doxifluridine, 835 Ceftibuten Hydrate, 657 Capsules, 598 Droxidopa, 842 for Syrup, 599 Ceftizoxime Sodium, 658, 2678 Emedastine Fumarate Extended-Cefalexin, 599 Ceftriaxone Sodium Hydrate, 660 release, 858 Capsules, 601 Cefuroxime Axetil, 662 Ethyl Icosapentate, Combination Granules. Cellacefate, 663, 2679 Flopropione, 939 for Syrup, 603 Cellulose Fluconazole, 941 Cefalotin Sodium, 604 , 2-hydroxypropyl ether, Hypromellose, 574 for Injection, 2912 Acetate Phthalate, 663 Indometacin, 1063 Cefatrizine Propylene Glycolate, 605 Microcrystalline, 664, 2914

Powdered, 667	Cinnamon	Codeine Phosphate
Celmoleukin (Genetical Recombina-	Bark, 1830	Hydrate, 746
tion), 668	Bark, Powdered, 1831	Powder, 1%, 746
Cetanol, 671	Oil, 1831	Powder, 10%, 747
Cetirizine Hydrochloride, 671	Cinoxacin, 709	Tablets, 748
Tablets, 672	Capsules, 710	Cod Liver Oil, 745
Cetotiamine Hydrochloride Hydrate,	Ciprofloxacin, 711	Codonopsis Root, 1835
673	Hydrochloride Hydrate, 712	Coix Seed, 1836
Cetraxate Hydrochloride, 674	Cisplatin, 714	Powdered, 1836
Chenodeoxycholic Acid, 675	Cistanche Herb, 1831	Colchicine, 749
Cherry Bark, 1826	Citicoline, 715	Colestimide, 750
Chloral Hydrate, 676	Citric Acid	Granules, 751
Chloramphenicol, 677, 2917	Anhydrous, 716	Tablets, 752
and Colistin Sodium Methanesul-	Hydrate, 717	Colistin
fonate Ophthalmic Solution,	Citrus Unshiu Peel, 1832	Sodium Methanesulfonate, 752,
2679	Clarithromycin, 717, 2680	2683
Palmitate, 677	for Syrup, 2918	Sulfate, 753
Sodium Succinate, 678, 2679	Tablets, 719	Compound
Chlordiazepoxide, 679	Clebopride Malate, 720	Acrinol and Zinc Oxide Oil, 377
Powder, 680	Clemastine Fumarate, 721	Diastase and Sodium Bicarbonate
Tablets, 681	Clematis Root, 1833	Powder, 783
Chlorhexidine	Clindamycin	Hycodenone Injection, 1347
Gluconate Solution, 682	Hydrochloride, 721	Iodine Glycerin, 1077
Hydrochloride, 683	Hydrochloride Capsules, 722	Methyl Salicylate Spirit, 1238
Chlorinated Lime, 683	Phosphate, 724	Oxycodone and Atropine Injection
Chlormadinone Acetate, 684	Phosphate Injection, 725	1348
Chlorobutanol, 685	Clinofibrate, 725	Oxycodone Injection, 1347
Chlorphenesin Carbamate, 685	Clobetasol Propionate, 726	Phellodendron Powder for
Tablets, 686	Clocapramine Hydrochloride Hydrate,	Cataplasm, 1940
Chlorpheniramine Maleate, 687	727	Rhubarb and Senna Powder, 195
Injection, 688	Clofedanol Hydrochloride, 728	Salicylic Acid Spirit, 1541
Powder, 689	Clofibrate, 728	Scopolia Extract and Diastase Pow
Tablets, 690	Capsules, 729	der, 1976
d-Chloropheniramine Maleate, 691	Clomifene Citrate, 730	Thianthol and Salicylic Acid Solu-
Chlorpromazine Hydrochloride, 692, 2918	Tablets, 731	tion, 1679
Injection, 692	Clomipramine Hydrochloride, 731	Vitamin B Powder, 1765, 2971
-	Tablets, 2681	Concentrated
Tablets, 693 Chlorpropamide, 694	Clonazepam, 732 Fine Granules, 733	Glycerin, 993 Glycerol, 993
Tablets, 695	Tablets, 733	Condurango, 1836
Cholecalciferol, 696	Clonidine Hydrochloride, 734	Fluidextract, 1837
Cholera Vaccine, 696	Cloperastine Hydrochloride, 735,	Coptis Rhizome, 1837
Cholesterol, 697, 2918	2919	Powdered, 1838
Chorionic Gonadotrophin, 997	Clopidogrel	Corn
for Injection, 999	Sulfate, 736	Oil, 1839
Chotosan Extract, 1827, 2751	Sulfate Tablets, 738	Starch, 1601
Chrysanthemum Flower, 1829, 2975	Clorazepate Dipotassium, 739	Cornus Fruit, 1839, 2752
Cibenzoline Succinate, 697	Capsules, 740	Cortisone Acetate, 754
Tablets, 698	Clotiazepam, 741	Corydalis Tuber, 1840
Ciclacillin, 699	Tablets, 2682	Powdered, 1841
Ciclosporin, 700	Clotrimazole, 742	Crataegus Fruit, 1842
A, 700	Clove, 1833	Creams
Cilastatin Sodium, 701	Oil, 1834	Absorptive, 755
Cilazapril	Powdered, 1833	Betamethasone Valerate and Gen-
Hydrate, 702	Cloxacillin Sodium Hydrate, 743,	tamicin Sulfate, 502
Tablets, 703	2683	Butenafine Hydrochloride, 537
Cilnidipine, 704	Cloxazolam, 744	Hydrophilic, 756
Tablets, 705	Cnidium	Ibuprofen Piconol, 1044
Cilostazol, 707	Monnieri Fruit, 1834	Ketoconazole, 1122
Tablets, 708	Rhizome, 1835	Lanoconazole, 2945
Cimetidine, 709	Rhizome, Powdered, 1835	Terbinafine Hydrochloride, 1663
Cimicifuga Rhizome, 1830	Cocaine Hydrochloride, 745	Cresol, 756
Cinchocaine Hydrochloride, 786	Coconut Oil, 1835	Solution, 756

drate, 782

1572

Diagnostic Sodium Citrate Solution,

Index

Solution, Saponated, 757 Diastase, 782 Powder, 1396 Croconazole Hydrochloride, 757 Sodium for Injection, 1398 and Sodium Bicarbonate Powder, Croscarmellose Sodium, 585 Tablets, 1397 Crospovidone, 758 and Sodium Bicarbonate Powder, Diphtheria Crude Glycyrrhiza Extract, 1864, Compound, 783 Antitoxin, Equine, Freeze-dried, 2756, 2976 Diazepam, 783 Crystal Violet, 1237 Tablets, 784 -Purified Pertussis-Tetanus Com-Curcuma Rhizome, 2753 Dibasic bined Vaccine, Absorbed, 816 -Tetanus Combined Toxoid, 816 Cyanamide, 759 Calcium Phosphate, Anhydrous, Cyanocobalamin, 760 557, 2911 -Tetanus Combined Toxoid, Ab-Injection, 761 Calcium Phosphate Hydrate, 558 sorbed, 816 Cyclopentolate Hydrochloride, 761 Sodium Phosphate Hydrate, 1585 Toxoid, 816 Cyclophosphamide Dibekacin Sulfate, 785 Toxoid for Adult Use, Absorbed, Hydrate, 762 Ophthalmic Solution, 785 816 Tablets, 762 Dibucaine Hydrochloride, 786 Dipyridamole, 817 Cycloserine, 763 Dichlorphenamide, 787 Disodium Edetate Hydrate, 1573 Disopyramide, 818 Tablets, 789 Cyperus Rhizome, 1842, 2753 Powdered, 1843, 2753 Diclofenac Sodium, 787 Distigmine Bromide, 818 Tablets, 819 Cyproheptadine Hydrochloride Hy-Suppositories, 2920 Diclofenamide, 787, 2684 drate, 764 Disulfiram, 820 Tablets, 789, 2684 L-Cysteine, 765 Dobutamine Hydrochloride, 820 Hydrochloride Hydrate, 766 Dicloxacillin Sodium Hydrate, 789 Docetaxel L-Cystine, 766 for Injection, 823 Diethylcarbamazine Citrate, 790 Tablets, 790 Cytarabine, 767 Hydrate, 821 Difenidol Hydrochloride, 791 Injection, 822 D Diflorasone Diacetate, 792 Dolichos Seed, 1848 Diflucortolone Valerate, 793 Domperidone, 824 Daiokanzoto Extract, 1843, 2754 Digenea, 1846 Donepezil Hydrochloride, 825 Daisaikoto Extract, 1844, 2754 Digitoxin, 794, 2684 Fine Granules, 826 Tablets, 795, 2684 Danazol, 768 Tablets, 827 Dantrolene Sodium Hydrate, 768 Digoxin, 796, 2685 Dopamine Hydrochloride, 829 Daunorubicin Hydrochloride, 769 Injection, 797 Injection, 829 Deferoxamine Mesilate, 771 Tablets, 798 Doripenem Dehydrated Alcohol, 896 Dihydrocodeine Phosphate, 800 for Injection, 2924 Powder, 1%, 800 Dehydrocholate Sodium Injection, Hydrate, 2921 Powder, 10%, 801 Dorzolamide Hydrochloride, 830 773 Ophthalmic Solution, 831 Dehydrocholic Acid, 772, 2920 Dihydroergotamine Mesilate, 802 Injection, 773 Purified, 772, 2920 Dihydroergotoxine Mesilate, 803 Doxapram Hydrochloride Hydrate, Dilazep Hydrochloride Hydrate, 805 832 Demethylchlortetracycline Hydrochlo-Diltiazem Hydrochloride, 806 Doxazosin Mesilate, 823 ride, 774, 2684 Extended-release Capsules, 807 Tablets, 834 Dilute Doxifluridine, 835 Antiformin, 436 Hydrochloric Acid, 1024 Capsules, 835 Iodine Glycerin, 1078 Iodine Tincture, 1076 Doxorubicin Hydrochloride, 836, Paraformaldehyde Paste, 1367 Diluted Opium Powder, 1924 2685 Phenol with Camphor, 1392 Dimemorfan Phosphate, 808 for Injection, 837 Sodium Hypochlorite Solution, 436 Dimenhydrinate, 809 Doxycycline Hydrochloride Hydrate, 838, 2686 Triozinc Paste, 1734 Tablets, 809 Dermatol, 512 Dimercaprol, 810 Tablets, 840 Deslanoside, 775 Injection, 811 Injection, 776 Dimorpholamine, 811 Aluminum Hydroxide Gel, 401 Dexamethasone, 776 Injection, 812 Aluminum Hydroxide Gel Fine Gran-Dextran Dinoprost, 812 ules, 402 Aluminum Potassium Sulfate, 403 40, 777, 2684 Dionin, 905 Sodium Carbonate, 1568 40 Injection, 778 Dioscorea Rhizome, 1847 70, 779 Powdered, 1847 Sodium Sulfite, 1592 Sulfate Sodium Sulfur 5, 780 Diphenhydramine, 813 Thyroid, 1684 Yeast, 1779 Sulfate Sodium Sulfur 18, 780 and Bromovalerylurea Powder, 814 Dextrin, 781 Hydrochloride, 814 Droperidol, 841 , Phenol and Zinc Oxide Liniment, Dextromethorphan Hydrobromide Hy-Droxidopa, 842

815

Tannate, 815

Diphenylhydantoin, 1396

Capsules, 842

Fine Granules,

Dydrogesterone,

Tablets, 845

 \mathbf{E}

Ebastine, 845 Orally Disintegrating Tablets, 846 Tablets, 847 Ecabet Sodium Granules, 849 Hydrate, 849 Ecarazine Hydrochloride, 1700 Ecothiopate Iodide, 851 Edaravone, 851 Injection, 852, 2687 Edrophonium Chloride, 853 Injection, 854 EDTA Sodium Hydrate, 1573 Elcatonin, 854 Eleutherococcus Senticosus Rhizome, 1848 Emedastine Fumarate, 857 Extended-release Capsules, Emorfazone, 858 Tablets, 859 Enalapril Maleate, 860 Tablets, 861 Enflurane, 863 Enoxacin Hydrate. 863 Entacapone, 2689 Tablets, 2691 Enviomycin Sulfate, 864 Epalrestat, 865, 2692 Tablets, 866 Eperisone Hydrochloride, 867 Ephedra Herb, 1849 Ephedrine Hydrochloride, 868 Injection, 869 Powder, 870 Powder, 10%, 870 Tablets, 870 Epimedium Herb, 1850 Epinephrine, 379 Injection, 379 Solution, 380 Epirizole, 871 Epirubicin Hydrochloride, 872, 2925 Eplerenone, 874 Tablets, 875 **Epoetin** Alfa (Genetical Recombination), Beta (Genetical Recombination), 879 Ergocalciferol, 882 Ergometrine Maleate, 883 Injection, 883 Tablets, 884 Ergotamine Tartrate, 884 Erythromycin, 885, 2692 Delayed-release Tablets,

Ethylsuccinate, 887

Lactobionate, 887

Stearate, 888

Estazolam, 889 Estradiol Benzoate, 889 Injection (Aqueous Suspension), 890 Estriol, 891, 2925 Injection (Aqueous Suspension), Tablets, 892 Etacrynic Acid, Tablets, 893 Ethacridine Lactate, 376 Ethambutol Hydrochloride, Ethanol, 895, 2693 Anhydrous, 896, 2693 for Disinfection, 897 Ethenzamide, 897 Ether, 898 Anesthetic, 898 Guaiacol Glyceryl, 1002 Polyoxyethylene Lauryl Alcohol, 1144 Ethinylestradiol, 899 Tablets, 899 Ethionamide, 900 Ethosuximide, 901 Ethoxybenzamide, 897 Ethyl Aminobenzoate, 901 Cysteine Hydrochloride, 902 L-Cysteine Hydrochloride, 902 Icosapentate, 903 Icosapentate Capsules, 904 Parahydroxybenzoate, 906 Ethylcellulose, 2926 Ethylenediamine, 903 Ethylmorphine Hydrochloride Hydrate, 905 Etidronate Disodium, 907 Tablets, 908 Etilefrine Hydrochloride, 909 Tablets, 909 Etizolam, 910, 2927 Fine Granules, 911 Tablets, 912 Etodolac, 913 Etoposide, 914 Eucalyptus Oil, 1850, 2975 Eucommia Bark, 1851 Euodia Fruit, 1851, 2755 Exsiccated Gypsum, 1868 Extracts Bakumondoto, 1805, 2747 Belladonna, 1810 Bofutsushosan, 1813, 2748, 2973 Boiogito, 1817, 2750, 2974 Chotosan, 1827, 2751 Crude Glycyrrhiza, 1864, 2756, 2976

Daiokanzoto, 1843, 2754

Glycyrrhiza, 1863, 2755, 2975

Goshajinkigan, 1865, 2758

Daisaikoto, 1844, 2754

Goreisan, 2756

Goshuyuto, 2976 Hachimijiogan, 1868, 2759 Hangekobokuto, 1871, 2761 Hangeshashinto, 1873, 2761 Hochuekkito, 1876, 2763, 2977 Juzentaihoto, 1887, 2766, 2978 Kakkonto, 1889, 2768 Kakkontokasenkvushin'i, 1892, 2771, 2978 Kamikihito, 1895, 2772, 2978 Kamishoyosan, 1898, 2773, 2979 Keishibukuryogan, 1901, 2775 Maoto, 1912, 2776 Mukoi-Daikenchuto, 1917, 2777 Nux Vomica, 1921 Orengedokuto, 1926, 2778 Otsujito, 1929, 2778, 2980 Rikkunshito, 1957, 2780 Ryokeijutsukanto, 1961, 2781 Saibokuto, 1963, 2783 Saikokeishito, 1965, 2785 Saireito, 1968, 2787 Scopolia, 1975 Shakuyakukanzoto, 1984, 2790 Shimbuto, 1985, 2791 Shosaikoto, 1988, 2791 Shoseiryuto, 1990, 2793 Tokakujokito, 1999, 2796 Tokishakuyakusan, 2002, 2797, Yokukansan, 2010, 2798, 2983

F

Famotidine, 915 for Injection, 917 Injection, 915 Powder, 918 Tablets, 919 Faropenem Sodium for Syrup, 920 Hydrate, 919 Tablets, 921 Felbinac, 923 Cataplasm, 923 Tape, 924 Felodipine, 2927 Tablets, 2928 Fenbufen, 924 Fennel, 1851 Oil, 1852 Powdered, 1852 Fentanyl Citrate, 925 Ferrous Sulfate Hydrate, 926 Fexofenadine Hydrochloride, 926 Tablets, 927 Filgrastim (Genetical Recombination), 929 Injection. 931 Fine Granules Cefaclor, 596 Cefcapene Pivoxil Hydrochloride, 613

Cefdinir, 617	Formalin, 961	Injection, 2933
Cefditoren Pivoxil, 619	Water, 961	Ointment, 2933
Cefixime, 2913	Formoterol Fumarate Hydrate, 962	Ophthalmic Solution, 983
Cefteram Pivoxil, 655	Forsythia Fruit, 1852	Gentian, 1856
Clonazepam, 733	Fosfomycin	and Sodium Bicarbonate Powder,
Donepezil Hydrochloride, 826	Calcium	1857
Dried Aluminum Hydroxide Gel,	Calcium for Syrup, 963	Powdered, 1856
402	Calcium Hydrate, 962, 2694	Geranium Herb, 1857
Droxidopa, 843	Sodium, 964, 2694	Powdered, 1857
Etizolam, 911	Sodium for Injection, 965	Ginger, 1857
Haloperidol, 1005	Fradiomycin Sulfate, 966	Powdered, 1858
Ifenprodil Tartrate, 1050	Freeze-dried	Processed, 1951
Irsogladine Maleate, 1090	BCG Vaccine (for Percutaneous	Ginseng, 1859
Levofloxacin, 1155	Use), 473	Powdered, 1860
Nifedipine, 1311	Botulism Antitoxin, Equine, 520	Glacial Acetic Acid, 365
Nifedipine Delayed-release, 1309	Diphtheria Antitoxin, Equine, 816	Glehnia Root and Rhizome, 1861
Pravastatin Sodium, 1446	Habu Antivenom, Equine, 1004	Glibenclamide, 984
Precipitated Calcium Carbonate,	Inactivated Tissue Culture Rabies	Gliclazide, 984
548	Vaccine, 1503	Glimepiride, 985
Probucol, 1461	Japanese Encephalitis Vaccine,	Tablets, 986
Quetiapine Fumarate, 1493	1111	Glucose, 988
Risperidone, 1522	Live Attenuated Measles Vaccine,	Hydrate, 2695
Sarpogrelate Hydrochloride, 1544	1196	Injection, 989, 2698
Tranilast, 1713	Live Attenuated Mumps Vaccine,	Purified, 2696
Troxipide, 1735	1278	L-Glutamic Acid, 989
Flavin Adenine Dinucleotide Sodium,	Live Attenuated Rubella Vaccine,	L-Glutamine, 990
932	1535	Glutathione, 991
Flavoxate Hydrochloride, 934	Mamushi Antivenom, Equine, 1191	Glycerin, 992
Flecainide Acetate, 934	Smallpox Vaccine, 1563	and Potash Solution, 994
Tablets, 935	Smallpox Vaccine Prepared in Cell	Concentrated, 993
Flomoxef Sodium, 936	Culture, 1563	Glycerol, 992
for Injection, 938	Tetanus Antitoxin, Equine, 1670,	Concentrated, 993
Flopropione, 939	2967	Glyceryl Monostearate, 995
Capsules, 939	Fritillaria Bulb, 1853	Glycine, 995
Fluconazole, 940	Fructose, 967	Glycyrrhiza, 1862
Capsules, 941	Injection, 968	Extract, 1863, 2755, 2975
Injection, 942	Fudosteine, 968	Extract, Crude, 1864, 2756, 2976
Flucytosine, 943	Tablets, 969	Powdered, 1863
Fludiazepam, 944	Furosemide, 970	Prepared, 1947
Fludrocortisone Acetate, 944	Injection, 971	Gonadorelin Acetate, 996
Fluidextracts	Tablets, 972	Goreisan Extract, 2756
Condurango, 1837	Fursultiamine Hydrochloride, 973	Goshajinkigan Extract, 1865, 2758
Platycodon, 1943, 2980		Goshuyuto Extract, 2976
Uva Ursi, 2008	\mathbf{G}	Gramicidin, 1001, 2698
Flunitrazepam, 945		Granules
Fluocinolone Acetonide, 946	Gabexate Mesilate, 974	Aciclovir, 370
Fluocinonide, 947	β -Galactosidase	Aldioxa, 387
Fluorescein Sodium, 948	(Aspergillus), 975	Azelastine Hydrochloride, 462
Fluorometholone, 949	(Penicillium), 975	Calcium Paraaminosalicylate, 556
Fluorouracil, 950	Gallium (67Ga) Citrate Injection, 976	Cefaclor Combination, 594
Fluoxymesterone, 950, 2694	Gambir, 1853	Cefalexin Combination, 602
Fluphenazine Enanthate, 951	Powdered, 1853	Colestimide, 751
Flurazepam Hydrochloride, 952	Gardenia Fruit, 1854, 2755	Ecabet Sodium, 849
Flurbiprofen, 953	Powdered, 1854	L-Isoleucine, L-Leucine and L-Valine,
Flutamide, 954	Gas Gangrene Antitoxin, Equine, 977	1096
Flutoprazepam, 955	Gastrodia Tuber, 1855, 2975	Minocycline Hydrochloride, 2949
Tablets, 955	Gatifloxacin	Montelukast Sodium, 2724
Fluvoxamine Maleate, 956	Hydrate, 2930	Pas-calcium, 556
Tablets, 958	Ophthalmic Solution, 2931	Polaprezinc, 2953
Foeniculated Ammonia Spirit, 1852	Gefarnate, 977	Ursodeoxycholic Acid, 1749
Folic Acid, 959, 2694	Gelatin, 978	Guaiacol Glyceryl Ether, 1002
Injection, 959	Purified, 980	Guaifenesin, 1002
Tablets 960	Gentamicin Sulfate 982 2695	Guanahenz Acetate 1003

Guanethidine Sulfate, 1003 Gypsum, 1868

H

Hachimijiogan Extract, 1868, 2759 Haloperidol, 1004, 2934 Fine Granules, 1005 Injection, 1006 Tablets, 1006 Halothane, 1007 Haloxazolam, 1008 Hangekobokuto Extract, 1871, 2761 Hangeshashinto Extract, 1873, 2761 Hedysarum Root, 1875 Hemp Fruit, 1876 Heparin Calcium, 1009, 2698 Sodium, 1013, 2699 Sodium Injection, 1017 L-Histidine, 1018 Hydrochloride Hydrate, 1019 Hochuekkito Extract, 1876, 2763, 2977 Homatropine Hydrobromide, 1020 Homochlorcyclizine Hydrochloride, 1020 Honey, 1879 Houttuynia Herb, 1880 Human Chorionic Gonadotrophin, 997 Chorionic Gonadotrophin for Injection, 999 Menopausal Gonadotrophin, 999 Normal Immunoglobulin, 1021 Hycoato Injection, 1348 Hydralazine Hydrochloride, 1021 for Injection, 1022 Powder, 1022 Tablets, 1022 Hydrochloric Acid, 1023 Dilute, 1024 Lemonade, 1024 Hydrochlorothiazide, 1024 Hydrocortisone, 1025, 2934 Acetate, 1026, 2701, 2934 and Diphenhydramine Ointment, 1027, 2934 Butyrate, 1028, 2701 Sodium Phosphate, 1028 Sodium Succinate, 1030 Succinate, 1031 Hydrocotarnine Hydrochloride Hydrate, 1031 Hydrogenated Oil, 1032 Hydrophilic Cream, 756 Petrolatum, 1387 Hydrous Lanolin, 1903 Hydroxocobalamin Acetate, 1033, 2701 Hydroxyethylcellulose, 2934 Hydroxypropylcellulose, 1033, 2702,

2937
Hydroxypropylmethylcellulose, 1038
Hydroxyzine
Hydrochloride, 1036
Pamoate, 1036
Hymecromone, 1037
Hypromellose, 1038, 2704, 2938
Acetate Succinate, 1040
Capsules, 574
Phthalate, 1041

I Ibudilast, 1042 Ibuprofen, 1043 Piconol, 1044 Piconol Cream, 1044 Piconol Ointment, 1045 Ichthammol, 1046 Idarubicin Hydrochloride, 1046 for Injection, 1047 Idoxuridine, 1048 Ophthalmic Solution, 1049 Ifenprodil Tartrate, 1050 Fine Granules, 1050 Tablets, 1051 Imidapril Hydrochloride, 1052 Tablets, 1053 Imipenem and Cilastatin Sodium for Injection, 1056 Hydrate, 1054 Imipramine Hydrochloride, 1057, 2940 Tablets, 1057, 2940 Immature Orange, 1880 Imperata Rhizome, 1880 Indapamide, 1058 Tablets, 1059 Indenolol Hydrochloride, 1060 Indigocarmine, 1061 Injection, 1062 Indium (111In) Chloride Injection, 1062 Indometacin, 1062 Capsules, 1063 Suppositories, 1064 Influenza HA Vaccine, 1065 Injection Acetylcholine Chloride for, 367 Aciclovir, 371 Aciclovir for, 371 Adrenaline, 379 Alendronate Sodium, 389 Alprostadil, 398 Amikacin Sulfate, 409 Amikacin Sulfate for, 409 Aminophylline, 410 Amphotericin B for, 427 Ampicillin Sodium and Sulbactam Sodium for, 433 Ampicillin Sodium for, Arbekacin Sulfate, 440

L-Arginine Hydrochloride, 443 Ascorbic Acid, 446 Atropine Sulfate, 456 Aztreonam for, 467 Benzylpenicillin Potassium for, 487 Calcium Chloride, 551 Carboplatin, 582 Cefalotin Sodium for, 2912 Cefazolin Sodium for, 608 Cefepime Dihydrochloride for, Cefmetazole Sodium for, 628 Cefoperazone Sodium and Sulbactam Sodium for, 632 Cefoperazone Sodium for, 2677 Cefotiam Hydrochloride for, 640 Cefozopran Hydrochloride for, 641 Ceftazidime for, 653 Chlorpheniramine Maleate, 688 Chlorpromazine Hydrochloride, Chorionic Gonadotrophin for, Clindamycin Phosphate, 725 Compound Hycodenone, 1347 Compound Oxycodone, 1347 Compound Oxycodone and Atropine, 1348 Cvanocobalamin, 761 Dehydrocholate Sodium, 773 Dehydrocholic Acid, 773 Deslanoside, 776 Dextran 40, 778 Digoxin, 797 Dimercaprol, 811 Dimorpholamine, 812 Diphenylhydantoin Sodium for, 1398 Docetaxel, 822 Docetaxel for, 823 Dopamine Hydrochloride, 829 Doripenem for, 2924 Doxorubicin Hydrochloride for, Edaravone, 852, 2687 Edrophonium Chloride, 854 Ephedrine Hydrochloride, 869 Epinephrine, 379 Ergometrine Maleate, 883 Estradiol Benzoate (Aqueous Suspension), 890 Estriol (Aqueous Suspension), 892 Famotidine, 915 Famotidine for, 917 Filgrastim (Genetical Recombination), 931 Flomoxef Sodium for, 938 Fluconazole, 942 Folic Acid, 959 Fosfomycin Sodium for, 965 Fructose, 968 Furosemide, 971 Gallium (67Ga) Citrate, 976 Gentamicin Sulfate, Glucose, 989, 2698

Haloperidol, 1006
Heparin Sodium, 1017
Human Chorionic Gonadotrophin
for, 999
Hycoato, 1348
Hydralazine Hydrochloride for,
1022
Idarubicin Hydrochloride for, 1047
Imipenem and Cilastatin Sodium for,
1056
Indigocarmine, 1062
Indium (111In) Chloride, 1062
Insulin Glargine (Genetical Recombi-
nation), 1067
Insulin Human (Genetical Recombi-
nation), 1070
Interferon Alfa (NAMALWA),
1074
Iodinated (131I) Human Serum Albu-
min, 1075 Iohexol, 1082, 2710
1011ex01, 1062, 2/10
Iopamidol, 1084
Isepamicin Sulfate, 1094
Isoniazid, 1100
Isotonic Sodium Chloride, 1571
Levallorphan Tartrate, 1152
Levofloxacin, 1156
Lidocaine, 1162
Lidocaine Hydrochloride, 1162
Lincomycin Hydrochloride, 1164
Magnesium Sulfate, 1190
D-Mannite, 1195 D-Mannitol, 1195
Meglumine Iotalamate, 1206
Meglumine Sodium Amidotrizoate,
1207
Mepivacaine Hydrochloride, 1213
Meropenem for, 1219
Metenolone Enanthate, 1222
Minocycline Hydrochloride for,
1256
Mitomycin C for, 1262
Mitomycin C for, 1262 Morphine and Atropine, 1273
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recom-
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recom-
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322 Norepinephrine, 1322
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322 Norepinephrine, 1322
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322 Norepinephrine, 1322 Norepinephrine Hydrochloride,
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322 Norepinephrine, 1322 Norepinephrine Hydrochloride, 1322 Operidine, 1385
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322 Norepinephrine, 1322 Norepinephrine Hydrochloride, 1322
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322 Norepinephrine, 1322 Norepinephrine Hydrochloride, 1322 Operidine, 1385 Opium Alkaloids and Atropine, 1338
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322 Norepinephrine, 1322 Norepinephrine Hydrochloride, 1322 Operidine, 1385 Opium Alkaloids and Atropine, 1338 Opium Alkaloids and Scopolamine,
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322 Norepinephrine, 1322 Norepinephrine Hydrochloride, 1322 Operidine, 1385 Opium Alkaloids and Atropine, 1338 Opium Alkaloids and Scopolamine, 1339
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322 Norepinephrine, 1322 Norepinephrine Hydrochloride, 1322 Operidine, 1385 Opium Alkaloids and Atropine, 1338 Opium Alkaloids and Scopolamine, 1339 Opium Alkaloids Hydrochlorides,
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322 Norepinephrine, 1322 Norepinephrine Hydrochloride, 1322 Operidine, 1385 Opium Alkaloids and Atropine, 1338 Opium Alkaloids and Scopolamine, 1339 Opium Alkaloids Hydrochlorides, 1337
Mitomycin C for, 1262 Morphine and Atropine, 1273 Morphine Hydrochloride, 1271 Nartograstim for (Genetical Recombination), 1293 Neostigmine Methylsulfate, 1297 Nicardipine Hydrochloride, 1299 Nicotinic Acid, 1307 Noradrenaline, 1322 Noradrenaline Hydrochloride, 1322 Norepinephrine, 1322 Norepinephrine Hydrochloride, 1322 Operidine, 1385 Opium Alkaloids and Atropine, 1338 Opium Alkaloids and Scopolamine, 1339 Opium Alkaloids Hydrochlorides,

Ozagrel Sodium for, 1357 Panipenem and Betamipron for, Papaverine Hydrochloride, 1364 Pazufloxacin Mesilate, 2728 Peplomycin Sulfate for, 1380 Pethidine Hydrochloride, 1385 Phenolsulfonphthalein, 1393 Phenytoin Sodium for, 1398 Piperacillin Sodium for, 1418 Prednisolone Sodium Succinate for, 1456 Procainamide Hydrochloride, 1463 Insulin Procaine Hydrochloride, 1465 Progesterone, 1470 Protamine Sulfate, 1484 Purified Sodium Hyaluronate, 1576 1065 Pyridoxine Hydrochloride, 1490 Reserpine, 1508 Riboflavin Phosphate, 1516 Riboflavin Sodium Phosphate, 1068 1516 Ritodrine Hydrochloride, 2954 Roxatidine Acetate Hydrochloride for, 1533 Sivelestat Sodium for, 1562 Sodium Bicarbonate, 1566 Sodium Chloride, 0.9%, 1571 Sodium Chloride, 10%, 1570 Sodium Chromate (51Cr), 1571 Sodium Citrate for Transfusion, 1572 1079 Sodium Iodohippurate (131I), 1580 Sodium Iotalamate, 1580 Sodium Pertechnetate (99mTc), 1585 Sodium Thiosulfate, 1593 Spectinomycin Hydrochloride for, 1598 Sterile Water for, in Containers, Streptomycin Sulfate for, 1609 Sulfobromophthalein Sodium, 1619 Sulpyrine, 1623 Suxamethonium Chloride, 1628, Suxamethonium Chloride for, 1628 Tazobactam and Piperacillin for, 1645 Teceleukin for (Genetical Recombination), 1652 2711 Testosterone Enanthate, 1668 Testosterone Propionate, 1669 Thallium (201Tl) Chloride, 1672 2940 Thiamine Chloride Hydrochloride, 1675 Thiamylal Sodium for, 1678 Thiopental Sodium for, 1681 Tobramycin, 1696 Tranexamic Acid, 1710 Vancomycin Hydrochloride for, Vasopressin, 1759, 2740 1096

Verapamil Hydrochloride, 2970 Vinblastine Sulfate for, 1763 Vitamin B₁ Hydrochloride, 1675 Vitamin B₂ Phosphate Ester, 1516 Vitamin B₆, 1490 Vitamin B₁₂, 761 Vitamin C, 446 Voriconazole for. Water for, 1774 Weak Opium Alkaloids and Scopolamine, 1340 Xylitol, 1778 Aspart (Genetical Recombination), Glargine (Genetical Recombination), Glargine (Genetical Recombination) Injection, 1067 Human (Genetical Recombination), Human (Genetical Recombination) Injection, 1070 Interferon Alfa (NAMALWA), 1071 Injection, 1074 Iodinated (131 I) Human Serum Albumin Injection, 1075 Iodine, 1075 Glycerin, Compound, 1077 Glycerin, Dental, 1078 , Salicylic Acid and Phenol Spirit, Tincture, 1076 Tincture, Dilute, 1076 Iodoform, 1080 Iohexol, 1080 Injection, 1082, 2710 Iopamidol, 1083 Injection, 1084 Iotalamic Acid, 1085 Iotroxic Acid, 1086 Ipecac, 1881 Powdered, 1881 Syrup, 1882 Ipratropium Bromide Hydrate, 1086 Ipriflavone, 1087 Tablets, 1088 Iproveratril Hydrochloride, 1761 Irbesartan, 1089 and Amlodipine Besilate Tablets, Tablets, 2710 Irinotecan Hydrochloride Hydrate, Irsogladine Maleate, 1090 Fine Granules, 1090 Tablets, 1091 Isepamicin Sulfate, 1093 Injection, 1094 Isoflurane, 1094 L-Isoleucine, 1096 , L-Leucine and L-Valine Granules,

Japanese Angelica Root, 1883, 2977 Angelica Root, Powdered, 1883, 2978 Encephalitis Vaccine, 1111 Gentian, 1884 Gentian, Powdered, 1884 Valerian, 1884 Valerian, Powdered, 1885 Zanthoxylum Peel, 1885 Zanthoxylum Peel, Powdered, 1886 Josamycin, 1112 Propionate, 1113 Tablets, 1113 Jujube, 1886, 2765 Seed, 1886, 2765 Juzentaihoto Extract, 1887, 2766, 2978

K

Kainic Acid
and Santonin Powder, 1115
Hydrate, 1114
Kakkonto Extract, 1889, 2768
Kakkontokasenkyushin'i Extract,
1892, 2771, 2978
Kallidinogenase, 1116
Kamikihito Extract, 1895, 2772, 2978
Kamishoyosan Extract, 1898, 2773,
2979
Kanamycin
Monosulfate, 1118
Sulfate, 1119
Kaolin, 1120

Keishibukuryogan Extract, 1901, 2775
Ketamine Hydrochloride, 1120
Ketoconazole, 1121
Cream, 1122
Lotion, 1123
Solution, 1123
Ketoprofen, 1124
Ketotifen Fumarate, 1125
Kitasamycin, 1126
Acetate, 1127
Tartrate, 1128
Koi, 1903

 \mathbf{L} Labetalol Hydrochloride, 1129 Tablets, 1130 Lactic Acid, 1131 L-Lactic Acid, 1131 Lactose, 1133 Hydrate, 1133, 2715 Anhydrous, 1132, 2714 Lactulose, 1134 Lafutidine, 1135 Tablets, 1136 Lanatoside C, 1137, 2715 Tablets, 1138, 2715 Lanoconazole, 2944 Cream, 2945 Cutaneous Solution, 2946 Ointment, 2946 Lanolin Hydrous, 1903 Purified, 1904 Lansoprazole, 1139 Delayed-release Capsules, 1140 Delayed-release Orally Disintegrating Tablets, 1141 Lard, 1905 Latamoxef Sodium, 1143 Lauromacrogol, 1144 Lemonades Hydrochloric Acid, 1024 Lenampicillin Hydrochloride, 1144 Lenograstim (Genetical Recombination), 1146 Leonurus Herb, 1905 L-Leucine, 1149 Leucomycin, 1126 Acetate, 1127 Tartrate, 1128 Leuprorelin Acetate, 1150 Levallorphan Tartrate, 1152 Injection, 1152

Levodopa, 1153

Hydrate, 1154

Injection, 1156

Tablets, 1157

Fine Granules, 1155

Ophthalmic Solution, 1156

Levomepromazine Maleate, 1159

Levofloxacin

Levothyroxine Sodium Hydrate, 1159 Tablets, 1160 Lidocaine, 1161 Hydrochloride Injection, 1162 Injection, 1162 Anhydrous Silicic Acid, 1552, 2955 Liquid Paraffin, 1365 Lilium Bulb, 1906 Limaprost Alfadex, 1162 Lincomycin Hydrochloride Hydrate, 1164 Injection, 1164 Lindera Root, 1906 Liniments Diphenhydramine, Phenol and Zinc Oxide, 815 Phenol and Zinc Oxide, 1392 Liothyronine Sodium, 1165 Tablets, 1166 Liquefied Phenol, 1391 Liquid Paraffin, 1365 Lisinopril Hydrate, 1167 Tablets, 1168 Lithium Carbonate, 1169 Lithospermum Root, 1907, 2979 Live Oral Poliomyelitis Vaccine, 1428 Lobenzarit Sodium, 1171

Lobenzarit Sodium, 1171
Longan Aril, 1907
Longgu, 1907
Powdered, 1908
Lonicera Leaf and Stem, 1908
Loquat Leaf, 1909
Lorazepam, 1171
Losartan Potassium, 1172
and Hydrochlorothiazide Tablets, 1174
Tablets, 1173

Lotions

Ketoconazole, 1123
Sulfur and Camphor, 1619
Tacalcitol, 1630
Low Substituted Hydroxypropylcellulose, 1035, 2703
Loxoprofen Sodium
Hydrate, 1177
Tablets, 1178
Lycium
Bark, 1909, 2979
Fruit, 1909
Lysine Hydrochloride, 1180
L-Lysine
Acetate, 1179

M

Hydrochloride, 1180

Lysozyme Hydrochloride, 1181

Macrogol 400, 1182 1500, 1182

Mercurochrome, 1216, 2718

4000, 1183 Solution, 1217, 2718 Hydrate, 1258 6000, 1183 Tablets, 1259 Meropenem 20000, 1184 Mitomycin C, 1261 for Injection, 1219 Hydrate, 1217 for Injection, 1262 Ointment, 1184 Magnesium Mesalazine, 2718 Mizoribine, 1262 Aluminosilicate, 2715 Extended-release Tablets, 2720 Tablets, 1263 Aluminometasilicate, 2716 Mestranol, 1220, 2947 Monobasic Calcium Phosphate Hy-Carbonate, 1185 Metenolone drate, 558 Oxide, 1186 Acetate, 1220 Monosodium Trichloroethyl Phos-Silicate, 1187 Enanthate, 1221 phate, 1725 Stearate, 1188 Enanthate Injection, 1222 Syrup, 1726 Sulfate Hydrate, 1189 Metformin Hydrochloride, 1222 Montelukast Sodium, 1264 Sulfate Injection, 1190 Tablets, 1223 Chewable Tablets, 1267 Sulfate Mixture, 1190 Methamphetamine Hydrochloride, Granules, 2724 Magnolia Tablets, 1269 1223 L-Methionine, 1224 Bark, 1910 Morphine Bark, Powdered, 1910 and Atropine Injection, 1273 Methotrexate, 1225 Flower, 1911 Capsules, 1225 Hydrochloride Hydrate, 1270 Mallotus Bark, 1911 Tablets, 2722 Hydrochloride Injection, 1271 Malt, 1912 Methoxsalen, 1227 Hydrochloride Tablets, 1272 Maltose Hydrate, 1190 Methyl Sulfate Hydrate, 1274 Manidipine Hydrochloride, 1192 Parahydroxybenzoate, 1235 Mosapride Citrate Tablets, 1193 Salicylate, 1238 Hydrate, 1275 Powder, 1276 D-Mannite Injection, 1195 Salicylate Spirit, Compound, 1238 Tablets, 1277 D-Mannitol, 1194, 2718 Methylbenactyzium Bromide, 1227 Injection, 1195 Methylcellulose, 1228, 2723, 2947 Moutan Bark, 1915 Powdered, 1916 Maoto Extract, 1912, 2776 Methyldopa Maprotiline Hydrochloride, 1196 Hydrate, 1229 Mukoi-Daikenchuto Extract, 1917, Meclofenoxate Hydrochloride, 1197 Tablets, 1230 2777 Mecobalamin, 1197 dl-Methylephedrine Hydrochloride, Mulberry Bark, 1918 Tablets, 1198 Mupirocin Calcium 1231 Medazepam, 1200 Powder, 1232 Hydrate, 1279 Medicinal Powder, 10%, 1232 Ointment, 1280 Carbon, 1200 Methylergometrine Maleate, 1233 Soap, 1201 Tablets, 1233 Medroxyprogesterone Acetate, 1202 Methylprednisolone, 1236, 2949 Mefenamic Acid, 1203 Succinate, 1236 Nabumetone, 1281 Methylrosanilinium Chloride, 1237 Mefloquine Hydrochloride, 1203 Tablets, 1282 Mefruside, 1204 Methyltestosterone, 1239 Nadolol, 1283 Tablets, 1240 Meticrane, 1241 Nafamostat Mesilate, 1283 Tablets, 1205 Meglumine, 1206 Naftopidil, 1284 Iotalamate Injection, 1206 Metildigoxin, 1242 Orally Disintegrating Tablets, 1285 Sodium Amidotrizoate Injection, Metoclopramide, 1243 Tablets, 1286 1207 Tablets, 1243 Nalidixic Acid, 1287 Melphalan, 1208 Metoprolol Tartrate, 1244 Naloxone Hydrochloride, 1288 Tablets, 1245 Menatetrenone, 1209 Naphazoline Mentha Metronidazole, 1246 and Chlorpheniramine Solution, Herb, 1914 1290 Tablets, 1246 Oil, 1915 Metyrapone, 1247 Hydrochloride, 1289 Water, 1915 Mexiletine Hydrochloride, 1248 Nitrate, 1289 dl-Menthol, 1210 Miconazole, 1249 Naproxen, 1291 l-Menthol, 1211 Nitrate, 1249 Narcotine, 1326 Mepenzolate Bromide, 1211 Microcrystalline Cellulose, 664, 2914 Hydrochloride, 1327 Mepirizole, 871 Micronomicin Sulfate, 1250 Nartograstim (Genetical Recombina-Midecamycin, 1251 Mepitiostane, 1212 tion), 1291 Acetate, 1252 Mepivacaine Hydrochloride, 1213 for Injection, 1293 Injection, 1213 Miglitol, 1253 Natamycin, 1404 Mequitazine, 1214 Nateglinide, 1294 Migrenin, 1254 Tablets, 1215 Minocycline Hydrochloride, 1255 Tablets, 1295 Merbromin, 1216 for Injection, 1256 Natural Aluminum Silicate, 404 Solution, 1217 Granules, 2949 Nelumbo Seed, 1919 Mercaptopurine Hydrate, 1215 Tablets, 1256 Neomycin Sulfate, 966

Mitiglinide Calcium

Neostigmine Methylsulfate, 1297

Injection, 1297	Ibuprofen Piconol, 1045	Oxycodone Hydrochloride Hydrate,
Nicardipine Hydrochloride, 1298	Lanoconazole, 2946	1346
Injection, 1299	Macrogol, 1184	Oxydol, 1349
Nicergoline, 1300	Mupirocin Calcium, 1280	Oxygen, 1350
Powder, 1301	Polyethylene Glycol, 1184	Oxymetholone, 1351
Tablets, 1302	Simple, 1993	Oxytetracycline Hydrochloride, 1351,
Niceritrol, 1303	Sulfur, Salicylic Acid and Thianthol,	2726
Nicomol, 1304	1620	Oxytocin, 1353
Tablets, 1304	Tacalcitol, 1631	Injection, 1355
Nicorandil, 1305	White, 1775	Oyster Shell, 1932
Nicotinamide, 1306	Zinc Oxide, 1784	Powdered, 1932
Nicotinic Acid, 1307	Olive Oil, 1923	Ozagrel Sodium, 1356
Injection, 1307	Olmesartan Medoxomil, 1329	for Injection, 1357
Nifedipine, 1308	Tablets, 1330	Injection, 1357
Delayed-release Fine Granules,	Olopatadine Hydrochloride, 1332	111gev1011, 100.
1309	Tablets, 1333	P
Extended-release Capsules, 1310	Omeprazole, 1334	•
Fine Granules, 1311	Delayed-release Tablets, 1335	Panax Japonicus Rhizome, 1932
Nilvadipine, 1312	Operidine, 1385	Powdered, 1933
Tablets, 1313	Injection, 1385	Pancreatin, 1358
Nitrazepam, 1314	Ophiopogon Root, 1923	Pancuronium Bromide, 1358
Nitrendipine, 1314	Ophthalmic Solution	
	Bromfenac Sodium, 2910	Panipenem, 1359
Tablets, 1315	•	and Betamipron for Injection, 1361
Nitrogen, 1316	Chloramphenicol and Colistin Sodi-	Pantethine, 1362
Nitroglycerin Tablets, 1317	um Methanesulfonate, 2679	Papaverine Hydrochloride, 1363
Nitrous Oxide, 1318	Dibekacin Sulfate, 785	Injection, 1364
Nizatidine, 1319	Dorzolamide Hydrochloride, 831	Paracetamol, 363
Capsules, 1320	Gatifloxacin, 2931	Paraffin, 1364
Noradrenaline, 1321, 2726	Gentamicin Sulfate, 983	Light Liquid, 1365
Hydrochloride Injection, 1322	Idoxuridine, 1049	Liquid, 1365
Injection, 1322	Levofloxacin, 1156	Paraformaldehyde, 1366
Norepinephrine, 1321	Pemirolast Potassium, 1373	Paste, Dental, 1367
Hydrochloride Injection, 1322	Purified Sodium Hyaluronate, 1577	Parnaparin Sodium, 1367
Injection, 1322	Silver Nitrate, 1557	Paroxetine Hydrochloride
Norethisterone, 1322	Tranilast, 1714	Hydrate, 1369
Norfloxacin, 1323	Zinc Sulfate, 1785	Tablets, 1371
Norgestrel, 1324	Ophthalmic Ointment	Pas-calcium
and Ethinylestradiol Tablets, 1324	Aciclovir, 372	Granules, 556
Nortriptyline Hydrochloride, 1326	Opium	Hydrate, 555
Tablets, 2951	Ipecac Powder, 1925	Paste
Noscapine, 1326	Powder, Diluted, 1924	Arsenical, 445
Hydrochloride Hydrate, 1327	Powdered, 1924	Paraformaldehyde, Dental, 1367
Notopterygium, 1919	Tincture, 1925	Triozinc, Dental, 1734
Nuphar Rhizome, 1919	Opium Alkaloids	Pazufloxacin Mesilate, 2726
Nutmeg, 1920	and Atropine Injection, 1338	Injection, 2728
Nux Vomica, 1920, 2979	and Scopolamine Injection, 1339	Peach Kernel, 1933
Extract, 1921	Hydrochlorides, 1336	Powdered, 1934
Extract Powder, 1922	Hydrochlorides Injection, 1337	Peanut Oil, 1935
Tincture, 1922	Orange	Pemirolast Potassium, 1372
Nystatin, 1328	Oil, 1925	for Syrup, 1374
	Peel Syrup, 1926	Ophthalmic Solution, 1373
0	Peel Tincture, 1926	Tablets, 1375
	Orciprenaline Sulfate, 1341	Penbutolol Sulfate, 1376
Ofloxacin, 1328	Orengedokuto Extract, 1926, 2778	Penicillin G Potassium, 486
Ointments	Oriental Bezoar, 1928, 2778, 2979	Pentazocine, 1376
Aciclovir, 372	Otsujito Extract, 1929, 2779, 2980	Pentobarbital Calcium, 1377, 2728
Aciclovir Ophthalmic, 372	Oxapium Iodide, 1342	Tablets, 2729
Acrinol and Zinc Oxide, 378	Oxaprozin, 1343	Pentoxyverine Citrate, 1378
Betamethasone Valerate and Gen-	Oxazolam, 1343	Peony Root, 1935
tamicin Sulfate, 504	Oxetacaine, 1344	Powdered, 1936
Gentamicin Sulfate, 2933	Oxethazaine, 1344	Peplomycin Sulfate, 1379
Hydrocortisone and Diphenhydra-	Oxprenolol Hydrochloride, 1345	for Injection, 1380

Oxybuprocaine Hydrochloride, 1345

Perilla Herb, 1937

mine, 1027, 2934

Darmhanagina 1201	Dhambata Hydrata 1410	Cadaina Phaenhata 10/ 746
Perphenazine, 1381	Phosphate Hydrate, 1419	Codeine Phosphate, 1%, 746
Maleate, 1383	Phosphate Tablets, 1419	Codeine Phosphate, 10%, 747
Maleate Tablets, 1383	Pirarubicin, 1420	Compound Diastase and Sodium Bi-
Tablets, 1382	Pirenoxine, 1421	carbonate, 783
Pethidine Hydrochloride, 1385	Pirenzepine Hydrochloride Hydrate,	Compound Phellodendron, for
Injection, 1385	1422	Cataplasm, 1940
Petroleum Benzin, 1387	Piroxicam, 1423	Compound Rhubarb and Senna,
Peucedanum Root, 1937	Pitavastatin Calcium	1957
Pharbitis Seed, 1938	Hydrate, 1424	Compound Scopolia Extract and Di-
Phellodendron	Tablets, 1425	astase, 1976
, Albumin Tannate and Bismuth Sub-	Pivmecillinam Hydrochloride, 1427	Compound Vitamin B, 1765, 2971
nitrate Powder, 1940	Tablets, 1428	Diastase and Sodium Bicarbonate,
Bark, 1938	Plantago	783
Bark, Powdered, 1939	Herb, 1942	Dihydrocodeine Phosphate, 1%,
Powder for Cataplasm, Compound,	Seed, 1942	800
1940	Platycodon	Dihydrocodeine Phosphate, 10%,
Phenazone, 437	Fluidextract, 1943, 2980	801
Phenethicillin Potassium, 1388	Root, 1942, 2779, 2980	Diluted Opium, 1924
Phenobarbital, 1389	Root, Powdered, 1943, 2980	Diphenhydramine and
Powder, 1390	Pogostemon Herb, 1944	Bromovalerylurea, 814
Powder, 10%, 1390	Polaprezinc, 2952	Diphenylhydantoin, 1396
Phenol, 1390	Granules, 2953	Ephedrine Hydrochloride, 870
and Zinc Oxide Liniment, 1392	Polyethylene Glycol	Ephedrine Hydrochloride, 10%,
for Disinfection, 1391	400, 1182	870
Liquefied, 1391	1500, 1182	Famotidine, 918
with Camphor, Dental, 1392	4000, 1183	Gentian and Sodium Bicarbonate,
Phenolated Water, 1392	6000, 1183	1857
for Disinfection, 1392	20000, 1184	Hydralazine Hydrochloride, 1022
Phenolsulfonphthalein, 1393	Ointment, 1184	Kainic Acid and Santonin, 1115
Injection, 1393	Polygala Root, 1944, 2981	dl-Methylephedrine Hydrochloride,
L-Phenylalanine, 1394	Powdered, 1944, 2981	1232
Phenylbutazone, 1395	Polygonatum Rhizome, 1945	dl-Methylephedrine Hydrochloride,
Phenylephrine Hydrochloride, 1395	Polygonum Root, 1945	10%, 1232
Phenytoin, 1396	Polymixin B Sulfate, 1429, 2730	Mosapride Citrate, 1276
Powder, 1396	Polyoxyethylene Lauryl Alcohol Ether,	Nicergoline, 1301
Sodium for Injection, 1398	1144	Opium Ipecac, 1925
Tablets, 1397	Polyoxyl 40 Stearate, 1430	Phellodendron, Albumin Tannate
Phytomenadione, 1398	Polyporus Sclerotium, 1946	and Bismuth Subnitrate, 1940
Phytonadione, 1398	Powdered, 1946	Phenobarbital, 1390
Picrasma Wood, 1941	Polysorbate 80, 1430, 2730	Phenobarbital, 10%, 1390
Powdered, 1941	Polyvidone, 1439	Phenytoin, 1396
Pilocarpine Hydrochloride, 1399	Polyvinylpyrrolidone, 1439	Reserpine, 1508
Tablets, 1400	Poria Sclerotium, 1946	Reserpine, 0.1%, 1508
Pilsicainide Hydrochloride	Powdered, 1947	Riboflavin, 1514
Capsules, 1402	Potash Soap, 1432	Scopolia Extract, 1975
Hydrate, 1402	Potassium	Scopolia Extract and Carbon, 1976
Pimaricin, 1404	Bromide, 1432	Scopolia Extract and Ethyl
Pimozide, 1405	Canrenoate, 1433	Aminobenzoate, 1976
Pindolol, 1406	Carbonate, 1433	Scopolia Extract, Papaverine and
Pinellia Tuber, 1941	Chloride, 1434	Ethyl Aminobenzoate, 1977,
Pioglitazone Hydrochloride, 1406	Clavulanate, 1434	2790
and Glimepiride Tablets, 1408,	Guaiacolsulfonate, 1436	Swertia and Sodium Bicarbonate,
2952	Hydroxide, 1436	1998, 2982
and Metformin Hydrochloride	Iodide, 1437	Thiamine Chloride Hydrochloride,
Tablets, 1411	Permanganate, 1438	1675
Tablets, 1407	Sulfate, 1438	Vitamin B ₁ Hydrochloride, 1675
Piperneillin	Potato Starch, 1602 Povidone, 1439, 2731	Vitamin B ₂ , 1514
Piperacillin		Vitamin C, 446
Hydrate, 1414 Sodium, 1416	-Iodine, 1441 Powder	Zinc Oxide Starch, 1785 Powdered
Sodium, 1416 Sodium for Injection, 1418	Ascorbic Acid, 446	Acacia, 1791
Piperazine	Chlordiazepoxide, 680	Acacia, 1791 Agar, 1793
Adipate, 1418	Chlorpheniramine Maleate, 689	Agai, 1793 Alisma Tuber, 1794, 2747
1 mipaic, 1710	omorphomiamme maleate, 009	111101110 1 UUUI, 1/77, 4/7/

Aloe, 1795	Tablets, 1448	Pueraria Root, 1952
Amomum Seed, 1797	Prazepam, 1450	Pullulan, 1487
Atractylodes Lancea Rhizome,	Tablets, 1450	Capsules, 574
1804	Prazosin Hydrochloride, 1451	Purified
Atractylodes Rhizome, 1805	Precipitated Calcium Carbonate, 548	Dehydrocholic Acid, 772, 2920
Calumba, 1822, 2974	Fine Granules, 548	Gelatin, 980
Capsicum, 1823	Tablets, 549	Glucose, 2696
Cellulose, 667	Prednisolone, 1452	Lanolin, 1904
Cinnamon Bark, 1831	Acetate, 1454	Shellac, 1551
Clove, 1833	Sodium Phosphate, 1455	Sodium Hyaluronate, 1575
Cnidium Rhizome, 1835	Sodium Succinate for Injection,	Sodium Hyaluronate Injection,
Coix Seed, 1836	1456	1576
Coptis Rhizome, 1838	Succinate, 1456	Sodium Hyaluronate Ophthalmic So-
Corydalis Tuber, 1841	Tablets, 1453	lution, 1577
Cyperus Rhizome, 1843, 2753	Prepared Glycyrrhiza, 1947	Water, 1773
Dioscorea Rhizome, 1847	Primidone, 1458	Water in Containers, 1774
Fennel, 1852	Probenecid, 1458	Pyrantel Pamoate, 1487
Gambir, 1853	Tablets, 1459	Pyrazinamide, 1488
Gardenia Fruit, 1854	Probucol, 1460	Pyridoval Phosphata Hydrota 2721
Gentian, 1856	Fine Granules, 1461 Tablets, 1462	Pyridoxal Phosphate Hydrate, 2731 Pyridoxine Hydrochloride, 1489
Geranium Herb, 1857 Ginger, 1858	Procainamide Hydrochloride, 1462	Injection, 1490
Ginseng, 1860	Injection, 1463	Pyroxylin, 1491
Glycyrrhiza, 1863	Tablets, 1463	Pyrrolnitrin, 1491
Ipecac, 1881	Procaine Hydrochloride, 1464	1 yrronntrini, 1491
Japanese Angelica Root, 1883,	Injection, 1465	Q
2978	Procarbazine Hydrochloride, 1466	¥
Japanese Gentian, 1884	Procaterol Hydrochloride Hydrate,	Quercus Bark, 1952
Japanese Valerian, 1885	1466	Quetiapine Fumarate, 1492
Japanese Zanthoxylum Peel, 1886	Processed	Fine Granules, 1493
Longgu, 1908	Aconite Root, 1948	Tablets, 1494
Magnolia Bark, 1910	Aconite Root, Powdered, 1949	Quick Lime, 554
Moutan Bark, 1916	Ginger, 1951	Quinapril Hydrochloride, 1496
Opium, 1924	Prochlorperazine Maleate, 1467	Tablets, 1497
Oyster Shell, 1932	Tablets, 1468	Quinidine Sulfate Hydrate, 1498
Panax Japonicus Rhizome, 1933	Progesterone, 1469	Quinine
Peach Kernel, 1934	Injection, 1470	Ethyl Carbonate, 1499
Peony Root, 1936	Proglumide, 1470	Hydrochloride Hydrate, 1500
Phellodendron Bark, 1939	L-Proline, 1471	Sulfate Hydrate, 1501
Picrasma Wood, 1941	Promethazine Hydrochloride, 1472	n
Platycodon Root, 1943, 2980	Propafenone Hydrochloride, 1473	R
Polygala Root, 1944, 2981	Tablets, 1474	D 1
Polyporus Sclerotium, 1946	Propantheline Bromide, 1475	Rabeprazole Sodium, 1502
Poria Sclerotium, 1947 Processed Aconite Root, 1949	Propiverine Hydrochloride, 1476 Tablets, 1477	Ranitidine Hydrochloride, 1503 Rape Seed Oil, 1953, 2780
Rhubarb, 1956	Propranolol Hydrochloride, 1478	Rebamipide, 1504
Rose Fruit, 1960	Tablets, 1479	Tablets, 1506
Scutellaria Root, 1979	Propylene Glycol, 1480	Red Ginseng, 1953
Senega, 1981	Propyl Parahydroxybenzoate, 1481	Rehmannia Root, 1954
Senna Leaf, 1982	Propylthiouracil, 1482	Reserpine, 1507
Smilax Rhizome, 1993	Tablets, 1482	Injection, 1508
Sophora Root, 1995	Propyphenazone, 1103	Powder, 1508
Sweet Hydrangea Leaf, 1996, 2796	Prostaglandin	Powder, 0.1%, 1508
Swertia Herb, 1997, 2982	E ₁ , 396	Tablets, 1509
Tragacanth, 2004	E ₁ α-Cyclodextrin Clathrate Com-	Retinol
Turmeric, 2006, 2798	pound, 400	Acetate, 1510
Pranlukast Hydrate, 1442	F_{2a} , 812	Palmitate, 1510
Pranoprofen, 1443	Protamine Sulfate, 1483	Rhubarb, 1955
Prasterone Sodium Sulfate Hydrate,	Injection, 1484	and Senna Powder, Compound,
1444	Prothionamide, 1484	1957
Pravastatin Sodium, 1444	Protirelin, 1485	Powdered, 1956
Fine Granules, 1446 Solution, 1447	Tartrate Hydrate, 1486	Ribavirin, 1511
SOUTION I/I/I/	Prunella Spike, 1951	Capsules, 1512

Index

Butenafine Hydrochloride, 538 Chlorhexidine Gluconate, 682

Compound Thianthol and Salicylic

Riboflavin, 1513	Schisandra Fruit, 1973	Benzoate, 1565
Butyrate, 1514	Schizonepeta Spike, 1973	Bicarbonate, 1566
Phosphate, 1515	Scopolamine	Bicarbonate and Bitter Tincture Mix-
Phosphate Injection, 1516	Butylbromide, 1547	ture, 1994
Powder, 1514	Hydrobromide Hydrate, 1548	Bicarbonate Injection, 1566
Sodium Phosphate, 1515	Scopolia	Bisulfite, 1567
Sodium Phosphate Injection, 1516	Extract, 1974	Borate, 1567
Ribostamycin Sulfate, 1517	Extract and Carbon Powder, 1976	Bromide, 1568
Rice Starch, 1603	Extract and Ethyl Aminobenzoate	Carbonate Hydrate, 1569
Rifampicin, 1518	Powder, 1976	Chloride, 1569
Capsules, 1519	Extract and Tannic Acid Supposito-	Chloride Injection, 0.9%, 1571
Rikkunshito Extract, 1957, 2780	ries, 1978	Chloride Injection, 10%, 1570
Ringer's Solution, 1520	Extract, Papaverine and Ethyl	Chromate (51Cr) Injection, 1571
Risperidone, 1521	Aminobenzoate Powder, 1977,	Citrate Hydrate, 1571
	2790	
Fine Granules, 1522		Citrate Injection for Transfusion,
Oral Solution, 1523	Extract Powder, 1975	1572
Tablets, 1524	Rhizome, 1973, 2982	Cromoglicate, 1573
Ritodrine Hydrochloride, 1525	Scutellaria Root, 1978	Fusidate, 1574, 2959
Injection, 2954	Powdered, 1979	Hyaluronate, Purified, 1575
Tablets, 1526	Senega, 1980	Hyaluronate Injection, Purified,
Rokitamycin, 1528, 2733	Powdered, 1981	1576
Tablets, 1529, 2733	Syrup, 1981	Hyaluronate Ophthalmic Solution,
Rose Fruit, 1959	Senna Leaf, 1981	Purified, 1577
Powdered, 1960	Powdered, 1982	Hydrogen Carbonate, 1566
Rosin, 1960	L-Serine, 1548	Hydrogen Sulfite, 1567
Roxatidine Acetate Hydrochloride,	Serrapeptase, 1549, 2734	Hydroxide, 1578
1530	Sesame, 1983	Iodide, 1579
Extended-release Capsules, 1531	Oil, 1984	Iodide (123I) Capsules, 1579
Extended-release Tablets, 1532	Sevoflurane, 1550	Iodide (131 I) Capsules, 1580
for Injection, 1533	Shakuyakukanzoto Extract, 1984,	Iodide (131I) Solution, 1580
Roxithromycin, 1534	2790	Iodohippurate (131I) Injection, 1580
Tablets, 2733	Shellac	Iotalamate Injection, 1580
Royal Jelly, 1960	Purified, 1551	L-Lactate Ringer's Solution, 1582
Ryokeijutsukanto Extract, 1961, 2781	White, 1552	
Ryokeljutsukanto Extract, 1901, 2701		L-Lactate Solution, 1581
	Shimbuto Extract, 1985, 2791	Lauryl Sulfate, 1584, 2734
S	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587
${f S}$	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection,
Saccharated Pepsin, 1535	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585
Saccharated Pepsin, 1535 Saccharin, 1535, 2955	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1994 Thiosulfate Hydrate, 1592
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1994 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1994 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate, 1593
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542 Salicylic Acid, 1539	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560 Sinomenium Stem and Rhizome,	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1994 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate, 1593 Valproate Extended-release Tablets
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542 Salicylic Acid, 1539 Adhesive Plaster, 1540	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560 Sinomenium Stem and Rhizome, 1993, 2982	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1994 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate, 1593 Valproate Extended-release Tablets A, 2961
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542 Salicylic Acid, 1539 Adhesive Plaster, 1540 Spirit, 1540	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560 Sinomenium Stem and Rhizome, 1993, 2982 Sitagliptin Phosphate	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1592 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate, 1593 Valproate Extended-release Tablets A, 2961 Valproate Extended-release Tablets
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542 Salicylic Acid, 1539 Adhesive Plaster, 1540 Spirit, 1540 Spirit, Compound, 1541	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560 Sinomenium Stem and Rhizome, 1993, 2982 Sitagliptin Phosphate Hydrate, 2955	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1592 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate, 1593 Valproate Extended-release Tablets A, 2961 Valproate Extended-release Tablets B, 2963
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542 Salicylic Acid, 1539 Adhesive Plaster, 1540 Spirit, 1540 Spirit, Compound, 1541 Salvia Miltiorrhiza Root, 1971	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560 Sinomenium Stem and Rhizome, 1993, 2982 Sitagliptin Phosphate Hydrate, 2955 Tablets, 2957	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1994 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate, 1593 Valproate Extended-release Tablets A, 2961 Valproate Extended-release Tablets B, 2963 Valproate Syrup, 1594
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542 Salicylic Acid, 1539 Adhesive Plaster, 1540 Spirit, 1540 Spirit, Compound, 1541 Salvia Miltiorrhiza Root, 1971 Santonin, 1542	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560 Sinomenium Stem and Rhizome, 1993, 2982 Sitagliptin Phosphate Hydrate, 2955 Tablets, 2957 Sivelestat Sodium	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1591 Sulfate Hydrate, 1994 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate, 1593 Valproate Extended-release Tablets A, 2961 Valproate Extended-release Tablets B, 2963 Valproate Syrup, 1594 Valproate Tablets, 1594
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542 Salicylic Acid, 1539 Adhesive Plaster, 1540 Spirit, 1540 Spirit, Compound, 1541 Salvia Miltiorrhiza Root, 1971 Santonin, 1542 Saponated Cresol Solution, 757	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560 Sinomenium Stem and Rhizome, 1993, 2982 Sitagliptin Phosphate Hydrate, 2955 Tablets, 2957 Sivelestat Sodium for Injection, 1562	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1592 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate, 1593 Valproate Extended-release Tablets A, 2961 Valproate Extended-release Tablets B, 2963 Valproate Syrup, 1594 Valproate Tablets, 1594 Solution
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542 Salicylic Acid, 1539 Adhesive Plaster, 1540 Spirit, 1540 Spirit, Compound, 1541 Salvia Miltiorrhiza Root, 1971 Santonin, 1542 Saponated Cresol Solution, 757 Saposhnikovia Root and Rhizome,	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560 Sinomenium Stem and Rhizome, 1993, 2982 Sitagliptin Phosphate Hydrate, 2955 Tablets, 2957 Sivelestat Sodium for Injection, 1562 Hydrate, 1561	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1994 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate Extended-release Tablets A, 2961 Valproate Extended-release Tablets B, 2963 Valproate Syrup, 1594 Valproate Tablets, 1594 Solution Adrenaline, 380
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542 Salicylic Acid, 1539 Adhesive Plaster, 1540 Spirit, 1540 Spirit, Compound, 1541 Salvia Miltiorrhiza Root, 1971 Santonin, 1542 Saponated Cresol Solution, 757 Saposhnikovia Root and Rhizome, 1971, 2789, 2981	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560 Sinomenium Stem and Rhizome, 1993, 2982 Sitagliptin Phosphate Hydrate, 2955 Tablets, 2957 Sivelestat Sodium for Injection, 1562 Hydrate, 1561 Slaked Lime, 553	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1994 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate, 1593 Valproate Extended-release Tablets A, 2961 Valproate Extended-release Tablets B, 2963 Valproate Syrup, 1594 Valproate Tablets, 1594 Solution Adrenaline, 380 Alum, 401
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542 Salicylic Acid, 1539 Adhesive Plaster, 1540 Spirit, 1540 Spirit, Compound, 1541 Salvia Miltiorrhiza Root, 1971 Santonin, 1542 Saponated Cresol Solution, 757 Saposhnikovia Root and Rhizome, 1971, 2789, 2981 Sappan Wood, 1972	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560 Sinomenium Stem and Rhizome, 1993, 2982 Sitagliptin Phosphate Hydrate, 2955 Tablets, 2957 Sivelestat Sodium for Injection, 1562 Hydrate, 1561 Slaked Lime, 553 Smilax Rhizome, 1993	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1592 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate, 1593 Valproate Extended-release Tablets A, 2961 Valproate Extended-release Tablets B, 2963 Valproate Syrup, 1594 Valproate Tablets, 1594 Solution Adrenaline, 380 Alum, 401 Benzalkonium Chloride, 480
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542 Salicylic Acid, 1539 Adhesive Plaster, 1540 Spirit, 1540 Spirit, Compound, 1541 Salvia Miltiorrhiza Root, 1971 Santonin, 1542 Saponated Cresol Solution, 757 Saposhnikovia Root and Rhizome, 1971, 2789, 2981 Sappan Wood, 1972 Sarpogrelate Hydrochloride, 1543	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560 Sinomenium Stem and Rhizome, 1993, 2982 Sitagliptin Phosphate Hydrate, 2955 Tablets, 2957 Sivelestat Sodium for Injection, 1562 Hydrate, 1561 Slaked Lime, 553 Smilax Rhizome, 1993 Powdered, 1993	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1994 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate, 1593 Valproate Extended-release Tablets A, 2961 Valproate Extended-release Tablets B, 2963 Valproate Syrup, 1594 Valproate Tablets, 1594 Solution Adrenaline, 380 Alum, 401 Benzalkonium Chloride, 480 Benzethonium Chloride, 482
Saccharated Pepsin, 1535 Saccharin, 1535, 2955 Sodium, 1536 Sodium Hydrate, 1536, 2734, 2955 Safflower, 1962, 2981 Saffron, 1963 Saibokuto Extract, 1963, 2783 Saikokeishito Extract, 1965, 2785 Saireito Extract, 1968, 2787 Salazosulfapyridine, 1537 Salbutamol Sulfate, 1538 Salicylated Alum Powder, 1542 Salicylic Acid, 1539 Adhesive Plaster, 1540 Spirit, 1540 Spirit, Compound, 1541 Salvia Miltiorrhiza Root, 1971 Santonin, 1542 Saponated Cresol Solution, 757 Saposhnikovia Root and Rhizome, 1971, 2789, 2981 Sappan Wood, 1972	Shimbuto Extract, 1985, 2791 Shosaikoto Extract, 1988, 2791 Shoseiryuto Extract, 1990, 2793 Silodosin, 1553 Tablets, 1555 Silver Nitrate, 1557 Nitrate Ophthalmic Solution, 1557 Protein, 1557 Protein Solution, 1558 Simple Ointment, 1993 Syrup, 1558 Simvastatin, 1559 Tablets, 1560 Sinomenium Stem and Rhizome, 1993, 2982 Sitagliptin Phosphate Hydrate, 2955 Tablets, 2957 Sivelestat Sodium for Injection, 1562 Hydrate, 1561 Slaked Lime, 553 Smilax Rhizome, 1993	Lauryl Sulfate, 1584, 2734 Metabisulfite, 1587 Pertechnetate (99mTc) Injection, 1585 Picosulfate Hydrate, 1585 Polystyrene Sulfonate, 1586 Pyrosulfite, 1587 Risedronate Hydrate, 1588 Risedronate Tablets, 1589 Salicylate, 1590 Starch Glycolate, 1591 Sulfate Hydrate, 1592 Thiosulfate Hydrate, 1592 Thiosulfate Injection, 1593 Valproate, 1593 Valproate Extended-release Tablets A, 2961 Valproate Extended-release Tablets B, 2963 Valproate Syrup, 1594 Valproate Tablets, 1594 Solution Adrenaline, 380 Alum, 401 Benzalkonium Chloride, 480

Saussurea Root, 1972

Acetate Hydrate, 1563 Aurothiomalate, 1564

Solution (continued)	Sucrose, 1611	Pemirolast Potassium for, 1374
Acid, 1679	Sulbactam Sodium, 1612, 2735	Senega, 1981
Cresol, 756	Sulbenicillin Sodium, 1613	Simple, 1558
Dental Sodium Hypochlorite, 436	Sulfadiazine Silver, 1614	Sodium Valproate, 1594
Diagnostic Sodium Citrate, 1572	Sulfafurazole, 1617	Tranilast for, 1715
Epinephrine, 380	Sulfamethizole, 1615	Triclofos Sodium, 1726
Glycerin and Potash, 994	Sulfamethoxazole, 1616	
Isotonic Salt, 1571	Sulfamonomethoxine Hydrate, 1616	T
Isotonic Sodium Chloride, 1571	Sulfasalazine, 1537	
Ketoconazole, 1123	Sulfisomezole, 1616	Tablets
Lanoconazole Cutaneous, 2946	Sulfisoxazole, 1617	Acemetacin, 362
Merbromin, 1217	Sulfobromophthalein Sodium, 1618	Acetylsalicylic Acid, 450
Mercurochrome, 1217, 2718	Injection, 1619	Aciclovir. 374
Naphazoline and Chlorpheniramine,		Ajmaline, 382
· · · · · · · · · · · · · · · · · · ·	Sulfur, 1619	•
1290	and Camphor Lotion, 1619	Alacepril, 383
Pravastatin Sodium, 1447	, Salicylic Acid and Thianthol Oint-	Aldioxa, 387
Ringer's, 1520	ment, 1620	Alendronate Sodium, 390
Risperidone Oral, 1523	Sulindae, 1620	Allopurinol, 392
Saponated Cresol, 757	Sulpiride, 1621	Alminoprofen, 394
Silver Protein, 1558	Capsules, 1621	Amiodarone Hydrochloride, 412
Sodium Iodide (131I), 1580	Tablets, 1622	Amitriptyline Hydrochloride, 414
Sodium L-Lactate, 1581	Sulpyrine	Amlexanox, 416
Sodium L-Lactate Ringer's, 1582	Hydrate, 1623	Amlodipine Besilate, 420
D-Sorbitol, 1597	Injection, 1623	Amlodipine Besilate Orally Disin-
Terbinafine Hydrochloride, 1664	Sultamicillin Tosilate	tegrating, 418
Tolnaftate, 1705	Hydrate, 1624, 2736	Amosulalol Hydrochloride, 422
Sophora Root, 1995	Tablets, 1625	Amphotericin B, 428, 2909
Powdered, 1995	Sultiame, 1626	Ascorbic Acid and Calcium Pan-
Sorbitan Sesquioleate, 1595	Suppositories	tothenate, 447
D-Sorbitol, 1596	Bisacodyl, 512	Aspirin, 450
Solution, 1597	Diclofenac Sodium, 2920	Atorvastatin Calcium, 454
Soybean Oil, 1995	Indometacin, 1064	Auranofin, 458
Spectinomycin Hydrochloride	Scopolia Extract and Tannic Acid,	Azathioprine, 460
for Injection, 1598	1978	Azelnidipine, 464
Hydrate, 1597	Suxamethonium Chloride	Azosemide, 2672
· ·		Baclofen, 470
Spiramycin Acetate, 1599, 2735	for Injection, 1628	
Spirit 1924	Hydrate, 1627	Benidipine Hydrochloride, 476
Capsicum and Salicylic Acid, 1824	Injection, 1628, 2737	Bepotastine Besilate, 489
Compound Methyl Salicylate, 1238	Sweet Hydrangea Leaf, 1996	Beraprost Sodium, 492
Compound Salicylic Acid, 1541	Powdered, 1996, 2796	Betahistine Mesilate, 496
Foeniculated Ammonia, 1852	Swertia	Betamethasone, 498
Iodine, Salicylic Acid and Phenol,	and Sodium Bicarbonate Powder,	Bezafibrate Extended-release, 508
1079	1998, 2982	Bisoprolol Fumarate, 515
Salicylic Acid, 1540	Herb, 1996, 2982	Brotizolam, 524
Spironolactone, 1600	Herb, Powdered, 1997, 2982	Bucillamine, 526
Tablets, 1600	Synthetic Aluminum Silicate, 405	Buformin Hydrochloride, 531
Spray	Syrup	Buformin Hydrochloride Delayed-
Butenafine Hydrochloride, 538	Aciclovir, 372	release, 529
Terbinafine Hydrochloride, 1665	Aciclovir for, 373	Cadralazine, 542
Starch	Amphotericin B, 428	Candesartan Cilexetil, 566
Corn, 1601	Cefadroxil for, 599	Candesartan Cilexetil and Amlodi-
Potato, 1602	Cefalexin for, 603	pine Besylate, 567
Rice, 1603	Cefatrizine Propylene Glycolate for,	Candesartan Cilexetil and Hydro-
Wheat, 1604	606	chlorothiazide, 570
Stearic Acid, 1606	Cefpodoxime Proxetil for, 646	L-Carbocisteine, 579
Stearyl Alcohol, 1607	Cefroxadine for, 649	Carvedilol, 590
Sterile	Clarithromycin for, 2918	Cefcapene Pivoxil Hydrochloride,
Purified Water in Containers, 1774	Faropenem Sodium for, 920	614
Water for Injection in Containers,	Fosfomycin Calcium for, 963	Cefditoren Pivoxil, 620
1775	Ipecac, 1882	Cefpodoxime Proxetil, 647
	Monosodium Trichloroethyl Phos-	Cefteram Pivoxil, 656
Streptomycin Sulfate, 1608		
for Injection, 1609	phate, 1726	Cetirizine Hydrochloride, 672
Sucralfate Hydrate, 1609	Orange Peel, 1926	Chlordiazepoxide, 681

Chlorphenesin Carbamate, 686	2940	Pemirolast Potassium, 1375
Chlorpheniramine Maleate, 690	Indapamide, 1059	Pentobarbital Calcium, 2729
Chlorpromazine Hydrochloride,	Ipriflavone, 1088	Perphenazine, 1382
693	Irbesartan, 2710	Perphenazine Maleate, 1383
Chlorpropamide, 695	Irbesartan and Amlodipine Besilate,	Phenytoin, 1397
Cibenzoline Succinate, 698	2711	Pilocarpine Hydrochloride, 1400
Cilazapril, 703	Irsogladine Maleate, 1091	Pioglitazone Hydrochloride, 1407
Cilnidipine, 705	Isoniazid, 1101	Pioglitazone Hydrochloride and
Cilostazol, 708	Isosorbide Dinitrate, 1105	Glimepiride, 1408, 2952
Clarithromycin, 719	Isosorbide Mononitrate, 1107	Pioglitazone Hydrochloride and Met-
Clomifene Citrate, 731	Isoxsuprine Hydrochloride, 1109	formin Hydrochloride, 1411
Clomipramine Hydrochloride, 2681	Josamycin, 1113	Piperazine Phosphate, 1419
Clonazepam, 733	Labetalol Hydrochloride, 1130	Pitavastatin Calcium, 1425
Clopidogrel Sulfate, 738	Lafutidine, 1136	Pivmecillinam Hydrochloride, 1428
Clotiazepam, 2682	Lanatoside C, 1138, 2715	Pravastatin Sodium, 1448
Codeine Phosphate, 748	Lansoprazole Delayed-release Orally	Prazepam, 1450
Colestimide, 752	Disintegrating, 1141	Precipitated Calcium Carbonate,
Cyclophosphamide, 762	Levofloxacin, 1157	549
Diazepam, 784	Levothyroxine Sodium, 1160	Prednisolone, 1453
	Liothyronine Sodium, 1166	
,		Probenecid, 1459
Diclofenamide, 789, 2784	Lisinopril, 1168	Probucol, 1462
Diethylcarbamazine Citrate, 790	Losartan Potassium, 1173	Procainamide Hydrochloride, 1463
Digitoxin, 795, 2684	Losartan Potassium and Hydrochlo-	Prochlorperazine Maleate, 1468
Digoxin, 798	rothiazide, 1174	Propafenone Hydrochloride, 1474
Dimenhydrinate, 809	Loxoprofen Sodium, 1178	Propiverine Hydrochloride, 1477
Diphenylhydantoin, 1397	Manidipine Hydrochloride, 1193	Propranolol Hydrochloride, 1479
Distigmine Bromide, 819	Mecobalamin, 1198	Propylthiouracil, 1482
Donepezil Hydrochloride, 827	Mefruside, 1205	Quetiapine Fumarate, 1494
Doxazosin Mesilate, 834	Mequitazine, 1215	Quinapril Hydrochloride, 1497
Doxycycline Hydrochloride, 840	Mesalazine Extended-release, 2720	Rebamipide, 1506
Dydrogesterone, 845	Metformin Hydrochloride, 1223	Reserpine, 1509
Ebastine, 847	Methotrexate, 2722	Risperidone, 1524
Ebastine Orally Disintegrating, 846	Methyldopa, 1230	Ritodrine Hydrochloride, 1526
Emorfazone, 859	Methylergometrine Maleate, 1233	Rokitamycin, 1529, 2733
Enalapril Maleate, 861	Methyltestosterone, 1240	Roxatidine Acetate Hydrochloride
Entacapone, 2691	Metoclopramide, 1243	Extended-release, 1532
Epalrestat, 866	Metoprolol Tartrate, 1245	Roxithromycin, 2733
Ephedrine Hydrochloride, 870	Metronidazole, 1246	Sarpogrelate Hydrochloride, 1545
Eplerenone, 875	Minocycline Hydrochloride, 1256	Silodosin, 1555
Ergometrine Maleate, 884	Mitiglinide Calcium, 1259	Simvastatin, 1560
Erythromycin Delayed-release, 886	Mizoribine, 1263	Sitagliptin Phosphate, 2957
Estriol, 892	Montelukast Sodium, 1269	Sodium Risedronate, 1589
Etacrynic Acid, 893	Montelukast Sodium Chewable,	Sodium Valproate, 1594
Ethinylestradiol, 899	1267	Sodium Valproate Extended-release,
Etidronate Disodium, 908	Morphine Hydrochloride, 1272	A, 2961
Etilefrine Hydrochloride, 909	Mosapride Citrate, 1277	Sodium Valproate Extended-release,
Etizolam, 912	Nabumetone, 1282	В, 2963
Famotidine, 919	Naftopidil, 1286	Spironolactone, 1600
Faropenem Sodium, 921	Naftopidil Orally Disintegrating,	Sulpiride, 1622
Felodipine, 2928	1285	Sultamicillin Tosilate, 1625
Fexofenadine Hydrochloride, 927	Nateglinide, 1295	Taltirelin, 1638
Flecainide Acetate, 935	Nicergoline, 1302	Taltirelin Orally Disintegrating,
Flutoprazepam, 955	Nicomol, 1304	1637
Fluvoxamine Maleate, 958	Nilvadipine, 1313	Tamsulosin Hydrochloride Exten-
Folic Acid, 960	Nitrendipine, 1315	ded-release, 1642
Fudosteine, 969	Nitroglycerin, 1317	Telmisartan, 1657
Furosemide, 972	Norgestrel and Ethinylestradiol,	Telmisartan and Hydrochlorothia-
Glimepiride, 986	1324	zide, 2964
Haloperidol, 1006	Nortriptyline Hydrochloride, 2951	Temocapril Hydrochloride, 1659
Hydralazine Hydrochloride, 1022	Olmesartan Medoxomil, 1330	Terbinafine Hydrochloride, 1665
Ifenprodil Tartrate, 1051	Olopatadine Hydrochloride, 1333	Thiamazole, 1673
Imidapril Hydrochloride, 1053	Omeprazole Delayed-release, 1335	Tiapride Hydrochloride, 1685
Iminramine Hydrochloride 1057	Paroxetine Hydrochloride 1371	Tiaramide Hydrochloride 1686

Tablets (continued)	Tablets, 1665	1700
Ticlopidine Hydrochloride, 1688	Terbutaline Sulfate, 1666	Tofisopam, 1701
Tipepidine Hibenzate, 1692, 2967	Testosterone	Tokakujokito Extract, 1999, 2796
Tolbutamide, 1703	Enanthate, 1667, 2967	Tokishakuyakusan Extract, 2002,
Tosufloxacin Tosilate, 1707	Enanthate Injection, 1668	2797, 2983
Tranexamic Acid, 1711	Propionate, 1668	Tolazamide, 1702, 2738
Trichlormethiazide, 1722	Propionate Injection, 1669	Tolbutamide, 1703
Trihexyphenidyl Hydrochloride,	Tetracaine Hydrochloride, 1670	Tablets, 1703
1728	Tetracycline Hydrochloride, 1670,	Tolnaftate, 1704
Trimetazidine Hydrochloride, 1731	2738	Solution, 1705
Troxipide, 1736	Thallium (²⁰¹ Tl) Chloride Injection,	Tolperisone Hydrochloride, 1705
Ursodeoxycholic Acid, 1750	1672	Tosufloxacin Tosilate
Valaciclovir Hydrochloride, 1753	Theophylline, 1672	Hydrate, 1706
Valsartan, 1756	Thiamazole, 1672	Tablets, 1707
Valsartan and Hydrochlorothiazide,	Tablets, 1673	Tragacanth, 2004
2968	Thiamine Chloride Hydrochloride,	Powdered, 2004
Verapamil Hydrochloride, 1761,	1674	Tramadol Hydrochloride, 2739
2741	Injection, 1675	Tranexamic Acid, 1708
Voglibose, 1767	Powder, 1675	Capsules, 1709
Voriconazole, 1770	Thiamine Nitrate, 1676	Injection, 1710
Warfarin Potassium, 1772	Thiamylal Sodium, 1677	Tablets, 1711
Zaltoprofen, 1780	for Injection, 1678	Tranilast, 1711
Zolpidem Tartrate, 1789	Thianthol, 1678	Capsules, 1712
Zonisamide, 2745	and Salicylic Acid Solution, Com-	Fine Granules, 1713
Facalcitol	pound, 1679	for Syrup, 1715
Hydrate, 1629	Thiopental Sodium, 1680	Ophthalmic Solution, 1714
Lotion, 1630	for Injection, 1681	Trapidil, 1716
Ointment, 1631	Thioridazine Hydrochloride, 1681	Trehalose Hydrate, 1717
Гасrolimus	L-Threonine, 1682	Trepibutone, 1718
Capsules, 1633	Thrombin, 1682, 2738	Tretoquinol Hydrochloride, 1733
Hydrate, 1632	Thymol, 1683	Triamcinolone, 1718
Γalampicillin Hydrochloride, 1633	Tiapride Hydrochloride, 1684	Acetonide, 1719, 2967
Γalc, 1634	Tablets, 1685	Triamterene, 1720
Taltirelin	Tiaramide Hydrochloride, 1685	Tribulus Fruit, 2005
Hydrate, 1636	Tablets, 1686	Trichlormethiazide, 1721
Orally Disintegrating Tablets, 1637	Ticlopidine Hydrochloride, 1687	Tablets, 1722
Tablets, 1638	Tablets, 1688	Trichomycin, 1724
Γamoxifen Citrate, 1640	Timepidium Bromide Hydrate, 1689	Trichosanthes Root, 2005, 2798
Γamsulosin Hydrochloride, 1640	Timolol Maleate, 1690	Triclofos Sodium, 1725
Extended-release Tablets, 1642	Tincture	Syrup, 1726
Γannic Acid, 1643	Bitter, 1813, 2973	Trientine Hydrochloride, 1726
Γartaric Acid, 1643	Capsicum, 1824	Capsules, 1727
Γaurine, 1643	Iodine, 1076	Trihexyphenidyl Hydrochloride, 1727
Γazobactam, 1644	Iodine, Dilute, 1076	Tablets, 1728
and Piperacillin for Injection, 1645	Nux Vomica, 1922	Trimebutine Maleate, 1729
Teceleukin	Opium, 1925	Trimetazidine Hydrochloride, 1730
for Injection (Genetical Recombina-	Orange Peel, 1926	Tablets, 1731
tion), 1652	Tinidazole, 1691	Trimetaguir al Hudra chlarida Hudrata
(Genetical Recombination), 1647	Tipepidine Hibenzate, 1691	Trimetoquinol Hydrochloride Hydrate
Гедаfur, 1653 Геісорlanin, 1654, 2737, 2964	Tablets, 1692, 2967 Titanium Oxide, 1693	1733 Tropicamide, 1734
Генсоргани, 1654, 2757, 2964 Гентіsartan, 1656	Tizanidine Hydrochloride, 1694	Tropicalinde, 1734 Troxipide, 1735
and Hydrochlorothiazide Tablets,	Toad Cake, 1998, 2983	Fine Granules, 1735
2964	Tobramycin, 1695, 2738	Tablets, 1736
Tablets, 1657	Injection, 1696	L-Tryptophan, 1737
Γemocapril Hydrochloride, 1658	Tocopherol, 1697	Tulobuterol, 1738
Tablets, 1659	Acetate, 1697	Hydrochloride, 1740
Feprenone, 1660	Calcium Succinate, 1698	Transdermal Tape, 1739
Capsules, 1661	Nicotinate, 1699	Turmeric, 2005, 2798
Ferbinafine Hydrochloride, 1662	dl - α -Tocopherol, 1697	Powdered, 2006, 2798
Cream, 1663	Acetate, 1697	Turpentine Oil, 2007
Solution, 1664	Nicotinate, 1699	L-Tyrosine, 1740
Spray, 1665	Todralazine Hydrochloride Hydrate,	2 2 3 2 0 0 1 1 1 0
· 1 · · · · · · · · · · · · · · · · · ·		

2968

and Hydrochlorothiazide Tablets,

Zonisamide, 2744 Tablets, 2745

\mathbf{U}	Tablets, 1756 Vancomycin Hydrochloride, 1757 for Injection, 1758	in Containers, Purified, 1774 in Containers, Sterile Purified, 1774
Ubenimex, 1741	Vasopressin Injection, 1759, 2740	Purified, 1773
Capsules, 1742	Verapamil Hydrochloride, 1761, 2741	Weak Opium Alkaloids and Scopola-
Ubidecarenone, 1743	Injection, 2970	mine Injection, 1340
Ulinastatin, 1744	Tablets, 1761, 2741	Weil's Disease and Akiyami Combined
	Vinblastine Sulfate, 1762, 2742	Vaccine, 1775
Uncaria Hook, 2007		•
Urapidil, 1746	for Injection, 1763	Wheat Starch, 1604 White
Urea, 1747	Vincristine Sulfate, 1764, 2742	
Urokinase, 1747	Vitamin A	Beeswax, 1808
Ursodeoxycholic Acid, 1748, 2968	Acetate, 1510	Ointment, 1775
Granules, 1749	Oil, 1765	Petrolatum, 1386
Tablets, 1750	Palmitate, 1510	Shellac, 1552
Ursodesoxycholic Acid, 1748	Vitamin B ₁	Soft Sugar, 1611
Uva Ursi Fluidextract, 2008	Hydrochloride, 1674	Whole Human Blood, 1775
${f v}$	Hydrochloride Injection, 1675	Wine, 1776
V	Hydrochloride Powder, 1675	Wood Creosote, 2008
	Nitrate, 1676	V
Vaccine	Vitamin B_2 , 1513	X
BCG, Freeze-dried, (for Percutane-	Phosphate Ester, 1515	
ous Use), 473	Phosphate Ester Injection, 1516	Xylitol, 1778
Cholera, 696	Powder, 1514	Injection, 1778
Diphtheria-Purified Pertussis-Teta-	Vitamin B ₆ , 1489	X 7
nus Combined, Adsorbed, 816	Injection, 1490	Y
Hepatitis B, Adsorbed, 1018	Vitamin B_{12} , 760	
Influenza HA, 1065	Injection, 761	Yellow
Inactivated Tissue Culture Rabies,	Vitamin C, 445	Beeswax, 1809
Freeze-dried, 1503	Injection, 446	Petrolatum, 1387
Japanese Encephalitis, 1111	Powder, 446	Yokukansan Extract, 2010, 2798,
Japanese Encephalitis, Freeze-dried,	Vitamin D_2 , 882	2983
1111	Vitamin D_3 , 696	7
Live Attenuated Measles, Freeze-d-	Vitamin E, 1697	${f Z}$
ried, 1196	Acetate, 1697	
Live Attenuated Mumps, Freeze-d-	Calcium Succinate, 1698	Zaltoprofen, 1780
ried, 1278	Nicotinate, 1699	Tablets, 1780
Live Attenuated Rubella, Freeze-d-	Vitamin K_1 , 1398	Zedoary, 2012, 2753
ried, 1535	Voglibose, 1766	Zidovudine, 1781
Live Oral Poliomyelitis, 1428	Tablets, 1767	Zinc
Purified Pertussis, Adsorbed, 1384	Voriconazole, 1769	Chloride, 1783
Smallpox, Freeze-dried, 1563	for Injection, 2743	Oxide, 1783
Smallpox, Freeze-dried, Prepared in	Tablets, 1770	Oxide Oil, 1784
Cell Culture, 1563	***	Oxide Ointment, 1784
Weil's Disease and Akiyami Com-	${f W}$	Oxide Starch Powder, 1785
bined, 1775		Sulfate Hydrate, 1785
Valaciclovir Hydrochloride, 1751	Warfarin Potassium, 1771	Sulfate Ophthalmic Solution, 1785
Tablets, 1753	Tablets, 1772	Zinostatin Stimalamer, 1786, 2744
L-Valine, 1754	Water, 1773	Zolpidem Tartrate, 1788
Valsartan, 1755	for Injection, 1774	Tablets, 1789

for Injection in Containers, Sterile,

1775

INDEX IN LATIN NAME

Fructus Pulveratus, 1823 Cardamomi Fructus, 1825 Carthami Flos, 1962 Achyranthis Radix, 1792 Caryophylli Flos, 1833 Aconiti Radix Processa, 1948 Flos Pulveratus, 1833 Radix Processa et Pulverata, 1949 Cassiae Semen, 1825 Adeps Lanae Purificatus, 1904 Catalpae Fructus, 1826 Suillus, 1905 Cera Alba, 1808 Agar, 1792 Carnauba, 1825 Flava, 1809 Pulveratum, 1793 Akebiae Caulis, 1793 Chrysanthemi Flos, 1829 Alismatis Tuber, 1793 Cimicifugae Rhizoma, 1830 Tuber Pulveratum, 1794 Cinnamomi Cortex, 1830 Aloe, 1794 Cortex Pulveratus, 1831 Pulverata, 1795 Cistanchis Herba, 1831 Alpiniae Fructus, 1812 Citri Unshiu Pericarpium, 1832 Officinari Rhizoma, 1796 Clematidis Radix, 1833 Cnidii Monnieris Fructus, 1834 Amomi Semen, 1797 Semen Pulveratum, 1797 Rhizoma, 1835 Rhizoma Pulveratum, 1835 Anemarrhenae Rhizoma, 1797 Angelicae Acutilobae Radix, 1883 Codonopsis Radix, 1835 Acutilobae Radix Pulverata, 1883 Coicis Semen, 1836 Dahuricae Radix, 1798 Semen Pulveratum, 1836 Apilac, 1960 Condurango Cortex, 1836 Araliae Cordatae Rhizoma, 1800 Coptidis Rhizoma, 1837 Arctii Fructus, 1821 Rhizoma Pulveratum, 1838 Arecae Semen, 1800 Corni Fructus, 1839 Armeniacae Semen, 1798 Corydalis Tuber, 1840 Artemisiae Capillaris Flos, 1800 Tuber Pulveratum, 1841 Folium, 1801 Crataegi Fructus, 1842 Asiasari Radix, 1801 Creosotum Ligni, 2008 Asparagi Radix, 1802 Crocus, 1963 Astragali Radix, 1803 Curcumae Longae Rhizoma, 2798 Atractylodis Lanceae Rhizoma, Rhizoma Purveratum, 2798 Lanceae Rhizoma Pulveratum, Curcumae Rhizoma, 2005, 2753 1804 Rhizoma Purveratum, 2006 Rhizoma, 1804 Cyperi Rhizoma, 1842 Rhizoma Pulveratum, 1805 Rhizoma Pulveratum, Aurantii Fructus Immaturus, 1880 Pericarpium, 1812

В

Belladonnae Radix, 1809 Benincasae Semen, 1811 Benzoinum, 1812 Bezoar Bovis, 1928 Bufonis Crustum, 1998 Bupleuri Radix, 1820

C

Calumbae Radix, 1821 Radix Pulverata, 1822 Cannabis Fructus, 1876 Capsici Fructus, 1822

D

Digenea, 1846 Dioscoreae Rhizoma, 1847 Rhizoma Pulveratum, 1847 Dolichi Semen, 1848

\mathbf{E}

Eleutherococci senticosi Rhizoma, 1848 Ephedrae Herba, 1849 Epimedii Herba, 1850 Eriobotryae Folium, 1909 Eucommiae Cortex, 1851 Euodiae Fructus, 1851

F

Fel Ursi, 1807 Foeniculi Fructus, 1851 Fructus Pulveratus, 1852 Forsythiae Fructus, 1852 Fossilia Ossis Mastodi, 1907 Ossis Mastodi Pulveratum, 1908 Fritillariae Bulbus, 1853 Fructus Hordei Germinatus, 1912

G

Gambir, 1853 Pulveratum, 1853 Gardeniae Fructus, 1854 Fructus Pulveratus, 1854 Gastrodiae Tuber, 1855 Gentianae Radix, 1856 Radix Pulverata, 1856 Scabrae Radix, 1884 Scabrae Radix Pulverata, 1884 Geranii Herba, 1857 Herba Pulverata, 1857 Ginseng Radix, 1859 Radix Pulverata, 1860 Radix Rubra, 1953 Glehniae Radix cum Rhizoma, 1861 Glycyrrhizae Radix, 1862 Radix Praeparata, 1947 Radix Pulverata, 1863 Gummi Arabicum, 1791 Arabicum Pulveratum, 1791 Gypsum Exsiccatum, 1868 Fibrosum, 1868

H

Hedysari Radix, 1875 Houttuyniae Herba, 1880 Hydrangeae Dulcis Folium, 1996 Dulcis Folium Pulveratum, 1996

I

Imperatae Rhizoma, 1880 Ipecacuanhae Radix, 1881 Radix Pulverata, 1881

K

Kasseki, 1796 Koi, 1903

\mathbf{L}

Leonuri Herba, 1905 Lilii Bulbus, 1906 Linderae Radix, 1906 Lithospermi Radix, 1907 Longan Arillus, 1907 Lonicerae Folium Cum Caulis, 1908 Lycii Cortex, 1909 Fructus, 1909

M

Magnoliae Cortex, 1910
Cortex Pulveratus, 1910
Flos, 1911
Malloti Cortex, 1911
Mel, 1879
Menthae Herba, 1914
Mori Cortex, 1918
Moutan Cortex, 1915
Cortex Pulveratus, 1916
Myristicae Semen, 1920

N

Nelumbis Semen, 1919 Notopterygii Rhizoma, 1919 Nupharis Rhizoma, 1919

Oleum Arachidis, 1935

0

Aurantii, 1925 Cacao, 1821 Camelliae, 1822 Caryophylli, 1834 Cinnamomi, 1831 Cocois, 1835 Eucalypti, 1850 Foeniculi, 1852 Maydis, 1839 Menthae Japonicae, 1915 Olivae, 1923 Rapae, 1953 Ricini, 1826 Sesami, 1984 Soiae, 1995 Terebinthinae, 2007 Ophiopogonis Radix, 1923 Opium Pulveratum, 1924 Oryzae Fructus, 1820 Ostreae Testa, 1932 Testa Pulverata, 1932

P

Paeoniae Radix, 1935 Radix Pulverata, 1936 Panacis Japonici Rhizoma, 1932 Japonici Rhizoma Pulveratum, 1933 Perillae Herba, 1937 Persicae Semen, 1933 Semen Pulveratum, 1934 Peucedani Radix, 1937 Pharbitidis Semen, 1938 Phellodendri Cortex, 1938 Cortex Pulveratus, 1939 Picrasmae Lignum, 1941 Lignum Pulveratum, 1941 Pinelliae Tuber, 1941 Plantaginis Herba, 1942 Semen, 1942 Platycodi Radix, 1942 Radix Pulverata, 1943 Pogostemoni Herba, 1944 Polygalae Radix, 1944 Radix Pulverata, 1944 Polygonati Rhizoma, 1945 Polygoni Multiflori Radix, 1945 Polyporus, 1946 Pulveratus, 1946 Poria, 1946 Pulveratum, 1947 Prunellae Spica, 1951 Pruni Cortex, 1826 Puerariae Radix, 1952

Q

Ouercus Cortex, 1952

R

Rehmanniae Radix, 1954 Resina Pini, 1960 Rhei Rhizoma, 1955 Rhizoma Pulveratum, 1956 Rosae Fructus, 1959 Fructus Pulveratus, 1960

S

Sal Mirabilis, 1994 Mirabilis Anhydricus, 1994 Salviae Miltiorrhizae Radix, 1971 Saposhnikoviae Radix, 1971 Sappan Lignum, 1972 Saussureae Radix, 1972 Schisandrae Fructus, 1973 Schizonepetae Spica, 1973 Scopoliae Rhizoma, 1973 Scutellariae Radix, 1978 Radix Pulverata, 1979 Senegae Radix, 1980 Radix Pulverata, 1981 Sennae Folium, 1981 Folium Pulveratum, Sesami Semen, 1983 Sevum Bovinum, 1808 Sinomeni Caulis et Rhizoma, 1993 Smilacis Rhizoma, 1993 Rhizoma Pulveratum, 1993 Sophorae Radix, 1995 Radix Pulverata, 1995 Strychni Semen, 1920 Swertiae Herba, 1996 Herba Pulverata, 1997

T

Tinctura Amara, 1813 Tragacantha, 2004 Pulverata, 2004 Tribuli Fructus, 2005 Trichosanthis Radix, 2005

U

Uncariae Uncis Cum Ramulus, 2007 Uvae Ursi Folium, 1807

V

Valerianae Fauriei Radix, 1884 Fauriei Radix Pulverata, 1885

\mathbf{Z}

Zanthoxyli Piperiti Pericarpium, 1885
Piperiti Pericarpium Pulveratum,
1886
Zedoariae Rhizoma, 2012
Zingiberis Rhizoma, 1857
Rhizoma Processum, 1951
Rhizoma Pulveratum, 1858
Ziziphi Fructus, 2765
Semen, 2765
Zizyphi Fructus, 1886
Semen, 1886

INDEX IN JAPANESE

ァ

亜鉛華デンプン 1785 亜鉛華軟膏 1784 アカメガシワ 1911 アクチノマイシン D 378 アクラルビシン塩酸塩 375 アクリノール・亜鉛華軟膏 378 アクリノール水和物 376 アクリノール・チンク油 377 アザチオプリン 459 アザチオプリン錠 460 **亜酸化窒素 1318** アシクロビル 369 アシクロビル顆粒 370 アシクロビル眼軟膏 372 アシクロビル錠 374 アシクロビルシロップ 372 アシクロビル注射液 371 アシクロビル軟膏 372 アジスロマイシン水和物 465 アジマリン 381 アジマリン錠 382 亜硝酸アミル 436 アスコルビン酸 445 アスコルビン酸・パントテン酸カルシ ウム錠 447 アスコルビン酸散 446 アスコルビン酸注射液 446 アズトレオナム 466 L-アスパラギン酸 449 アスピリン 449 アスピリンアルミニウム 450 アスピリン錠 450 アスポキシシリン水和物 451 アセグルタミドアルミニウム 359, 2671 アセタゾラミド 364 アセチルシステイン 368 アセトアミノフェン 363 アセトヘキサミド 366 アセブトロール塩酸塩 359 アセメタシン 361 アセメタシンカプセル 361 アセメタシン錠 362 アゼラスチン塩酸塩 461 アゼラスチン塩酸塩顆粒 462 アゼルニジピン 463 アゼルニジピン錠 464 アセンヤク 1853 アセンヤク末 1853 アゾセミド 2672 アゾセミド錠 2672 アテノロール 452 アトルバスタチンカルシウム錠 454 アトルバスタチンカルシウム水和物 453 アドレナリン 379 アドレナリン液 380 アドレナリン注射液 379 アトロピン硫酸塩水和物 456 アトロピン硫酸塩注射液 456 亜ヒ酸パスタ 445 アプリンジン塩酸塩 437 アプリンジン塩酸塩カプセル 438 アフロクアロン 380 アヘン・トコン散 1925 アヘンアルカロイド・アトロピン注射 液 1338 アヘンアルカロイド・スコポラミン注 射液 1339 アヘンアルカロイド塩酸塩 1336 アヘンアルカロイド塩酸塩注射液 1337 アヘン散 1924 アヘンチンキ 1925 アヘン末 1924 アマチャ 1996 アマチャ末 1996, 2796 アマンタジン塩酸塩 405 アミオダロン塩酸塩 411 アミオダロン塩酸塩錠 412 アミカシン硫酸塩 408 アミカシン硫酸塩注射液 409 アミドトリゾ酸 407 アミドトリゾ酸ナトリウムメグルミン 注射液 1207 アミトリプチリン塩酸塩 414 アミトリプチリン塩酸塩錠 414 アミノ安息香酸エチル 901 アミノフィリン水和物 410 アミノフィリン注射液 410 アムホテリシンB 426 アムホテリシンB錠 428,2909 アムホテリシン B シロップ 428 アムロジピンベシル酸塩 417 アムロジピンベシル酸塩口腔内崩壊錠 418 アムロジピンベシル酸塩錠 420 アモキサピン 424 アモキシシリンカプセル 425 アモキシシリン水和物 424, 2671 アモスラロール塩酸塩 421 アモスラロール塩酸塩錠 422

アモバルビタール 421

アラセプリル 382

L-アラニン 384

アラセプリル錠 383

アラビアゴム 1791

アラビアゴム末 1791

アリメマジン酒石酸塩 391

亜硫酸水素ナトリウム 1567 アルガトロバン水和物 441 L-アルギニン 442 L-アルギニン塩酸塩 443 L-アルギニン塩酸塩注射液 443 アルジオキサ 386 アルジオキサ顆粒 387 アルジオキサ錠 387 アルプラゾラム 395 アルプレノロール塩酸塩 396 アルプロスタジル 396 アルプロスタジル アルファデクス アルプロスタジル注射液 398 アルベカシン硫酸塩 439 アルベカシン硫酸塩注射液 440 アルミノプロフェン 393 アルミノプロフェン錠 394 アレンドロン酸ナトリウム錠 390 アレンドロン酸ナトリウム水和物 アレンドロン酸ナトリウム注射液 389 アロエ 1794 アロエ末 1795 アロチノロール塩酸塩 444 アロプリノール 392 アロプリノール錠 392 安息香酸 482 安息香酸ナトリウム 1565 安息香酸ナトリウムカフェイン 544 安息香酸ベンジル 484 アンソッコウ 1812 アンチピリン 437 アンピシリン水和物 430, 2671 アンピシリンナトリウム 431 アンピロキシカム 434 アンピロキシカムカプセル 435 アンベノニウム塩化物 406 アンモニア・ウイキョウ精 1852 アンモニア水 420 アンレキサノクス 415 アンレキサノクス錠 416

1

イオウ 1619 イオウ・カンフルローション 1619 イオウ・サリチル酸・チアントール軟膏 1620 イオタラム酸 1085 イオタラム酸ナトリウム注射液 1580 イオタラム酸メグルミン注射液 1206 イオトロクス酸 1086 イオパミドール 1083 イオパミドール注射液 1084

イオヘキソール 1080 イオヘキソール注射液 1082,2710 イクタモール 1046 イコサペント酸エチル 903 イコサペント酸エチルカプセル 904 イセパマイシン硫酸塩 1093 イセパマイシン硫酸塩注射液 1094 イソクスプリン塩酸塩 1108 イソクスプリン塩酸塩錠 1109 イソソルビド 1103 イソニアジド 1099 イソニアジド錠 1101 イソニアジド注射液 1100 イソフェンインスリン ヒト(遺伝子 組換え)水性懸濁注射液 2705 イソフルラン 1094 *l*-イソプレナリン塩酸塩 1102 イソプロパノール 1102 イソプロピルアンチピリン 1103 イソマル水和物 1098, 2942 L-イソロイシン 1096 イソロイシン・ロイシン・バリン顆粒 1096 イダルビシン塩酸塩 1046 一硝酸イソソルビド錠 1107 70%一硝酸イソソルビド乳糖末 1105, 2714 イドクスウリジン 1048 イドクスウリジン点眼液 1049 イトラコナゾール 1110 イフェンプロジル酒石酸塩 1050 イフェンプロジル酒石酸塩細粒 1050 イフェンプロジル酒石酸塩錠 1051 イブジラスト 1042 イブプロフェン 1043 イブプロフェンピコノール 1044 イブプロフェンピコノールクリーム イブプロフェンピコノール軟膏 1045 イプラトロピウム臭化物水和物 1086 イプリフラボン 1087 イプリフラボン錠 1088 イミダプリル塩酸塩 1052 イミダプリル塩酸塩錠 1053 イミプラミン塩酸塩 1057, 2940 イミプラミン塩酸塩錠 1057, 2940 イミペネム水和物 1054 イリノテカン塩酸塩水和物 2940 イルソグラジンマレイン酸塩 1090 イルソグラジンマレイン酸塩細粒 1090 イルソグラジンマレイン酸塩錠 1091 イルベサルタン 1089 イルベサルタン・アムロジピンベシル 酸塩錠 2711 イルベサルタン錠 2710 イレイセン 1833 インジゴカルミン 1061 インジゴカルミン注射液 1062 インスリン アスパルト(遺伝子組換 え) 2708 インスリン グラルギン(遺伝子組換

え) 1065

インスリン グラルギン(遺伝子組換 え)注射液 1067 インスリン ヒト(遺伝子組換え) 1068 インスリン ヒト(遺伝子組換え)注射 液 1070 インダパミド 1058 インダパミド錠 1059 インターフェロン アルファ (NAMALWA) 1071 インターフェロン アルファ (NAMALWA)注射液 1074 インチンコウ 1800, 2747, 2973 インデノロール塩酸塩 1060 インドメタシン 1062 インドメタシンカプセル 1063 インドメタシン坐剤 1064 インフルエンザ HA ワクチン 1065 インヨウカク 1850

ウ

ウイキョウ 1851 ウイキョウ末 1852 ウイキョウ油 1852 ウコン 2005, 2798 ウコン末 2006, 2798 ウベニメクス 1741 ウベニメクスカプセル 1742 ウヤク 1906 ウラピジル 1746 ウリナスタチン 1744 ウルソデオキシコール酸 1748, 2968 ウルソデオキシコール酸顆粒 1749 ウルソデオキシコール酸錠 1750 ウロキナーゼ 1747 ウワウルシ 1807 ウワウルシ流エキス 2008

エ

エイジツ 1959 エイジツ末 1960 エカベトナトリウム顆粒 849 エカベトナトリウム水和物 849 液状フェノール 1391 エコチオパートヨウ化物 851 エスタゾラム 889 エストラジオール安息香酸エステル エストラジオール安息香酸エステル水 性懸濁注射液 890 エストリオール 891, 2925 エストリオール錠 892 エストリオール水性懸濁注射液 892 エタクリン酸 893 エタクリン酸錠 893 エタノール 895, 2693 エダラボン 851 エダラボン注射液 852, 2687 エタンブトール塩酸塩 894 エチオナミド 900 エチゾラム 910, 2927

エチゾラム細粒 911 エチゾラム錠 912 エチドロン酸二ナトリウム 907 エチドロン酸二ナトリウム錠 908 エチニルエストラジオール 899 エチニルエストラジオール錠 899 L-エチルシステイン塩酸塩 902 エチルセルロース 2926 エチルモルヒネ塩酸塩水和物 905 エチレフリン塩酸塩 909 エチレフリン塩酸塩錠 909 エチレンジアミン 903 エデト酸カルシウムナトリウム水和物 561 エデト酸ナトリウム水和物 1573 エーテル 898 エテンザミド 897 エトスクシミド 901 エトドラク 913 エトポシド 914 エドロホニウム塩化物 853 エドロホニウム塩化物注射液 854 エナラプリルマレイン酸塩 860 エナラプリルマレイン酸塩錠 861 エノキサシン水和物 863 エバスチン 845 エバスチン口腔内崩壊錠 846 エバスチン錠 847 エパルレスタット 865, 2692 エパルレスタット錠 866 エピリゾール 871 エピルビシン塩酸塩 872, 2925 エフェドリン塩酸塩 868 エフェドリン塩酸塩散10% 870 エフェドリン塩酸塩錠 870 エフェドリン塩酸塩注射液 869 エプレレノン 874 エプレレノン錠 875 エペリゾン塩酸塩 867 エポエチン アルファ(遺伝子組換え) エポエチン ベータ(遺伝子組換え) 879 エメダスチンフマル酸塩 857 エメダスチンフマル酸塩徐放カプセル 858 エモルファゾン 858 エモルファゾン錠 859 エリスロマイシン 885, 2692 エリスロマイシンエチルコハク酸エス テル 887 エリスロマイシンステアリン酸塩 888 エリスロマイシン腸溶錠 886 エリスロマイシンラクトビオン酸塩 887 エルカトニン 854 エルゴカルシフェロール 882 エルゴタミン酒石酸塩 884

エルゴメトリンマレイン酸塩 883

エルゴメトリンマレイン酸塩錠 884

エルゴメトリンマレイン酸塩注射液

883

塩化亜鉛 1783 塩化インジウム(111In)注射液 1062 塩化カリウム 1434 塩化カルシウム水和物 550 塩化カルシウム注射液 551 塩化タリウム(201Tl)注射液 1672 塩化ナトリウム 1569 10%塩化ナトリウム注射液 1570 エンゴサク 1840 エンゴサク末 1841 塩酸 1023 塩酸リモナーデ 1024 エンタカポン 2689 エンタカポン錠 2691 エンビオマイシン硫酸塩 864 エンフルラン 863

ォ

オウギ 1803 オウゴン 1978 オウゴン末 1979 黄色ワセリン 1387 オウセイ 1945 オウバク 1938 オウバク・タンナルビン・ビスマス散 1940 オウバク末 1939 オウヒ 1826 オウレン 1837 黄連解毒湯エキス 1926, 2778 オウレン末 1838 オキサゾラム 1343 オキサピウムヨウ化物 1342 オキサプロジン 1343 オキシコドン塩酸塩水和物 1346 オキシテトラサイクリン塩酸塩 1351, 2726 オキシトシン 1353 オキシトシン注射液 1355 オキシドール 1349 オキシブプロカイン塩酸塩 1345 オキシメトロン 1351 オキセサゼイン 1344 オクスプレノロール塩酸塩 1345 オザグレルナトリウム 1356 オザグレルナトリウム注射液 1357 乙字湯エキス 1929, 2779, 2980 オフロキサシン 1328 オメプラゾール 1334 オメプラゾール腸溶錠 1335 オーラノフィン 457 オーラノフィン錠 458 オリブ油 1923 オルシプレナリン硫酸塩 1341 オルメサルタン メドキソミル 1329 オルメサルタン メドキソミル錠 1330 オレンジ油 1925 オロパタジン塩酸塩 1332 オロパタジン塩酸塩錠 1333

オンジ 1944, 2981

オンジ末 1944, 2981

力

カイニン酸・サントニン散 1115 カイニン酸水和物 1114 ガイヨウ 1801 カオリン 1120 カカオ脂 1821 加香ヒマシ油 1826 カゴソウ 1951 カシュウ 1945 ガジュツ 2012, 2753 加水ラノリン 1903 ガスえそウマ抗毒素 977 ガチフロキサシン水和物 2930 ガチフロキサシン点眼液 2931 カッコウ 1944 カッコン 1952 葛根湯エキス 1889,2768 葛根湯加川芎辛夷エキス 1892, 2771, 2978 カッセキ 1796 過テクネチウム酸ナトリウム(99mTc) 注射液 1585 果糖 967 果糖注射液 968 カドララジン 541 カドララジン錠 542 カナマイシン一硫酸塩 1118 カナマイシン硫酸塩 1119 カノコソウ 1884 カノコソウ末 1885 カフェイン水和物 543 カプセル 574 カプトプリル 575 ガベキサートメシル酸塩 974 過マンガン酸カリウム 1438 加味帰脾湯エキス 1895, 2772, 2978 加味逍遙散エキス 1898, 2773, 2979 カモスタットメシル酸塩 562 **β**-ガラクトシダーゼ(アスペルギルス) **β**-ガラクトシダーゼ(ペニシリウム) 975 カリジノゲナーゼ 1116 カリ石ケン 1432 カルシトニン サケ 545 カルテオロール塩酸塩 587 カルナウバロウ 1825 カルバゾクロムスルホン酸ナトリウム 水和物 576 カルバマゼピン 575 カルビドパ水和物 577 カルベジロール 589 カルベジロール錠 590 L-カルボシステイン 578 L-カルボシステイン錠 579 カルボプラチン 580 カルボプラチン注射液 582 カルメロース 583 カルメロースカルシウム 583

カルメロースナトリウム 584

カルモナムナトリウム 587

カルモフール 586 カロコン 2005, 2798 カンキョウ 1951 カンゾウ 1862 乾燥亜硫酸ナトリウム 1592 カンゾウエキス 1863, 2755, 2975 乾燥甲状腺 1684 乾燥酵母 1779 乾燥細胞培養痘そうワクチン 1563 乾燥ジフテリアウマ抗毒素 816 乾燥弱毒生おたふくかぜワクチン 1278 乾燥弱毒生風しんワクチン 1535 乾燥弱毒生麻しんワクチン 1196 乾燥水酸化アルミニウムゲル 401 乾燥水酸化アルミニウムゲル細粒 402 カンゾウ粗エキス 1864, 2756, 2976 乾燥組織培養不活化狂犬病ワクチン 1503 乾燥炭酸ナトリウム 1568 乾燥痘そうワクチン 1563 乾燥日本脳炎ワクチン 1111 乾燥破傷風ウマ抗毒素 1670, 2967 乾燥はぶウマ抗毒素 1004 乾燥 BCG ワクチン 473 乾燥ボツリヌスウマ抗毒素 520 カンゾウ末 1863 乾燥まむしウマ抗毒素 1191 乾燥硫酸アルミニウムカリウム 403 カンデサルタン シレキセチル 565 カンデサルタン シレキセチル・アム ロジピンベシル酸塩錠 567 カンデサルタン シレキセチル・ヒド ロクロロチアジド錠 570 カンデサルタン シレキセチル錠 566 カンテン 1792 カンテン末 1793 含糖ペプシン 1535 d-カンフル 563 dl-カンフル 564 肝油 745 カンレノ酸カリウム 1433

丰

希塩酸 1024 キキョウ 1942, 2779, 2980 キキョウ末 1943, 2980 キキョウ流エキス 1943, 2980 キクカ 1829, 2975 キササゲ 1826 キジツ 1880 キシリトール 1778 キシリトール注射液 1778 キタサマイシン 1126 キタサマイシン酢酸エステル 1127 キタサマイシン酒石酸塩 1128 キナプリル塩酸塩 1496 キナプリル塩酸塩錠 1497 キニジン硫酸塩水和物 1498 キニーネエチル炭酸エステル 1499 キニーネ塩酸塩水和物 1500 キニーネ硫酸塩水和物 1501 牛脂 1808 吸水クリーム 755 キョウカツ 1919 キョウニン 1798 キョウニン水 1799 希ヨードチンキ 1076 金チオリンゴ酸ナトリウム 1564

ク

グアイフェネシン 1002 グアナベンズ酢酸塩 1003 グアネチジン硫酸塩 1003 グアヤコールスルホン酸カリウム 1436 クエチアピンフマル酸塩 1492 クエチアピンフマル酸塩細粒 1493 クエチアピンフマル酸塩錠 1494 クエン酸ガリウム(67Ga)注射液 976 クエン酸水和物 717 クエン酸ナトリウム水和物 1571 クコシ 1909 クジン 1995 クジン末 1995 苦味重曹水 1994 苦味チンキ 1813, 2973 クラブラン酸カリウム 1434 グラミシジン 1001, 2698 クラリスロマイシン 717, 2680 クラリスロマイシン錠 719 グリクラジド 984 グリシン 995 グリセリン 992 グリセリンカリ液 994 クリノフィブラート 725 グリベンクラミド 984 グリメピリド 985 グリメピリド錠 986 クリンダマイシン塩酸塩 721 クリンダマイシン塩酸塩カプセル クリンダマイシンリン酸エステル クリンダマイシンリン酸エステル注射 液 725 グルコン酸カルシウム水和物 552 グルタチオン 991 L-グルタミン 990 L-グルタミン酸 989 クレゾール 756 クレゾール水 756 クレゾール石ケン液 757 クレボプリドリンゴ酸塩 720 クレマスチンフマル酸塩 721 クロカプラミン塩酸塩水和物 727 クロキサシリンナトリウム水和物 743, 2683 クロキサゾラム 744 クロコナゾール塩酸塩 757

クロスカルメロースナトリウム 585

クロスポビドン 758

クロチアゼパム 741 クロチアゼパム錠 2682 クロトリマゾール 742 クロナゼパム 732 クロナゼパム細粒 733 クロナゼパム錠 733 クロニジン塩酸塩 734 クロピドグレル硫酸塩 736 クロピドグレル硫酸塩錠 738 クロフィブラート 728 クロフィブラートカプセル 729 クロフェダノール塩酸塩 728 クロベタゾールプロピオン酸エステル 726 クロペラスチン塩酸塩 735, 2919 クロミフェンクエン酸塩 730 クロミフェンクエン酸塩錠 731 クロミプラミン塩酸塩 731 クロミプラミン塩酸塩錠 2681 クロム酸ナトリウム(51Cr)注射液 1571 クロモグリク酸ナトリウム 1573 クロラゼプ酸二カリウム 739 クロラゼプ酸二カリウムカプセル 740 クロラムフェニコール 677, 2917 クロラムフェニコールコハク酸エステ ルナトリウム 678, 2679 クロラムフェニコール・コリスチンメ タンスルホン酸ナトリウム点眼液 2679 クロラムフェニコールパルミチン酸エ ステル 677 クロルジアゼポキシド 679 クロルジアゼポキシド散 680 クロルジアゼポキシド錠 681 クロルフェニラミンマレイン酸塩 687 d-クロルフェニラミンマレイン酸塩 691 クロルフェニラミンマレイン酸塩散 クロルフェニラミンマレイン酸塩錠 690 クロルフェニラミンマレイン酸塩注射 液 688 クロルフェネシンカルバミン酸エステ ル 685 クロルフェネシンカルバミン酸エステ ル錠 686 クロルプロパミド 694 クロルプロパミド錠 695 クロルプロマジン塩酸塩 692,2918 クロルプロマジン塩酸塩錠 693 クロルプロマジン塩酸塩注射液 692 クロルヘキシジン塩酸塩 683 クロルヘキシジングルコン酸塩液 682 クロルマジノン酢酸エステル 684

クロロブタノール 685

ケ

ケイガイ 1973 経口生ポリオワクチン 1428 ケイ酸アルミン酸マグネシウム 2715 ケイ酸マグネシウム 1187 軽質無水ケイ酸 1552, 2955 軽質流動パラフィン 1365 桂枝茯苓丸エキス 1901, 2775 ケイヒ 1830 ケイヒ末 1831 ケイヒ油 1831 ケタミン塩酸塩 1120 結晶セルロース 664, 2914 ケツメイシ 1825, 2974 ケトコナゾール 1121 ケトコナゾール液 1123 ケトコナゾールクリーム 1122 ケトコナゾールローション 1123 ケトチフェンフマル酸塩 1125 ケトプロフェン 1124 ケノデオキシコール酸 675 ゲファルナート 977 ケンゴシ **1938** ゲンタマイシン硫酸塩 982,2695 ゲンタマイシン硫酸塩注射液 2933 ゲンタマイシン硫酸塩点眼液 983 ゲンタマイシン硫酸塩軟膏 2933 ゲンチアナ 1856 ゲンチアナ・重曹散 1857 ゲンチアナ末 1856 ゲンノショウコ 1857 ゲンノショウコ末 1857

\Box

コウイ 1903 コウカ 1962, 2981 硬化油 1032 コウジン 1953 合成ケイ酸アルミニウム 405 コウブシ 1842, 2753 コウブシ末 1843, 2753 コウベイ 1820 コウボク **1910** コウボク末 1910 ゴオウ 1928, 2778, 2979 コカイン塩酸塩 745 ゴシツ 1792 牛車腎気丸エキス 1865, 2758 ゴシュユ 1851, 2755 呉茱萸湯エキス 2976 コデインリン酸塩散1% 746 コデインリン酸塩散10% 747 コデインリン酸塩錠 748 コデインリン酸塩水和物 746 ゴナドレリン酢酸塩 996 ゴボウシ 1821 ゴマ 1983 ゴマ油 1984 ゴミシ 1973 コムギデンプン 1604

コメデンプン 1603 コリスチンメタンスルホン酸ナトリウ サンヤク末 1847 스 752, 2683 コリスチン硫酸塩 753 コルチゾン酢酸エステル 754 コルヒチン 749 五苓散エキス 2756 コレカルシフェロール 696 コレスチミド 750 コレスチミド顆粒 751 コレスチミド錠 752 コレステロール 697, 2918 コレラワクチン 696 コロンボ 1821, 2974 コロンボ末 1822, 2974 コンズランゴ 1836 コンズランゴ流エキス 1837

サ

サイクロセリン 763 サイコ 1820 柴胡桂枝湯エキス 1965, 2785 サイシン 1801 柴朴湯エキス 1963, 2783 **柴苓湯エキス 1968, 2787** 酢酸 365 酢酸ナトリウム水和物 1563 サッカリン 1535, 2955 サッカリンナトリウム水和物 1536、 2734, 2955 サフラン 1963 サラシ粉 683 サラシミツロウ 1808 サラゾスルファピリジン 1537 サリチル・ミョウバン散 1542 サリチル酸 1539 サリチル酸精 1540 サリチル酸ナトリウム 1590 サリチル酸絆創膏 1540 サリチル酸メチル 1238 ザルトプロフェン 1780 ザルトプロフェン錠 1780 サルブタモール硫酸塩 1538 サルポグレラート塩酸塩 1543 サルポグレラート塩酸塩細粒 1544 サルポグレラート塩酸塩錠 1545 酸化亜鉛 1783 酸化カルシウム 554 酸化チタン 1693 酸化マグネシウム 1186 サンキライ 1993 サンキライ末 1993 サンザシ 1842 三酸化二ヒ素 445 サンシシ 1854, 2755 サンシシ末 1854 サンシュユ 1839, 2752 サンショウ 1885 サンショウ末 1886

酸素 1350

サントニン 1542

サンソウニン 1886, 2765

サンヤク 1847

ジアスターゼ 782 ジアスターゼ・重曹散 783 ジアゼパム 783 ジアゼパム錠 784 シアナミド 759 シアノコバラミン 760 シアノコバラミン注射液 761 ジエチルカルバマジンクエン酸塩 790 ジエチルカルバマジンクエン酸塩錠 790 ジオウ 1954 歯科用アンチホルミン 436 歯科用トリオジンクパスタ 1734 歯科用パラホルムパスタ 1367 歯科用フェノール・カンフル 1392 歯科用ヨード・グリセリン 1078 ジギトキシン 794, 2684 ジギトキシン錠 795, 2684 シクラシリン 699 ジクロキサシリンナトリウム水和物 789 シクロスポリン 700 ジクロフェナクナトリウム 787 ジクロフェナクナトリウム坐剤 2920 シベレスタットナトリウム水和物 ジクロフェナミド 787, 2684 ジクロフェナミド錠 789, 2684 シクロペントラート塩酸塩 761 シクロホスファミド錠 762 シクロホスファミド水和物 762 シゴカ 1848 ジゴキシン 796, 2685 ジゴキシン錠 798 ジゴキシン注射液 797 ジコッピ 1909, 2979 シコン 1907, 2979 次硝酸ビスマス 513 ジスチグミン臭化物 818 ジスチグミン臭化物錠 819 L-シスチン 766 L-システイン 765 L-システイン塩酸塩水和物 766 シスプラチン 714 ジスルフィラム 820 ジソピラミド 818 シタグリプチンリン酸塩錠 2957 シタグリプチンリン酸塩水和物 2955 臭化カリウム 1432 シタラビン 767 シチコリン 715 シツリシ 2005 ジドブジン 1781 ジドロゲステロン 844 ジドロゲステロン錠 845 シノキサシン 709 シノキサシンカプセル 710 ジノスタチン スチマラマー 1786,

2744

ジノプロスト 812

ジヒドロエルゴタミンメシル酸塩 ジヒドロエルゴトキシンメシル酸塩 803 ジヒドロコデインリン酸塩 800 ジヒドロコデインリン酸塩散1% 800 ジヒドロコデインリン酸塩散10% 801 ジピリダモール 817 ジフェニドール塩酸塩 791 ジフェンヒドラミン 813 ジフェンヒドラミン・バレリル尿素散 814 ジフェンヒドラミン・フェノール・亜鉛 華リニメント 815 ジフェンヒドラミン塩酸塩 814 ジブカイン塩酸塩 786 ジフテリアトキソイド 816 ジフテリア破傷風混合トキソイド 816 ジフルコルトロン吉草酸エステル 793 シプロフロキサシン 711 シプロフロキサシン塩酸塩水和物 712 シプロヘプタジン塩酸塩水和物 764 ジフロラゾン酢酸エステル 792 ジベカシン硫酸塩 785 ジベカシン硫酸塩点眼液 785 1561 シベンゾリンコハク酸塩 697 シベンゾリンコハク酸塩錠 698 シメチジン 709 ジメモルファンリン酸塩 808 ジメルカプロール 810 ジメルカプロール注射液 811 ジメンヒドリナート 809 ジメンヒドリナート錠 809 次没食子酸ビスマス 512 ジモルホラミン 811 ジモルホラミン注射液 812 シャカンゾウ 1947 弱アヘンアルカロイド・スコポラミン 注射液 1340 シャクヤク 1935 芍薬甘草湯エキス 1984, 2790 シャクヤク末 1936 ジャショウシ 1834 シャゼンシ 1942 シャゼンソウ 1942 臭化ナトリウム 1568 十全大補湯エキス 1887, 2766, 2978 ジュウヤク 1880 シュクシャ 1797, 2973 シュクシャ末 1797 酒石酸 1643 ショウキョウ 1857 ショウキョウ末 1858 小柴胡湯エキス 1988, 2791 硝酸イソソルビド 1104 硝酸イソソルビド錠 1105

硝酸銀 1557 硝酸銀点眼液 1557 常水 1773 ショウズク 1825 小青竜湯エキス 1990, 2793 焼セッコウ 1868 消毒用エタノール 897 消毒用フェノール 1391 消毒用フェノール水 1392 ショウマ 1830 ジョサマイシン 1112 ジョサマイシン錠 1113 ジョサマイシンプロピオン酸エステル 1113 シラザプリル錠 703 シラザプリル水和物 702 シラスタチンナトリウム 701 ジラゼプ塩酸塩水和物 805 ジルチアゼム塩酸塩 806 ジルチアゼム塩酸塩徐放カプセル シルニジピン 704 シルニジピン錠 705 シロスタゾール 707 シロスタゾール錠 **708** シロップ用アシクロビル 373 シロップ用クラリスロマイシン 2918 シロップ用セファトリジンプロピレン グリコール 606 シロップ用セファドロキシル 599 シロップ用セファレキシン 603 シロップ用セフポドキシム プロキセ チル 646 シロップ用セフロキサジン 649 シロップ用トラニラスト 1715 シロップ用ファロペネムナトリウム 920 シロップ用ペミロラストカリウム 1374 シロップ用ホスホマイシンカルシウム シロドシン 1553 シロドシン錠 1555 シンイ 1911

ス

診断用クエン酸ナトリウム液 1572

シンギ 1875

親水クリーム 756

親水ワセリン 1387

シンバスタチン 1559

シンバスタチン錠 1560

真武湯エキス 1985, 2791

水酸化カリウム 1436 水酸化カルシウム 553 水酸化ナトリウム 1578 スキサメトニウム塩化物水和物 1627 スキサメトニウム塩化物注射液 1628, 2737 スクラルファート水和物 1609 スコポラミン臭化水素酸塩水和物 1548 ステアリルアルコール 1607 ステアリン酸 1606 ステアリン酸カルシウム 562 ステアリン酸ポリオキシル40 1430 ステアリン酸マグネシウム 1188 ストレプトマイシン硫酸塩 1608 スピラマイシン酢酸エステル 1599, 2735 スピロノラクトン 1600 スピロノラクトン錠 1600 スペクチノマイシン塩酸塩水和物 1597 スリンダク 1620 スルタミシリントシル酸塩錠 1625 スルタミシリントシル酸塩水和物 1624, 2736 スルチアム 1626 スルバクタムナトリウム 1612,2735 スルピリド 1621 スルピリドカプセル 1621 スルピリド錠 1622 スルピリン水和物 1623 スルピリン注射液 1623 スルファジアジン銀 1614 スルファメチゾール 1615 スルファメトキサゾール 1616 スルファモノメトキシン水和物 1616 スルフイソキサゾール 1617 スルベニシリンナトリウム 1613 スルホブロモフタレインナトリウム 1618 スルホブロモフタレインナトリウム注 射液 1619

t

成人用沈降ジフテリアトキソイド 816 精製水 1773 精製水(容器入り) 1774 精製ゼラチン 980 精製セラック 1551 精製デヒドロコール酸 772, 2920 精製白糖 1611 精製ヒアルロン酸ナトリウム 1575 精製ヒアルロン酸ナトリウム注射液 1576 精製ヒアルロン酸ナトリウム点眼液 1577 精製ブドウ糖 2696 精製ラノリン 1904 生理食塩液 1571 石油ベンジン 1387 セタノール 671 セチリジン塩酸塩 671 セチリジン塩酸塩錠 672 セッコウ 1868 セトチアミン塩酸塩水和物 673 セトラキサート塩酸塩 674 セネガ 1980 セネガシロップ 1981 セネガ末 1981 セファクロル 592

セファクロルカプセル 593 セファクロル細粒 596 セファクロル複合顆粒 594 セファゾリンナトリウム 607 セファゾリンナトリウム水和物 609 セファトリジンプロピレングリコール 605 セファドロキシル 597 セファドロキシルカプセル 598 セファレキシン 599 セファレキシンカプセル 601 セファレキシン複合顆粒 602 セファロチンナトリウム 604 セフィキシムカプセル 624 セフィキシム細粒 2913 セフィキシム水和物 623,2676 セフェピム塩酸塩水和物 621 セフォジジムナトリウム 629 セフォゾプラン塩酸塩 640 セフォタキシムナトリウム 634 セフォチアム塩酸塩 639 セフォチアム ヘキセチル塩酸塩 637 セフォテタン 635 セフォペラゾンナトリウム 631 セフカペン ピボキシル塩酸塩細粒 613 セフカペン ピボキシル塩酸塩錠 614 セフカペン ピボキシル塩酸塩水和物 612 セフジトレン ピボキシル 618 セフジトレン ピボキシル細粒 619 セフジトレン ピボキシル錠 620 セフジニル 615 セフジニルカプセル 617 セフジニル細粒 617 セフスロジンナトリウム 650 セフタジジム水和物 652 セフチゾキシムナトリウム 658,2678 セフチブテン水和物 657 セフテラム ピボキシル 654 セフテラム ピボキシル細粒 655 セフテラム ピボキシル錠 656 セフトリアキソンナトリウム水和物 660 セフピラミドナトリウム 642 セフピロム硫酸塩 643 セフブペラゾンナトリウム 610 セフポドキシム プロキセチル 644 セフポドキシム プロキセチル錠 647 セフミノクスナトリウム水和物 629 セフメタゾールナトリウム 627 セフメノキシム塩酸塩 625 セフロキサジン水和物 648 セフロキシム アキセチル 662 セボフルラン 1550 セラセフェート 663, 2679 ゼラチン 978

セラペプターゼ

L-セリン 1548

1549, 2734

セルモロイキン(遺伝子組換え) 668

センキュウ 1835 センキュウ末 1835 ゼンコ 1937 センコツ 1919 センソ 1998, 2983 センナ 1981 センナ末 1982 センブリ 1996, 2982 センブリ・重曹散 1998, 2982 センブリ末 1997, 2982

ソ

ソウジュツ 1803 ソウジュツ末 1804 ソウハクヒ 1918 ゾニサミド 2744 ゾニサミド錠 2745 ソボク 1972 ソヨウ 1937 ソルビタンセスキオレイン酸エステル 1595 ゾルピデム酒石酸塩 1788 ゾルピデム酒石酸塩錠 1789 D-ソルビトール 1596 D-ソルビトール液 1597

夕

ダイオウ 1955 大黄甘草湯エキス 1843, 2754 ダイオウ末 1956 大柴胡湯エキス 1844, 2754 ダイズ油 1995 タイソウ 1886, 2765 ダウノルビシン塩酸塩 769 タウリン 1643 タカルシトール水和物 1629 タカルシトール軟膏 1631 タカルシトールローション 1630 タクシャ 1793 タクシャ末 1794, 2747 タクロリムスカプセル 1633 タクロリムス水和物 1632 タゾバクタム 1644 ダナゾール 768 タムスロシン塩酸塩 1640 タムスロシン塩酸塩徐放錠 1642 タモキシフェンクエン酸塩 1640 タランピシリン塩酸塩 1633 タルク 1634 タルチレリン口腔内崩壊錠 1637 タルチレリン錠 1638 タルチレリン水和物 1636 炭酸カリウム 1433 炭酸水素ナトリウム 1566 炭酸水素ナトリウム注射液 1566 炭酸ナトリウム水和物 1569 炭酸マグネシウム 1185 炭酸リチウム 1169 単シロップ 1558 タンジン 1971 ダントロレンナトリウム水和物 768 単軟膏 1993 タンニン酸 1643 タンニン酸アルブミン 386 タンニン酸ジフェンヒドラミン 815 タンニン酸ベルベリン 494

チ

チアプリド塩酸塩 1684 チアプリド塩酸塩錠 1685 チアマゾール 1672 チアマゾール錠 1673 チアミラールナトリウム 1677 チアミン塩化物塩酸塩 1674 チアミン塩化物塩酸塩散 1675 チアミン塩化物塩酸塩注射液 1675 チアミン硝化物 1676 チアラミド塩酸塩 1685 チアラミド塩酸塩錠 1686 チアントール 1678 チオペンタールナトリウム 1680 チオリダジン塩酸塩 1681 チオ硫酸ナトリウム水和物 1592 チオ硫酸ナトリウム注射液 1593 チクセツニンジン 1932 チクセツニンジン末 1933 チクロピジン塩酸塩 1687 チクロピジン塩酸塩錠 1688 チザニジン塩酸塩 1694 窒素 1316 チニダゾール 1691 チペピジンヒベンズ酸塩 1691 チペピジンヒベンズ酸塩錠 1692, 2967 チメピジウム臭化物水和物 1689 チモ 1797 チモール 1683 チモロールマレイン酸塩 1690 注射用アシクロビル 371 注射用アズトレオナム 467 注射用アセチルコリン塩化物 367 注射用アミカシン硫酸塩 409 注射用アムホテリシンB 427 注射用アンピシリンナトリウム 432 注射用アンピシリンナトリウム・スル バクタムナトリウム 433 注射用イダルビシン塩酸塩 1047 注射用イミペネム・シラスタチンナト リウム 1056 注射用オザグレルナトリウム 1357 注射用シベレスタットナトリウム 1562 注射用水 1774 注射用水(容器入り) 1775 注射用スキサメトニウム塩化物 1628 注射用ストレプトマイシン硫酸塩 1609 注射用スペクチノマイシン塩酸塩 1598 注射用セファゾリンナトリウム 608 注射用セファロチンナトリウム 2912

注射用セフェピム塩酸塩 622

注射用セフォゾプラン塩酸塩 641

注射用セフォチアム塩酸塩 640 注射用セフォペラゾンナトリウム 2677 注射用セフォペラゾンナトリウム・ス ルバクタムナトリウム 632 注射用セフタジジム 653 注射用セフメタゾールナトリウム 628 注射用タゾバクタム・ピペラシリン 注射用チアミラールナトリウム 1678 注射用チオペンタールナトリウム 1681 注射用テセロイキン(遺伝子組換え) 1652 注射用ドキソルビシン塩酸塩 837 注射用ドセタキセル 823 注射用ドリペネム 2924 注射用ナルトグラスチム(遺伝子組換 え) 1293 注射用パニペネム・ベタミプロン 1361 注射用バンコマイシン塩酸塩 1758 注射用ヒト絨毛性性腺刺激ホルモン 999 注射用ヒドララジン塩酸塩 1022 注射用ピペラシリンナトリウム 1418 注射用ビンブラスチン硫酸塩 1763 注射用ファモチジン 917 注射用フェニトインナトリウム 1398 注射用プレドニゾロンコハク酸エステ ルナトリウム 1456 注射用フロモキセフナトリウム 938 注射用ペプロマイシン硫酸塩 1380 注射用ベンジルペニシリンカリウム 487 注射用ホスホマイシンナトリウム 965 注射用ボリコナゾール 2743 注射用マイトマイシン C 1262 注射用ミノサイクリン塩酸塩 1256 注射用メロペネム 1219 注射用ロキサチジン酢酸エステル塩酸 塩 1533 チョウジ 1833 チョウジ末 1833 チョウジ油 1834 チョウトウコウ 2007 釣藤散エキス 1827, 2751 チョレイ 1946 チョレイ末 1946 L-チロシン 1740 チンク油 1784 沈降ジフテリア破傷風混合トキソイド 816 沈降精製百日せきジフテリア破傷風混 合ワクチン 816 沈降精製百日せきワクチン 1384 沈降炭酸カルシウム 548 沈降炭酸カルシウム細粒 548 沈降炭酸カルシウム錠 549 沈降破傷風トキソイド 1670

沈降はぶトキソイド 1004, 2934

沈降 B 型肝炎ワクチン 1018 チンピ 1832

ツ

ツバキ油 1822 ツロブテロール 1738 ツロブテロール塩酸塩 1740 ツロブテロール経皮吸収型テープ 1739

テ

テイコプラニン 1654, 2737, 2964 低置換度ヒドロキシプロピルセルロー ス 1035, 2703 テオフィリン 1672 テガフール 1653 デキサメタゾン 776 デキストラン40 777, 2684 デキストラン40注射液 778 デキストラン70 779 デキストラン硫酸エステルナトリウム イオウ5 780 デキストラン硫酸エステルナトリウム イオウ18 780 デキストリン 781 デキストロメトルファン臭化水素酸塩 水和物 782 テストステロンエナント酸エステル 1667, 2967 テストステロンエナント酸エステル注 射液 1668 テストステロンプロピオン酸エステル 1668 テストステロンプロピオン酸エステル 注射液 1669 デスラノシド 775 デスラノシド注射液 776 テセロイキン(遺伝子組換え) 1647 テトラカイン塩酸塩 1670 テトラサイクリン塩酸塩 1670,2738 デヒドロコール酸 772,2920 デヒドロコール酸注射液 773 デフェロキサミンメシル酸塩 771 テプレノン 1660 テプレノンカプセル 1661 デメチルクロルテトラサイクリン塩酸 塩 774, 2684 テモカプリル塩酸塩 1658 テモカプリル塩酸塩錠 1659 テルビナフィン塩酸塩 1662 テルビナフィン塩酸塩液 1664 テルビナフィン塩酸塩クリーム 1663 テルビナフィン塩酸塩錠 1665 テルビナフィン塩酸塩スプレー 1665 テルブタリン硫酸塩 1666 テルミサルタン 1656 テルミサルタン錠 1657 テルミサルタン・ヒドロクロロチアジ 下錠 2964 テレビン油 2007

天然ケイ酸アルミニウム 404

デンプングリコール酸ナトリウム 1591 テンマ 1855, 2975 テンモンドウ 1802

卜 桃核承気湯エキス 1999, 2796 トウガシ 1811 トウガラシ 1822 トウガラシ・サリチル酸精 1824 トウガラシチンキ 1824 トウガラシ末 1823 トウキ 1883, 2977 当帰芍薬散エキス 2002, 2797, 2983 トウキ末 1883, 2978 トウジン 1835 トウニン 1933 トウニン末 1934 トウヒ 1812 トウヒシロップ 1926 トウヒチンキ 1926 トウモロコシデンプン 1601 トウモロコシ油 1839 ドキサゾシンメシル酸塩 833 ドキサゾシンメシル酸塩錠 834 ドキサプラム塩酸塩水和物 832 ドキシサイクリン塩酸塩錠 840 ドキシサイクリン塩酸塩水和物 838, 2686 ドキシフルリジン 835 ドキシフルリジンカプセル 835 ドキソルビシン塩酸塩 836,2685 ドクカツ 1800 トコフェロール 1697 トコフェロールコハク酸エステルカル シウム 1698 トコフェロール酢酸エステル 1697 トコフェロールニコチン酸エステル 1699 トコン 1881 トコンシロップ 1882 トコン末 1881 トスフロキサシントシル酸塩錠 1707 トスフロキサシントシル酸塩水和物 1706 ドセタキセル水和物 821 ドセタキセル注射液 822 トチュウ 1851 トドララジン塩酸塩水和物 1700 ドネペジル塩酸塩 825 ドネペジル塩酸塩細粒 826 ドネペジル塩酸塩錠 827 ドパミン塩酸塩 829 ドパミン塩酸塩注射液 829 トフィソパム 1701 ドブタミン塩酸塩 820 トブラマイシン 1695, 2738 トブラマイシン注射液 1696 トラガント 2004

トラガント末 2004

トラニラスト 1711

トラザミド 1702, 2738

トラニラストカプセル 1712 トラニラスト細粒 1713 トラニラスト点眼液 1714 トラネキサム酸 1708 トラネキサム酸カプセル 1709 トラネキサム酸錠 1711 トラネキサム酸注射液 1710 トラピジル 1716 トラマドール塩酸塩 2739 トリアムシノロン 1718 トリアムシノロンアセトニド 1719, 2967 トリアムテレン 1720 トリエンチン塩酸塩 1726 トリエンチン塩酸塩カプセル 1727 トリクロホスナトリウム 1725 トリクロホスナトリウムシロップ 1726 トリクロルメチアジド 1721 トリクロルメチアジド錠 1722 トリコマイシン 1724 L-トリプトファン 1737 トリヘキシフェニジル塩酸塩 1727 トリヘキシフェニジル塩酸塩錠 1728 ドリペネム水和物 2921 トリメタジオン 1732 トリメタジジン塩酸塩 1730 トリメタジジン塩酸塩錠 1731 トリメトキノール塩酸塩水和物 1733 トリメブチンマレイン酸塩 1729 ドルゾラミド塩酸塩 830 ドルゾラミド塩酸塩点眼液 831 トルナフタート 1704 トルナフタート液 1705 トルブタミド 1703 トルブタミド錠 1703 トルペリゾン塩酸塩 1705 L-トレオニン 1682 トレハロース水和物 1717 トレピブトン 1718 ドロキシドパ 842 ドロキシドパカプセル 842 ドロキシドパ細粒 843 トロキシピド 1735 トロキシピド細粒 1735 トロキシピド錠 1736 トロピカミド 1734 ドロペリドール 841 トロンビン 1682, 2738 豚脂 1905 ドンペリドン 824

ナ

ナイスタチン 1328 ナタネ油 1953, 2780 ナテグリニド 1294 ナテグリニド錠 1295 ナドロール 1283 ナファゾリン・クロルフェニラミン液 1290 ナファゾリン塩酸塩 1289 ナファゾリン硝酸塩 1289 ナファモスタットメシル酸塩 1283 ナフトピジル 1284 ナフトピジルロ腔内崩壊錠 1285 ナフトピジル錠 1286 ナブメトン 1281 ナブメトン錠 1282 ナプロキセン 1291 ナリジクス酸 1287 ナルトグラスチム(遺伝子組換え) 1291 ナロキソン塩酸塩 1288

=

ニガキ 1941 ニガキ末 1941 ニカルジピン塩酸塩 1298 ニカルジピン塩酸塩注射液 1299 ニクジュヨウ 1831 ニクズク 1920 ニコチン酸 1307 ニコチン酸アミド 1306 ニコチン酸注射液 1307 ニコモール 1304 ニコモール錠 1304 ニコランジル 1305 ニザチジン 1319 ニザチジンカプセル 1320 二酸化炭素 580 ニセリトロール 1303 ニセルゴリン 1300 ニセルゴリン散 1301 ニセルゴリン錠 1302 二相性イソフェンインスリン ヒト (遺伝子組換え)水性懸濁注射液 2706 ニトラゼパム 1314 ニトレンジピン 1314 ニトレンジピン錠 1315 ニトログリセリン錠 1317 ニフェジピン 1308 ニフェジピン細粒 1311 ニフェジピン徐放カプセル 1310 ニフェジピン腸溶細粒 1309 日本脳炎ワクチン 1111 乳酸 1131 L-乳酸 1131 乳酸カルシウム水和物 553 L-乳酸ナトリウム液 1581 L-乳酸ナトリウムリンゲル液 1582 乳糖水和物 1133, 2715 尿素 1747 ニルバジピン 1312 ニルバジピン錠 1313 ニンジン 1859 ニンジン末 1860

ネ

ニンドウ 1908

ネオスチグミンメチル硫酸塩 1297 ネオスチグミンメチル硫酸塩注射液 1297

ノ

濃グリセリン 993
濃ベンザルコニウム塩化物液50 480
ノスカピン 1326
ノスカピン塩酸塩水和物 1327
ノルアドレナリン 1321, 2726
ノルアドレナリン注射液 1322
ノルエチステロン 1322
ノルゲストレル 1324
ノルゲストレル・エチニルエストラジオール錠 1324
ノルトリプチリン塩酸塩 1326
ノルトリプチリン塩酸塩錠 2951
ノルフロキサシン 1323

ハ

バイモ 1853 バカンピシリン塩酸塩 バクガ 1912 白色セラック 1552 白色軟膏 1775 白色ワセリン 1386 白糖 1611 バクモンドウ 1923 麦門冬湯エキス 1805, 2747 バクロフェン 470 バクロフェン錠 470 バシトラシン 469, 2674 パズフロキサシンメシル酸塩 2726 パズフロキサシンメシル酸塩注射液 2728 バソプレシン注射液 1759, 2740 八味地黄丸エキス 1868, 2759 ハチミツ 1879 ハッカ 1914 ハッカ水 1915 ハッカ油 1915 パップ用複方オウバク散 1940 パニペネム 1359 パパベリン塩酸塩 1363 パパベリン塩酸塩注射液 1364 ハマボウフウ 1861 バメタン硫酸塩 471 パラアミノサリチル酸カルシウム顆粒 556 パラアミノサリチル酸カルシウム水和 物 555 パラオキシ安息香酸エチル 906 パラオキシ安息香酸ブチル 540 パラオキシ安息香酸プロピル 1481 パラオキシ安息香酸メチル 1235 バラシクロビル塩酸塩 1751 バラシクロビル塩酸塩錠 1753 パラフィン 1364 パラホルムアルデヒド 1366 L-バリン 1754 バルサルタン 1755 バルサルタン錠 1756 バルサルタン・ヒドロクロロチアジド

錠 2968

パルナパリンナトリウム 1367 バルビタール 472 バルプロ酸ナトリウム 1593 バルプロ酸ナトリウム錠 1594 バルプロ酸ナトリウム徐放錠 A 2961 バルプロ酸ナトリウム徐放錠 B 2963 バルプロ酸ナトリウムシロップ 1594 バレイショデンプン 1602 ハロキサゾラム 1008 パロキセチン塩酸塩錠 1371 パロキセチン塩酸塩水和物 1369 ハロタン 1007 ハロペリドール 1004, 2934 ハロペリドール細粒 1005 ハロペリドール錠 1006 ハロペリドール注射液 1006 パンクレアチン 1358 パンクロニウム臭化物 1358 ハンゲ 1941 半夏厚朴湯エキス 1871, 2761 半夏瀉心湯エキス 1873, 2761 バンコマイシン塩酸塩 1757 パンテチン 1362 パントテン酸カルシウム 554

匕 ピオグリタゾン塩酸塩 1406 ピオグリタゾン塩酸塩・グリメピリド 錠 1408, 2952 ピオグリタゾン塩酸塩・メトホルミン 塩酸塩錠 1411 ピオグリタゾン塩酸塩錠 1407 ビオチン 510 ピコスルファートナトリウム水和物 1585 ビサコジル 511 ビサコジル坐剤 512 L-ヒスチジン 1018 L-ヒスチジン塩酸塩水和物 1019 ビソプロロールフマル酸塩 514 ビソプロロールフマル酸塩錠 515 ピタバスタチンカルシウム錠 1425 ピタバスタチンカルシウム水和物 1424 ビタミン A 油 1765 ヒト下垂体性性腺刺激ホルモン 999 ヒト絨毛性性腺刺激ホルモン 997 人全血液 1775 人免疫グロブリン 1021 ヒドララジン塩酸塩 1021 ヒドララジン塩酸塩散 1022 ヒドララジン塩酸塩錠 1022 ヒドロキシエチルセルロース ヒドロキシジン塩酸塩 1036 ヒドロキシジンパモ酸塩 1036 ヒドロキシプロピルセルロース 1033, 2702, 2937 ヒドロキソコバラミン酢酸塩 1033, 2701

ヒドロクロロチアジド 1024

ヒドロコタルニン塩酸塩水和物

ヒドロコルチゾン 1025, 2934

1031

ヒドロコルチゾン・ジフェンヒドラミ ン軟膏 1027, 2934 ヒドロコルチゾンコハク酸エステル ヒドロコルチゾンコハク酸エステルナ トリウム 1030 ヒドロコルチゾン酢酸エステル 1026, 2701, 2934 ヒドロコルチゾン酪酸エステル 1028, 2701 ヒドロコルチゾンリン酸エステルナト リウム 1028 ピブメシリナム塩酸塩 1427 ピブメシリナム塩酸塩錠 1428 ヒプロメロース 1038, 2704, 2938 ヒプロメロースカプセル 574 ヒプロメロース酢酸エステルコハク酸 エステル 1040 ヒプロメロースフタル酸エステル 1041 ピペミド酸水和物 1414 ピペラシリン水和物 1414 ピペラシリンナトリウム 1416 ピペラジンアジピン酸塩 1418 ピペラジンリン酸塩錠 1419 ピペラジンリン酸塩水和物 1419 ビペリデン塩酸塩 510 ビホナゾール 509 ヒマシ油 1826 ピマリシン 1404 ヒメクロモン 1037 ピモジド 1405 ビャクゴウ 1906 ビャクシ 1798 ビャクジュツ 1804 ビャクジュツ末 1805 氷酢酸 365 ピラジナミド 1488 ピラルビシン 1420 ピランテルパモ酸塩 1487 ピリドキサールリン酸エステル水和物 ピリドキシン塩酸塩 1489 ピリドキシン塩酸塩注射液 1490 ピリドスチグミン臭化物 1489 ピルシカイニド塩酸塩カプセル 1402 ピルシカイニド塩酸塩水和物 1402 ピレノキシン 1421 ピレンゼピン塩酸塩水和物 1422 ピロ亜硫酸ナトリウム 1587 ピロカルピン塩酸塩 1399 ピロカルピン塩酸塩錠 1400 ピロキシカム 1423 ピロキシリン 1491 ピロールニトリン 1491 ビワヨウ 1909 ビンクリスチン硫酸塩 1764, 2742 ピンドロール 1406 ビンブラスチン硫酸塩 1762,2742 ビンロウジ 1800

フ

ファモチジン 915 ファモチジン散 918 ファモチジン錠 919 ファモチジン注射液 915 ファロペネムナトリウム錠 921 ファロペネムナトリウム水和物 919 フィトナジオン 1398 フィルグラスチム(遺伝子組換え) 929 フィルグラスチム(遺伝子組換え)注射 液 931 フェキソフェナジン塩酸塩 926 フェキソフェナジン塩酸塩錠 927 フェニトイン 1396 フェニトイン散 1396 フェニトイン錠 1397 L-フェニルアラニン 1394 フェニルブタゾン 1395 フェニレフリン塩酸塩 1395 フェネチシリンカリウム 1388 フェノバルビタール 1389 フェノバルビタール散10% 1390 フェノール 1390 フェノール・亜鉛華リニメント 1392 フェノール水 1392 フェノールスルホンフタレイン 1393 フェノールスルホンフタレイン注射液 1393 フェルビナク 923 フェルビナクテープ 924 フェルビナクパップ 923 フェロジピン 2927 フェロジピン錠 2928 フェンタニルクエン酸塩 925 フェンブフェン 924 複方アクリノール・チンク油 377 複方オキシコドン・アトロピン注射液 1348 複方オキシコドン注射液 1347 複方サリチル酸精 1541 複方サリチル酸メチル精 1238 複方ジアスターゼ・重曹散 783 複方ダイオウ・センナ散 1957 複方チアントール・サリチル酸液 1679 複方ビタミンB散 1765,2971 複方ヨード・グリセリン 1077 複方ロートエキス・ジアスターゼ散 1976 ブクモロール塩酸塩 527 ブクリョウ 1946 ブクリョウ末 1947 ブシ 1948 フシジン酸ナトリウム **1574, 2959** ブシ末 1949 ブシラミン 525 ブシラミン錠 526 ブスルファン 536 ブチルスコポラミン臭化物 1547

ブテナフィン塩酸塩 536

ブテナフィン塩酸塩液 538 ブテナフィン塩酸塩クリーム 537 ブテナフィン塩酸塩スプレー 538 ブドウ酒 1776 ブドウ糖 988 ブドウ糖水和物 2695 ブドウ糖注射液 989,2698 フドステイン 968 フドステイン錠 969 ブトロピウム臭化物 539 ブナゾシン塩酸塩 532 ブピバカイン塩酸塩水和物 ブフェトロール塩酸塩 528 ブプラノロール塩酸塩 534 ブプレノルフィン塩酸塩 535 ブホルミン塩酸塩 529 ブホルミン塩酸塩錠 531 ブホルミン塩酸塩腸溶錠 529 ブメタニド 531 フラジオマイシン硫酸塩 966 プラステロン硫酸エステルナトリウム 水和物 1444 プラゼパム 1450 プラゼパム錠 1450 プラゾシン塩酸塩 1451 プラノプロフェン 1443 プラバスタチンナトリウム 1444 プラバスタチンナトリウム液 1447 プラバスタチンナトリウム細粒 1446 プラバスタチンナトリウム錠 1448 フラビンアデニンジヌクレオチドナト リウム 932 フラボキサート塩酸塩 934 プランルカスト水和物 1442 プリミドン 1458 フルオキシメステロン 950, 2694 フルオシノニド 947 フルオシノロンアセトニド 946 フルオレセインナトリウム 948 フルオロウラシル 950 フルオロメトロン 949 フルコナゾール 940 フルコナゾールカプセル 941 フルコナゾール注射液 942 フルジアゼパム 944 フルシトシン 943 フルスルチアミン塩酸塩 973 フルタミド 954 フルトプラゼパム 955 フルトプラゼパム錠 955 フルドロコルチゾン酢酸エステル 944 フルニトラゼパム 945 フルフェナジンエナント酸エステル 951 フルボキサミンマレイン酸塩 956 フルボキサミンマレイン酸塩錠 958 フルラゼパム塩酸塩 952 プルラン 1487 プルランカプセル 574 フルルビプロフェン 953

ブレオマイシン塩酸塩 516

ブレオマイシン硫酸塩 518

フレカイニド酢酸塩 934 フレカイニド酢酸塩錠 935 プレドニゾロン 1452 プレドニゾロンコハク酸エステル プレドニゾロン酢酸エステル 1454 プレドニゾロン錠 1453 プレドニゾロンリン酸エステルナトリ ウム 1455 プロカインアミド塩酸塩 1462 プロカインアミド塩酸塩錠 1463 プロカインアミド塩酸塩注射液 1463 プロカイン塩酸塩 1464 プロカイン塩酸塩注射液 1465 プロカテロール塩酸塩水和物 1466 プロカルバジン塩酸塩 1466 プログルミド 1470 プロクロルペラジンマレイン酸塩 1467 プロクロルペラジンマレイン酸塩錠 1468 プロゲステロン 1469 プロゲステロン注射液 1470 フロセミド 970 フロセミド錠 972 フロセミド注射液 971 プロタミン硫酸塩 1483 プロタミン硫酸塩注射液 1484 プロチオナミド 1484 ブロチゾラム 523 ブロチゾラム錠 524 プロチレリン 1485 プロチレリン酒石酸塩水和物 1486 プロテイン銀 1557 プロテイン銀液 1558 プロパフェノン塩酸塩 1473 プロパフェノン塩酸塩錠 1474 プロパンテリン臭化物 1475 プロピベリン塩酸塩 1476 プロピベリン塩酸塩錠 1477 プロピルチオウラシル 1482 プロピルチオウラシル錠 1482 プロピレングリコール 1480 プロブコール 1460 プロブコール細粒 1461 プロブコール錠 1462 プロプラノロール塩酸塩 1478 プロプラノロール塩酸塩錠 1479 フロプロピオン 939 フロプロピオンカプセル 939 プロベネシド 1458 プロベネシド錠 1459 ブロマゼパム 520 ブロムフェナクナトリウム水和物 2909 ブロムフェナクナトリウム点眼液 2910 ブロムヘキシン塩酸塩 521 プロメタジン塩酸塩 1472 フロモキセフナトリウム 936 ブロモクリプチンメシル酸塩 522 ブロモバレリル尿素 523

L-プロリン 1471

粉末セルロース 667

^

ベカナマイシン硫酸塩 474 ベクロメタゾンプロピオン酸エステル 473, 2909 ベザフィブラート 507 ベザフィブラート徐放錠 508 ベタキソロール塩酸塩 506 ベタネコール塩化物 507 ベタヒスチンメシル酸塩 495 ベタヒスチンメシル酸塩錠 496 ベタミプロン 505 ベタメタゾン 497 ベタメタゾン吉草酸エステル 502 ベタメタゾン吉草酸エステル・ゲンタ マイシン硫酸塩クリーム 502 ベタメタゾン吉草酸エステル・ゲンタ マイシン硫酸塩軟膏 504 ベタメタゾンジプロピオン酸エステル 499, 2909 ベタメタゾン錠 498 ベタメタゾンリン酸エステルナトリウ لم **500** ペチジン塩酸塩 1385 ペチジン塩酸塩注射液 1385 ベニジピン塩酸塩 475 ベニジピン塩酸塩錠 476 ヘパリンカルシウム 1009, 2698 ヘパリンナトリウム 1013, 2699 ヘパリンナトリウム注射液 1017 ペプロマイシン硫酸塩 1379 ベポタスチンベシル酸塩 488 ベポタスチンベシル酸塩錠 489 ペミロラストカリウム 1372 ペミロラストカリウム錠 1375 ペミロラストカリウム点眼液 1373 ベラドンナエキス 1810 ベラドンナコン 1809, 2973 ベラドンナ総アルカロイド 1810 ベラパミル塩酸塩 1761,2741 ベラパミル塩酸塩錠 1761,2741 ベラパミル塩酸塩注射液 2970 ベラプロストナトリウム 490 ベラプロストナトリウム錠 492 ペルフェナジン 1381 ペルフェナジン錠 1382 ペルフェナジンマレイン酸塩 1383 ペルフェナジンマレイン酸塩錠 1383 ベルベリン塩化物水和物 493 ベンザルコニウム塩化物 479 ベンザルコニウム塩化物液 480 ベンジルアルコール 483 ベンジルペニシリンカリウム 486, 2674 ベンジルペニシリンベンザチン水和物 485 ヘンズ 1848

ベンズブロマロン 481

ベンゼトニウム塩化物 481

ベンセラジド塩酸塩 478

ベンゼトニウム塩化物液 482

ペンタゾシン 1376 ペントキシベリンクエン酸塩 1378 ベントナイト 478 ペントバルビタールカルシウム 1377, 2728 ペントバルビタールカルシウム錠 2729 ペンプトロール硫酸塩 1376

木

ボウイ 1993, 2982 防已黄耆湯エキス 1817, 2750, 2974 ボウコン 1880 ホウ酸 520 ホウ砂 1567 ボウショウ 1994 抱水クロラール 676 ボウフウ 1971, 2789, 2981 防風通聖散エキス 1813, 2748, 2973 ボクソク 1952 ボグリボース 1766 ボグリボース錠 1767 ホスホマイシンカルシウム水和物 962, 2694 ホスホマイシンナトリウム 964, 2694 ボタンピ 1915 ボタンピ末 1916 補中益気湯エキス 1876, 2763, 2977 ポビドン 1439, 2731 ポビドンヨード 1441 ホマトロピン臭化水素酸塩 1020 ホミカ 1920, 2979 ホミカエキス 1921 ホミカエキス散 1922 ホミカチンキ 1922 ホモクロルシクリジン塩酸塩 1020 ポラプレジンク 2952 ポラプレジンク顆粒 2953 ボリコナゾール 1769 ボリコナゾール錠 1770 ポリスチレンスルホン酸カルシウム 559 ポリスチレンスルホン酸ナトリウム 1586 ポリソルベート80 1430, 2730 ホリナートカルシウム 551 ポリミキシン B 硫酸塩 1429, 2730 ホルマリン 961 ホルマリン水 961 ホルモテロールフマル酸塩水和物 962 ボレイ 1932 ボレイ末 1932

マ

マイトマイシン C 1261 マオウ 1849 麻黄湯エキス 1912, 2776 マーキュロクロム 1216, 2718 マーキュロクロム液 1217, 2718 マクリ 1846 マクロゴール400 1182 マクロゴール1500 1182 マクロゴール4000 1183 マクロゴール6000 1183 マクロゴール20000 1184 マクロゴール軟膏 1184 マシニン 1876 麻酔用エーテル 898 マニジピン塩酸塩 1192 マニジピン塩酸塩錠 1193 マプロチリン塩酸塩 1196 マルトース水和物 1190 D-マンニトール 1194, 2718 D-マンニトール注射液 1195

Ξ

ミグリトール 1253 ミグレニン 1254 ミクロノマイシン硫酸塩 1250 ミコナゾール 1249 ミコナゾール硝酸塩 1249 ミゾリビン 1262 ミゾリビン錠 1263 ミチグリニドカルシウム錠 1259 ミチグリニドカルシウム水和物 1258 ミツロウ 1809 ミデカマイシン 1251 ミデカマイシン酢酸エステル 1252 ミノサイクリン塩酸塩 1255 ミノサイクリン塩酸塩顆粒 2949 ミノサイクリン塩酸塩錠 1256 ミョウバン水 401

Δ

無コウイ大建中湯エキス 1917, 2777 無水アンピシリン 429 無水エタノール 896, 2693 無水カフェイン 543 無水クエン酸 716 無水乳糖 1132, 2714 無水ボウショウ 1994 無水リン酸水素カルシウム 557, 2911 ムピロシンカルシウム水和物 1279 ムピロシンカルシウム軟膏 1280

X

メキシレチン塩酸塩 1248 メキタジン 1214 メキタジン錠 1215 メグルミン 1206 メクロフェノキサート塩酸塩 1197 メコバラミン 1197 メコバラミン錠 1198 メサラジン 2718 メサラジン徐放錠 2720 メストラノール 1220, 2947 メタケイ酸アルミン酸マグネシウム 2716 メダゼパム 1200 メタンフェタミン塩酸塩 1223 L-メチオニン 1224 メチクラン 1241 メチラポン 1247 dl-メチルエフェドリン塩酸塩 1231 dl-メチルエフェドリン塩酸塩散10% 1232 メチルエルゴメトリンマレイン酸塩 1233 メチルエルゴメトリンマレイン酸塩錠 1233 メチルジゴキシン 1242 メチルセルロース 1228, 2723, 2947 メチルテストステロン 1239 メチルテストステロン錠 1240 メチルドパ錠 1230 メチルドパ水和物 1229 メチルプレドニゾロン 1236, 2949 メチルプレドニゾロンコハク酸エステ ル 1236 メチルベナクチジウム臭化物 1227 メチルロザニリン塩化物 1237 滅菌精製水(容器入り) 1774 メテノロンエナント酸エステル 1221 メテノロンエナント酸エステル注射液 1222 メテノロン酢酸エステル 1220 メトキサレン 1227 メトクロプラミド 1243 メトクロプラミド錠 1243 メトトレキサート 1225 メトトレキサート錠 2722 メトトレキサートカプセル 1225 メトプロロール酒石酸塩 1244 メトプロロール酒石酸塩錠 1245 メトホルミン塩酸塩 1222 メトホルミン塩酸塩錠 1223 メドロキシプロゲステロン酢酸エステ ル 1202 メトロニダゾール 1246 メトロニダゾール錠 1246 メナテトレノン 1209 メピチオスタン 1212 メピバカイン塩酸塩 1213 メピバカイン塩酸塩注射液 1213 メフェナム酸 1203 メフルシド 1204 メフルシド錠 1205 メフロキン塩酸塩 1203 メペンゾラート臭化物 1211 メルカプトプリン水和物 1215

Ŧ

メルファラン 1208

dl-メントール 1210

l-メントール 1211

メロペネム水和物 1217

木クレオソート 2008 モクツウ 1793 モサプリドクエン酸塩散 1276 モサプリドクエン酸塩錠 1277 モサプリドクエン酸塩水和物 1275 モッコウ 1972 モノステアリン酸アルミニウム 402 モノステアリン酸グリセリン 995 モルヒネ・アトロピン注射液 1273 モルヒネ塩酸塩錠 1272 モルヒネ塩酸塩水和物 1270 モルヒネ塩酸塩注射液 1271 モルヒネ硫酸塩水和物 1274 モンテルカストナトリウム 1264 モンテルカストナトリウム顆粒 2724 モンテルカストナトリウム錠 1269 モンテルカストナトリウムチュアブル 錠 1267

ヤ

ヤクチ 1812 ヤクモソウ 1905 薬用石ケン 1201 薬用炭 1200 ヤシ油 1835

ユ

ユウタン 1807 ユーカリ油 1850, 2975 輸血用クエン酸ナトリウム注射液 1572 ユビデカレノン 1743

\exists

ヨウ化カリウム 1437 ヨウ化ナトリウム 1579 ヨウ化ナトリウム(¹²³I)カプセル 1579 ヨウ化ナトリウム(131I)液 1580 ヨウ化ナトリウム(131I)カプセル 1580 ヨウ化人血清アルブミン(131I)注射液 1075 ヨウ化ヒプル酸ナトリウム(131I)注射 液 1580 葉酸 959, 2694 葉酸錠 960 葉酸注射液 959 ョウ素 1075 ヨクイニン 1836 ヨクイニン末 1836 抑肝散エキス 2010, 2798, 2983 ヨード・サリチル酸・フェノール精 1079 ヨードチンキ 1076 ヨードホルム 1080

ラ

ラウリル硫酸ナトリウム 1584, 2734 ラウロマクロゴール 1144 ラクツロース 1134 ラタモキセフナトリウム 1143 ラッカセイ油 1935 ラナトシド C 1137, 2715 ラナトシド C 錠 1138, 2715 ラニチジン塩酸塩 1503 ラノコナゾール 2944 ラノコナゾール外用液 2946 ラノコナゾールクリーム 2945 ラノコナゾール軟膏 2946 ラフチジン 1135 ラフチジン錠 1136 ラベタロール塩酸塩 1129 ラベタロール塩酸塩錠 1130 ラベプラゾールナトリウム 1502 ランソプラゾール 1139 ランソプラゾール腸溶カプセル 1140 ランソプラゾール腸溶性口腔内崩壊錠 1141

IJ

リオチロニンナトリウム 1165 リオチロニンナトリウム錠 1166 リシノプリル錠 1168 リシノプリル水和物 1167 L-リシン塩酸塩 1180 L-リシン酢酸塩 1179 リスペリドン 1521 リスペリドン細粒 1522 リスペリドン錠 1524 リスペリドン内服液 1523 リセドロン酸ナトリウム錠 1589 リセドロン酸ナトリウム水和物 1588 リゾチーム塩酸塩 1181 六君子湯エキス 1957, 2780 リドカイン 1161 リドカイン注射液 1162 リトドリン塩酸塩 1525 リトドリン塩酸塩錠 1526 リトドリン塩酸塩注射液 2954 リバビリン 1511 リバビリンカプセル 1512 リファンピシン 1518 リファンピシンカプセル 1519 リボスタマイシン硫酸塩 1517 リボフラビン 1513 リボフラビン散 1514 リボフラビン酪酸エステル 1514 リボフラビンリン酸エステルナトリウ لا 1515 1515 リボフラビンリン酸エステルナトリウ ム注射液 1516

リマプロスト アルファデクス 1162 レボフロキサシン水和物 1154 リュウガンニク 1907 リュウコツ 1907 リュウコツ末 1908 硫酸亜鉛水和物 1785 硫酸亜鉛点眼液 1785 硫酸アルミニウムカリウム水和物 403 硫酸カリウム 1438 硫酸鉄水和物 926 硫酸バリウム 473 硫酸マグネシウム水 1190 硫酸マグネシウム水和物 1189 硫酸マグネシウム注射液 1190 リュウタン 1884 リュウタン末 1884 流動パラフィン 1365 リュープロレリン酢酸塩 1150 リョウキョウ 1796 苓桂朮甘湯エキス 1961,2781 リンゲル液 1520 リンコマイシン塩酸塩水和物 1164 リンコマイシン塩酸塩注射液 1164 リン酸水素カルシウム水和物 558 リン酸水素ナトリウム水和物 1585 リン酸二水素カルシウム水和物 558

レセルピン 1507 レセルピン散0.1% 1508 レセルピン錠 1509 レセルピン注射液 1508 レチノール酢酸エステル 1510 レチノールパルミチン酸エステル 1510 レナンピシリン塩酸塩 1144 レノグラスチム(遺伝子組換え) 1146 レバミピド 1504 レバミピド錠 1506 レバロルファン酒石酸塩 1152 レバロルファン酒石酸塩注射液 1152 レボチロキシンナトリウム錠 1160 レボチロキシンナトリウム水和物 1159 レボドパ 1153 レボフロキサシン細粒 1155 レボフロキサシン錠 1157

レボフロキサシン注射液 1156 レボフロキサシン点眼液 1156 レボホリナートカルシウム水和物 2675 レボメプロマジンマレイン酸塩 1159 レンギョウ 1852 レンニク 1919

П

ロキサチジン酢酸エステル塩酸塩

L-ロイシン 1149

1530 ロキサチジン酢酸エステル塩酸塩徐放 カプセル 1531 ロキサチジン酢酸エステル塩酸塩徐放 錝 1532 ロキシスロマイシン 1534 ロキシスロマイシン錠 2733 ロキソプロフェンナトリウム錠 1178 ロキソプロフェンナトリウム水和物 1177 ロキタマイシン 1528, 2733 ロキタマイシン錠 1529, 2733 ロサルタンカリウム 1172 ロサルタンカリウム・ヒドロクロロチ アジド錠 1174 ロサルタンカリウム錠 1173 ロジン 1960 ロートエキス 1974 ロートエキス・アネスタミン散 1976 ロートエキス・カーボン散 1976 ロートエキス・タンニン坐剤 1978 ロートエキス・パパベリン・アネスタミ ン散 1977, 2790 ロートエキス散 1975 ロートコン 1973, 2982 ロベンザリットナトリウム 1171 ローヤルゼリー 1960 ロラゼパム 1171

ワ

ワイル病秋やみ混合ワクチン 1775 ワルファリンカリウム 1771 ワルファリンカリウム錠 1772